


Quadstar 32-bit

Application Program for
 Quadrupole Mass Spectrometers QMS 422 and QMS 200

Documentation



Publisher

PFEIFFER VACUUM
Berliner Str. 43
D-35614 Asslar
GERMANY

Website

<http://www.pfeiffer-vacuum.de>

Phone number

+49-(0)6441-802-0

Fax number

+49-(0)6441-802-202

Trademarks

Channeltron[®] is a registered trademark of Burle Industries.

Windows[®], Microsoft[®] and Excel[®] are registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Copyright

© INFICON AG 2003

Subject to technical alterations without notice.

No part of these Instructions may be reproduced in any way (photocopies, microfilms or any other reproduction processes) nor may it be manipulated with electronic systems, duplicated or distributed without written permission from INFICON AG. Offenders are liable to pay damages.

The original software delivered with the system and updated, state of the art versions of the system software are intended for use with the Quadstar 32-bit only. The software may not be used for purposes other than those intended nor is it permitted to make copies of the software. In particular, it is strictly forbidden to give copies of the software to other people.

The use of trade names, brand names, trademarks, etc. in these Instructions does not entitle third parties to consider these names to be unprotected and to use them freely. This is in accordance with the meaning of the laws and acts covering brand names and trademarks.

Table of contents

1.	Introduction	
1.1	About this manual.	17
1.1.1	What is its structure?	17
1.1.2	How do I work with it?	17
1.2	About Quadstar 32-bit	18
1.2.1	Product description	18
1.2.2	The individual programs	19
1.2.3	Quadstar 32-bit main features	19
1.2.4	Quadstar 32-bit «Quantitative Analysis» features.	20
1.2.5	Quadstar 32-bit «Process Control».	21
1.2.6	Scope of delivery	21
1.2.7	Computer system requirements	22
1.3	Abbreviations and terms	23
1.3.1	Hardware modules and components	23
1.3.2	Quadstar 32-bit specific	24
1.3.3	General	24
1.3.4	Conventions	25
1.4	Installation	26
1.4.1	Using the SETUP installation program	26
1.4.2	Installation/update on a hard disk	27
1.4.3	Connecting the mass spectrometer to the computer	28
1.4.4	Configuring Quadstar 32-bit	29
1.4.4.1	Communication with the mass spectrometer	30
1.4.4.2	Configuring the mass spectrometer	31
1.4.4.3	Changing the directories.	33
1.5	Getting started	34
1.5.1	Installation	34
1.5.2	Choose operating mode «Simulation»	34
1.6	Air analysis	36
1.6.1	Scan Analog overview measurement	36
1.6.2	Save measure data.	37
1.6.3	Analyzing saved measure data.	38
1.7	Measuring time dependent processes	41
1.7.1	Basic measuring techniques	41
1.7.1.1	The MID measurement.	41
1.7.1.2	How the MCD measurement works	43
1.7.1.3	Preparatory measurements	44
1.7.2	Preparing the measurements	44
1.7.2.1	Mass Scale Calibration	44
1.7.2.2	Gas Specific Calibration	44
1.7.2.3	Analysis file	45
1.7.3	Measurements	45
1.7.3.1	Calibrating the mass scale	45

1.7.3.2	Measure the zero gas	46
1.7.3.3	Calibrate the mass spectrometer	46
1.7.3.4	Measure and store concentrations	46
1.7.4	Analysis	48
1.7.4.1	Display modes	48
1.7.4.2	Display additional data as a table	50
1.7.4.3	Print data	51
1.7.5	Auxiliary tools	53
1.7.5.1	Documentation	53
1.7.5.2	Filter	54
1.7.5.3	Interpretation support	54
1.8	Measuring real gas components	55
1.9	Hints concerning the work with Quadstar 32-bit	56
1.9.1	Spectrum interpretation	56
1.9.2	Starting the Quadstar 32-bit programs	57
1.9.2.1	Open a program with a certain file	57
1.9.2.2	Automatic start of a program	58
1.9.3	Transferring screens via the clipboard	58
1.9.3.1	Copying an entire screen	58
1.9.3.2	Copying an active window	58
1.9.4	Inserting screens from the clipboard	59
1.9.5	Save or print Quadstar 32-bit screens	59

2. Measure

2.1	General information	61
2.1.1	Measuring with Quadstar 32-bit	61
2.1.1.1	Measurement types	61
2.1.1.2	Measurement methods	62
2.1.1.3	FIR Filter	62
2.1.1.4	Fast MID/MCD measurements	62
2.1.2	Choosing the measuring range	64
2.1.3	Status bar	65
2.1.3.1	Message text output during the measurement	65
2.1.4	Configuring the graphic display range	65
2.1.4.1	Upper limit of display (Disp. F.S.R.)	66
2.1.4.2	Number of decades to display (Disp. Decades)	66
2.1.4.3	Display range in sequencer operation	67
2.1.4.4	Y-axis lettering (MID/MCD)	68
2.1.5	Documenting the measurement data (File Info)	69
2.1.6	Multiplex operation with several control units	70
2.1.7	Calibrating the measurement equipment	71
2.1.8	Measurement of analog input signals	72
2.1.9	Total pressure monitoring	73
2.2	Measurement main menu	74
2.3	File menu	75
2.3.1	About Measure	75
2.3.2	Exit	75
2.4	Scan menu	76
2.4.1	Setup submenu	76
2.4.2	Scan Analog	77
2.4.2.1	Scan Analog measurement menu	78
2.4.2.2	File submenu	78
2.4.2.3	Operation submenu	79

2.4.2.4	Parameters submenu	79
2.4.2.5	Setup parameters	80
2.4.2.6	Channel parameters	81
2.4.3	Scan Bargraph	82
2.4.3.1	File submenu	83
2.4.3.2	Operation submenu	84
2.4.3.3	Parameters submenu	84
2.4.3.4	Setup parameters	85
2.4.3.5	Channel parameters	86
2.5	MID measurement type.	87
2.5.1	Setup submenu	87
2.5.2	MID measurement menu	88
2.5.2.1	File submenu	88
2.5.2.2	Operation submenu	89
2.5.2.3	Display submenu	89
2.5.2.4	Parameters submenu	90
2.5.3	Offset alignment	90
2.5.4	Background subtraction (Zero Gas)	91
2.5.5	Temporary save	91
2.5.6	MID Table	92
2.5.6.1	Parameter setup	92
2.5.6.2	Data display	93
2.5.7	MID Bargraph	96
2.5.7.1	Parameter setup	96
2.5.7.2	Data display	98
2.5.8	MID Versus Time	98
2.5.8.1	Parameter setup	99
2.5.8.2	Data display	100
2.5.9	MID last measured values	101
2.5.9.1	File submenu	101
2.5.9.2	Display submenu	101
2.5.9.3	Parameters submenu	102
2.6	MCD measurement type.	103
2.6.1	Setup submenu	104
2.6.2	MCD measurement menu	105
2.6.2.1	File submenu	105
2.6.2.2	Operation submenu	106
2.6.2.3	Display submenu	106
2.6.2.4	Parameters submenu	107
2.6.3	Data display	107
2.6.4	Offset alignment	108
2.6.5	Background subtraction	108
2.6.6	Operation with the QMU 112	108
2.6.7	MCD Table	108
2.6.7.1	Parameter setup	109
2.6.7.2	Data display	110
2.6.8	MCD Bargraph	112
2.6.8.1	Parameter setup	113
2.6.8.2	Data display	115
2.6.9	MCD Versus Time.	115
2.6.9.1	Parameter setup	116
2.6.9.2	Data display	118
2.6.10	MCD last measured values.	118
2.6.10.1	File submenu	118
2.6.10.2	Display submenu	119
2.6.10.3	Parameters submenu	119
2.7	Measuring with the sequencer (Sequence)	120

2.7.1	Setup submenu	120
2.7.2	Execute submenu	122
2.7.2.1	Testing a sequence	122
2.7.2.2	Status display	122
2.7.2.3	Display of measurement data	123
2.7.3	Controlling the display via sequencer instructions	125
2.7.3.1	Measurement data of the MID() instruction	125
2.7.3.2	Scan measurement data	127
2.7.3.3	Measurement data of the MCD() instruction	129
2.7.3.4	Measurement data of the GSC() instruction (gas specific calibration)	131
2.7.3.5	Measurement data of the ZeroGas() instruction (zero gas measurement)	131
2.7.3.6	Measurement data of the MSC() instruction (mass scale calibration)	132
2.7.4	User-controlled sequence	133
2.8	Calibration	134
2.8.1	Mass scale calibration (Mass Scale)	134
2.8.1.1	File submenu	136
2.8.1.2	Operation submenu	136
2.8.1.3	Parameters submenu	136
2.8.1.4	Setup	137
2.8.1.5	Channel	137
2.8.2	Background measurement (Zero Gas)	137
2.8.2.1	File submenu	139
2.8.2.2	Operation submenu	139
2.8.2.3	Parameters submenu	139
2.8.3	QMS Offset	140
2.8.3.1	QMS Offset Calibration Setup	141
2.8.4	Gas specific sensitivity	142
2.8.4.1	File submenu	143
2.8.4.2	Operation submenu	144
2.8.4.3	Parameters submenu	144
2.9	Manual Parameters	145
2.9.1	DI/DO Manual for QMS 422/421	145
2.9.2	DI/DO Manual for QMS 200	146
2.10	Communication menu	147
2.11	Setup	148

3. Parameter Setup

3.1	File	153
3.1.1	System Information	153
3.1.2	About Parameter Setup	154
3.1.3	Print Settings	154
3.1.4	Exit	155
3.2	Measure Parameters	156
3.2.1	Selecting the parameter file	156
3.2.2	Storing the parameters	156
3.2.3	Create or copy a parameter file	157
3.2.4	The channel parameter editor	158
3.2.4.1	Measuring parameters overview	159
3.2.4.1.1	'Main' group	159
3.2.4.1.2	'Detector' group	160
3.2.4.1.3	'Mass' group	160
3.2.4.1.4	'Amplifier' group	161

3.2.4.1.5	'Output' group	161
3.2.4.1.6	'Trip' group	162
3.2.4.1.7	'Display' group	162
3.2.4.1.8	'<<Rest>>' group	162
3.2.4.2	Modification of parameters	162
3.2.4.3	Configurations (Parameter groups/register cards)	164
3.2.5	Channel Parameter Editor menu	168
3.2.5.1	File submenu	169
3.2.5.2	Edit submenu	170
3.2.5.3	Display submenu	171
3.2.5.3.1	Configuration	172
3.2.5.3.2	Font selection	173
3.2.5.3.3	Measure Data	174
3.2.5.4	Functions submenu	174
3.2.5.4.1	Adjust function	175
3.2.5.4.2	Peak Search function	178
3.2.5.5	Editing the MCD parameter files (Quantitative Analysis)	180
3.2.5.5.1	Analysis matrix	180
3.2.5.5.2	Expansion of the analysis matrix	181
3.2.5.5.3	MCD Measure Parameter menu	181
3.2.5.6	Activation of analog input channels (MID/MCD)	186
3.2.5.7	Switching functions (Trip)	186
3.2.5.8	Outputting the switching function state to a DO channel	187
3.2.5.8.1	Behavior of the switching functions	188
3.3	Sequence programming	190
3.3.1	Quadstar 32-bit Sequence editor	190
3.3.1.1	Programming	191
3.3.1.2	Inserting a new Sequencer line (Choose Line)	191
3.3.1.3	Compiling/storing the sequences	192
3.3.1.3.1	Sequence Editor	192
3.3.1.4	Parameter limits (array limits)	197
3.3.1.4.1	Variable types	197
3.3.1.5	Parameter Files	198
3.4	Calibration Parameters	199
3.4.1	Mass Scale	199
3.4.2	Gas-specific sensitivity (Quantitative Analysis)	200
3.4.2.1	Calibration matrix	201
3.4.2.2	Extending the calibration matrix	202
3.4.3	Gas Specific Calibration menu	202
3.4.3.1	File submenu	203
3.4.3.2	Edit submenu	203
3.4.3.2.1	Channel Parameters	204
3.4.3.3	Options submenu	204
3.4.3.3.1	Defining the internal standard	205
3.4.3.3.2	Setting the external normalization	205
3.5	Manual parameters	207
3.5.1	DI/DO Manual (for QMS 422/421)	207
3.5.2	DI/DO Manual (for QMS 200)	208
3.6	Communication parameters	209
3.6.1	Overview	209
3.6.2	Communication via RS232C	210
3.6.2.1	Protocol/transmission reliability	210
3.6.2.2	Configuration of multiple serial interfaces under Windows	210
3.6.2.3	RS232C normal	210
3.6.2.4	Communication via ArcNet	211
3.6.2.4.1	Installation of the ArcNet board	212

3.6.2.4.2	Multiplex operation	212
3.6.2.4.3	Installation of the hub	212
3.6.2.4.4	Data collection possibilities	212
3.6.2.4.5	Unit-specific files	213
3.6.3	Setting the interface parameters	214
3.6.3.1	Selecting the communication type	214
3.6.3.2	Setting the parameters for the selected communication type	214
3.6.3.3	QMS communication type	214
3.6.3.4	QMS communication parameters	215
3.6.3.4.1	Communication via RS232C Normal	215
3.6.3.4.2	Communication via ArcNet	216
3.6.3.5	Connect	221
3.6.3.6	Disconnect	222
3.7	Setup parameters	223
3.7.1	SEM/Emission Control	223
3.7.1.1	Normal Operation (without QMU 112)	224
3.7.1.2	Operation with QMU 112	225
3.7.2	General setup parameters	226
3.7.2.1	Simulation mode	228
3.7.3	Printer	229
3.7.3.1	Configuring the printer	229
3.7.3.2	Set up the printer page (Page Setup)	230
3.7.3.3	Select the printer font	230
3.7.3.4	Print test output	231
3.7.3.5	Log output	231
3.7.4	ASCII Conversion	232
3.8	Instrument Configuration	234
3.8.1	QMS Configuration	234
3.8.1.1	First time operation, configuration changes	236
3.8.2	Digital inputs and outputs	237
3.8.2.1	QMS 422/421	237
3.8.2.2	QMS 200	237
3.8.2.3	DI/DO identification	237
3.8.3	Analog input characteristic curve configuration	238
3.8.3.1	Overview	238
3.8.3.2	Specifying an AI sensor	239
3.8.4	AO Characteristic Curve Configuration	240
3.8.4.1	Overview	240
3.8.4.2	Specifying an analog output	240
3.8.4.3	Programming the characteristic	240
3.8.4.4	Enter the characteristic curves (AI)	242
3.8.5	Screen Colors	245
3.8.5.1	Changing the colors of the individual screen elements	245

4. Display Saved Values

4.1	File	250
4.1.1	About Display Saved Values	250
4.2	Scan measurement data	251
4.2.1	Scan Analog data	251
4.2.1.1	File submenu	254
4.2.1.2	Display submenu	255
4.2.1.3	Select submenu	256
4.2.1.4	Setup submenu	257
4.2.1.5	Function submenu	260

4.2.1.5.1	Magnify	261
4.2.1.5.2	Cursor	263
4.2.1.5.3	Graphic editor	267
4.2.1.6	Special submenu	269
4.2.1.6.1	Filter	269
4.2.1.7	Info submenu	270
4.2.2	Scan Bargraph Cycles	271
4.2.2.1	File submenu	272
4.2.2.2	Mode submenu	273
4.2.2.2.1	TIC display mode (Total Ion Current)	274
4.2.2.2.2	IVT display mode (Intensity Versus Time)	274
4.2.2.3	Display submenu	275
4.2.2.3.1	Cycle	276
4.2.2.4	Select submenu	277
4.2.2.4.1	Data	277
4.2.2.5	Setup the display	279
4.2.2.5.1	X-Axis	280
4.2.2.5.2	Y Axis	280
4.2.2.6	Function submenu	281
4.2.2.6.1	Magnify a Screen Section	282
4.2.2.6.2	Determine measurement data graphically (Cursor)	282
4.2.2.6.3	Display format of [Function] > [Table]	283
4.2.2.6.4	Display format of [Function] > [Spectrum]	284
4.2.2.6.5	Graphic editor	284
4.2.2.7	Special submenu	284
4.2.2.8	Info submenu	284
4.3	Process	286
4.3.1	Cycle/Reference measurement data	286
4.3.1.1	File submenu	288
4.3.1.2	Display setup	289
4.3.1.2.1	Cycle	289
4.3.1.2.2	Statistical Functions	289
4.3.1.2.3	Integration (Integral)	291
4.3.1.2.4	Bargraph	293
4.3.1.2.5	Intensity Versus Time (IVT)	294
4.3.1.2.6	Graphic Data	294
4.3.1.3	Select data to be displayed	294
4.3.1.3.1	Select Blocks	295
4.3.1.3.2	Select Data	296
4.3.1.4	Setup the display	298
4.3.1.4.1	X-Axis	298
4.3.1.5	Function submenu	299
4.3.1.5.1	Magnifying a Screen Section	300
4.3.1.5.2	Determine measurement data graphically (Cursor)	301
4.3.1.5.3	Graphic editor	303
4.3.1.6	Special submenu	303
4.3.1.6.1	Filter	303
4.3.1.7	Info submenu	304
4.4	Auxiliary	304
4.4.1	Mass scale calibration	304
4.4.1.1	File submenu	305
4.4.1.2	Edit submenu	306
4.4.1.2.1	Masses	306
4.4.1.3	Info submenu	308
4.4.1.3.1	File	309
4.4.2	Background measurement (Zero Gas)	309
4.4.2.1	File submenu	310

4.4.2.2	Edit submenu	310
4.4.2.2.1	Values	310
4.4.2.3	Info submenu	314
4.4.2.3.1	File	314
4.4.3	Error History (Stored Error Messages)	315

5. Tune Up

5.1	File	318
5.1.1	About Tune Up	318
5.1.2	Exit	318
5.2	Tune	319
5.2.1	Ion Source	319
5.2.1.1	Parameter setup for the QMA 200/QMA 4XX	320
5.2.1.1.1	QMA 200 (Prisma™)	320
5.2.1.1.2	QMA 4XX	320
5.2.1.2	Modifying parameters	321
5.2.1.2.1	Channel parameters	321
5.2.1.2.2	Filament	322
5.2.1.2.3	RF supply (QMS 200)	322
5.2.1.2.4	Ion source voltages	322
5.2.1.2.5	Filament Protection Extern (QMA 200/QMA 4XX)	323
5.2.1.3	QMA 125	323
5.2.1.3.1	Standard editor (Edit Channel Par)	326
5.2.1.3.2	Copy/Exch. IS-Set	327
5.2.1.3.3	Copy/Exch. Channel	328
5.2.1.3.4	Default Ion Source Set	329
5.2.1.3.5	Optimizing the parameters (Measure Data)	329
5.2.2	Peak Shape Optimizer	333
5.2.2.1	Channel parameters	333
5.2.2.2	Ion source voltages	333
5.2.2.3	Display width	334
5.2.2.4	Mean value	334
5.2.2.5	Start measurement (Peak Shape)	335
5.2.3	QMS 200 Tune Data	337
5.2.4	QMS 200 Tune Mass Scale	338
5.3	Manual parameters	341
5.3.1	Digital I/O	341
5.4	Communication Parameters	343
5.5	Setup Parameters	344
5.5.1	Degas Control	344
5.5.1.1	Operation with QMA 200/QMA 4XX	344
5.5.1.2	Operation with QMA 125	345
5.6	PPM Tune Up	349
5.7	PPM 422 Plasma Process Monitor	349
5.7.1	Functional diagram Ion-optical potentials	350
5.7.2	Parameter overview	350
5.8	File	352
5.8.1	File menu	352
5.8.1.1	About Tune Up	352
5.9	Tune	353
5.9.1	PPM 422	353

5.9.1.1	Parameter setup	353
5.9.1.1.1	Parameter setup of an E-Scan	354
5.9.1.1.2	Parameter setup of an M-Scan	355
5.9.1.1.3	Modifying Parameters	356
5.9.1.2	Copy/Exch. IS-Set.	359
5.9.1.3	Copy/Exch. Channel	360
5.9.1.4	Default IS-Set	361
5.9.1.5	Single/Common Incr/Decr.	362
5.9.1.6	Relations (Ion energy and Ve, CORR and SHIFT)	363
5.9.1.7	Measurement parameters for determining the energy zero	363
5.9.1.7.1	Energy scan in Sampling Mode Neutrals	364
5.9.1.7.2	Effect of the correction value CORR.	365
5.9.1.8	Measurement parameters for positive ions.	366
5.9.1.8.1	Energy distribution of positive ions	367
5.9.1.8.2	Effect of SHIFT	368
5.9.2	Optimizing the parameters (Measure Data)	369
5.9.2.1	Measure Data menu	371
5.9.2.1.1	File submenu	371
5.9.2.1.2	Operation submenu	372
5.9.2.1.3	Function submenu (crosshairs).	372
5.9.2.1.4	Setup submenu.	373
5.9.2.1.5	Display submenu	374
5.10	Communication Parameters	375
5.11	Setup Parameters	376
5.11.1	Turning the SEM high voltage and emission on/off	376

6. Service

6.1	File	378
6.1.1	About Service	378
6.2	Measure	379
6.2.1	Leak Test	379
6.2.1.1	Setup.	380
6.2.1.2	Channel.	381
6.2.2	Total Pressure.	382
6.2.2.1	QMS 422/421: PE 420 and/or PI 420	383
6.2.2.2	QMS 200: TPR 250 or PKR 250/251	383
6.2.2.3	Setup.	383
6.3	Manual	386
6.3.1	Digital I/O (DI/DO).	386
6.3.2	Analog I/O (AI/AO)	387
6.4	Communication Parameters	389
6.5	Setup Parameters	389

7. Utility

7.1	File menu	392
7.1.1	About Utility.	392
7.2	Spectra library menu	393
7.2.1	Open Library	394
7.2.2	Displaying a spectrum.	395

7.2.2.1	Spectrum submenu	396
7.2.2.1.1	Load a comparison spectrum (Select/Load S2)	396
7.2.2.1.2	Display spectrum as a table (Display Numeric)	397
7.2.2.1.3	Overlap S2 on S1	398
7.2.2.1.4	Store spectrum (S1 to Library)	398
7.2.2.1.5	Subtract S1 - S2	399
7.2.2.1.6	Subtraction List	400
7.2.2.1.7	Undo last Subtraction	400
7.2.2.1.8	Normalize S1	400
7.2.2.2	Copying a spectrum	401
7.2.2.3	Removing a spectrum	402
7.2.2.4	Searching a spectrum	402
7.2.2.5	Go to a spectrum	403
7.2.2.6	Sorting the library	403
7.2.2.7	Entering a new spectrum	403
7.2.2.8	Editing a spectrum	404
7.2.2.9	Printing the spectra	405
7.2.3	Export Library to Text File	406
7.3	Installation menu	408

8. Process Control

8.1	Introduction	409
8.2	Programming	410
8.2.1	Input via the keyboard	410
8.2.2	Input via dialog editor	410
8.3	Variables	411
8.3.1	Local integer variables i[.]	412
8.3.2	Global integer variables gi[.]	412
8.3.3	Local float variables f[.]	412
8.3.4	Global float variables gf[.]	413
8.3.5	Global field of float variables gfa[.][.].	413
8.3.6	Global character string gs[.]	413
8.4	Sequencer instructions	414
8.4.1	General information	414
8.4.1.1	Syntax of a sequencer instruction	414
8.4.1.2	Default values	415
8.4.1.3	File names with indices	416
8.4.1.4	String values	417
8.4.1.5	Syntax of a sequence	418
8.4.2	Control structures	419
8.4.2.1	If-Else instructions	419
8.4.2.2	Loop instructions	420
8.4.2.2.1	Loop()	420
8.4.2.2.2	While instructions	421
8.4.3	Log output	422
8.4.4	Data exchange with other programs (DDE)	424
8.4.4.1	Passive data exchange	424
8.4.4.2	Active data exchange	424
8.5	Auxiliary instructions	426
8.5.1	Definition of blocks	426
8.5.1.1	Begin instruction block	426
8.5.1.2	End instruction block	426
8.5.1.3	Entering a blank line	427

8.5.2	Else instruction	427
8.5.3	Comments	428
8.5.3.1	Comments comprising multiple lines	428
8.5.3.1.1	/* - Begin comment	428
8.5.3.1.2	*/ - End comment	428
8.5.3.2	// - Single line comment	428
8.6	Example of a sequence	430
8.6.1	Sequence AIRDEMO	430
8.6.2	User-controlled analysis (AIRDEMO2)	431
8.6.2.1	Quadstar 32-bit dialog box	431
8.6.2.2	Sequence AIRDEMO2	431
8.6.2.2.1	Description of AIRDEMO2	434
8.7	Sequencer instructions	439
8.7.1	Calculate() - Perform a calculation	439
8.7.2	DDEClose() - Close DDE Service	441
8.7.3	DDEExecute() - Call an instruction in an application by DDE	441
8.7.4	DDEInit() - Initialize DDE Service	442
8.7.5	DDEPoke() - Send data by DDE	442
8.7.6	DDERequest() - Request data by DDE	443
8.7.7	Delay() - Time delay	444
8.7.8	Dialog() - Displaying a dialog box	445
8.7.9	LoopDialog() - Loop with Dialog window	448
8.7.10	DispEventWindow() - Place and scale the Event window	448
8.7.11	DispWindow() - Place and scale the window	449
8.7.12	Error() - Outputting an error message	450
8.7.13	Execute() - Execute a program	451
8.7.14	Exit() - Quitting the sequencer	451
8.7.15	ExtractString() - Extracting characters from a string	452
8.7.16	FileCopy() - Copy a file	453
8.7.17	FileDelete() - Delete a file	453
8.7.18	FileMove() - Move a file	454
8.7.19	FileRename() - Rename a file	454
8.7.20	FormatDate() - Format date	455
8.7.21	FormatTime() - Format time	456
8.7.22	GetDate() - Get current date	457
8.7.23	GetTime() - Get current time	457
8.7.24	GetErrorCode() - Read error and status register	459
8.7.25	GetSysInfo() - Get System Information	462
8.7.26	GSC() - Gas-specific calibration	463
8.7.27	IfBit() - Testing a bit of an integer variable	464
8.7.28	IfDI() - Testing a digital input	465
8.7.29	IfDO() - Testing a digital output	466
8.7.30	IfMeaVal() - Testing a measured value	467
8.7.31	IfString() - Compare two text strings	469
8.7.32	IfVar() - Testing an array or field variable	469
8.7.33	IniFileRead() - Read a value from an INI file	470
8.7.34	IniFileWrite() - Write a value to an INI file	471
8.7.35	InitDisp() - Initialize the display	473
8.7.36	Loop() - Loop	473
8.7.37	LoopTable() - Loop Table	474
8.7.38	MCD() - MCD measurement	475
8.7.39	Message() - Message output	477
8.7.40	MID() - MID measurement	478
8.7.41	MSC() - Mass scale calibration	480
8.7.42	OffsetCal() - QMS Offset calibration	481
8.7.43	Process() - Data processing	482

8.7.44	QMSConnect() - Switch on the communication to the mass spectrometer	485
8.7.45	QMSDisconnect() - Switch off the communication to the mass spectrometer	485
8.7.46	ScanAnalog() - Scan Analog	485
8.7.47	ScanBar() - Scan Bargraph	487
8.7.48	SeqExit() - Defining a sequence for sequence cancellation	488
8.7.49	Sequence() - Starting a (sub-)sequence	489
8.7.50	SetAO() - Setting an analog output channel	490
8.7.51	SetDefault() - Set default values	491
8.7.52	SetDO() - Set a digital output channel	491
8.7.53	SetFileInfo() - Changing the measure data documentation settings	493
8.7.54	SetISPar() - Set Ion Source Parameters of the current unit	494
8.7.55	SetPar() - Setting Quadstar 32-bit parameters	496
8.7.56	SetParFile() - Set value(s) in a channel parameter file	499
8.7.57	SetQMS() - Setting the QMS parameter(s)	500
8.7.58	SetScale() - Define a scale	502
8.7.59	SetString() - Setting a character string	503
8.7.60	SetUnit() - Selecting the mass spectrometer in multiplex mode	504
8.7.61	SetVar() - Setup of variables	505
8.7.62	Sound() - Tone output	506
8.7.63	TextClear() - Clear all text	507
8.7.64	TextClose() - Close text window	507
8.7.65	TextInit() - Open a text window	507
8.7.66	TextSet() - Write a text to the text window	508
8.7.67	WhileDI() - While loop for monitoring a digital input	510
8.7.68	WhileDO() - While loop for monitoring a digital output	512
8.7.69	WhileMeaVal() - While loop for monitoring a measured value	513
8.7.70	WhileString() - While loop for monitoring a string	515
8.7.71	WhileVar() - While loop for monitoring a variable	517
8.7.72	ZeroGas() - Zero gas measurement	519
8.8	Process Monitor PPM 422	521
8.8.1	PPMGetCorr() - Get correction value CORR	521
8.8.2	PPMMID() - PPM 422 MID measurement	521
8.8.3	PPMNeutrals() - Neutral particle analysis	523
8.8.4	PPMScanEnergy() - PPM 422 Energy Scan	523
8.8.5	PPMScanMassAnalog() - PPM 422 Scan Mass Analog	525
8.8.6	PPMSetEnergy() - Entering the ion energy	526
8.8.7	PPMSetEPolarity() - Polarity of the energy voltage	527
8.8.8	PPMSetIS() - Entering the ion source parameters	527
8.8.9	PPMSetPSM() - Entering the PPM 422 Plasma Sampling Mode (PSM)	530

9. Access Rights Control

9.1	Login	532
9.1.1	Login to a Quadstar 32-bit program	532
9.1.1.1	Forgotten user password	532
9.1.1.2	Logging into ACCESSC	533
9.1.1.3	Forgotten ACCESSC password	533
9.1.2	Safeguarding the password protection	533
9.2	Control of the access privileges	534
9.2.1	Selecting the Quadstar 32-bit program	534
9.2.2	Defining the access privileges	534
9.2.3	Switch the password protection on/off	535

9.2.4	Main menu	535
9.2.4.1	File menu	536
9.2.4.1.1	About Access Rights Control	536
9.2.4.2	User menu	536
9.2.4.2.1	Change name	537
9.2.4.2.2	Insert	537
9.2.4.2.3	Delete	538
9.2.4.2.4	Copy	538
9.2.4.2.5	Define the work environment (Path Setup)	538
9.2.4.2.6	Set Predefined Access	539
9.2.4.2.7	Delete user Password	539
9.2.5	Protection menu	539
9.2.5.1	Administrator password (ACCESSC Login Setup...)	539

10. File Handling, Data Formats, Error Messages, Troubleshooting

10.1	File Manager	541
10.2	Path and file names	544
10.3	Quadstar 32-bit files	544
10.3.1	Overview	544
10.4	File formats	549
10.5	Data export ASCII	549
10.5.1	Scan Analog data	549
10.5.2	Scan Bargraph data	550
10.5.3	MID/MCD Data	552
10.6	Error messages/warnings	553
10.6.1	General	553
10.6.2	Quadstar 32-bit Error Messages/Warnings	553
10.7	Troubleshooting	573

11. Residual Gas Analysis Program for Prisma™

11.1	Overview of the operating modes	576
11.1.1	Main selection box	576
11.1.1.1	Scan	577
11.1.1.2	Currents	577
11.1.1.3	Concentr.	578
11.1.1.4	Part. Press.	578
11.1.1.5	Auxiliaries	578
11.2	Working with the RGA program	580
11.2.1	Start of the RGA program	580
11.2.2	Introduction	581
11.2.3	Switching on the emission	582
11.2.4	Display of an overview spectrum	583
11.2.5	Switching on the Channeltron	584
11.2.5.1	Calibrating the Channeltron	585
11.3	Main Selection Box	586
11.3.1	Terminate RGA program [EXIT SEQUENCE]	586
11.3.2	Back to previous screen	587

11.3.3	Measurement programs	587
11.3.3.1	Scan measurements	587
11.3.3.1.1	Enter measure and display parameters	587
11.3.3.1.2	Storage of measured data	588
11.3.3.1.3	Analog	589
11.3.3.1.4	Bargraph	590
11.3.3.1.5	Stair	591
11.3.3.2	Currents	591
11.3.3.2.1	Table	591
11.3.3.2.2	Bargraph	592
11.3.3.2.3	Versus Time	593
11.3.3.3	Concentr.	593
11.3.3.3.1	Table	594
11.3.3.3.2	Bargraph	594
11.3.3.3.3	Versus Time	595
11.3.3.4	Part. Press. (Partial Pressures)	595
11.3.3.4.1	Table	595
11.3.3.4.2	Bargraph	596
11.3.3.4.3	Versus Time	597
11.3.3.5	Auxiliaries	597
11.3.3.5.1	Calibrate the mass scale	598
11.3.3.5.2	Determine the residual gas background	598
11.3.3.5.3	Perform a helium leak test	601
11.3.3.5.4	Change the mass resolution	602
11.3.3.5.5	Calibrate the SEM voltage	602
11.3.3.5.6	Calibrate the amplifier offset	602
11.3.3.5.7	Determine the mass spectrometer sensitivity	603

12. Sputter Process Monitor Program

12.1	Summary of the different operation modes	606
12.1.1	The Measurement Selection Box	606
12.2	General functions of the SPM 200 sequence	608
12.2.1	Starting the SPM 200 sequence	608
12.2.2	Leaving the SPM 200 sequence	608
12.2.3	Jumping backwards within the SPM 200 sequence	609
12.2.4	Switching on the filament and the emission	609
12.2.5	Switching on the Channeltron Voltage	611
12.2.6	Calibration of the Channeltron Voltage	611
12.2.7	Display of an Overview Spectrum	612
12.2.8	Saving of measurement data	613
12.2.8.1	Saving data from Analog Scans	613
12.2.8.2	Assigning data file names	614
12.2.8.3	File Manager	614
12.3	Impurity determinations	616
12.3.1	Selection of different Display Modes	616
12.3.1.1	Adapting the concentration scale for Versus Time displays	617
12.4	Residual Gas Measurements	619
12.4.1	Selection of different Display Modes	620
12.4.1.1	Adapting the pressure scale for Residual Gas measurements	621

1. Introduction


1.1 About this manual


1.1.1 What is its structure?


In this chapter you'll find the product description, an installation guide and an introduction to the work with Quadstar 32-bit.

The chapters 2 to 7 contain detailed descriptions of the individual programs (Measure, Service etc.).



Chapter 8. Process Control,  409 is the reference guide for sequence-programming.

Chapter 9. Access Rights Control,  531 shows you how to setup permission of access to Quadstar 32-bit for different users.

Chapter 10. File Handling, Data Formats, Error Messages, Troubleshooting,  541 describes data formats, error messages and troubleshooting in Quadstar 32-bit.

Chapter 11. Residual Gas Analysis Program for Prisma™,  575 explains the Residual Gas Analysis program for QMS 200 (Prisma™).

1.1.2 How do I work with it?

- Use it mainly as a reference book. It doesn't make too much sense to read it from the beginning to the end.
- If you are definitely not familiar with Quadstar 32-bit, you should study thoroughly the Chapter 1.5 Getting started,  34; the basic processes are briefly explained there.
- If you have questions concerning a certain part of the program while you are working, search for this part in the table of contents of the program and look it up there.
- If you have questions you can not associate to a certain program part, then find matching keywords and search for them in the Index. Perhaps you have to try several different keywords until you find the reference you are looking for.
- If Quadstar 32-bit doesn't work the way it should, you might find help in Chapter 10.7 Troubleshooting,  573. Often encountered problems and possible solutions are described there.

1.2 About Quadstar 32-bit

1.2.1 Product description

Quadstar 32-bit is a powerful, modular application software for the quadrupole mass spectrometers QMS 422/421 and QMS 200 Prisma™. Quadstar 32-bit (standard) is capable of performing qualitative and quantitative analyses and it is user programmable.

Quadstar 32-bit simplifies the operation of the mass spectrometer, improves the reproducibility of the measurements, and allows all relevant parameters and data to be stored. In addition, the consistency of certain parameters is checked and warning messages are output in case of error.

Quadstar 32-bit allows you to store measurement data, to read them back in, and to represent them in various formats (table, bargraph, intensity/concentration versus time).

Quadstar 32-bit has two possibilities to communicate with the mass spectrometer:

- Serial via RS232C, available on all PC systems as a standard feature. The required connection cable is provided. Only one QMS connectable, medium transmission speed.
- Using fiber optics and ArcNet (has to be ordered separately). Several QMSs may be operated at the same time at very high transmission speed.

Main features

The Main features of the Quadstar 32-bit package are:

- Ease of operation through graphic user interface
- Powerful measurement programs for:
 - Qualitative analysis: Multiple Ion Detection (MID), Scan Analog, Scan Bargraph
 - Quantitative analysis: Multiple Concentration Detection (MCD)
 - Optimized analysis of the stored measurement data
 - 3-dimensional representation of Scan Analog spectra
 - Representation of multiple channel measurements as table, bargraph or Versus Time
 - Documentation of the measurements is supported
- Measuring simultaneously with several control units via the ArcNet interface
- Freely programmable, automatic measurement sequences
- Online data interchange with other programs by DDE
- Comprehensive service functions
- Online parameter setup of the mass spectrometer
- Post-measurement storage of data (temporary save)
- Threshold monitoring with alarm functions
- Automatic electrometer range selection
- Programmable analog and digital inputs and outputs
- Support for special QMS functions such as Chopper Lock In, Ion Counter etc.

1.2.2 The individual programs

Quadstar 32-bit currently consists of the following main programs:

Measure

Measure is the actual measure program. It provides all kinds of measurements and storage of the acquired data. From here, sequences are started as well.

Dispsav

Dispsav (Display Saved Values) is the analyzing program. Stored data can be represented in different ways, magnified and processed.

Parset

Parset (Parameter Setup) is the setup program. All kinds of settings for measurements, QMS, communication and the working environment can be done here. You can program sequences here.

Tune Up

Tune Up is the tuning program. You can tune the ion source, adjust the RF generator and optimize the peak shape or degas the filament.

Service

Service is the service program. You can test for leaks, measure total pressure and operate analog or digital inputs and outputs manually.

Utility

Utility is an additional program. It contains functions for the spectra library.

Accessc

Accessc is the administration program. You can allow or prohibit access to Quadstar 32-bit components for individual users.

1.2.3 Quadstar 32-bit main features

Quadstar 32-bit provides qualitative analysis of measurement data. It also contains various auxiliary and service programs such as leak test, service scan, total pressure measurement, etc. The main features of Quadstar 32-bit are:

Powerful measurement programs

- Multiple Ion Detection (MID) with simultaneous representation of ion currents and external analog signals. Also threshold monitoring, integration, auto ranging, etc.
- Scan Analog with online zero gas or reference subtraction, mean value calculation across several cycles if desired
- Scan Bargraph for fast peak detection and display across a certain mass range

Storage functions

- Online storage of an unlimited number of data files
- Post-measurement storage of data through Temporary Save feature
- Documenting measurement data by storing parameters and remarks concerning the measurement

Optimized analysis of the stored measurement data

- Magnify function for representing selected measurement data
- Library support for Scan Bargraph measurements
- Simplified interpretation by means of the measure data documentation
- Integral and statistics functions
- Creation of ASCII files for data export (e.g. to EXCEL, etc.)

Auxiliary functions for Service

- Leak test with audible signal
- Total pressure monitoring
- Digital I/O manual mode
- Password protection with access privileges
- Calibration of the mass scale
- Easily understandable parameter setup
- Standard parameter editors for the different measurement types
- Editor for entering/modifying spectra in the spectra library
- Editor for characteristic curves of the analog inputs and outputs
- Graphic editor to process spectra
- Screen color editor

1.2.4**Quadstar 32-bit «Quantitative Analysis» features**

Quantitative Analysis enhances Quadstar 32-bit with quantitative analysis through Multiple Concentration Detection (MCD measurement).

The calculation of the concentrations is based on a matrix comprising max. 64 x 64 elements. The mathematical solution of the resulting equations is based on the process of elimination according to Givens, and also allows super-determinate equation systems.

The calibration factors are stored in a library and can be entered via the keyboard, or be determined through a special measurement method. The matrices access the individual elements in the calibration factor library. Recalibration of a component thus updates all matrices containing this component.

The main features of the Quadstar 32-bit Quantitative Analysis are:

- Analysis of a gas mixture by Multiple Concentration Detection (MCD measurement) with simultaneous display of concentrations, ion currents and external analog signals. Furthermore threshold monitoring, integration, auto ranging, etc.
- Calibration with internal and external normalization

1.2.5 Quadstar 32-bit «Process Control»

Using the sequencer that's supplied with Quadstar 32-bit you can program measurement sequences and control gas inlet systems automatically. Sophisticated, application-specific measurement programs can be created with the editor. Data may be exchanged online with other programs by DDE (Dynamic Data Exchange).

The sequencer also supports online calculations. By entering mathematical formulas you can calculate ion current ratios, standard deviations etc. from the measurement data that have been stored in a buffer during the execution of a sequence. The Calculation Module supports not only basic arithmetic functions, but also functions (Log(), Exp(), Sqr(), Abs(), etc.).

With the Process Control Module, the following main functions can be programmed in any order:

- MID (Multiple Ion Detection), Scan Bargraph and Scan Analog measurements
- Storage of measure data and their documentation
- Programmable calculations with measure data
- Individual control of a process sequence through programmable dialog boxes
- Integrated gas inlet control
- Nested operation loops with conditional branching
- Wait loops
- Zero gas measurement for background analysis (Zero Gas)
- Selection of the display mode (Table, Bargraph, Versus Time)
- Measurement log output on the printer
- Concentration measurements at gas mixtures
- Programmable calculations with concentrations
- Gas specific calibration

1.2.6 Scope of delivery

Software

Quadstar 32-bit comprises various compressed files that are supplied on CD-ROM. When the software is installed with the supplied setup program, the required directories are automatically created and the decompressed files stored in them.

When the software is installed using the given default settings, on drive C the following directories are created:

<code>c:\qs32bit\</code>	Program and setup directory
<code>c:\qs32bit\auxil\</code>	Auxiliary file directory
<code>c:\qs32bit\dat\</code>	Measurement data file directory
<code>c:\qs32bit\par\</code>	Measurement parameter file directory

For each measurement type, parameters and measurement data files are supplied. These are intended as sample files for the measurement and analysis programs. The following demo files are available:

<i>airdemo.xxx</i>	Parameter/measured data files for Simulation mode
<i>airfar.xxx</i>	Parameter/measured data files for Faraday detector
<i>airsem.xxx</i>	Parameter/measured data files for SEM/CH-TRON detector

Quadstar 32-bit distinguishes the different file types by their file extension. An overview of those files is given in Chapter 10.3 Quadstar 32-bit files, 544.

Cables/ connectors

- RS232C connection cable for the data exchange between the mass spectrometer and the computer
- Additional accessories (ArcNet board, optic fiber etc.) if ordered

1.2.7 Computer system requirements

Quadstar 32-bit can be run on any system that satisfies the following requirements:

CPU	133 MHz or higher Pentium-compatible CPU
Disk drive	2GB hard disk with a minimum of 650MB of free space; CD-ROM drive
Monitor	800 x 600 or higher resolution, 1024x768 recommended
RAM	64 MB minimum; 256MB recommended
Interfaces	RS232C (Com1-Com9 supported); PCI or PCMCIA for Arcnet NIC
Operating System	Windows® 2000 or later; Windows® 2000 professional, english recommended

1.3 Abbreviations and terms

1.3.1 Hardware modules and components

Term	Significance
AI 421	Analog input module, 16 channels for QMS 422/421
AO/IC 421	Ion counter including 12 analog output channels for QMS 422/421
CD-SEM	Secondary electron multiplier with high voltage conversion dynode
CH-TRON	Channeltron multiplier
CLA 400	Chopper lock-in amplifier for QMS/QMI 421
CP 400	Ion counter preamplifier QMS 422/421
CS 421	Front panel of the QMS 421
DI 420	Digital input module, 32 channels for QMS 422/421
DO 420	Digital output module, 32 channels QMS 422/421
EP 422	Electrometer preamplifier for QMS 422
EP 112	Electrometer preamplifier for QMS 422/421
FIR	Finite Impulse Response filter
HV 420/421	High voltage supply for the SEM of the QMS/QMI 421
HV 200	High voltage supply for the CH-TRON of the QMS
IS 420	Ion source supply for QMS 422/421
OH 421	Optical HUB for multiplex operation using ArcNet
PE 420	Inverted magnetron (cold cathode) gauge module for QMS 422/421
PI 420	Pirani gauge module for QMS 422/421
PKR 250/251	Compact full range gauge for QMS 200
QC 421	Quadrupole controller of the QMS/QMI 421
QC	Generic term for quadrupole controller
QMA 4XX	Analyzer of the 400 family
QMA 200	Analyzer of the 200 family
QMA 125	Analyzer of the 125 family
QME 125	Quadrupole electronics for the 125 family
QMH 400/410	RF generator of the 400 family
QMS 422	Quadrupole mass spectrometer 422
QMS 421	Quadrupole mass spectrometer 421

Term	Significance
QMS 200	Quadrupole mass spectrometer 200 (Prisma™)
QMS	Generic term for quadrupole mass spectrometers, applies to QMS 422/421, QMI 421, QMS 421C, QMS 200, etc
QMU 112	Gauge head selector
RF 20x	RF generator of the QMS 200, mass range x00
SEM	Secondary electron multiplier
TPR 250	Compact Pirani gauge for QMS 200

1.3.2 Quadstar 32-bit specific

Term	Significance
GSC	Gas Specific Calibration
IVT	Intensity Versus Time
MCD	Multiple Concentration Detection
MID	Multiple Ion Detection
MSC	Mass Scale Calibration
PCM	Process Control Module
QAM	Quantitative Analysis Module
SIM	Simulation Mode
TIC	Total Ion Current

1.3.3 General

Term	Significance
ArcNet	LAN with «Token Passing» protocol
ASCII	American Standard Code for Information Interchange
Baudrate	Data signaling rate in words/sec
F.S.R.	Full Scale Range. The highest displayable value of a scale
Hardware	Physical equipment/modules such as QMS 422, QME 125, etc.
LAN	Local Area Network
LED	Light Emitting Diode
Node	Every device in the network that has its own address
RAM	Random access memory (read/write)

Term	Significance
ROM	Read-only memory
RS232C	Serial standard interface
SECS	Semi Equipment Communications Standards
Windows®	Microsoft® Windows® graphic user interface

1.3.4 Conventions

Select [xxx] > [yyy] > [zzz] means:

In the current menu select the menu item [xxx]. In the subsequently appearing field choose [yyy]. And there choose the entry [zzz].

1.4 Installation

1.4.1 Using the SETUP installation program

Quadstar 32-bit comprises various compressed files that are supplied on CDROM and requires approximately 10 MBytes of hard disk space. SETUP creates the required directories automatically and stores the decompressed files there.

The main directory (default: QS32BIT) contains all Quadstar 32-bit program files. Special subdirectories are created for the measurement data, parameter and auxiliary files.

NOTE:

The subdirectory names can later be changed, but note that the directory names must be entered under [Setup] > [General] in the Parameter Setup program.

When an update installation is performed, SETUP searches for the default directory QS32BIT on all disks. If this directory can be found, it is shown together with the subdirectories in the corresponding dialog boxes and does not have to be re-entered.

SETUP supports:

- New installation of Quadstar 32-bit
- Update installation of an old version to a new version of Quadstar 32-bit

1.4.2 Installation/update on a hard disk

Install Quadstar 32-bit by performing the steps outlined below.

- 1 Start Windows®. If Windows® is already running, close all other applications.
- 2 Insert the Quadstar 32-bit installation CD-ROM
 - 2.1 Installation should start automatically after a few seconds
 - 2.2 If it does not start up after 30 seconds, the autostart feature may be disabled. In this case, use Windows® Explorer® to navigate to the CD-ROM drive and start "SETUP.EXE".
- 3 Select the desired language (e.g. English) by clicking on the corresponding command button. The online documentation will be installed in the selected language.
- 4 Enter the main directory for Quadstar 32-bit. If in the given directory an older Quadstar 32-bit installation is found, then enter the desired installation type (New or Update)

NOTE:

If you have an older installation of Quadstar 32-bit and you want to leave it untouched, please see Section «Multiple installations», 27.

- 5 Continue with the installation procedure and enter the requested parameters or accept the given standard entries

The SETUP program creates a Quadstar 32-bit entry in the "programs" folder in the start menu (if it does not exist yet):

Nonstandard directory

Setup installs Quadstar 32-bit in the desired directory. If this is not the default directory «QS32BIT», sequences supplied for the Process Control Module must be recompiled:

- 1 Start the Parameter Setup program in Quadstar 32-bit program group
- 2 Open [Sequence] > [Editor], choose a sequence and confirm by [OK]
- 3 Open [File] > [Compile All] to activate the Sequence Compiler dialog box and recompile all sequences by [Compile all]
 - After the recompilation, Quadstar 32-bit reports that all sequences have been compiled successfully
- 4 Click on [Close] and leave the Parameter Setup program by [File] > [Close] and finally [File] > [Exit]

Multiple installations

For multiple installations, the characteristics of the entry in start menu must be changed before the second installation.

- 1 Right-click the start menu entry (Quadstar 32-bit)
- 2 Chose [rename]
- 3 Enter a a new name for the entry, for example "Quadstar 32-bit #1"
- 4 Use another directory than the standard (qs32bit) for the second installation

Software key

Quadstar 32-bit is protected against unauthorized use by a key which contains your individual customer number. This key is stored in the software.

To order expansion modules or software updates, please use the order card at the beginning of this manual and do not forget to specify your key number on it. The key number is printed on the CD-ROM.

1.4.3 Connecting the mass spectrometer to the computer

Connection cables

Quadstar 32-bit is supplied with a set of connection cables for serial operation. Which cables are to be used depends on the type of connector on the computer.

NOTE:

If you use ArcNet/fiber optics, you can skip this section. Installation support for ArcNet can be found in Chapter 3.6 Communication parameters, 209.

Connection sockets

For connecting the serial communication cable to the mass spectrometer, two types of sockets are normally available, depending on the type of interface being used. The pin assignment of these two connectors (DB9 male, DB25 male) is as follows:

Pin	Significance
1	Carrier Detect (CD)
2	Receive Data (RXD)
3	Transmit Data (TXD)
4	Data Terminal Ready (DTR)
5	Signal Ground (GND)
6	Data Set Ready (DSR)
7	Request To Send (RTS)
8	Clear To Send (CTS)
9	Ring Indicator (RI)

Tab. 1-1

DB9 connector (male)

Pin	Significance
1	Chassis ground
2	Transmit Data (TXD)
3	Receive Data (RXD)
4	Request To Send (RTS)
5	Clear To Send (CTS)
6	Data Set Ready (DSR)
7	Signal Ground (GND)
8	Carrier Detect (CD)

Pin	Significance
20	Data Terminal Ready (DTR)
22	Ring indicator (RI)

Tab. 1-2

DB25 connector (male)

RS232C null modem cable (DB9 female to DB9 female)

DB9 connector (female)	1	2	3	4	5	6/8	7	9
DB9 connector (female)	7	3	2	6/8	9	4	1	5

Extension cable (DB9 male to DB9 female)

DB9 connector (female)	1	2	3	4	5	6	7	8	9
DB9 connector (male)	1	2	3	4	5	6	7	8	9

Installing the connection cable

The following combinations are feasible:

Computer connector	Connection	Connector on QC
DB9connector (male)	RS232C cable	DB9 connector (male)
DB9 connector (male)	RS232C cable + extension cable	DB9 connector (male)

Connect the QMS with the supplied cable to your computer and note the selected interface COM1, COM2, etc.

If you use an additional extension, pins 2 and 3 (DB-25 female) of this cable must not be crossed. If you have problems with an extension cable, a reverser may have to be used.

If several interfaces are used under Windows®, problems can occur. Their solution is described in Chapter 3.6.2 Communication via RS232C, ¶ 210.

For additional information please see Chapter 3.6 Communication parameters, ¶ 209.

1.4.4 Configuring Quadstar 32-bit

If you have just installed Quadstar 32-bit, the following must be configured:

- Communication from PC to mass spectrometer
- The mass spectrometer’s characteristics

If the system has already been configured, you may skip these two sections.

1.4.4.1 Communication with the mass spectrometer

- 1 Start the Parset program in the program group Quadstar 32-bit
- 2 Open [Comm] > [QMS Comm. Type]




Fig. 1-1

Quadstar 32-bit shows the following dialog box:




Fig. 1-2

At the time of installation Quadstar 32-bit assigns the port for communication with the mass spectrometer to the serial (RS232C Normal) interface COM 1 by default. If your mouse is connected to COM 1, the interface to the mass spectrometer must be assigned to COM 2.

NOTE:

If more than 2 serial interfaces are operated concurrently, problems can occur.

If a «PS/2 mouse» or «USB mouse» is used, this problems do not occur if the serial lines are not already occupied by other applications.

For additional information about communication (e.g. by ArcNet fiber optic) please refer to Chapter 3.6 Communication parameters, 209.

If you want to communicate via RS232C Normal, please proceed as follows:

- 3 Select RS232C Normal (if not selected already)
- 4 Choose [Comm] > [QMS Comm. Parameters]. Quadstar 32-bit shows the following dialog box:




Fig. 1-3

- 5 Select the desired Baud Rate (Default = 19200)
- 6 In the «Connection Line» list box select the desired interface and press [OK]
- 7 Start the mass spectrometer and wait approximately 5 seconds until the equipment is ready
- 8 Choose [Comm] > [Connect]

The QMS is automatically adapted to the selected baud rate. When the system is ready for operation, the LED simulated in the lower left corner of the active window changes to green.

Additional information can be found in the description of the Parameter setup program under Communication via RS232C.

Error Message

If the connection cannot be established, this is signaled with the error message «QMS off, not available or wrong node ID configured». Check the following possible sources of trouble:

- Is the mass spectrometer switched on and was it initialized when [Connect] was activated?
- Is the connection cable plugged into the mass spectrometer and into the computer?
- If another cable is used, check whether or not the TXD and RXD lines (pins 2 and 3 of the DB25 female connector) are crossed
- If an extension cable is used, make sure that the TXD and RXD lines (pins 2 and 3 of the DB-25 socket) are not crossed again
- Is the interface specified under QMS Comm. Parameters the one that's actually used?

1.4.4.2 Configuring the mass spectrometer

Declare the hardware you are using (analyzer, ion source, extensions etc.):

- 1 Start the Parameter setup program and choose [Comm] > [Connect]
- 2 Select [Config] > [QMS]. Quadstar 32-bit displays one of the following dialog boxes, according to the QMS in use:

Fig. 1-4

Fig. 1-5

It is recommended to connect the mass spectrometer and to activate the interface via [Comm] > [Connect] to setup the configuration, although it's not mandatory. This way, Quadstar 32-bit is able to read back part of the configuration from the mass spectrometer. For the QMS 422/421 this is particularly useful with respect to the optional Hardware Boards (AI 421, DO 420, DI 420, etc.) and for the QMS 200 with respect to the available analyzers, mass ranges (RF generators) etc.

- 3 Please enter the following data:
 - MASS-R - The mass range of your analyzer
 - QMA - The type of analyzer being used
 - Detector - The type of collector being used
 - IS-Type - The type of ion source being used
 - Options - The installed options
 - Protection (SEM+FIL) - Select INTERN if you are not using an external vacuum control system. See Chapter 5.2.1.2.5 Filament Protection Extern (QMA 200/QMA 4XX), 323.

Confirm these settings by [OK], so they will be stored automatically

1.4.4.3 Changing the directories

The names of the QS32BIT directories created by Quadstar 32-bit can be changed after the software has been installed.

To declare e.g. a different data file directory, proceed as follows:

- 1 Start the Parameter Setup program
- 2 Select [Setup] > [General...]
- 3 Position the cursor at the «Data File Path» input field and change the directory (and/or drive) by entering the name directly in this field, or via the [Select Path] command button.

NOTE:

The directory or path names are to be changed by Quadstar 32-bit only (as described above). If you change them e.g. with Windows® Explorer®, Quadstar 32-bit can't find its files any more and thus works with standard setup.

1.5 Getting started

1.5.1 Installation

Make sure that Quadstar 32-bit and the mass spectrometer (QMS) are installed and set up properly. See Chapter 1.4 Installation, 26.

1.5.2 Choose operating mode «Simulation»

For the following exercises we choose the operating mode «Simulation». This way we can work independently from the real vacuum conditions, because the QMS delivers simulated measure data.

- 1 Open the parameter setup program (Parset) and choose [Setup] > [General...]
- 2 In the now appearing dialog box set the «Simulation Mode» to «Intern» and confirm by pressing the [OK] button




Fig. 1-6

- 3 Activate the communication to the mass spectrometer by [Comm] > [Connect]




Fig. 1-7

Now you see on the left of the status bar of Parset the green LED (green = communication is active), then a field for text display and on the right the symbol for the emission (filament gray = emission is off) and the abbreviation SIM (for Simulation Mode). If you work with ArcNet communication type, then the name of the connected mass spectrometer is displayed in the last field on the right.



Fig. 1-8

- 4 Leave Parset by [File] > [Exit]

1.6 Air analysis

The measuring task is to analyze the composition of air. For that purpose we do not use a real measurement under vacuum conditions but the simulation data stored in the QMS instead. These data represent a typical measurement of air with a Faraday detector. The procedures are shown for a Prisma™ mass spectrometer, but of course they are nearly the same with other mass spectrometers.

1.6.1 Scan Analog overview measurement

To get an overview of the gases being present, the measuring type Scan Analog (ion current as a function of m/e) is suitable.

- 1 Start the program Measure and choose [Scan] > [Analog]. The following Open dialog appears:




Fig. 1-9

- 2 Click on «sca1f.sap» (DO NOT double click!)
- 3 Now click [Copy] to copy this file to another one
- 4 In the «Copy Files» dialog, in the field «To», change the name to *myairsim* and confirm with [OK]




Fig. 1-10

- 5 Confirm in the Open dialog the choice of *myairsim.sap* by [OK]. The measurement is started and shows, after a while, the following:




Fig. 1-11

You can clearly see the characteristic peaks of air on the mass numbers 14, 16, 18, 28, 32, 40 and 44 m/e.

1.6.2 Save measure data

Now, store some measuring cycles (passes):

- 1 Choose [File] > [Save Cycle Data] and enter the file name *airsim1.sac* under «File Name». The measured data will be stored in the file *airsim1.sac* in the path *c:\qs32bit\dat*.




Fig. 1-12

- 2 If you press [OK], you will be asked for the number of cycles to be stored. Enter 5 and press [OK].




Fig. 1-13

You can enter some information on the measurement here. Press "Save measure data + Info" to start the measurement.




Fig. 1-14

The QMS now performs the same measurement five times and stores the results. This takes about 4 minutes. In the text field of the status bar the measure/save progress is displayed continuously («Saving Cycle ... of ...»). When the text disappears, the storage is finished, but the measurement continues. Now you can stop the measurement by [File] > [Close] and leave Measure by [File] > [Exit].

1.6.3 Analyzing saved measure data

You can analyze stored measure data very easily and thoroughly with the program Dispsav (Display saved values).

- 1 Open the Dispsav program and choose [Scan] > [Analog Data]. In the now appearing Open dialog choose the file *airsim1.sac* and confirm the choice by [OK]. Now you see the following display of measure data:




Fig. 1-15

If you set the cursor to a certain spot on the graph, you see in the status bar the corresponding X- and Y-values and the information C:1 (displayed cycle: 1). This representation is not ideal because the peaks are small and their intensities can't be determined accurately. Therefore we enlarge the area we are interested in:

- 2 Choose [Function] > [Magnify]. Set the cursor to the upper left corner of the area to be enlarged (e.g. 10 amu, 0.85 E-8 A) and press the left mouse button. Drag the cursor (while still holding the left mouse button pressed) to the lower right corner of the area (e.g. 35 amu, 0.00 E-7 A) and press the right mouse button additionally. Now release both mouse buttons.




Fig. 1-16

The peaks can now be seen much better and their intensities may be determined more precisely. You can use the magnifier function repeatedly to enlarge certain details to a great extent.

- The last magnification step can be undone by [Special] > [Undo Magnify]
- To get back to the original representation, choose [Special] > [Set Default]

Until now, you have only worked with the first one of the five stored cycles. To see the changes in the spectrum from cycle to cycle, you can display the other cycles as well:

- 3 Choose [Setup] < [3 Dimensional], so the cycles will be displayed one after each other in the same diagram, i.e. 3 Dimensional:




Fig. 1-17

Of course, in our simulated measurement the values do not drift over the time. In a real measurement you could see such effects in this display mode. Try the settings and functions of Dispsav (e.g. [Setup] > [Graphic]) to become familiar with the features of the analyzing tool.

- If you should encounter a stage where the display is not sensible any more, you get back to the original representation by [Special] > [Set Default]
- Exit the program by [File] > [Close] and finally [File] > [Exit]

1.7 Measuring time dependent processes

The Scan Analog measurement gives you information like «The gas consists of ..., ..., ..., most of it is ... The amount of ... is very small». This measuring type has some disadvantages thus. The quantitative composition of the gas in % and ppm can only be figured out by a lot of hard, manual calculation work. The information about the changes of the gas composition is even in 3 dimensional representation very hard to extract and after all the amount of measure data gets very large when long term measurements are performed.

The solution for this problems is measuring with MID (Multiple Ion Detection) or MCD (Multiple Concentration Detection) measurement type. This way, only the mass numbers of interest and their peak heights are determined and analyzed. For the MID measurement, those intensities are already the result, while in MCD measurement type the intensities are used to calculate concentrations in %, ppm etc.

1.7.1 Basic measuring techniques

1.7.1.1 The MID measurement

First let's have a look at the MID measurement. In this measuring type you declare the m/e ratios for which the intensities (ion currents) are to be measured. Each of these single measurements is represented by a separate channel of the QMS, on which all parameters can be set.

Example: The m/e ratios (we call them «masses» from now on) 14, 16, 28 and 32 are to be measured.

- 1 Start the Parset program
- 2 Choose [Measure] > [MID]
- 3 Enter *airsim2.mip* in the file name field of the open dialog and confirm by [OK]
- 4 Quadstar 32-bit reports that this file does not exist (yet) and asks whether to create it or not. Confirm the creation by [Yes].
- 5 The channel parameter editor appears:




Fig. 1-18

- 6 Switch on the channels 0...3, according to the four masses of which we want to measure the intensities
- 7 Click in column CH-0, «State» row, on the field with the text OFF with the right mouse button
- 8 In the now appearing list that contains all possible values for this parameter (channel state: ENABLE/SKIP/OFF meaning active/currently inactive/off), click with the left mouse button on ENABLE
- 9 Switch on the channels 1, 2 and 3 in the same way. Now we have four active channels we can assign the four masses to.
- 10 Click with the right mouse button on the field CH-0, «Mass» row
- 11 Change the value from 0.00 to 14.00 and click with the left mouse button on the gray button on the left of the new value
- 12 Enter the masses 16, 28 and 32 in the other three channels in the same way
- 13 Click at the bottom of the editor with the left mouse button on the register card «Display». The editor should now display the following:




Fig. 1-19

The other parameters have already been set to reasonable values by the basic initialization and need not to be changed for the planned measuring task.

- 14 Close the file and save the changes by [OK]
- 15 Close the Parset program in order to get back to the Windows® Desktop
- 16 Now start the Measure program and choose [MID] > [Versus Time]. Perform the measurement with the parameter file *airsim2.mip*
- 17 Now you should modify some settings under [Parameters] > [Setup]
 - 17.1 Y-Scale: Set to «log» in order to see small signals as well
 - 17.2 Cycles: Enter 10, in order to get visible results quickly
 - 17.3 Leave the rest of the settings in their standard conditions
 - 17.4 Confirm the changes by [OK]. After a few seconds, the following will be displayed:




Fig. 1-20

You see the intensities of the masses versus the measuring cycles. At the right border you can see what mass is represented by which color.

- 18 Now click on the mass number 14.00 with the left mouse button and read its intensity in the status bar. Read the intensity of the mass 32.00 as well. The intensities are almost equal, therefore the lines coincide. The displayed intensities are the ones of the last cycle.
- 19 Try the settings of [Parameters] > [Setup] and [Parameters] > [Channel]. You can still change the channel parameters under [Parameters] > [Channel].
- 20 Then close the measure program by [File] > [Close] and eventually [File] > [Exit]

1.7.1.2 How the MCD measurement works

The concentrations of the gases are determined by multiplying every intensity measured on a certain mass with a calibration factor (related to the gas to be detected) and solving the resulting equation system.

In other words: An MID measurement is performed and the parts of the gas are calculated by multiplying the intensities with previously determined factors. Those calibration

factors can not be declared universally, but they may be found and stored in a library by performing a GSC (Gas Specific Calibration) measurement using a calibration gas. The solution method provided with Quadstar 32-bit allows determination of gases as well, that occur on identical or partly identical masses in the spectrum.

1.7.1.3 Preparatory measurements

Because the gas components are always represented by integer mass numbers, the easiest way would be to measure the intensities on integer masses. The measure equipment always causes small deviations between the nominal mass and the real peak maximum. Therefore, the (correct) maxima of the peaks are normally located a bit beside the theoretical number (e.g. 13.95 instead of 14.00). To perform the measurement with the highest possible precision, we use the measure type MSC (Mass Scale Calibration) to determine the exact positions of the peak maxima.

Mass spectrometer background

A possible cause for faulty measurements is the so-called «mass spectrometer background». It consists of gases that are present in the vacuum system which do not come from the sample to be measured. Because the background parts of a well-conditioned measuring equipment are nearly constant, they may easily be recorded by the Zero Gas measurement and eliminated in the actual analysis. The intensities recorded by the Zero Gas measurement can be subtracted from all subsequent MID/MCD measurements. This is selectable under [Parameters] > [Setup] > [Zero Gas Subtraction]. This way, the result consists only of the real parts of the analyzed gas.

1.7.2 Preparing the measurements

First we prepare all measurements, i.e. we create the required parameter files. Open the Parset program in order to do so.

1.7.2.1 Mass Scale Calibration

- 1 Choose [Calibration] > [Mass Scale]
- 2 Copy the file *adjustms.msp* to *myadjust.msp* by [Copy] and open it. This file is already prepared for measurements in simulation mode.

NOTE:

For real measurements the mass numbers must match the ones that really occur in the gas to be measured.

1.7.2.2 Gas Specific Calibration

- 1 Choose [Calibration] > [Gas Specific Sensitivity]
- 2 Copy the file *airdemo.gcp* to *mysense.gcp* by [Copy] and open it. This file is already prepared for measurements of air and does not require any modifications.

NOTE:

For real measurements, the gases that occur in the calibration gas bottles, their

concentrations and the occurring mass numbers must be entered. For every calibration gas mixture a new file must be created.

1.7.2.3 Analysis file

- 1 Choose [Measure] > [MCD]
- 2 Copy the file *mcdmedf.mcp* to *mymcd.mcp* by [Copy] and open it. This file is already prepared for measurements of air and does not require any modifications.

NOTE:

For real measurements, the gas components of the gas to be analyzed and the mass numbers to measure the intensities on must be entered.

- 3 Choose [Edit] > [Channel Par]. Make sure that the display of all measure channels is switched to OFF in the register card Display, «Disp. Mode» row. We do not want to display any ion currents but only the concentrations calculated from them by using the calibration factors
- 4 Finally you can define the range of concentrations to be displayed by [Edit] > [Concentration Par]. Choose Disp.F.S.R.=100% and DispDecades=4.
- 5 Then close all files and confirm by [OK] that they shall be saved
- 6 Eventually exit Parset

1.7.3 Measurements

Now perform the actual measurements. For that purpose start the Measure program. Let's begin with the auxiliary measurements.

1.7.3.1 Calibrating the mass scale

- 1 Choose [Calibration] > [Mass Scale] and open the file *myadjust.msp*. Quadstar 32-bit performs the mass scale calibration and displays the newly found peak maxima. If in later measurements mass numbers are used that have not been calibrated, their positions will be interpolated from the calibrated values.

Nbr	Old Mass	New Mass	Intensity	State
0	2.00	2.00	5.718E-10A	ok
1	18.00	17.97	1.231E-09A	ok
2	28.00	28.00	9.980E-09A	ok
3	32.00	32.00	8.298E-09A	ok
4	44.00	44.00	6.266E-10A	ok

Fig. 1-21

- 2 Exit the window by [File] > [Close]
- 3 To let the mass scale calibration take effect, in the MID or MCD measurement under [Parameters] > [Setup] the option «Mass Scale Correction» must be checked

1.7.3.2 Measure the zero gas

While you are working with simulated measure data, the Zero Gas measurement must not be performed. The simulated data are exactly the same for the analysis and the Zero Gas measurement. If you would perform a Zero Gas measurement here, then in a later analysis those values would be subtracted from the measured data and the results would be ion currents that are all zero.

NOTE:

In real measurements too, the Zero Gas measurement may only be done when the analysis chamber is on minimal pressure or a real zero gas has been let in. Otherwise severe measurement errors may occur.

You will get more information about the Zero Gas measurement in Chapter 2.8.2 Background measurement (Zero Gas), [137](#).

- 1 Exit the window by [Close]

1.7.3.3 Calibrate the mass spectrometer

- 1 Choose [Calibration] > [Gas Specific sensitivity] and open the file *mysense.mcp*

Quadstar 32-bit carries out the gas specific calibration and shows the resulting gas components, the ion currents and the calculated calibration factors. You will get more information about the meaning of the calibration factors in Chapter 2.8 Calibration, [134](#).

Nbr	Mass	Component	Concentration	Intensity	Calib.Factor
0	14.00	N2	78.1400%	8.304E-09A	8.320E-01
1	28.00	N2	78.1400%	9.981E-09A	* 1.000E+00
2	29.00	N2	78.1400%	3.925E-10A	3.932E-02
3	32.00	O2	20.9000%	8.303E-09A	3.110E+00
4	40.00	Ar	9300.0ppm	1.652E-09A	1.391E+01
5	44.00	CO2	300.00ppm	6.109E-10A	1.594E+02

Fig. 1-22

1.7.3.4 Measure and store concentrations

Now measure the composition of the gas. In order to store measuring conditions and comments additionally to the measure data, check the option «Save File Info».

- 1 Choose [Setup] > [File Info]
- 2 Check «Measurement Parameters» and «Save File Info»
- 3 Confirm the choice by [OK]




Fig. 1-23

- 4 Choose [MCD] > [Versus Time] and open *mymcd.mcp*. The concentrations calculated from the measure data are displayed versus the number of cycles. You will get more information about the calculation of concentrations in Chapter 2.6 MCD measurement type, 103.
- 5 Try the settings here as well and note the influence of parameter changes




Fig. 1-24

- 6 To store data, choose [File] > [Save Cycle Data]. Enter *myconc.mdc* as the file name and confirm by [OK].
- 7 In the following dialog enter 30 for the number of cycles to be saved and confirm by [OK].
- 8 Then the File Info dialog appears, where you can comment and document your measurements. Enter notes etc. and confirm by [Save Measure Data + Info]




Fig. 1-25

- The measurement is started and 30 cycles are stored. This takes about 3 minutes. When the measurement is done, i.e. no text is displayed any more in the text field of the status bar, you can close the measurement and exit Measure.

1.7.4 Analysis

1.7.4.1 Display modes

- Start Dispsav and choose [Process] > [Cycles]
- Load the file *myconc.mdc*. The following will be displayed:




Fig. 1-26

You see all measured ion currents and concentrations. At the upper left corner, just below the menu bar, you see a list box that provides change over of the axis. At the

moment, the axis is set to concentrations in % ('Concentr. [%]'), which causes the concentration names on the right to be marked by a short vertical line.

- 3 Change the axis to ion currents ('Ion Current [A]'), so the axis changes and the mass names are marked as valid
- 4 You may magnify these stored data as well
 - 4.1 Choose [Function] > [Magnify]
 - 4.2 Set the cursor to the upper left corner of the area to be enlarged (e.g. Cycle 0, 4.0[E+01%]) and press the left mouse button
 - 4.3 Drag the cursor to the lower right corner of the desired area (e.g. Cycle 25, 0.5[E+01%]) while still holding the left mouse button pressed and press the right mouse button additionally
 - 4.4 Now release both mouse buttons
- 5 From the stored data as well, you are only interested in the concentrations and therefore you want to hide the ion currents. Choose [Select] > [Blocks]. The following dialog appears:




Fig. 1-27

Under «Available Blocks» you see the existing blocks, i.e. concentrations and ion currents. Under «Selected Blocks» you see the selected ones, which are at the moment all available blocks.

- 6 Under «Selected Blocks» select the entry «Ion Current» and press [Remove]. Now, only the data block concentrations is selected.
- 7 Confirm this choice by [OK]. Only the concentrations are displayed now.

Apart from the display Versus Time shown above (Intensity (Concentration) Versus Time, IVT) there are also table and bargraph display modes.

- 8 Change the cursor function from magnification to selection by choosing [Function] > [Cursor]
- 9 Press the left mouse button and drag the cursor (which now appears as crosshairs) across the measure display. The exact x-and y-position is displayed in the status bar.
- 10 Now drag the cursor to Cycle number 10 while still holding the left mouse button pressed. Press the right mouse button additionally and release both. You see the data of measure cycle 10 represented as a table. Individual measure values can be changed manually by [Change].




Fig. 1-28

- 11 Now choose [Bargraph], so the measure data of the chosen cycle are displayed as a bargraph




Fig. 1-29

- 12 Return to the display type «Intensity Versus Time» by choosing [Display] > [IVT]

1.7.4.2 Display additional data as a table

Quadstar 32-bit calculates from all displayed process data (MID, MCD, calculated values etc.) statistical and integral data that can be displayed as a table. If you have chosen a detail display by Magnify, the statistical data will be calculated only for the current detail area.

- [Display] > [Statistic] provides statistical analysis of the cycles. For every item minimum, maximum, mean value, absolute and relative standard deviation are calculated and displayed.
- [Display] > [Integral] shows the integrals of all items across the time

1.7.4.3 Print data

To get a hardcopy of the measured data there are two possibilities:

[File] > [Print]

Use this to print the displayed measure data as a table. In our case it doesn't make sense to print all data, so we only want the concentrations of the components Ar, N₂, O₂ and CO₂ of the cycles 1...5.

- 1 Choose [Setup] > [X Axis]




Fig. 1-30

- 2 Enter Start Cycle 1 and End Cycle 5
- 3 Set the Axis Type to «Time Absolute» to have the x-axis labeled for time
- 4 Then press [OK]
- 5 Now choose [Select] > [Data] and check whether the chosen components are those you want to print
- 6 If so, then close the Select Data dialog and choose [File] > [Print]. If the printer is ready, the table will now be printed.

Key combination CTRL+SHIFT+F12

Use this to print the currently displayed window of Quadstar 32-bit or the complete screen. This way you can easily get screenshots of measurements, settings etc. To make such screenshots even richer, there is a graphic editor in Dispsav:

- 7 Choose [Function] > [Graphic Editor]. The following toolbar will appear:



Fig. 1-31

- 8 Click on the tool [T] to enter a text. Then click on the spot in the display, where the text shall begin. In the Text Editor dialog enter the text (e.g. «Simulated Air Measurement») and confirm by [OK].
- 9 Now we draw a frame around the title: Click on the tool [FRAME] and set the cursor to the upper left corner of the frame to be drawn. Drag the cursor to the lower right corner of the desired frame while holding the left mouse button pressed. Release left mouse button. Your measure image should now display the following:




Fig. 1-32

- 10 Exit the graphic editor with the tool [EXIT] and confirm that the changes are to be saved
- 11 You should not change the window size of the picture to be printed any more, because the size and position of the title wouldn't be matched to the new window size
- 12 By [Display] > [Graphic Data] you can switch the display of the just created additional graphic data on or off
- 13 Make sure the printer is ready. Then press the key combination CTRL+SHIFT+F12, so the Print Screen dialog appears.
- 14 Choose the active window as the printing region and 1.0 as the zoom factor
- 15 Confirm the settings by [OK]




Fig. 1-33

The picture is now printed in the same manner as it is displayed on the screen.

1.7.5 Auxiliary tools

1.7.5.1 Documentation

With the measure data you have stored a documentation consisting of measure parameters and notes

- 1 You can recall this documentation by [Info] > [File]:




Fig. 1-34

- 2 Press [Measurement Parameters] to see the measurement parameters or [About File] to get the general file information. The entries of «Additional Information» and «Notes» can still be modified.
- 3 Exit the File Info dialog by [OK]. Use [Cancel], if you don't want to save changes.

1.7.5.2 Filter

If your measure data contain a lot of noise (of course, our simulation data contain very low noise), you can smoothen them now by using the filter of Dispsav.

- 1 Choose [Special] > [Filter...], so the following dialog appears:




Fig. 1-35

- 2 First get a new file by the [File Manager] that will hold the filtered data. Otherwise the original data would be overwritten.
- 3 Press [Start]. The raw data are filtered and stored to the new file.

1.7.5.3 Interpretation support

In Quadstar 32-bit there is a table of key fragments ([Help] > [Key Fragment Ions]), that supports the interpretation of spectra. See Chapter 1.9.1 Spectrum interpretation, 56.

1.8 Measuring real gas components

Now you have learned the basic procedures and possibilities of Quadstar 32-bit by using simulated data. If you change over to real measurements with the mass spectrometer, you must check the following:

- 1 Start Parset and choose [Comm] > [Connect]
- 2 Switch off the simulation mode by setting the Simulation Mode to Off. See Chapter 1.5.2 Choose operating mode «Simulation», 34.

NOTE:

Otherwise you will never get real measure data, but always the simulated ones that are completely unrelated to your vacuum system.

- 3 Switch on the emission, SEM if necessary. If you are sure that the pressure in the vacuum chamber is low enough, press the key combination CTRL + S.
- 4 Check Emission (and SEM, if you want to use the SEM). In the status bar of Parset the emission symbol (and SEM Symbol if necessary) changes its color from gray to yellow.

Now you should be able to perform your own measurements. You will find detailed instructions about that and about subjects that haven't been explained here (e.g. ion source tuning) in the following chapters or in the (equipment-)specific operating manuals.

1.9 Hints concerning the work with Quadstar 32-bit

1.9.1 Spectrum interpretation

Sometimes you may want to know what a peak on a certain mass means, if you didn't expect it e.g. There is a «Key Fragment Ions» table that shows mass number, key fragment, possible parent molecule and additional mass numbers. It may help you to find the desired information.

- 1 You can access the table through the main menu of any Quadstar 32-bit module by [Help] > [Key Fragment Ions] or by the key combination Ctrl + L

The example below shows the corresponding entry of the Dispsav module:



Fig. 1-36

- 2 First have a look at the main mass number (e.g. 12) to determine the key fragment
- 3 Then compare the additional mass numbers with the ones in your spectrum
- 4 Try to find the parent molecule that matches your spectrum most accurately (e.g. if there is an additional peak on mass number 13 and none on 28, the molecule is most likely to be C_xH_y and not CO or CO₂)

Key fragment ions			
<i>Mass Number (m/e)</i>	<i>Key fragments</i>	<i>Probable Parent Molecule(s)</i>	<i>Additional Mass Number (m/e)</i>
1	H+	H ₂	2
	H+	H ₂ O	18, 17, 16
	H+	C _x H _y	12, 13, 14, 26, 27 etc.
2	H ₂ ⁺	H ₂	1
	He ⁺⁺	He	4
4	He ⁺	He	2
6	C ⁺⁺	CO	12, 28, 29
	C ⁺⁺	CO ₂	12, 28, 44
	C ⁺⁺	C _x H _y	12, 13, 14, 26, 27 etc.
7	N ⁺⁺	N ₂	14, 28, 29
8	O ⁺⁺	O ₂	16, 32, 34
	O ⁺⁺	H ₂ O	16, 17, 18
12	C ⁺	CO	28, 29
	C ⁺	CO ₂	28, 29, 44
	C ⁺	C _x H _y	13, 14, 26, 27 etc.

Fig. 1-37

Of course there are much more sophisticated methods to analyze a spectrum, but in many cases the Key Fragment Ions table is sufficient.

1.9.2 Starting the Quadstar 32-bit programs

The programs are started from the Windows® start menu. Click the Windows® START button, select the “programs” folder, select the entry “Quadstar 32bit” and finally click to the corresponding icon of the program needed. Note that changes made in Parameter Setup program such as Communication, Setup Configuration etc. are not going to take effect in other programs until they are restarted.

1.9.2.1 Open a program with a certain file

If you intend to open a Quadstar 32-bit program always with a certain file, you may declare its name in the properties of the shortcut in the Windows® start menu.

Example: You want to see the measure data of *airdemo.sac* immediately, if you open the DISPSAV-program.

- 1 Right click the shortcut to “DISPSAV” and select “Properties”
- 2 Enter the name of the data file (*airdemo.sac*) in the field “Target” behind the program’s name (in this case *DispSav.exe*). Separate them with a single space.




Fig. 1-38

Of course, you can copy the shortcuts and give them different descriptions and files, e.g. in order to start different measurements directly.

For detailed information about managing the start menu and shortcuts please refer to the Windows® documentation or online-help.

1.9.2.2 Automatic start of a program

For automatically starting a measurement, sequence or data display when the PC is turned on, an icon that is modified as follows can be inserted into the *Startup* entry of the start menu.

Example: When the PC is booted, you want the measure data of *airdemo.sac* to be displayed automatically with the Dispsav program:




Fig. 1-39

- 1 Copy the icon «Dispsav» into the folder «Startup»
- 2 Change the description and the command line under «Properties» as shown above

When you turn on the PC and start the Windows® OS, the Dispsav program will be started automatically, showing the measure data of *airdemo.sac*.

For detailed information about the startup folder, please refer to the Windows® documentation or online-help.

For additional information refer to the Windows® Program Manager description under «Starting Applications from the Startup Group».

1.9.3 Transferring screens via the clipboard

Windows® uses a simple method for transferring screens via the clipboard. The clipboard is a temporary memory that is always available as long as Windows® is active. The screens (dialog boxes, etc.) can be copied from Quadstar 32-bit, pasted into the clipboard, and then transferred from there into other application programs.

1.9.3.1 Copying an entire screen

The contents of an entire screen can be copied into the clipboard at any time:

- 1 Make sure that the required information is shown on the screen
- 2 Press the PRINT SCREEN key

A «snapshot» of the screen is transferred to the clipboard.

1.9.3.2 Copying an active window

You may as well copy the active window only to the clipboard:

- 1 Make sure that the required information is shown on the screen
- 2 Press ALT + PRINT SCREEN key

A «snapshot» of the active window is transferred to the clipboard.

1.9.4 Inserting screens from the clipboard

The screens in the clipboard can be inserted (using SHIFT + INSERT) at any time into other programs such as Windows® Word®, Paint®, etc.

Additional information can be found in the Windows® help files.

1.9.5 Save or print Quadstar 32-bit screens

With Quadstar 32-bit you can print or save the active window or the entire screen as a picture

- 1 Call the Save/Print Screen dialog box in any Quadstar 32-bit program by pressing the key combination CTRL + SHIFT + F12 (or CTRL + SHIFT + ALT + F2). The following dialog box appears:




Fig. 1-40

Control elements

Region

Selection of the region to be captured.

- Active Window: The active window will be captured
- Entire Screen: The entire screen will be captured

Bitmap to...

Device to save or print the captured picture.

- Printer: Print out the picture on the active printer
- File (.BMP): Save the picture as a file using the Windows® bitmap format

Printer Option

- Zoom Factor: Magnification factor for the picture to be printed

2. Measure

Measure is used for capturing or displaying and storing ion currents, analog input signals, ion counter signals, etc. Concentrations can be calculated from the raw data via calibration matrices. With the sequencer, automatic measurement and control sequences may be programmed.

When the Measure program is started, Quadstar 32-bit attempts to establish a connection to the quadrupole mass spectrometer based on the selected transmission mode (only for RS232C or ArcNet). If successful, the LED symbolically displayed in the lower left corner of the active window changes to green. If the QMS is not connected or not switched on, a corresponding error message is output.

After the fault has been remedied, the connection setup can be reinitiated by [Comm] > [Connect].

2.1 General information

2.1.1 Measuring with Quadstar 32-bit

2.1.1.1 Measurement types

In Quadstar 32-bit there are different types of measurements:

Scan Analog

Continuous measurement of the ion current as a function of the mass. Representation as a curve across the chosen mass range.

Scan Bargraph

Search for peaks across the chosen mass range. The found peaks are represented as vertical lines over their masses.

MID

Multiple Ion Detection. On one or several predefined masses the ion current is measured. Various representation types.

MCD

Multiple Concentration Detection. On one or several predefined masses the ion current is measured and concentrations are calculated with them. Various representation types.

2.1.1.2 Measurement methods

Quadstar 32-bit provides different measurement methods for all measurement types (Scan Analog, Scan Bargraph, MID and MCD). The measurement method may be defined in advance ([Scan] > [Setup] > [Analog], [MID] > [Setup] > [Table] etc.), or later within the measurement ([Parameters] > [Setup]).

Single

The measurement is performed once and then stopped. By [Operation] > [Go] or when a new measurement is loaded it can be restarted.

Repeat

The measurement is continuously repeated. The speed depends mainly on the Speed or Dwell-time specified in the channels. Repeat measures continuously without interruption.

Timer

The measurement is restarted in the intervals defined under Timer Step.

Dynamic

Is a special measurement mode and only available for MID/MCD measurements. See Chapter 2.1.1.4 Fast MID/MCD measurements, 62.

2.1.1.3 FIR Filter

In Scan Analog and Scan Bargraph mode, a FIR (Finite Impulse Response) filter can be activated by setting the mode of SCAN-N or PEAK-L to SCAN-F or PEAK-F in the channel parameter editor (in group «Mass»).

2.1.1.4 Fast MID/MCD measurements

In Multiple Ion Detection (MID) or Multiple Concentration Detection (MCD) mode, Quadstar 32-bit provides different possibilities of data acquisition ([MID/MCD] > [Setup] > [Versus Time]):

- In method Single, Repeat and Timer:
 - Measure and Display
 - Display after Measurement
- The method Dynamic




Fig. 2-1

The actual measurement method (Single, Repeat, Timer, Dynamic) is selected in the upper left corner in the group «Measurement». The display mode (Measure and Display, Display after Measurement) is selected in the group «Measure/Display Mode» in the lower right corner.

Measure and Display

If fast measurements are to be performed over longer periods of time, «Measure and Display» is the right mode to use. Cycle by cycle is measured, transmitted to the PC and displayed.

Fast measurement with concurrent saving of data

When the save function is activated (in MID or MCD mode, [File] > [Save Cycle Data]), the following dialog box is displayed:




Fig. 2-2

To achieve a higher speed, unmark the checkbox «Display» to switch off the display. The PC doesn't consume computing power for the display then and it stores at maximum speed instead. If you have a fast graphics cards the effect of «Display» is small. Of course, the channel parameters have to be set accurately (short Speed/Dwell times). Please note that higher speed means lower measuring accuracy.

NOTE:

Close all Windows® applications that are not needed at the moment.

If the currently stored measurement cycle should be displayed continuously in the status bar, activate «Display» and switch the display function for all channels in the channel parameter editor off (in group «Display»).

Display after Measurement

Very fast short time measurements can be achieved using the «Display after Measurement» mode. The measurement is performed with the chosen speed. The data is stored in the internal buffer of the QMS and when the job is done, the data is transferred to the PC. The disadvantage is that the number of buffered measurement data is limited by the capacity of the internal buffer in the control unit. Therefore it is possible that the desired number of cycles is not achievable.

NOTE:

Please take into account that the «Display after Measurement» mode is always faster than «Measure and Display».

Dynamic Measurement

This measurement method is a combination of the two above-mentioned. It provides best measure speed for slow or irregular communication lines. The measure data is stored in the buffer of the QMS and concurrently transferred from there to the PC as fast as possible. The transmission speed can be optimized by the Performance Test. See Chapter 2.5 MID measurement type, ¶ 87 and Chapter 2.6 MCD measurement type, ¶ 103. If the buffer overflows, the measurement is stopped, the buffer flushed and the measurement started again. In that case, Quadstar 32-bit reports «Measure data drop-out, measurement too fast», because there was a period of time when no measurement was on.

2.1.2 Choosing the measuring range

In Quadstar 32-bit, the measuring range can be adapted manually or automatically. Open the channel parameter editor and click on the register card «Amplifier». If you choose [Amp. Mode] > [FIX], then you can enter the measuring range under [Amp. Range]. By [Amp. Mode] > [AUTO] or [AUTO D] the QMS can be instructed to find the correct measuring range automatically:

Load-Ch.:01	CH-0	CH-1	CH-2
State	ENABLE	ENABLE	ENABLE
Det. Type	FARADAY	FARADAY	FARADAY
Mass	1.01	2.01	14.03
Amp. Mode	FIX	FIX	FIX
Amp. Range	E-05	E-05	E-05
Range - L	----	----	----
Pause-Cal.	0.1	0.1	0.1
Offset	ON	ON	ON

Fig. 2-3

AUTO

The measuring range is chosen automatically from all available ranges.

AUTO D

The measuring range is chosen automatically from all ranges between the highest (most insensitive) and Range-L.

2.1.3 Status bar

Information concerning the measurement in progress as well as the status of the emission, the SEM and the simulation are displayed on the status bar. In ArcNet mode the name of the QMS is also shown.








Fig. 2-4

Control elements

Symbol	Significance
	Symbol for emission. Orange light when active.
	Symbol for SEM operation. Orange light when active.
	Simulation ON SIM on/off. See Chapter 3.7.2 General setup parameters, 226
	Name of the selected QMS in operation via ArcNet. The name is entered in the Parset program under [Comm] > [Unit Identifications]. In order to make use of the QMS, it has to be activated under [Comm] > [Unit Identifications] by the checkbox «Enable». Within a sequence, you can switch over between the activated QMS units by the instruction SetUnit().

2.1.3.1 Message text output during the measurement

In the text field of the status bar messages are output. In manual mode the output is done automatically. In sequencer mode any text can be output using the sequencer command Message().

Example: The output of the message text «Saving cycle 1 of 31 to sca1310.sca» originates from a manually started Scan Analog measurement with data storage.



Fig. 2-5

2.1.4 Configuring the graphic display range

The graphic display range can be set up in channel parameter editor, register card «Display». Both of the parameters [Disp. F.S.R.] (Display Full Scale Range, upper limit of display) and [Disp. Decades] (number of decades to display) may be selected manually or automatically (AUTO).




Fig. 2-6

2.1.4.1 Upper limit of display (Disp. F.S.R.)

AUTO

Quadstar 32-bit chooses the upper display limit automatically according to the other channel parameters.

Manual

The channel parameter editor shows all possible upper limits. Click with the right mouse button on the field [Disp. F.S.R.]. The choice may be restricted by other channel parameters.

2.1.4.2 Number of decades to display (Disp. Decades)

AUTO

Quadstar 32-bit chooses the number of decades to display automatically according to the other channel parameters which determine the maximum number of displayable decades. Particularly the resolution of the measured signal is a limiting factor. Example: 10 bits of resolution equals a range of 0... 1024, i.e. three decades.

Manual

The channel parameter editor shows all possible numbers of decades to display. Click with the right mouse button on the field [Disp. Decades]. This choice may be restricted by other channel parameters as well.

If the parameters Disp. F.S.R. and Disp. Decades are not the same for all channels, the display range is selected in such a way that the measurement signals of all channels can be displayed in the same chart.

The displayed measure signal may be distorted if e.g. Amp.Range is not set adequately.

Example: Measured signal, 16 bits resolution, highest value 5E-14 A, Amp. Range 1E-9 A. By reducing the Disp. F.S.R. manually to 1E-13 A, the displayed signal looks stair-like. This effect is caused by the internal quantization of the mass spectrometer.

Therefore, make sure the manually selected parameter values match each other, or use the AUTO mode.

Example: Scan Analog (First Mass = 10, Width = 40, Speed = 1 s/amu) with the following settings:

- Amp Mode: Auto (Amplifier register card of the Channel Parameter Editor)
- Decades: 5
- Disp. F.S.R.: 1.0E-8




Fig. 2-7

2.1.4.3 Display range in sequencer operation

For displaying all possible values (e.g. values calculated by the Process() instruction), the range for sequences is selectable separately under [Sequence] > [Setup] > [Display Range]. These settings are valid if nothing else is specified by the SetScale() instruction.




Fig. 2-8

Control elements

Display Ranges

- Maximum: Largest decade to be displayed
- Decades: Number of decades to be displayed

2.1.4.4 Y-axis lettering (MID/MCD)




Fig. 2-9

The Y-axis lettering can be changed between display in % (concentrations), Amperes (for ion currents) and other units of measure (for example, for analog input measurement values, etc.)

Total pressure values

- QMS 422/421: The total pressure values are marked with [Pirani] and/or [Cold Cathode]
- QMS 200: The total pressure values are marked with [TPR 250] or [PKR 251]

Marking the affiliation

The affiliation with the Y-axis is marked by identifying the names of the components, the names of the AI channels, or the mass numbers of the ion currents (in the table) with a vertical line.

2.1.5 Documenting the measurement data (File Info)

To make measure data evident and reproducible, you can store a documentation with it in Quadstar 32-bit. This may include measurement parameters, notes, title etc. In the main menu under [Setup] > [File Info], check the option «Save File Info», so the following dialog appears when you store measure data:

Nbr	Name	Value
0	temperature	23°C
1	total pressure	2.8 E-8 mbar
2	humidity	59%

Fig. 2-10

Control elements

Parameters/Data to save

Here you can check which parameters and files you want to store together with the measurement data.

- Measurement Parameters: Measure and channel parameters
- Ion Source Parameters: Current settings of the ion source
- State Parameters: Status information
- Zero Gas Data: Zero gas data that was used for this measurement

Additional Information

In these fields you can enter up to 20 user-defined parameter sets (a name and a value each). The previously used entries appear here by default.

- [Delete All Names]: Delete all user-defined parameter names
- [Delete All Values]: Delete all user-defined parameter values

Notes

Notes related to the current measurement (up to 1000 characters of text).

Title

Enter a title for the file (shows up in the file manager).

[Save Measure Data + Info]

Save measurement data and documentation.

[Save Measure Data Only]

Save only the measurement data (without documentation).

[Cancel]

Cancel data storage.

If data are appended (Append) to an existing data file, but these data have been obtained with a different parameter setup than the original data, the documentation is always related to the first measured data.


NOTE:

In measure mode (menu [File] > [File Info Notes]) or when recalling data (Display Saved Values under [Info] > [File]) you can add additional notes to report special events or make clear under which conditions the individual sections of the measurement data file have been created.

2.1.6 Multiplex operation with several control units

If Quadstar 32-bit is run with the optical HUB via the ArcNet-interface, simultaneous operation of several control units type QMS 422/421 or QMS 200 is supported.

The configuration is defined in the Parameter Setup program. For each control unit declared in the network, a separate path is set up in which the equipment-specific parameters and files are stored. After they have been configured, the individual control units can be logged on or off the network.

For detailed information on configuring, activating the individual control units etc. see Chapter 3.6.2.4 Communication via ArcNet,  211.

2.1.7 Calibrating the measurement equipment

For simple overview measurements with little accuracy, no calibration is necessary. If you want to carry out high precision measurements, Quadstar 32-bit supplies methods (depending on the QMS) to compensate different effects that otherwise would cause deviations. The following shows an overview of these effects and how they are handled in Quadstar 32-bit. The corresponding steps and further information are found in Chapter 2.8 Calibration, [134](#).

Mass scale

The nominal mass scale doesn't match the real masses exactly. These deviations are mostly small, nonlinear and unsystematic. By means of selected known masses, those deviations are measured and corrected. See Chapter 2.8.1 Mass scale calibration (Mass Scale), [134](#).

Background

In the vacuum chamber there's some pollution, also referred to as residual gas background. Those gas contributions can be subtracted from the measured data so that only the contributions from the gas to be measured remain. See Chapter 2.8.2 Background measurement (Zero Gas), [137](#).

Measure amplifier

The amplifier in the QMS that amplifies the (very small) ion currents to be measured causes different offsets in the different measure ranges. These offsets can be measured and subtracted from the actual measure data (only QMS 422/200). Since the offsets drift over longer periods of time, this calibration should be repeated periodically. See Chapter 2.8.3 QMS Offset, [140](#).

Detector

Besides the ions to be measured, the detector is hit as well by accidentally appearing ions. This causes another kind of offset which can be measured on a mass where there is no peak (channel parameter «Zero Mass», default: 5.5 amu) and then be subtracted (only QMS 422/421, in measure modes MID/MCD). This correction is performed for each channel before the measurement, if «Zero Mass» is not set to OFF. While the correction is in progress, the message «Channel x: Zero Equalize in progress» is displayed in the status bar. See Chapter 2.5.3 Offset alignment, [90](#) and Chapter 2.6.4 Offset alignment, [108](#).

Sensitivity

For concentration measurements (MCD mode), the gas specific sensitivity of the QMS must be determined. For that purpose, the relative sensitivity of the mass spectrometer is determined on several masses by using a special calibration gas. See Chapter 2.8.4 Gas specific sensitivity, [142](#).

2.1.8 Measurement of analog input signals

In MID and MCD mode Quadstar 32-bit supports the capture of measurement data via analog inputs. The number of channels depends on the available hardware:

- QMS 422/421: 16 analog input channels of the AI 421 board, which can be installed as an option in the QMS 421/422, are available. Measurement range: -10.24...+10.24 volts.
- QMS 200: The QMS 200 has two analog inputs as a standard feature. The number of inputs cannot be expanded. Measurement range -5.12...+5.12 volts.

In the channel parameter editor the analog input channels can be assigned to the individual QMS channels under «Det. Type» in the «Main Param.» group. The selected analog input channel is subsequently used as a signal source for the QMS channel.

The characteristic curve (AI curve) of the connected sensor can be assigned in the Parameter Setup program to each AI channel that is used.

Default values of the AI characteristic curve

Quadstar 32-bit is shipped with all available channels (QMS 422/421: 0... 15, QMS 200: 0, 1) programmed with a linear characteristic for positive values in millivolts. Thus, an input voltage of 0.100 V is interpreted as 100, 5.000 V means 5000 etc.




Fig. 2-11




Fig. 2-12

For additional information on analog inputs please see Chapter 3.8.3 Analog input characteristic curve configuration, 238.

2.1.9 Total pressure monitoring

In MID and MCD mode (manual and sequencer operation), Quadstar 32-bit also supports the measurement of the total pressure if the corresponding hardware options (PI 420 and/or PE 420) are installed in the QMS 422/421, or if a TPR 250 or a PKR 250/251 is connected to the QMS 200.

The total pressure gauge heads can be selected as the signal source in the Detector group of the channel parameter editor. The measured data can be displayed and stored.

Turning the PE 420/TPR 250 option ON/OFF

Whereas the PI 420/TPR 250 hardware is always active, the PE 420/PKR 250 gauge head must be turned on and off. Turning on and off is done by Quadstar 32-bit, but the procedure is different for operation with or without sequencer. The following should be noted:

- The PE 420 (QMS 422/421) or PKR 250/251 (QMS 200) hardware option is not continually ON because the plasma in the cell may interfere with the gas composition
- Quadstar 32-bit only turns on the high voltage if in the current MID/MCD parameter file for at least one channel under Detector the PE 420 or PKR 251 option has been chosen as signal source
- The high voltage is turned off again when the corresponding measurement type (MID, MCD, Sequencer) is terminated
- Please note that it may take some time before the gauge head ignites and indicates the correct pressure. See instruction manual of the gauge head.

Operation without sequencer

In manual operation, i.e. when an MID or MCD measurement is started directly via the menu, Quadstar 32-bit tests the parameter file and activates the PE 420/TPR 250 hardware if PE or PKR/TPR has been selected as the detector in at least one of the activated channels. After the measurement the PE 420/TPR 250 hardware is switched off again.

Operation with sequencer

In sequencer mode, the PE 420 hardware should be switched on by using the keyword TotGauge=on available in the SetQMS() commands.

When the sequence is loaded, a sequence test is performed to check whether or not the PE 420/TPR 250 hardware is needed in an MID or MCD measurement and, if so, whether the SetQMS(TotGauge=on) command is used in the sequence. If this command is missing, the sequence test displays a warning that the PE 420/TPR 250 hardware will only be switched on immediately before the measurement.

As with manual operation, the PE 420/TPR 250 hardware is switched off when the sequencer is terminated.

QMS 200 (Prisma™) using External Protection

If you run a QMS 200 Prisma™ in conjunction with a PKR 250/251 and you have selected in Parset program under [Config] > [QMS] > [Protection (SEM+FIL)] the option EXTERN or EXT-PROT, then the PKR 250/251 is used for pressure monitoring to protect filament and SEM. Therefore it stays switched on continuously. Switching it off by sequencer command would cause a malfunction.

2.2 Measurement main menu




Fig. 2-13

Menus

File

Display the Quadstar 32-bit Measure program version and terminate the program. See Chapter 2.3 File menu, 75.

Scan

Measurement program for Scan Analog and Scan Bargraph mode. See Chapter 2.4 Scan menu, 76.

MID

Measurement program for Multiple Ion Detection (MID) measurement type. See Chapter 2.5 MID measurement type, 87.

MCD

Measurement program for Multiple Concentration Detection (MCD) measurement type. See Chapter 2.6 MCD measurement type, 103.

Sequence

Starting a sequence. See Chapter 2.7 Measuring with the sequencer (Sequence), 120.

Calibration

Calibration programs. See Chapter 2.8 Calibration, 134.

Manual

Operation/display of the digital input and output channels (DI/DO). See Chapter 2.9 Manual Parameters, 145.

Comm

Switch the communication to the QMS on/off. See Chapter 2.10 Communication menu, 147.

Setup

Switch the emission and the SEM on/off. See Chapter 2.11 Setup, 148.

Help

Display the «Key Fragment Ions» table and help contents. See Chapter 1.9.1 Spectrum interpretation, 56.

NOTE:

The measurement programs are dependent on the parameter files created in the Parameter Setup Program. If any configuration changes are made ([Parameter Setup] > [Config]), Measure must be restarted before the changed parameters can be used.

2.3 File menu




Fig. 2-14

Submenus

About Measure

Information on the Measure program. See Chapter 2.3.1 About Measure, 75.

Exit

Quit the Measure program. See Chapter 2.3.2 Exit, 75.

2.3.1 About Measure

Under this menu entry the version number of the Quadstar 32-bit Measure program is displayed:




Fig. 2-15

2.3.2 Exit

Quit the Quadstar 32-bit Measure program.

2.4 Scan menu

Use this menu to select the scan measurement type.




Fig. 2-16

Submenus

Analog

Start Scan Analog measurement. See Chapter 2.4.2 Scan Analog, [177](#).

Bargraph

Start Scan Bargraph measurement. See Chapter 2.4.3 Scan Bargraph, [82](#).

2.4.1 Setup submenu

Setup activates a submenu to setup the scan measurement:

- Analog: Input of control parameters for the Scan Analog measurement type. See Chapter 2.4.2.5 Setup parameters, [80](#).
- Bargraph: Input of control parameters for the Scan Bargraph measurement type. See Chapter 2.4.3.4 Setup parameters, [85](#).

NOTE:

The measurement can only be started if the communication with the QMS is activated. See Chapter 2.10 Communication menu, [147](#).

2.4.2 Scan Analog

Open the Quadstar 32-bit file manager by choosing [Scan] > [Analog]. It offers the already existing measurement parameter files for Scan Analog (e.g. *airdemo.sap*) to choose from.

Fig. 2-17, 77 shows the measurement of three analog spectra. The status text «Number of mean value cycles» appears when the entry «Mean Value of Scans» is marked in the «Scan Analog Setup» dialog box.




Fig. 2-17

NOTE:

Multichannel operation, i.e. measurement of different scans within the same cycle(s), can be activated in the channel parameter editor by setting several channels to State=Enable. By [Parameters] > [Setup] the display can be set up.

2.4.2.1 Scan Analog measurement menu

The acquisition and display of the Scan Analog measurement data is controlled via the following menu:



Fig. 2-18

Submenus

File

Handling of parameter and measurement data files and return to the main menu of the Measure program. See Chapter 2.4.2.2 File submenu, 78.

Operation

Start/Stop the measurement. See Chapter 2.4.2.3 Operation submenu, 79.

Parameters

Parameter handling for control (display) and measurement parameters. See Chapter 2.4.2.4 Parameters submenu, 79.

2.4.2.2 File submenu




Fig. 2-19

Load Parameter

Load a different measurement parameter file and start the measurement.

Save Cycle Data

Store the next cycle(s) as a Scan Analog measurement data file. Cycle data can be stored only if «Temporary Save» has been turned off in the «Setup» dialog box.

After having chosen the measure data file in the file manager, the «Cycle Save» dialog box is displayed in which the number of cycles can be entered which are to be measured and stored:




Fig. 2-20

If an existing file containing compatible measurement data is selected, the new data are alternatively appended (Append) or the file is overwritten (Reset).

Save Reference

Store the next cycle as a Scan Analog reference data file.

Load Reference

Load a reference file to compare with the current measure data.

Close Reference

Close the reference file.

File Info Notes

Add notes and remarks to the documentation and store them.

Close

Quit the Scan Analog measurement type. A cycle that has been stopped will not be saved.

NOTE:

Do not modify the setup or channel parameters during the storage and don't change the size of the measure window. Otherwise the measurement is started anew.

2.4.2.3**Operation submenu**

Fig. 2-21

Go

Start the measurement.

Halt

Stop the measurement.

2.4.2.4**Parameters submenu**

Fig. 2-22

Setup

Enter various parameters that control the type of measurement as well as the display of the measured data.

Channel

Modify the current measurement parameters (channel parameters).

2.4.2.5 Setup parameters

Choose [Parameters] > [Setup] to enter the setup parameters:




Fig. 2-23

Control elements

Measurement

- Single: Only one Scan Analog measurement (1 Cycle) is performed
- Repeat: Repetitive measurement of the activated channels
- Timer: The measurement is periodically restarted

See Chapter 2.1.1.2 Measurement methods, 62.

Y-Scale

- Lin: Linear display of measurement data
- Log: Logarithmic display of measurement data

Scan Arrangement

Arrangement of the scans if several scans are to be displayed.

- Vertical: The individual scans are arranged one below the other
- Horizontal: The individual scans are arranged side by side

Display Mode

Display mode if a reference is loaded. In this case, the current measure values (Current), the reference values (Reference) and the difference between both (Difference) can be displayed.

- Multi Chart: Current, Reference and Difference are each displayed in an individual diagram with individual axes
- Single Chart: Current, Reference and Difference are displayed in the same diagram with the same axes
- Current: Display current measure values. The color of Current is determined by the channel parameters under [Parameters] > [Channel] > [Display] > [Color].
- Reference: Display values of the reference measurement. The color of Reference correspond to the ones of Bar Reference and Bar Last/Difference found in Parset program under [Config] > [Screen Color].

- **Difference:** Display difference between Current and Reference. The color of Difference correspond to the ones of Bar Reference and Bar Last/Difference found in Par-set program under [Config] > [Screen Color].
- **Decades:** Number of decades to display the difference

NOTE:

Reference and Difference are only available if a reference measurement has been loaded before by [File] > [Load Reference].

Save Data As

Determine the data that is to be saved:

- **Difference:** Save difference between Current and Reference.
- **Current:** Save current measure values (Current).

Temporary Save

Save the measurement data in a temporary file. When the measurement is stopped, a dialog box for entering a file name is displayed. Exactly the values displayed in the measurement window are stored, even a started (incomplete) cycle and a remaining part of the preceding cycle.

When this option is chosen, it is not possible to save data by [File] > [Save Cycle Data].

Mean Value of Scan

Performs continuous calculation of mean values during the measurement.

Timer Step

Timer interval for the measurement method Timer.

2.4.2.6 Channel parameters

Choose [Parameters] > [Channel] to change the QMS channel parameters.

Operating modes

For a Scan Analog measurement you can select the following operating modes in the channel parameter editor. Click [Parameters] > [Channel] > [Mass Mode].

- **SCAN-N:** Scan Analog measurement. The intensities are recorded in the form of a mass scan.
- **SCAN-F:** Scan Analog measurement with FIR filter. The intensity values are additionally subjected to an FIR (Finite Impulse Response) filter algorithm. This filter largely eliminates superimposed noise from the raw electrometer signal. In this way also very small peaks can be detected against the background.

Integral scan

For measuring e.g. the total pressure, an integral scan can be produced. To activate it, set the Resolution parameter in the group «Mass» to OFF.

Interface load

In Scan Analog the interface load is the highest because up to 64 measured values per amu are transmitted. At a scan speed of < 0.5 s/amu this can lead to a relatively high load on the interface and the computer. If the measurement data are not stored but only displayed graphically, a lower data rate should normally be adequate. In Quadstar 32-bit it is possible to reduce the data rate via the Steps parameter (group «mass»), e.g. from 64 to 16 measured values/amu (QMS 422/421 only). Please note that if concurrent stor-

age is active, only the reduced number of measured values/amu will be stored. Steps can be changed in the same group (register card Mass) as Speed.

NOTE:

Due to fast CPUs and graphics cards in modern computers, the effect of the interface load is generally not that high anymore.

2.4.3 Scan Bargraph

In Scan Bargraph mode, only the maximum peak intensities and the corresponding mass numbers are collected. This is done by the peak detection (Peak-L and Peak-F) integrated in the QMS.

Open the Quadstar 32-bit file manager by choosing [Scan] > [Bargraph]. It offers the already existing measurement parameter files for Scan Bargraph to choose from. Fig. 2-24, 82 shows the measurement of three bargraph spectra.




Fig. 2-24

NOTE:

Multichannel operation, i.e. measurement of different scans within the same cycle(s), can be activated in the channel parameter editor by setting several channels to State=Enable. By [Parameters] > [Setup] the display can be set up.

The acquisition and display of the San Bargraph measurement data is controlled via the following menu:



Fig. 2-25

Submenus

File

Handling of parameter and measurement data files and return to the main menu of the Measure program. See Chapter 2.4.3.1 File submenu, 83.

Operation

Start/Stop the measurement. See Chapter 2.4.3.2 Operation submenu, 84.

Parameters

Parameter handling for control (display) and measurement parameters. See Chapter 2.4.3.3 Parameters submenu, 84.

2.4.3.1**File submenu**

Fig. 2-26

Load Parameter

Load a different measurement parameter file and start the measurement.

Save Cycle Data

Store the next cycle(s) as a Scan Bargraph measurement data file. Cycle data can be stored only if «Temporary Save» has been turned off in the «Setup» dialog box. Save Cycle Data is disabled in Multiple Scan, i.e. if more than one channel is switched on ([State] > [Enable]) in the channel parameter editor.

After having chosen the measure data file in the file manager, the «Cycle Save» dialog box is displayed in which the number of cycles can be entered which are to be measured and stored:




Fig. 2-27

If an existing file containing compatible measurement data is selected, the new data are alternatively appended (Append) or the file is overwritten (Reset).

Save Reference

Store the next cycle as a Scan Bargraph reference data file.

Load Reference

Load a reference file to compare with the current measure data.

Close Reference

Close the reference file.

File Info Notes

Add notes and remarks to the documentation and store them.

Close

Quit the Scan Bargraph measurement type. A cycle that has been stopped will not be saved.

NOTE:

Do not modify the setup or channel parameters during the storage and don't change the size of the measure window. Otherwise the measurement is started anew.

2.4.3.2 Operation submenu



Fig. 2-28

Go

Start the measurement.

Halt

Stop the measurement.

2.4.3.3 Parameters submenu



Fig. 2-29

Setup

Enter various parameters that control the type of measurement as well as the display of the measured data.

Channel

Modify the current measurement parameters (channel parameters).

2.4.3.4 Setup parameters

Choose [Parameters] > [Setup] to enter the setup parameters:




Fig. 2-30

Control elements

Measurement

- Single: Only one Scan Bargraph measurement (1 Cycle) is performed
- Repeat: Repetitive measurement of the activated channels
- Timer: The measurement is periodically restarted

See Chapter 2.1.1.2 Measurement methods, 62.

Y-Scale

- Lin: Linear display of measurement data
- Log: Logarithmic display of measurement data

Scan Arrangement

Arrangement of the scans if several scans are to be displayed.

- Vertical: The individual scans are arranged one below the other
- Horizontal: The individual scans are arranged side by side

NOTE:

No data can be stored with several scans.

Display Mode

Display mode if a reference is loaded. In this case, the current measure values (Current), the reference values (Reference) and the difference between both (Difference) can be displayed.

- Multi Chart: Current, Reference and Difference are each displayed in an individual diagram with individual axes
- Single Chart: Current, Reference and Difference are displayed in the same diagram with the same axes
- Current: Display current measure values. The color of Current is determined by the channel parameters under [Parameters] > [Channel] > [Display] > [Color].

- Reference: Display values of the reference measurement. The color of Reference correspond to the ones of Bar Reference and Bar Last/Difference found in Parset program under [Config] > [Screen Color].
- Difference: Display difference between Current and Reference. The difference is always displayed as rounded mass numbers. The color of Difference correspond to the ones of Bar Reference and Bar Last/Difference found in Parset program under [Config] > [Screen Color].
- Decades: Number of decades to display the difference

NOTE:

Reference and Difference are only available if a reference measurement has been loaded before by [File] > [Load Reference].

Save Data As

Determine the data that is to be saved:

- Difference: Save difference between Current and Reference.
- Current: Save current measure values (Current).

Timer Step

Timer interval for the measurement method Timer.

2.4.3.5 Channel parameters

Choose [Parameters] > [Channel] to change the QMS channel parameters.

Operating modes

For a Scan Bargraph measurement you can select the following operation modes in the channel parameter editor. Click [Parameters] > [Channel] > [Mass Mode].

- PEAK-L: Scan Bargraph measurement with threshold definition. Only values greater or equal to the threshold are displayed.
- PEAK-F: Scan Bargraph measurement with threshold and FIR filter. In the QC the raw measurement signal is subjected to a FIR (Finite Impulse Response) filter algorithm. This filter largely eliminates superimposed noise from the raw measure signal. In this way also very small peaks can be detected against the background.
- STAIR: Record a mass spectrum with integer mass jumps across the mass range defined by the QMS channel parameters First and Width.

2.5 MID measurement type

In MID (Multiple Ion Detection) mode, one measurement on one defined mass is done per channel. The measured intensities can be displayed as a table (Table), as a bargraph (Bargraph), or as a function of time (Versus Time). If you change the display mode while a measurement is running, it will be restarted.

MID Menu

Use this menu to select the MID measurement type.




Fig. 2-31

Submenus

Table...

Start the measurement, display measured data as a table. See Chapter 2.5.6 MID Table, 92.

Bargraph...

Start the measurement, display measured data as a bargraph. See Chapter 2.5.7 MID Bargraph, 96.

Versus Time...

Continuous display of the measured data as a function of time, measured cycles, etc. See Chapter 2.5.8 MID Versus Time, 98.

- To start a measurement, a parameter file must be specified. Quadstar 32-bit responds by displaying the file manager in order to choose the appropriate file.
- If display of reference data is chosen in MID Table setup ([Display Option] > [Reference]) or in MID Bargraph setup ([Display Reference]), the Quadstar 32-bit file manager is called again in order to choose a reference file (*.MDR). Versus Time does not support reference display.

2.5.1 Setup submenu

Setup activates a submenu to setup the MID measurement:

- Table: Input of control parameters for the MID Table measuring mode. See Chapter 2.5.6.1 Parameter setup, 92.
- Bargraph: Input of control parameters for the MID Bargraph measuring mode. See Chapter 2.5.7.1 Parameter setup, 96.
- Versus Time: Input of control parameters for the MID Versus Time measuring mode. See Chapter 2.5.8.1 Parameter setup, 99.

Last as Table/Last as Bargraph

With Last as Table/Last as Bargraph the measurement data of the last measure cycle can be recalled. See Chapter 2.5.9 MID last measured values, 101.

NOTE:

Last as Table and Last as Bargraph can only be chosen, if values of the last measurement are available.

2.5.2 MID measurement menu

Recording and display of the MID measurement data is controlled via the following menu:



Fig. 2-32

File

Handling of parameter and measurement data files and return to the main menu of the Measure program.

Operation

Start/stop a measurement, optimize transmission speed.

Display

Choose the display mode (Table, Bargraph or Versus Time).

Parameters

Handling of control parameters (display) and measurement parameters.

2.5.2.1 File submenu




Fig. 2-33

Load Parameter

Load a different measurement parameter file and start the measurement.

Save Cycle Data

Store the next measurement cycle(s) as MID data file.

After having chosen the measure data file in the file manager, the Cycle Save dialog box is displayed.

If an existing file with compatible data is selected, the new data can alternatively be appended (Append), or the file is overwritten (Reset).




Fig. 2-34

Control elements

- Number of cycles to save: Enter the number of cycles, you want to measure and store.
- Display: For fast measurements the display of the measurement data can be turned off. See Chapter 2.1.1.4 Fast MID/MCD measurements, 62.

Save Reference

Store the last measurement (1 cycle) as a MID reference data file.

Load Reference

Load reference data, compare to the current data.

File Info Notes

Add notes and remarks to the documentation and store them.

Close

Quit the MID measurement.

NOTE:

Do not modify the setup or channel parameters during the storage and don't change the size of the measure window, 'cause otherwise the measurement is started anew.

2.5.2.2

Operation submenu



Fig. 2-35

Go/Halt

Start/stop the measurement.

2.5.2.3

Display submenu

Choose the display mode. When you change the display mode while a measurement is running, the measurement will be restarted.



Fig. 2-36

Table

Display data as a table.

Bargraph

Display data as a bargraph.

Versus Time

Display data Versus Time.

2.5.2.4 Parameters submenu



Fig. 2-37

Setup

Enter various parameters that control the type of measurement as well as the display of the measured data.

- MID Table - See Chapter 2.5.6 MID Table, 92
- MID Bargraph - See Chapter 2.5.7 MID Bargraph, 96
- MID Versus Time - See Chapter 2.5.8 MID Versus Time, 98

Channel

Modify the measurement parameters (channel parameters).

2.5.3 Offset alignment

For MID or MCD measurements, the mass dependent (detector-) offset of the QMS 422/421 can be compensated by an automatic correction measurement. On each channel, whose parameter Zero Mass contains a mass number, a measurement is carried out on this zero mass. The intensity received thereby is the offset that will be subtracted in the subsequent MID (or MCD) measurement from the measured value of the corresponding channel.

Procedure

- 1 Open the MID or MCD parameter file that you intend to use for the measurement and click on the register card Mass of the channel parameter editor.
- 2 Enter for each channel under Zero Mass the mass number, on which the offset is to be measured (or OFF, if the offset of this channel is disregarded). Choose a zero mass, where no peak appears (Default: 5.5 amu).
- 3 Close and save the parameter file and start the MID or MCD measurement. Now you see for each channel, for which a Zero Mass has been defined, the message «Channel x: Zero Equalize in progress» shown in the status bar; after that, the

MID or MCD measurement is started. This correction measurement (Zero Equalize) is automatically redone, if the MID/MCD measurement is restarted or a channel parameter has been changed.

NOTE:

The automatic offset alignment can be disabled in measurements performed with the sequencer (MID(), MCD(), GSC() and ZeroGas() instruction) by the SetPar() command (ZeroEq=off). In manual operation, this automatic offset alignment can only be switched off by setting the Zero Mass of all channels to OFF.

2.5.4 Background subtraction (Zero Gas)

Please note that the background should be measured before each measurement. See Chapter 2.8.2 Background measurement (Zero Gas), 137. The results of the last background measurement are stored and subtracted from the measured intensities during each measurement, if under [Parameters] > [Setup] the Option Zero Gas Subtraction is checked. This can lead to incorrect measurement data if the background has changed.

Example

If previously you had measured «air» as the background (Zero Gas) and subsequently want to measure air, the measured values (intensities) you obtain are only the difference from the currently measured values to the stored values.

NOTE:

Make sure that the background is not measured under incorrect conditions, e.g. with an opened inlet valve. In this case the intensities stored as backgrounds would be nearly equal to the actual measurement data, so the result of the subtraction would be some small differences that are hardly related to the measured data.

If the background cannot be measured:

- Check that in the Display Saved Values program, Auxiliary menu, under Zero Gas, the stored background measurements are still valid.
- Delete invalid entries (and store again) in the File menu.

NOTE:

In measurements with the sequencer (MID(), MCD() command), the background subtraction can be turned off by the SetPar() command (ZeroSub=off).

2.5.5 Temporary save

In MID and MCD mode Quadstar 32-bit allows temporary storage of measurement data. Under [Parameters] > [Setup], check the option Temporary Save, so the measurement data are written to a temporary file. The number of stored cycles is displayed in the status bar («XX cycles saved to temporary file»). After the desired number of cycles has been measured (or if you leave the measurement), Quadstar 32-bit asks whether to save the displayed measurement data permanently or not.



Fig. 2-38

2.5.6 MID Table

MID Table represents the measured ion currents or the measured values of an analog input channel (AI) as a table. Depending on the selected display mode the results are normalized to the largest value (=100%) or shown together with comparative data (reference, last, etc.).

2.5.6.1 Parameter setup

Under [Parameters] > [Setup] you can enter the control parameters:




Fig. 2-39

Control elements

Measurement

- Single: Perform only one measurement (= 1 cycle) across all activated channels.
- Repeat: Repetitive measurement of the activated channels.
- Timer: The measurement is periodically restarted.
- Dynamic: The measurement data is temporarily stored in the QMS buffer. See Chapter 2.1.1.2 Measurement methods, 62.

Relative Values

Display the ion currents.

- None: Indicate the intensities in amperes.

- Sum=100%: Normalize the displayed data to «sum of all intensities = 100%».
- Max=100%: Normalize the displayed data to «highest intensity = 100%».

Display Option

- None: No additional display options.
- Integral: Additional display of the integral of all measurements performed since the start.
- Last: Additional display of the measure data of the last cycle.
- Reference: Additional display of the measure data of a reference file.
- Multiple column: Display the measure data in several columns (scrolling). See Section «Example 2», 95.

Mass scale correction

Mass scale correction on/off. See Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

Zero Gas Subtraction

Zero gas subtraction on/off. See Chapter 2.8.2 Background measurement (Zero Gas), 137.

Alarm Display

Display the alarm status for all channels with activated threshold monitoring.

Alarm Sound

Enable the alarm buzzer.

Alarm Error Message

Display an alarm additionally as an error.

If the Error History is enabled, incoming threshold monitoring alarms are logged in the Error History file. Configurable in the Parameter Setup program under [Setup] > [General].

Temporary Save

Store the measurement data in a temporary file.

When the measurement is stopped, a dialog box is displayed for entering the file name. See Chapter 2.5.5 Temporary save, 91.

Mult. Column Width

Column width in pixels in Multiple Column display mode.

Digits for Float Values

Number of digits after the decimal point for float values.

Timer Step

Timer interval for the measurement method Timer.

2.5.6.2

Data display

Example 1

Display Option: Reference

In this mode the current measurement data are compared to the measurement data of a reference file. When the measurement is started, Quadstar 32-bit asks for the name of the reference file.

The screenshot shows a software window titled "[M1] QUADSTAR 32-bit Measurement - [MID < airdemo.mip >]". The window has a menu bar with "File", "Operation", "Display", and "Parameters". Below the menu bar is a table with the following columns: Nbr, Type, Ident, Cur Val, Rel Val, Ref Val, Delta, and Alarm. The table contains 7 rows of data. At the bottom of the window, there are several icons and labels: a green circle, a yellow triangle, a red square, "SIM", and "QMS 200".

Nbr	Type	Ident	Cur Val	Rel Val	Ref Val	Delta	Alarm
0	Cur.	14.00	8.316E-06A	25.558%	8.316E-06A	-0.007%	
1	Cur.	16.00	2.487E-06A	7.645%	2.487E-06A	-0.002%	
2	Cur.	18.00	1.249E-06A	3.838%	1.249E-06A	-0.002%	
3	Cur.	28.00	9.871E-06A	30.338%	9.877E-06A	-0.059%	
4	Cur.	32.00	8.316E-06A	25.560%	8.316E-06A	-0.000%	
5	Cur.	40.00	1.668E-06A	5.125%	1.668E-06A	0.000%	
6	Cur.	44.00	6.299E-07A	1.936%	6.300E-07A	-0.001%	

Fig. 2-40

Specify the desired number of decimal places under [Parameters] > [Setup] > [Digits for Float ValuesF].

Displayed Columns

Number and type of the displayed columns depend on the options selected under [Parameters] > [Setup] > [Display Option]:

Nbr

Current number of the displayed measured value.

Type

Data type of the displayed measured value.

- Cur.: Ion current.
- Cnt. R: Ion counter content.
- AI: Measured value of an AI channel.
- PI: Pirani measurement signal(QMS 422/421).
- PE: Cold cathode gauge signal.(QMS 422/421).
- TPR: Measured value of a Compact Pirani Gauge TPR 250(QMS 200).
- PKR: Measured value of a Compact Full Range Gauge PKR 250/251 (QMS 200).

Ident

Data type identification.

- XX.XX: Mass number of the ion current (Cur.) or the counter value (Cnt.)
- 'nnnnn': Logical name of the (AI) channel.
- PIR0/PIR1: Pirani 0, Pirani 1 (PI)(QMS 422/421).
- PEN: Cold cathode gauge (PE)(QMS 422/421).
- TPR: Compact Pirani Gauge (QMS 200).
- PKR: Compact Full Range Gauge (QMS 200).

QMU

QMU 112 Channel number on which the specified mass has been measured.

Cur Val

Current measured value.

Rel Val

Relative measured value (sum = 100% or max = 100%).

Integral

Integrated measurement signal since the start of the measurement.

Last Val

Measured value of the last cycle.

Ref Val

Reference value.

Delta

Deviation Cur Val - Last Val in %.

Alarm

Switching functions.

Threshold monitoring

The threshold status is displayed graphically in the last column under alarm.

See Chapter 3.2.5.7 Switching functions (Trip), 186.

Example 2

Display Option: Multiple Column

The measure data is displayed in several columns (width selectable under Setup), where each column is one measured cycle.

Measurement Number: 17				MID airdemo.mip 12:47:28	MID airdemo.mip 12:47:27	MID airdemo.mip 12:47:27	MID airdemo.mip 12:47:28
Nbr	Type	Ident	Unit				
0	Cur.	14.00	A	8.305E-09	8.305E-09	8.304E-09	8.304E-09
1	Cur.	15.97	A	2.471E-09	2.471E-09	2.472E-09	2.470E-09
2	Cur.	17.97	A	1.231E-09	1.231E-09	1.231E-09	1.231E-09
3	Cur.	28.00	A	9.980E-09	9.980E-09	9.981E-09	9.979E-09
4	Cur.	32.00	A	8.306E-09	8.305E-09	8.306E-09	8.305E-09
5	Cur.	40.00	A	1.653E-09	1.654E-09	1.654E-09	1.654E-09
6	Cur.	44.00	A	6.108E-10	6.108E-10	6.106E-10	6.108E-10

Fig. 2-41

The individual columns are overwritten cyclically from the left to the right:

Entering Parameter setup

For setting the display mode choose [Parameters] > [Setup] (or double click the left mouse button) and for changing the QMS channel parameters choose [Parameters] > [Channel] (or double click the right mouse button).

2.5.7 MID Bargraph

In MID Bargraph the measurement is started with the parameter file selected with the Quadstar 32-bit file manager and the measurement data are displayed as a bar chart.




Fig. 2-42

In the comments field, the data of the bar marked with the mouse are displayed.

2.5.7.1 Parameter setup

Under [Parameters] > [Setup] you can enter the control parameters:




Fig. 2-43

Control elements**Measurement**

- Single: Perform only one measurement (= 1 cycle) across all activated channels.
 - Repeat: Repetitive measurement of the activated channels.
 - Timer: The measurement is periodically restarted.
 - Dynamic: The measurement data is temporarily stored in the QMS buffer.
- See Chapter 2.1.1.2 Measurement methods, 62.

Y-Scale

- Lin: Linear display of measurement data.
- Log: Logarithmic display of measurement data.

Color Mode

- Single Color: The measurement data are displayed in the colors defined in the Parameter Setup program under [Config] > [Screen Color].
- Multiple Color: The measurement data are displayed in the colors that are defined in the channel parameter editor for the corresponding channels. The bars of a channel (normal, last, reference) can only be distinguished by their patterns.

Mass scale correction

Mass scale correction on/off. See Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

Zero Gas Subtraction

Zero gas subtraction on/off. See Chapter 2.8.2 Background measurement (Zero Gas), 137.

Alarm Display

Display the alarm status for all channels with activated threshold monitoring.

Alarm Sound

Enable the alarm buzzer.

Alarm Error Message

Display an alarm additionally as an error.

If the Error History is enabled, incoming threshold monitoring alarms are logged in the Error History file. Configurable in the Parameter Setup program under [Setup] > [General].

Temporary Save

Store the measurement data in a temporary file.

When the measurement is stopped, a dialog box is displayed for entering the file name. See Chapter 2.5.5 Temporary save, 91.

Display Last

Display the measured data of the last cycle. The bar is shown on the right of the current measured value.

Display Reference

Display a reference measurement. The bar is shown on the right of the current measured value and marked with the pattern chosen by Bar Reference Display. The reference file is selected by the Quadstar 32-bit file manager when the measurement is started.

X-Axis Lettering

X-axis lettering on/off.

Y-Grid

Display a Y-grid on the measurement screen.

Bar Normal Display

Pattern of the bar for the current measured value.

Bar Last Display

Pattern of the bar for the measured value of the last cycle.

Bar Reference Display

Pattern of the bar for the reference measurement.

Bar Width

Bar width in pixels.

Timer Step

Timer interval for the measurement method Timer.

2.5.7.2 Data display

Parameter setup

For setting the display mode choose [Parameters] > [Setup] (or double click the left mouse button) and for changing the QMS channel parameters choose [Parameters] > [Channel] (or double click the right mouse button).

Operation with the QMU 112

When operating with the QMU 112, the QMU channel is shown in square brackets behind the mass number. For example «28 [7]» means mass 28, measure on the QMU channel 7.

Threshold monitoring

The threshold status is displayed graphically below the X-axis. See Chapter 3.2.5.7 Switching functions (Trip), 186.

2.5.8 MID Versus Time

In MID Versus Time the measurement is started with the parameter file selected in the file manager, and the measurement data are displayed as a function of time or number of cycles:




Fig. 2-44

2.5.8.1 Parameter setup

Under [Parameters] > [Setup] you can enter the control parameters:




Fig. 2-45

Control elements

Measurement

- Single: Perform only one measurement (= 1 cycle) across all activated channels.
 - Repeat: Repetitive measurement of the activated channels.
 - Timer: The measurement is periodically restarted.
 - Dynamic: The measurement data is temporarily stored in the QMS buffer.
- See Chapter 2.1.1.2 Measurement methods, 62.

X-Scale

- Cycle: Assign the cycle number to the X-axis label.
- Time relative: Assign the relative time since the start to the X-axis label.
- Time absolute: Assign the absolute time (time of day) of the measurement to the X-axis label.

Y-Scale

- Lin: Linear display of measurement data.
- Log: Logarithmic display of measurement data.

Line Type

- Solid: The measured values are connected by a solid line.
- Dot: The measured values are shown as individual dots.

Mass scale correction

Mass scale correction on/off. See Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

Zero Gas Subtraction

Zero gas subtraction on/off. See Chapter 2.8.2 Background measurement (Zero Gas), 137.

Alarm Display

Display the alarm status for all channels with activated threshold monitoring.

Alarm Sound

Enable the alarm buzzer.

Alarm Error Message

Display an alarm additionally as an error.

If the Error History is enabled, incoming threshold monitoring alarms are logged in the Error History file. Configurable in the Parameter Setup program under [Setup] > [General].

Temporary Save

Store the measurement data in a temporary file.

When the measurement is stopped, a dialog box is displayed for entering the file name. See Chapter 2.5.5 Temporary save, 91.

Help table

On the right, beside the measure data, the measured values are displayed additionally to the measure data identifications.

X-Grid

Display an X-grid.

Y-Grid

Display a Y-grid.

Marker

Mark the measured values with symbols.

Roll Axis

Each time the end of the X-axis is reached, scroll the screen to the left by 10%.

Cycles

Define the number of cycles that are displayed in the measurement data window before the old data are overwritten. For the «Display after Measure» method, Cycles is the number of cycles to be measured.

Time-Scale

Select the period of time on the X-axis.

Timer Step

Timer interval for the measurement method Timer.

Measure/Display Mode

- Measure and Display: The values are measured and displayed cycle by cycle.
- Display after Meas.: The values are displayed after the measurement.

2.5.8.2**Data display****Threshold monitoring**

Changes of the threshold status are indicated with small arrows above and below the measured values.

The arrows above the measured values signalize threshold state 1 (=A), the arrows below the measured values signalize threshold state 2 (=B). An upward pointing arrow

indicates that the threshold monitoring has responded. A downward pointing arrow indicates that the threshold monitoring has returned to normal condition. See Chapter 3.2.5.7 Switching functions (Trip), 186.

Parameter setup

For setting the display mode choose [Parameters] > [Setup] (or double click the left mouse button) and for changing the QMS channel parameters choose [Parameters] > [Channel] (or double click the right mouse button).

2.5.9 MID last measured values

[MID] > [Last as Table] and [MID] > [Last as Bargraph display] the last measured cycle.



Fig. 2-46

File

File handling and return to the main menu of the Measure program.

Display

Select the display mode (Table or Bargraph).

Parameters

Control parameter (display) and measurement parameter handling.

2.5.9.1 File submenu



Fig. 2-47

Load Reference

Load a reference data file.

Save Reference

Store the last measurement (1 cycle) as a MID reference data file.

Close

Quit the MID Last Values mode.

2.5.9.2 Display submenu



Fig. 2-48

Table

Display as a table.

Bargraph

Display as a bargraph.

2.5.9.3

Parameters submenu



Fig. 2-49

Setup

Setup the display of the measured data.

2.6 MCD measurement type

The MCD (Multiple Concentration Detection) measurement type provides concentration measurements. The measured ion currents are input to a solution matrix and the individual concentrations of the components in the gas are determined via calibration factors.

Calibration factor library

The Quantitative Analysis of Quadstar 32-bit works with a calibration factor library. The size of this library is not fixed, i.e. it grows with the input of each new component. For each component that is referenced in an application or calibration matrix, this library contains the calibration factors and the corresponding mass numbers. The number of calibration factors is limited to 24 per component.

The application matrix does not contain the actual calibration factors but only «pointers» to the location in the calibration factor library. The advantage of this method is, that the latest calibration factors are always used when the application matrix is loaded.

Matrix solution method

The matrix offers the possibility to define a calibration factor for every gas component and mass number that occur in the measurement task (calibration factors may also be 0). The resulting equation system must then be solved by a suitable process.

For solving the equation system, Quadstar 32-bit uses the method of elimination according to Givens. In this method the linear equation system is transformed to an equivalent system with a triangular matrix as a matrix of coefficients. Such equation systems are easy to solve by means of backward insertion. The process of elimination according to Givens allows also super-determinate equations, i.e. it may contain more equations than «unknowns» (more measured intensities than wanted concentrations). For solving such equations Givens employs the method of least squares. In this way inaccuracies can be minimized.

Additional information on the algorithm for solving linear equation systems used in the Quantitative Analysis of Quadstar 32-bit can be found in the books «COMPUTER-MATHEMATIK» and «COMPUTERMATHEMATIK Lösungen» by Prof. Dr. Walter Gander.

Normalization

The calculated values must be normalized (related to a reference quantity). If all gas components are known, normalization modes 'Maximum=100%' or 'Sum=100%' can be used, otherwise 'Extern' must be used. See Section «Setting the external normalization», § 205.

Literature

- COMPUTERMATHEMATIK by Prof. Dr. Walter Gander
Birkhäuser Verlag, Basel/Boston/Stuttgart 1985
ISBN 3-7643-1688-8 (Language: German)
- COMPUTERMATHEMATIK Lösungen by Prof. Dr. Walter Gander
Birkhäuser Verlag, Basel/Boston/Stuttgart 1986
ISBN 3-7643-1802-3 (Language: German)

MCD Menu

In MCD measurements (Multiple Concentration Detection) the individual measurement channels are set up similarly to MID measurements. The results can be represented as a table (Table), as a bar chart (Bargraph), or as a function of time (Versus Time).




Fig. 2-50

Submenus

Table...

Start the measurement, display measured data as a table. See Chapter 2.6.7 MCD Table, [108](#).

Bargraph...

Start the measurement, display measured data as a bargraph. See Chapter 2.6.8 MCD Bargraph, [112](#).

Versus Time...

Start the measurement, continuous display of measured data as a function of time, measured cycles, etc. See Chapter 2.6.9 MCD Versus Time, [115](#).

- To start a measurement, a parameter file must be specified. Quadstar 32-bit responds by displaying the file manager in order to choose the appropriate file.
- If display of reference data is chosen in MCD Table setup [Display Option] > [Reference] or in MID Bargraph setup [Display Reference], the Quadstar 32-bit file manager is called again in order to choose a reference file (*.MDR). Versus Time does not support reference display.

2.6.1

Setup submenu

Setup activates a submenu to setup the MCD measurement:

- Table: Control parameter input (measurement method, display mode, etc) for the MCD Table measuring mode. See Chapter 2.6.7.1 Parameter setup, [109](#).
- Bargraph: Control parameter input (measurement method, display mode, etc) for the MCD Bargraph measuring mode. See Chapter 2.6.8.1 Parameter setup, [113](#).
- Versus Time: Control parameter input (measurement method, display mode, etc.) for the MCD Versus Time measuring mode. See Chapter 2.6.9.1 Parameter setup, [116](#).

Last as Table/Last as Bargraph

With [Last as Table] and [Last as Bargraph], measurement data of the last measured cycle can be recalled. See Chapter 2.6.10 MCD last measured values, [118](#).

NOTE:

[Last as Table] and [Last as Bargraph] can only be chosen, if measured values of a last measurement are available.

2.6.2 MCD measurement menu

The recording and display of MCD measurement data is controlled via the following menu:



Fig. 2-51

File

Handling of parameter and measurement data files and return to the main menu of the Measure program.

Operation

Start/stop a measurement, optimize transmission speed.

Display

Choose the display mode (Table, Bargraph or Versus Time).

Parameters

Handling of control parameters (display) and measurement parameters.

2.6.2.1 File submenu




Fig. 2-52

Load Parameter

Load a different measurement parameter file and start the measurement.

Save Cycle Data

Store the next several measurement data cycles as MCD data file.

After having chosen the measure data file in the file manager, the Cycle Save dialog box is displayed.

If an existing file with compatible data is selected, the new data can alternatively be appended (Append), or the file is overwritten (Reset).




Fig. 2-53

Control elements

- Number of cycles to save: Enter the number of cycles, you want to measure and store.
- Display: For fast measurements the display of the measurement data can be turned off. See Chapter 2.1.1.4 Fast MID/MCD measurements, 62.

Save Reference

Store the last measurement (1 cycle) as a MCD reference data file.

Load Reference

Load reference data, compare to the current data.

File Info Notes

Add notes and remarks to the documentation and store them.

Close

Quit the MCD measurement.

NOTE:

Do not modify the setup or channel parameters during the storage and don't change the size of the measure window. Otherwise the measurement is started anew.

2.6.2.2

Operation submenu



Fig. 2-54

Go/Halt

Start/stop the measurement.

2.6.2.3

Display submenu

Choose the display mode. When you change the display mode while a measurement is running, the measurement will be restarted.



Fig. 2-55

Table

Display data as a table.

Bargraph

Display data as a bargraph.

Versus Time

Display data Versus Time.

2.6.2.4 Parameters submenu



Fig. 2-56

Setup

Enter various parameters that control the type of measurement as well as the display of the measured data.

- MCD Table - See Chapter 2.6.7.1 Parameter setup, 109
- MCD Bargraph - See Chapter 2.6.8.1 Parameter setup, 113
- MCD Versus Time - See Chapter 2.6.9.1 Parameter setup, 116

Channel

Modify the measurement parameters (channel parameters).

Concentration

Setup the trip parameters (threshold monitoring), color assignment and concentration display.

2.6.3 Data display


How the data is displayed, depends on the display mode (Table, Bargraph or Versus Time).

Which data is displayed, depends on the parameter setup:


- Concentrations are displayed, if under [Parameters] > [Concentration] the Disp.Mode of the corresponding concentrations is set to ON.
- Ion currents for concentration calculation are displayed, if under [Parameters] > [Channel] on the register card Display the Disp.Mode of the corresponding channels is set to ON.
- Incidentally measured ion currents (which are not used for concentration calculation) are displayed, if under [Parameters] > [Channel], after the channels being in [State]

> [Matrix], more channels are switched on ([State] > [Enable]) and their Disp.Mode is set to ON.

2.6.4 Offset alignment

Quadstar 32-bit performs an offset alignment (Zero Equalize) before the MCD measurement, if the 'Zero Mass' of the measure channels are not switched OFF. See Chapter 2.5.3 Offset alignment,  90.

2.6.5 Background subtraction

Before every measurement the background should be determined, in order to subtract it from the measure data. Otherwise the Zero Gas Subtraction should be turned off. See Chapter 2.5.4 Background subtraction (Zero Gas),  91.

2.6.6 Operation with the QMU 112

MCD measurements don't provide full support for the QMU 112. The calibration factors are only valid for the QMU channel on which the calibration was performed. For correct analysis the measurements are to be performed exclusively on this channel.

2.6.7 MCD Table

MCD Table starts the measurement with the paramere file selected via the Quadstar 32-bit file manager and displays the calculated concentrations as a table.

2.6.7.1 Parameter setup

Under [Parameters] > [Setup] you can enter the control parameters:




Fig. 2-57

Control elements

Measurement

- Single: Perform only one measurement (= 1 cycle) across all activated channels.
 - Repeat: Repetitive measurement of the activated channels.
 - Timer: The measurement is periodically restarted.
 - Dynamic: The measurement data is temporarily stored in the QMS buffer.
- See Chapter 2.1.1.2 Measurement methods, 62.

Normalization Mode

Display of the concentrations.

- None: Indicate the intensities in amperes.
- Sum=100%: Normalize the displayed data to «sum of all intensities = 100%».
- Max=100%: Normalize the displayed data to «highest intensity = 100%».
- Extern: External normalization. See Section «Setting the external normalization», 205.

Display Option

- None: No additional display options.
- Integral: Additional display of the integral of all measurements performed since the start.
- Last: Additional display of the measure data of the last cycle.
- Reference: Additional display of the measure data of a reference file.
- Multiple column: Display the measure data in several columns (scrolling).

Mass scale correction

Mass scale correction on/off. See Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

Zero Gas Subtraction

Zero gas subtraction on/off. See Chapter 2.8.2 Background measurement (Zero Gas), 137.

Alarm Display

Display the alarm status for all channels with activated threshold monitoring.

Alarm Sound

Enable the alarm buzzer.

Alarm Error Message

Display an alarm additionally as an error.

If the Error History is enabled, incoming threshold monitoring alarms are logged in the Error History file. Configurable in the Parameter Setup program under [Setup] > [General].

Temporary Save

Store the measurement data in a temporary file.

When the measurement is stopped, a dialog box is displayed for entering the file name. See Chapter 2.5.5 Temporary save, 91.

Neg. Conc. Suppression

Suppress negative concentrations.

Mult. Column Width

Column width in pixels in Multiple Column display mode.

Digits for Float Values

Number of digits after the decimal point for float values.

Timer Step

Timer interval for the measurement method Timer.

2.6.7.2

Data display

Example 1:

Display Option: Last

In this mode the current measurement data are compared to those of the last cycle.

Nbr	Type	Ident	Cur Val	Last Val	Delta	Alarm
0	Conc.	Nitrogen	86.380%	86.680%	-0.347%	
1	Conc.	CarbonDiox	-185.682ppm	-185.410ppm	0.147%	
2	Conc.	Argon	50.479ppm	59.351ppm	-17.575%	
3	Conc.	Oxygen	13.633%	13.333%	2.204%	
4	Cur.	28.00	-1.044E-09A	-1.043E-09A	0.174%	
5	Cur.	32.00	-6.373E-10A	-6.200E-10A	2.711%	
6	Cur.	18.00	-3.055E-11A	-2.247E-11A	26.425%	

Fig. 2-58

Displayed Columns

Number and type of the displayed columns depend on the options selected under [Parameters] > [Setup] > [Display Option]:

Nbr

Current number of the displayed measured value.

Type

Data type of the displayed measured value.

- Conc.: Concentration.
- Cur.: Ion current.
- Cnt. R: Ion counter content.
- AI: Measured value of an AI channel.
- PI: Pirani measurement signal(QMS 422/421).
- PE: Cold cathode gauge signal.(QMS 422/421).
- TPR: Measured value of a Compact Pirani Gauge TPR 250(QMS 200).
- PKR: Measured value of a Compact Full Range Gauge PKR 250/251 (QMS 200).

Ident

Data type identification.

- XX.XX: Mass number of the ion current (Cur.), the counter value (Cnt.) or the name of the component.
- 'nnnnn': Logical name of the (AI) channel.
- PIR0/PIR1: Pirani 0, Pirani 1 (PI)(QMS 422/421).
- PEN: Cold cathode gauge (PE)(QMS 422/421).
- TPR: Compact Pirani Gauge (QMS 200).
- PKR: Compact Full Range Gauge (QMS 200).

Cur Val

Current measured value.

Last Val

Measured value of the last cycle.

Delta

Deviation Cur Val - Last Val in %.

Alarm

Switching functions.

Threshold monitoring

The threshold status is displayed graphically in the last column under alarm.

See Chapter 3.2.5.7 Switching functions (Trip), 186.

Example 2

Display Option: Multiple Column

The measure data is displayed in several columns (width selectable under Setup), where each column is one measured cycle.

Measurement Number: 3				MCD airdemo.mcp 13:51:16	MCD airdemo.mcp 13:51:19	MCD airdemo.mcp 13:51:22	MCD airdemo.mcp 13:51:11
Nbr	Type	Ident	Unit				
0	Conc.	Nitrogen		82.269%	82.269%	82.269%	82.269%
1	Conc.	CarbonDiox		80.148ppm	80.152ppm	80.149ppm	80.148ppm
2	Conc.	Argon		192.498ppm	192.500ppm	192.499ppm	192.499ppm
3	Conc.	Oxygen		17.704%	17.704%	17.704%	17.704%
4	Cur.	28.00	A	9.997E-06	9.997E-06	9.997E-06	9.997E-06
5	Cur.	32.00	A	8.318E-06	8.318E-06	8.318E-06	8.318E-06
6	Cur.	18.00	A	1.249E-06	1.249E-06	1.249E-06	1.249E-06

Fig. 2-59

The individual columns are overwritten cyclically from the left to the right:

**Entering
Parameter setup**

For setting the display mode choose [Parameters] > [Setup] (or double click the left mouse button) and for changing the QMS channel parameters choose [Parameters] > [Channel] (or double click the right mouse button).

2.6.8 MCD Bargraph

In MCD Bargraph the measurement is started with the parameter file selected with the Quadstar 32-bit file manager and the measured data is displayed as a bar chart.




Fig. 2-60

In the comments field, the data of the bar marked with the mouse is displayed.

2.6.8.1 Parameter setup

Under [Parameters] > [Setup] you can enter the control parameters:




Fig. 2-61

Control elements

Measurement

- Single: Perform only one measurement (= 1 cycle) across all activated channels.
 - Repeat: Repetitive measurement of the activated channels.
 - Timer: The measurement is periodically restarted.
 - Dynamic: The measurement data is temporarily stored in the QMS buffer.
- See Chapter 2.1.1.2 Measurement methods, 62.

Y-Scale

- Lin: Linear display of measurement data.
- Log: Logarithmic display of measurement data.

Color Mode


- Single Color: The measurement data are displayed in the colors defined in the Parameter Setup program under [Config] > [Screen Color].
- Multiple Color: The measurement data are displayed in the colors that are defined in the channel parameter editor for the corresponding channels. The bars of a channel (normal, last, reference) can only be distinguished by their patterns.

Normalization Mode


Display of the concentrations.

- Sum=100%: Normalize the displayed data to «sum of all concentrations = 100%».
- Max=100%: Normalize the displayed data to «highest concentration = 100%».
- Extern: External normalization. See Section «Setting the external normalization», 205.

Mass scale correction

Mass scale correction on/off. See Chapter 2.8.1 Mass scale calibration (Mass Scale),  134.

Zero Gas Subtraction

Zero gas subtraction on/off. See Chapter 2.8.2 Background measurement (Zero Gas),  137.

Alarm Display

Display the alarm status for all channels with activated threshold monitoring.

Alarm Sound

Enable the alarm buzzer.


Alarm Error Message

Display an alarm additionally as an error.

If the Error History is enabled, incoming threshold monitoring alarms are logged in the Error History file. Configurable in the Parameter Setup program under [Setup] > [General].

Temporary Save

Store the measurement data in a temporary file.

When the measurement is stopped, a dialog box is displayed for entering the file name. See Chapter 2.5.5 Temporary save,  91.

Neg. Conc. Suppr.

Suppress negative concentrations.

Display Last

Display the measured data of the last cycle. The bar is shown on the right of the current measured value.

Display Reference

Display a reference measurement. The bargraph is shown on the right of the current measured value and marked with the pattern chosen by Bar Reference Display. The reference file is selected by the Quadstar 32-bit file manager when the measurement is started.

X-Axis Lettering

X-axis lettering on/off.

Y-Grid

Display a Y-grid on the measurement screen.

Bar Normal Display

Pattern of the bar for the current measured value.

Bar Last Display

Pattern of the bar for the measured value of the last cycle.

Bar Reference Display

Pattern of the bar for the reference measurement.

Bar Width

Bar width in pixels.

Timer Step

Timer interval for the measurement method Timer.

2.6.8.2 Data display

Parameter setup For setting the display mode choose [Parameters] > [Setup] (or double click the left mouse button) and for changing the QMS channel parameters choose [Parameters] > [Channel] (or double click the right mouse button).

Threshold monitoring The threshold status is displayed graphically below the X-axis. See Chapter 3.2.5.7 Switching functions (Trip), 186.

2.6.9 MCD Versus Time

In MCD Versus Time the measurement is started with the parameter file selected in the file manager, and the measurement data are displayed as a function of time or number of cycles:




Fig. 2-62

2.6.9.1 Parameter setup

Under [Parameters] > [Setup] you can enter the control parameters:




Fig. 2-63

Control elements

Measurement

- Single: Perform only one measurement (= 1 cycle) across all activated channels.
 - Repeat: Repetitive measurement of the activated channels.
 - Timer: The measurement is periodically restarted.
 - Dynamic: The measurement data is temporarily stored in the QMS buffer.
- See Chapter 2.1.1.2 Measurement methods, 62.

X-Scale

- Cycle: Assign the cycle number to the X-axis label.
- Time relative: Assign the relative time since the start to the X-axis label.
- Time absolute: Assign the absolute time (time of day) of the measurement to the X-axis label.

Y-Scale

- Lin: Linear display of measurement data.
- Log: Logarithmic display of measurement data.

Line Type

- Solid: The measured values are connected by a solid line.
- Dot: The measured values are shown as individual dots.

Normalization Mode

Indication of the concentrations.

- Sum=100%: Normalize the displayed data to «sum of all concentrations = 100%».
- Max=100%: Normalize the displayed data to «highest concentrations = 100%».
- Extern: External normalization. See Section «Setting the external normalization», 205.

Mass scale correction

Mass scale correction on/off. See Chapter 2.8.1 Mass scale calibration (Mass Scale),
 134.

Zero Gas Subtraction

Zero gas subtraction on/off. See Chapter 2.8.2 Background measurement (Zero Gas),
 137.

Alarm Display

Display the alarm status for all channels with activated threshold monitoring.

Alarm Sound

Enable the alarm buzzer.

Alarm Error Message

Display an alarm additionally as an error.

If the Error History is enabled, incoming threshold monitoring alarms are logged in the Error History file. Configurable in the Parameter Setup program under [Setup] > [General].

Temporary Save

Store the measurement data in a temporary file.

When the measurement is stopped, a dialog box is displayed for entering the file name. See Chapter 2.5.5 Temporary save, 91.

Help table

On the right, beside the measure data, the measured values are displayed additionally to the measure data identifications.

Neg. Conc. Suppr.

Suppress negative concentrations.

X-Grid

Display an X-grid.

Y-Grid

Display a Y-grid.

Marker

Mark the measured values with symbols.

Roll Axis

Each time the end of the X-axis is reached, scroll the screen to the left by 10%.

Cycles

Define the number of cycles that are displayed in the measurement data window before the old data are overwritten. For the «Display after Measurement» method, Cycles is the number of cycles to be measured.

Time-Scale

Select the period of time on the X-axis.

Timer Step

Timer interval for the measurement method Timer.

Measure/Display Mode

- Measure and Display: The values are measured and displayed cycle by cycle.
- Display after Measurement: The values are displayed after the measurement.

2.6.9.2 Data display

Threshold monitoring

Changes of the threshold status are indicated with small arrows above and below the measured values.

The arrows above the measured values signalize threshold state 1 (=A), the arrows below the measured values signalize threshold state 2 (=B). An upward pointing arrow indicates that the threshold monitoring has responded. A downward pointing arrow indicates that the threshold monitoring has returned to normal condition. See Chapter 3.2.5.7 Switching functions (Trip), 186.

Parameter setup

For setting the display mode choose [Parameters] > [Setup] (or double click the left mouse button) and for changing the QMS channel parameters choose [Parameters] > [Channel] (or double click the right mouse button).

2.6.10 MCD last measured values

[MCD] > [Last as Table] and [MCD] > [Last as Bargraph] display the last measured cycle.



Fig. 2-64

File

File handling and return to the main menu of the Measure program.

Display

Select the display mode (Table or Bargraph).

Parameters

Control parameter (display) and measurement parameter handling.

2.6.10.1 File submenu



Fig. 2-65

Load Reference

Load a reference data file.

Save Reference

Store the last measurement (1 cycle) as a MCD reference data file.

Close

Quit the MCD Last Values mode.

2.6.10.2 Display submenu



Fig. 2-66

Table

Display as a table.

Bargraph

Display as a bargraph.

2.6.10.3 Parameters submenu



Fig. 2-67

Setup

Setup the display of the measured data.

2.7 Measuring with the sequencer (Sequence)

All measurements, that can be done manually, can be programmed in a sequence as well. By using decisions, user calls etc., measurement procedures can be programmed, that are adapted to a specific task. This way, dull routine jobs may be automated, and sophisticated measurement procedures can be controlled by less skilled users as well. The sequencer works with the same parameter files as the manual measurements do, and the display of measurement data is the same as well.

While a sequence is running, the only menu entry left is [File] > [Close]. By this entry, the sequence can be cancelled. If an exit sequence has been defined by the SeqExit() instruction, important closing work can be done such as switching off the emission, opening/closing valves, etc.

NOTE:

More information about sequences, programming, instructions etc. can be found in Chapter 8. Process Control, 409.

Sequence menu




Fig. 2-68

Submenus

Execute

Start a sequence.

Setup

Setup the display (digital in-/output, Lin/Log, etc.).

2.7.1 Setup submenu

Under [Sequence] > [Setup] you can set up the display:

QMS 422/421




Fig. 2-69

QMS 200




Fig. 2-70

Lines

Number of lines to display

- 16: 16 lines (for screens with low resolution)
- 32: 32 lines (for screens with high resolution)

DI/DO Identification

Channel name display on/off. See Chapter 2.7.2.2 Status display, 122.

Y-Scale

- Lin: Linear display of measurement data.
- Log: Logarithmic display of measurement data.
- [Display range]: Display range for float and integer values, if no scale has been defined by the SetScale() instruction:
 - Maximum: Largest decade to be displayed.
 - Decades: Number of decades to be displayed.




Fig. 2-71

Selection of DI/DO Boards

- #: Number 0...9 of the entry.
- Unit Name: Logical name of the control unit (in Multiplex mode).
- DI/DO Board: Select the DI/DO board and its channels to be displayed:
 - QMS 422/421: 0...15, 16...31 for [Lines] : [16] or 0...31 for [Lines] : [32]
 - QMS 200: 0...1

NOTE:

Please note that the selected DI/DO channels are displayed only if the display option Disp has not been set to OFF in the sequencer instructions (SetDO(), While-DI() etc.).

2.7.2 Execute submenu

Starting a sequence

By [Sequence] > [Execute] the file manager is called. Choose the desired sequence and start it by [OK].

2.7.2.1 Testing a sequence

Before the sequence is executed, Quadstar 32-bit checks it for programming errors. During the test, the name of the sequence and the current line number are displayed in the status bar (in the lower left corner). In the following example, the sequence *demoair.sqe* could not be found:




Fig. 2-72

However, a complete error test is not possible, if e.g. file names with consecutive numbering (e.g. DEMO00.MIP, DEMO01.MIP, ... DEMOxx.MIP) are used in the sequence. In this case, only a warning is output, that files are used, which could not be tested.

NOTE:

Sequences of this type should be tested by copying and simplifying them. Remove all lines that contain no file access of this type, and set the loop counters for repetitive measures to 1. Now you can start the sequence to check path and file names.

2.7.2.2 Status display

The status of the digital in-/outputs and the loop state is only displayed, if in the corresponding instruction (SetDO(), WhileDI(), WhileMeaVal(), WhileVar() etc.) the display has been switched on by the keyword Disp=on.

NOTE:

Please note that when Disp=on is set, time is required for building up the screen, so the switch-on (e.g. due to a digital input) may be slightly delayed.

Example:

```
Loop( i[0]=0;8 )
  SetDO( i[0]=on;1, Disp=on )
WhileVar( i[5]=0, Disp=on, BreakTime=30 )
Begin
  SetDO( 16=on;5, Disp=on )
  SetDO( 16=off;1, Disp=on )
  Loop( i[0]=0;8 )
```

```
SetDO( i[0]=off;1, Disp=on )
End
```

The example above is a gas inlet control. Start the sequence, so the following is displayed:




Fig. 2-73

DI#X/DO#X

State of the digital in-/output channels. The name of the channel is displayed only if its display has been switched on under [Sequence] > [Setup] by enabling DI/DO Identification. See Chapter 2.7.1 Setup submenu, 120.

While Loop State

Loop state.

NOTE:

The display area is limited. Only as many DI/DO channels are shown, as space is available, even if under [Sequence] > [Setup] more channels have been selected.

2.7.2.3 Display of measurement data

The display of measurement data is the same as in manual mode. If several types of display are possible, the desired type may be chosen by the keyword Disp.

The MID() instruction e.g. has the following possibilities:

- off: Display is switched off.
- tab: Display the measured values as a table.
- bar: Display the measured values as a bargraph.
- vt: Display the measured values as a function of time.
- default: Display according to the mode set by the SetDefault() instruction (off/on, where on=tab).

By the SetPar() instruction, you can do the settings, that are done under [Parameters] > [Setup] in manual mode, for example:

- AlarmDisp: Graphic indication of the threshold monitoring.
- AlarmError: Alarm error messages.
- AlarmSound: Audible alarm signal.
- BarWidth: Bargraph bar width.
- ColorMode: Color display mode.

- Cycles: Number of cycles.
- DispLast: Display Last Values for the bargraph.
- DispOpt: Display option for table output.
- LineType: Line type.
- Marker: Marker display.
- MSC: Mass scale correction.
- NorMode: Normalization mode.
- RelVal: Display of relative values.
- RollAxis: Axis roll mode.
- Scale: Axis mode.
- ScanArran: Display mode, if more than one Scan Analog or Scan Bargraph spectrum is to be displayed.
- TimeScale: Time scale extent.
- XLetter: X-axis lettering.
- XRaster: Display a X-grid on the measurement screen.
- XScale: X-scale display.
- YRaster: Display a Y-grid on the measurement screen.
- ZeroEq: Zero equalization.
- ZeroSub: Zero gas subtraction.

Example:

In the Parameter Setup program, open the sequence editor by [Sequence] > [Editor] and choose a sequence. Set the mouse cursor to a blank line, type SetPar and press the right mouse button. The following dialog box appears that contains all possible parameters of the SetPar() instruction:




Fig. 2-74

If you check the box on the left of a parameter (e.g. x AlarmDisp), then you can click on the button and enter its value (e.g. on).

2.7.3 Controlling the display via sequencer instructions

2.7.3.1 Measurement data of the MID() instruction

The following sequence performs the same MID measurement three times and displays the results in different ways:

```
SetPar(
BarWidth=25,ColorMode=single,Cycles=50,DispOpt=last,LineType=so
lid, Marker =0 )
Loop( i[0]=0;50 )
Begin
  MID( Par="d:\qs32bit\par\airdemo.mip", Disp=tab, SaveGfa=0 )
  ..MID( Par="d:\qs32bit\par\airdemo.mip", Disp=bar, SaveGfa=0 )
  ..MID( Par="d:\qs32bit\par\airdemo.mip", Disp=vt, SaveGfa=0 )
  SetString(gs[0] = «MID measure cycle @, intensity of channel 0
= @";i[0])
  Message( Text=gs[0];gfa[0][0] )
End
```

- The line containing the SetPar() instruction sets the parameters for the display in mode Table (tab), Bargraph (bar) and Versus Time (vt).
- The lines containing the MID() instructions control the measurement of the data and the display in the different modes Table (tab), Bargraph (bar) and Versus Time (vt).
- The SetString() instruction loads a text into the string table.
- The Message() instruction outputs the string on the status bar (on the bottom left).

On the following pages you can see the measurement results of this sequence.

Specialties of the MID/MCD measurements

In sequencer mode as well, the background should be taken into consideration. If you have done a correct background measurement, switch on the background subtraction by the SetPar() instruction (ZeroSub=on); otherwise, switch it off by (ZeroSub=off). See Chapter 2.8.2 Background measurement (Zero Gas), 137.

In the same way as in manual mode, an offset alignment (Zero Equalize) is performed before every MID/MCD measurement, if you switch it on by the SetPar() instruction (ZeroEq=on). The alignment is not done, if you switch it off by the SetPar() instruction (ZeroEq=off). See Chapter 2.5.3 Offset alignment, 90.

MID Table measurement data

The sequencer line

```
MID (Par="d:\qs32bit\par\airdemo.mip", Disp=tab, SaveGfa=0 )
```

stores the measured values in array 0 of the global field (gfa[0][i]) and produces the following table:




Fig. 2-75

MID Bargraph measurement data

The sequencer line

```
MID (Par="d:\qs32bit\run\par\airdemo.mip", Disp=bar, SaveGfa=0 )
```

stores the measured values in array 0 of the global field (gfa[0][i]) and produces the following bar chart:




Fig. 2-76

MID Versus Time measurement data

Here the measured data are displayed as a function of time. The extent of the X-axis can be specified at will, but it should be adapted to the measure loop. This display mode is only useful if the MID() instruction is used within a loop.

The sequencer line

```
MID (Par="d:\qs32bit\run\par\airdemo.mip", Disp=vt, SaveGfa=0 )
```

stores the measured values in array 0 of the global field (gfa[0][i]) and produces the following picture:




Fig. 2-77

2.7.3.2 Scan measurement data

Scan Analog and Scan Bargraph display the measurement data of one or several scans. The number of displayed scans depends on the number of activated channels (State=Enable) in the selected parameter file. If several scans are available, the display mode can be chosen by the SetPar() instruction by the keyword ScanArran:

- Vert: The individual spectra are displayed vertically, i.e. above each other, or
- Horiz: The individual spectra are displayed horizontally, i.e. side by side.

NOTE:

If the SetPar() parameter ScanArran is not set explicitly, it defaults to the last setting chosen manually under [Parameters] > [Setup] > [Scan Arrangement].

Measurement data of the ScanAnalog() instruction

The following sequencer lines execute the three Scan Analog measurements of the *airdemo.sap* parameter file 25 times, while the current cycle number is shown in the status bar:

```

....
Loop( i[5]=1;25 )
Begin
  Message( Text="Scan Measure, Scan Analog # @";i[5] )
  ScanAnalog( Par="airdemo.sap", Disp=on )
End
....

```

This sequence produces the following picture:




Fig. 2-78

Measurement data of the ScanBar() instruction

The following sequencer lines execute the three Scan Bargraph measurements of the airdemo.sbp parameter file 25 times, while the current cycle number is shown in the status bar:

```

Loop( i[5]=1;25 )
Begin
  Message( Text="Scan Measure, Scan Bargraph # @";i[5] )
  ScanBar( Par="airdemo.sbp", Disp=on, Prot=off )
End
....
    
```

This sequence produces the following picture:




Fig. 2-79

2.7.3.3 Measurement data of the MCD() instruction

The following sequence executes the same MCD measurement three times and shows the results in three different ways:


```
SetPar( AlarmDisp=on, BarWidth=25, ColorMode=single, Cycles=50,
DispOpt=last,                                     LineType=solid,
Marker=on, Scale=log, YRaster=off )
Loop( i[0]=0;50 )
Begin
  MCD( Par="d:\qs32bit\par\airdemo.mcp", Disp=tab, SaveGfa=0 )
//MCD( Par="d:\qs32bit\par\airdemo.mcp", Disp=bar, SaveGfa=0 )
//MCD( Par="d:\qs32bit\par\airdemo.mcp", Disp=vt, SaveGfa=0 )
  SetString( gs[0] = "MCD measure cycle @, Concentration N2 = @
%" ;i[0] )
  Message( Text=gs[0];gfa[0][0] )
End
```

MCD Table measurement data

The sequencer line

```
MCD ( Par="d:\qs32bit\par\airdemo.mcp", Disp=tab, SaveGfa=0 )
```

stores the measured value in array 0 of the global field (gfa[0][i]) and produces the following table:



Nbr	Type	Ident	Cur Val	Last Val	Delta	Alarm
0	Conc.	Nitrogen	86.880%	87.103%	-0.257%	<input type="checkbox"/>
1	Conc.	CarbonDiox	-177.783ppm	-160.270ppm	9.851%	<input checked="" type="checkbox"/>
2	Conc.	Argon	39.482ppm	74.661ppm	-89.102%	<input checked="" type="checkbox"/>
3	Conc.	Oxygen	13.134%	12.906%	1.739%	<input checked="" type="checkbox"/>
4	Cur.	28.00	-1.035E-09A	-1.101E-09A	-6.324%	<input type="checkbox"/>
5	Cur.	32.00	-6.052E-10A	-6.307E-10A	-4.208%	<input type="checkbox"/>
6	Cur.	18.00	1.249E-06A	1.249E-06A	-0.000%	<input type="checkbox"/>

MCD measure cycle 1, Concentration N2 = 86.8799 %

Fig. 2-80

MCD Bargraph measurement data

The sequencer line

```
MCD ( Par="d:\qs32bit\par\airdemo.mcp", Disp=bar, SaveGfa=0 )
```

stores the measured values in array 0 of the global field (gfa[0][i]) and produces the following bar chart:




Fig. 2-81

MCD Versus Time measurement data

Here the measured data are displayed as a function of time. The extent of the X-axis can be specified at will, but it should be adapted to the measure loop. This display mode is only useful if the MCD() instruction is used within a loop.

The sequencer line

```
MCD ( Par="d:\qs32bit\par\airdemo.mcp", Disp=vt, SaveGfa=0 )
```

stores the measured values in array 0 of the global field (gfa[0][i]) and produces the following picture:





Fig. 2-82

2.7.3.4 Measurement data of the GSC() instruction (gas specific calibration)

The sequencer line

```
GSC( Par="d:\qs32bit\run\par\airdemo.gcp", Disp=on, Prot=off )
```

produces the following screen:



Nbr	Mass	Component	Concentration	Intensity	Calib.Factor
0	14.00	Nitrogen	78.0850%	-9.732E-10A	1.143E+00
1	28.00	Nitrogen	78.0850%	-8.513E-10A	* 1.000E+00
2	29.00	Nitrogen	78.0850%	4.103E-07A	-4.819E+02
3	32.00	Oxygen	20.9500%	-8.813E-10A	3.859E+00
4	40.00	Argon	9300.0ppm	-2.592E-11A	2.557E+00
5	44.00	CarbonDiox	300.00ppm	1.559E-10A	-4.766E+02

Fig. 2-83

The table shows the component names, the mass numbers, the concentrations, and the determined calibration factors. The line marked with the symbol «*» is the Internal Standard.

NOTE:


More information about this measurement is found in Chapter 2.8.4 Gas specific sensitivity, 142.

2.7.3.5 Measurement data of the ZeroGas() instruction (zero gas measurement)

The sequencer line

```
ZeroGas( Par="d:\qs32bit\run\par\airdemo.mip", Disp=on )
```

produces the following table:



Nbr	Mass	Type	Intensity
0	14.00	Cur.	8.316E-06A
1	16.00	Cur.	2.487E-06A
2	18.00	Cur.	1.249E-06A
3	28.00	Cur.	9.877E-06A
4	32.00	Cur.	8.316E-06A
5	40.00	Cur.	1.668E-06A
6	44.00	Cur.	6.300E-07A
7	****	AI	*****
8	****	PI	*****

Fig. 2-84

The measurement data are stored in Quadstar 32-bit and are subtracted in all subsequent measurements from the measured intensities, if the background subtraction is switched on.

NOTE:


More information about this measurement is found in Chapter 2.8.2 Background measurement (Zero Gas), 137.

2.7.3.6 Measurement data of the MSC() instruction (mass scale calibration)

The sequencer line

```
MSC(Par="d:\qs32bit\run\par\airdemo.msp", CalMode=fine,
SaveMode=reset, Disp=on)
```

produces the following table:



Nbr	Old Mass	New Mass	Intensity	State
0	14.00	14.00	8.317E-06A	ok
1	16.00	16.00	2.489E-06A	ok
2	18.00	17.98	1.249E-06A	ok
3	28.00	28.00	9.997E-06A	ok
4	32.00	32.00	8.317E-06A	ok
5	40.00	40.00	1.669E-06A	ok
6	44.00	44.00	6.300E-07A	ok

Fig. 2-85

Please note the significance of the keyword SaveMode.

- reset: The previous data of the mass scale calibration are deleted and replaced by the new data.
- append: The new mass scale calibration data are inserted into the previous data. Already existing basic points are overwritten and new ones are appended.

Calibration of the mass scale in CalMode=fine should only be performed, if it was previously possible to perform at least one (successful) calibration in coarse mode.

NOTE:

More information about this measurement is found in Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

2.7.4 User-controlled sequence

Quadstar 32-bit offers the possibility of programming a user-controlled sequence by the Dialog() instruction. The Dialog() instruction produces a Quadstar 32-bit dialog box that enables the operator to influence the execution of a sequence. For example, the dialog box could contain two pushbuttons that are labeled with «0... 100 amu» and «100...200 amu». The operator would then choose, whether he wants to measure in the range of 0... 100 amu or 100...200 amu.

See Chapter 8.6.2 User-controlled analysis (AIRDEMO2),  431.

2.8 Calibration

Different calibrations are necessary before a gas analysis is performed for the first time:

- Mass scale calibration
- Background measurement
- QMS Offset calibration (QMS 422/200)
- Gas specific sensitivity calibration.

The calibrations can be performed manually or in the sequencer.

Calibration menu



Fig. 2-86

Mass Scale

Mass scale calibration.

Zero Gas

Background measurement. See Chapter 2.8.2 Background measurement (Zero Gas), 137.

Gas Specific Sensitivity

Gas specific sensitivity calibration. See Chapter 2.8.4 Gas specific sensitivity, 142.

QMS Offset

Determine the QMS offset. See Chapter 2.8.3 QMS Offset, 140.

2.8.1 Mass scale calibration (Mass Scale)

Why calibrate the mass scale?

For exact MID or MCD measurements it is vital to measure at the peak maximum. For that purpose it is necessary to calibrate the mass scale. In this way, deviations of the mass scale within the range of ± 0.5 amu from the nominal value can be corrected.

Larger deviations must be corrected before:

- QMS 422/421: Perform a rough alignment on the RF generator of the spectrometer.
- QMS 200: Rough alignment in the Tune Up program by [Tune] > [QMS 200 Tune Mass Scale].

Under [Calibration] > [Mass Scale] the peak maxima of a few defined masses are accurately determined and the positions of the remaining masses are interpolated or extrapolated linearly (for extrapolation only the first and the last basic point are used).

The calibration data that is obtained this way can be checked or modified in the Display Saved Values program under [Auxiliary] > [Mass Scale Calibration].

Application: In MCD measurements only an integer value (14, 16, 28, etc.) is inserted for the mass number and the correction value is calculated from the calibration data. This method ensures that when the mass scale is recalibrated, the mass numbers of previously prepared parameter files are also updated.

Procedure

- 1 Let a suitable calibration gas into the vacuum chamber. If no specific calibration gas is available, you can work with air instead.
- 2 Choose [Calibration] > [Mass Scale] and open the matching file (e.g. *airdemo.msp* if you work with air).
- 3 In the now appearing setup dialog box choose [Calibration Mode] > [Coarse]. If you want to replace the previous mass scale calibration data, then choose [Save Mode] > [Reset], otherwise [Append]. Confirm the choice by [OK], so the calibration is started:

Nbr	Old Mass	New Mass	Intensity	State
0	14.00	14.00	8.315E-07A	ok
1	16.00	16.00	2.486E-07A	ok
2	18.00	17.98	1.247E-07A	ok
3	28.00	28.00	9.997E-06A	ok
4	32.00	32.00	8.315E-07A	ok
5	40.00	40.00	1.667E-07A	ok
6	44.00	43.98	6.300E-07A	ok

Fig. 2-87

The calibration produces the table above. The State column shows the status of the individual channels:

State column

- ok: The calibration has been successfully performed, the exact mass numbers are shown in the column New Mass.
 - no peak found: No peak was found within the search range.
 - mass too low: The entered mass number is too far to the left of the peak.
 - mass too high: The entered mass number is too far to the right of the peak.
 - int. below threshold: Intensity lower than the level defined by the Threshold parameter.
 - int. too high: Intensity equal to or higher than full scale range.
 - int. drop insufficient: Intensity is too constant to determine a maximum.
 - peak too narrow: The peak is too narrow to determine a maximum.
- 4 If all mass numbers have been successfully calibrated, perform the calibration in [Calibration Mode] > [Fine] as well; the search range is smaller and the accuracy even higher.

If many errors have been reported, you should check the peaks and the measure parameters by a Scan Analog measurement (measure range, peak position, peak shape etc.). Choose [Parameters] > [Channel] or double-click with the right mouse button on the table, so the channel parameter editor appears for optimizing the measurement parameters.

NOTE:

The parameter optimization is explained in the Parameter Setup program. See Chapter 3.4.1 Mass Scale, 199

In the Display Saved Values program (under [Auxiliary] > [Mass Scale Calibration]) the data of the mass scale calibration can be checked and corrected if necessary.

Mass Scale Calibration menu




Fig. 2-88

Submenus

File

File handling and return to the main menu of the Measure program.

Operation

Repeat the measurement.

Parameters

Parameter handling for control and measurement parameters.

2.8.1.1

File submenu



Fig. 2-89

Load Parameter

Load a parameter file.

Close

Quit the mass scale calibration.

2.8.1.2

Operation submenu



Fig. 2-90

Remeasure

Repeat the measurement.

2.8.1.3

Parameters submenu



Fig. 2-91

Setup

Enter control parameters.

Channel

Enter channel parameters.


2.8.1.4**Setup**

Fig. 2-92

Calibration Mode

- Coarse: The QMS works with a fast, coarse search. The search range is wide and the result fairly exact.
- Fine: The peak top is searched within a reduced range, the result is very exact.

Save Mode

- Reset: The previous data of the mass scale calibration will be replaced by the new data.
- Append: The new data of the mass scale calibration are inserted into the previous. Already existing basic points are replaced and new ones are inserted.

2.8.1.5**Channel**

Change or optimize channel parameters.

2.8.2**Background measurement (Zero Gas)****Why background measurement ?**

The residual gas background in the analysis chamber can falsify the results if it is not handled properly. The partial pressures and the composition of the residual gas in the analysis chamber are essentially determined by:

- type of vacuum pump (turbomolecular pump, diffusion pump, etc.), and
- history of the system.

Zero Gas enables you to determine the residual gas background that is found in any analysis chamber. The ion currents determined this way can be subtracted in all subsequent measurements. The residual gas is measured either with a closed gas inlet or with the aid of a zero gas (ultra-pure gas that is admitted into the measurement chamber at the same pressure that is subsequently used for analysis). Through the admission of a zero gas for residual gas determination, adsorption and desorption effects are better taken into consideration than if the chamber is pumped down to the ultimate vacuum.

Mass numbers on which peaks of the zero gas occur must not be included in the measurement parameter list.

NOTE:

The values of a background measurement are stored and subtracted from the intensities measured in subsequent measurements (MID, MCD or GSC measurements), if the Zero Gas Subtraction is switched on.

The stored data of the background measurement can be checked and corrected if necessary in the Display Saved Values program under [Auxiliary] > [Zero Gas].

Zero Gas measurement with SEM

The calibration and measurement operations are greatly simplified if the Common SEM Voltage option is used. This ensures that the selected SEM voltage is the same for the current background measurement and all subsequent analyses.

NOTE:

The Common SEM Voltage is entered under [Setup] > [SEM/Emission Control] (or by Ctrl+S).

Assign the Common SEM Voltage to a channel: Click (in the channel parameter editor) on the SEM Voltage of the corresponding channel with the right mouse button. Then click with the right mouse button on either of the arrows [up] or [down]. The Voltage is now displayed between brackets (e.g. <<2500>>); this indicates that it is always set to the current Common SEM voltage.

Procedure

- 1 Let a suitable zero gas into the vacuum chamber or pump it down to the ultimate pressure.
- 2 Choose [Calibration] > [Zero Gas] and open the matching parameter file; this is normally the MID or MCD file that you use later for the analysis (e.g. *airdemo.mip* if you analyze air).
- 3 If you want to replace the previous zero gas data, then choose [Save Mode] > [Reset], otherwise [Append]. Confirm the choice by [OK], so the zero gas measurement is started. On all mass numbers of the current parameter file the intensities are measured and displayed:

Nbr	Mass	Type	Intensity
0	14.00	Cur.	8.315E-06A
1	16.00	Cur.	2.487E-06A
2	17.98	Cur.	1.249E-06A
3	28.00	Cur.	9.877E-06A
4	32.00	Cur.	8.316E-06A
5	40.00	Cur.	1.668E-06A
6	43.98	Cur.	6.300E-07A
7	****	AI	*****
8	****	PI	*****

Zero Gas measurement done !

Fig. 2-93

The measurement data are displayed as a table. In the background measurement only the Detector types FARADAY, SEM, CH-TRON and ION COUNT are taken into consideration; for other types (AI, PI, PE, TPR or PKR) no measurement data are collected:

Type

Data type of the displayed measurement value.

- Cur.: Ion current (in Ampere) for Detector Type: FARADAY, SEM or CH-TRON.
- Cnt.R.: Counter reading (in cps) of the ion counter for Detector Type: ION COUNT.

Zero Gas menu



Fig. 2-94

Submenus

File

File handling and return to the main menu of the Measure program.

Operation

Repeat the measurement.

Parameters

Parameter handling for control and measurement parameters.

2.8.2.1

File submenu



Fig. 2-95

Load Parameter

Load a parameter file and start the measurement.

Close

Quit the background measurement.

2.8.2.2

Operation submenu



Fig. 2-96

Remeasure

Repeat the measurement.

2.8.2.3

Parameters submenu



Fig. 2-97

Channel

Modify channel parameters.

2.8.3 QMS Offset

The QMS 422/200 determines all necessary correction values to eliminate offsets of the measure amplifier under different conditions. For that, the ion current is measured on the mass specified under [Mass Parameters] > [Zero Mass] (default=5.5 amu). The offset values are stored in the QMS.

Procedure

- 1 Choose [Calibration] > [QMS Offset]. The current values are read from the QMS and displayed as 'Previous Offset Values'.

NOTE:

Make sure that the measurement window is large enough to display all columns and rows. Otherwise, the display may be mutilated.

- 2 If you want to remeasure the offset values (e.g. if all values are zero, i.e. not determined yet), then choose [Operation] > [Remeasure]; the offset values of all measuring ranges are now measured and displayed. During the measurement, which takes some seconds, the status bar message «QMS Offset Calibration in progress !» is displayed. After measurement the offset values are displayed as 'New Offset Values', and the end of the measurement is confirmed with the status bar message «QMS Offset Calibration done !».

QMS 422

Range	Previous Offset Values SEM	New Offset Values SEM
1E-05	-0.0001E-05	-0.0001E-05
1E-06	0.0001E-06	0.0001E-06
1E-07	-0.0001E-07	-0.0001E-07
1E-08	0.0001E-08	0.0001E-08
1E-09	-0.0001E-09	-0.0001E-09
1E-10	0.0001E-10	0.0001E-10
1E-11	-0.0001E-11	-0.0001E-11
1E-12	0.0001E-12	0.0001E-12

Fig. 2-98

QMS 200

Range	Previous Offset Values			New Offset Values		
	Fast	Normal	Slow	Fast	Normal	Slow
1E-05	0.9930E-05	-0.0217E-05	-0.4483E-05	0.0000E-05	0.0000E-05	0.0000E-05
1E-06	-0.4985E-06	-0.7568E-06	-0.6332E-06	-0.0004E-06	-0.0004E-06	-0.0008E-06
1E-07	-0.2181E-07	-0.0215E-07	0.0014E-07	-0.0110E-07	-0.0110E-07	-0.0120E-07
1E-08	-0.2504E-08	-1.0142E-08	-0.8446E-08	0.0000E-08	0.0000E-08	0.0000E-08
1E-09	-0.1767E-09	-0.5310E-09	0.8250E-09	-0.0005E-09	-0.0005E-09	-0.0008E-09
1E-10	-0.0670E-10	0.1360E-10	0.0215E-10	-0.0110E-10	-0.0110E-10	-0.0120E-10
1E-11	-0.2119E-11	-0.7436E-11	0.8910E-11	-0.1168E-11	-0.0005E-11	-0.0008E-11
1E-12	****	0.2129E-12	0.8088E-12	****	-0.0110E-12	-0.0120E-12

Fig. 2-99

The values are displayed only for check or service purposes. Since they are stored only in the QMS and not in Quadstar 32-bit, it's not possible to display them in the Display Saved Values program.

Fast, Normal, Slow

QMS 200 needs different offset values for the different measurement speeds (Fast, Normal, Slow), whereas QMS 422 has one value for all speeds.

QMS Offset calibration menu

The following functions are available through the QMS Offset Calibration menu:

File

- Close: Terminate the QMS Offset Calibration function.

Operation

- Remeasure: Restart the offset measurement.
- Clear Offset Values: Set all Offset values to zero (only for QMS422 available).

Parameters

- Setup: Call the QMS Offset Calibration Setup dialog box.

2.8.3.1**QMS Offset Calibration Setup****QMS 422**

Fig. 2-100

QMS 200

Fig. 2-101

Control elements**Detector Parameters**

- Type: Type of detector (FARADAY, SEM or CH-TRON).
- SEM Voltage: Voltage of the SEM or CH-TRON.

Mass Parameters

- Zero Mass: Mass number on which the offset is to be measured.
- Resolution: Resolution.

Output Parameters

- Monitor: Service signal at the monitor output of the QMS on/off.

2.8.4 Gas specific sensitivity

The objective of the quantitative gas analysis (MCD measurement) is to determine the concentrations of the individual gas components from the mass spectrum. The measured ion currents are input to a solution matrix and the individual concentrations of the components in the gas to be analyzed are determined via calibration factors.

For calculating the gas concentrations from ion currents, the mass spectrometer sensitivity for the individual gas components must be known. Those relative mass spectrometer sensitivities are determined by the following measurement and stored as calibration factors.

Calibration factor library

The Quadstar 32-bit Quantitative Analysis works with a calibration factor library. The size of this library is open, i.e. it grows with the input of each new component. For each component that is referenced in an application or calibration matrix, the calibration factors and the corresponding mass numbers must be stored in this library. The number of calibration factors is limited to 24 per component.

The calibration factors are stored in a calibration factor library. They are defined as follows:

$$\text{Concentration} = \frac{\text{Intensity}}{\text{CalibrationFactor}}$$

which means that the calibration factors are proportional to the intensities.

Calibration factors can be determined for each fragment ion of a component. If it is certain that on a given mass number only the desired component exists, the concentration can also be determined by measuring only this single mass number of the component. But the accuracy with which the concentration is determined usually increases with the number of measured mass numbers.

NOTE:

Only gas mixtures are to be used for calibration, whose components do not overlap in the mass spectrum.

Since all calibration factors are given as relative values, a fixed calibration factor (internal standard) must be defined. Normally the calibration factor for N₂ (m/e = 28) is set to 1 (default setting), and all other calibration factors relate to this value. In this way a consistent calibration factor library can be built.

Procedure**NOTE:**

Do not forget to measure the background before the calibration, and calibrate the mass scale before if necessary.

- 1 Let a suitable calibration gas into the vacuum chamber.
- 2 Choose [Calibration] > [Gas Specific Sensitivity] and open the parameter file that belongs to the calibration gas in use (e.g. *airdemo.gcp* if you use air for calibration), so the gas specific sensitivity calibration is started:

Nbr	Mass	Component	Concentration	Intensity	Calib.Factor
0	14.00	Nitrogen	78.0850%	8.303E-09A	8.320E-01
1	28.00	Nitrogen	78.0850%	9.980E-09A	* 1.000E+00
2	29.00	Nitrogen	78.0850%	4.083E-10A	4.091E-02
3	32.00	Oxygen	20.9500%	8.304E-09A	3.101E+00
4	40.00	Argon	9300.0ppm	1.652E-09A	1.390E+01
5	44.00	CarbonDiox	300.00ppm	6.274E-10A	1.636E+02

Fig. 2-102

The calculated calibration factors are displayed as a table. The Internal Standard (in this example N₂/28amu) is marked by an * sign. The analysis matrix does not contain the actual calibration factors but only «pointers» to the location in the calibration factor library. The advantage of this method is that always the latest calibration factors are used when the application matrix is loaded. Already existing matrices are «automatically» updated this way.

Menu for Gas Specific Calibration



Fig. 2-103

Submenus**File**

File handling and return to the main menu of the Measure program.

Operation

Restart the measurement.

Parameters

Measurement parameter handling.

2.8.4.1**File submenu**

Fig. 2-104

Load Parameter

Load a parameter file and start the measurement.

Close

Quit the gas specific calibration.

2.8.4.2 Operation submenu



Fig. 2-105

Remeasure

Start the measurement again.

2.8.4.3 Parameters submenu



Fig. 2-106

Channel

Modify/optimize channel parameters.

NOTE:

The optimization of the parameters is described in the Parameter Setup program under [Calibration] > [Gas Specific Sensitivity].

2.9 Manual Parameters

Manual menu

Choose Manual for displaying the state of digital input channels and for manually switching digital output channels on or off. The communication to the mass spectrometer must be enabled for that.



Fig. 2-107

For turning a digital output channel on or off, double-click on the corresponding command button with the left mouse button. The channel state is displayed by the colored dot on the command button: green = on, grey = off.

2.9.1 DI/DO Manual for QMS 422/421




Fig. 2-108

Display left

Select the DI/DO group to be displayed on the left side of the screen.

Display right

Select the DI/DO group to be displayed on the right side of the screen.

A group cannot be displayed at the left and the right side simultaneously.

2.9.2 DI/DO Manual for QMS 200




Fig. 2-109

Only the two standard digital input and output channels are available:

2.10 Communication menu

Under Communication you can turn on or off the interface to the QMS. If several QMS's are available, the desired unit can be chosen.

Comm Menu




Fig. 2-110

Connect

Establish the connection to the QMS.

Disconnect

Cancel the connection to the QMS.

SetUnit

Select a control unit.

See Section «Selecting a control unit in the ArcNet network (Set Unit)», 218.

NOTE:

For detailed information see Chapter 3.6 Communication parameters, 209.

2.11 Setup

In the Setup menu you can turn the emission and the SEM high voltage on/off and setup the documentation of measure data (file info):

Setup menu



Fig. 2-111

Submenus

SEM/Emission Control

Turn the emission and the SEM high voltage on/off. Here you can enter the generally usable Common SEM Voltage as well.

The menu entries SEM and Emission Control are only accessible when the QMS-communication is enabled.

You can also call the dialog box directly by <Ctrl+S>.

File Info

Set up the documentation of measure data. You may here preset the parameter groups to save (e.g. useful for sequencer operation) and determine generally whether to save a documentation or not.




Fig. 2-112

Control elements

Parameters/Data to save

Check the groups here you want to save as a documentation:

- Measurement Parameters: Channel parameters valid at measure time.
- Ion Source Parameters: Ion source parameters valid at measure time.
- State Parameters: Settings that may be important for the interpretation of the measured results.
- Zero Gas Data: Zero gas data valid at measure time.

Additional Information

You may define here up to 20 sets of parameters, each consisting of a Name and a Value.

- [Delete All Names]: Delete self defined parameter names.
- [Delete All Values]: Delete self defined parameter values.

Save File Info

Check here whether to save a measure data documentation at all or not.

[OK]

Adopt the setup and leave the dialog box.

[Cancel]

Drop the setup and leave the dialog box.

3. Parameter Setup

In Parameter Setup program you can enter and check different settings:

- Get general system information.
- Create parameter files for the different measurement types and calibrations.
- Program sequences.
- Enter control parameters of Quadstar 32-bit.
- Setup communication with the QMS.
- Configure the system.

Quadstar 32-bit does not automatically initiate communication with the QMS when the Parameter Setup program is started. This communication must be established by [Comm] > [Connect].

Main menu



Fig. 3-1

Submenus

File

Display the Quadstar 32-bit software version, recall information about the software key and the PC system. General settings for printing. Terminate the Parameter Setup program.

Measure

Create measurement parameters for the Scan Analog, Scan Bargraph, MID and MCD measurement types. See Chapter 3.2 Measure Parameters, 156.

Sequence

Create/edit a sequence. See Chapter 3.3 Sequence programming, 190.

Calibration

Create calibration parameter files. See Chapter 3.4 Calibration Parameters, 199.

Manual

Operation/display of the digital input and output channels (DI/DO). See Chapter 3.5 Manual parameters, 207.

Comm

Select the interface to the mass spectrometer (RS232C, ArcNet) and set it up. See Chapter 3.6 Communication parameters, 209.

Setup

Set up mass spectrometer, printer and ASCII-export. See Chapter 3.7 Setup parameters, 223.

Config

Enter the mass spectrometer configuration (e.g. detector type, mass range, options such as DO 420, PKR 250, etc.), digital and analog inputs and outputs and the screen colors. See Chapter 3.8 Instrument Configuration, [5](#) 234.

Help

Display the 'Key Fragment Ions' table and help contents. See Chapter 1.9 Hints concerning the work with Quadstar 32-bit, [5](#) 56.

3.1 File

File menu




Fig. 3-2

System Info

Display information on the installed hardware and software.

About Parameter Setup

Information on the Parameter Setup program.

Print Settings

General settings for printing.

Exit

Quit the Parameter Setup program.

3.1.1 System Information




Fig. 3-3

Software Key Info/Option File Info

Provides information on the QMS key (if available), the key number, and the enabled options of your Quadstar 32-bit software.

Software Key Info

- Key Number: Identification of the QMS key (e.g.: 1234).

Option File Info

- Key Number: Number of this Quadstar 32-bit software.
- Control Option: Process Control Module installed/not installed.
- Calculation Option: Quantitative Analysis Module installed/not installed.

Free Disk Space

Displays the available storage space on the selected disk. Free Disk Space can be used before a measurement to check whether or not sufficient space is available on the disk drive selected in [Setup] > [General] under Data File Path.

3.1.2 About Parameter Setup

Use this menu option to determine the version of the Quadstar 32-bit Parameter Setup program:




Fig. 3-4

3.1.3 Print Settings

In this dialog box you can choose the settings to be printed:




Fig. 3-5

Global Settings

Settings that apply to all connected units:

- Communication Parameters: Settings of the communication line PC <=> QMS.
- General Setup: Quadstar 32-bit settings (Paths, Error History, Alarm etc.)
- ASCII Conversion Parameters: ASCII-Conversion setup.
- Screen Colors: Colors used for the screen display.

Unit Specific Settings

Settings that apply to a single unit:

- DI Identifications: Named digital inputs.

- DO Identifications: Named digital outputs.
- AI Characteristic Curve: Analog inputs with characteristic curves.
- AO Characteristic Curve: Analog outputs with characteristic curves.
- of [(Unit Specification)]: Units, for which these settings are to be printed:
 - Current Unit: Currently used unit.
 - All Enabled Units: All enabled units of the net.
 - All Units: All declared units.

3.1.4

Exit

Quit the Quadstar 32-bit Parameter Setup program.

3.2 Measure Parameters

Measure menu

In the Measure menu of the Parameter Setup program you can choose the type of measurement and afterwards set up the measurement channels.




Fig. 3-6

The following measurement types are available:

Submenus

Scan Analog

Analog mass spectrum beginning with a start mass (FIRST MASS) across a selectable mass range (WIDTH) on one or more channels.

Scan Bargraph

Peak search, beginning with a start mass (FIRST MASS) across a selectable mass range (WIDTH) on one or more channels. Only the intensities of the found peaks with their corresponding mass numbers are recorded.

MID

Multiple Ion Detection, a localized measurement on one mass number per channel. To get the intensities of several masses, several channels must be used.

MCD

Multiple Concentration Detection, an MID-measurement with subsequent calculation of the concentrations. See Chapter 3.2.5.5 Editing the MCD parameter files (Quantitative Analysis), 180.

3.2.1 Selecting the parameter file

- Select the desired measuring type (Scan Analog, Scan Bargraph, MID, MCD).
- the File Manager appears that shows the existing parameter files. Choose one of them, make a copy of one, or create a new one. See Chapter 10.1 File Manager, 541.

3.2.2 Storing the parameters

When you exit the editor, Quadstar 32-bit stores all parameters in the opened file that has already been assigned a name, and closes it. Each file contains all parameters required for the corresponding measuring type. The file size depends on the number of activated channels (total max. 64 channels). A channel is activated, if its parameter State is set to ENABLE. The different measurement types create parameter files that have different file extensions. For an explanation of file extensions see Chapter 10.3 Quadstar 32-bit files, 544.

3.2.3 Create or copy a parameter file

There are two methods for creating a new parameter file. An existing file can be copied (a) and then modified, or a new file can be created (b) that is initialized with default parameter values. In the latter case all parameters must be adapted, in the first one only a few.

- 1 Choose the desired measurement type. The Quadstar 32-bit file manager shows up:




Fig. 3-7

- 2 a) Copy an existing file to a new one: Select a file and click [Copy].
b) Enter a new file name; in this case, Quadstar 32-bit asks whether to create a new file with this name or not.
- 3 Press the [OK] command button.

3.2.4 The channel parameter editor

When you have opened a parameter file, the channel parameter editor appears in the standard configuration shown below:




Fig. 3-8

- 1 menu
- 2 cursor position
- 3 QMS channel number
- 4 group "Main"
- 5 group "Detector"
- 6 no available Parameter
- 7 Register cards
- 8 set cursor to activ channel
- 9 navigate through register cards
- 10 first channel which is displayed at the active scroll bar position
- 11 channel selection scroll bar

The number of displayed channels depends on the screen resolution, the font and the window size.

The parameters shown in the upper half of the picture (State, Det. Type and Mass) belong to the Main group. This group is always displayed and has no register card.

The other parameters are subdivided into groups which are displayed as register cards. Those groups are by default arranged as follows:

Detector:

Detector properties (e.g. SEM voltage)

Mass

Parameters for a mass scan.

Amplifier

Signal amplifier properties (measure range, offset etc.).

Output

Programming of the analog measurement output.

Trip

Programming of the switching functions (for MID/MCD measurements).

Display

Color and display range of the measured values.

NOTE:

You can define your own groups. These can contain any parameters, for example, the main parameters for a measuring type (MID, Scan Bargraph, etc.). See Chapter 3.2.5.3.1 Configuration, 172.

To changeover the display to a different group (register card), click on the name of the desired group. To display the channel parameters of the Mass group click on the Mass register card:

Load-Ch.:00	CH-0	CH-1	CH-2	CH-3	CH-4	CH-5
State	ENABLE	ENABLE	ENABLE	ENABLE	ENABLE	ENAE
Det. Type	FARADAY	FARADAY	FARADAY	FARADAY	FARADAY	FARA
Mass	14.00	16.00	18.00	28.00	32.00	40.00
Dwell	0.2s	50ms	1s	1s	1s	1s
Resolution	25	25	25	25	25	25
Zero Mass	5.50	5.50	5.50	5.50	5.50	5.50

Navigation bar: Detector \ **Mass** / Amplifier / Output / Trip / Display / First-Ch.:00

Fig. 3-9

Parameters, that have no influence in a channel (due to the selection of higher ranking parameters), are marked by a dash only. Parameters, that are not available (e.g. due to the configuration), are not displayed at all.

3.2.4.1 Measuring parameters overview

Below you will find a list of the most important measuring parameters (also called 'channel parameters') and their meanings. Depending on the mass spectrometer, the measurement type, the hardware configuration etc. the appropriate parameters are displayed. The subdivision into groups (register cards) corresponds to the default setup.

3.2.4.1.1 'Main' group

This group contains the most important parameters, they are displayed constantly:

State

Channel On/Skip/Off (Enable/Skip/Off).

- Enable: The channel is measured and the parameters are stored.
- Skip: The channel is not measured, but its parameters remain stored.
- Off: The channel is not measured and its parameters are deleted.

Det. Type

Detector type that is used for the measurement.

Mass Mode

Type of the mass scan.

- Scan-N: Normal analog scan.
- Scan-F: Analog scan using the FIR filter.
- Stair: Mass spectrum with integer mass jumps.
- Peak-L: Search for peaks, that reach at least a specified intensity.
- Peak-F: Search for peaks, that reach at least a specified intensity; with FIR Filter.

Mass

Mass on which the measurement is performed.

First Mass

Start mass, at which the measurement begins.

3.2.4.1.2 'Detector' group

This group contains the detector parameters:

SEM Voltage

Operating voltage of the SEM.

3.2.4.1.3 'Mass' group

This group contains the parameters, that control the mass scan.

Speed

Measuring time per amu for scan measurements.

Dwell

Dwell time for the measurement on this channel.

Width

Mass range from the start mass onward.

Resolution

Mass resolution (smaller number = better separation of the peaks).

Steps

Number of measured points per amu for scan measurements.

Threshold

Intensity that a peak must reach at least to get displayed.

Zero Mass

Mass on which the electrometer offset is measured.

3.2.4.1.4 'Amplifier' group

This group contains the parameters that control the signal amplifier.

Amp. Mode

Amplifier operation mode (Fix/Auto/Auto-D).

- Fix: The measuring range is unchanging and has to be set.
- Auto: The measuring range is variable, it is determined by the mass spectrometer.
- Auto-D: The measuring range is variable, it is determined by the mass spectrometer, but a lower limit can be set that is not undercut.

Amp. Range

Measuring range of the amplifier.

Calibration

Fine tuning of the amplifier gain.

Range-L

Lower limit of the automatically selectable amplifier ranges.

Filter

Time constant of the built-in analog filter.

Pause-Cal.

Amount of time between the measurement of two channels, during which no values are measured; is used to wait until the amplifier has settled.

Offset

Offset correction value of the amplifier.

Gain

Amplifier gain for detector type Extern.

CP-Level

Threshold of the ion counter preamplifier.

3.2.4.1.5 'Output' group

This group contains the parameters that control the analog output of measured values:

AO-Channel

Number of the analog output, on which the measured values of the channel are output.

AO-Mode

Output the measured values in a linear or logarithmic style on the analog output.

Monitor

Output the measured values in a linear or logarithmic style on the Monitor output.

AO-Decades

Number of decades that are output on the analog output.

AO-Range

Range of measured values that are output on the analog output.

3.2.4.1.6 'Trip' group

This group contains the parameters of the threshold monitoring function:

Trip Type

Type of the threshold monitoring function.

Level A

Threshold A.

DO# A

Digital output, on which the threshold state A or the result of the threshold monitoring function is output.

Level B

Threshold B.

DO# B

Digital output, on which the threshold state B or the result of the threshold monitoring function is output.

3.2.4.1.7 'Display' group

This group contains the parameters of the display control:

Disp. Mode

Display measured values On/Off.

Color

Color, in which the measured values of this channel are to be displayed.

Disp. F.S.R.

Full Scale Range, Upper limit of the display.

Disp. Decades

Number of decades.

3.2.4.1.8 '<<Rest>>' group

This group contains the residual parameters:

Baseline Shift

Value by which the display is shifted up or down, relative to the baseline.

3.2.4.2 Modification of parameters

Click on the desired parameter with the right mouse button, so an input field appears. This can be (depending on the parameter) a context menu (Speed, Range, etc.) or an editing field (Mass, Width, etc.).

Input of parameters via the context menu

No direct input is allowed. In the displayed context menu (and depending on the installed hardware) the alternatives available for the selected parameter are displayed. In the following example these are the parameters Faraday and AI.

With the right mouse button click on the detector type (Faraday) under CH-1:

Load-Ch.:01	CH-0	CH-1	CH-2
State	ENABLE	ENABLE	ENABLE
Det. Type	FARADAY	FARADAY	FARADAY
Mass	14.00	16.00	18.00
Dwell	50ms	50ms	50ms
Resolution	50	50	50

Fig. 3-10

The arrow behind AI (=Analog Input) means that for this parameter a subgroup with additional parameters (here: analog input channels) is available.

These can be called on the screen by clicking on AI with the left mouse button:

Load-Ch.:01	CH-0	CH-1	CH-2
State	ENABLE	ENABLE	ENABLE
Det. Type	FARADAY	FARADAY	FARADAY
Mass	14.00	16.00	18.00
Dwell	50ms	50ms	50ms
Resolution	50	50	50

Fig. 3-11

If logical names are assigned to the analog and digital input and output channels, these are shown in brackets behind the channel number:

FARADAY	
AI	0 temp.
	1 press.

Fig. 3-12

Numeric input via the editing field

In contrast to the context menu, clicking on Mass with the right mouse button opens an editing field in which direct input is possible.

Load-Ch.:01	CH-0	CH-1	CH-2
State	ENABLE	ENABLE	ENABLE
Det. Type	FARADAY	FARADAY	FARADAY
Mass	14.00	16.00	18.00
Dwell	50ms	50ms	50ms
Resolution	50	50	50

Fig. 3-13

The parameter value can be typed or changed by the [UP] or [DOWN] arrow buttons. Accept the change by clicking the button on the left side of the editing field or reject it by ESC.

The [UP] or [DOWN] buttons have the following function:

- Left mouse button + [UP] or [DOWN]: current value ± small steps

- Right mouse button + [UP] or [DOWN]: current value = Default value
- Right + left mouse button + [UP] or [DOWN]: current value \pm large steps
For this input, the right mouse button must be pressed while the cursor is still located in the numeric input field. Then while holding the right mouse button shift the cursor to the [UP] or [DOWN] button and click on it with the left mouse button while still holding the right mouse button.

Special cases

For the following parameters a special effect can be achieved with the right mouse button:

Parameter	Condition	Display	Effect
SEM Voltage	Detector = SEM	<<xxxx>>	SEM Voltage = Common SEM
AO-Channel	---	OFF	No AO channel active
DO# A/DO# B	MID measurement	OFF	No DO channel active
Zero Mass	MID measurement	OFF	Offset calibration disabled
Disp. F.S.R.	---	Auto	Select the display range automatically

Fig. 3-14

3.2.4.3 Configurations (Parameter groups/register cards)

Separate configurations exist for each application of the editor (MID, ScanAnalog, Peak Search parameters, etc.) with differently arranged register cards (Groups). The group names remain the same, but they may contain different parameters. You can change this groups or create new groups, if you want to:

Configuring the editor




Fig. 3-15

- 1 Choose configuration
- 2 Current register card
- 3 Register cards of this configuration
- 4 Parameters not available for the current hardware configuration
- 5 All available Parameters
- 6 Parameters of the current register card
- 7 Buttons either for moving parameters within "In Group" or for transferring parameters between "In Group" and possible

In this dialog box you can define new editor configurations. Each configuration can contain up to 10 register cards that you can put together individually.

NOTE:

Please note that parameters, that are not available in the current hardware configuration, are set between brackets in the Possible (and In Group) fields (e.g. QMU Channel). They can be transferred from the Possible field to the In Group field, but they are not displayed there. See Chapter 3.2.5.3.1 Configuration, 172)

Main Parameters

The State, Det. Type and Mass parameters shown in the upper part of the editor are put by default in the << Main Param. >> register card in the << Default >> configuration. See Chapter 3.2.5.3.1 Configuration, 172).

This group is always displayed. It can be modified and expanded as well.

State

Channel state (on/skip/off).




Fig. 3-16

The operating states ENABLE, SKIP and OFF have the following meaning:

- **ENABLE:** The channel is activated and is taken into consideration in a measurement. The parameters entered in the parameter editor are stored.
- **SKIP:** The channel is not executed in the measurement (skipped). Its parameters are preserved and are also stored. The channel can be reactivated by ENABLE.
- **OFF:** The channel is deleted from the parameter file after saving. The parameters must be entered again when the channel is reactivated.

Det. Type

Chosen detector (FARADAY, SEM, ION COUNT, etc.). The detector type must be defined under [Config] > [QMS...]

Mass Mode

Scan mode (PEAK-N, SCAN-N, STAIR), FIR-Filter (PEAK-F, SCAN-F). The modes that can be selected depend on the chosen measuring type (Scan Analog, Scan Bargraph, MID, MCD).

Mass

Mass number on which the intensity is to be measured.

Configuration example

You want to configure the channel parameter editor in such a way, that it contains only the most important parameters for a Scan Analog measurement.

- 1 Choose [Measure] > [Scan Analog...]:




Fig. 3-17

- Open any file in the now appearing File Manager (e.g. AIRDEMO.SAP).

Create a new configuration (called ‘QMS 200 Scan Analog’)

- Choose [Display] > [Configuration...]
- Click on the [NEW] command button in the Defined Configurations field, enter the name of the configuration to be created (here: ‘QMS 200 Scan Analog’) and confirm your input with [OK]:




Fig. 3-18

NOTE:

The configurations shown in the Defined Configuration list box are only available for Scan Analog. Consequently the new configuration QMS 200 Scan Analog is only available when Scan Analog parameter files are edited.

Setup the group “<< Main Param. >>”

- In the Groups field select the (empty) group name << Main Param. >>, select the parameters State, Det. Type, Mass Mode and Color; press [RIGHT] each time. This transfers the corresponding parameters into the field In Group, i.e. they are now found in the group << Main Param. >>.




Fig. 3-19

Create the register card “Scan”

- Click on the [New] key in the Groups field, enter the group name (‘Scan’) and confirm the entry by [OK].



Fig. 3-20

Adapt the “Scan” register card

- 7 In the Groups field select the new group name Scan.
- 8 Select the desired parameters under Parameters in the Possible field
=> Please note that when the <CTRL> and <SHIFT> keys are pressed, several parameters can be selected and shifted jointly (in both directions).
- 9 With the [RIGHT] button shift these into the In Group field:




Fig. 3-21

- 10 Sort the parameters as desired by [UP] and [DOWN]
- 11 quit the configuration editor by [OK].

The editor now shows the following picture:




Fig. 3-22

3.2.5 Channel Parameter Editor menu

The editor provides the following menu and shows the name of the opened file (here: 'airdemo.sap' in the title bar.



Fig. 3-23

Submenus

File

Call the File Manager, print the parameters and return to the main menu of the program. See Chapter 3.2.5.1 File submenu, 169.

Edit

Insert and delete channels. Sorting and copying functions. See Chapter 3.2.5.2 Edit submenu, 170.

Display

Set up the display options in the parameter editor. See Chapter 3.2.5.3 Display submenu, 171.

Functions

Adjust and Peak Search functions. See Chapter 3.2.5.4 Functions submenu, 174.

3.2.5.1 File submenu




Fig. 3-24

New

Create a new file.

Open...

Call the File Manager for opening a parameter file.

Save

Store the active parameter file under its current name.

Save As...

Store the active parameter file under a different name and/or in a different directory.

Print...

Print the active channel parameter file as a table.

Close

Return to the Parameter Setup main menu.

Print...

Print opens a dialog box in which additional information for the printout can be specified.




Fig. 3-25

print QMS configuration

Print the current QMS configuration on the first page.

column headers on all pages

Print channel headers on all subsequent vertically arranged pages (Enabled)/only on the top page (disabled).

row headers on all pages

Print the corresponding parameter name on all subsequent horizontally arranged pages (enabled)/only on the first page (disabled).

The last two options should be disabled for printouts on which the pages are subsequently pasted together.

The font and the paper orientation is taken from the current printer setup and the printer font dialog box. See Chapter 3.7.3.1 Configuring the printer, 229, Chapter 3.7.3.2 Set up the printer page (Page Setup), 230 and Chapter 3.7.3.3 Select the printer font, 230.


3.2.5.2**Edit submenu**

Fig. 3-26

Insert Channel

Insert a new measurement channel; shift all channels (including the selected channel) one position to the right, starting at the selected channel.

Remove Channel

Remove (delete) the selected channel; shift all channels one position to the left, starting at the selected channel.

Compress Channels

Compress all channels that are not in the 'OFF' State.

Copy/Exchange...

Copy individual parameters or entire channels.

Sort Mass

Sort the channels by mass numbers (First/Mass).

Sort Range

Sort the channels by their amplifier ranges.

Copy/Exchange

Copy opens a dialog box in which you can copy individual parameters or all parameters of a channel:




Fig. 3-27

Mode Selection

Select the desired action:

- Copy: Copy the object chosen under Object Selection to one or to all channels.
- Exchange: Swap all parameters of the channel selected under Source with those of the channel selected under Destination. The Exchange function is only accessible if [Object Selection] > [Channel] has been chosen.

Object Selection

Select the object to be copied or exchanged:

- Channel Copy or exchange all parameters of the channel.
- Parameter Copy or exchange only the selected parameter.

The channels are selected under Source and Dest:

Source

Enter the channel from which the parameters are to be copied or exchanged.

Dest

Enter the destination channel to which the parameters are to be copied or exchanged. To specify ALL (copying/exchanging all channels) click on [UP] or [DOWN] with the right mouse button.

Special cases

The following special cases apply to the copying of channel parameters:

- Threshold Parameters: are only copied to channels for which the same Detector has been selected under Det. Type.
- Range: If a CLA400 (chopper lock-in amplifier) option is installed in the QMS, Parameters are only copied to channels with identical Lock-In state. Without CLA400, there is no restriction.

3.2.5.3

Display submenu



Fig. 3-28

Configuration

Create new configurations and change existing configurations of the selected measurement type (Scan Analog, MID, etc.).

Font

Choose the font type, style and size to be used for the parameter editor.

Measure Data

Start the service measurement (only if the communication to the QMS is switched on). See Chapter 3.2.5.3.3 Measure Data, 174.

3.2.5.3.1 Configuration




Fig. 3-29

Control elements Defined Configuration

List of available configurations.

- [Save As]: Save the current settings and register cards as a new configuration:




Fig. 3-30

- [New]: Create a new configuration.




Fig. 3-31

- [Delete]: Delete the currently selected configuration:



Fig. 3-32

Groups

Groups shows the current register cards (parameter groups).

- [UP]/[DOWN]: With these keys you can change the order of the register cards.
- [Rename]: Rename a register card.
- [New]: Create a new register card.
- [Delete]: Delete the selected register card. The <<Main Param.>> group name cannot be deleted.

Parameters

Setup the register card.

- Possible: Lists the parameters that can be assigned to the register card selected under Groups. => Those parameters which are not available in the current hardware configuration are shown in brackets. They can be transferred to the In Group field, but they will not be displayed by the editor.
- In Group: Lists the parameters that have already been assigned to the register card selected under Groups. These parameters are no longer listed under Possible.
- [UP]/[DOWN]: With these keys you can change the order of the parameters within a register card.
- [LEFT]: Remove the selected parameter from this group.
- [RIGHT]: Move the selected parameter to this group.

3.2.5.3.2 Font selection

Here you can choose the font type, the font style and the font size of the channel parameter editor.

The selected font is shown in the Sample box:




Fig. 3-33

3.2.5.3.3 Measure Data

If a QMS is connected and the communication is on, Measure Data starts a service measurement with the parameters in the currently selected channel.

The display responds to all parameter changes; the type of display depends on the previously selected measuring type (Bargraph (MID), Scan Bargraph or Scan Analog).

If no measurement data are displayed, this may (among others), be caused by one of these reasons:

- Emission OFF
- SEM Voltage OFF
- Range or SEM Voltage setting too low
- Wrong detector (Det. Type parameter in << Main Parameter >> group)
- Channel State OFF or SKIP
- Inlet system closed, that is, the working pressure is too low (recommended: 2.0E-06 mbar).

3.2.5.4 Functions submenu

This menu is only available for MID/MCD parameter files. Adjust and Peak Search are only enabled if the interface to the mass spectrometer is switched on:



Fig. 3-34

Adjust

Adjust the mass number entered under First/Mass to the peak maximum.


Peak Search

Search the highest peaks (selectable number) and store their mass numbers in the channels.

3.2.5.4.1 Adjust function

The Adjust function is used for optimizing the mass numbers entered under Mass to the peak maximum. In Adjust the Range setting can also be tested and optimized.

Adjust opens a window containing an empty table for the adjust data:



Chan	Mass	New Mass	Intensity	State	Comment
0	13.98	-	-	-	-
1	15.98	-	-	-	-
2	17.98	-	-	-	-
3	27.98	-	-	-	-
4	31.98	-	-	-	-
5	39.98	-	-	-	-
6	43.98	-	-	-	-

Fig. 3-35

The Adjust function is controlled via the following main menu:

File

Quit the adjust function.

Adjust

Start and stop the adjust function.

Options

Select the adjust mode.

Edit

Edit the adjust parameters.

File submenu



Fig. 3-36

Close

Return to the main menu of the Channel Parameter Editor. The optimized mass numbers are transferred into the corresponding channels of the editor when an Adjust has been performed successfully.

Adjust submenu



Fig. 3-37

Start

Start the Adjust function.

Stop

Stop the running Adjust function.

The result is displayed as a table:

 A screenshot of a window titled 'Adjust < airdemo.mip >'. The window contains a table with the following data:

Chan	Mass	New Mass	Intensity	State	Comment
0	13.98	14.02	8.3069E-06	000000	ok
1	15.98	16.02	2.4853E-06	000000	ok
2	17.98	18.02	1.2469E-06	000000	ok
3	27.98	28.02	9.9844E-06	000000	ok
4	31.98	32.02	8.3069E-06	000000	ok
5	39.98	40.02	1.6666E-06	000000	ok
6	43.98	44.02	6.2813E-07	000000	ok

Fig. 3-38

The Comment column shows the status of the individual channels:

Comment column display

- ok: The Adjust function was successfully completed and the mass number in the Mass column has been corrected to the value under New Mass.
- no peak found: No peak was found within the search range.
- mass too low: The Adjust function was not completed successfully because the entered mass number is too far to the left of the peak.
- mass too high: The Adjust function was not completed successfully because the entered mass number is too far to the right of the peak.
- intensity too low: Intensity lower than the level specified by the Threshold parameter.
- intensity too high: Intensity greater than or equal to full scale range.

Under [Edit] > [Adjust Parameter], the Mass, Range, and Threshold parameters can be modified.

Options submenu

This menu controls the Adjust function:



Fig. 3-39

Coarse

The mass spectrometer performs a fast, coarse search across the specified range, starting at the mass number entered under Mass.

Fine

Fine search of the peak top within a reduced search range but with greater accuracy, starting with the mass number entered under Mass or the mass number found by Coarse.

All Channels

Adjust all channels.

Single Channels

Adjust only the selected channel.

NOTE:

The time required for completing the Adjust depends on the number of channels and the selected Dwell measurement time. If the adjust function is cancelled before it has been completed, only those mass numbers are accepted for which the adjustment was completed without an error message.

Edit submenu

Fig. 3-40

Adjust Parameter

Modify the parameters Mass, Range and ThresholdAdjust:




Fig. 3-41

File

Close Quit the dialog box

NOTE:

The parameters above are stored in the corresponding channels of the MID/MCD parameter file.

3.2.5.4.2 Peak Search function

Peak Search helps you to very quickly setup the parameters of an MID measurement containing the highest peaks. A chosen mass range is searched for the highest peaks and the mass numbers with the optimized settings are stored in the channels of the current parameter file. Using this file, you can start an MID measurement that shows all values of the found peaks. Peak Search is only available for MID and Mass Scale Calibration parameter files.

[Functions] > [Peak Search] opens a window with an empty table for the peak search data:




Fig. 3-42

Peak Search has the following main menu:

File

Quit the peak search.

Search

Start and stop the search function.

Options

Sort the channels by mass numbers or intensities

Edit

Edit the Peak Search parameters.

Under Number of Peaks, specify the number of peaks that shall be found. Then start the peak search by [Search] > [Start]. In order to obtain accurate results, it may be useful to modify some parameters under [Edit] > [Search Parameters] (recommended: 'Threshold' value of 3 ... 10%, 'Resolution' value of 1 ... 10).

File submenu




Fig. 3-43

Close

Quit the Peak Search function.

Search submenu




Fig. 3-44

Start

Start the Peak Search function.

Stop

Stop the running Peak Search function.

After the start, the peaks found are shown in a table. If the spectrum within the search range contains more than the desired number of peaks, the largest peaks are listed:




Fig. 3-45

Peak Search performs as many scans with increased gain each time, until the desired number of peaks has been found.

Options submenu




Fig. 3-46

Sort Mass

Sort the found peaks by mass number.

Sort Int.

Sort the found peaks by intensity.

Edit submenu



Fig. 3-47

Search Parameters...

Enter the parameters (search criteria) for the Peak Search function:

Search Parameters




Fig. 3-48

3.2.5.5 Editing the MCD parameter files (Quantitative Analysis)

The MCD measuring type enables you to measure concentrations. For that, an analysis matrix is used, that cooperates with the calibration factor library (if you change a calibration factor in the analysis matrix, it will be changed in the calibration factor library as well and vice versa). By the following editor you can edit such an analysis matrix.

3.2.5.5.1 Analysis matrix

Choose [Measure] > [MCD] and select a file, so the following dialog box appears. It shows the analysis matrix of the measurement, i.e. the gas components and their calibration factors for the individual mass numbers:




Fig. 3-49

In this example the calibration factor for CO₂ on mass 28 has not been taken into consideration, because it has only a minimal influence in the analysis of air.

Comp.

Name of the component.

Mass

Mass number.

- The maximum size of the matrix is limited to 64 mass numbers and 64 components.
- In an analysis each component is defined by the name and the mass number(s). Due to overlaps it is possible that a mass number has intensity portions of several components.
- If a calibration factor is missing, a “?” is displayed at its place.

3.2.5.5.2 Expansion of the analysis matrix

Additional components and mass numbers can be added to the displayed analysis matrix:

- by inserting spectra from the Spectra Library ([Edit] > [Insert from Spectra Library])
- by manual insertion (by [Edit] > [Insert Mass] or [Edit] > [Insert Component])
- by appending (and subsequent sorting)

When the last line or column has been used, an additional cell is appended automatically. Enter the component names and the mass numbers that are relevant to the analysis:

- The cursor can be positioned in any field.
- An entry is made when the cursor is moved or the <TAB> key is pressed; <ENTER> has no function.

When a component or mass number is removed, the line or the column is deleted.

When a new mass number or a new component is entered, the existing calibration factors are displayed. Undefined calibration factors are displayed as “?”. By [Options] > [Set Calib. Factors Zero] you can set all these undefined factors to zero.

NOTE:

In an analysis matrix, no undefined calibration factors are allowed. All factors must contain a positive value or zero.

If a mass number contains concentration portions of different components, a calibration factor of each component must exist for accurate analysis of this mass number. Missing calibration factors can be determined by the gas-specific calibration, see Chapter 2.8.4 Gas specific sensitivity, 142, or they can be entered manually here.

3.2.5.5.3 MCD Measure Parameter menu



Fig. 3-50

MCD measurement parameters are edited via the following main menu:

File

File handling, print analysis matrix, quit the editor.

Edit

Edit and sort the mass numbers and the component names.

Options

Set unknown calibration factors to 0.


File submenu

Fig. 3-51

Open...

Call the File manager.

Save

Store the current file.

Save As...

Store the current file under a different name or in a different directory.

Print

Print the current analysis matrix.

Close

Return to the Parameter Setup main menu.


Edit submenu

Fig. 3-52

Channel Par.

Enter the channel parameters. See Section «Channel Parameters», 183.

Concentration Par.

Enter the threshold values (concentrations) and display parameters. See Section «Concentration Parameters», 184.

Insert Mass

Extend the matrix by inserting a mass number (column) at the cursor position.

Delete Mass

Delete a mass number from the matrix.

Sort Mass

Sort the mass numbers in ascending order.

Insert Component

Extend the matrix by inserting a component (line) at the cursor position.

Delete Component

Delete a component from the matrix.

Sort Component

Sort the components alphabetically.

Insert from Spectra Library

Insert components and mass numbers from the spectra library. See Section «Insert Spectra into the matrix», 185.

Options submenu

Fig. 3-53

Set Calib. Factors Zero

Set all undefined calibration factors (marked with “?” or “*”) to zero.

Channel Parameters

Choose [Edit] > [Channel Par.] to bring up the editor:




Fig. 3-54

- Input or modification of mass numbers is inhibited for the masses ([State] > [Matrix]) that have been entered in the analysis matrix editor.
- After the channels marked with [State] > [Matrix], additional channels can be activated. These channels are additionally measured during the analysis and can be saved as well. Additionally measured intensities are ignored in the calculation of the concentrations.
- For optimizing the channel parameters use [Functions] > [Adjust].

NOTE:

In the Display group under Disp. Mode you can determine whether in addition to the calculated concentrations also the intensities (raw data) are to be displayed.

Concentration Parameters

In this editor you can enter under [Group] > [Trip] the thresholds of the concentrations, and under [Group] > [Display] the display parameters (color, etc.) for the individual concentration values.

Choose [Edit] > [Concentration Par.] to bring up the MCD concentration editor:




Fig. 3-55

File

- Print: Print the parameters of the groups Display and Trip.
- Close: Return to the MCD parameter editor.

Edit

- Copy/Exchange Supports parameter copying and swapping in the Trip and Display group:
 - If the cursor is positioned on a channel, that is, on COMPONENT, only the entire group can be copied.
 - If the cursor is positioned on one of the parameters at the bottom of the editor, you can switch between this parameter (Parameter) and all parameters (Channel) in the Copy/Exchange dialog box under Object Selection.

Display

Display parameter group.

- Color: Choose the color.
- Disp Mode: Display concentration on/off.
- Disp. F.S.R.: Full scale range.
- Disp Decades: Number of decades to be displayed.

Trip

Trip parameter group for programming the switching functions. See Chapter 3.2.5.7 Switching functions (Trip), 186.

Insert Spectra into the matrix

[Edit] > [Insert from Spectra Library] opens the Spectra Library File Manager:




Fig. 3-56

Choose the desired spectrum (have a look at it by Spectrum, if necessary) and then press the [Export] button. The following dialog box appears:




Fig. 3-57

Here you can choose the mass numbers that are to be inserted into the analysis matrix in conjunction with the molecular formula. You can mark one or several mass numbers and move them between the two groups Available and Selected by the arrows [LEFT] and [RIGHT].

[More]

Insert the current spectrum into the analysis matrix. Return to the Spectra Library to choose more spectra.

[OK]

Insert the current spectrum into the analysis matrix and return to the matrix editor.

[Cancel]

Cancel the current action and return to the matrix editor.

The spectra inserted this way do not contain valid calibration factors. The spots in the matrix, where these calibration factors belong to, are marked by '*'.

3.2.5.6 Activation of analog input channels (MID/MCD)

A measurement channel can be used as an analog input by setting its Det. Type to AI.




Fig. 3-58

If AI has been set under Det. Type, Quadstar 32-bit automatically shows the channel number (QMS 422/421: 0 ... 15, QMS 200: 0,1) of the analog input channel in the input field behind AI.

3.2.5.7 Switching functions (Trip)

In MID- and MCD measurements you can program the switching functions A (measure value $A < \text{threshold } A$?) and B (measure value $B > \text{threshold } B$?). The status (OFF/ON) can be displayed on the screen. You can combine the switching functions by logical conditions and control a digital output by the result. The parameters below are displayed only when a function type has been selected under Type:

Trip			
Type	ABS A,B	DO# A	8
Level A	0.0100ppm	Level B	0.0800ppm

Fig. 3-59

Type

Choose the type of function for the two switching functions A and B.

- OFF: Switching functions A and B not active.
- ABS A: Switching function A is turned ON, if the measured intensity is below the threshold entered under Level A. (Switching function B not active).
- ABS B: Switching function B is turned ON, if the measured intensity is above the threshold entered under Level B. (Switching function A not active).
- ABS A,B: Switching function A and B are both active (See ABS A and ABS B).

- ABS A AND B: If A AND B are both switched ON, the result is ON.
- ABS A OR B: If A OR B OR BOTH are switched ON, the result is ON.
- ABS A XOR B: If EITHER A OR B are switched ON, the result is ON.

Level A and Level B can also be used as threshold values for the Hysteresis switching function. For this switching function the entered Level B must be greater than Level A:

- HYS HI OFF: The result will be OFF when the signal value is above Level B and it changes to ON when the signal value drops below Level A.
- HYS HI ON: The result will be ON when the signal value is above Level B and it changes to OFF when the signal value drops below Level A.

Threshold values

- Level A: Threshold value of switching function A.
- Level B: Threshold value of switching function B.

DO Channel

- DO# A: Digital output channel of switching function A.
- DO# B: Digital output channel of switching function B.

3.2.5.8 Outputting the switching function state to a DO channel

The channel number is chosen under DO# A or DO# B.

QMS 422/421

The input range depends on the number of installed DO 420 boards (0 ... 31 for one board, 0...63 for two boards, and 0...95 for 3 boards; declared under [Config] > [QMS]).

QMS 200

Two digital outputs are available as standard equipment (REL_1 and REL_2).

ABS A, ABS B, ABS A,B

In these modes DO# A and DO# B correspond to the states of the switching functions A and B.

In logical combinations the following applies:

ABS A OR B, ABS A XOR B

DO# A and DO# B are inhibited.

ABS A AND B

Under DO# B the channel number of DO# A is adopted.

HYS HI OFF

Only DO# A is enabled.

HYS HI ON

Only DO# B is enabled.

3.2.5.8.1 Behavior of the switching functions

Level B > Level A




Fig. 3-60

Level A > Level B




Fig. 3-61

Level A > Level B "Hysteresis"

When HYS HI OFF or HYS HI ON is chosen, the threshold value of switching function B must be above the value of switching function A:





Fig. 3-62


Graphic symbols

The result of the threshold monitoring is displayed with graphic symbols. The colors of the following symbols are fixed and have the following meaning:

Green

 Switching function OFF

Red

 Switching function ON

For the MID Bargraph and MID Table measurements the following symbols are used:




Fig. 3-63

In MID Versus Time mode the monitor trip is signaled by an arrow. The color of the symbols can be selected under [Config] > [Screen Color].

When the two switching functions are logically combined, the symbol of switching function A is used:




Fig. 3-64

Display the switching functions

The display of the switching functions can be activated during a measurement (in the Measure program) under [Parameters] > [Setup] > [Alarm Display].

3.3 Sequence programming

The Quadstar 32-bit Process Control enables you to create and execute measurement programs (sequences). These programs can be created with the Sequence Editor and started in the Measure program. In the following description this is referred to as the Sequencer.

The Sequencer executes a previously created instruction sequence instead of directly starting a measurement. A sequence may call and start another sequence (subroutine, e.g. measurement, valve control etc.) or start additional programs.

The main functions of the Sequencer are:


- Start and control measurements (Scan Analog, Scan Bargraph, MID, MCD, etc.)
- Control valves or auxiliary equipment
- Control data storage
- Exchange data with other programs
- Perform calculations automatically
- Display of measurement data and calculated values
- Display of self defined messages

NOTE:

This section describes only the principal operating elements of the Sequence Editor. For an overview of all instructions and a complete description of them, please refer to Chapter 8. Process Control, 409.

3.3.1 Quadstar 32-bit Sequence editor

[Sequence] > [Editor] opens the File Manager in order to choose or create a sequence file. This file is opened and its content (the sequence) is displayed.



```

Sequence Editor < mid.seq >
File Edit Options
|Message( Text="Okay, let's go" )
Loop( i[0]=0;500 )
Begin
  Sequence( Par="c:\qs32bit\par\clearall.seq" )
  Sequence( Par="c:\qs32bit\par\inlges20.seq" )
  Loop( i[0]=0;3 )
  Begin
    MID( Par="c:\qs32bit\par\spr.mip", Disp=tab, SaveGfa=0 )
    SetString( gs[0] = "MID measurement @ intensity channel 0 = @";i[0] )
    Message( Text=gs[0],gfa[0][0] )
  End

  Loop( i[0]=0;3 )
  Begin
    MID( Par="c:\qs32bit\par\spr.mip", Disp=bar, SaveGfa=0 )
    SetString( gs[0] = "MID measurement @ intensity channel 0 = @";i[0] )
  End
End

```

Fig. 3-65

The Quadstar 32-bit Sequencer is designed as a command interpreter: the sequence is a series of instructions (e.g. MID(), ScanAnalog(), etc.) that are executed one after the other. An instruction consists of its name followed by a parameter list in parentheses:

- The specification of certain parameters is mandatory, whereas other parameters are optional.
- If an optional parameter is not specified explicitly, the Sequencer automatically takes the default value (for example the parameter “prot=off” does not have to be specified).

3.3.1.1 Programming

There are two ways of programming the individual sequencer lines:

- Input via the keyboard.
- Compose instructions by using the dialog editor.

Input via the keyboard

This method is intended for experienced users and requires familiarity with the command syntax. The text is entered directly at the cursor position of the selected Sequencer line. A new line is inserted with the Enter key. The [Del], [←] and [Ins] keys have the usual text editor functions.

Input via dialog editor

Set the cursor to an instruction (e.g. MID) and choose [Edit] > [Dialog] or click with the right mouse button. The dialog editor containing all parameters that are valid for the selected instruction is displayed. If the cursor is positioned e.g. on a MID Sequencer line, the following dialog box is displayed:




Fig. 3-66

You may now, guided by the dialog editor, compose a MID instruction without knowing its syntax.

3.3.1.2 Inserting a new Sequencer line (Choose Line)

Set the cursor to an empty line and choose [Edit] > [Dialog] or click with the right mouse button. The Choose Line dialog box appears that offers all possible instructions:




Fig. 3-67

3.3.1.3 Compiling/storing the sequences

When the sequence is compiled, the instructions are translated into a Quadstar 32-bit specific, executable code that is stored separately. Errors are flagged by the compiler and the cursor is parked at the corresponding position (of a line) where the first error has been detected.

Sequence menu




Fig. 3-68

Editor...

Load and edit a sequence file.

Parameter Files...

Display/edit the parameter files used in a sequence. See Chapter 3.3.1.5 Parameter Files, 198.

3.3.1.3.1 Sequence Editor



Fig. 3-69

The sequence editor is operated via the following main menu:

File

File handling, compile and print a sequence.

Edit

Edit a sequence.

Options

Display control, font selection, storage size etc.

File submenu




Fig. 3-70

Open...

Call the File Manager to open a file.

Save

Store the current sequence.

Save As...

Store the current sequence under a different name or in a different directory.

Print...

Print the sequence or a part of it. To configure the printer use [Setup] > [Printer].

Save & Compile

Compile the sequence and store it.

Compile All...

Translate all sequences in the selected path.

Close

Return to the Parameter Setup main menu.

Print

Print prints the current sequence or a part of it:




Fig. 3-71

Print Range

- All Lines: Print all Lines.
- Selected Lines: Print the selected lines (only possible, if lines are marked in the editor).
- Lines: Print lines starting at From up to To.

Print Options

- Print Line Numbers Print the numbers of the lines as well.

Compile All

Compile All translates all sequences found in the selected path.




Fig. 3-72

[Compile All]

Translates all sequences.

[Continue]

- Translates all sequences beginning at the selected sequence, in the illustrated example from the sequence airdemo.seq onward.
- Translates the next sequence if an error has occurred in a sequence.

NOTE:

Please note that a subsequence that is called by another sequence while the sequencer is running, may only be modified but not compiled in the sequence editor.

Edit submenu




Fig. 3-73

Insert Line

Insert a new line between existing lines.

Dialog...

Call the dialog editor.

If the cursor is positioned on a blank line, the Choose Line dialog box (See Section «Choose Line», 195) is opened for inserting a new line. If the cursor is positioned on a line containing a sequencer instruction, the dialog editor with the matching dialog box is displayed directly.

NOTE:

Instead of choosing [Edit] > [Dialog] you can click with the right mouse button.

To mark text (for the Copy, Cut and Delete commands) drag the cursor with pressed left mouse button from the start to the end of the desired text section, or use the cursor buttons in combination with Shift, Ctrl, Home, End, Page Up, Page Down.

Undo

“Reverse” the last action such as input/deletion of text, etc.

Cut

Delete the text marked with the cursor and save it on the clipboard. This text can be inserted in a different place with the Paste function.

Copy

Copy the text marked with the cursor to the clipboard.

Paste

Insert the clipboard content at the current cursor position.

Delete

Delete the marked text.

Goto Line...

Position the cursor on a specific line.

This function can be used for positioning the cursor on a line number flagged by the compiler as “containing errors”. The line numbers may be printed with the lines.

Find...

Search text.

Replace...

Replace text.

Choose Line

When [Edit] > [Insert Line] is selected, the Choose Line dialog is displayed:




Fig. 3-74

The Controls are described in Chapter 8.5.3 Comments, 428.

The sequencer instructions are found in the list fields Conditions, Loops and Instructions.

Options submenu



Fig. 3-75

Display...

Enter the display parameters, select font type etc.

Parameter Limits...

Define the size of the local (i[], f[]) and the global (gi[], gf[], gfa[[]], gs[]) storage arrays of the sequencer.

Display




Fig. 3-76

Tabulator Space

Indentation size. The indentation corresponds to twice the width of the letter 'B' and is valid for the printout as well (=> using [File] > [Print]).

[Choose Font]

With the Choose Font command you can select the font, the type face and size.




Fig. 3-77

3.3.1.4 Parameter limits (array limits)

3.3.1.4.1 Variable types

The Sequencer works with different variables:




Fig. 3-78

There are the following types of variables:

Integer

Integer number, range -1'000'000...+1'000'000.

Float

Floating point variables, range -1E+37 ... +1E+37, resolution 15 digits.

String

Text, stored as a character string.

The variables are organized as arrays (one dimensional tables), except the global float field gfa[][] that is a two dimensional table. Each element of the table is a variable. The size of these tables is specific to the sequence and can be defined in the Parameter Limits dialog box by specifying the maximal index.

There are local and global variables:

Local variables i[], f[]

Local variables are valid only within the sequence being currently processed. In every subsequence new local variables are created that are usable independently from the ones of the calling sequence.

i - Local integer variable array.

f - Local float variable array.

Global variables gi[], gf[], gfa[][], gs[]

Global variables are valid for all main- and subsequences. A single set of them is created, so they can be used e.g. for exchanging data between main- and subsequences.

gi -Global integer variable array.

gf - Global float variable array.

gfa - Global float variable field of the dimension $N \times 64$ for intermediate storage of measured values, calculated values, etc.

gs - Global array of character strings for defining text, error messages, file names etc. The maximum length of a character string is 100 characters.

Memory requirements

The last column shows the memory requirements. Make sure that the size of the array is matched to the application. Dimensioning large arrays may be simple and convenient, but wastes memory resources, which on a PC with little memory left can cause problems. This should be particularly avoided if nested sequences are used, because each sequence sets up its own local variables in memory.

Note that the maximum index is entered for dimensioning the array. When you enter, for example, 24 for Maximal Index, array elements 0 ... 24 (=25 variables) can be used.

NOTE:

For a detailed description of the variables please refer to Chapter 8. Process Control, 409.

3.3.1.5 Parameter Files




Fig. 3-79

Parameter Files displays the parameter files (Parameter Files) and subsequences (Sequences) used in a sequence. They can be selected and edited:

[Edit Seq.]

Edit sequences

[Edit Par.]

Edit parameter files.

3.4 Calibration Parameters

Calibration menu




Fig. 3-80

Mass Scale...

Set the mass scale calibration parameters.

Gas Specific Sensitivity...

Set the parameters for a gas-specific measurement for determining calibration factors.

3.4.1 Mass Scale

By using the mass scale calibration, deviations of the mass scale in the range of ± 0.5 amu from the nominal mass can be compensated. For that, a calibration gas with a known spectrum is required; on its mass numbers the deviations are measured.

NOTE:

For additional information on the mass scale calibration please refer to Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

Choose [Calibration] > [Mass Scale] and open a parameter file:




Fig. 3-81

Input/optimize the parameters

To optimize the parameters, the communication to the mass spectrometer must be switched on and the calibration gas must be let in. Then enter the mass numbers of the calibration gas.

- If the mass numbers are known, they may be entered manually.
- [Functions] > [Peak Search] automatically searches for the highest peaks and the corresponding mass numbers and puts them in the parameter file.

The mass numbers entered this way should be precalibrated using [Functions] > [Adjust]. If errors are reported, optimize the corresponding channels:

- 1 In the Mass Scale Calibration Parameter editor select the desired channel.
- 2 Choose [Display] > [Measure Data]. Quadstar 32-bit starts a single-channel measurement that responds immediately to any parameter changes.




Fig. 3-82

- 3 Optimize the following parameters as a function of the error messages produced by Adjust:
 - Mass: Enter the desired mass number and shift it until it has reached the peak maximum (as shown by the bargraph).
 - SEM Voltage: Set the SEM voltage (if a SEM is used) to a value that is appropriate for your SEM.
 - Amp. Range: Optimize the Amp. Range (only in Amp. Mode FIX), so that the bargraph shows approx. 10 ... 90% full scale. Please note that with increasing gain the electrometer becomes slower so that the Dwell speed must also be adjusted.
 - Threshold: Set the threshold (in % full scale) to half the peak intensity.

Hint: With the drag pointer you can observe the maximum deflection.

If errors still occur, you should check the Resolution. The easiest way is to program a Scan Analog ([Measure] > [Scan Analog]). Then choose [Display] > [Measure Data] and optimize the Resolution according to the shown spectrum.

NOTE:

The resolution should be set in such a way that the smallest value between two adjacent peaks is approx. 10% of the intensity of the smaller peak.

Afterwards transfer the resolution into the channels of the parameter file.

3.4.2 Gas-specific sensitivity (Quantitative Analysis)

In order to measure concentrations, calibration factors have to be determined. These describe the sensitivity of the mass spectrometer for a certain gas. Under Gas Specific Sensitivity you can determine such calibration factors.

3.4.2.1 Calibration matrix

When [Calibration] > [Gas Specific Sensitivity] is activated, a parameter file (calibration matrix) must be selected or created via the File Manager. This file is subsequently opened and the dialog box for entering the calibration gas specific parameters is displayed.

Each component of a calibration gas measurement is defined by its name, the mass number and the concentration. Several calibration factors of a component can be determined simultaneously:

#	Component	Conc.	Mass							
			14	28	29	32	40	44		
0	Argon	9300.0ppm	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	
1	CarbonDiox	300.00ppm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	
2	Nitrogen	78.0850%	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
3	Oxygen	20.9500%	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
4										

Fig. 3-83

Int. Std

Component selected as the internal standard.

Ext. Norm

Component selected as the reference point for external normalization.

Component

Component name.

Conc.

Concentration of the component. The concentration is entered in % or in ppm but if it is <1% it is always displayed in ppm.

Mass

Mass number.

3.4.2.2 Extending the calibration matrix

The displayed calibration matrix can be extended with additional components and mass numbers:

- by inserting spectra from the Spectra Library ([Edit] > [Insert from Spectra Library]) or
- by inserting manually ([Edit] > [Insert Mass/Insert Component]) or
- by appending with subsequent sorting (via [Edit] > [Sort Mass/Sort Component]).

When the last line or column has been filled, an additional line or column is automatically appended. Enter the various component names and mass numbers (that exist in the calibration gas) and the mass numbers for which the calibration factors are to be determined:

- Position the cursor in the desired field.
- An entry is accepted when the cursor is moved or when the <TAB> key is pressed; <ENTER> has no function.

When you remove a component or a mass number, the line or the column is deleted.

Concentrations can be entered in % (without unit of measure, e.g.: 20 or 1.5) or in ppm (with unit of measure, e.g.: 350 p, 350 ppm); 200 p, 200 ppm and 0.02 therefore is always the same concentration, shown in different representation each time.

Subsequently, define the positions on which the calibration should take place. No overlapping of gases is allowed, i.e. a mass can be assigned to only one component.

At the end you have to tell the program which component is to be used as the standard. See Chapter 3.4.3.3.1 Defining the internal standard, 205.

3.4.3 Gas Specific Calibration menu



Fig. 3-84

This function is handled via the following main menu:

File

File handling, print calibration parameters, quit the editor.

Edit

Edit and sort the mass numbers and the components.

Options

Sort the entered parameters, set the internal standards.

3.4.3.1 File submenu




Fig. 3-85

Open...

Call the File Manager.

Save

Store the current file.

Save As...

Store the current file under a different name or in a different directory.

Print

Print the current calibration parameters.

Close

Return to the Parameter Setup main menu.

3.4.3.2 Edit submenu




Fig. 3-86

Channel Par.

Enter the channel parameters. The editor can also be called by double clicking the right mouse button if the cursor is located outside the input fields and not on a scroll bar.

Insert Mass

Extend the matrix by inserting a mass number (column) at the current cursor position.

Delete Mass

Remove a mass number from the matrix.

Sort Mass

Sort the mass numbers in ascending order.

Insert Component

Insert a component (row) into the matrix at the current cursor position.

Delete Component

Remove a component from the matrix.

Sort Component

Sort the components alphabetically.

Insert from Spectra Library

Insert components and mass numbers from the spectra library. See Chapter Insert Spectra into the matrix, 185.

3.4.3.2.1 Channel Parameters

The mass numbers are adopted from the calibration matrix and the corresponding channels are marked by [State] > [Matrix]. Under [Display] > [Measure Data] you can use the Adjust-function to optimize the peak heights.

Choose [Edit] > [Channel Par.] to call the channel parameter editor:




Fig. 3-87

3.4.3.3 Options submenu



Fig. 3-88

Set Internal Standard

Define the internal standard.

Set External Normalization

Define the external normalization.

Remove External Normalization

Disable the external normalization.

3.4.3.3.1 Defining the internal standard

Quadstar 32-bit must be told which component is to be used as the internal Standard (i.e. which calibration factor is to be considered as fixed). It is displayed in the upper left corner of the dialog box. Because the relative sensitivities of the mass spectrometer are determined, pressure variations etc. have no influence on the results of an analysis, if the internal standard is used.

NOTE:

The Set Internal Standard command is enabled only if the selected mass number has been assigned to a component.

The calibration factor of the Internal Standard must already be defined at the time of the calibration. It establishes the relationship between the new and the previously determined calibration factors.

3.4.3.3.2 Setting the external normalization

In certain cases it may be useful to work with external normalization, e.g.:

- If only one component is to be measured in a gas mixture.
- If a gas mixture contains unknown components.

Through the external normalization the absolute sensitivity of the mass spectrometer is calculated based on the selected standard. The mass number selected with the cursor and the associated component become the reference point for the external normalization. The same component as for the internal standard or a different one can be chosen. The chosen component is shown in the upper left corner of the dialog box.

NOTE:

When performing the measurement, make sure the pressure conditions are the same as during the calibration, because the External normalization cannot eliminate pressure variations. Changed ion source parameters, changed SEM voltage etc. influence the measuring results as well.

The normalization factor NF is calculated as follows:

$$NF = \frac{I(g,m)}{S(g,m) \cdot C(g)}$$

In an analysis with external normalization the calculated concentration is normalized as follows:

$$C(g) = \frac{C'(g)}{NF}$$

$$C'(g) = \frac{I(g,m)}{S(g,m)}$$

NF = Normalization factor

m = Index for mass number

g = Index for gas component

I(g,m) = Ion current of the gas component g on mass number m

S(g,m) = Calibration factor of the gas component g on mass number m

C(g) = Concentration of the gas component g (%)

C'(g) = Unnormalized concentration of the gas component g (I(g,m)/S(g,m))U

NOTE:

For more information on how to calibrate the mass spectrometer please refer to Chapter 2.8.4 Gas specific sensitivity, 142.

3.5 Manual parameters

Manual menu

By Manual you can display the state of digital input channels and switch on or off digital output channels manually, provided that the communication to the mass spectrometer is switched on.



Fig. 3-89

For turning the individual digital output channels on or off, double-click the corresponding push button with the left mouse button. The switching state is indicated by the colored dot on the push button: green = on, grey = off.

3.5.1 DI/DO Manual (for QMS 422/421)




Fig. 3-90

Display left

Select the DI/DO group to be displayed on the left side of the screen.

Display right

Select the DI/DO group to be displayed on the right side of the screen.

The same group cannot be displayed at the left and the right side simultaneously.

3.5.2 DI/DO Manual (for QMS 200)

Only the two standard digital input and output channels are available:




Fig. 3-91

The digital input and output channels of the QMS 200 have the following designation in the hardware description:

DI 00

EXT_PROT (Connection TPR/PKR, Pin No. 1 <=> 2)

DI 01

EXT_RUN (Connection Control, Pin No. 10 <=> 15)

DO 00

REL_1 (Connection Control, Pin No. 11 <=> 12)

DO 01

REL_2 (Connection Control, Pin No. 13 <=> 14)

3.6 Communication parameters

3.6.1 Overview

Transmission types

Quadstar 32-bit supports several data transmission types (called 'communication types' in the following text) to the mass spectrometer (QMS):

Transmission via RS232C Normal

Transmission via ArcNet (duplex fibre optic)

only with ArcNet board in the PC.

Communication menu

In the menu Comm you can choose the communication type and set up its parameters, or switch the communication to the QMS on or off:




Fig. 3-92

QMS Comm. Type...

Select the type of communication with the QMS. See Chapter 3.6.3.1 Selecting the communication type, 214.

QMS Comm. Parameters...

Enter the communication parameters for the selected communication type. See Chapter 3.6.3.2 Setting the parameters for the selected communication type, 214.

Unit Identifications...

Set the Destination Node Identification (DID) for communication via ArcNet. See Section «Setting the node address in the QMS», 217.

Connect

Establish the connection to the QMS. See Chapter 3.6.3.5 Connect, 221.

Disconnect

Clear the connection to the QMS. See Chapter 3.6.3.6 Disconnect, 222.

SetUnit

Select a control unit activated in the Unit Identifications dialog box under Enable for measurement. See Section «Setting the node address in the QMS», 217.

3.6.2 Communication via RS232C

3.6.2.1 Protocol/transmission reliability

The communication is based on data packet (frame) transmission and is monitored for data loss and corruption. The communication protocol conforms to the "SECS-I" standard (Semi Equipment Communication Standard I).

A data transmission request transmitted from the sender to the responder must be acknowledged by the responder before the sender may start to transmit the actual message. After the message (data, parameters, etc.) has been transmitted, a checksum is forwarded and a positive or negative acknowledgment is awaited. If a positive acknowledgment is received, the transmission is completed, if a negative acknowledgment is received (error in the checksum, etc.), the transmission is repeated up to N times (N = retries for data block transmit).

NOTE:

For detailed information such as protocol, data frame formatting, etc., please consult the User's guide for the QMS interface.

3.6.2.2 Configuration of multiple serial interfaces under Windows

If multiple serial interfaces are operated, problems can occur. The following restrictions should be noted:

- The COM1 and COM3 interfaces as well as the COM2 and COM 4 interfaces cannot be used simultaneously.
- If a mouse is connected to COM1, COM3 does not function (and vice versa). The same applies analogously to COM2 and COM4. Valid combinations are COM1 and COM2 or COM3 and COM4.

3.6.2.3 RS232C normal

The normal communication between Quadstar 32-bit and mass spectrometer is based on a serial RS232C interface. RS232C is not electrically insulated from the equipment and consequently susceptible to noise. The QMS is automatically adapted to the programmed baud rate.

Connection cable

The reliable transmission distance depends on the noise field intensity in the environment and the selected baud rate. The following values should not be exceeded:

Baud rate	Max. cable length (shielded)
1200	90 m
2400	30 m
4800	30 m

Baud rate	Max. cable length (shielded)
9600	15 m
19200	7.5 m

The send and receive lines must be crossed (null modem cable). Quadstar 32-bit operates with a software handshake, i.e. the control lines for the hardware handshake (DSR, DTS, RTS, CTS) are not used.

The pin assignment of the DB9 connector on the QC is described in the User's guide of the mass spectrometers.

Pin	Significance
1	Carrier Detect (CD)
2	Receive Data (RXD)
3	Transmit Data (TXD)
4	Data Terminal Ready (DTR)
5	Signal Ground (GND)
6	Data Set Ready (DSR)
7	Request To Send (RTS)
8	Clear To Send (CTS)
9	Ring Indicator (RI)

Tab. 3-1

DB9 connector (male, on QC)

Make sure that the send and receive lines are crossed only once in the cable. The supplied cable conforms to this requirement. If problems occur, install a Reverser.

NOTE:

To protect the line drivers, the computer and the QMS should be connected to the same phase of the power line.

Only shielded cables should be used.

To prevent ground loops, the shield of the cable should only be grounded at the QMS.

3.6.2.4 Communication via ArcNet

For high data throughput as well as multiplex operation, Quadstar 32-bit supports the ArcNet interface of the mass spectrometer via Duplex Fibre Optic lines.

NOTE:

For communication via the ArcNet interface, a special board OPA 200 or PCMCIA (PC-Card) with a matched optical interface is required in the PC. If several mass spectrometers are to be connected to the same PC, an optical hub OHA 200 is required. For additional information please contact our service office.

3.6.2.4.1 Installation of the ArcNet board

For installing and configuring the ArcNet board please consult the instructions of the board manufacturer. The ArcNet parameters are preset to the following values:

- Source identification (SID): [hex] A0

For choosing the different parameters call the ArcNet-Parameter dialog box. See Chapter 3.6.3.4.2 Communication via ArcNet, 216.

NOTE:

For the OPA 200 the values in Quadstar 32-bit must match the ones that are really set on the ArcNet board. To change values on the ArcNet board, perhaps jumpers must be changed. Please refer to the user's guide of the ArcNet board.

Verify that the memory area D000-D07F (HEX) is not being used by another BUS board. If it is although, assign the ArcNet board to another area.

If you have done a hardware search with built in ArcNet board by accident, delete the ArcNet board completely from the operating system. Quadstar 32-bit uses the shared Arcnet board driver included on the installation cd.

3.6.2.4.2 Multiplex operation

With the ArcNet interface, Quadstar 32-bit is able to communicate with several mass spectrometers via an optical hub which is available as an option. If you choose [Comm] > [Connect], all units that were set to Enable in the Unit Identifications dialog box are addressed. See Section «Setting the node address in the QMS», 217.

3.6.2.4.3 Installation of the hub

The hub is an optical star distributor and can be used for setting up a network comprising several control units. For larger networks, several hubs can be used. For the installation procedure, please refer to the user's guide of the hub.

3.6.2.4.4 Data collection possibilities

The type of data collection depends on the configuration of the Quadstar 32-bit:

- If the sequencer is used, multiplex operation with several control units QMS 422/421 and QMS 200 is possible in conjunction with the SetUnit() instruction.
- In operation without a sequencer the desired control unit is selected manually through a corresponding dialog box:




Fig. 3-93

- The Set Unit dialog box is normally called via the Comm menu.
- This dialog box is automatically displayed if several control units are initialized and equipment-specific parameters or calibration data are accessed. See Section «Setting the node address in the QMS», 217.

3.6.2.4.5 Unit-specific files

The unit specific files (e.g. required for calibration) are stored in different subdirectories associated with the corresponding control unit. The directory is identified via the logical number of the control unit.

Example:

For a description of this dialog box see Section «Setting the node address in the QMS», 217.




Fig. 3-94

For a standard installation of Quadstar 32-bit on drive C and input of the Unit Identifications as shown on the left, the directories A001 ... A005 are created by Quadstar 32-bit:

- C:\q_s32bit\AUXI Auxiliary directory for unit 000
- C:\q_s32bit\A001 Auxiliary directory for unit 001
- C:\q_s32bit\A002 Auxiliary directory for unit 002
- C:\q_s32bit\A003 Auxiliary directory for unit 003

C:\qs32bit\A004 Auxiliary directory for unit 004

C:\qs32bit\A005 Auxiliary directory for unit 005

NOTE:

The number of the directories A001 ... A005 corresponds to the Unit number and not the Node-ID.

Equipment-specific files (for parameters and data) are assigned to each control unit. Always use the dialog box Unit Identifications for changes, do not modify those files or directories manually.

3.6.3 Setting the interface parameters

3.6.3.1 Selecting the communication type

- 1 From the Parameter Setup program choose [Comm] > [QMS Comm. Type].
- 2 Select the desired type (RS232C, ArcNet, ArcNet Server) via the displayed option command buttons.

3.6.3.2 Setting the parameters for the selected communication type

- 1 From the Parameter Setup program choose [Comm] > [QMS Comm. Parameters].
- 2 Enter the desired transmission parameters.

3.6.3.3 QMS communication type

In the QMS Communication Type dialog box you can select the desired communication type via the option command buttons:




Fig. 3-95

RS232C normal

Communication via a serial RS232C interface.

ArcNet

Communication of one Quadstar 32-bit application with one or several control units via ArcNet.

ArcNet Server

Communication of one or several Quadstar 32-bit applications with several control units via ArcNet. Two independent Measure programs can be run on two control units e.g., but it's not allowed to let two programs communicate with the same control unit.

3.6.3.4 QMS communication parameters

For each of the possible communication types a separate dialog box with separate storage of the transmission parameters is available.

For each communication type the standard setting can be retrieved via the [Default] command button. The entered parameters are saved only if you quit the dialog box by choosing the [OK] button.

3.6.3.4.1 Communication via RS232C Normal




Fig. 3-96

The baud rate to be selected depends on the transmission distance and the level of the noise fields in the environment. The QMS is automatically adapted to the selected baud rate.

Connection related parameters

Baud Rate

Speed in bit/sec. at which the data are transmitted between the PC and the mass spectrometer.

Connection Line

Select the serial PC interface.

Protocol related parameters

The communication protocol is derived from the block transfer protocol specified in SEMI E4-87 (SEMI = Semiconductor Equipment and Material International). For the communication of Quadstar 32-bit with the mass spectrometer the protocol timeouts T1, T2 and T3 have fixed settings:

Transmit protocol timeout T2

Maximum time (in ms) the sender waits for:

- A response to a data transmission request.

- A positive or negative acknowledgement from the responder after the checksum has been transmitted.
- Maximum time (in ms) the responder waits for incoming data after an <EOT> has been sent (in response to <ENQ>).

Reply timeout T3

Maximum waiting time (in ms) after a data transmission request by the mass spectrometer.

Inter-character timeout T1

Maximum waiting time (in ms) for the next byte during the reception of the individual data bits within a frame.

The number of retries is recorded in two separate counters, i.e. a distinction is made between the attempt to establish a connection and the frame transmission retries after receipt of a <NAK>.

Retries for handshake

Number of attempts to establish a connection to the QMS after a transmit protocol timeout.

Retries for data block transmit

Number of times Quadstar 32-bit has tried to transmit a frame after receipt of a negative acknowledgment.

Control Elements

[Default]

Restore the communication parameters to their standard factory settings.

3.6.3.4.2

Communication via ArcNet




Fig. 3-97

Protocol related parameters

Transmit protocol timeout T2

Maximum waiting time (in ms) for the confirmation of a correct data transmission to the QMS.

Reply timeout T3

Maximum waiting time (in ms) after a data transmission request from the QMS.

Retries for data block transfer

Number of transmit retries after a Transmit Timeout before a Transmit Timeout error is signaled.

Retries for data block receive

Number of receive retries after a Receive Timeout before a Receive Timeout error is signalled.

Retries after failed transmission

Number of times Quadstar 32-bit attempts to retransmit a frame after a negative acknowledgment has been received.

ArcNet related parameters**Source identification (SID)**

Station number 01 ... FF of the sender (PC).

The default values accessible via the [Default] command button correspond to the standard configuration normally used.

Control elements**[Default]**

Restore the communication parameters to their standard factory settings.

The ArcNet related parameters are restored to the following default values:

Source identification (SID)

A0h

Retries for data block transfer

3

Retries for data block receive

3

Retries after failed transmission

6

NOTE:

The [Default] command button may change settings, that have to be modified on the mass spectrometer (node address DID).

Setting the node address in the QMS

The node address 01...FF (Destination identification DID) of each QMS can be entered in the following dialog box. This address must be programmed as a unique identification number on every QMS:




Fig. 3-98

Control elements

Identifications

- Unit: Internal logical number of the control unit
- Name: Logical name of the control unit.
- Node-ID: Node address (Node) of the control unit on the ArcNet network.
- Enable: Activate the unit in the net.

[Set Node Address at 'Q_Unit1']

This button is used to program a new node address (HEX 01 ... FF) in the selected mass spectrometer, in this case in the unit Q_Unit_5.

NOTE:

This option is only available if the cursor is located in an input field (Name, Node-ID or Enable) and if the communication is switched off.

Setting the node addresses of several control units

To put the units into operation, a node address must be programmed to each control unit in the network by the button Set Node Address at "Unit Name".

NOTE:

Only the mass spectrometer to be programmed right now may be connected to the ArcNet board in the PC. Remove the connectors of all other units of the optical hub or switch off the power supply of all other units or connect the control unit with a direct line to the PC.

Selecting a control unit in the ArcNet network (Set Unit)

In Set Unit you can choose (with a double click of the left mouse button) one of the QMS units activated in the Unit Identifications dialog box under Enable:




Fig. 3-99

Parallel multiplex operation of several mass spectrometers

The operation mode ArcNet (Communication Type: ArcNet) provides the fastest measuring and data acquisition, but allows only peer to peer communication (one Quadstar 32-bit application to one control unit). It's not possible to operate another Quadstar 32-bit application at the same time on ArcNet. In contrast, operation mode ArcNet Server (Communication Type: ArcNet Server) allows to operate several Quadstar 32-bit applications at the same time on ArcNet, but the transmission gets the slower the more applications are run simultaneously.




Fig. 3-100

For example, you can tune the ion source of unit 1 (TuneUp) while you run a Scan Analog measurement on unit 2 (Measure). Every application communicates independently from the other with the units assigned to it. Although, it is never allowed to access a unit simultaneously by different applications. You must not run a Scan Analog measurement (Measure) and a Leak Test (Service) on unit 2 at the same time for example.

Example

You intend to use two mass spectrometers named 'MS1' and 'MS2'. On MS1, a Scan Analog measurement shall be performed while MS2 is doing a MID measurement. In order to do so, proceed as follows:

- 1 Select all units you want to use (in Parset program under [Comm] > [Unit Identifications]), i.e. here MS1 and MS2. Confirm the choice by pressing the [OK] button and leave Parset.




Fig. 3-101

- 2 Start a Measure program and select unit 'MS1' under [Comm] > [Set Unit]. Connect the unit by [Comm] > [Connect] and start the Scan Analog measurement:




Fig. 3-102

- 3 Start a second Measure program and select unit 'MS2' under [Comm] > [Set Unit]. Connect the unit by [Comm] > [Connect] and start the MID measurement:




Fig. 3-103

The second (new) Measure window conceals the first one, so place the mouse cursor on the title bar and drag it aside while holding the left mousebutton pressed. Now you see the two simultaneously running measurements on MS1 and MS2. For better distinction, the title bars of the measure windows are consecutively numbered ([1], [2], ..). In the lower right corner of the status bar you see the current unit name (e.g. 'MS1').

Naturally the parallel multiplexing works with the sequencer (Process Control) as well, but the used units are selected by the instruction SetUnit() then. The windows can be positioned by the instruction DispWindow(). To avoid interference between different sequences, make sure the groups of units used by the sequences are strictly separated.

Example

You want to run two sequences. The first works with the units 'SPM1' and 'SPM2' and the second uses 'MS1' and 'MS2'. So you mark all units you intend to use (under 'Unit Identifications'):




Fig. 3-104

Sequence 1	Sequence 2	
... SetUnit (Unit="SPM 1") ... SetUnit (Unit="SPM 2") SetUnit (Unit="MS 1") ... SetUnit (Unit="MS 2") ...	Right The groups of units for sequence 1 and sequence 2 are completely separated.
... SetUnit (Unit="SPM 1") ... SetUnit (Unit="SPM 2") SetUnit (Unit="MS 1") ... SetUnit (Unit="SPM 2") ... SetUnit (Unit="MS 2") ...	Wrong ! The unit 'SPM2' appears in both sequences and will presumably cause a malfunction.

3.6.3.5 Connect

The communication to the QMS is switched on by Connect.

Connect establishes the connection based on the communication type RS232C Normal, ArcNet or ArcNet Server selected under QMS Communication Type.

Automatic adaptation of the baud rate on the QMS for RS232C Normal

In RS232C normal mode the Quadstar 32-bit communication module tests the baud rate setting and changes it to the value set with the QMS Comm. Parameter.




Fig. 3-105

- During the baud rate search the transmission speed is displayed.
- To terminate the process choose [Cancel].

When the communication link is ready for operation, the LED simulated in the lower left corner of the active window changes to green:




Fig. 3-106

- 1 LED for communication state
- 2 Message field

3.6.3.6 Disconnect

Choose Disconnect to switch off an activated communication link.

3.7 Setup parameters

Setup menu




Fig. 3-107

SEM/Emission Control...

Display the dialog box of the SEM/Emission Control for:

- Switching the SEM voltage on/off.
- Entering the Common SEM voltage.
- Switching the emission on/off.

General...

Display the General Setup dialog box for general settings. See Chapter 3.7.2 General setup parameters, 226.

Printer...

Configure the printer. See Chapter 3.7.3.1 Configuring the printer, 229.

Page Setup...

Setup the printer page (margins). See Chapter 3.7.3.2 Set up the printer page (Page Setup), 230.

Printer Fonts...

Select printer font. See Chapter 3.7.3.3 Select the printer font, 230.

Print Testoutput

Print a specimen. See Chapter 3.7.3.4 Print test output, 231.

Sequence Printing...

Change over the printing mode for sequencer operation. See Chapter 3.7.3.5 Log output, 231.

ASCII Conversion...

Parameter selection for binary to ASCII conversion of measurement data. See Chapter 3.7.4 ASCII Conversion, 232.

3.7.1 SEM/Emission Control

SEM/Emission Control supports:

- Turning the Emission on/off
- Turning the SEM high voltage on/off and entering the Common SEM Voltage, if a SEM is used.

NOTE:

This dialog box can also be called by entering <Ctrl+S>.

SEM and Emission can only be switched ON/OFF if the communication to the QMS is established.

The type of displayed dialog box depends on the selected operating mode (normal or with a QMU 112):

3.7.1.1 Normal Operation (without QMU 112)




Fig. 3-108

Control elements

SEM

The SEM status is displayed symbolically in the status bar: when the SEM voltage is ON, the symbol changes to orange:



SEM high voltage OFF (grey).



SEM high voltage ON (orange).

Common SEM Voltage

Enter the SEM voltage in common for all channels.

NOTE:

The common voltage is set in the channel parameter editor in Detector group under the SEM Voltage parameter by clicking on either arrow with the right mouse button. The voltage is now marked by << xxxx >> on the display.

Emission

The status of the ion source emission is displayed symbolically in the status bar: when the emission is on, the symbol changes to Orange:



Emission switched off (grey).



Emission switched on (orange).



Emission switched on, but FAULT (red).

3.7.1.2 Operation with QMU 112

If [Option] > [QMU 112] has been selected under [Config] > [QMS], the dialog box changes.

The connected QME 125 units must be tagged under QME 125 # and the emission for each QME 125 must be switched on individually.






Fig. 3-109

SEM

The SEM status is shown jointly in the status bar for all activated QME 125 units; when the SEM voltage is on, the symbol changes to Orange:

-  SEM high voltage switched off (grey).
-  SEM high voltage switched on (orange).

Common SEM

Input of the common SEM for all channels of each QME 125.

NOTE:

The entered voltage is activated and displayed symbolically as << xxxx >> by the channel parameter editor in the Detector group under SEM Voltage when you click the right mouse button.

Emission

The emission status is shown in the dialog box with ON/OFF and also displayed symbolically in the status bar as a LED:



Fig. 3-110

The symbol on the far left corresponds to QME 125 #0, the one on the far right to QME 125 #7. The symbol colors have the following meaning:

- QME 125 not selected (light grey).
- QME 125 selected, emission switched off (dark green).
- QME 125 selected, emission switched on (light green).
- QME 125 selected, emission switched on, but Fault (red).

3.7.2 General setup parameters

Various general settings:




Fig. 3-111

Control elements

Unit of Total Pressure

Select unit of measure as mbar, Pa(scal), Torr for total pressure measurement

NOTE:

For total pressure measurement, additional hardware is required:

QMS 422/421: The total pressure is measured by using extension boards; one of each type can be installed in the QMS:

- PI 420: Evaluation electronics for two Pirani gauges.
- PE 420: Evaluation electronics for one cold cathode gauge.

QMS 200: The total pressure is measured by a sensor (with built-in evaluation electronics) that can be connected to the QMS 200:

- TPR 250: Compact Pirani Gauge
- PKR 250/251: Compact Full Range Gauge

Disk Memory Space Warning Thres. [MB]

Warning threshold for low remaining disk space.

A warning is output to the operator if after the data save less than the specified amount of storage capacity remains on the selected disk drive.

Enter the warning threshold numerically (0.01 ... 655.35 MB) or by clicking on the [UP] or [DOWN] arrows on the right side of the list box. To turn the threshold monitoring OFF, click on either arrow with the right mouse button.

Execution File Path

Directory of the Quadstar 32-bit program files (.EXE, .DLL) (has strictly informative character and cannot be changed).

Auxiliary File Path

Directory of the auxiliary parameter files. Can be entered directly or via the Select Path command button.

Parameter File Path

Directory of the QMS measurement parameter files. Can be entered directly or via the Select Path command button.

Data File Path

Directory of the measurement data files. Can be entered directly or via the Select Path command button.

User Name

Enter the operator name via the ACCESSC program (password protection). See Chapter 9. Access Rights Control, 531.

Alarm Sound Definition

Enter the duration and frequency of the alarm tone.

- Frequency: Tone frequency - input range 1 ... 20,000 Hz.
- Duration [ms]: Warning tone duration - input range: 10 ... 60,000 ms.
- [Test Sound]: Test the tone.

Simulation Mode

Switch the Simulation Mode in the QMS on/off. See Chapter 3.7.2.1 Simulation mode, 228.

- OFF: Simulation switched off.
- Internal: Simulation via QC internal measurement path.
- External: Simulation via QC external connection.

Error History

Turn the saving of the error messages on/off.

This function is used for recording errors that have occurred during a measurement. If the programmed maximum number of recorded error messages is exceeded, the oldest messages are overwritten.

Max. number of Errors in Error History

Maximum number of error messages in the error history.

Digits for Float Values

Number of decimal positions for float values.

This parameter can also be modified in the setup dialog boxes of MID Table and MCD Table.

[Select Path]

Select the drive and directory under which the corresponding files are to be stored or searched.

The [Select Path] command button is only accessible if the cursor is located in the General Setup dialog box on one of the three input boxes such as Auxiliary File Path, Parameter File Path or Data File Path.

3.7.2.1 Simulation mode

Simulation Mode provides a simulated measurement mode without vacuum, i.e. with the ion source switched off. For this purpose the air spectrum covering the masses 1 ... 63 is stored in the Quadrupole Controller (QC). To allow "testing" of higher mass numbers, the spectrum is repeated after mass 64. The parameters of the Amplifier group are not relevant in the simulation.

The simulation works in Scan Analog, Scan Bargraph, MID and MCD mode. For operation with Quadstar 32-bit, two simulation modes are available:




Fig. 3-112

Internal simulation

When [Simulation Mode] > [Intern] is selected, the D/A converter signal is transmitted via an analog MUX and the filter to the amplifier. This operating mode requires no external connection on the rear panel of the QC. However, it can only be used for Detector Type FARADAY, SEM or CH-TRON.

External simulation

With [Simulation Mode] > [Extern], noise can be added to the signal for test purposes.

QMS 422/421:

The output signal of the D/A converter can be tapped on the rear panel of the QC 421/422 at the QME/QMH connector, and be returned to the corresponding input, depending on the type of collector.

- If Detector Type FARADAY/SEM is selected, the electrometer input of the QME/QMH-connector must be used for the signal return.
- If Detector Type ION COUNT is selected, the signal input of the IC 421 (Ion counter preamplifier) installed on the QC must be used for the signal return.

QMS 200:

The output signal of the D/A converter can be tapped at the Control connector and be returned via the same connector.

Status bar information

The simulation state is displayed in the status bar: SIM is displayed when simulation is switched on:

 Simulation switched on.

SIM

Simulation switched on.

NOTE:

Always switch the simulation mode to OFF when not intentionally in use! It might be very irritating to get simulated measure data instead of real data.

Behavior of the digital outputs

Please note that for safety reasons the digital outputs are disabled in simulation mode.

3.7.3 Printer

Quadstar 32-bit supports the output of measurement data, error messages and parameters in text form as well as the Print Screen function.

3.7.3.1 Configuring the printer

Only the installed printers are displayed. The printer identified with the color bar is the currently active standard printer:




Fig. 3-113

In this dialog box you can define the (standard) printer and its configuration.

The printer ports must be configured through the Windows Control Panel. Information on how to install a printer can be found in the Windows Manual or Online Help

Control elements**Install**

Under Install you can call a dialog box in which you can define the standard settings of the selected printer.

3.7.3.2 Set up the printer page (Page Setup)

In this dialog box you can enter the margins of the printer page in cm or inches (left, right, top, bottom). The values shown below correspond to the default values:




Fig. 3-114

Control elements

Margins

Page margins in cm or inches.

- Left: Left margin.
- Right: Right margin.
- Top: Top margin.
- Bottom: Bottom margin.

[Default]

Restores the standard settings.

3.7.3.3 Select the printer font

In this dialog box you can select the type of printer font, the style and the points for outputting Quadstar 32-bit documents. The selected font is displayed in the Sample field.




Fig. 3-115

NOTE:

To get quick and adequate printouts, we recommend to use a Printer font instead of a True Type font; moreover, it should not be larger than 12 points.

3.7.3.4 Print test output

Choose this option to create a test output using the selected font. This option can also be used for testing the correct functioning of the connected printer.

Print Testoutput produces the following:

```
PRINTER TEST Selected Configuration with NORMAL
TYPEFACE: Arial
SIZE 10
PRINTER: HP LaserJet4/4M
PORT: LPT1:
```

```
PRINTER TEST: Selected Configuration with BOLD
TYPEFACE: Arial
SIZE 10
PRINTER: HP LaserJet4/4M
PORT: LPT1:
```

3.7.3.5 Log output

Quadstar 32-bit supports the output of a buffered line-by-line log (Line Printing) instead of the full page output that's normally performed by Windows (Page Printing).

In the following dialog box you can choose the printing mode (page-by-page/line-by-line) and the printer if necessary. Choose [Setup] > [Sequence Printing]:




Fig. 3-116

Printing Mode

- Page Printing: The sequencer log is output on the standard printer. A page is only printed when it is full.
- Line Printing: The sequencer log is output line-by-line on the selected printer. The printer specified here does not have to be the same as the standard printer. You can use several printers this way (e.g. an inkjet and a matrix printer).

NOTE:

Laser or inkjet printers do not output a page unless it is full. An incomplete page can normally be output by manual form feed.

For more information about the log output in sequencer operation see Chapter 8. Process Control, 409.

Line Printing Parameter

- Printer Port: Port, to which the line printer to be used is connected (LPT1 ... LPT3).

3.7.4 ASCII Conversion

The Display Saved Values program of Quadstar 32-bit is able to convert binary measurement data to ASCII data so that they can be exported to other applications such as EXCEL, etc.

Control characters, time format and options for this export can be specified in the following dialog box:




Fig. 3-117

Column Separator

Separator for dividing columns

- Char 1: Character 1 (Default: Horizontal tabulator 'HT')
- Char 2: Character 2 (Default: OFF)

Line Separator

Separator for dividing lines

- Char 1: Character 1 (Default: Carriage return 'CR')
- Char 2: Character 2 (Default: OFF)

Fraction of Seconds

Outputting fraction of seconds

- Char: Separating character (Default: Colon ':')
- Digits: Number of digits (Default: 2)

Conversion Options

Select auxiliary information

- Write Header: Create complete file header containing auxiliary information and column labeling (Default selected)
- Write Column Title: Label columns only

NOTE:

The option **Write Column Title** is only available if **Write Header** is not checked already. If no option is selected, an ASCII-file is created that contains only the measure data as columns of numbers. This is useful if the data is automatically processed, because the data is always located at the same file position.

[Default]

Set the above values to their defaults (e.g. for Excel).

NOTE:

Make sure that in the Windows Control Panel/International/Number Format, the 1000 Separator is set to the quotation mark (') in order to avoid problems with the exponents under Microsoft® EXCEL®.

3.8 Instrument Configuration

Configuration menu




Fig. 3-118

QMS...

Define the hardware configuration of the QMS. See Chapter 3.8.1 QMS Configuration, [234](#).

DI/DO Identification...

Define the names of analog and digital I/O channels. See Chapter 3.8.2 Digital inputs and outputs, [237](#).

AI Characteristic Curve...

Specify the sensors, enter the characteristic curve of the analog input channels. See Chapter 3.8.3 Analog input characteristic curve configuration, [238](#).

AO Characteristic Curve...

Enter the characteristic curve of the analog output channels. See Chapter 3.8.4 AO Characteristic Curve Configuration, [240](#).

Screen Color...

Select the Quadstar 32-bit specific desktop colors. See Chapter 3.8.5 Screen Colors, [245](#).

3.8.1 QMS Configuration

The different parameter editors of Quadstar 32-bit read the QMS hardware configuration table to determine which options are installed. Only those options appear in the list boxes of the channel parameter editors which are feasible according to the [Config] > [QMS...] menu.




Fig. 3-119




Fig. 3-120

The displayed picture depends upon whether or not the interface to the QMS is active.

- If the interface is switched on, the configuration (the hardware installed in the QMS) is read back and displayed under Hardware Boards.
- If the interface is switched off, the existing components can be switched on or off manually via the control fields displayed under Hardware Boards.

The following information items are read in by the QMS if the interface is switched on but they cannot be modified:

- QMS EPROM Version: Software version of the main processor on the QC 422/421/200.
- DSP EPROM Version: Software version of the digital signal processor on the QC 422/421/200.

In the following list boxes the installed hardware can be specified (except if only one entry is possible, then there is no list box but a text instead, which can't be modified):

QMS

Define the mass spectrometer type. (QMS422, QMS421, QMS200)

QMA

Define the analyzer type. (QMA 125, QMA 4xx, QMS 200).

MASS-R

Mass range of the analyzer.

Detector

Declare the detector type: FARADAY, SEM, CD-SEM, CH-TRON, depending on the installed analyzer (QMA 125, QMA 4xx or QMS 200).

IS-Type

Define the ion source type. The possible selections depend on the QMS.

Option

Declare the hardware extensions installed in the QMS.

The possible choice (QMU 112, CLA 421, CP 400) depends on the installed analyzer.

Protection (Sem+Fil)

Control the filament and the SEM supply (INTERN, EXTERN or EXT-PROT).

- INTERN: SEM and Emission are switched on and off manually (=>Ctrl+S).
- EXTERN: SEM and Emission are controlled by an external system (e.g. a pressure gauge or a switch).
- EXT-PROT: SEM and Emission are protected by an external system (e.g. a pressure gauge or a switch), but otherwise switched on and off manually.

For further information on this topic see Chapter 5.2.1.2.5 Filament Protection Extern (QMA 200/QMA 4XX), [§ 323](#).

HV External Supply

Operate the QMS with an external SEM supply.

With this switch it is possible to operate the QMS in SEM mode even if the HV 420 or HV 421 boards are not installed.

Emission External

Operate the QMS with external Emission control.

With this switch on, the emission is neither controlled nor monitored by Quadstar 32-bit. Control must be supplied by a separate system.

3.8.1.1**First time operation, configuration changes**

In the parameter editors only those parameters appear which are valid for the installed hardware. For this reason the hardware configuration must be updated in the QMS XXX Configuration dialog boxes when Quadstar 32-bit is put into service for the first time or whenever system changes have been made.

Example:

In the Type list box under [Group] > [Detector] of the channel parameter editor only those ion collectors are offered for selection, which are actually feasible according to [Config] > [QMS] > [Detector].

NOTE:

After a change in the QMS XXX configuration the old measurement parameter files are not updated automatically.

3.8.2 Digital inputs and outputs

3.8.2.1 QMS 422/421

In the QMS 422/421 the number of available channels depends on the hardware configuration, that is, the number of installed supplementary boards (Options: DI#0, DI#1, DO#0, DO#1, DO#2):

- DO 420: Board 0 ... Board 2 (DO 00 ... DO 31, DO 32 ... DO 63, DO 64 ... DO 95)
- DI 420: Board 0 ... Board 1 (DI 00 ... DI 31, DI 32 ... DI 63)

Up to three DO 420 boards and two DI 420 boards can be installed in the QMS 422/421 in any slot sequence. Quadstar 32-bit enables the corresponding DO channels.

3.8.2.2 QMS 200

The QMS 200 features two digital inputs and two digital outputs as standard equipment, but no additional options can be installed:

- two digital inputs: DI 00 ('ext. Protection') and DI 01 ('ext. Run')
- two digital outputs: DO 00 (REL_1) and DO 01 (REL_2)

3.8.2.3 DI/DO identification

With the DI/DO Identification editor, logical names (max. 8 characters) can be assigned to the individual lines of the digital inputs/outputs:




Fig. 3-121

Identifications

Logical name of the DI/DO channels

I/O Boards

Declaration of the DI/DO board

- DI#0, DI#1: board 0 ... board
- DO#0 ... DO#2: board 0 ... board 2




Fig. 3-122

DI 00

Logical name DI 0 ('Ext. Protection')

DI 01

Logical name DI 1 ('Ext. Run')

DO 00

Logical name DO 0 ('REL_1')

DO 01

Logical name DO 1 ('REL_2')

Procedure

- 1 Choose [Config] > [DI/DO Identifications].
- 2 Select the desired hardware option by activating the command buttons of the I/O Boards group.
- 3 In the Identifications list box select the desired channel by means of the vertical scroll bar, position the cursor, and enter the names.

NOTE:

If several identifications are entered, move the cursor by TAB key or mouse.

3.8.3 Analog input characteristic curve configuration

3.8.3.1 Overview

In MID/MCD mode Quadstar 32-bit supports the data acquisition via analog inputs. The number of channels depends on the installed hardware:

QMS 422/421

The AI 421 board, which can be installed in the QMS as an Option, offers 16 analog input channels.

Resolution 5 mV/lsb

Measurement range -10.24 ... +10.24 Volt

QMS 200

The QMS 200 features two analog inputs. This configuration cannot be extended.

Resolution	156.25 $\mu\text{V}/\text{lsb}$
Measurement range	-5.12 V ... +5.12 V

The measurement data acquisition via the analog inputs is based on the following two possibilities:

- An Analog Input channel is assigned via Det.Type = AI to a QMS channel. During each measurement cycle the analog input channels are measured and the measurement data (analogously to the intensities) made available for further processing (saving, displaying, etc.).
- Only in sequencer operation mode: The analog inputs are read in by the QMS as channel-independent parameters at any time. This method is used if a sensor of a secondary function is to be measured or monitored.

3.8.3.2 Specifying an AI sensor

For operation with Quadstar 32-bit, the AI characteristic curve of the connected sensor must be assigned to the used AI channel; Quadstar 32-bit works with this curve. The sensor type (AI Type) and unit of measure (Unit) can be specified as well as the sensor characteristic.

If no Characteristic Curve is specified, the factory supplied, linear Characteristic Curve (0 ... 10'000 mV \Leftrightarrow 0 ... 10'000 Units for QMS 422/421, or 0 ... 5'119 mV \Leftrightarrow 0 ... 5'119 Units for QMS 200) is used.

Programming an AI characteristic

When programming the characteristic curve of a sensor, the curve must be entered by means of basic points or formula. A linear characteristic can be fixed with two set points; for programming a non-linear characteristic, up to 32 basic points are possible. Intermediate values are calculated through linear interpolation.

NOTE:

The characteristic curve should only be defined across the expected range of values.

Input ranges

AI 421

Input voltage: -10240 ... 10238 mV

Value range: -1.000E+37 ... 1.000E+38

QMS 200

Input voltage: -5120 ... 5119 mV

Value range: -1.000E+37 ... 1.000E+38

3.8.4 AO Characteristic Curve Configuration

3.8.4.1 Overview

In Process Control (sequencer, SetAO() command) Quadstar 32-bit can output any analog values, that are assigned to voltages by the characteristic curve. The number of available channels depends on the installed hardware:

QMS 422/421

The AO 421 board, which can be installed in the QC 421 as an Option, offers 12 analog output channels.

Resolution	5 mV/lb
Output range	-10.24 ... +10.24 Volt

QMS 200

The QMS 200 features four analog output as standard equipment. This configuration cannot be extended.

Resolution	156.25 μ V/lb
Output range	-5.12 V ... +5.12 V

3.8.4.2 Specifying an analog output

For operation with the sequencer, an output characteristic (AO characteristic) must be assigned to the used AO channel. The sensor type (AO Type) and unit of output (Unit) can be specified as well.

If no Characteristic Curve is specified, the factory supplied, linear Characteristic Curve (0 ... 10'000mV \Leftrightarrow 0 ... 10'000 Units for QMS 422/421, or 0 ... 5'119mV \Leftrightarrow 0 ... 5'119 Units for QMS 200) is used.

When operated manually, the AO characteristic curve is not used. The value to be output is then scaled in a linear or logarithmic mode to the output voltage range of the AO, according to the setting made under AO-Mode (Register card Output of the channel parameter editor).

3.8.4.3 Programming the characteristic

The programming of the analog outputs (AO's) is similar to the programming of the analog inputs (AI's). The curve can be defined by means of basic points or formula. A linear characteristic can be fixed with two set points; for programming a non-linear characteristic, up to 32 basic points are possible. Intermediate values are calculated through linear interpolation.

NOTE:

Please note that the programmed output voltage (entered under mV) shall not exceed -10240 ... +10240 mV in the case of the AO 421 board, or -5120 ... +5119 mV for the QMS 200.

Example**Specifying an analog input (AI)****Procedure**

- 1 In the Parameter Setup program choose [Config] > [AI Characteristic Curve].
- 2 In the AI Channel text field, select the desired channel by entering its number directly, or use the [UP]/[DOWN] buttons.
- 3 In the AI-Type text field enter the type (pressure, temperature etc.) and in the Unit field enter the unit of measure (mbar, °C etc.).
- 4 The type of sensor characteristic determines how to continue:
 - The sensor has a linear characteristic:
 - Enter start- and endpoint in the Basic points field. Press the TAB key after each entry or move the cursor to the next input field with the mouse.
 - The characteristic of the sensor is non-linear but given by an algebraic expression:
 - Call the Formula Editor by pressing the Formula button. See Section «[Formula]», 244.
 - Enter the formula for your measurement as a function of the X-value in mV (variable X).
 - Enter the minimal and maximal voltage for X in mV and press [OK]; the basic points are now calculated in such a way that the largest relative deviation of the values between the points is minimized.
 - The characteristic of the sensor is non-linear and not given by an algebraic expression that contains only allowed elements as described in Section «[Formula]», 244:
 - Under Basic points enter the desired basic points. Press the TAB key after each entry or move the cursor to the next input field with the mouse.

Expanding an existing table of basic points

Make sure not to exceed the maximum of 32 basic points.

- 1 In the AI Channel text field select the desired channel.
- 2 Enter the new basic points at the end of the Basic points table.
- 3 Activate the Sort command button for resorting the basic points.

Deleting individual basic points

- 1 In the AI Channel text field select the desired channel.
- 2 Position the cursor in the Basic points table at the basic point to be removed.
- 3 Delete the basic point by pressing the Delete command button.

Displaying the sensor characteristic

- 1 Select the desired channel in the AI Channel text field.
- 2 Click on the [Diagram] command button for graphic representation.
- 3 Under Option choose the desired display mode: Linear or Logarithmic.

Copying/exchanging AI channel parameters

- 1 Press the [Copy] command button.
- 2 Choose Copy for copying one channel to another or Exchange to exchange two channels.
- 3 Under Source and Dest enter the corresponding channel numbers and press the [OK] button.

3.8.4.4 Enter the characteristic curves (AI)

In this dialog box you can program the characteristics of the individual AI channels.




Fig. 3-123

Control elements**AI Channel**

Number of analog input channel.

AI Type

Type of AI channel (sensor type, max. 12 characters).

Unit

Unit of measure applicable to the AI channel (max. 12 characters).

Basic points

Basic points table that defines the AI-characteristic.

[Sort]

Sort the table of basic points.

[Copy]

Copy all parameters of the selected AI channel to another AI channel or exchange their characteristics.




Fig. 3-124

Copy control elements

Mode Selection

- Copy: Copy one channel to another.
- Exchange: Exchange two channels.

Source

Channel from which will be copied or that is to be exchanged.

Destination

Channel to which will be copied or that is to be exchanged with Source.

[Delete]

Delete the table entry selected by the cursor.

[Diagram]

Choose the Diagram command button to show the programmed sensor characteristic:




Fig. 3-125

Under Options you can choose the display mode:

Linear

Linear scale.

Logarithmic

Logarithmic scale. No negative values can be displayed on the Y-axis if [Options] > [Logarithmic] is chosen

X-Grid

Display the vertical grid for the X-axis.

Y-Grid

Display the horizontal grid for the Y-axis.

[Formula]

Formula enables you to enter the sensor characteristic by an algebraic formula:




Fig. 3-126

Enter Formula...

In this field, enter the formula that characterizes the sensor. Only the capital X is allowed as variable.

Minimal X

Lower limit for X in mV.

Maximal X

Upper limit for X in mV.

The formula describes the Y-value (any measure unit) as a function of X (voltage in mV). The following operations are allowed:

- Basic operations: + - */
- Calculation of e to the power of p: $\exp(p)$
- Calculation of 10 to the power of p: $\exp(p)$
- Calculation of the logarithm, base e: $\lg(n)$
- Calculation of the logarithm, base 10: $\log()$
- Squaring: $\text{sqr}()$
- Calculation of the square root: $\text{srt}()$
- Determination of the absolute value: $\text{abs}()$
- Expressions in brackets (...)

Blanks are not allowed; they will cause a syntax error to be reported.

Example

You want to measure with your pressure gauge in a range of 5E-8 to 5E-3 mbar. Its characteristic is given by the expression ' $p[\text{mbar}] = 10^{(U[\text{V}] * 1.667 - 10.5)}$ '.

- 1 'Translate' the expression to the form that Quadstar 32-bit expects. Note that the voltage must be given in mV. The formula therefore reads as follows:

' $\exp(X/1000 * 1.667 - 10.5)$ '

- 2 Determine the limits for minimal and maximal X as restrictively as possible (i.e. calculate the voltage for highest and lowest expected pressure) and enter them (in mV) as Minimal X and Maximal X. In our example, voltages between 2 and 5 Volts are expected:

Minimal X: 2000
Maximal X: 5000

Actually, you could forget this and simply set the largest possible range (-10 ... 10 Volts). But you would lose precision by doing so: the 32 basic points would be spread over the whole (too large) range.

In some cases the formula might not even be defined across the whole range (e.g. division by zero!). Formula and limits must always be set in a manner that results of the calculation (intermediate results too) remain in the range of $-1E+37$... $+1E+38$.

- 3 Start the calculation by [OK]. The program will return to the previous dialog box and the calculated points will occur in the table 'Basic Points'. Check the curve by using 'Diagram' and modify limits and formula if necessary or delete and add manually calculated points.

3.8.5 Screen Colors

[Config] > [Screen Color] is used for defining the Quadstar 32-bit specific desktop colors. With the aid of three test screens (Versus Time - Graphic, bar chart and table) colors can be selected individually for the elements of graphically displayed measurement data. The colors of the individual screen elements can be stored as a set, i.e. as color palettes, and reactivated via the Color Palette list box. The selected colors affect only the measurement data pictures of Quadstar 32-bit.

Quadstar 32-bit is supplied with a standard color palette. You can either adopt this palette or change the colors and store the palette under a different name.

Preparation of a suitable color palette is particularly important when outputting measurement data graphics on a printer. Colors that appear as white on a monochrome display should not be used. If a suitable color palette for the printout is defined, all visible elements will be drawn.

NOTE:

Quadstar 32-bit can display up to 192 colors, provided that the PC supports this.

3.8.5.1 Changing the colors of the individual screen elements

This change relates to the color palette selected under Color Palette. There are two ways of selecting the individual screen elements: indirect selection via the Screen Element list box, or direct selection by clicking on the desired screen element in one of the three test pictures.

The following screen elements are available in the Screen Element list box:

- Alarm Off: Alarm deactivated.
- Alarm On: Alarm activated.
- Axis: Horizontal and vertical axes.
- Axis Text: Numeric labeling of the axes.
- Axis Unit: Text Name of the vertical axis.
- Background: Background of the test pictures
- Bar Current: Currently measured value (Bargraph).
- Bar Last/Difference: Measured value of the last cycle (Bargraph)/difference (Analog).
- Bar Reference: Reference value (Bargraph/Analog).
- Cursor: Line cursor.

- Curve: Versus Time line.
- Grid: Grid lines.
- Table Frame: Frame for Multiple Table display.
- Table Head: Header line of the table.
- Table Line Separator: Line separator for Multiple Table display.
- Table Marker: Mark the current cycles for Multiple Tables display.
- Table Text: Contents of the table.
- Text Lightened: Emphasized text.
- Text: Normal Text.
- Trip Level A: Threshold A.
- Trip Level B: Threshold B

Selecting a color palette

Select the colors via the Color Palette field. You can either modify an existing color palette or create a new palette and store it under a different name.

- 1 Choose [Config] > [Screen Color].
- 2 Open the Color Palette list box and click on the desired color palette.

Indirect selection through the list box

- 1 Choose [Config] > [Screen Color].
- 2 Open the Screen Element list box and click on the desired screen element. (you can access the next element without opening the list box by pressing the up/down arrow key).
- 3 Assign a different color to the selected screen element by clicking on the desired color in the Colors field

Direct selection in the test picture

- 1 Choose [Config] > [Screen Color].
- 2 Click on the desired element in one of the test pictures. The selection is confirmed in the Screen Element list box.
- 3 Assign a different color to the selected screen element by clicking on the desired color in the Colors field.

Storing a color palette

- 1 Choose Save. Quadstar 32-bit responds with the Save Color Table dialog box.
- 2 Change the displayed name if you want to store a new color palette.
- 3 Activate the [OK] command button.

Screen Colors dialog box




Fig. 3-127

The Screen Colors dialog box is subdivided into the following Groups and control elements:

Control elements

Color Palette

Choose a color palette.

Screen Element

Choose a screen element via the list box.

Colors

Color range of the installed video card.

Test screens

- Versus Time graphic.
- Bar chart graphic.
- Table output.

NOTE:

The entered values are stored automatically when the Screen Color editor is closed with the [OK] command.

[Alternating Colors]

Choose the colors to be used alternatingly (used in TuneUp and Sequencer instructions, where no fixed colors can be set).

- #0 ... #63: The alternatingly used colors, maximum 64.
- [Set Default Colors]: Choose the standard colors (suitable for a bright background).

[Save]

Store the current state in the selected color palette.

[Delete]

Delete the selected color palette.

4. Display Saved Values

Display Saved Values is used for displaying and analyzing stored measurement data. Such files may contain intensities, concentrations, calculated values and measured values from analog inputs. With Display Saved Values you can select these blocks and elements individually.

Main menu

Choose the type of measure data to work on:



Fig. 4-1

Submenus

File

Display the version number of the Quadstar 32-bit Display Saved Values program and exit the program.

Scan

Load and evaluate the scan measurement data.

- Scan Analog measurement mode. See Chapter 4.2.1 Scan Analog data, 251.
- Scan Bargraph measurement mode. See Chapter 4.2.2 Scan Bargraph Cycles, 271.

Process

Load and evaluate MID (Multiple Ion Detection) or MCD (Multiple Concentration Detection) measurement data, or data that has been created by the sequencer instructions Process.

Auxiliary

Display zero gas data, mass scale calibration and error history.

Help

Display the 'Key Fragment Ions' table and help contents. See Chapter 1.9 Hints concerning the work with Quadstar 32-bit, 56.

4.1 File

File menu



Fig. 4-2

About Display Saved Values...

Info on the Display Saved Values program.

Exit

Quit the Display Saved Values program.

4.1.1 About Display Saved Values

Read out the version number of the Quadstar 32-bit Display Saved Values program:




Fig. 4-3

4.2 Scan measurement data

Scan menu



Fig. 4-4

Analog Data...

Display Scan Analog measurement data.

Analog Reference...

Display Scan Analog reference measurement data

Bargraph Cycles...

Display Scan Bargraph measurement data. See Chapter 4.2.2 Scan Bargraph Cycles, [271](#).

4.2.1 Scan Analog data

Analog Data are data, that contain one or more cycles. Analog Reference are data, that contain just one cycle and serve as reference data.

In the subsequent Sections we distinguish between 2-dimensional and 3-dimensional representation. The display mode can be changed over in the Setup menu.

- [3 Dimensional] checked: 3-dimensional representation.
- [3 Dimensional] not checked: 2-dimensional representation.

2-Dimensional representation

In this display mode only the measurement data of one cycle is displayed. The file can contain the measurement data of one or several cycles.




Fig. 4-5

The status bar shows the coordinates of the mouse. The function Magnify or Cursor can be assigned to the mouse. See Chapter 4.2.1.5.1 Magnify, 261 and Chapter 4.2.1.5.2 Cursor, 263.

3-Dimensional representation

With this display mode, Scan Analog measurement data that have been measured and stored across several cycles can be represented in a 3-dimensional diagram.

Depending on the choice of X-Axis Type under [Setup] > [Graphic], the X-axis can show, for example, the mass numbers...




Fig. 4-6

...the measured cycles...




Fig. 4-7

...or also the relative and the absolute time.

The status bar shows the coordinates of the mouse.

NOTE:

Please note that 3-dimensional representation is only feasible, if several cycles of a scan have been stored in the measurement data file.

Analog data menu



Fig. 4-8

File

Load, save or ASCII-convert a file. Return to the main menu of the evaluation program. See Chapter 4.2.1.1 File submenu, 254.

Display

Display selected intensities of a scan numerically, show/hide added graphic elements. See Chapter 4.2.1.2 Display submenu, 255.

Select

Select the scans that are to be displayed. See Chapter 4.2.1.3 Select submenu, 256.

Setup

Set up display range and options. See Chapter 4.2.1.4 Setup submenu, 257.

Function

Set the cursor function ('magnify'/'show position'), call the graphic editor.
Chapter 4.2.1.5 Function submenu, 260.

Special

Show data using standard settings, call the filter function, and reset magnification. See Chapter 4.2.1.6 Special submenu, 269.

Info

Display measure data documentation. See Chapter 4.2.1.7 Info submenu, 270.

4.2.1.1

File submenu



Fig. 4-9

Open Data...

Open a measurement data file (*.sac).

Open Reference...

Open a reference data file (*.sar).

Load Reference...

Load and display a reference file for the currently displayed data file.

Close Reference...

Close the reference file for the currently displayed data file.

Save Difference As...

Save the difference between the measure data and the reference data.

Reload Cycles

Reload newly measured cycles of the currently displayed data file.

Reload Setup...

Reload newly measured cycles automatically:




Fig. 4-10

- Number of Cycles: Number of cycles to reload (of the last measured ones).
- Automatic Reload: Reload automatically on/off.
- Repeat Time[s]: Time between automatic reloads.

Save

Store the active measurement data file under its current name.

Save As...

Store the current measurement data file under a different name and path

Save Selected Cycles As...

Store the cycles selected under [Setup] > [Graphic] as a new file.

Convert to ASCII...

Convert the current Scan Analog file to an ASCII file.

Convert to Data...

Convert a Scan Analog reference data file to a Scan Analog measurement data file.

Convert to Reference...

Convert a Scan Analog measurement data file to a Scan Analog reference data file.

Close

Return to the main menu of Display Saved Values.

4.2.1.2 Display submenu

Display the intensity of a selectable mass number, show or hide added graphic elements:



Fig. 4-11

Value

Numerically display or change the intensity of a selected mass number in the dialog box Scan Analog Display. See Chapter 4.2.1.5.2 Cursor, 263.

Graphic Data

Display the graphic elements that have been added to the measurement data picture by the Graphic Editor. See Chapter 4.2.1.5.3 Graphic editor, 267.

4.2.1.3 Select submenu



Fig. 4-12

Scans/Additional Data

Choose those scans or additional data that shall actually be displayed. Opens the dialog box Select Scans/Additional Data:




Fig. 4-13

Control elements

Scans

- Scan Number: Choose the scan.
- Display: Show or hide the chosen scan.
- Display All: Show all scans.
- Comment: Short comment (will be displayed together with this scan).

Scan Data

- Scan Name: Mass Axis Name
- Scan Unit: Mass Axis Unit
- Start Value: First value of Mass Axis
- End Value: Last value of Mass Axis
- Value Name: Intensity Axis Name
- Value Unit: Intensity Axis Unit
- Comment: Short Comment (will be displayed together with this scan)

Additional Data

Additional measurement data (if present)

- Name: Designation of the additional measurement data.
- Unit: Unit of measure of the additional measurement data.
- Display: Display additional measurement data on/off.

4.2.1.4 Setup submenu




Fig. 4-14

3 Dimensional

Display several Scan Analog spectra as a 3-dimensional diagram. This option is only available if the file contains the measurement data of more than one cycle.

Graphic...

Setup the display.

Fill Curve

Fill the envelope with the Fill Color.

Curve Color...

Enter the Curve Color and the Fill Color:




Fig. 4-15

- Curve Color: Select the color of the envelope.
- Fill Color: Select the color to fill the envelope with.

Arrangement

Arrange scans, if several of them are stored in this file:

- Vertical: Display scan below previous.
- Horizontal: Display scan beside previous.

The following menu entries are only relevant for 3-dimensional representation of measurement data:

Cursor

Show the cursor projection lines.

- Cross: Normal cursor.
- Cross & Projection: Projection of the cursor line to the X-, Y- and Z-axis.

Draw Box

The graphic is drawn in a box.

Align Cursor

The cursor jumps from cycle to cycle.

Graphic...

The type of dialog box that is displayed depends on the setting of [Setup] > [3 Dimensional]: enabled or disabled.


2-Dimensional Representation

Fig. 4-16

NOTE:

The dialog box can also be called by double-clicking the left mouse button inside the displayed measurement data.

Scan

- Scan Number: Choose the scan, if more than one scan is stored in the file.
- Start Value: Enter the starting value at which the graph should start.
- End Value: Enter the last value to be displayed.

Grid

- X-Axis: Display the vertical grid for the X-axis.
- Y-Axis: Display the horizontal grid for the Y-axis.

Cycle

- Number: Select the cycle, if the file contains several cycles.

Value

Name (Ion current, Ion Counts, ...) and unit (A, cps, ...) of measure data.

- Highest Value: End value of Y-axis.
- Lowest Value: Start value of Y-axis.

Scale

- Lin: Display the measurement data on a linear scale.
- Log: Display the measurement data on a logarithmic scale.

NOTE:

Number-parameters (e.g. Start Value/End Value) can be restored to the original values by clicking on the or arrows with the right mouse button

3-Dimensional representation




Fig. 4-17

NOTE:

This dialog box can also be called up by double clicking the left mouse button inside the displayed measurement data.

Cycle

Select the cycles to be displayed

- Start Cycle: First cycle to be displayed
- End Cycle: Last cycle to be displayed

X-Axis Type

The following data types may be shown on the X-axis:

- Cycle: The X-axis is labeled with the cycle numbers.
- Scan: The X-axis is labeled with scan values.
- Data: The X-axis is labeled with additional measurement data.
- Time Relative: The X-axis is labeled with the relative time since the start.
- Time Absolute: The X-axis is labeled with the absolute time (time of day).

Z-Axis Type

The following data types may be shown on the Z-axis:

- Cycle: The Z-axis is labeled with the cycle numbers.
- Scan: The Z-axis is labeled with the scan values.
- Data: The Z-axis is labeled with additional measurement data.

Scan

Name (Mass,...) and unit (amu,...) of the scan.

- Start Value: First cycle to be included in the displayed graph
- End Value: Last cycle to be included in the displayed graph.

Data

Additional data with name and unit of measure.

- Start Value: Starting value of the display.

- End Value: Ending value of the display.

Display

- Angle between X- and Z-axis: Enter the angle between the X- and the Z-axis
- Length of Y-axis in % of X-axis: Length of the Y-axis in % of the X-axis.
- Length of Z-axis in % of X-axis: Length of the Z-axis in % of the X-axis.

Grid

- X-Axis: Display the vertical grid for the X-axis.
- Y-Axis: Display the horizontal grid for the Y-axis.
- Z-Axis: Display the horizontal grid for the Z-Axis.

Data Reduction

Reduction of the displayed points by a factor of 2, 4 or 8.

Value

Name (Ion Current, Ion Counts ,...) and unit (A, cps, ...) of measure data.

- Highest Value: End value of Y-axis
- Lowest Value: Start value of Y-axis.

Scale

- Lin: Display the measurement data on a linear scale.
- Log: Display the measurement data on a logarithmic scale.

NOTE:

Number-parameters (e.g. Highest Value/Lowest Value) can be restored to the original values by clicking on the or arrows with the right mouse button

4.2.1.5

Function submenu



Fig. 4-18

To select the cursor function type or call the graphic editor, open the Function menu:

Magnify

Switch to block marking for the magnifier function.

Cursor

Switch to crosshairs.

Graphic Editor

Call the graphic editor.

NOTE:

Changeover from Magnify to Cursor and vice versa is also possible with a click of the right mouse button anywhere within the measurement data screen.

4.2.1.5.1 Magnify

Magnifying a screen section

If [Function] > [Magnify] is chosen, a section of the displayed screen can be marked and magnified:

2-Dimensional representation

Hold down the left mouse button and drag the mouse to mark a rectangle:




Fig. 4-19

If you click the right mouse button while still holding the left mouse button pressed, the selected rectangle is magnified on the display.




Fig. 4-20

- This magnify function can be nested several times.
- The last picture magnification can be reversed using [Special] > [Undo Magnify].
- The graphic can be restored to its original size by [Special] > [Set Default].

3-Dimensional representation

Hold down the left mouse button and drag the mouse to mark a rectangular solid:




Fig. 4-21

If you click the right mouse button while still holding the left mouse button pressed, the selected part is magnified on the display.




Fig. 4-22

In this way it is possible e.g., to view the peaks that have been concealed by the peak of mass 20 before:

4.2.1.5.2 Cursor

Determine measurement data graphically

If [Function] > [Cursor] is selected, a measured value within the graphic display can be marked by the crosshairs. Its values for X and Y are displayed in the status bar.

2-Dimensional representation

Press the left mouse button within the graphic display, so the cursor appears as crosshairs and the values for X and Y are displayed in the status bar:




Fig. 4-23

If you click the right mouse button while still holding the left mouse button pressed, the mass selected with the vertical line of the crosshairs (X: 28.00) and the corresponding intensity are displayed numerically in a dialog box.




Fig. 4-24

The mass number can be selected directly, the intensity is modified by Change:

Control elements

Scan Number

Number of the selected scan (if more than one scan is available).

Cycle

Number of the cycle (if more than one cycle is available).

Cycle Time

Relative start time of the cycle.

Date

Date and time when the scan was started.

Mass

Selected Mass number.

Ion Current/Countrate

Intensity of the selected mass number.

[Change]

Manual elimination of freak values. See Section «Manual elimination of freak values»,
 302.




Fig. 4-25

Faulty measurements (“freak values”) can also be corrected automatically by [Special] > [Filter].

3-Dimensional representation

Press the left mouse button within the displayed graphic, so the cursor is displayed as crosshairs. The values of the current cursor position appear in the status bar:




Fig. 4-26

About cursor operation:

- The cursor always moves to the scan of the first cycle. If the Ctrl key is pressed simultaneously, the cursor can be shifted on the Z-axis.
- If [Setup] > [Align Cursor] is active and the Ctrl key is pressed, the cursor jumps from cycle to cycle on the Z-axis.

- When a peak top is marked within a group of peaks, it may be difficult to find the correct peak in the correct cycle. By using [Setup] > [Cursor] > [Cross & Projection] this can be done much easier since the cursor is projected to the different axes then.

If you click the right mouse button while still holding the left mouse button pressed, the mass selected with the vertical line of the crosshairs (for example, mass 16.00) and the corresponding intensity are displayed numerically in a dialog box:




Fig. 4-27

Control elements The selected mass number can be modified directly, the intensity is modified by Change.

Scan Number

Number of the selected scan (if more than one scan is available)

Cycle

Selected cycle.

Cycle Time

Relativ start time of the cycle.

Date

Date and time when the scan was started.

Mass

Selected mass number.

Ion Current/Countrate

Intensity of the selected mass number.

[Change]

Manual elimination of freak values. See Section «Manual elimination of freak values», 302.

Additional Data

Additional data (one value per cycle) stored with this cycle.

4.2.1.5.3 Graphic editor

You can add text and drawings to your measurement pictures, for example, to prepare them for printout or as an aid for interpretation. The dialog box shown below contains ten functions, each of which is identified by a symbol:



Fig. 4-28



Insert text

Position the insert marker where the text should start and click the left mouse button. The following text input dialog box is displayed:




Fig. 4-29

In addition you can specify the desired Font and the background mode (Background mode: normal overprints existing picture elements, transparent only superposes the characters), Delete the entered text, if desired, and accept ([OK]) or cancel ([Cancel]) the input.



Draw line

Position the crosshairs where the line should start. While keeping the left mouse button pressed drag the mouse to the point where the line should end and then release the mouse button.



Draw rectangle

Position the crosshairs at one corner of the rectangle. While keeping the left mouse button pressed drag the mouse to the diagonally opposite corner and then release the mouse button.



Draw solid bar

Position the crosshairs at one corner of the solid bar. While keeping the left mouse button pressed drag the mouse to the diagonally opposite corner and then release the mouse button.



Choose options

Click on the graphic element to be processed or choose one of the above functions. If you now click on the Options button and depending on the graphic element, a dialog box for Draw lines is displayed, for example as the one shown below:




Fig. 4-30

You can now adjust the options (for example, color and style) of the graphic element to your requirements and accept the setting with [OK] or quit with [Cancel].



Resize element

Resize element: Click on the graphic element to be processed which is then displayed with eight handles. While holding down the left mouse button drag one of the handles (larger/smaller) or move the element itself. Release the mouse button.



Bring to front

Click on the graphic element to be processed which is then displayed with eight handles. When you click on this command the marked element will be brought forward by one layer relative to the next higher graphic element, but not relative to the measurement data.



Bring to back

Send to back: Click on the graphic element to be processed which is then displayed with eight handles. When you click on this command the marked element will be sent one layer back relative to the next lower graphic element, but not relative to the measurement data.



Delete

Click on the graphic element to be processed which is then displayed with eight handles. Click on the [Delete] command to remove the marked graphic element. If no graphic element is marked, all graphic elements are removed. To prevent inadvertent removal the operation has to be conformed with [OK] or [Cancel] in the subsequent dialog box.



Quit

Click on this command button to quit the graphic editor.

NOTE:

Do not change the size of the measurement picture after you have processed it with the graphic editor. This would have the result that the graphic elements would be shifted relative to the measurement curve and consequently be positioned in the wrong place.

4.2.1.6 Special submenu



Fig. 4-31

Undo Magnify

Reverse the changes of the last magnification.

Set Default

Restore all settings for measurement data display to their standard values.

Filter...

Automatic elimination of freak values.

4.2.1.6.1 Filter

Choose filter for automatic elimination of freak values:




Fig. 4-32

Filter Type

Select the desired filter type.

- **Moving Median:** Determines a moving median. For this purpose a certain number (selectable via Rank) of the last measured values is sorted by their magnitude and the one in the middle is used as the filter output (this way, extremely high or low values do not affect the results at all, since they are never taken into consideration).
- **FIR QMS...:** The individual measured values are subjected to the same FIR filter algorithm (Finite Impulse Response) that is implemented in the QMS.

Moving Median

- **Rank:** Determines the number of measured values taken into consideration for the filter type Moving Median (number of measured values = $2 * \text{Rank} + 1$).

Output File Name

Name of the file in which the results are to be written. The name can be entered by activating the File Manager command button.

NOTE:

Record the filtered data on a new file so that the original data will be preserved in any case.

4.2.1.7 Info submenu

Display measure data documentation



Fig. 4-33

File...

Display measure data documentation, if there is one. Otherwise a new documentation is created that consists of additional info, notes and title only:




Fig. 4-34

Control elements**[Measurement Parameters]**

Display the measurement or channel parameters in effect at the time of the measurement.

[Ion Source Parameters]

Display the ion source parameters in effect at the time of the measurement.

Additional Information

Display the user-defined parameter sets, each comprising Name and Value. These names and values can be added or modified later.

- [Delete All Names]: Delete all names of the user-defined parameters.
- [Delete All Values]: Delete all values of the user-defined parameters.

Notes

User notes/appendices to the measurements can be added directly by positioning the insertion mark within the notice field and entering the text on the keyboard.

Title

File title (shows up in the file manager).

[OK]

Accept and save the file information.

[Cancel]

Cancel the input, leave the old information unchanged.

NOTE:

It is important that you quit this dialog box with [Cancel] if you have inadvertently changed something. (for example, name of additional information deleted). In the subsequent dialog box (“...save current changes?”) enter [NO], you can then recall the original file information from [Info] > [File].

[About File]

Display certain specifications of the measurement data file.




Fig. 4-35

Filename

Name of the selected measurement data file.

Software Version

Quadstar 32-bit version used for creating this file.

Date

File creation date.

Time

Time of day the file has been created.

Author

Operator name.

4.2.2 Scan Bargraph Cycles

Menu

Fig. 4-36

Submenus**File**

Load, save, ASCII-convert or print a file. Return to the main menu.

Mode

Select the display mode:

- Total Ion Current (TIC)
- Intensity Versus Time (IVT)

Display

Display all mass numbers and peak intensities of a scan in numeric and graphic form, show/hide added graphic elements. See Chapter 4.2.2.3 Display submenu, 275.

Select

Select the mass numbers for which graphic representation of the peak intensities is desired. See Chapter 4.2.2.4 Select submenu, 277.

Setup

Set up the X-axis and Y-axis parameters. See Chapter 4.2.2.5 Setup the display, 279.

Function

Define the cursor type or function and start the graphic editor. See Chapter 4.2.2.6 Function submenu, 281.

Special

Show data using standard settings or reset magnification. See Chapter 4.2.2.7 Special submenu, 284.

Info

Display measure data documentation or supplementary information of the opened measurement data file. See Chapter 4.2.2.8 Info submenu, 284.


4.2.2.1**File submenu**

Fig. 4-37

Open...

Call the Quadstar 32-bit file manager for opening a measurement data file.

Reload Cycles

Reload newly measured cycles of the currently displayed data file.

Reload Setup...

Reload newly measured cycles automatically:




Fig. 4-38

- Number of Cycles: Number of cycles to reload (of the last measured ones).
- Automatic Reload: Reload automatically on/off.
- Repeat Time[s]: Time between automatic reloads.

Save

Store the active measurement data file under its current name.

Save As...

Store the active measurement data file under a different file name and/or path.

Print

Print the selected masses and their peaks of the selected cycles.

Convert to ASCII...

Convert the current Scan Bargraph file to an ASCII file.

NOTE:

The format of the ASCII file is described in more detail in Chapter 10.5 Data export ASCII, 549.

Close

Return to the Display Saved Values main menu.

4.2.2.2

Mode submenu



Fig. 4-39

Scan Bargraph measurement data can be represented in two formats:

TIC

Total Ion Current: All measured intensities per cycle are added. The results of this addition are displayed as a function of time.

IVT

Intensity Versus Time: The peak intensities of the chosen masses are displayed as a function of time.

4.2.2.2.1 TIC display mode (Total Ion Current)

This window shows the Total Ion Current [A] (calculated by adding all measured peak intensities per cycle) across a selectable time range. The time range, i.e. the number of cycles to be displayed, can be selected under [Setup] > [X-Axis].




Fig. 4-40

When displayed for the first time, the screen is automatically stretched in the Y-Axis; the display range is reduced to the range, in which the sum of all ion currents (TIC) varies. This display area can be changed by [Setup] > [Y Axis] or by double clicking on the Y-Axis.

The measured cycles are displayed on the X-Axis. The display area can be changed to Time Relative or Time Absolute in the Setup menu under X-axis (or by double clicking the X-axis label).

4.2.2.2.2 IVT display mode (Intensity Versus Time)

The picture shows the peak intensities of the masses chosen under Select Data, across a selectable time range. This time range can be defined under [Setup] > [X-Axis] by specifying the number of cycles to be displayed.




Fig. 4-41

4.2.2.3 Display submenu



Fig. 4-42

Cycle...

Numerically display the Peak intensities of a selectable scan/cycle.

Graphic Data

Display the supplementary graphic elements that have been added to the measurement data picture with the Graphic Editor.

4.2.2.3.1 Cycle




Fig. 4-43

Cycle

Number of the cycle, of which the spectrum is to be displayed.

Cycle Time

Relative measurement data time (time, since the measurement has started).

Date

Date and absolute time of the measurement.

TIC

Total Ion Current of the selected spectrum.

[Spectrum]

Graphic representation of the measurement data.

[Library]

Call the spectra library. You may compare the current spectrum with spectra of the library and save spectra. (=> Command button Spectrum, Menu [Spectrum] > [S1 to Library])

NOTE:

All peaks of a spectrum are listed, even if under [Select] > [Data] not all peaks have been selected.


[Spectrum]

Fig. 4-44

Setup

- Y-Axis: Setup the display range of the Y-Axis.
- X-Axis: Setup the mass range to display.

You may access these setup boxes as well by double clicking on the corresponding axis.

Print

- Cycle: Print the peaks of the displayed spectrum.

Library

- Manual: Call the spectrum library.

4.2.2.4**Select submenu**

Fig. 4-45

Data...

Dialog box for selecting the mass numbers whose peak intensities are to be displayed in IVT mode.

NOTE:

This dialog box can also be called up by double-clicking the left mouse button in the measurement data display.

4.2.2.4.1**Data**

In the Select Data dialog box you can select the mass numbers of which the intensities are to be displayed in IVT mode. The "intensity" shown in graphic mode is of a static nature and represents the frequency in percent (Freq. [%]), with which the peak was found in the selected cycles.

Remove, selected color etc. always relate to the currently marked peak. To select the color, click on the desired field in the color palette.

By Select, several masses can be selected simultaneously:




Fig. 4-46

Mass

Mass number of the selected peak.

Freq. [%]

Frequency in % at which the peak was measured in the selected number of cycles, where 100% means that the peak was found in every cycle, 50% means that the peak was found only in half the cycles.

Color

Peak color.

Marker

Marker for monochrome display. To get this marker on the menu, activate the Marker option under [Setup] > [Y-Axis].

[Remove]

Remove the marked mass number from the display table.

[Remove All]

Remove all mass numbers from the display table.

[Select]

Choose peaks of a range.

[Select]




Fig. 4-47

Control elements**From Mass**

Starting mass number.

To Mass

Ending mass number.

Threshold

Threshold value in % that must be exceeded for the peak to be selected.

NOTE:

A new peak can be selected via **Select** or via the displayed spectrum: hold down the left mouse button and click simultaneously on the desired peak with the right mouse button.

4.2.2.5**Setup the display**

Fig. 4-48

X-Axis...

Determine the display range of the X-axis.

Y-Axis...

Determine the display range of the Y-Axis.

NOTE:

The dialog boxes for setting up the X-axis and Y-Axis can also be activated by double clicking the lettering of the corresponding axes.

Number-parameters (e.g. Start Cycle/End Cycle) can be restored to the original values by clicking on the [UP] or [DOWN] arrows with the right mouse button

4.2.2.5.1 X-Axis




Fig. 4-49

Start Cycle

First cycle to be displayed.

End Cycle

Last cycle to be displayed.

X-Scale

The X-axis can be lettered according to the following types:

- Cycle: Cycle numbers.
- Time Relative: Relative time since the start.
- Time Absolute: Absolute time (time of day).

X-Grid

Display the vertical grid for the X-axis.

4.2.2.5.2 Y Axis




Fig. 4-50

Highest Value

End value for lettering the Y-Axis.

Lowest Value

Start value for lettering the Y-Axis.

Y-Scale

- Lin: Display the measurement data on a linear scale.

- Log: Display the measurement data on a logarithmic scale.

Line Type

- Solid: Display the measurement data with a solid line.
- Dot: Display the components:
 - Dots: when Marker is switched off.
 - Symbols: when Marker is switched on.

Y-Grid

Display the horizontal grid for the Y-Axis.

Marker

Display the symbols.

In [Mode] > [TIC], you can choose the Type of the marker symbol:




Fig. 4-51

4.2.2.6

Function submenu




Fig. 4-52

You can select the type of cursor function via the Function menu:

Magnify

Magnify a screen section by marking a block with the cursor.

Cursor

Display the scan measurement data numerically.

NOTE:

Switching from Magnify to Cursor and vice versa is also possible directly with a click of the right mouse button within the picture of the measurement data.

Table

Display the measurement data as a table. See Chapter 4.3.1.2.1 Cycle, 289.

Spectrum

Display the measurement data graphically as a spectrum. See Chapter 4.2.2.6.4 Display format of [Function] > [Spectrum], 284.

Graphic Editor

Call the graphic editor. See Chapter 4.2.2.6.5 Graphic editor, 284.

4.2.2.6.1 Magnify a Screen Section

If [Function] > [Magnify] is chosen, an area of the measurement data screen can be marked that is to be magnified:




Fig. 4-53

The cursor position is shown on the status bar. C: XX represents the cycle number; the position on the Y-Axis [A] can be read under Y. Now press the right mouse button while still holding the left mouse button pressed to magnify the marked area.




Fig. 4-54

This magnifier function can be nested several times.

The graphic can be restored to its original size by [Special] > [Set Default].

4.2.2.6.2 Determine measurement data graphically (Cursor)

If [Function] > [Cursor] is active, you can select a cycle with the left mouse button.




Fig. 4-55

When the right mouse button is pressed, while still holding the left mouse button pressed, the measurement data are displayed as a table (if [Function] > [Table] is set) or as a spectrum (if [Function] > [Spectrum] is set).

4.2.2.6.3 Display format of [Function] > [Table]




Fig. 4-56

This table is the same as the one that is called by [Display] > [Cycle]. See Chapter 4.2.2.3.1 Cycle, 276.

4.2.2.6.4 Display format of [Function] > [Spectrum]




Fig. 4-57

4.2.2.6.5 Graphic editor

You can add text and drawings to your measurement pictures, for example, to prepare them for printout or as an aid for interpretation. See Chapter 4.2.1.5.3 Graphic editor, 267.

4.2.2.7 Special submenu



Fig. 4-58

Undo Magnify

Undo the last magnification.

Set Default


Restore all measurement data display settings to their standard values. All components and measurement data blocks are displayed across the stored range (all cycles).

4.2.2.8 Info submenu



Fig. 4-59

File...

Display measure data documentation, if there is one. Otherwise a new documentation is created that consists of additional info, notes and title only. See Chapter 4.2.1.7 Info submenu,  270.

4.3 Process

Process menu



Fig. 4-60

Cycles...

MID or MCD measure data that contain one or more cycles.

Reference...

MID or MCD reference measure data that contain just one cycle.

4.3.1 Cycle/Reference measurement data

Cycle data are measurement data that contain one or more cycles of the Multiple Ion Detection (MID) or the Multiple Concentration Detection (MCD) measurement type. Reference data are actually the same as cycle data, except that they consist of just one cycle. By the sequencer instruction Process, Cycle- or Reference- data may be created out of random data records.




Fig. 4-61

- The Y-axis lettering always relates to the selected measurement data block (in the above example the name of the selected block is: Ion Current). The names (here: mass numbers) of the components that belong to this block are marked by a short vertical line in front of them.
- If several blocks have been selected under [Select] > [Blocks], the Y-axis display can be assigned to a different block via the list box displayed in the upper left corner:




Fig. 4-62

- The X-axis is lettered with the cycle numbers.

Data block

The data are summarized and displayed (separated by type and origin) as blocks. Several blocks may be displayed simultaneously.

QMU 112

If the data have been collected by a QMU 112, the QMU-channel is shown between brackets behind the mass number. For example 28 [7] signifies Mass 28, measured via QMU channel 7.

Menu



Fig. 4-63

Submenus

File

Load, reload, save, print or ASCII-convert a file; return to the main menu of the Dispsav program. See Chapter 4.3.1.1 File submenu, 287.

Display

Display the selected components of a block. See Chapter 4.3.1.2 Display setup, 288.

Select

Select the measurement data groups (blocks) and the components that are to be displayed. See Chapter 4.3.1.3 Select data to be displayed, 294.

Setup

Set up the display of the X-axis and Y-axis. See Chapter 4.3.1.4 Setup the display, 298.

Function

Define the cursor function, call the graphic editor. See Chapter 4.3.1.5 Function submenu, 299.

Special

Restore the standard display mode, call the filter function. See Chapter 4.3.1.6 Special submenu, 303.

Info

Display measure data documentation. See Chapter 4.3.1.7 Info submenu, 304.

4.3.1.1 File submenu




Fig. 4-64

Open...

Call the File Manager for opening a measurement data file.

Reload Cycles

Reload newly measured cycles of the currently displayed data file.

Reload Setup...

Reload newly measured cycles automatically:




Fig. 4-65

- Number of Cycles: Number of cycles to reload (of the last measured ones).
- Automatic Reload: Reload automatically on/off.
- Repeat Time[s]: Time between automatic reloads.

Save

Store the current measurement data file.

Save As...

Store the current measurement data file under a different file name and/or path.

Print

Print the displayed cycles as a table.

Convert to ASCII...

Store the displayed measure data as an ASCII file under a user-assigned name.

NOTE:

The ASCII file format is described in Chapter 10.5 Data export ASCII, 549.

Close

Return to the Display Saved Values main menu.

4.3.1.2 Display setup



Fig. 4-66

Cycle...

Display the Cycle Data Display dialog box.

Statistic...

Display statistical information.

Integral...

Integration, i.e. summing of the displayed measured values across time.

Bargraph

Display an individual cycle as a bargraph.

IVT

Display several cycles in Intensity Versus Time mode.

Graphic Data

Display supplementary graphic elements.

4.3.1.2.1 Cycle

Display the Cycle Data Display dialog box with the measured values of the chosen data block. See Chapter 4.2.2.6.2 Determine measurement data graphically (Cursor), 282.

4.3.1.2.2 Statistical Functions

Statistic provides statistical information on the stability, reproducibility etc. of the measured data. The statistical calculations are performed block by block for the data records from Start Cycle to End Cycle :

Name	Minimum	Maximum	Mean	STD ABS	STD REL
14.00	1.776E-09	4.541E-09	3.275E-09	1.027E-09	3.137E-01
16.00	6.637E-10	1.578E-09	1.157E-09	3.389E-10	2.928E-01
18.00	5.195E-10	8.718E-10	7.039E-10	1.299E-10	1.846E-01
28.00	2.388E-08	6.238E-08	4.472E-08	1.431E-08	3.199E-01
32.00	5.444E-09	1.443E-08	1.029E-08	3.349E-09	3.255E-01
40.00	3.593E-10	8.350E-10	6.169E-10	1.792E-10	2.905E-01
44.00	1.083E-10	1.457E-10	1.284E-10	1.371E-11	1.067E-01

Fig. 4-67

The dialog box displays the statistical information of the Ion Current block. The number of cycles (Start Cycle and End Cycle) can be changed by [Setup] > [X Axis].

To display the other activated data blocks, press [Next] or [Previous]:

[Previous]

Display the statistical information of the previous block.

[Next]

Display the statistical information of the next block.

Block related information

The Title, the mass Unit and the Comment can be entered in the corresponding input fields of the Select menu under Data. They always relate to the selected data block. See Chapter 4.3.1.3.2 Select Data, 296: Title, Unit and Comment parameters.

Calculating statistical data

Name

Designation of the data (concentrations Ar, CO₂, ...; intensities 28, 32, ...; AI measured values, etc.)

Minimum

Smallest measured value stored between Start Cycle and End Cycle:

$$\text{Minimum} = \text{Minimum}(x_1, x_2, x_3, \dots, x_n)$$

$x_1 \dots x_n$ = Measured values 1...n (beginning at Start Cycle, ending at End Cycle)

Maximum

Largest measured value stored between Start Cycle and End Cycle:

$$\text{Maximum} = \text{Maximum}(x_1, x_2, x_3, \dots, x_n)$$

$x_1 \dots x_n$ = Measured values 1...n (beginning at Start Cycle, ending at End Cycle)

Mean

Mean value of the measured values stored between Start Cycle and End Cycle:

$$\text{Mean} = \left| \frac{\sum x_i}{n} \right|$$

x_i = Measured values i

n = Number of measured values

STD ABS

Absolute standard deviation of the measured values stored between Start Cycle and End Cycle:

$$\text{STDABS} = \sqrt{\left| \frac{\sum x_i^2 - n \cdot \text{Mean}^2}{n - 1} \right|}$$

x_i = Measured values i

n = Number of measured values

STD REL

Relative standard deviation of the measured values stored between Start Cycle and End Cycle:

$$\text{STDREL} = \frac{\text{STDABS}}{\text{Mean}}$$

4.3.1.2.3 Integration (Integral)

Integral integrates (from Start Cycle to End Cycle) across the time all selected signals of the chosen data block:




Fig. 4-68

Zero Based shows the integral over the x axis:




Fig. 4-69

Curve Based shows the integral over the line between the first and the last displayed measure value:




Fig. 4-70

The number of cycles (Start Cycle and End Cycle) can be changed under [Setup] > [X Axis].

To display other activated data blocks, press the [Next] or [Previous] command button:

[Previous]

Previous block.

[Next]

Next block.

The displayed unit of measure (in this picture As) is composed of Ampere (A, defined under Unit) and seconds (s).

Block related information

The Title, the mass Unit and the Comment can be entered in the corresponding input fields of the Select menu under Data. They always relate to the selected data block. See Chapter 4.3.1.3.2 Select Data, 296: Title, Unit and Comment parameters.

4.3.1.2.4 Bargraph

Bargraph displays an individual cycle as a bar chart. If there are several cycles, one can be selected under [Display] > [Cycle] (in this example: Cycle No. 183):




Fig. 4-71

4.3.1.2.5 Intensity Versus Time (IVT)

IVT displays the cycles selected under Select Data across a selectable time range. This display mode can only be chosen if more than one cycle is available, otherwise only Bar-graph display mode is possible.




Fig. 4-72

4.3.1.2.6 Graphic Data

Display the supplementary graphic elements that have been added to the measurement data picture with the Graphic Editor.

4.3.1.3 Select data to be displayed



Fig. 4-73

Choose here the blocks and components to be displayed:

Blocks...

Select the blocks available for display from the opened measurement data file.

Data...

Select the block components to be displayed.

4.3.1.3.1 Select Blocks

Quadstar 32-bit stores the different measurement data of a cycle measurement (intensities, analog input values, etc.) in separate blocks within the data file. Before the evaluation of the measurement data, the available blocks can be selected by Select Blocks:




Fig. 4-74

Control elements

Filename

Name of the file.

Available Blocks

Title of the available data blocks.

Selected Blocks

Measurement data blocks that have already been selected.

Title

Title of the measurement data block marked under Available Blocks or Selected Blocks.

Comment

Comments field of the measurement data block marked under Available Blocks or Selected Blocks.

[Select]

Add the block selected under Available Blocks to Selected Blocks. This function can also be activated by double-clicking the block to be selected in the Available Blocks field.

[Remove]

Remove the block selected under Selected Blocks.

[Select Data]

Call the Select Data dialog box. This function can also be activated by double-clicking the corresponding block in the Selected Blocks field.

4.3.1.3.2 Select Data

In the Select Data dialog box you can choose individual components of the selected data block. In Select Data you can also set display properties such as colors, markers, etc.:




Fig. 4-75

NOTE:

This dialog box can also be called up by double-clicking with the left mouse button on a mass number in the measurement data window.

Control elements

[Previous]

Previous measurement data block.

[Next]

Next measurement data block.

[Select Blocks]

Branch to the Select Blocks dialog box.

The following input fields apply to a measurement data block:

Title

Title of the measurement data block.

Comment

Comment for the measurement data block.

Unit

Unit of measure (A, etc.) for the measurement data block.

The following fields relate to the individual components of the displayed measurement data block:

Name

Input fields for displaying/modifying the names of the individual components.

Display

Control field for switching the measurement data display of the individual components on/off.

Color

Color in which the current component will be displayed. Click on the color bar to choose a different color.

Marker

Marker symbol for the current component.

The following values have strictly informative character and cannot be changed.

Level A

Threshold of switching function A.

Level B

Threshold of switching function B.

Thr.-Type

Switching function type.

The following parameter groups are used for display control:

Line Type

- Solid: Display the components with a solid line.
- Dot: Display the components with:
 - Dots, if Display Marker is switched off.
 - With the symbols selected under Marker, if Display Marker is on.

Y-Scale

- Lin: Display the measurement data on a linear scale.
- Log: Display the measurement data on a logarithmic scale.

Range

- Low. Value: Start value for the Y-axis labeling.
- High. Value: End value for the Y-axis labeling.

NOTE:

Number parameters (e.g. Start Cycle/End Cycle) can be restored to the original values by clicking on the [UP] or [DOWN] arrows with the right mouse button.

Display Thresh.

Display the switching thresholds.

Display Marker

Display the symbols selected under Marker.

Display Grid

Display the horizontal Y-axis grid.

4.3.1.4 Setup the display



Fig. 4-76

X-Axis...

Enter the parameters for selecting the display range applicable to the X-axis.

Y-Axis...

Enter the parameters for selecting the display range applicable to the Y-axis. [Setup] > [Y-Axis] opens the same dialog box as [Select] > [Data]. See Chapter 4.3.1.3.2 Select Data, 296.

NOTE:

The dialog boxes for setting up the X-axis and Y-Axis can also be activated by double clicking the lettering of the corresponding axes.

Number parameters (e.g. Start Cycle/End Cycle) can be restored to the original values by clicking on the [UP] or [DOWN] arrows with the right mouse button.

4.3.1.4.1 X-Axis




Fig. 4-77

Display Range

The display range of the X-axis relates to the measured Cycles:

- Start Cycle: First cycle to be displayed.
- End Cycle: Last cycle to be displayed.

X-Grid

Display the vertical grid for the X-axis.

Axis Type

The X-axis label can be chosen as follows:

- Cycle: Label the X-axis with the cycle number.
- Time Relative: Label the X-axis with the relative time since the measurement started.
- Time Absolute: Label the X-axis with the absolute time (time of day).
- Data: The measured values of one component are used as the X-Axis:

Data Axis	
Datablock	Ion Current
Dataname	14.00
Lowest Value	9.2302E-11
Highest Value	6.3072E-08
Scale	<input checked="" type="radio"/> Lin <input type="radio"/> Log

Fig. 4-78

- Datablock: Data block for the X-Axis
- Dataname: Component of the data block.
- Lowest Value: Lowest X-value to be displayed.
- Highest Value: Highest X-value to be displayed.
- Scale: Linear (Lin) or logarithmic (Log) X-axis scale.

4.3.1.5 Function submenu



Fig. 4-79

You can select the type of cursor function via the Function menu or call the graphic editor:

Magnify

Magnify a screen section.

Cursor

Display crosshairs to determine measure values.

Graphic Editor

Call graphic editor.

NOTE:

Switching from Magnify to Cursor and vice versa is also possible with a click of the right mouse button in the measurement data window.

4.3.1.5.1 Magnifying a Screen Section

With Magnify you can define an area of the measurement data screen with the left mouse button.




Fig. 4-80

To magnify the marked area, hold down the left mouse button, then click the right mouse button:




Fig. 4-81

- The magnifier function can be nested several times.
- To reverse the last magnification, choose [Special] > [Undo Magnify].
- The graphic can be restored to its original size by [Special] > [Set Default]:

4.3.1.5.2 Determine measurement data graphically (Cursor)

If [Function] > [Cursor] is selected, a measured value within the graphic display can be marked by the crosshairs. Its values for X and Y are displayed in the status bar.




Fig. 4-82

Hold down the left mouse button and position the crosshairs on the spot in the measure data graphic, of which you want to know the exact values. They are displayed in the status bar then. In mode Bargraph, a horizontal line appears instead of the crosshairs, which serves as cursor. The cycle selected by the cursor can be displayed as a table by holding the left mouse button down and simultaneously clicking the right mouse button:




Fig. 4-83

Measurement data block and cycle can be changed within this dialog box:

[Previous]

Display the measurement data of the previous block.

[Next]

Display the measurement data of the next block.

Cycle

Number of the selected cycle.

Cycle Time

Measurement time of the cycle (since the measurement started).

Date

Current date and absolute measurement time of the cycle.

[Change]

Manual elimination of freak values. See Section «[Change]», 302.

[Bargraph]

Display the bargraph of the current cycle. See Chapter 4.3.1.2.4 Bargraph, 292.

Displaying the threshold values**L/L, L/H, H/L, H/H**

If switching functions were activated during the measurement, an L (Low) or H (High) is shown behind the exponent of Value for documenting the switching functions A and B.

Operation with QMU 112

If the measure data have been recorded by a QMU 112, the QMU channel is displayed in square brackets after the mass number. 28 [7] means thus, mass number 28, measured on QMU channel 7.

Manual elimination of freak values**[Change]**

Measured values can be corrected manually:




Fig. 4-84

The measured value to be changed is selected with the cursor, modified in the Value field and stored by [Save].

Measurement errors (“freak values”) can also be corrected automatically by [Special] > [Filter].

4.3.1.5.3 Graphic editor

You can add text and drawings to your measurement pictures, for example, to prepare them for printout or as an aid for interpretation. See Chapter 4.2.1.5.3 Graphic editor, 267.

4.3.1.6 Special submenu



Fig. 4-85

Undo Magnify

Cancel the changes of the last magnify function.

Set Default

Restore all display settings to their basic values. All components and measurement data blocks are displayed across the stored range (all cycles).

Filter...

Automatic elimination of freak values.

4.3.1.6.1 Filter

Choose filter for automatic elimination of freak values:

Moving Median determines a moving median. For this purpose a certain number (selectable via Rank) of the last measured values is sorted by their magnitude and the one in the middle is used as the filter output:




Fig. 4-86

Moving Median

- Rank: Determines the number of measured values taken into consideration for the filter type Moving Median (number of measured values = Rank * 2 + 1).

Output File Name

- Name of the file in which the results are to be written. The name can be entered by activating the [File Manager] command button.

NOTE:

Write the filtered data on a new file so that the original data will be preserved in any case.

4.3.1.7 Info submenu



Fig. 4-87

File...

Display measure data documentation, if there is one. Otherwise a new documentation is created that consists of additional info, notes and title only. See Chapter 4.2.1.7 Info submenu, 270.

4.4 Auxiliary

Menu

Fig. 4-88

Mass Scale Calibration...

Display the measurement data (characteristic curve) of the mass scale calibration.

Zero Gas...

Display the measurement data of the background measurement.

Error History...

Display the stored error messages and warnings.

4.4.1 Mass scale calibration

Choose Mass Scale Calibration to display the calibration curve of the current analyzer:




Fig. 4-89

The exact peak position is determined for each defined mass. In the graphic, the deviation of the peak maximum from the nominal mass is shown. Deviations of up to ± 0.5 amu can be corrected. For larger deviations, a coarse alignment of the HF unit must be performed first.

NOTE:

More information about the mass scale calibration can be found in Chapter 2.8.1 Mass scale calibration (Mass Scale), 134.

Mass Scale Calibration menu



Fig. 4-90

File

Store/close the basic point file of the calibration curve.

Edit

Edit the mass numbers.

Info

Retrieve supplementary information about the file.

4.4.1.1

File submenu



Fig. 4-91

Save

Store the current calibration curve.

Close

Return to the Display Saved Values main menu.

4.4.1.2 Edit submenu



Fig. 4-92

4.4.1.2.1 Masses

With Edit Masses you can define mass numbers through manual input:




Fig. 4-93

NOTE:

This dialog box can also be called up by double-clicking the left mouse button in the measurement data window.

Defined Masses

- Number of Masses: Number of stored mass numbers.
- First Nominal Mass: Smallest stored nominal mass number.
- Last Nominal Mass: Largest stored nominal mass number.

Valid Range

Validity range of the correction curve. Prevents errors in case of an unfavorably defined mass range. Mass numbers outside this range are not extrapolated.

- Start Mass: Lower validity limit of the calibration curve.
- End Mass: Upper validity limit of the calibration curve.

In the above window (calibration gas = air) these two limits are calculated as follows:

- Start Mass = $14 - (44 - 14)/2 = -1 \Rightarrow 0$
- End Mass = $44 + (44 - 14)/2 = 59$

X-Grid

Display the vertical grid for the X-axis.

Y-Grid

Display the horizontal grid for the Y-axis.

[Add/Replace]

A mass number that is marked in the List of Masses box is transferred to the edit box by Add/Replace where it can be modified. A new mass number can also be entered in this field. The changed mass number is adopted, when Add/Replace is activated.

You can enter the mass number directly or change it with the or buttons. A click of the left mouse button changes the fraction of the mass by ± 0.01 , and if the right mouse button is pressed simultaneously, the step width of the change increases to ± 1.00 .

[Remove]

Remove the selected mass number.

**Inserting a
nominal mass
with the cursor**

Press the left mouse button in the calibration curve picture. A line cursor appears and the selected mass number is displayed numerically in the status bar.




Fig. 4-94

Please note that only integer mass numbers can be selected on the horizontal axis. The fraction (max: ± 0.49 amu) must be changed with the cursor position in the vertical direction.

When the right mouse button is pressed simultaneously the selected mass number is accepted and the Edit Masses dialog box is displayed:




Fig. 4-95

The inserted mass number can now be edited. Afterwards, press Add/Replace again.

NOTE:

The change is stored when [File] > [Save] is chosen.




Fig. 4-96

If you press [OK], the calibration curve is corrected:

4.4.1.3 Info submenu



Fig. 4-97

File...

Retrieve supplementary information about the file.

4.4.1.3.1 File




Fig. 4-98

Filename

Name of the file.

Software Version

Quadstar 32-bit version under which this file has been created.

Date

File creation date.

Time

Time of day the file has been created.

Author

Operator name.

4.4.2 Background measurement (Zero Gas)

[Auxiliary] > [Zero Gas] opens the file containing the measurement data of the background measurement, that belongs to the current analyzer. The intensities and mass numbers are displayed as a bar chart.




Fig. 4-99

The mass numbers are displayed below the X-axis and at the right edge.

NOTE:

More information about the background measurement can be found in Chapter 2.8.2 Background measurement (Zero Gas), 137.

Zero Gas menu

Fig. 4-100

File

Store/close the file.

Edit

Edit the mass numbers and intensities.

Info

Retrieve supplementary information about the file.

4.4.2.1**File submenu**

Fig. 4-101

Save

Store the mass numbers and their intensities.

Close

Return to the main menu of Display Saved Values.

4.4.2.2**Edit submenu**

Fig. 4-102

4.4.2.2.1**Values**

Values shows the intensities and the corresponding mass numbers of the chosen signal source:




Fig. 4-103

NOTE:

This dialog box can also be called up by double-clicking the left mouse button in the measurement data window.

Last modification

Date/time of the last modification.

Detector

Select the signal source. (For each selectable collector (Faraday .. Extern) a separate set of peaks is stored).

Defined Masses

Number and limits of the mass numbers.

- Number of Masses: Number of stored mass numbers.
- First Nominal Mass: Mass number of the first peak.
- Last Nominal Mass: Mass number of the last peak.

Defined Values

Limits of the peak intensities.

- Lowest: Lowest intensity of the stored peaks.
- Highest: Highest intensity of the stored peaks.

X-Grid

Display the vertical grid for the X-axis.

Y-Grid

Display the horizontal grid for the Y-axis.

Values

Edit values.

- Mass: Mass number.
- Value: Intensity.

The mass number and its intensity displayed in the Values dialog box under Mass and Value can be changed or deleted. Additional (new) peaks can also be entered:

[Add/Replace]

Click on the desired peak to edit it in the Mass and Value input fields. In these fields new values can also be entered numerically. To transfer the new or changed values to the list, press Add/Remove.

[Remove]

Remove removes the marked peak from the list.

Inserting a new peak with the cursor

Click on the spectrum with the left mouse button. Crosshairs are displayed in the measure picture and the cursor position is displayed numerically in the status bar:




Fig. 4-104

Now press the right mouse button additionally, so the selected mass is accepted and the Edit Masses dialog field is displayed:




Fig. 4-105

Press Add/Replace to include the new peak in the list of peaks:




Fig. 4-106

You can now edit the inserted peak. Press Add/Replace when you have finished.

NOTE:

The peak is not stored until [File] > [Save] is chosen.

Press [OK] to correct the displayed spectrum:




Fig. 4-107

4.4.2.3 Info submenu



Fig. 4-108

File

Retrieve supplementary information about the file.

4.4.2.3.1 File




Fig. 4-109

Filename

Name of the file.

Software

Version Quadstar 32-bit version used for creating this file.

Date

File creation date.

Time

Time of day the file has been stored.

Author
Operator name.

4.4.3 Error History (Stored Error Messages)

Display the stored error messages and warnings. This is only possible if the Error History control field has been activated in the Parameter Setup program under [Setup] > [General]. If the same error occurs several times consecutively (in the window below: 'MS 1': RF error), the error is not shown each time, but from the second time onwards a number is appended to the message. This number indicates the repetitions of the error (in the example: 'MS 1': RF error(15) there were 15 repetitions of the error). This way, the first and last occurrence of a repeated error is documented.




Fig. 4-110

Print...

Print the listed error messages. First you are asked for the message numbers to be printed. The printer can be configured in the Parameter Setup program under [Config] > [Printer].

Reset

Delete the error history file.

Close

Close the error history window.

5. Tune Up

Quadstar 32-bit Tune Up is used to set the QMS 422/421 or QMS 200 into an operational state. For example, the ion source can be tuned or the peak shape may be optimized graphically. The parameters available depend on the basic setting. See Chapter 3.8.1 QMS Configuration, [p. 234](#).

When the Tune Up program is started, Quadstar 32-bit attempts to establish a connection to the QMS. When this connection has been successfully established, the LED symbolically displayed in the lower left corner changes to green. If the QMS is not connected or not switched on, an error message is output. After the fault has been remedied, the connection setup can be reinitiated by [Comm] > [Connect].

Main menu



Fig. 5-1

File

Display version number of the Quadstar 32-bit Tune Up program and terminate it.

Tune

Tune the ion source, optimize peak shape, print ion source settings. In addition (QMS200): adjust mass scale, show operating conditions. See Chapter 5.2 Tune, [p. 319](#).

Manual

Operation/display of the digital input/output channels (DI/DO). See Chapter 5.3 Manual parameters, [p. 341](#).

Comm

Switch the interface to the QMS on/off, select active QMS. See Chapter 5.4 Communication Parameters, [p. 343](#).

Setup

Switch the emission and the SEM on/off, degas the filament. See Chapter 5.5 Setup Parameters, [p. 344](#).

Help

Display the 'Key Fragment Ions' table and help contents. See Chapter 1.9 Hints concerning the work with Quadstar 32-bit, [p. 56](#).

5.1 File

File Menu



Fig. 5-2

5.1.1 About Tune Up

In this menu entry the version number of your Quadstar 32-bit Tune Up program is shown:



Fig. 5-3

5.1.2 Exit

Terminate the Quadstar 32-bit Tune Up program.

5.2 Tune

Tune menu




Fig. 5-4

Ion Source...

Parameter setup for the different ion source types:

- Ion source parameter setup when a QMA 200 analyzer is used. See Chapter 5.2.1.1.1 QMA 200 (Prisma™), 320.
- Ion source parameter setup when a QMA 4XX analyzer is used. See Chapter 5.2.1.1.2 QMA 4XX, 320.
- Optimize ion source settings when a QMA 125 analyzer is used. See Chapter 5.2.1.3 QMA 125, 323.
- Print ion source settings.

Peak Shape...

Optimize peak shape graphically. See Chapter 5.2.2 Peak Shape Optimizer, 333.

QMS 200 Tune Data...

Graphic and numeric display of QMS 200 operating conditions. See Chapter 5.2.3 QMS 200 Tune Data, 337.

QMS 200 Tune Mass Scale...

Calibration of the QMS 200 mass scale. See Chapter 5.2.4 QMS 200 Tune Mass Scale, 338.

QMS Offset...

Offset calibration of the measure amplifier. See Chapter 2.8.3.1 QMS Offset Calibration Setup, 141.

5.2.1 Ion Source

Depending on the selected hardware configuration Tune Ion Source supports:

- Parameter setup of the Ion Source
- Optimizing the ion source settings by tuning it while a measurement is performed
- Display of the emission current (QMA 125)
- Printing the ion source parameters

The type of dialog box that is presented depends on the installed analyzer and RF generator. If the parameters or the appearance of the displayed dialog box do not agree with the installed hardware, you must first define the correct QMS configuration by calling up the [Config] > [QMS] menu of the Parameter Setup program.

5.2.1.1 Parameter setup for the QMA 200/QMA 4XX

In this dialog box you can enter the ion source parameters Ion Source, IS-Emission and IS-Voltages. To optimize these parameters choose the [Display] > [Measure Data] menu where you will find a test measurement in Sample, Scan Bargraph and Scan Analog mode.

5.2.1.1.1 QMA 200 (Prisma™)




Fig. 5-5

5.2.1.1.2 QMA 4XX




Fig. 5-6

5.2.1.2 Modifying parameters

All ion source and channel parameters can be modified while the measurement is in progress. The parameters can be entered via the input field located below the Tune Ion Source menu bar. This field shows the parameter that has been selected by the cursor.

In the dialog boxes above, the cursor is positioned in the emission current input field:



Fig. 5-7

To enter a changed value

- Press the <ENTER> key
- Click on the small command box on the left of the input field or
- Press the Tab key.

5.2.1.2.1 Channel parameters

- 1 First choose a channel and switch it on (Enable):

Channel

Channel number (QMS channel 0...63) and state (Enable/Skip/Off)

- 2 Then enter the required parameters in the corresponding groups:

Detector

Detector parameter group.

Mass

Mass parameter group.

Amplifier

Amplifier parameter group.

NOTE:

For More information on the channel parameters see Chapter 3.2.4.1 Measuring parameters overview, 159 or in the User's Guide of the corresponding QMS.

Ion source parameters (Ion Source)

QMA 200:

The ion source type can not be chosen for the QMA 200, but it must have been set correctly under [Config] > [QMS]. It is displayed here in the field IS Type.

Filament

Select the desired filament (Fil 1, Fil 2 or Fil 1+2). Fil 1+2 means here, that both filaments are switched on simultaneously.

IS-Set #

Load an IS set (Set0 ... Set3). There are four IS sets with ion source voltages that contain one record each for all filament configurations (Fil 1, Fil 2 and FIL 1+2).

QMA 4XX:**Type**

Selection of special ion source types that are feasible according to the configuration (Parameter Setup program, [Config] > [QMS...] > [IS-Type]).

Filament #

Select the desired filament (Fil 1, Fil 2 or Fil 1+2). The number of available filaments depends on the ion source type.

Set for Fil 1

Assign an IS set (0 ... 3) to Filament 1.

Set for Fil 2

Assign an IS set (0 ... 3) to Filament 2.

5.2.1.2.2 Filament**IS-Emission**

These parameters relate to the filament selected under Ion Source.

- Current: Enter the desired emission current (0.00 mA corresponds to OFF).
- Protect: Enter the maximum overload current (Protection).

5.2.1.2.3 RF supply (QMS 200)**RF Polarity**

QMA 200: Allows the changeover of the RF polarity on the mass filter rods.

5.2.1.2.4 Ion source voltages**IS-Voltages**

Four independent IS sets are available for each ion source type:

- QMA 4XX: The IS-Voltages relate to the IS set that has been chosen under Set for Fil x.
- QMA 200: The IS Voltages relate to the IS set selected under IS-Set #. In this set a separate ion source voltage group (IS-Voltages) exists for each filament configuration (Fil 1, Fil 2 and Fil 1+2).

NOTE:

For detailed information on these parameters and their optimization please consult the User's guide of the QMS, or the installed ion source.

5.2.1.2.5 Filament Protection Extern (QMA 200/QMA 4XX)

QMA 200

On the QMA 200 (Prisma™) the SEM and filament can automatically be switched on when evacuated and off when flooded, if a total pressure sensor (TPR 250/PKR 250/251) is connected.

Fil. Prot. Thresh.

On/Off thresholds for the measured total pressure.

- ON below: Lower threshold to be attained before the filament and the SEM voltage are switched on in automatic operation [Protection (SEM+FIL)] > [EXTERN].
- OFF above: Upper threshold at which the filament and the SEM voltage are switched off in automatic operation [Protection(SEM+FIL)] > [EXTERN] and in [EXT-PROT] mode.

In operation under [Protection(SEM+FIL)] > [EXT-PROT] the filament and the SEM voltage must be switched on manually.

The function (INTERN/EXTERN/EXT-PROT) is chosen in the Parameter Setup program under [Config] > [QMS...] in the Protection(SEM+FIL) field.

QMA 4XX

When a QMA 4xx is used, no thresholds can be entered. The protection mechanisms in [SEM+FIL] > [EXT-PROT] and [SEM+FIL] > [EXTERN] mode relate to external switching contacts.

5.2.1.3 QMA 125

When operating with a QMA 125 (without QMU 112) a test measurement may be started in the same way as with QMA4XX/200 ([Display] > [Measure Data]). The actual ion source parameters have to be adjusted at the QME itself. The parameters marked by a '*' (*IS Type and *Emiss) are only shown in the editor and cannot be modified via Quadstar 32-bit:

Operation without QMU 112




Fig. 5-8

Operation with QMU 112




Fig. 5-9

Emission

- *Emiss/Current: Emission current.
- Filament #: Filament (Fil1 or Fil2). The number of available filaments depends on the type of ion source.

QMU Channel

Choose the QMU 125 channel. Only the QME 125 units activated in the Parameter Setup program ([Config] > [QMS]) can be selected.

NOTE:

If Filament Select Remote is selected on the individual QME 125 units in configurations with a QMU 112, all filaments are switched jointly.

*IS Type/Ion Source Type

Ion source type in use.

Tune Ion Source menu



Fig. 5-10

File

Save or print the ion source parameters.

Display

Display the measurement data. See Chapter 5.2.1.3.5 Optimizing the parameters (Measure Data), 329.

Option

Edit/copy channel parameters and ion source parameters.

File submenu



Fig. 5-11

Save

Save the ion source parameters.

Print Ion Source Settings

Print the ion source settings.

Close

Close the parameter file and return to the Tune Up main menu.

Display submenu

Fig. 5-12

Measure Data

Perform a test measurement. See Chapter 5.2.1.3.5 Optimizing the parameters (Measure Data), 329.

Option submenu

Fig. 5-13

Edit Channel Par...

Edit the channel parameters. See Chapter 5.2.1.3.1 Standard editor (Edit Channel Par), 326.

Copy/Exch. IS-Set...

Copy/exchange IS parameter sets. See Chapter 5.2.1.3.2 Copy/Exch. IS-Set, 327.

Copy/Exch. Channel...

Copy/exchange channel parameters. See Chapter 5.2.1.3.3 Copy/Exch. Channel, 328.

Default IS-Set...

Set the current ion source parameters to default values. See Chapter 5.2.1.3.4 Default Ion Source Set, 329.

5.2.1.3.1 Standard editor (Edit Channel Par)

Set up the QMS channel parameters with the standard editor of Quadstar 32-bit:

Load-Ch.:00	CH-0	CH-1	CH-2	CH-3	CH-4	CH-5	
State	ENABLE	ENABLE	ENABLE	OFF	OFF	OFF	OFF
Det. Type	FARADAY	FARADAY	FARADAY	----	----	----	----
Mass Mode	SCAN-N	SCAN-F	SAMPLE	----	----	----	----
First Mass	0.00	0.00	----	----	----	----	----
Mass	----	----	0.00	----	----	----	----
Speed	0.2s	1s	----	----	----	----	----
Dwell	----	----	1s	----	----	----	----
Width	30	20	----	----	----	----	----
Resolution	10	25	50	----	----	----	----
Threshold	----	----	----	----	----	----	----
Steps	64	64	----	----	----	----	----

Fig. 5-14

The Ion Source Channel Parameter editor basically corresponds to the standard editor described in Chapter 3.2 Measure Parameters, 156. It fulfills the same functions as the input fields described in Chapter 5.2.1.1 Parameter setup for the QMA 200/QMA 4XX, 320 and Chapter 5.2.1.3 QMA 125, 323, but in addition to that, other parameters may be changed as well (e.g. Display which controls the display range). Because several channels can be displayed simultaneously, it provides greater clarity (e.g. for copying functions).

NOTE:

More information is found in Chapter 3.2 Measure Parameters, 156.

Selecting the measuring mode

Click on the Mass Mode field with the right mouse button of the desired channel, so the context menu for selecting the desired measuring mode appears:




Fig. 5-15

5.2.1.3.2 Copy/Exch. IS-Set

Copy or exchange complete ion source parameter sets (IS-Sets):




Fig. 5-16

Mode Selection

Desired action:

- Copy: Copy an IS set to another IS set or to all other IS sets of the current ion source.
- Exchange: Exchange all parameters of the Source IS set with those of the Dest. IS set of the current ion source.

Source

IS set, whose parameters are to be copied or exchanged.

Dest.

Target IS set to copy the parameters to or exchange with. ALL (for copying to all IS sets) can be specified by clicking on the [UP] or [DOWN] arrow with the right mouse button.

5.2.1.3.3 Copy/Exch. Channel

Copy or exchange individual parameters or all parameters of the selected channels:




Fig. 5-17

Mode Selection

Desired action:

- Copy: Copy the object chosen under Object Selection to a different channel or to all channels.
- Exchange: Exchange all parameters of the Source channel with those of the Dest. channel. Exchange is only possible if [Object Selection] > [Channel] is active.

Object Selection

Choose the object to be copied or exchanged:

- Channel: All parameters of the channel are copied or exchanged.
- Parameter: Only the current parameter (marked before Copy/Exch. Channel was called) is copied or exchanged. Parameter cannot be chosen if no parameter is selected.

Source

Channel whose parameters are to be copied or exchanged.

Dest

Target channel to copy the parameters to or exchange with. ALL (for copying to all channels) can be specified by clicking on the [UP] or [DOWN] arrow with the right mouse button.

Special cases

The following special cases must be taken into consideration when copying parameters:

Range

QMS 421: Copies only into those channels that have the same Lock-In status (only meaningful, if the optional Chopper Lock In Amplifier CLA 421 is in use).

Threshold

Copies only into those destination channels for which the same Detector has been selected under Det. Type.

5.2.1.3.4 Default Ion Source Set

Default IS set loads default values for the ion source parameters of the current IS-set.




Fig. 5-18

Default Selection

- Factory: Load factory default values for the current filament.
- User Fil #X: Load or save a self defined Default ion source set (User Fil #X) by Load or Save.

Mode Selection

- Load: Load factory default or self defined IS-set.
- Save: Store the current IS set as User Default (for that, [Default Selection] > [User Fil #X] must be selected).

5.2.1.3.5 Optimizing the parameters (Measure Data)

[Display] > [Measure Data] starts a measurement on the selected channel, provided the communication with the QMS is switched on.

Optimize the ion source by watching the influence of parameter changes on the measurement. To check whether there is an offset, the measurement scale reaches from 100% positive to 10% negative (of full scale range). If there is an offset, QMS 422/200 can be calibrated by [Tune] > [QMS Offset], the QMS 421 must be adjusted by the channel parameter Offset.

Scan Analog measurement data




Fig. 5-19

The measured values of the different cycles are drawn in different colors in the same picture. The display field is only cleared if:

- the size of the Measure Data window is changed,
- the First, Width or Range parameters are changed, or
- if this is explicitly chosen by the user through [Refresh] > [Display].

NOTE:

For printing pictures in Quadstar 32-bit (by Ctrl+Shift+F12), stop the measurement by [Operation] > [Halt] in order to prevent errors of the interface.

Sample Measurement data




Fig. 5-20

In this measurement, the largest measured intensity is marked with a drag pointer. The status of the measurement (Go or Halt) is shown behind the mass number.

NOTE:

For printing this picture you should select [Operation] > [Halt] in order to prevent error messages (Time Out, etc.) from the interface.

Measure Data menu



Fig. 5-21

File

Store the reference data (only Scan Analog)

Operation

Start/stop the measurement. Its status is shown in the measurement window.

Setup

Scale the Y-axis, change alternating colors.

Display

Refresh the display, show the reference on/off

File submenu



Fig. 5-22

Save Reference

Store the measurement data of the next scan as a reference (only in Scan Analog mode SCAN-N und SCAN-F).

Close

Close the measurement window.

Operation submenu



Fig. 5-23

Go

Start the measurement.

Halt

Stop the measurement.

Setup submenu



Fig. 5-24

Axis...

Choose between linear and logarithmic Y-axis scale. This box can also be accessed by double-clicking the left mouse button.




Fig. 5-25

- Y-Scale
 - Lin: Display the measurement data on a linear Y-axis.
 - Log: Display the measurement data on a logarithmic Y-axis.

Colors...

Dialog box to choose the alternating colors for a repeated scan analog.




Fig. 5-26

In the fields #0 ... #63, enter the colors to draw the cycles with (alternatingly). The standard color table, that contains 8 colors, produces scans with good contrast that are easily distinguishable in most cases. It can be recalled by [Set Default Colors], but the current settings of the Alternating Colors for the current color palette will be lost then. See Chapter 3.8.5 Screen Colors, 245.

Display submenu

This menu is only accessible in Scan Analog measuring mode (SCAN-N and SCAN-F):



Fig. 5-27

Refresh

Clear the display.

Reference

Display a reference scan.

Displaying reference data (Reference)

In SCAN-N and SCAN-F measuring mode you can also display the measurement data of a reference scan. When the Reference option (Display menu) is active, a reference scan is automatically stored

- When the first scan is performed, or
- If either the First, Width or Range parameter changes.

You may store a reference scan manually as well by [File] > [Save Reference]. This may be necessary for example, to make the option Reference accessible again, if it has been grayed.

5.2.2 Peak Shape Optimizer

[Tune] > [Peak Shape] allows:

- to shape peaks by varying the ion source voltages.
- to compare the shapes of two peaks close to their maximums.
- to determine mass number and relative intensity precisely.

The displayed parameters depend on the hardware configuration. If the parameters or the appearance of the displayed dialog box do not agree with the installed hardware, you must first define the correct QMS configuration by calling the [Config] > [QMS] menu of the Parameter Setup program.




Fig. 5-28

5.2.2.1 Channel parameters

In group Channels you can set up the two measure channels you want to use. These are 'virtual' channels, they do not interfere with the ones you normally use on the QMS. You may set all channels to ENABLE, but only the ones shown as Peak A and Peak B are really measured. If you often switch between different measurements, just enter the corresponding channel numbers for Peak A and Peak B to do so. If you wish to measure Peak A only (here: channel 0), just set Peak B to State 'SKIP'.

5.2.2.2 Ion source voltages

Under IS-Voltages you find those voltages of the used ion source, that can be set up by software. They are channel independent and provide peak shape conditioning and optimization.

NOTE:

The ion source parameters entered here do apply to the Peak Shape Optimizer

only. As soon as you quit it, the settings shown under [Tune] > [Ion Source] will take effect again.

5.2.2.3 Display width

Under [Display] > [Width] you can enter the mass range to display. This value is a minimum that is rounded up to the next multiple of the current stepwidth (QMS 422/421: inverse value of the channel parameter Steps; QMS 200: 1/32).

5.2.2.4 Mean value

[Mean Value] > [Scans] allows you to choose the number of scans to be performed per channel. Then the mean value of these scans is built in order to achieve higher resolution and precision.

Peak Shape Optimizer menu



Fig. 5-29

File

Save measurement parameters, quit Peak Shape Optimizer.

Display

Perform and display measurement (See Chapter 5.2.2.5 Start measurement (Peak Shape), 335).

Option

Call standard editor, set measurement parameters to default values.

File submenu



Fig. 5-30

Save

Save all parameters (channel, ion source, display and mean value parameters).

Close

Quit the editor, return to the Tune Up main menu.

Display submenu



Fig. 5-31

Peak Shape...

Perform and display measurement using the entered parameters. See Chapter 5.2.2.5 Start measurement (Peak Shape), 335.

Option submenu



Fig. 5-32

Edit Channel Par...

Call the standard editor to edit channel parameters that are not present in the Peak Shape Optimizer Editor. You should use this option only if you are sure that you really want to change those values that are set to default and are not displayed usually (e.g. Steps).

Set Default

Set all parameters to default values.

5.2.2.5 Start measurement (Peak Shape)

Start measurement by [Display] > [Peak Shape]. The QMS first measures Peak A (from Mass -0.5 to Mass+0.5) as many times as entered under [Mean Value] > [Scans]. After that, the peak maximum is scaled and displayed across the mass range specified by [Display] > [Width]. Then Peak B is measured the same way. The status (Op: Go/Halt/Done) and progress (Progress: XX %) of the measurement are displayed at the bottom of the measure window.

The x axis is lettered with mass numbers, the y axis with relative intensities in %.




Fig. 5-33

Should the shown curves seem unsensible or 'untrustworthy', please check the used parameters (esp. Amp.Range and Mass) on a scan analog measurement. Particularly a slight Overrange, although indicated during the measurement instead of Progress, might not obviously show up in the finished measurement.

Peak Shape menu



Fig. 5-34

File

Close measurement window.

Operation

Start/stop measurement. The status of the measurement is shown at the bottom of the measurement window.

Setup

Setup display.

File submenu



Fig. 5-35

Close

Close measurement window.

Operation submenu



Fig. 5-36

Go

Start measurement.

Halt

Stop measurement.

Setup submenu



Fig. 5-37

Axis...

Setup display (lin/log, grid).




Fig. 5-38

- Y-Scale: linear (Lin) or logarithmic (Log) display of the y axis.

- X-Grid: show or hide the mass scale grid.
- Y-Grid: show or hide the intensity scale grid.

5.2.3 QMS 200 Tune Data

QMS 200 Tune Data displays certain operating parameters of the QMS 200 Prisma™ and provides tuning of the RF 20X RF generator:




Fig. 5-39

The individual fields have the following meaning:

Emission Current

Actual value of the emission current.

Filament Current

Actual value of the filament current.

Tune Voltage

Tuning voltage of the RF generator. The RF generator is tuned at the QMS 200 itself.

Sem Voltage

Actual value of the SEM voltage.

Temperature

Temperature inside the QMS 200. Normal: 30 ... 50 °C/maximum: <80 °C. If the temperature rises too high: change the filter mat of the QMS 200 (See User's Guide of the QMS 200), check the fan if necessary.

Reset Arrows (Display menu)

The drag pointers for minimum and maximum values can be reset to the current value by [Display] > [Reset Arrows].

Ion Source (Parameters menu)

Ion Source opens the same dialog box as [Tune] > [Ion Source] in the main menu does.

Setup (Parameters menu)

[Parameters] > [Setup] (or double-clicking the mouse within the window) opens the dialog box shown to choose the colors for the individual display bars:




Fig. 5-40

Click on the desired color above the horizontal scroll bar or change the color number in the corresponding edit field.

5.2.4 QMS 200 Tune Mass Scale

The mass scale of the QMS 200 can be calibrated with one peak in the upper and one in the lower region of the mass scale. If no special calibration gas is available, the peaks of air can be used (for example, for measurements within the mass range 0 ... 50):




Fig. 5-41

A measurement can be stopped by [Operation] > [Halt] and restarted by [Operation] > [Go].

Calibrating the mass scale

- Manually

- With the keys [<<], [<], [>] and [>>] the peak tops can be shifted slowly or fast in either direction of the mass scale within the given limits (Offset: -767 ... +767, Slope: -394 ... +394).



Fig. 5-42

- The values can also be entered directly in the fields Offset and Slope and the result can be seen on the display.
- Using the mouse
 - 1 Press the left mouse button and hold it.
 - 2 Position the cursor on the peak top.
 - 3 Press the right mouse button and hold it.
 - 4 Drag the cursor to the correct position while still holding both mouse buttons pressed.
 - 5 Release the right mouse button.
 - 6 Release the left mouse button.

Parameter setup

In the Parameters menu you can change the setup of the display and the measurement channels:



Fig. 5-43

Setup

[Parameters] > [Setup] opens a dialog box to choose between linear and logarithmic Y-axis scale. This dialog box can also be activated by a double-click of the left mouse button.




Fig. 5-44

Y-Scale

- Lin: Display the measurement data on a linear Y-axis.
- Log: Display the measurement data on a logarithmic Y-axis.

Channel

[Parameters] > [Channel] enables you to change certain channel parameters that are relevant for the scans. This dialog box can also be called by double-clicking the right mouse button:

QMS 200 Tune Mass Scale Parameter

Detector Parameters

Type: CH-TRON

SEM voltage: 1400

Mass Parameters

Mass Mode: SCAN-F

First Mass: 0.00 41.00

Width: 6 6

Speed: 1 s

Resolution: 50

Display Parameters

Disp. Range: E-07 E-09

Disp. Decades: 3 2

OK Cancel

Fig. 5-45

5.3 Manual parameters

5.3.1 Digital I/O

By [Manual] you can display the state of digital input channels and switch on or off digital output channels manually, provided that the communication to the mass spectrometer is switched on.



Fig. 5-46

For turning the individual digital output channels on or off, double-click the corresponding push button with the left mouse button. The switching state is indicated by the colored dot on the push button: green = on, grey = off.

DI/DO Manual for QMS 422/421




Fig. 5-47

Display left

Select the DI/DO group to be displayed on the left side of the screen.

Display right

Select the DI/DO group to be displayed on the right side of the screen.

A group cannot be displayed at the left and the right side simultaneously.

DI/DO Manual for QMS 200

Only the two standard digital input and output channels are available:




Fig. 5-48

The digital input and output channels of the QMS 200 have the following designation in the hardware description:

DI 00 : EXT_PROT (Connection TPR/PKR, Pin No. 1 <=> 2)

DI 01 : EXT_RUN (Connection Control, Pin No. 10 <=> 15)

DO 00 : REL_1 (Connection Control, Pin No. 11 <=> 12)

DO 01 : REL_2 (Connection Control, Pin No. 13 <=> 14)

5.4 Communication Parameters



Fig. 5-49

Connect

Establish connection to the QMS.

Disconnect

Cancel the connection to the QMS.

SetUnit...

Select current measuring device. See Section «Setting the node addresses of several control units», 218.

NOTE:

For detailed information see Chapter 3.6 Communication parameters, 209.

5.5 Setup Parameters

Setup menu



Fig. 5-50

SEM/Emission Control...

Open the dialog box SEM/Emission Control to:

- Switch the SEM high voltage on/off.
- Input the Common SEM voltage.
- Switch the Emission on/off.

See Chapter 3.7 Setup parameters, 223.

Degas Control...

Degas the filament

Operation with QMA 200/QMA 4XX, see Chapter 5.5.1.1 Operation with QMA 200/QMA 4XX, 344, Operation with QMA 125, see Chapter 5.5.1.2 Operation with QMA 125, 345.

5.5.1 Degas Control

If the communication to the QMS is switched on, you can degas the current filament of the ion source (Fil 1 or Fil 2).

5.5.1.1 Operation with QMA 200/QMA 4XX

For the QMA series 400 and the QMA 200 (Prisma™), the Filament #, degas time, emission current and the maximum filament current can be programmed:

QMA 200




Fig. 5-51

QMA 4XX




Fig. 5-52

[Start]

Start degassing.

[Stop]

Stop degassing.

Degas

- Filament #: Filament to be degassed.
- Current: Emission current for degassing.
- Protect: Maximum filament current for degassing.
- Time: Time in minutes during which the filament is to be degassed; when this time is over, the degassing stops automatically. If you prefer to start and stop manually rather than state a time, click on [UP] or [DOWN] with the right mouse button; now in the Time field, man. (manual) is displayed.

Degasing

Elapsed time since the start, during which the ion source has already been degassed.

5.5.1.2 Operation with QMA 125

When a QMA 125 is used, only the filament can be selected:

Operation without QMU 112




Fig. 5-53

Operation with QMU 112




Fig. 5-54

[Start]

Start degassing.

[Stop]

Stop degassing.

QMU Channel

QMU 125 channel, whose ion source is to be degassed.

Degas

- Filament #: Filament to be degassed.
- Time: Time in minutes during which the filament is to be degassed; when this time is over, the degassing stops automatically. If you prefer to start and stop manually rather than state a time, click on [UP] or [DOWN] with the rightleft mouse button; now in the Time field, man. (manual) is displayed.

Current

Emission current during degassing (display only; Emission and Protection must be set up at the spectrometer itself).

Degasing

Elapsed time since the start, during which the ion source has already been degassed.

5.6 PPM Tune Up

Quadstar 32-bit PPM Tune Up is used for entering and adjusting/optimizing various parameters of the energy filter or the mass filter (Ion mode), and the ion source (Neutral mode) of the PPM 422 Plasma Process Monitor. You can manually perform measurements in Tune Up as well.

When the Tune Up program is started, Quadstar 32-bit PPM attempts to establish a connection to the PPM 422 Plasma Process Monitor. If this attempt is successful, the LED symbolically indicated in the lower left corner of the active window changes to Green. If a QMS 422 is not connected or not switched on, a corresponding error message is displayed.

After the fault has been remedied, a new attempt to connect can be made by [Comm] > [Connect].

5.7 PPM 422 Plasma Process Monitor

The PPM 422 Plasma Process Monitor is a differentially pumped mass spectrometer with an energy filter between the ion source and the mass filter. It is used in plasma processes for analyzing ions and neutral particles (molecules, radicals) which are extracted (for example, from a sputter plasma) by means of an entrance orifice. The analysis is not limited to the mass of the particles, it also includes their energy.

NOTE:

For a detailed description please refer to the user's guide "PPM 422 Plasma Process Monitor".

The potentials required for operating the various function groups of the PPM 422 are illustrated in Fig. 5-55, 350. Depending on the operating mode the potentials are connected differently. The detailed function of the individual potentials is described in the PPM 422 User's guide. The Quadstar 32-bit PPM Tune Up module described in this chapter allows access to all relevant system parameters. These are the potentials for operating the transfer lens/ion source, the energy filter, and the mass filter (including detector), as well as the filament current and emission current of the ion source. Tab. 5-1, 351 provides an overview of the parameters accessible in Quadstar 32-bit PPM Tune Up. These can be optimized for various operating modes while the measurement is running.

5.7.1 Functional diagram Ion-optical potentials




Fig. 5-55

For operating the PPM 422 Plasma Process Monitor, two QMS 422's are needed: "QMS-E" or "E-Chassis" for the energy analysis and "QMS-M" or "M-Chassis" for the mass analysis. Depending on the selected operating mode the ions are detected through one of the QMS chassis. A separate operating voltage for the extraction hood, if required, must be supplied by the user. (See PPM 422 user's guide).

The QMS-E generates the voltage V_e for the energy scan. Its IS 420, designated as "IS-E", essentially supplies the energy analyzer and the mass filter. In Neutrals mode, V_e defines the potential of the IS-E to ground. In Ions mode the IS-E and IS-M are both connected to the V_e potential.

5.7.2 Parameter overview

Name	Operation mode	Voltage	Default value	Function
CATH C	N	V2-M	70 V	El. energy Coarse
IONREF	N	V1-M	-	Ions' Reference
ENTRY	I'	V9-M	30 V	Entrance lens
ITRO	I'	V5-M	200 V	Ion lens
EXTR	N	V5-M	4 V	Extraction Coarse
EXTR	N	V5-M and div.	-	Extraction Fine
FOC 1	I' N	V8-M	60 V	Lens 1
CENTR	I' N'	V2-E	41 V	Inner cylinder
MIRR	I' N'	V8-E	17 V	Outer cylinder
INFL	I N	V9-E	26 V	Inflection/deflection

Name	Operation mode	Voltage	Default value	Function
FOC 2	I ' N '	V3-E	7 V	Lens 2
CASE	I ' N '	V5-E	28 V	QMA tube
F.AXIS	I ' N '	V4-E	4 V	Field axis
DFLI	I ' N '	V6-E	300 V	Inner deflection
DFLO	I ' N '	V7-E	145 V	Outer deflection
Current	N	-	1.0 mA	Emission current
Protect	N	-	3 A	Max. filament current
ORIF	I N	AO 422	-	Extr. hood
Ve/Energy	I ' N '	E-Scan	-	Ion energy

Tab. 5-1

Parameter Overview

Operation modes: N = Neutrals, I = Ions

The apostrophe in the operating mode designation means that this value can be changed for hardware adjustment (or: intentional shift with the SHIFT function) of the energy scale.

Main menu

The menu with its parameter groups deviates from the standard Quadstar 32-bit Tune Up software and is specially fitted to the PPM 422 Plasma Process Monitor:



Fig. 5-56

File

Display the version number of the Quadstar 32-bit PPM Tune Up program and exit the program. See Chapter 5.8.1.1 About Tune Up, 352.

Tune

Input of ion source parameters, perform measurement. See Chapter 5.9 Tune, 353.

Comm

Turn the interface to the PPM 422 on/off. See Chapter 5.10 Communication Parameters, 375.

Setup

Turn the emission and the SEM on/off. See Chapter 5.11 Setup Parameters, 376.

Help

Display the 'Key Fragment Ions' table and help contents. See Chapter 1.9.1 Spectrum interpretation, 56.

5.8 File

5.8.1 File menu



Fig. 5-57

5.8.1.1 About Tune Up

Under this menu entry you can identify the version number of the Quadstar 32-bit PPM Tune Up program.



Fig. 5-58

Exit

Quit the Quadstar 32-bit PPM Tune Up program.

5.9 Tune

Tune menu

Choose the Tune menu to adjust and optimize the transfer lens, the ion source and the energy- or mass filter. You may perform measurements with it as well.



Fig. 5-59

PPM 422...

Adjust and optimize the PPM 422 specific parameters, perform measurements manually.

NOTE:

[Tune] > [PPM 422...] is only accessible when the communication is switched ON.

5.9.1 PPM 422

Choose [Tune] > [PPM 422] for entering and adjusting/optimizing the parameters of the energy- or mass filter, the ion transfer optics (Ions mode) and the ion source (Neutrals mode), in order to prepare the PPM 422 Plasma Process Monitor for operation.

The following parameters can be optimized:

- Operating modes Ions+/Ions- Voltages of the mass/energy filter and of the ion transfer optics
- Operating mode Neutrals Voltages of the mass/energy filter and the ion source

5.9.1.1 Parameter setup

In the following dialog boxes you can enter the mass/energy parameter (M- or E- parameters) as well as the ion source parameters Ion Source, IS-Emission and IS-Voltages:

- For selecting the PPM 422 operating mode the groups Scan, Sampling Mode, Ions and Neutrals are available.
- For adjusting/optimizing these parameters, an on-line measurement in different scan modes is available under [Display] > [Measure Data] menu.

5.9.1.1.1 Parameter setup of an E-Scan




Fig. 5-60




Fig. 5-61

5.9.1.1.2 Parameter setup of an M-Scan




Fig. 5-62




Fig. 5-63

5.9.1.1.3 Modifying Parameters

The displayed parameters can be modified while a measurement is in progress, by clicking on the desired parameter. It then appears in the edit field displayed below the menu line of Tune PPM 422 and can be modified there.

The change takes effect, when you

- Press the Enter key,
- Click on the small [*] command button on the left of the input field, or
- Press the Tab key.

In the dialog box below, the cursor is positioned at the Mode parameter of the group Scan.




Fig. 5-64

The parameters are subdivided into groups:

Channel

Enter the channel number (QMS 422 channel 0 ... 63)

Detector

Parameter group Detector

Mass/Energy

Parameter group Mass/Energy

Amplifier

Parameter group Amplifier

NOTE:

For detailed information on the parameters please refer to the QMS 422 user's guide (Parameters, sorted alphabetically).

Please note that the parameter designations are modified by the prefixes "E-" and "M-" (for Energy/Mass).

IS-Voltages

IS Voltages relates to the IS sets that have been selected under E-Set and M-Set. Sixteen IS sets are available for each QMS 422. The voltage display corresponds to the selected sets of both QMS 422.

Special

Special parameters

- **SHIFT:** SHIFT shifts the data display (only in Energy scan mode) in such a way, that negative values can be displayed as well and therefore the distribution near zero can be seen (adjustment range 0 ... +10 Volt). The energy calibration remains unchanged. SHIFT has only a temporary effect and is only useful when a measurement is made.
- **CORR:** Software implemented correction of the residual deviation between the energy scale and the true energy. The possible hardware adjustments are not fine

enough for this purpose. CORR is determined separately for each E-Set and saved together with the set. (Adjustment range -2.0 ... +2.0 Volt).

NOTE:

For detailed information, please refer to the user's guide of the PPM 422 Plasma Process Monitor, Interrelation between ion energy and Ve.

- ORIF: Extraction hood potential. The extraction hood potential can be output via an analog output channel and the corresponding AO characteristic curve. It serves as a control voltage for an external power supply. ORIF has only an effect, if the corresponding AO characteristic curve has been defined and an external power supply is connected.

Ion Source

- Filament #: Select the desired filament (Fil1 or Fil2) for the Neutrals mode.
- E-Set: Load an IS-E- Set (0 ... 15) for the QMS E.
- M-Set: Load an IS-M-Set (0 ... 15) for the QMS M.

The QMS E and the QMS M each have sixteen separate IS-sets.

IS-Emission

These parameters relate to the filament selected under Ion Source (Fil1 or Fil2).

- Current: Emission current in mA (0.00 means OFF).
- Protect: Maximum filament current in A (Protection).

Scan

These parameters determine the operating mode of the measurement. They are stored together with the other channel parameters.

- Mass: Mass scan. The parameter Fixed E determines the fixed energy value in the QMS-E that is used for the mass scan.
- Energy: Energy scan. The parameter Fixed Mass determines the fixed mass value in the QMS M that is used for the energy scan.

NOTE:

In a changeover between MASS and ENERGY all parameters of a channel are preserved, only the designations change. A channel whose parameters were originally set for M-Scan would now be switched to E-Scan. A large part of the parameters would no longer be applicable. For example "Fixed Ve 5.4" would become "Fixed Mass 5.4". For changing the scan mode, it is therefore more convenient to select a different channel.

It is recommended to reserve the channels for the corresponding scan mode based on their numbers, for example odd numbers for M-scan. In the supplied software Channel 0 is set for E-Scan and Channel 1 for M-Scan.

Sampling Mode

Select the operating mode (plasma sampling mode PSM).

- Ions +: Measurement of positive ions.
- Ions -: Measurement of negative ions.
- Neutrals: Measurement of neutral species by means of electron impact ionization in the ion source.

NOTE:

In a change of the Plasma Sampling Mode (PSM) from "Ions" to "Neutrals" (and vice versa), the original parameter values remain stored. When you change back

to the previous operation mode, the stored values are reinserted. So it is not necessary to reserve specific M-Sets for the different PSM's.

Ions

Select the polarity (in mode Ions):

- Ve+: Positive polarity of the energy scan voltage.
- Ve-: Negative polarity of the energy scan voltage.

Neutrals

Select the mode for the detection of Neutrals:

- Coarse: Normal operating mode for measuring Neutrals. The minimum step width of the extraction voltage EXTR is 2 Volts.
- Fine: Measurement of Neutrals with accurately adjustable extraction voltage EXTR (minimum step width is 0.02 Volts).

Tune PPM 422 menu



Fig. 5-65

File

Store the ion source parameters, leave the editor.

Display

Display measurement data. (See Chapter 5.9.2 Optimizing the parameters (Measure Data), 369)

Option

Copy the channel and IS set parameters, set default parameters, single- or common increment/decrement of potentials.

File submenu



Fig. 5-66

Save

Store the ion source parameters.

Print Ion Source Settings

Print the current ion source settings.

Close

Close the parameter file and return to the Tune Up main menu.

Display submenu



Fig. 5-67

Measure Data

Start the measurement. (See Chapter 5.9.2 Optimizing the parameters (Measure Data), 369)

Option submenu




Fig. 5-68

Copy/Exch. IS-Set...

Copy/exchange IS parameter sets.

Copy/Exch. Channel...

Copy/exchange channel parameters.

Default IS-Set

Reset the ion source parameters to default values.

Single Incr/Decr

Increment/decrement the potentials of the IS-Voltages group individually.

Common Incr/Decr

Increment/decrement the potentials of the IS-Voltages group together, that are marked with ****....****. This simplifies the coarse correction of the energy scale without losing a previous optimization.

5.9.1.2

Copy/Exch. IS-Set

Copy/Exch IS-Set opens a dialog box for selecting an E-set or M-set. This dialog box is skipped when the cursor is positioned within the Ion Source group on E-Set or M-Set.




Fig. 5-69

Choose [Copy/Exch. E-Set] or [Copy/Exch. M-Set] to open a dialog box for copying or exchanging the selected E-Sets or M-Sets:




Fig. 5-70

Mode Selection

Select the desired action:

- Copy: Copy an IS set to a different set or to the other 15 sets of the selected QMS E or QMS M.
- Exchange: Exchange all parameters of the IS set selected under Source with those in the IS set selected under Dest. (within the selected QMS-E or QMS-M).

The desired IS set can be selected under Source and Dest:

Source

IS set whose parameters are to be copied.

Dest

Target IS set to which the parameters are to be copied. For copying to all IS sets, click on the [UP] or [DOWN] arrow with the right mouse button to choose ALL.

5.9.1.3 Copy/Exch. Channel

Choose Copy/Exch. Channel to open a dialog box for copying individual channel parameters or for copying all parameters of the current channel:




Fig. 5-71

Mode Selection

Select the desired action:

- Copy: Copy the object (Channel or Parameter) selected under Object Selection to a different channel or to all channels.
- Exchange: Exchange all parameters of the channel selected under Source with those of the channel selected under Dest.. You can choose Exchange only in conjunction with [Object Selection] > [Channel].

Object Selection

Choose the object to be copied:

- Channel: Copy or exchange whole channels (all parameters).
- Parameter: The selected parameter is copied (exchange is impossible). Parameter cannot be chosen if no parameter has been selected.

The channels are selected under Source and Dest:

Source

Enter the channel whose parameters are to be copied.

Dest

Enter the destination channel to which the parameters are to be copied. For copying to all channels click on the [UP] or [DOWN] arrow with the right mouse button to choose ALL.

Special cases

In the parameter copying functions the following Special cases must be taken into consideration:

[Scan Mode] > [Mass/Energy]

Cannot be copied as a parameter.

Threshold

It is only possible to copy into channels in which the same Detector has been selected.

5.9.1.4

Default IS-Set

[Option] > [Default IS-Set] loads standard values for the E-Set and the M-Set.

NOTE:

[Option] > [Default IS-Set] deletes and overwrites the current E-Set and M-Set without request for confirmation. Do not choose this option unless you are sure, you don't need the parameter values of the chosen IS-Sets any more.

5.9.1.5 Single/Common Incr/Decr

Single Incr/Decr

By choosing Single Incr/Decr the Ion source potentials of the IS-Voltages group can be modified individually.

IS-Voltages [V]	
CATH C	—
IONREF	—
ENTRY	14.50
ITRO	420
EXTR	—
FOC 1	120
CENTR	44.0
MIRR	24
INFL	22.50
FOC 2	14.00
CASE	26
F.AXIS	14.00
DFLI	336
DFLO	145

Fig. 5-72

Common Inc/Decr

If Common Incr/Decr is chosen, a modification of a single potential influences all potentials of the IS-Voltages group that are marked by “** ... **”.

The step width for all those potentials is ± 2.0 V then.

IS-Voltages [V]	
CATH C	—
IONREF	—
ENTRY	** 14.50
ITRO	** 420 **
EXTR	—
FOC 1	** 120 **
CENTR	** 44.0 **
MIRR	** 24 **
INFL	22.50
FOC 2	** 14.00
CASE	** 26 **
F.AXIS	** 14.00
DFLI	** 336 **
DFLO	** 145 **

Fig. 5-73

Difference between Common Incr/Decr and SHIFT

With Common Incr/Decr all potentials marked with “** ... **” can be shifted in steps of ± 2 Volt. This simplifies the coarse correction of the energy scale because a previous optimization is not destroyed. The changed potentials are stored in the same way as in Single Incr/Decr.

SHIFT shifts the energy scale into the negative, without affecting the energy calibration. The scale and the measurement curve are shifted jointly, so that distributions close to zero can be watched. For that purpose, an intentional misadjustment (all potentials marked with “** ... **”) in the range of 0 ... +10 Volt is created and the resulting deviation between the energy scale and V_e scale is corrected. Shift is not stored. In measuring sequences, Shift can be calculated and set automatically based on the selected parameters.

NOTE:

For detailed information, please refer to the User's guide of the Plasma Process Monitor PPM 422, Interrelation between ion energy and V_e .

5.9.1.6 Relations (Ion energy and V_e , CORR and SHIFT)

The detailed interrelation is described in the User's guide of the PPM 422 Plasma Process Monitor. The following Examples provide a brief overview:

5.9.1.7 Measurement parameters for determining the energy zero

Set-up of the Quadstar 32-bit PPM Tune Up measurement parameters for determining the energy zero in the Sampling Mode Neutrals, Scan Mode Energy and Fixed Mass 40 (Argon peak):




Fig. 5-74




Fig. 5-75

5.9.1.7.1 Energy scan in Sampling Mode Neutrals

Fixed mass = 40 (Argon peak), CORR = 0, SHIFT = 0.

In this measure mode you can determine the CORR Value graphically: choose [Function] > [CORR] from the measure window menu, move the crosshairs to energy zero while having the left mouse button pressed, and press the right mouse button additionally. The crosshairs move only within ± 2.0 eV. In the now appearing Set CORR Value dialog box you can adopt the determined value by pressing [OK].




Fig. 5-76

5.9.1.7.2 Effect of the correction value CORR

CORR can be entered by [Operation] > [Set CORR] or by the mouse into the dialog box *Set CORR Voltage* shown below. If you choose [OK], the value is adopted:




Fig. 5-77

CORR = 1.4 V shifts the measurement curve relative to the energy scale by this amount:




Fig. 5-78

Since CORR changes neither the V_e scale nor the IS voltages, the measurement curve does not shift relative to the window frame (First V_e always meets the starting point of the horizontal scale of the window).

5.9.1.8 Measurement parameters for positive ions

The following picture shows the measurement parameters for measuring positive ions produced in the argon sputter plasma (Sampling Mode Ions +, Scan Mode Energy and Fixed Mass = 40 [argon peak]):




Fig. 5-79

5.9.1.8.1 Energy distribution of positive ions

The following picture shows the energy distribution of positive ions produced in the argon sputter plasma. The different curves originate from the optimization process:




Fig. 5-80

5.9.1.8.2 Effect of SHIFT

SHIFT = 2Volts influences the groups IS-Voltages and Special as follows:

SHIFT = 0 Volt		SHIFT = 2 Volt	
IS-Voltages [V]		IS-Voltages [V]	
CATH C	—	CATH C	—
IONREF	—	IONREF	—
ENTRY	30.00	ENTRY	32.00
ITRO	300	ITRO	302
EXTR	—	EXTR	—
FOC 1	80	FOC 1	82
CENTR	38.0	CENTR	40.0
MIRR	16	MIRR	18
INFL	25.75	INFL	25.75
FOC 2	5.00	FOC 2	7.00
CASE	98	CASE	100
F.AXIS	9.00	F.AXIS	11.00
DFLI	298	DFLI	300
DFLO	118	DFLO	120
Special		Special	
SHIFT	0	SHIFT	2
CORR	1.4	CORR	1.4
ORIF	0.0	ORIF	0.0

Fig. 5-81

The next picture shows the same measurement data (energy distribution of positive ions produced in the argon sputter plasma) after SHIFT = 2 Volts has been entered.




Fig. 5-82

The measurement curve does not shift relative to the energy scale. Both are shifted the same way relative to the window frame (or relative to the V_e scale that is not explicitly shown). By using the [Function] > [Cursor] - Mode (See Chapter 5.9.2.1.3 Function sub-

menu (crosshairs), [Fig. 372](#)) you can determine the measured values directly within the graphic: move the crosshairs to the desired location while having the left mouse button pressed and press the right mouse button additionally. The Values are now shown in the status bar of the measure window: (X: 0.83, Y: 2.380E+02).

The same measurement with identical parameters in sample mode (E-Mode SAMPLE) on the maximum of the energy curve (fixed E: 0.83 eV) produces the following result:




Fig. 5-83

The field on the right in the status bar shows the measured value numerically (2.380E+02).

“Fixed Energy” setting

The energy value cannot be entered directly in the Tune Up menu. Fixed Ve is entered instead (See Section «Scan», [Fig. 357](#)). Fixed E is calculated and displayed by taking CORR and SHIFT into consideration (see [Fig. 5-83](#), [Fig. 369](#)). The desired energy value can be set by modifying Fixed Ve with the UP/DOWN keys.

5.9.2 Optimizing the parameters (Measure Data)

By choosing [Display] > [Measure Data] you can open a measure window and start a measurement on the currently used Channel.

You can perform measurements manually and optimize the parameters by means of the measurement. See also example in Chapter 5.9.1.8.1 Energy distribution of positive ions, [Fig. 367](#).

Scan Analog measurement data




Fig. 5-84

In the display, the old measurement data are displayed together with the new measurement data in a different color. The display field is not cleared before each new measurement. The displayed scans are only erased if:

- the size of the Measure Data window is changed,
- either the First, Width or Range parameter is changed, or
- if the user explicitly requests this by choosing [Display] > [Refresh].

NOTE:

For printing pictures in Quadstar 32-bit (by Ctrl+Shift+F12), stop the measurement by [Operation] > [Halt] in order to prevent errors of the interface.

Displaying reference measurement data

In measurement modes SCAN-N and SCAN-F also the measurement data of a reference scan can be displayed. If the option Reference (menu Display) is checked, a reference scan is stored automatically

- if the first scan after Measure Data has been chosen in the menu Display, or
- when either the First, Width or Range parameter has been changed.

You may store a reference scan manually as well, by [File] > [Save Reference]. This is necessary for example, to make the option Reference accessible again, if it is grayed.

Sample measurement data




Fig. 5-85

In this measurement the highest measured intensity is flagged with a drag pointer. The measurement Operation status (Go or Halt) and the measured value (Val: 9.864E+05) are displayed in the status bar.

5.9.2.1 Measure Data menu



Fig. 5-86

File

Store reference data, quit measure window.

Operation

Start/stop the measurement, set CORR.

Function

Setup the function mode of the crosshairs.

Setup

Scale Y-axis, set its display range, choose colors.

Display

Refresh the display, display reference data.

5.9.2.1.1 File submenu



Fig. 5-87

Save Reference

Store the measurement data of the next scan as reference data. This command is only accessible in Scan Analog mode (SCAN-N and SCAN-F).

Close

Quit the measurement window.

5.9.2.1.2 Operation submenu



Fig. 5-88

Go

Start the measurement.

Halt

Stop the measurement.

Set CORR...

Enter the CORR parameter for the selected E-set.

Set CORR... parameter input

Choose Set CORR to enter the CORR parameter. CORR performs a software shift of the energy scale relative to the V_e scale. See Chapter 5.9.1.7.2 Effect of the correction value CORR, [365](#).




Fig. 5-89

NOTE:

This dialog box can as well be invoked in mode Scan Energy, Function: CORR by mouse (See Chapter 5.9.1.7.2 Effect of the correction value CORR, [365](#)).

5.9.2.1.3 Function submenu (crosshairs)



Fig. 5-90

CORR

The crosshairs are used to determine and input the CORR-value, restricted to X-values within ± 2.0 . The X- and Y- values are not displayed in the status bar. This mode is only accessible in scan mode Energy, sampling mode Neutrals.

Cursor

The crosshairs are used to determine the X- and Y- values of a certain point in the measured scan. The values are displayed in the status bar.

You bring up the crosshairs by holding the left mouse button pressed and moving the mouse cursor into the graphic. The crosshairs now follow the mouse cursor and are quite exactly positionable on a certain point of the graphic. With a click of the right mouse button the crosshairs are released and the last displayed values in the status bar are frozen (e.g. useful for printing) or the CORR-dialog box appears.

If CORR and Cursor are both possible, you can switch between these two function modes by a click of the right mouse button.

NOTE:

The X- and Y-values are displayed with constant number of significant digits. The real resolution may though be far lower, depending on the graphic; it is generally determined by the ratio of pixelwidth to window width/height.

5.9.2.1.4 Setup submenu



Fig. 5-91

Axis...

Display the Setup dialog box to scale the Y-axis (lin/log). This setup dialog box may as well be activated with a double-click of the left mouse button.




Fig. 5-92

- Lin: Display measurement data on a linear scale.
- Log: Display measurement data on a logarithmic scale.

Display Range...

Setup the display range of the Y-axis.




Fig. 5-93

- Disp. F.S.R.: Upper limit for the values to be displayed.
- Disp. Decades: Number of decades to display.

Colors...

Choose the colors to draw the scans with.




Fig. 5-94

In the fields #0 ... #63, enter the colors to draw the scans with (alternatingly). The standard color table, that contains eight colors, produces scans with good contrast that are easily distinguishable in most cases. It can be recalled by [Set Default Colors], but the current settings of the Alternating Colors for the current color palette will be lost then. See Chapter 3.8.5 Screen Colors, 245.

5.9.2.1.5 Display submenu

This menu is only available in Scan Analog Mode (SCAN-N and SCAN-F):



Fig. 5-95

Refresh

Erase the graphic that's already drawn, get new measure data and draw it.

Reference

Show reference scan on/off.

5.10 Communication Parameters

Under Communication the interface to the PPM 422 can be turned on or off:




Fig. 5-96

Connect

Establish connection to the PPM 422.

Disconnect

Cancel the connection to the PPM 422.

⚠ CAUTION	
	<p>The Plasma Process Monitor PPM 422 works only with communication type ArcNet. The communication type RS232C is reserved for service purpose and does not provide a correct, safe operation of the Plasma Process Monitor.</p> <p>For safety reasons, no electrical conductive connection between PC and PPM is allowed. ArcNet and fibre optics provide this, while RS232C cables don't!</p>

NOTE:

For detailed information please refer to **Chapter 3.6 Communication parameters**, **209**.

5.11 Setup Parameters

Setup menu



Fig. 5-97

Choose Setup to turn the emission and the SEM high voltage on or off (and to enter the Common SEM Voltage).

5.11.1 Turning the SEM high voltage and emission on/off

SEM/Emission Control...

Display the SEM/Emission Control dialog box in which the emission and the SEM high voltage can be turned on/off and the Common SEM Voltage can be entered. See Chapter 3.7 Setup parameters, [p. 223](#).

NOTE:

This dialog box can also be called directly by pressing <Ctrl+S>.

SEM and Emission can only be switched on or off, when the communication with the PPM 422 is enabled.

6. Service

Quadstar 32-bit Service provides test measurements (in the Measure menu) and manual operation of analog and digital input and output channels (in the Manual menu). The available options depend on the quadrupole mass spectrometer (QMS):

QMS 422/421

DO 420, DI 420, AI 421, AO 421, PE 420 and PI 420

QMS 200

PKR 250/251 or TPR 250 (total pressure gauges). The QMS 200 is equipped with two analog inputs and four analog outputs, two digital inputs and two digital outputs as a standard feature. The only possible options are PKR 250/251 or TPR 250.

When the program is started, Service automatically activates the communication with the QMS and reads in the hardware configuration. The individual functions of Measure and Manual are enabled only if the required hardware is present.

Main menu



Fig. 6-1

File

Display the version number of the Quadstar 32-bit Service program and terminate the program. See Chapter 6.1.1 About Service, 378.

Measure

Start the service programs

- Leak Test - See Chapter 6.2.1 Leak Test, 379.
- Total Pressure - See Chapter 6.2.2 Total Pressure, 382.

Manual

Operate/display the DI/DO and AI/AO channels. See Chapter 6.3 Manual, 386.

Comm

Switch the interface to the QMS on/off, choose the measuring device. See Chapter 6.4 Communication Parameters, 389.

Setup

Switch the emission and the SEM on/off. See Chapter 6.5 Setup Parameters, 389.

Help

Display the 'Key Fragment Ions' table and help contents. See Chapter 1.9.1 Spectrum interpretation, 56.

6.1 File



Fig. 6-2

6.1.1 About Service

Display the version number of the Quadstar 32-bit Service program:



Fig. 6-3

Exit
Quit the Quadstar 32-bit Service program.

6.2 Measure



Fig. 6-4

The measurements can only be started if the communication link to the QMS is switched on:

Leak Test

Find leaks in the vacuum system.

Total Pressure

Measure and display the total pressure.

6.2.1 Leak Test

Leak Test provides assistance in the search for a leak in the vacuum system. A special test gas (helium, for example) can be used for locating the leak. The measurement data are displayed Versus Time (on the left) and as a bar (on the right):




Fig. 6-5

File

- Close: Close file, terminate measurement.

Parameters

- Setup: Control the measurement data display. This dialog box can also be activated by a double-click of the left mouse button.
- Channel: Channel parameters (Range, Mass etc.). This dialog box can also be activated by a double-click of the right mouse button.

Alarm display

The alarm threshold set in the Channel Parameters dialog box under Trip Level B is indicated with the symbol 'B>' in the bar chart. It can be switched on/off by the enable/disable Alarm Display control field under [Leak Test] > [Setup].

6.2.1.1

Setup




Fig. 6-6

Y-Scale

- Lin: Display the measurement data on a linear scale.
- Log: Display the measurement data on a logarithmic scale.

Time Scale

Time scale on the X-axis.

X-Grid

X-grid in the measurement window on/off.

Y-Grid

Y-grid in the measurement window on/off.

Alarm Display

Indication of the alarm threshold and the alarm state on/off.

Alarm Sound

Alarm buzzer on/off.

Colors

Assign colors to the two alarm states. The colors can be entered numerically or by a mouse click on the color palette.

- Alarm On: Color of the bar and the alarm symbol if there is an alarm.
- Alarm Off: Color of the bar and the alarm symbol if there is no alarm.

6.2.1.2

Channel

Parameter	Value
Mass	14.00
Resolution	25
Amp. Range	E-09
Amp. Mode	FIX
Amp. Offset	ON
Dwell	0.1s
Trip Level B	8.00E-10
Detector	CH-TRON
SEM Voltage	1250
Disp. F.S.R.	AUTO
Disp. Decades	2

Buttons: Default, Adjust, OK, Cancel

Fig. 6-7

Command buttons**[Default]**

Set all parameters to default values (test gas = helium)

[Adjust]

Optimize the selected mass number to the peak maximum.

Parameters

Number and type of displayed parameters depend on the options selected in the Parameter Setup Program under [Config] > [QMS...]. The following parameters are possible:

Mass

Mass number on which the specific measurement is conducted.

Resolution

Resolution of the quadrupole mass filter.

Amp. Range

Electrometer measurement range.

Amp. Mode

Amplifier operation mode (automatic/manual).

Amp. Offset

Amplifier offset.

Dwell

Measurement time.

Trip Level B

Alarm threshold B.

Detector

Detector type.

SEM Voltage

SEM or CH-TRON high voltage.

Activate the Common SEM Voltage (<< XXXX>>) by clicking on the [UP] or [DOWN] arrow with the right mouse button.

DO# B

Digital output channel for the alarm status message.

Disp. F.S.R.

Upper limit for display.

Disp. Decades

Number of decades to display when a logarithmic scale is used.

**Special parameters
(QMS 422/421)**

When a QMS 422/421 is used and additional options are installed, then there are special parameters to operate those options:

**Ion counter CP
400**

CP Level

Response threshold of the ion counter preamplifier.

**Multiplex
operation with
QMU 112**

QMU Channel

Channel selection in multiplex operation with QMU 112.

**Operation with
chopper lock-in
amplifier CLA 421**

Lock In

Switch the CLA 421 Lock In amplifier on/off. OFF means normal electrometer mode.

Phase

Phase shift (0° or 180°) of the Chopper Lock In amplifier.

6.2.2**Total Pressure**

Total Pressure shows the total pressure in the vacuum system, provided the required hardware options (measurement gauge, additional plug-in board for QMS 422/421) are installed.

6.2.2.1 QMS 422/421: PE 420 and/or PI 420




Fig. 6-8

6.2.2.2 QMS 200: TPR 250 or PKR 250/251




Fig. 6-9

Parameters

- Setup: Control the measurement data display. This dialog box can also be activated by a double click of the left mouse button.

6.2.2.3 Setup

[Parameters] > [Setup] controls the display of measurement data. This dialog box (Total Pressure Setup) can also be activated by double-clicking the left mouse button:


QMS 422/421

Fig. 6-10


QMS 200

Fig. 6-11

Unit of Total Pressure

Unit of measure (mbar, Pa(scal), Torr) for total pressure measurement with the PI 420 and/or PE 420 when the QMS 422/421 is used, or a PKR 250/251 or TPR 250 when a QMS 200 is used.

Time Scale

Time scale on the X-axis

X-Grid

X-grid on the measurement screen on/off.

Y-Grid

Y-grid on the measurement screen on/off.

Marker

Mark the measured values with a symbol. The symbol is selected automatically.

Colors

Assign colors to the different sensors. The color can be entered numerically or by a mouse click on the color palette.

- Cold cathode: Sensor (Cold cathode Gauge) of the PE 420 for QMS 422/421.
- Pirani #1: Sensor 1 (Pirani Gauge #1) of the PI 420 for QMS 422/421.
- Pirani #2: Sensor 2 (Pirani Gauge #2) of the PI 420 for QMS 422/421.
- PKR 250/251: Compact Full Range Gauge for QMS 200.
- TPR 250: Compact Pirani Gauge for QMS 200.

6.3 Manual

Choose Manual to display analog or digital input channels and manually change the output of analog or digital output channels. For that, the communication with the mass spectrometer must be switched on.



Fig. 6-12

6.3.1 Digital I/O (DI/DO)

For turning the individual digital output channels on or off, double-click the corresponding push button with the left mouse button. The switching state is indicated by the colored dot on the push button: green = on, grey = off.

QMS 422/421




Fig. 6-13

Display left

Select the DI/DO group to be displayed on the left side of the screen.

Display right

Select the DI/DO group to be displayed on the right side of the screen.

A group cannot be displayed at the left and the right side simultaneously.

QMS 200

Only the two standard digital input and output channels are available:




Fig. 6-14

The digital input and output channels of the QMS 200 have the following designation in the hardware description:

- DI 00: EXT_PROT (Connection TPR/PKR, Pin No. 1 <=> 2)
- DI 01: EXT_RUN (Connection Control, Pin No. 10 <=> 15)
- DO 00: REL_1 (Connection Control, Pin No. 11 <=> 12)
- DO 01: REL_2 (Connection Control, Pin No. 13 <=> 14)

6.3.2 Analog I/O (AI/AO)

The input and output voltages of the analog channels are displayed numerically and graphically. Please note, that the AI channels might not display 0.0 V, if their input is open (floating). The single channels can be configured and given names in the Parameter Setup program under [Config] > [AI Characteristic Curve] and [Config] > [AO Characteristic Curve].

NOTE:

The values shown here are the real input and output voltages in Volts, not the physical values (such as “ °C “ for temperatures) assigned to them by the Characteristic Curve.

QMS 422/421




Fig. 6-15

QMS 200

Only the two analog input and four output channels provided as a standard feature are available:




Fig. 6-16

6.4 Communication Parameters



Fig. 6-17

Connect

Establish connection to the QMS.

Disconnect

Cancel connection to the QMS.

SetUnit...

Select current measuring device. See Section «Setting the node addresses of several control units», 218.

NOTE:

For detailed information see Chapter 3.6 Communication parameters, 209.

6.5 Setup Parameters



Fig. 6-18

SEM/Emission Control...

Open the dialog box SEM/Emission Control to:

- Switch the SEM high voltage on/off.
- Input the Common SEM voltage.
- Switch the Emission on/off.

See Chapter 3.7 Setup parameters, 223.

7. Utility

Quadstar 32-bit Utility contains the spectra library and an installation utility.

The spectra library allows to store, load and edit spectra as well as to compare them with other spectra. You can export the library to a text file. The installation utility is used for installing custom-specific Quadstar 32-bit extensions.

Main menu



Fig. 7-1

File

Display the version number of Quadstar 32-bit Utility and quit the program.

Spectra Library

Call the Spectra Library File Manager or export the library to a text file. See Chapter 7.2 Spectra library menu, [393](#)

Installation

Install custom-specific extensions. See Chapter 7.3 Installation menu, [408](#).

Help

Display the 'Key Fragment Ions' table and help contents. See Chapter 1.9.1 Spectrum interpretation, [56](#).

7.1 File menu



Fig. 7-2

7.1.1 About Utility



Fig. 7-3

Exit

Quit the Quadstar 32-bit Utility program.

7.2 Spectra library menu

The spectra library supplied with Quadstar 32-bit contains 72 spectra and can be expanded to up to 10,000 spectra. New spectra can be entered manually ([New] or [Edit] command button), or by measuring a spectrum under Scan Bargraph (Quadstar 32-bit Measure program).

The spectra library features the following functions:

- Displaying stored spectra.
- Modifying stored spectra.
- Inserting new spectra.
- Subtracting spectra.
- Removing spectra.

The following information is stored for each spectrum:

- Title of the component.
- Formula of the component.
- Spectrum of the component (max. 32 peaks normalized, the largest peak corresponding to 100%).

Spectra Library menu



Fig. 7-4

Open Library...

Call the Spectra Library File Manager to work with one or more spectra of the library. See Chapter 7.2.1 Open Library, 394.

Export Library to Text File...

Export the complete Quadstar 32-bit Spectra Library to a text file that is readable e.g. by the NIST mass spectral search program. See Chapter 7.2.3 Export Library to Text File, 406.

7.2.1 Open Library




Fig. 7-5

Library State

State of the library (sorted/unsorted).

Stored Spectra

Number of spectra stored in the spectra library.

Please note that any change in the stored spectra or the input of a new spectrum sets the Library State to UNSORTED. In order to work efficiently with the library, you should Sort it.

The following functions are available for modifying the spectra library:

[Spectrum]

Display the selected spectrum - See Chapter 7.2.2 Displaying a spectrum, 395.

[Copy Spec]

Copy the selected spectrum - See Chapter 7.2.2.2 Copying a spectrum, 401.

[Remove]

Remove a spectrum - See Chapter 7.2.2.3 Removing a spectrum, 402.

[Search]

Search a spectrum by the known formula - See Chapter 7.2.2.4 Searching a spectrum, 402.

[Goto]

Select a spectrum by the known spectrum number - See Chapter 7.2.2.5 Go to a spectrum, 403.

[Sort]

Sort the library - See Chapter 7.2.2.6 Sorting the library, 403.

[New]

Enter a new spectrum - See Chapter 7.2.2.7 Entering a new spectrum, 403.

[Edit]

Modify the selected spectrum - See Chapter 7.2.2.8 Editing a spectrum, 404.

[Print]

Print spectra - See Chapter 7.2.2.9 Printing the spectra, 405.

7.2.2 Displaying a spectrum

[Spectrum]

Open the window SPECTRUM Comparison.

The selected spectrum (in our example: Argon) is displayed as S1 in the upper half of the window:




Fig. 7-6

In the lower half of the screen you can load a comparison spectrum by [Spectrum] > [Select/Load S2].

SPECTRUM Comparison Menu



Fig. 7-7

Spectrum

Load and compare spectra.

Subtraction

Subtract and normalize spectra.

Parameters

Set up the display (lin/log, X/Y grid on/off)

7.2.2.1 Spectrum submenu




Fig. 7-8

Select/Load S2...

Select the spectrum to be displayed under S2.

Display Numeric...

Display the spectrum as a table.

Overlap S2 on S1

Load and compare spectra.

S1 to Library...

Store the spectrum displayed under S1 in the spectra library.

Close

Close the Spectrum Comparison window.

7.2.2.1.1 Load a comparison spectrum (Select/Load S2)




Fig. 7-9

[Spectrum]

Return to Spectrum Comparison, display the selected spectrum under S2. Double clicking the spectrum does the same.

[Search]

Search a spectrum by its known formula - See Chapter 7.2.2.4 Searching a spectrum, 402.

[Goto]

Select a spectrum by its known number. See Chapter 7.2.2.5 Go to a spectrum, 403.

The selected spectrum is shown in the lower half of the window as S2:




Fig. 7-10

7.2.2.1.2 Display spectrum as a table (Display Numeric)

You can display the spectra as a table:




Fig. 7-11

7.2.2.1.3 Overlap S2 on S1

Overlap spectrum S2 to spectrum S1:




Fig. 7-12

Both spectra are displayed in the upper half of the screen, differences are easy to see.

7.2.2.1.4 Store spectrum (S1 to Library)

S1 to Library stores the spectrum displayed as S1 in the upper half of the window to the library. Before it is stored definitely, you can check it by the (automatically displayed) Edit window and make last changes if necessary. See Chapter 7.2.2.8 Editing a spectrum, 404.

Subtraction submenu



Fig. 7-13

Subtract S1-S2...

Subtract spectrum S2 from spectrum S1. The result is displayed as S1.

Subtraction List...

Display the subtractions that have already been performed.

Undo last Subtraction...

Reverse the last subtraction.

Normalize S1...

Normalize the spectrum displayed under S1 (largest peak = 100%).

7.2.2.1.5 Subtract S1 - S2

Before the subtraction you must determine the mass number to be eliminated:




Fig. 7-14

Mass to eliminate

S2 is normalized in such a way, that the chosen mass number (in our example: mass 18) is eliminated in spectrum S1. Only a mass number, that appears in both spectra, can be chosen. The ratio of their intensities can be checked in advance by Overlap S2 on S1. See Chapter 7.2.2.1.3 Overlap S2 on S1, 398.

Ratio

Calculated normalization factor.

[Subtract]

Start the subtraction.

The result of the subtraction is displayed in the upper half of the display as S1.




Fig. 7-15

NOTE:

In the logarithmic representation, the upper spectrum is always displayed in the range of -100% ... +100%.

- The last subtraction can be reversed by [Subtraction] > [Undo last Subtraction].
- Before another subtraction is performed, the spectrum should possibly be normalized again if the largest peak (=100%) has been changed. See Chapter 7.2.2.1.8 Normalize S1, 400.

7.2.2.1.6 Subtraction List

Subtraction List shows a list containing the subtractions that have been performed so far.




Fig. 7-16

7.2.2.1.7 Undo last Subtraction




Fig. 7-17

With Undo last Subtraction you can reverse the previous subtraction. In our example, spectrum nbr. 48 (H2O: water) has been subtracted from S1 and mass 18 (Mass to eliminate = 18) has been eliminated.

7.2.2.1.8 Normalize S1

Normalize the spectrum so that the intensity of the largest peak corresponds to 100%:




Fig. 7-18

Parameters submenu



Fig. 7-19

Setup opens a dialog box for controlling the display




Fig. 7-20

Y-Scale

- Linear: Display the measurement data on a linear scale.
- Logarithmic: Display the measurement data on a logarithmic scale.

X-Grid

Grid for the X-axis on/off.

Y-Grid

Grid for the Y-axis on/off.

Colors

- Spec S1: Color of the spectrum S1.
- Spec S2: Color of the spectrum S2.

7.2.2.2 Copying a spectrum

[Copy Spec]

Copy spectrum S1 to another, already existing spectrum.




Fig. 7-21

NOTE:

The target spectrum, whose number is entered, is replaced by S1. The original spectrum (in our example: nbr. 48) will be lost.

7.2.2.3 Removing a spectrum

[Remove]

Remove a spectrum from the library by deleting it.



Fig. 7-22

7.2.2.4 Searching a spectrum

[Search]

Search a spectrum by its known formula.




Fig. 7-23

You can search a spectrum by:

- A known formula (such as Air, H₂O, N₂, O₂, etc.)
- The starting characters of the formula name and '*' (such as H₂*, SO₂*, C₂Cl*, C₂F*, etc.)

The search function does not distinguish between upper and lower case letters. The search works only with a sorted library.

7.2.2.5 Go to a spectrum

[Goto]

Go to a spectrum, whose spectra number is known.




Fig. 7-24

NOTE:

By "Print Spectra" you can create a list of all spectra, which provides the required numbers quickly. See Chapter 7.2.2.9 Printing the spectra, 405.

7.2.2.6 Sorting the library

[Sort]

Sort the spectra library.



Fig. 7-25

The formulas are sorted alphabetically, the peaks by their intensities.

NOTE:

Sort the library, when you have entered all changes. Depending on the number of stored spectra, the sort may take some time.

7.2.2.7 Entering a new spectrum

[New]

Enter a new spectrum - See also Chapter 7.2.2.8 Editing a spectrum, 404).

7.2.2.8 Editing a spectrum

[Edit]

Edit a spectrum.

In the Edit <spectrum> dialog box you can edit an existing spectrum or enter a new one:




Fig. 7-26

[Insert]

Add to the current spectrum S1 another peak (mass number, intensity). The intensity must be entered in normalized form, i.e. in %. The maximum of peaks that can be present in a spectrum is 32.

[Spectrum]

Display the spectrum.

[Remove]

Remove a peak from the current spectrum S1:




Fig. 7-27

[Norm]

Normalize the spectrum (largest peak = 100%).

[Sort Mass]

Sort the spectrum by ascending mass numbers.

[Sort Intern]

Sort the spectrum by descending intensities.

[Spectrum]

Display the changed spectrum S1:




Fig. 7-28

Options submenu




Fig. 7-29

Linear

Display the measurement data on a linear scale.

Logarithmic

Display the measurement data on a logarithmic scale.

X-Grid

Grid for the X-axis on/off.

Y-Grid

Grid for the Y-axis on/off.

7.2.2.9 Printing the spectra

[Print]

Print one or several spectra as a table:




Fig. 7-30

from Nbr

First spectrum to be printed.

to Nbr

Last spectrum to be printed.

Options

- Print all information: Print all information of the selected spectra.
- Print Formula and Titles: Print only the number, the formula and the title of the selected spectra.

7.2.3 Export Library to Text File

[Spectra Library] > [Export Library to Text File] allows to export the complete Quadstar 32-bit Spectra Library to a text file. This text file can be read e.g. by the NIST Mass Spectral Search program that uses the NIST Mass Spectral Library. This way you can compare any spectrum stored in Quadstar 32-bit to any spectrum of an external library such as the NIST library.




Fig. 7-31

Choose path and...

Path and name of the text file to be created. The default path is the Quadstar 32-bit data path and the default name is 'qspeclib.msp'.

[File Manager]

Choose the file using the file manager instead of entering it manually.

Function State

Status of the export function.

[Start]

Start export.

[Cancel]

Leave dialogbox.

The resulting text file consists of the following elements:

- Name of the component, normally capitals (e.g. 'ACETIC ACID')
- Formula of the component (e.g. 'C2H4O2')
- Number of peaks for this component (e.g. '11')

- Masses and relative intensities of this component (e.g. '43 999; 45 935; 60 576; 29 156; 42 145; 28 56; 44 49; 31 45; 41 45; 46 14; 61 12; '). The mass number of the highest peak (here: 43 amu) is followed by its intensity (999), a semi-colon and a blank. Then the next mass number (45), intensity (935), semi-colon and blank follows etc. The highest intensity is always 999, all other intensities refer to this one. Intensities lower than 0.05% of the highest intensity are disregarded, those between 0.05% and 0.15% are displayed as '1' etc.

This structure is repeated for every component.

7.3 Installation menu

Install custom-specific extensions such as special sequences, parameter sets etc.:




Fig. 7-32

Current Installation

Current pathnames of your installation (they can be changed in Parameter Setup program under [Setup] > [General]).

Installation Disk

Information concerning the running installation (will only occur when the installation is started).

Source Drive

Source drive, from which the extensions are to be read.

[Install]

Start installation.

[Cancel]


Leave dialogbox.

For such installations proceed as follows:

- 1 Make sure the above declared paths for parameters, data etc. match the ones you desire.
- 2 Choose the drive to read from in the field Source Drive (in this example: Drive a:, which contains the floppy disk that includes the extensions to install).
- 3 Choose [Install] to install or [Cancel] to leave without taking action; once the installation process is started, it cannot be cancelled any more.

8. Process Control

NOTE:

If the Plasma Process Monitor PPM 422 is used, please see Chapter 8.8 Process Monitor PPM 422,  521.

8.1 Introduction

The Quadstar 32-bit sequencer is designed as a command interpreter. There are two methods for programming the individual sequencer lines (in program Parset, under [Sequence] > [Editor]):

- indirectly via dialog boxes ([Edit] > [Dialog]).
- through direct text input in the text editor.

A sequence essentially comprises a series of instructions (e.g. **MID()**, **ScanAnalog()**, etc.). An instruction comprises the name, followed by an opening and closing parenthesis. These two parentheses delimit the parameter input field.

NOTE:

Certain parameters are mandatory while others can be selected as options.

If an optional parameter is not specified explicitly, the sequencer inserts the default value (for example, setting “prot=off” is not necessary).

8.2 Programming

8.2.1 Input via the keyboard

This input method is intended for experienced users and requires thorough familiarity with the command syntax. The text is entered directly at the position in the sequence line selected with the cursor. Enter inserts a new line, Del , Backspace and Ins have the same functions normally found in text input programs.

8.2.2 Input via dialog editor

Set the cursor to an instruction (e.g. MID) and select [Edit] > [Dialog] or click with the right mouse button. The dialog editor containing all possible parameters of this instruction will show up. You may as well select [Edit] > [Dialog] on an empty line; in the now appearing window, choose the desired instruction and click [OK]. The dialog editor will be displayed then as well.

NOTE:

During the familiarization phase you should use the dialog editor. The advantage is, that incorrect entries are avoided to a great extent.

8.3 Variables

The Quadstar 32-bit sequencer works with the following types of variables:

Integer

Integer variables, range -1'000'000 ... +1'000'000.

Float

Floating point variables, range -1E+37 ... +1E+37, resolution 15 decimal places.

String

Text, stored as character strings.

Integer- as well as Float- variables are internally handled with higher accuracy than outlined above. But the input and output arguments of Quadstar 32-bit instructions must be within those ranges.

The variables are represented as arrays and, depending on the type, have either local or global validity.

These arrays have fixed names and are defined as follows:

i[]	Local integer variable array.
f[]	Local float variable array.
gi[]	Global integer variable array.
gf[]	Global float variable array.
gfa[]	Global float variable field size $N \times 64$ for storing measured values, calculated values etc.
gs[]	Global array of character strings for defining messages, error messages and file names. The length of a character string must not exceed 100 characters.

The size of these arrays can be entered in the Parameter Setup program under [Sequence] > [Editor] > [Options] > [Parameter Limits] and is designated as the configurable maximum in the description of the individual sequencer instructions.

NOTE:

The first group, the local arrays i[] and f[], are specific to the sequence. These variables are known only to the current sequence and the declared array size applies only to this sequence.

The second group (gi[], gf[], gfa[] and gs[]) are global arrays. The declared size applies to all sequences and the arrays are also known to the sequences called by Sequence() or SeqExit().

gs[] is a global array in which character strings such as filenames, messages etc. can be stored. The different array types can be loaded with the **SetString()** instruction.

The index of an array must be either an integer number or an element of the array **i[]** or **gi[]**; further nesting of indices is not allowed.

Examples

Right:

```
f[2], gfa[gi[12]], gfa[i[24]], gs[gi[5]], etc.
```

Wrong:

```
f[i[i[2]]], gfa[gi[i[3]]], gfa[f[2]], etc.
```

8.3.1 Local integer variables i[.]**Description:**

Local array of integer variables (known only to the current sequence)

Syntax:

```
i[xx]
```

Parameters:

xx - Index ($0 \leq xx \leq$ configurable maximum)

Example(s):

```
i[5], i[0], i[i[7]], ...
```

8.3.2 Global integer variables gi[.]**Description:**

Global array of integer variables (known also to sequences called by the **Sequence()** and **SeqExit()** instructions).

Syntax:

```
gi[xx]
```

Parameters:

xx - Index ($0 \leq xx \leq$ configurable maximum)

Example(s):

```
gi[22], gi[10], gi[i[12]], ...
```

8.3.3 Local float variables f[.]**Description:**

Local Array of float variables (known only to the current sequence).

Syntax:

```
f[xx]
```

Parameters:

xx - Index ($0 \leq xx \leq$ configurable maximum)

Example(s):

```
f[43], f[1], f[gi[9]], ...
```

8.3.4 Global float variables gf[.]

Description:

Global array of float variables (known also to sequences called by the **Sequence()** and **SeqExit()** instructions).

Syntax:

```
gf[xx]
```

Parameters:

xx - Index ($0 \leq xx \leq$ configurable maximum)

Example(s):

```
gf[12], gf[90], gf[i[0]], ...
```

8.3.5 Global field of float variables gfa[.][.]

Description:

Global field of float variables (known also to sequences called by the **Sequence()** and **SeqExit()** instructions).

Syntax:

```
gfa[xx][yy]
```

Parameters:

xx - Array vector number ($0 \leq xx \leq$ configurable maximum)

yy - Array index number ($0 \leq yy \leq 63$)

Example(s):

```
gfa[0][5], gfa[i[8]][gi[28]]
```

8.3.6 Global character string gs[.]

Description:

Global array of character strings (known also to sequences called by the **Sequence()** and **SeqExit()** instructions). The maximum length of a character string is 100 characters.

Syntax:

```
gs[xx]
```

Parameters:

xx - Index ($0 \leq xx \leq$ configurable maximum)

Example(s):

```
gs[0], gs[i[4]], gs[gi[2]]
```

8.4 Sequencer instructions

8.4.1 General information

8.4.1.1 Syntax of a sequencer instruction

Per line only one instruction is allowed. The instructions all have the same format. They consist of the name and a list of parameters (between parentheses):

```
Instruction(parameter list)
```

Example:

The instruction **Delay()** has one mandatory (*tim*) and two optional parameters (*unit* and *d*). The first optional parameter (*unit*) relates to the mandatory *tim* parameter (can be recognized by the separator ';'). By default, *unit* is set to **s** (seconds):

```
Delay(Time=tim[ ;unit][ ,Disp=d])
```

Delay Designates the name of this instruction.

(...) The name is followed by a parameter list between parentheses.

Time Time is a keyword. Since it's not located between parentheses, it is mandatory. Keywords are emphasized with bold letters and are always followed by a '=' character.

= Assign a value to the parameter

tim *tim* is the parameter value assigned to the keyword **Time**.

[] Brackets identify optional parameters. Optional parameters need to be declared only if they differ from the default value. These brackets have only symbolic value and are not actually entered.

{ } These braces identify configuration-dependent parameters. They need to be declared only if this is required by the configuration (hardware, communication type etc.). These braces have only symbolic value and are not actually entered.

Example: Define the Unit (QMS) if Quadstar 32-bit runs in multiplex mode

```
SetDO(var{ ;unit }=state[ ;tim[ ;unit ]][ ,Disp=d])
```

_ The underscore indicates that the command does not end here and that it is continued on the next line. It has symbolic value only and is not actually entered. Its purpose is to make the syntax and code lines more readable in this documentation.

Example:

```
DDERequest(Service=ServName , Topic=TopName , _  
           Item=ItemName , Data=Dat)
```

;
Is located inside the brackets and serves as a separator between the individual values if a parameter consists of multiple specifications such as **Time=5;min**. The character ';' indicates that **min** relates to **5**.

unit *unit* is an optional parameter. Only the specific character sequences **ms**, **s**, **min**, **h**, may be entered for this parameter.

Disp **Disp** is the second optional parameter. It is used to switch the display on or off and can be set to **[default]**, **off** or **on** with '='.

, When a keyword is entered manually, it must be preceded by a comma.

Examples:

```
Delay(Time=10)
Delay(Time=10;min)
Delay(Time=10;min,Disp=on)
Delay(Time=10;min,Disp=default)
```

8.4.1.2 Default values

In the description of the instruction parameters, the default values (they don't have to be declared explicitly) are listed between brackets. An example is the default value "seconds" in an instruction containing a time specification. The Delay instruction e.g. has the following syntax:

```
Delay(Time=tim[ ;unit][ ,Disp=d])
```

The parameters ;*unit* and ,**Disp**=*d* in brackets are optional. If they are not declared, the compiler sets them to default values:

- Parameters without a keyword (such as *unit* etc.) are initialized to a default value that cannot be changed.
- Parameters with a keyword (such as **Disp**, **Prot** etc.) are initialized to a default value that can be altered by the **SetDefault()** instruction.

In the sequencer line:

```
Delay(Time=10;min)
```

the time specification 10 is interpreted as minutes because of the subsequent (optional) parameter min.

The following initialization table applies to the **SetDefault()** instruction:

Keyword	Designation	Initialization value
Disp	Display mode	tab
Prot	Protocol mode	off

The default value is loaded with the initialization value when the main sequence is started and is valid also for the called sub-sequences. This way it is possible, for example, to switch the printer log on or off from a central location, i.e. it does not have to be specified in each instruction.

If you are programming by [Edit] > [Dialog], you can choose the values from list boxes. Optional parameters are accessed by pushbuttons (e.g. [Disp]). Without an explicit declaration, all optional parameters are set to their default values and are not displayed. If you use the dialog editor to program a **Delay()** instruction e.g.




Fig. 8-1

and you click on [Options] > [Disp] > [default],




Fig. 8-2

this means 'Disp not specified'; ",Disp=default" will not be displayed therefore:

```
Delay(Time=10;min)
```

8.4.1.3 File names with indices

In all sequencer instructions that contain file or program names, an index may be appended to those names. This enables you e.g. to

- save measure data to files, whose names are consecutively 'numbered' (indiced).
- perform a series of different measurements by a single loop, if the corresponding parameter files are indexed.
- execute different indexed batch programs comfortably (by using one loop).

The file index is a number in the range $0 \leq Index \leq 1'000'000$, that is appended to the file's or program's name. Please note that this new name including the index number may not exceed 255 characters.

The actual file name is built using the following rules:

- File name with index contains less than 8 characters:
Zeros are inserted before the index until the length equals 8 characters.
- File name contains less than 8 characters, and the file name with index contains more than 8 characters:
Index is simply appended to the file name.
- File name contains 8 or more characters:
Index is extended to 5 digits with leading zeros and then appended to the file name.

Examples:

```
'Name="air";5'  
results in "air00005"
```



```
'Name="c:\par\air";17'
```

results in "c:\par\air00017"

```
'Name="c:\par\airdemo";5'
```

results in "c:\par\airdemo5"

```
'Name="c:\par\airdemo";125'
```

results in "c:\par\airdemo125"

```
'Name="c:\par\newairdemo";125'
```

results in "c:\par\newairdemo00125"

You may set the index by:

cc Numeric value (integer)

i[xx] Element xx of the local integer array

gi[xx] Element xx of the global integer array

Application examples:

```
// Start 'batprog1', then 'batprog5'
Execute(Prog="batprog";1)
Execute(Prog="batprog";5)
// Measure 10 AnalogScans by using the parameter files 'airmeas0'
// to 'airmeas9' and save the data in 'scan0000' to 'scan0009'
Loop(i[0]=0;9)
    ScanAnalog(Par="c:\par\airmeas";i[0], SaveCyc="scan";i[0])
```

8.4.1.4 String values

In all sequencer instructions that contain strings, those strings may be extended by so called string values: any data object (text, integer number, float number; direct input or element of an array) is being defined as the string value and inserted at a certain position of the original string. The position to insert the string value is marked by an '@' character. The first encountered '@' character in the string is replaced by the string value. If the string doesn't contain an '@' character, the string value is appended to the original string.

Examples:

The string value is defined as the integer number '422':

```
'Name="q@name";422'
```

results in the text "q422name".

```
'Name="qname";422'
```

results in the text "qname422".

Now the string value is 3.1416 and is stored in the global float array gf[0]; it is to be inserted in a text:

```
'Text="Pi is @ approx.";gf[0]'
```

results in the text "Pi is 3.1416 approx."

Finally the name of a unit (“Prisma5”) that’s stored in the global string array gs[11] shall be appended to a message:

```
'Message="Alarm at ";gs[11]'
```

results in the text “Alarm at Prisma5”.

You may set string values by:

cc	Direct input of a numeric value (integer or float).
i[xx]	Element xx of the local integer array.
gi[xx]	Element xx of the global integer array.
f[xx]	Element xx of the local float array.
gf[xx]	Element xx of the global float array.
gfa[xx][yy]	Element xx/yy of the global float field.
gs[xx]	Element xx of the global character string array.

8.4.1.5 Syntax of a sequence

The compiler indents the lines belonging to one instruction by one horizontal tab stop in order to keep the program readable. If an instruction is followed by a block, the latter must be marked by Begin and End:

```
// Test digital input 12 (channel 12) for off
IfDI(12=off)
Begin
  // Wait 5 seconds
  Delay(Time=5)
  // Switch the digital output 8 (channel 8) on
  SetDo(8=on)
End
```

Begin and End themselves are not indented, but only the instructions belonging to the block. Wrong indentations are detected by the compiler and corrected. In the case of an If-Else instruction this can lead to ambiguities. See Chapter 8.4.2.1 If-Else instructions, 419.

NOTE:

The width of the indentation is twice the width of the letter ‘B’. It can be changed in the Parameter Setup program under [Sequence] > [Editor] > [Options] > [Display].

8.4.2 Control structures

8.4.2.1 If-Else instructions

For conditional instructions (**IfDI()**, **IfString()** etc.) the same **If-Else** syntax is always used, where the **Else** portion is optional; the following syntax applies:

```
Ifxxx(condition)
    statement1
```

or

```
Ifxxx(condition)
    statement1
Else
    statement2
```

xxx = instruction type like **DI**, **DO**, **Emi**, etc.

If *condition* is true, *statement1* is executed. Otherwise *statement2* is executed if an **Else** part exists.

statement1 and *statement2* represent individual instructions such as:

```
IfDI(12=off) // Test digital input channel 12 for off state
    Delay(Time=5) // Wait 5 seconds
Else
    Delay(Time=10) // Wait 10 seconds
```

or a block of instructions (See Chapter 8.5.1 Definition of blocks, 426):

```
IfDI(12=off) // Test digital input (channel 12) for off
Begin
    Delay(Time=5) // Wait 5 seconds
    SetDo(8=on) // Switch digital output 8 (channel 8) to on
End
Else
Begin
    Delay(Time=10) // Wait 10 seconds
    SetDo(12=on) // Switch digital output 12 (channel 12) on
End
```

Since the **Else** part of an **If-Else** instruction is optional, a missing **Else** part in nested **If** instructions may cause ambiguity. The compiler associates any **Else** part to the last **If** for which no **Else** part exists. For example in:

```
IfDI(12=off) // Test digital Input channel 12 for off state
    IfDi(14=on) // Test digital Input channel 14 for on state
        Delay(Time=1;min) // Wait 1 minute
    Else
        Delay(Time=10;min) // Wait 10 minutes
```

the **Else** part belongs to the inner **If**, as is implied by the indentation. But if the **Else** part belongs to the outer **If**, **Begin** and **End** statements must be inserted to force correct association:

```

IfDI(12=off) // Test digital input channel 12 for off state
Begin
  IfDi(14=on) // Test digital input channel 14 for on state
  Delay(Time=1;min) // Wait 1 minute
End
Else
  Delay(Time=10;min) // Wait 10 minutes

```

Always set instruction blocks within nested **If-Else** structures between **Begin** and **End** to avoid errors. Indentations are not recognized by the compiler as associations. They are corrected according to the compiler rules, which might cause logical mistakes that are hard to find.

8.4.2.2 Loop instructions

8.4.2.2.1 Loop()

The following syntax is used for defining a **Loop()**:

```

Loop(expr1; expr2; expr3)
  statement

```

- expr1* Is the initialization of the loop variable. The loop variable must be a valid element of the local (**i[xx]**) or the global (**gi[xx]**) integer array. The index *xx* must be within the range of $0 \leq xx \leq \text{Max}$ (Max = configurable maximum).
- expr2* Is the exit condition. When *expr1* reaches *expr2*, the loop is exited.
- expr3* Is the increment that is added to *expr1* after each loop execution before *expr2* is checked; *expr3* is a positive or negative integer within the range $-1'000'000 \leq \text{expr3} \leq +1'000'000$.
- statement* represents an individual instruction such as SetDO() or a block of instructions. See Chapter 8.5.1 Definition of blocks, 426.

NOTE:

The configurable maximum for the loop variable can be set in the Parameter Setup program under [Sequence] > [Editor] > [Options] > [Parameter Limits].

Examples:

```

// Perform loop with i[8]=0 to i[8]=30, Step = 2
// Set digital output 0,2,4,6, ... 30 on
Loop(i[8]=0;30;2)
  SetDO(i[8]=on)

// Perform loop with i[8]=12 to i[8]=0, Step = -1
// Set digital output 12,11,10 ... 0 on
// Wait 10 seconds after each DO setting
Loop(i[8]=12;0;-1)
Begin
  SetDO(i[8]=on)
  Delay(Time=10)

```

End

8.4.2.2.2 While instructions

There are two types of While instructions:

- **WhileDI()** and **WhileDO()** only accept equality comparison to on or off:

```
// Wait while DI channel 12 remains on
WhileDI(12=on)
```

```
// Wait while DO channel 3 remains off
WhileDO(3=off)
```

- **WhileMeaVal()** and **WhileVar()** accept several comparison operators, that is =, <, >, <=, >= and <>.

Additional loop termination conditions can be specified for all **While** loops. The termination criteria '*Count*' and '*Time*' defined with the two keywords **Break** and **BreakTime** can be used for exiting a loop if the main criteria is never reached due to malfunctions or unexpected process conditions, etc.

Break

Break sets a counter that is decremented by one during each loop pass.

BreakTime

BreakTime defines a time after which the loop is exited if no other condition has responded.

If one of these two conditions becomes true, i.e. if the test variable var of **Break** has reached zero or if the time specified with **BreakTime** has run down, the **Break** sequence defined as an option under the keyword **BreakSeq** is started.

Examples:

```
// Wait until the intensity of channel 12 >= 9E-10
WhileMeaVal(12<9E-10,par="d:\qs\par\myfile",Mode=integral)
  Delay(Time=60)
```

If several instructions shall be performed within the **While** loop, they must be defined as a block using **Begin** and **End** (See Chapter 8.5.1 Definition of blocks, 426):

```
// Perform a MID measurement and a calculation until the content
// of element 21 of the global float array >= 45.5
//
WhileVar(gf[21]<45.5)
  Begin
    MID(Par="midpar",SaveGfa=0)
    Calculate(gf[21]=gfa[0][0]*srt(2)/10)
  End

// Purge the gas inlet until element 2 of the local integer array
// has reached a value of 3.
// After 60 passes or max. 60 seconds the While loop shall be
// cancelled and the sequence "c:\qs\par\seq" be executed.
//
WhileVar(i[2]=3,Disp=on,Break=60,BreakTime=60;s,_)
```

```

        BreakSeq=" c:\qs\par\seq" )
Begin
  SetDO( 5=off )
  Delay( Time=5 )
  SetDO( 6=on )
  Delay( Time=1 )
  SetDO( 6=off )
  SetDO( 5=on )
  IfMeaVal( 3>1E-5, par="testpar" )
    SetVar( i[ 2]=0 )
End

```

8.4.3 Log output

In sequencer mode, Quadstar 32-bit allows buffered output of a line-by-line log instead of the page-by-page printing that's normally performed by Windows®. Changeover: See Chapter 3.7.3.5 Log output, 231.

Quadstar 32-bit uses a physical printer port for line-by-line printing. A printer driver has to be installed and connected to a port like shown in the following screenshot:




Fig. 8-3

For detailed information about installing printer drivers, please refer to the Windows® documentation or online-help.

If the printer fails for any reason (paper out, paper jam, power fail etc.), the sequence continues and the lines that could not be printed are stored in a temporary file. After a while the following message (which does not prevent the program execution) appears:




Fig. 8-4

Get the printer ready and switch it online, then choose [Yes] to print the missing lines. If the program has already been terminated, the lines can although be printed by using the Quadstar 32-bit PrintManager:



PrtExe.exe

Fig. 8-5

Double click on the icon (shown above) in the program group Quadstar 32-bit; the Print-Manager (which normally only shows up as a symbol while a sequence is running) displays the following message:




Fig. 8-6

Choose [No] to prevent the stored lines from being deleted. Click again on the icon , then this window will show up:




Fig. 8-7

Get the printer ready, switch it online and choose [File] > [Print/Delete Data Files] C:\TEMP\PRTEXE.DAT. The missing lines will now be printed, presumed they were written to the buffer file.

8.4.4 Data exchange with other programs (DDE)

The Quadstar 32-bit sequencer can exchange data with other programs. For that, the feature DDE (Dynamic Data Exchange) supplied by Windows® is used. It allows two running programs to pass data to each other. The DDE interface of Quadstar 32-bit passes data elements always as strings.

Quadstar 32-bit can act in two ways as a DDE supplier:

8.4.4.1 Passive data exchange

The sequence places its global variables (exclusively those) to other programs' disposal. They are allowed to read and write these variables. The sequence thus supplies a Service, that may be used by other programs.

The service is initialized by **DDEInit(...)**.

Example:

You want to perform a MID measurement of helium 100-times and display the result in an EXCEL sheet.

Quadstar 32-bit sequence:

```
DDEInit (Service="quadserv").
Loop (i[0]=1;100)
  MID (Par="c:\qs32bit\par\helium4.mip", SaveGfa=0)
DDEClose()
```

Entry in the field of the Excel sheet to show the result:

```
=quadserv|Data!'gfa[0][0]'
```

This is the far reference formula of Excel. The entries mean:

quadserv:	Name of the Quadstar 32-bit service, any name possible
	Vertical line
Data	Topic name, in Quadstar 32-bit always 'Data'
!	Exclamation mark
gfa[0][0]	Item, in Quadstar 32-bit always a global variable

You notice, that the value in the Excel sheet changes only, if you tell Excel to update it.

8.4.4.2 Active data exchange

The sequence sends data (global or local variables) to other applications (**DDEPoke(...)**) or requests data from them (**DDERequest(...)**) by its own. It uses the Service that must be supplied by the other application.

Example:

You want to perform a MID measurement of helium 100-times and display the result in an Excel sheet. This time it's not Excel that asks for the values, but Quadstar 32-bit itself who writes them (actively) in row 1, column 1 of the Excel sheet 'Sheet1'. At the end, Quadstar 32-bit asks for the value of the field in Sheet1, row 1, column 2 and saves it in **gi[0]**.

Quadstar 32-bit sequence:

```
Loop(i[0]=1;100)
Begin
  MID(Par="c:\qs32bit\par\helium4.mip", SaveGfa=0)
  DDEPoke(Service="Excel", Topic="Sheet1", Item="r1c1", _
    Data=gfa[0][0])
End
DDERequest(Service="Excel", Topic="Sheet1", Item="r1c2", _
  Data=gi[0])
```

The current values are continuously refreshed in the Excel sheet.

Other properties:

- You may let two or more Quadstar 32-bit or other applications exchange data with each other, but each service must have its own name.
- The properties of the service of an application must be known, if you want to use it (e.g. Excel, see example above).
- A Quadstar 32-bit sequence can provide a service (passively, **DDEInit()**) and send or request data (actively, **DDEPoke()/DDERequest()**) at the same time.
- The sequence can make a program execute a certain instruction by **DDEExecute(...)**. There are programs that allow execution of internal instructions (e.g. macros etc.) by DDE call. Quadstar 32-bit instructions cannot be executed directly by DDE calls.

8.5 Auxiliary instructions

8.5.1 Definition of blocks

The expressions **Begin** and **End** are used to combine declarations and statements into a block. Such a block is handled by the compiler like a single instruction.

For example, you can use **Begin** and **End** after functions such as **IfDI()**, **Loop()**, **WhileDI()** etc. to execute a set of instructions depending on certain conditions.

8.5.1.1 Begin instruction block

Description:

Mark the begin of an instruction block.

Syntax:

Begin

Parameters:

None.

Example(s):

```
Loop(i[5]=1;10)
Begin
  SetDO(5=on)
  Delay(Time=1;min)
End
```

The instructions between **Begin** and **End** are performed 10 times.

8.5.1.2 End instruction block

Description:

Mark the end of an instruction block.

Syntax:

End

Parameters:

None.

Example(s):

```
Loop(i[5]=1;10)
Begin
  SetDO(5=on)
  Delay(Time=1;min)
End
```

The instructions between **Begin** and **End** are performed 10 times.

8.5.1.3 Entering a blank line

Description:

Blank line. You can enter blank lines to keep the program readable; they do not influence the program execution.

Syntax:

Parameters:

None.

Example(s):

```
SetDO(12=on)

Loop(i[3]=1;5)
Begin
  SetDO(14=on)
  Delay(Time=1)
  SetDO(14=off)
  Delay(Time=1)
End

MID(Par="sample")
```

To give the program a structure and thus improve its readability, a blank line has been inserted before and after the loop and its instruction block.

8.5.2 Else instruction

If the condition of an **If** instruction is not true, the instruction following the **Else** statement is performed. The **Else** can also be followed by a block of instructions between **Begin** and **End** statements. If **Else** is true, the entire block of instructions will be processed.

Description:

Covers the "Else" situation of an If instruction.

Syntax:

Else

Parameters:

None.

Example(s):

```
IfVar(gi[1]=45)
  SetDO(5=on)
Else
  SetDO(5=off)
```

If **gi[1]** has a value of 45, open valve 5, else close valve 5.

```
IfVar(f[3]>1E-6)
  SetDO(i[2]=on)
Else
  SetDO("valve3"=off)
```

If **f[3]** is greater than 1E-6, set Digital Output **i[2]** on, else close “valve3”.

8.5.3 Comments

8.5.3.1 Comments comprising multiple lines

All lines between the begin comment (`/*`) and end comment (`*/`) symbols are interpreted as comment. Comment symbols may be nested, i.e. within a comment block there may be another comment block.

8.5.3.1.1 `/*` - Begin comment

Description:
Begin comment

Syntax:
`/*`

Example(s):
See Chapter 8.5.3.1.2 `*/` - End comment, 428

8.5.3.1.2 `*/` - End comment

Description:
End comment. Lines between the begin comment (`/*`) and end comment (`*/`) symbols are treated as comments.

Syntax:
`*/`

Example(s):
`/* Begin of comment`
`SetDO("Turbo"=on)`
`MID(Par="measure")`
`*/End of comment`

These two program lines are temporarily not needed but should later be reactivated. They can be disabled by placing them between comment symbols.

8.5.3.2 `//` - Single line comment

Description:
Single line comment. If `//` is set at the beginning of a line, the line is regarded as a comment.

Syntax:
`//`

Example(s):

```
// SetDO("Turbo"=on)
```

8.6 Example of a sequence

8.6.1 Sequence AIRDEMO

```

-----
01 // Test Sequence <AIRDEMO.SEQ>
02
03 Message(Text="Hello, are you ready ?", Confirm=on)
04 SetPar(BarWidth=15, ColorMode=single, Cycles=50, _
        DispOpt=last, LineType=solid, Marker=on, YRaster=on)
05
06 // GSC(Par="airdemo.gcp", Disp=on, Prot=off)
07 MSC(Par="airdemo.msp", CalMode=coarse, Disp=on, _
        SaveMode=reset)
08 // ZeroGas(Par="airdemo.mip", Disp=on)
09
10 Message(Text="Calibration done, continue ?", Confirm=on)
11
12 Loop(i[1]=0;1)
13 Begin
14     Sequence(Par="didodem1.seq")
15     Sequence(Par="didodem2.seq")
16
17     Loop(i[0]=0;20)
18     Begin
19         MID(Par="airdemo.mip", Disp=tab, SaveGfa=0)
20         SetString(gs[0] = "MID measure cycle @, intensity of_
                channel 0 = @ A";i[0])
21 // MCD(Par="airdemo.mcp", Disp=tab, SaveGfa=0)
22 // SetString(gs[0] = "MCD measure cycle @, concentration N2_
                = @ %";i[0])
23     Message(Text=gs[0];gfa[0][0])
24     End
25
26     Loop(i[0]=0;20)
27     Begin
28         MID(Par="airdemo.mip", Disp=bar, SaveGfa=0)
29         SetString(gs[0] = "MID measure cycle @, intensity of_
                channel 0 = @ A";i[0])
30 // MCD(Par="airdemo.mcp", Disp=bar, SaveGfa=0)
31 // SetString(gs[0] = "MCD measure cycle @, concentration N2_
                = @ %";i[0])
32     Message(Text=gs[0];gfa[0][0])
33     End
34
35     Loop(i[0]=0;75)
36     Begin
37         MID(Par="airdemo.mip", Disp=vt, SaveGfa=0)
38         SetString(gs[0] = "MID measure cycle @, intensity of_
                channel 0 = @ A";i[0])
39 // MCD(Par="airdemo.mcp", Disp=vt, SaveGfa=0)
40 // SetString(gs[0] = "MCD measure cycle @, concentration N2_
                = @ %";i[0])
41     Message(Text=gs[0];gfa[0][0])
42     End
43 End

```

NOTE:

The MCD measurements have been flagged with the comments symbol //.

8.6.2 User-controlled analysis (AIRDEMO2)

8.6.2.1 Quadstar 32-bit dialog box

Using the **Dialog()** instruction, you can program a user-controlled analysis.

The **Dialog()** instruction produces a Quadstar 32-bit dialog box that enables the operator to control the execution of an analysis. This control is supplied by Buttons; a return value can be assigned to each of the command buttons.

The dialog boxes produced with the **Dialog()** instruction offer the following functions:

- Labeled command buttons (keyword **Button**).
- Text output within the dialog box (keyword **Text**).
- Input of numeric values and texts (keyword **Edit**).
- Choose an entry from a list field (keyword **Combo**).
- Input a file name or choose it by using the file manager (keyword **Filename**).

Within the **Dialog()** instruction a return variable (keyword **ReturnVar**) is defined, that contains the return value:

- the assigned value **rval**, when one of the command Buttons is pressed.
- the value **EnterVal**, when the Enter key is pressed.
- the value **EscapeVal**, when the Esc key is pressed.

8.6.2.2 Sequence AIRDEMO2

```

01 // Test Sequence <AIRDEMO2.SEQ>
02 // Test Sequence <AIRDEMO2.SEQ>
03
04 Message(Text="Hello, are you ready ?", Confirm=on)
05
06 // Set Quadstar 32-bit parameter for sequence measure data_
    display
07 SetPar(AlarmDisp=on, BarWidth=25, ColorMode=single, _
    Cycles=50, DispOpt=last, LineType=solid, Marker=on, _
    Scale=log, YRaster=off)
08
09 // Mass Calibration
10 Dialog(Box="Mass Calibration Save Mode ";200;200;370;140, _
    ReturnVar=i[0], EnterVal=0, EscapeVal=-1, _
    Button="&Append";20;20;150;40;1, _
    Button="&Reset";190;20;150;40;2, _
    Button="&No Mass Calibration";50;80;260;40;3)
11
12 // Mass Calibration Append

```

```

13 IfVar(i[0]=1)
14 Begin
15   MSC(Par="airdemo.msp", CalMode=coarse, Disp=on)
16   Message(Text="Calibration done, continue ?", Confirm=on)
17 End
18
19 // Mass Calibration Reset
20 IfVar(i[0]=2)
21 Begin
22   MSC(Par="airdemo.msp", CalMode=coarse, Disp=on,
        SaveMode=reset)
23   Message(Text="Calibration done, continue ?", Confirm=on)
24 End
25
26 WhileVar(i[0]>=0, Disp=off)
27 Begin
28 // Select Measure Mode (MID, Scan Analog or Scan Bargraph)
    Dialog(Box="Selection Menu"; 200; 200; 460; 420,
           ReturnVar=i[0], EnterVal=0, EscapeVal=-1,
           Text="Quadstar 32-bit Gas Inlet_
Control";90;25;280;20,
           Button("&Open Inlet";20;60;200;60;1,
           Button("&Close Inlet";240;60;200;60;2,
           Text="Quadstar 32-bit Measurement Mode";80;145;300;20,
           Button("&Multiple Ion Detection (MID)";70;180;320;60;3,
           Button("Scan &Analog";20;260;200;60;6,
           Button("Scan &Bargraph";240;260;200;60;7,
           Button("C&ancel";40;340;380;60;-1)
29
30 // Select MID Measure Mode (Table, Bargraph or Versus Time)
31 IfVar(i[0]=3)
32   Dialog(Box="MID Measure Mode Select";200;200;240;200,
           ReturnVar=i[0], EnterVal=0, EscapeVal=-1,
           Button="MID &Table";20;20;200;40;3,
           Button="MID &Bargraph";20;80;200;40;4,
           Button="MID &Intensity vs. Time";20;140;200;40;5)
33
34 // Enter number of Cycles
35 IfVar(i[0]>=3)
36   Dialog(Box="Cycles";200;200;390;60, ReturnVar=i[1],
           EnterVal=0, EscapeVal=-1,
           Text="Enter number of cycles 0...500 :";20;20;250;20,
           Edit=280;15;100;30;i[6])
37
38 IfVar(i[0]=0)
39   Message(Text="Please choose a button !")
40
41
42 // Open Inlet
43 IfVar(i[0]=1)
44   Sequence(Par="didodem1.seq")
45
46 // Close Inlet
47 IfVar(i[0]=2)
48   Sequence(Par="didodem2.seq")
49
50
51 // MID Table

```



```

52   IfVar(i[0]=3)
53   Begin
54     Loop(i[5]=1;i[6])
55     Begin
56       MID(Par="airdemo.mip", Disp=tab, Prot=off, SaveGfa=0)
57       SetString(gs[0] = "MID measure cycle @, _
                    ion current of channel 0 = @ A";i[5])
58       Message(Text=gs[0];gfa[0][0])
59     End
60   End
61
62 // MID Bargraph
63   IfVar(i[0]=4)
64   Begin
65     Loop(i[5]=1;i[6])
66     Begin
67       MID(Par="airdemo.mip", Disp=bar, Prot=off, SaveGfa=0)
68       SetString(gs[0] = "MID measure cycle @, _
                    ion current of channel 0 = @ A";i[5])
69       Message(Text=gs[0];gfa[0][0])
70     End
71   End
72
73 // MID Versus Time
74   IfVar(i[0]=5)
75   Begin
76     Loop(i[5]=1;i[6])
77     Begin
78       MID(Par="airdemo.mip", Disp=vt, Prot=off, SaveGfa=0)
79       SetString(gs[0] = "MID measure cycle @, _
                    ion current of channel 0 = @ A";i[5])
80       Message(Text=gs[0];gfa[0][0])
81     End
82   End
83
84
85 // Scan Analog
86   IfVar(i[0]=6)
87   Begin
88     Loop(i[5]=1;i[6])
89     Begin
90       Message(Text="Scan Measure, Scan Analog # @";i[5])
91       ScanAnalog(Par="airdemo.sap", Disp=on)
92     End
93   End
94
95
96 // Scan Bargraph
97   IfVar(i[0]=7)
98   Begin
99     Loop(i[5]=1;5)
100    Begin
101      Message(Text="Scan Measure, Scan Bargraph # @";i[5])
102      ScanBar(Par="airdemo.sbp", Disp=on, Prot=off)
103    End
104  End
105
106 End
-----

```

8.6.2.2.1 Description of AIRDEMO2

The example AIRDEMO2.SEQ demonstrates partially nested **Dialog()** instructions. In all Quadstar 32-bit dialog boxes of AIRDEMO2, the variable **i[0]** is used as the return value.

With the first dialog box (Mass Calibration Save Mode), the user is requested to perform a mass calibration. A choice between the Append and Reset mode is offered.

```
10 Dialog(Box="Mass Calibration Save Mode ";200;200;370;140, _
    ReturnVar=i[0], EnterVal=0, EscapeVal=-1, _
    Button="&Append";20;20;150;40;1, _
    Button="&Reset";190;20;150;40;2, _
    Button="&No Mass Calibration";50;80;260;40;3)
```

This **Dialog()** instruction produces the following Quadstar 32-bit dialog box:




Fig. 8-8

The return variable can be loaded with the following values:

- -1: Esc key has been pressed, cancel the sequence.
- 0: Enter key has been pressed
- 1: Mass calibration Append
- 2: Mass calibration Reset
- 3: No mass calibration

The input is evaluated immediately:

```
12 // Mass Calibration Append
13 IfVar(i[0]=1)
14 Begin
15   MSC(Par="airdemo.msp", CalMode=coarse, Disp=on)
16   Message(Text="Calibration done, continue ?", Confirm=on)
17 End
18
19 // Mass Calibration Reset
20 IfVar(i[0]=2)
21 Begin
22   MSC(Par="airdemo.msp", CalMode=coarse, Disp=on, _
    SaveMode=reset)
23   Message(Text="Calibration done, continue ?", Confirm=on)
24 End
```

An analysis can be controlled via the return variable by systematically using **IfVar()** and **WhileVar()** instructions. For canceling via the Esc key ($i[0] = -1$) the **WhileVar()** instruction can, for example, be used as follows:

```
WhileVar(i[0]>=0, Disp=off)
```

```

Begin
// Select Measure Mode (MID, Scan Analog or Scan Bargraph)
Dialog(Box="Selection Menu";200;200;460;420, ReturnVar=i[0],_
EnterVal=0, EscapeVal=-1, Button="Open ...
...
...
...
End

```

In the individual dialog boxes, the variable can be loaded with the following values:

Value	Quadstar 32-bit dialog box "Selection Menu"	Quadstar 32-bit dialog box "MID Mode Selection Menu"
-1	Esc key has been pressed	Esc key has been pressed
0	Enter key has been pressed	Enter key has been pressed
1	Open gas inlet	----
2	Close gas inlet	----
3	Multiple Ion Detection MID	MID Table
4	----	MID Bargraph
5	----	MID Versus Time
6	Scan Analog	----
7	Scan Bargraph	----

First the Quadstar 32-bit dialog box Selection Menu is called. The array element i[0] contains a return value of 0 (Enter key), -1 (Esc key or [Cancel] command button) depending on the selected command button 1, 2, 3, 6 or 7:

```

28 // Select Measure Mode (MID, Scan Analog or Scan Bargraph)
Dialog(Box="Selection Menu"; 200; 200; 460; 420, _
ReturnVar=i[0], EnterVal=0, EscapeVal=-1, _
Text="Quadstar 32-bit Gas Inlet_
Control";90;25;280;20, _
Button="&Open Inlet";20;60;200;60;1, _
Button="&Close Inlet";240;60;200;60;2, _
Text="Quadstar 32-bit Measurement Mode";80;145;300;20, _
Button="&Multiple Ion Detection (MID)";70;180;320;60;3, _
Button="Scan &Analog";20;260;200;60;6, _
Button="Scan &Bargraph";240;260;200;60;7, _
Button="C&ancel";40;340;380;60;-1)

```

This **Dialog()** instruction produces the following Quadstar 32-bit dialog box:




Fig. 8-9

If Multiple Ion Detection (MID) is chosen, **i[0]** contains a value of 3. Through the decision **IfVar (i[0]=3)** the next Quadstar 32-bit dialog box MID Measure Mode Select is offered for choosing the MID measurement mode.

The return value is again assigned to the variable **i[0]**, i.e. **i[0]** contains a return value of 3, 4 or 5 (or 0, or -1), depending on the selected command button:

```

30 // Select MID Measure Mode (Table, Bargraph or Versus Time)
31 IfVar(i[0]=3)
32   Dialog(Box="MID Measure Mode Select";200;200;240;200, _
         ReturnVar=i[0], EnterVal=0, EscapeVal=-1, _
         Button="MID &Table";20;20;200;40;3, _
         Button="MID &Bargraph";20;80;200;40;4, _
         Button="MID &Intensity vs. Time";20;140;200;40;5)

```

This **Dialog()** instruction produces the following Quadstar 32-bit dialog box:




Fig. 8-10

In the last Quadstar 32-bit dialog box, the operator is requested to enter the number of measurement cycles. The entered value (number of cycles) is returned via the variable **i[6]**:

```

34 // Enter number of Cycles

```

```

35   IfVar(i[0]>=3)
36   Dialog(Box="Cycles";200;200;390;60, ReturnVar=i[1],_
        EnterVal=0, EscapeVal=-1,_
        Text="Enter number of cycles 0...500 :";20;20;250;20,_
        Edit=280;15;100;30;i[6])

```

This **Dialog()** instruction produces the following Quadstar 32-bit dialog box:



Fig. 8-11

After this input the desired measurement is started. The variable **i[5]** is used as the loop counter.

The following example illustrates the course of the sequence if e.g. Intensity vs Time has been selected:

```

...
...
26  WhileVar(i[0]>=0, Disp=off)
27  Begin
...
...
73  // MID Versus Time
74  IfVar(i[0]=5)
75  Begin
76  Loop(i[5]=1;i[6])
77  Begin
78  MID(Par="airdemo.mip", Disp=vt, Prot=off, SaveGfa=0)
79  SetString(gs[0] = "MID measure cycle @, _
        ion current of channel 0 = @ A";i[5])
80  Message(Text=gs[0];gfa[0][0])
81  End
82  End
...
...
106 End

```

The execution of this loop produces the following measurement data window:




Fig. 8-12

8.7 Sequencer instructions

8.7.1 Calculate() - Perform a calculation

Description:

Perform a calculation.

Calculate() determines the result of the formula defined under formula and assigns it to the variable var. The content of this variable can later be stored in a file via the **Process()** instruction and/or displayed on the screen.

Through the options Name and Unit, a name and/or unit of measure can be assigned to the variable (var). The name is subsequently used for identifying the block:

- when the data are processed by the **Process()** instruction.
- when stored data are shown in the Display Saved Values program.

An individual scale may be assigned to the variable as well.

NOTE:

Of the data calculated with the Calculate() instruction and stored in the gfa[][] data field, only those are taken into consideration in the Process() instruction to which names have been assigned. All other data of the Calculate() instruction are regarded as intermediate results.

Syntax:

```
Calculate( [ScaleID=scID, ] [Name=name[;strval], ] [Unit=unit, ]_
          var=formula)
```

Parameters:

ScaleID	Keyword for scale identification.
<i>scID</i>	Number of the scale (defined by the instruction 'SetScale'). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (0 ... 199). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Name	Keyword for variable name.
<i>name</i>	Name of the variable (max. 12 characters) <ul style="list-style-type: none"> • gs[xx] - Element xx of the global string array. • "xxxxxx" - Direct input of the variable name.
<i>strval</i>	String value - See Chapter 8.4.1.4 String values, 417
Unit	Keyword for unit of measure.
<i>unit</i>	unit of measure that fits the result (max. 12 characters) <ul style="list-style-type: none"> • gs[xx] - Element xx of the global string array. • "xxxxxx" - Direct input of the unit of measure.

<i>var</i>	<p>Variable to assign the result of the calculation to</p> <ul style="list-style-type: none"> • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.
<i>formula</i>	<p>Input of the calculation formula. This formula can comprise up to 128 characters. Nesting of function calls is permitted (e.g. <code>srt(abs(gf[3]))</code>). All float and integer array element types can be used:</p> <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer or float). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field. <p>In an index specification only simple nesting is permitted. In addition to the direct specification (e.g. i[12], f[10]) the index itself may be an element of an integer array, i.e. also expressions such as i[gi[10]], gi[gi[15]], etc. may be entered.</p> <p>The following functions and operations are possible:</p> <ul style="list-style-type: none"> • +, -, *, / - Basic operations. • () - Parenthesized expressions. • exn(x) - Calculation of e to the power of x • exp(x) - Calculation of 10 to the power of x. • lgn() - Calculation of the logarithm, base e. • log() - Calculation of the logarithm, base 10. • sqr() - Squaring. • srt() - Calculation of the square root. • abs() - Determination of the absolute value. <p>These functions may be nested, i.e. other functions may be specified within parentheses.</p> <p>Example:</p> <pre>log(srt(gfa[12][i[4]]*12.5)-3.2).</pre> <p>Input of blanks is admissible but they are removed during the syntax check of the formula.</p> <p>Integer fractions are rounded arithmetically:</p> <pre>0.49 => 0 0.51 => 1.</pre>

Example(s):

```
Calculate(i[3]=(i[4]+gi[5]-2)*3)
```

Simple calculation with basic operations and assignment of the result to the variable **i[3]**.

```
Calculate(Name="Temperature";i[8],Unit="°C", gf[24]=(f[4]/i[0]_
+log(srt(gfa[12][i[4]]*12.5)-3.2)+sqr(f[8]))/2)
```

A more complicated calculation with functions and nested expressions. The result is assigned to element 24 of the global float array **gf[.]**. This variable (**gf[24]**) is given the name TemperatureX before and the unit °C is assigned to it, where X is the content of **i[8]** (e.g. Temperature9).

8.7.2 DDEClose() - Close DDE Service

Description:

Close the DDE service that has been initialized in the sequence beforehand.

Syntax:

```
DDEClose()
```

Parameters:

None

Example:

See Chapter 8.4.4 Data exchange with other programs (DDE), 424.

8.7.3 DDEExecute() - Call an instruction in an application by DDE

Description:

Call an instruction in another program by DDE.

See Chapter 8.4.4 Data exchange with other programs (DDE), 424.

Syntax:

```
DDEExecute(Service=ServName, Topic=TopicName, Data=DatString)
```

Parameters:

Service	Keyword for service name.
<i>ServName</i>	Name of the service. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the name. • gs[xx] - Element <i>xx</i> of the global string array, if the name of the service is stored in gs[].
Topic	Keyword for topic name.
<i>TopicName</i>	Name of the topic. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the name. • gs[xx] - Element <i>xx</i> of the global string array, if the name of the topic is stored in gs[].
Data	Keyword for data.
<i>DatString</i>	Data that are additionally passed to the application as a string. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the data. • gs[xx] - Element <i>xx</i> of the global string array, if the data are stored in gs[].

8.7.4 DDEInit() - Initialize DDE Service

Description:

Initialize DDE Service.

Syntax:

```
DDEInit( Service=ServName )
```

Parameters:

Service	Keyword for service name.
ServName	Name of the service. <ul style="list-style-type: none"> • "nnn" - Direct input of the name. • gs[xx] - Element xx of the global string array, if the name of the service is stored in gs[].

Example:

See Chapter 8.4.4 Data exchange with other programs (DDE), 424.

8.7.5 DDEPoke() - Send data by DDE

Description:

Send data by DDE.

Data that exist as global or local variables are passed to another running program.

Syntax:

```
DDEPoke( Service=ServName, Topic=TopicName, Item=ItemName, _
        Data=Dat )
```

Parameters:

Service	Keyword for service name.
<i>ServName</i>	Name of the service. <ul style="list-style-type: none"> • "nnn" - Direct input of the name. • gs[xx] - Element xx of the global string array, if the name of the service is stored in gs[].
Topic	Keyword for topic name.
<i>TopicName</i>	Name of the topic. <ul style="list-style-type: none"> • "nnn" - Direct input of the name. • gs[xx] - Element xx of the global string array, if the name of the topic is stored in gs[].
Item	Keyword for item name.
<i>ItemName</i>	Name of the item. <ul style="list-style-type: none"> • "nnn" - Direct input of the name. • gs[xx] - Element xx of the global string array, if the name of the item is stored in gs[].

Data	Keyword for data.
Dat	Data that is passed to the program. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the data. • i[<i>xx</i>] - Element <i>xx</i> of the local integer array. • gi[<i>xx</i>] - Element <i>xx</i> of the global integer array. • f[<i>xx</i>] - Element <i>xx</i> of the local float array. • gf[<i>xx</i>] - Element <i>xx</i> of the global float array. • gfa[<i>xx</i>][<i>yy</i>] - Element <i>xx/yy</i> of the global float field. • gs[<i>xx</i>] - Element <i>xx</i> of the global string array, if the data are stored in gs[].

Example:

See Chapter 8.4.4 Data exchange with other programs (DDE), 424.

8.7.6 DDERequest() - Request data by DDE

Description:

Request data by DDE.

Request data from another program and store it in global or local variables.

Syntax:

```
DDERequest( Service=ServName, Topic=TopicName, Item=ItemName, _
            Data=Dat )
```

Parameters:

Service	Keyword for service name.
<i>ServName</i>	Name of the service. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the name. • gs[<i>xx</i>] Element <i>xx</i> of the global string array, if the name of the service is stored in gs[].
Topic	Keyword for topic name.
<i>TopicName</i>	Name of the topic. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the name. • gs[<i>xx</i>] Element <i>xx</i> of the global string array, if the name of the topic is stored in gs[].
Item	Keyword for item name.
<i>ItemName</i>	Name of the item. <ul style="list-style-type: none"> • “<i>nnn</i>” - Direct input of the name. • gs[<i>xx</i>] Element <i>xx</i> of the global string array, if the name of the topic is stored in gs[].

Data	Keyword for data.
Dat	Variable to assign the data to that is passed by the program. <ul style="list-style-type: none"> • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field. • gs[xx] - Element xx of the global string array, if the data are stored in gs[].

Example:

See Chapter 8.4.4 Data exchange with other programs (DDE), 424.

8.7.7 Delay() - Time delay

Description:

Delay the execution of the sequence by a certain amount of time.

The time is specified by the value *tim*. If unit is not specified, the time is interpreted as seconds.

Syntax:

```
Delay(Time=tim[;unit][,Disp=d])
```

Parameters:

Time	Keyword for time specification.
<i>tim</i>	Time specification (1 ... 1'000'000). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
unit	Unit of time <ul style="list-style-type: none"> • ms - in milliseconds. • [s] - in seconds (default). • min - in minutes. • h - in hours. • d - in days.
Disp	Keyword for display mode. If Disp=on is set, the remaining delay time will be displayed in the status bar.
<i>d</i>	Display mode <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p>

Example(s):

```
Delay(Time=i[77])
```

Delay the sequence by the number of seconds stored in *i[77]*.

```
Delay(Time=10;min,Disp=on)
```

Delay the sequence by 10 minutes and display the remaining time.

8.7.8 Dialog() - Displaying a dialog box

Description:

The **Dialog()** instruction displays a programmer-defined dialog window when the sequence is executed (or when the Test button is pressed). A background picture may be displayed if necessary. This way, the operator can control the sequence.

To interpret the input, a variable can be defined as ReturnVar which returns

- the value of **EnterVal** if the dialog box is closed by pressing 'Enter'
- the value of **EscapeVal** if the dialog box is closed by pressing 'Esc'
- the corresponding **ReturnVal** if any command button has been pressed.

NOTE:

The elements labeled with the keywords **Text**, **Button**, **Edit**, **Combo** and **Filename** can be listed included as often as desired in any order, but the line length cannot exceed 1000 characters.

Syntax:

```
Dialog(Box=dtext;dx;dy;dxw;dyw,ReturnVar=rvar,EnterVal=ent,_,_
      EscapeVal=esc_
      [,Text=tttext;tx;ty;txw;tyw]_
      [,Button=bttext;bx;by;bxw;byw;rval]_
      [,Edit=ex;ey;exw;eyw;eret]_
      [,Combo=cfirst;cx;cy;cxw;cyw;cret]_
      [,Filename=ftype;fx;fy;fxw;fyw;fret]_
      [,Picture=BkGndPict[,XPos=x][,YPos=y][,Height=h]_
      [,Width=w]]
      ...
)
```

Parameters:

Box	Keyword for dialog box.
<i>dtext</i>	Name of dialog box <ul style="list-style-type: none"> • "ttttttt" Text for name • gs[xx] Element xx of the global character string array if the window name is stored in gs[].
<i>dx, dy</i>	Position X and Y of the dialog box, relative to the coordinate origin, i.e. the upper left corner of the screen.
<i>dxw, dyw</i>	Extent of the dialog box in the X and Y direction. <ul style="list-style-type: none"> • ccc Numeric value 0...screen resolution
ReturnVar	Keyword for defining the return variable.
<i>rvar</i>	Definition of the array element. <ul style="list-style-type: none"> • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
EnterVal	Keyword for defining the return value for the Enter key.

<i>ent</i>	value for 'Enter' <ul style="list-style-type: none"> • <i>cc</i> - Numeric value (integer, -1'000'000 ... +1'000'000) • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
EscapeVal	Keyword for defining the return value for the Esc key.
<i>esc</i>	value for 'Esc' <ul style="list-style-type: none"> • <i>cc</i> - Numeric value (integer, -1'000'000 ... +1'000'000) • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
Text	Keyword for text.
<i>ttext</i>	Text to display <ul style="list-style-type: none"> • "ttttttt" Text • gs[xx] - Element <i>xx</i> of the global character string array if the text is stored in gs[].
<i>tx, ty</i>	Position X and Y of the text, relative to the upper left corner of the dialog box.
<i>txw, tyw</i>	Extent of the text field in the X and Y direction. <ul style="list-style-type: none"> • <i>ccc</i> - Numeric value 0 ... screen resolution.
Button	Keyword for a command button.
<i>btext</i>	Name of the command button <ul style="list-style-type: none"> • "&ttttttt" - Direct input of the name => The & marks the character to be underscored (shortcut). • gs[xx] - Element <i>xx</i> of the global character string array if the name of the command button is stored in gs[].
<i>bx, by</i>	Position X and Y of the command button, relative to the upper left corner of the dialog box.
<i>bxw, byw</i>	Extent of the command button in the X and Y direction. <ul style="list-style-type: none"> • <i>ccc</i> - Numeric value 0 ... screen resolution
<i>rval</i>	Value that is returned when the command button is pressed. <ul style="list-style-type: none"> • <i>cc</i> - Numeric value (integer, -1'000'000 ... +1'000'000) • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
Edit	Keyword for input (numeric value or text)
<i>ex, ey</i>	Position X and Y of the command button, relative to the upper left corner of the dialog box.
<i>exw, eyw</i>	Extent of the command button in the X and Y direction. <ul style="list-style-type: none"> • <i>ccc</i> - Numeric value 0 ... screen resolution.
<i>eret</i>	Return value (numeric value or text) <ul style="list-style-type: none"> • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] -Element <i>xx</i> of the global integer array. • f[xx] - Element <i>xx</i> of the local float array. • gf[xx] - Element <i>xx</i> of the global float array. • gs[xx] - Element <i>xx</i> of the global character string array in which the text is to be stored.
Combo	Keyword for list field (Combobox)

<i>cfirst</i>	String-Position <i>n</i> of the first list entry. The list begins at <i>gs[n]</i> and ends where the first empty gs[] is found. <ul style="list-style-type: none"> • <i>cc</i> - Numeric value • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
<i>cx, cy</i>	Position X and Y of the list field, relative to the upper left corner of the dialog box.
<i>cxw, cyw</i>	Extent of the list field in the X and Y direction. <ul style="list-style-type: none"> • <i>ccc</i> - Numeric value 0 ... screen resolution.
<i>cret</i>	Return value (numeric value or text) <ul style="list-style-type: none"> • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array. • f[xx] - Element <i>xx</i> of the local float array. • gf[xx] - Element <i>xx</i> of the global float array. • gs[xx] - Element <i>xx</i> of the global character string array in which the text is to be stored. The return value is either the number of the list entry (0, 1, .. for i[], gi[], f[], gf[]), or the chosen string (for gs[]).
Filename	Keyword for entering a file name.
<i>ftype</i>	Definition of the file type through the file extension.
<i>fx, fy</i>	Position X and Y of the command button, relative to the coordinate origin, that is, the upper left corner of the dialog box.
<i>fxw, fyw</i>	Command button size in the X and Y direction. <ul style="list-style-type: none"> • <i>ccc</i> - Numeric value 0 ... screen resolution
<i>fret</i>	Text return value (file name) <ul style="list-style-type: none"> • gs[xx] - Element <i>xx</i> of the global character array to be stored in the file name.
Picture	Keyword for background graphic.
<i>BkGndPict</i>	Name and path of the background graphic <ul style="list-style-type: none"> • "Pict.bmp" - Direct input of a name. • gs[xx] - Element <i>xx</i> of the global character string array.
XPos, YPos, Height, Width	Keywords for position and size of the background graphic. The position of the background graphic relative to the upper left corner of the dialog window can be defined by the values XPos and YPos, the size by Height and Width.
<i>x, y, h, w</i>	X- and Y- position of the graphic, and height and width. <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.

For the elements **Edit**, **Combo** and **Filename** you can set default values by preloading each of their return variables with the corresponding value. These default values are loaded and displayed when the dialog box is created.

Example(s):

See Chapter 8.6.2.2 Sequence AIRDEMO2, 431

```
Dialog(Box="Select Measure Mode";200;200;240;200, _
      ReturnVar= i[1], EnterVal=0, _
```

```

        EscapeVal=-1, Button="&Bargraph";20;20;200; 40;1, _
        Button="&Table";20;80;200;40;2, _
        Button="&Intensity Versus Time";20;140;200;40;3)
IfVar(i[1]=1)
Begin
    Loop(i[2]=1;20)
        MID(Par="d:\qs32bit\par\airdemo.mip",Disp=vt,Prot=off)
    End
...

```

8.7.9 LoopDialog() - Loop with Dialog window

Description:

LoopDialog() is a loop that starts with a dialog window. A block of sequencer instructions follows (delimited by **Begin** and **End**), that's executed repeatedly. This execution can be interrupted, the dialog window will be shown then again.

The LoopDialog loop is terminated when the return value RetVar corresponds to the value of EscapeVal.

If the key combination <Shift><F4> is pressed within the loop, the dialog window shows up immediately. This way the operator can e.g. interrupt a measurement and select a new action in the dialog (the second entry in the File menu (text=title of the dialog box) does the same as <Shift><F4>).

Syntax:

```
LoopDialog(...)
```

Parameters:

See Chapter 8.7.8 Dialog() - Displaying a dialog box, 445.

Example(s):

```

LoopDialog(Box="Selection";200;200;240;200, _
        ReturnVar=i[0], EnterVal=0, EscapeVal=-1, _
        Button="Scan Analog";20;70;200;30,1, _
...
        Button="MCD Bargraph";20;270;200;30,5)
Begin
    ...
    ...
End

```

Display a dialog box. When the <ESC> key is pressed within the dialog box, the program branches to the sequencer instruction that follows the End instruction.

8.7.10 DispEventWindow() - Place and scale the Event window

Description:

Place the event window somewhere on the screen and enlarge or reduce its scale.

By using this instruction you can place the event window of the running sequence in a chosen size at a chosen position on your screen. This instruction is equal to the instruction DispWindow, except that it controls the event window instead of the measure window.

Syntax:

```
DispEventWindow( . . . )
```

Parameters:

See Chapter 8.7.11 DispWindow() - Place and scale the window, 449.

Example(s):

See Chapter 8.7.11 DispWindow() - Place and scale the window, 449.

8.7.11**DispWindow() - Place and scale the window****Description:**

Place the sequencer's measure window somewhere on the screen and enlarge or reduce its scale.

By using this instruction you can place the measure window of the running sequence in a chosen size at a chosen position on your screen. This is especially helpful, if there are several sequences running at the same time.

Syntax:

```
DispWindow(Disp=dispmode [ , Xpos=x ] [ , Ypos=y ] [ , Width=w ]
           [ , Height=h ] )_
```

Parameter:

Disp Keyword for the display mode of the window (normal or iconic).

dispmode Text

- normal - Normal display.
- symbol - Window reduced to symbol (iconic size).
- top - Window is always visible.

x, y, w, h Horizontal and vertical coordinates that indicate the position of the upper left corner of the measure window relative to the upper left corner of the screen, where *x* is the offset to the right and *y* is the offset downwards. *w* is the width and *h* is the height of the measure window. The values *x, y, w, h* are meant as units of the screen resolution and are corrected automatically if they are out of range. The following entries are possible:

- *cc* - Numeric value (positive integer)
- *i[xx]* - Element *xx* of the local integer array.
- *gi[xx]* - Element *xx* of the global integer array.

Example(s):

Your PC works with a graphics resolution of 800 × 600 Pixels (X × Y).

```
DispWindow(Disp=normal, Xpos=200, Ypos=150, Width=400, _
           Height=300)
```

Place the window on the center of the screen and scale it to cover a quarter of the total screen area:

```
DispWindow(Disp=symbol)
```

Reduce the window to iconic size (symbol). The position is determined by Windows® in this case:

8.7.12 Error() - Outputting an error message

Description:

Output a user-specific error message.

NOTE:

Error messages are normally output on the screen.

The optional parameter Save provides storage in the error history file.

The optional parameter Prot provides "On-line" logging of the error messages on the connected printer.

Syntax:

```
Error(Text=text[;strval][,Save=s][,Prot=p][,Confirm=c])
```

Parameters:

Text Keyword for entering the text (error message) to be output.

text Text

- "Text" - Direct input of Text.
- **gs[xx]** - Element xx of the global character string array.

strval String value - See Chapter 8.4.1.4 String values, 417

Save Keyword for save mode.

The error messages can be stored in the error history file. The default setting for the save mode is **on**.

NOTE:

This will work only, if in Parameter Setup program [Setup] > [General] > [Error History] is switched on.

s Save mode

- **off** - Save function off.
- **[on]** - Save function on (default).

Prot Keyword for log mode.

p Log mode.

- **[default]** - Default log mode.
- **off** - Log function off.
- **on** - Printer log on.

See Chapter 8.4.3 Log output, 422.

If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.

Confirm Keyword for the confirmation mode. When an error message is output, user confirmation can be requested with Confirm.

c Confirmation mode

- **[off]** - Sequence proceeds without confirmation.
- **on** - Sequence is stopped until the [OK] button is pressed.

Example(s):

```
Error(Text="Calibration of bottle 5 in progress")
```

Output of an error.

```
Error(Text=gs[12], Save=on, Prot=on, Confirm=on)
```

Display an error, log it on the printer and store it in the error history, then wait until the operator confirms continuation of the sequence.

8.7.13 Execute() - Execute a program

Description:

Execute any program. You can e.g. start a second Measure, a batch file or an application program such as a text editor or similar.

Syntax:

```
Execute(Prog=name[;ind][,Par=p])
```

Parameters:

Prog	Keyword for program name and path.
<i>name</i>	(drive:\path\) program name <ul style="list-style-type: none"> • "pppppp" Program name (e.g. "c:\qs\myprog"). • gs[xx] Element xx of the global string array, if the program name is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Par	Keyword for parameter transfer to the program.
<i>p</i>	Any parameters may be transferred to the program; e.g. the name of a parameter file. <ul style="list-style-type: none"> • [off] no parameter is transferred. • "pppppp" transfer a string (e.g. parameterfile name "c:\qs\myfile"). • gs[xx] Element xx of the global string array, if the parameter to be transferred is stored in gs[].

Example(s):

```
Execute(Prog="myprog.exe")
```

The program 'myprog.exe' is started from the current directory.

```
Execute(Prog="c:\auxis\batsav.bat";3,Par=gs[0])
```

The batch file 'batsav03.bat' is started from the path 'c:\auxis' and the content of the global string gs[0] is transferred as a parameter to it.

8.7.14 Exit() - Quitting the sequencer

Description:

Execute any **SeqExit()** instructions and terminate the sequence. (= > See also Chapter 8.7.48 SeqExit() - Defining a sequence for sequence cancellation, 488). Terminate or don't terminate the program Measure that executes the sequence.

Syntax:

```
Exit(Terminate=OnOff)
```

Parameters:

Terminate Keyword 'terminate measure program automatically Yes/No'.

OnOff

- **[off]** Exit sequence, don't terminate measure program.
- **on** Exit sequence, terminate measure program.

8.7.15**ExtractString() - Extracting characters from a string****Description:**

Extract a part of a character string and assign it to another variable or string.

ExtractString() is especially useful to separate the different informations of a string that has been assembled and transmitted by an external application. As separation criterion, a number of characters or a special character can be used. In the latter case, the separation character itself is deleted after the separation.

Syntax:

```
ExtractString( Type=sep[ ,Number=nbr][ ,Character=char] , _
              Source=src ,Dest=dst )
```

Parameters:

Type Type of separation criterion.

sep

- Number - Get the first <number> characters.
- Character - Get all characters until a certain character is found.

Number Number of characters to extract.

nbr

- *cc* - Numeric value (integer, 0,....,<stringlength>)
- **i[xx]** - Element xx of the local integer array.
- **gi[xx]** - Element xx of the global integer array.

Character Character that indicates the end of the string part to be extracted.

char

- "c" - Direct input of the separation character.
- **gs[xx]** - Element xx of the global character string array.

Source String to extract parts of.

src **gs[xx]** - Element xx of the global character string array.

Destination Variable to store the extracted part of the string.

dst

- **i[xx]** - Element xx of the local integer array.
- **gi[xx]** - Element xx of the global integer array.
- **f[xx]** - Element xx of the local float array.
- **gf[xx]** - Element xx of the global float array.
- **gfa[xx][yy]** - Element xx/yy of the global float field.
- **gs[xx]** - Element xx of the global string array.

Example(s):

You want to collect the pressure in mbar and the temperature in °C. The values are measured by an external application, assembled as a string and sent to Quadstar 32-bit by DDE where they are stored in gs[0]. The string now looks like this:

"1.5E-07mbar ;25°C"

First you separate the pressure value and store it in a local float variable; in this example, the number consists of 7 characters, so the instruction reads:

```
ExtractString(Type=number, Number=7, Source=gs[0],_
             Destination=f[0])
```

The value "1.5E-07" is now found in f[0], gs[0] still holds "mbar ;25°C". The information "mbar" is not needed any more, so you cut off everything until ";":

```
ExtractString(Type=character, Character=";", Source=gs[0],_
             Destination=gs[1])
```

gs[0] now holds only "25°C". Because the temperature value can consist of 1, 2 or 3 characters, you use "°" as separation character. This way it is guaranteed that the number is stored correctly in i[0] in any case:

```
ExtractString(Type=character, Character="°", Source=gs[0],_
             Destination=i[0])
```

8.7.16 FileCopy() - Copy a file

Description:

Copy an existing file to a new or existing file.

NOTE:

The function does not request any confirmation. An already existing destination file will be replaced.

See Chapter 8.7.19 FileRename() - Rename a file,  454 for examples.

Syntax:

```
FileCopy(SrcFile=StrgSrc, DestFile=StrgDst)
```

Parameters:

- **SrcFile** - Keyword for source file (name and path).
- **DestFile** - Keyword for destination file (name and path).

where:

- *StrgSrc, StrgDst* - Elements of the global string array (**gs[xx]**) or direct text inputs.

8.7.17 FileDelete() - Delete a file

Description:

Delete an existing file.

NOTE:

The function does not request any confirmation for deleting the file.

See Chapter 8.7.19 FileRename() - Rename a file,  454 for examples.

Syntax:

```
FileDelete(File=StrgFile)
```

Parameters:

- **File** - Keyword for file (name and path).

where:

- *StrgFile* - Element of the global string array (**gs[xx]**) or a direct text input.

8.7.18 FileMove() - Move a file

Description:

Move an existing file to a different path.

NOTE:

The function does not request any confirmation. An already existing file in the destination path will be replaced.

See Chapter 8.7.19 FileRename() - Rename a file, 454 for examples.

Syntax:

```
FileMove(SrcFile=StrgSrc, DestFile=StrgDst)
```

Parameters:

- **SrcFile** - Keyword for source file (name and path).
- **DestFile** - Keyword for destination file (name and path).

where:

- *StrgSrc, StrgDst* - Elements of the global string array (**gs[xx]**) or direct text inputs.

8.7.19 FileRename() - Rename a file

Description:

Rename an existing file.

NOTE:

The function does not request any confirmation. An already existing destination file will be replaced.

Syntax:

```
FileRename(SrcFile=StrgSrc, DestFile=StrgDst)
```

Parameters:

- **SrcFile** - Keyword for source file (name and path).
- **DestFile** - Keyword for destination file (name).

where:

- *StrgSrc, StrgDst* - Elements of the global string array (**gs[xx]**) or direct text inputs.

Example:

You are repeating a measurement 20 times. The latest completed run is to be saved in the temporary path as 'tempold' each. After that it shall be stored in the regular data path under a different name.

```
// Repeat the following program 20 times
Loop(i[0]=1;20)
Begin
```

```

// MID measurement using 'airsim2', displayed as a bargraph,
// store in 'temp' path as 'tempfile'
MID(Par="c:\qs32bit\par\airsim2", Disp=bar, _
    SaveCyc="c:\temp\tempfile";reset)
// Report the nbr. of the current measurement to the status bar.
Message(Text="Measurement running, Nbr. ";i[0])
// Delete the file 'tempold', if already existing
FileDelete(File="c:\temp\tempold.mdc")
// Copy current measurement file 'tempfile' to 'tempold'
FileCopy(SrcFile="c:\temp\tempfile.mdc", _
    DestFile="c:\temp\tempold.mdc")
// Rename the current measurement file 'tempfile' to_
'midcyc...',
// i.e. append the loop counter value to the file name
FileRename(OldFile="c:\temp\tempfile.mdc",
NewFile="midcyc.mdc";i[0])
// Move the file 'midcyc...' from the temporary path to the_
data path
FileMove(SrcFile="c:\temp\midcyc.mdc";i[0], _
    DestFile="c:\qs32bit\dat\midcyc.mdc";i[0])
End

```

8.7.20


FormatDate() - Format date

Description:

Format system time as a date.

The absolute system time (a large integer number that indicates the seconds passed since a standardized zero time) can be converted to a date by using FormatDate().

NOTE:

See Chapter 8.7.21 FormatTime() - Format time,  456 for examples.

Syntax:

```
FormatDate(Abs=AbsInt, Date=Datestring, Day=DInt, _
    Month=MoInt, Year=YInt)
```

Parameters:

- **Abs** - Keyword to get the absolute system time in seconds (integer) from AbsInt.
- **Date** - Keyword to store the date string (Windows® format) to Datestring.
- **Day, Month, Year** - Keywords to store the integer values for **day, month, year** to DInt, MoInt, YInt.

where:

- **Datestring** - Element of the global string array (**gs[xx]**).
- **DInt, MoInt, YInt, AbsInt** - Elements of the global (**gi[xx]**) or the local (**i[xx]**) integer array.

8.7.21 FormatTime() - Format time

Description:

Format system time as a daytime.

The absolute system time (a large integer number that indicates the seconds passed since a standardized zero time) can be converted to a time by using FormatTime().

Syntax:

```
FormatTime(Abs=AbsInt, Time=Datestring, Hour=HInt, Min=MiInt, _
           Sec=SInt)
```

Parameters:

- **Abs** - Keyword to get the absolute system time in seconds (integer) from AbsInt.
- **Time** - Keyword to store the time string (Windows® format) to Timestring
- **Hour, Min, Sec** - Keywords to store the integer values for hours, minutes, seconds to HInt, MiInt, SInt.

where:

- **Timestring** - Element of the global string array (**gs[xx]**).
- **HInt, MiInt, SInt, AbsInt** - Elements of the global (**gi[xx]**) or the local (**ii[xx]**) integer array.

Example:

Your aim is to measure MID data repeatedly for at least one hour. Because you can't use incomplete spectra, you check whether the end time is already reached before starting a new cycle. Finally you want to know how much time the measurement really took and at what date and time it stopped:

```
// Initialize i[0] to current time
GetTime(Abs=i[0])

// Repeat MID measurement for at least 1 hour
WhileVar(i[0]>0, BreakTime=1;h)
    MID(...)

// Get Time when loop was left and calculate running time
GetTime(Abs=i[1])
Calculate(i[2] = i[1]-i[0])

// Display running time
SetString(gs[0] = "The Loop ran for @ seconds")
SetString(gs[0] = gs[0];i[2])
Message(Text=gs[0], Confirm=on)
// Display stop date and time
FormatDate(Abs=i[1], Date=gs[1])
FormatTime(Abs=i[1], Time=gs[2])
SetString(gs[0] = "The Loop was stopped the @ at @ ")
SetString(gs[0] = gs[0];gs[1])
SetString(gs[0] = gs[0];gs[2])
Message(Text=gs[0], Confirm=on)
```


8.7.22 GetDate() - Get current date

Description:

Get current date.

The current date is loaded from the system clock and may then be stored in different ways (Windows® string, integer values for year/month/day, absolute system time in seconds) in variables.

NOTE:

See Chapter 8.7.23 `GetTime()` - Get current time, 457 and Chapter 8.7.21 `FormatTime()` - Format time, 456 for examples.

Syntax:

```
GetDate([Date=Datestring][,Day=DInt][,Month=MoInt][,Year=YInt]_
        [,Abs=AbsInt])
```

Parameters:

- **Date** - Keyword to load the date string (Windows® format) and assign it to *Datestring*.
- **Day, Month, Year** - Keywords to load the integer values for *day, month, year* and assign them to *DInt, MoInt, YInt*.
- **Abs** - Keyword to load the absolute system time in seconds (integer) and assign it to *AbsInt*.

where:

- **Datestring** - Element of the global string array (**gs[xx]**).
- *DInt, MoInt, YInt, AbsInt* - Elements of the global (**gi[xx]**) or the local (**i[xx]**) integer array.

8.7.23 GetTime() - Get current time

Description:

Get current time.

The current time is loaded from the system clock and may then be stored in different ways (Windows® string, integer values for hours/minutes/seconds, absolute system time in seconds) in variables.

NOTE:

See also Chapter 8.7.21 `FormatTime()` - Format time, 456 for examples.

Syntax:

```
GetTime(Time=Timestring, Hour=HInt, Min=MiInt, Sec=SInt,
        Abs=AbsInt)
```

Parameters:

- **Time** - Keyword to load the time string (Windows® format) and assign it to *Timestring*.
- **Hour, Min, Sec** - Keywords to load the integer values for hours, minutes, seconds and assign them to *HInt, MiInt, SInt*.
- **Abs** - Keyword to load the absolute system time in seconds (integer) and assign it to *AbsInt*.

where:

- **Timestring** - Element of the global string array (**gs[xx]**).

- *HInt, MilInt, SInt, AbsInt* - Elements of the global (**gi[xx]**) or the local (**i[xx]**) integer array.

Example(s):

You want to create a file within a sequence. The file name shall consist of the year, month, day and hour of its creation as two digit numbers. You decide to build this name by storing its parts in the correct order in the local integer array, starting with *i[0]*. The result string is then stored in *gs[0]* and you make sure that all numbers (day, hour etc.) contain two digits each:

```
// Get current date and time
GetDate(Day = i[2], Month = i[1], Year = i[0])
GetTime(Hour = i[3])

// Initialize result string (4 alphanumeric pieces) and aux string
// (leading zero)
SetString(gs[0] = "####")
SetString(gs[1] = "0@")

// if i[0] (year) is single digit, add a leading 0
IfVar(i[0] < 10)
    SetString(gs[0] = gs[0]; gs[1])
//replace first '@' of result string by current year
SetString(gs[0] = gs[0]; i[0])

IfVar(i[1] < 10)
    SetString(gs[0] = gs[0]; gs[1])
SetString(gs[0] = gs[0]; i[1])

IfVar(i[2] < 10)
    SetString(gs[0] = gs[0]; gs[1])
SetString(gs[0] = gs[0]; i[2])

IfVar(i[3] < 10)
    SetString(gs[0] = gs[0]; gs[1])
SetString(gs[0] = gs[0]; i[3])

// Display the result string
Message(Text = gs[0], Confirm = on)
```

The string (e.g. '95061916' for June 19, 1995 at 4 p.m.) is now stored in *gs[0]* and may be used as a file name. The problem may be solved shorter but less obviously the following way:

```
...
// Get current date and time
GetDate(Day = i[2], Month = i[1], Year = i[0])
GetTime(Hour = i[3])

// Initialize result string (4 alphanumeric pieces) and aux string
// (leading zero)
SetString(gs[0] = "####")
SetString(gs[1] = "0@")

// check local integers 0..3 (Year..Hour) for single digits and
// build the result string
Loop(i[4] = 0;3)
    Begin
```

```

// if single digit, add a leading 0
IfVar(i[i[4]] < 10)
  SetString(gs[0] = gs[0]; gs[1])
// replace i[4]th '@' of result string by matching number
SetString(gs[0] = gs[0]; i[i[4]])
End

```

8.7.24

GetErrorCode() - Read error and status register

Description:

Read the error/status frames and store them in selectable variables. The codes are briefly explained here. For a detailed description refer to the User's guide of the QMS 422/421, or QMS 200.

NOTE:

- The **IfBit()** instruction helps you to evaluate the state of the flags (separate bits in a variable) contained in Error-, Warning- etc. variables.
- The Error-flags are cleared after being read and are only set again (after approx. 1 second) if the error occurs still/again.

In the following tables

- **X** means: valid for this equipment
- **-** means: not valid for this equipment

Error Variable

Bit No	Meaning when 1/0	Message No	QMS 422/421	QMS 200	QMU 112
0	QMS controller stack overflow/-	17	X	X	-
1	QMS reset error/-	18	X	X	-
2	QMS controller watchdog error/-	19	X	X	-
3	QMS controller EPROM defect/-	20	X	X	-
4	QMS controller NOVRAM defect/-	21	X	X	-
5	QMS controller dual port RAM defect/-	22	X	X	-
6	QMS controller buffer RAM defect/-	23	X	X	-
7	DSP EPROM defect/-	24	X	X	-
8	DSP dual port RAM defect/-	25	X	X	-
9	Monitor - DAC error/-	26	X	X	-
10	Resolution - DAC error/-	27	X	X	-
11	AO - DAC error/-	28	X	X	-
12	ADC error/-	29	X	X	-
13	QMS temperature underrange/-	30	-	X	-
14	QMS temperature overrange/-	31	-	X	-

Bit No	Meaning when 1/0	Mes- sage No	QMS 422/ 421	QMS 200	QMU 112
15	RF-ID error/-	32	-	X	-
16	QMS controller idle error/-	33	X	X	-
17	SEM error/-	34	X	X	-
18	CD error/-	35	X	X	-
19	Ion source error/-	36	X	X	-
20	Filament 1 defect/-	37	X	X	-
21	Filament 2 defect/-	38	X	X	-
22	Emission error/-	39	X	X	-
23	CAN error/-	40	X	X	-
24	Parameter lost in NOVRAM/-	41	X	X	-
25	Communication error controller <=> console/-	42	X	-	-
26	Communication error controller <=> DSP/-	43	X	X	-
27	Communication error on LAN (ArcNet)/-	44	X	X	-
28	Communication error on RS-232-C/-	45	X	X	-
29	RF error/-	46	X	X	-
30	Integrator: intensity < underground/-	47	X	X	-

Warning Variable

Bit No	Meaning when 1/0	Mes- sage No	QMS 422/ 421	QMS 200	QMU 112
0	No High Voltage supply/-	17	X	-	-
1	No Ion Source supply/-	18	X	-	-
2	Cannot Degas Fil 1+2/-	19	X	-	-
3	No Auto (Fil 1+2) on Degas/-	20	X	-	-
4	Only Filament 1 available	21	X	-	-
5	Buffer failure	22	X	-	-

Status Variable

Bit No	Meaning when 1/0	Message No	QMS 422/421	QMS 200	QMU 112
0	Cycle run/halt	-	X	X	-
1	Multi-/mono- channel	-	X	X	-
2	Emission on/off	-	X	X	-
3	SEM Supply on/off	-	X	X	-
4	Wait for external trigger/-	-	X	X	-
5	Settling run/halt	-	X	-	-
6	Integrator underground run/halt	-	X	-	-
7	Display emission current/electrometer value	-	X	X	-
8	Degas on/off	-	X	X	-
9	Adjust on/off	-	X	X	-
10	Adjust run/-	-	X	X	-
14	Ringbuffer empty/not empty	-	X	X	-
15	Ringbuffer overflow/-	-	X	X	-

QMUEmi

Bit No	Meaning when 1/0	Message No	QMS 422/421	QMS 200	QMU 112
0	QMU 0: Emission failure/ok	-	-	-	X
1	QMU 1: Emission failure/ok	-	-	-	X
..
7	QMU 7: Emission failure/ok	-	-	-	X

Syntax:

```
GetErrorCode([Error=err][,Warning=warn][,Status=sta]_
             [,QmuEmi=emi])
```

Parameters:

- Error** Keyword for error code.
- err* Variable to be filled with the error code.
- **i[xx]** - Local integer variable
 - **gi[xx]** - Global integer variable
- Warning** Keyword for warning code.

<i>warn</i>	Variable to be filled with the warning code. <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable
Status	Keyword for status code.
<i>sta</i>	Variable to be filled with the status code <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable
QmuEmi	Keyword for QMU emission status code.
<i>emi</i>	Variable to be filled with the QMU emission status code <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable

Example(s):

```
GetErrorCode(Error=i[0], Warning=i[1], Status=i[2])
Loop(i[3]=0;7)
  IfBit(i[0];i[3] = set)
    Error(Text="Error, Bit @ in error register: Set !";i[4])
Loop(i[3]=0;7)
  IfBit(i[1];i[3] = set)
    Error(Text="Warning, Bit @ in warning register: Set !";i[4])
...
```

Read the QMS error and status word and print out the flags that are set.

8.7.25**GetSysInfo() - Get System Information****Description:**

Get information about general system properties like free disc space, horizontal and vertical screen resolution. The screen resolution information is useful for using instructions like DispWindow independently of the graphic system, while the free disc space information can prevent a storage overflow.

Syntax:

```
GetSysInfo([DiscSpace=space ;drive][,ResHorz=hRes]_
           [,ResVert=vRes])
```

Parameters:

DiscSpace	Keyword for disc space (in KiloBytes).
<i>space</i>	Variable to hold the value of the free disc space in kBytes <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable
<i>drive</i>	Disc drive (e.g. A, C, ...) <ul style="list-style-type: none"> • "b" - Letter that identifies the disc drive. • gs[xx] - Element xx of the global string array, if the letter is stored in gs[].
ResHorz	Keyword for horizontal screen resolution.

<i>hRes</i>	Variable to hold the value of the horizontal resolution <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable
ResVert	Keyword for vertical screen resolution.
<i>vRes</i>	Variable to hold the value of the vertical resolution <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable

Example(s):

```
GetSysInfo(DiscSpace=i[0];"c", ResHorz=i[1], ResVert=i[2])
```

Store the number of free Kilobytes of hard drive C: in i[0], store horizontal and vertical screen resolution in i[1] and i[2].

```
GetSysInfo(DiscSpace=i[0];"c")
WhileVar(i[0] > 1000)
Begin
  MID(Par="c:\qs32bit\par\helium4.mip",_
      SaveCyc= "c:\qs32bit\dat\helium_c.mdc")
  GetSysInfo(DiscSpace=i[0];"c")
End
Message(Text="Drive C is full, switched over to drive D")
WhileVar(i[0] > 0)
Begin
  MID(Par="c:\qs32bit\par\helium4.mip",_
      SaveCyc= "d:\spare\helium_d.mdc")
End
...
```

Carry out a helium MID measurement and store the results on drive C, until the residual capacity of drive C falls below 1000 kBytes. Then store further data on disc drive D and report this in the status bar.

8.7.26**GSC() - Gas-specific calibration****Description:**

Perform a gas-specific calibration. To increase precision, a mean value calculation is available.

Syntax:

```
GSC(Par=file[;ind][,Disp=d][,Prot=p][,Mean=nbr])
```

Parameters:

Par	Keyword for parameter file. No file extension should be used because it is added automatically by Quadstar 32-bit.
<i>file</i>	(Drive:\Path\)\ Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416

Disp Keyword for display mode.

d Display mode

- **[default]** - Default display mode.
- **off** - Display off.
- **on** - Display on.

Without any explicit declaration, the default display mode is on. The default settings can be altered by **SetDefault()**.

Prot Keyword for log mode.

p Log mode.

- **[default]** - Default log mode.
- **off** - Log function off.
- **on** - Printer log on.

See Chapter 8.4.3 Log output, 422.

If nothing else is specified, the default log mode is off. It can be altered with the **SetDefault()** instruction.

Mean Keyword for mean value calculation.

nbr Number of measurements used for mean value calculation (1-10000)

- *cc* - Direct input of a numeric value (integer)
- **i[xx]** - Local integer variable
- **gi[xx]** - Global integer variable

Example(s):

```
GSC(Par="gcal", Mean=10)
```

Perform a gas-specific calibration using the "gcal" parameter file. Take the mean value out of 10 measurements.

```
...
Loop(i[0]=0;3)
Begin
    // Flush Inlet
    Sequence(Par="flush")
    // Open Inlet
    Sequence(Par="open";i[0])
    // Calibrate
    GSC(Par="gcal";i[0])
    // Close Inlet
    Sequence(Par="close";i[0])
End
...
```

Perform several gas-specific calibrations. The files "gcal0" ... "gcal3" with one inlet sequence each ("open0" ... "open3") are used as parameter files.

8.7.27

IfBit() - Testing a bit of an integer variable

Description:

Test an individual bit of an integer variable.

This instruction is used in conjunction with the **GetErrorCode()** instruction:

=> **GetErrorCode()** reads the Error and the Status register.

=> **IfBit()** analyzes the individual error and status flags represented by bits.

Syntax:

IfBit(*var*; *bit=state*)

Parameters:

<i>var</i>	Variable containing the bit to be tested <ul style="list-style-type: none"> • i[xx] - Local integer variable • gi[xx] - Global integer variable
<i>bit</i>	Number of the bit to be tested (0 ... 31) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer) • i[xx] - Local integer variable • gi[xx] - Global integer variable
<i>state</i>	Comparison state <ul style="list-style-type: none"> • clear - Test for bit not set • set - Test for bit set

Example(s):

Read the QMS 422/421 or QMS 200 error code and write it into *i[0]*. If bit 3 is set, perform a function:

```
GetErrorCode(Status=i[0])
IfBit(i[0];3=set)
Begin
  ...
End
```

8.7.28

IfDI() - Testing a digital input

Description:

If instruction for monitoring a digital input of the control unit, where the number of available channels on

- QMS 422/421 depends on the number of installed DI 420 boards
- QMS 200 is limited to 2 (DI 0: Extern Protection, DI 1: Extern Start)

The instruction that follows **IfDI()** is executed only if the condition (*var*{;unit}=*state*) is true.

Specify the digital input of the control unit as follows:

- enter the channel number or
- enter the logical name of the input channel.

=> The logical names are specified in the Parameter Setup program under [Config] > [DI/DO Identifications].

Syntax:

IfDI(*var*{;unit}=*state*)

Parameters:

- var* number of the digital input channel.
QMS 422/421: 0 ... 31(63)/QMS 200: 0 ... 1
- “zzzzz” - Logical name of the digital input.
 - *cc* - Direct input of a numeric value (integer).
 - *i[xx]* Element *xx* of the local integer array.
 - *gi[xx]* - Element *xx* of the global integer array.
- unit* Identification of the control unit, only available in multiplex mode.
- “zzzzz” - Logical name of the control unit.
 - *cc* - Direct input of the internal logical number.
- => Only those control units are available whose status has been set to Enable in the Parameter Setup program under [Comm] > [Node Identification]. This entry does not influence the unit selected with **SetUnit()**.
- state* DI comparison status
- **off** - Test DI channel status = off.
 - **on** - Test DI channel status = on.

Example(s):

```
IfDI(i[23]=on)
  Delay(time=5)
```

If the digital input whose number is stored in *i[23]* is in the on state, a five second pause is counted down.

```
IfDI(12; "Q_Unit_1"=off)
  Delay(Time=5)
```

Multiplex mode: If digital input 12 of unit “Q_Unit_1” is off, a pause of five seconds is counted down.

```
IfDI("inlet5"=off)
Begin
  Delay(Time=5)
  SetDO(10=on)
End
```

If the digital input “inlet5” is in the off state, a 5 second pause is counted down and digital output 10 is subsequently set.

8.7.29**IfDO() - Testing a digital output****Description:**

If instruction for monitoring a digital output of the control unit, where the number of available channels on:

- QMS 422/421 depends on the number of installed DO 420 boards
- QMS 200 is limited to 2 (DO 0, DO 1)

The instruction that follows *IfDO()* is executed only if the condition (*var*{;unit} =*state*) is true.

The digital output of the control unit is specified as follows:

- enter the channel number or
- enter the logical name of the output channel.

=> The logical names are specified in the Parameter Setup program under [Config] > [DI/DO Identifications].

Syntax:

```
IfDO(var{;unit}=state)
```

Parameters:

<i>var</i>	number of the digital output channel. QMS 422/421: 0 ... 31(63, 95)/QMS 200: 0 ... 1 <ul style="list-style-type: none"> • “zzzzz” - Logical name of the digital output. • <i>cc</i> - Direct input of a numeric value (integer). • <i>i</i>[<i>xx</i>] - Element <i>xx</i> of the local integer array. • <i>gi</i>[<i>xx</i>] - Element <i>xx</i> of the global integer array.
<i>unit</i>	Identification of the control unit, only available in multiplex mode. <ul style="list-style-type: none"> • “zzzzz” - Logical name of the control unit. • <i>cc</i> - Direct input of the internal logical number. <p>=> Only those control units are available whose status has been set to Enable in the Parameter Setup program under [Comm] > [Node Identification]. This entry does not influence the unit selected with SetUnit().</p>
<i>state</i>	DO comparison status <ul style="list-style-type: none"> • off - Test DO channel status = off. • on - Test DO channel status = on.

Example(s):

```
IfDO(gi[8]=on)
  Delay(time=5)
```

If the digital output whose number is stored in *gi*[8] is in the on state, a five second pause is counted down.

```
IfDO(63;"Q_Unit_1"=on)
  Delay(Time=5)
```

Multiplex mode: If the digital output 63 of the unit *Q_Unit_1* is on, a five second pause is counted down.

```
IfDO("cal gas 7"=off)
Begin
  Delay(time=5)
  SetDo(10=on)
End
```

If the digital output “cal gas 7” is in the off state, a five second pause is counted down and digital output 10 is subsequently set.

8.7.30

IfMeaVal() - Testing a measured value

Description:

If instruction for monitoring a measured value on one of the 64 measurement channels of the mass spectrometer. An MID measurement is started invisibly using the given

parameter file, the measured value is recorded and the measurement is stopped again. The instruction following `IfMeaVal()` is executed only if the condition (`var cond val`) is true.

Syntax:

```
IfMeaVal(var cond val,Par=file[,ind])
```

Parameters:

<i>var</i>	Measure channel number (0 ... 63). <ul style="list-style-type: none"> • cc Direct input of a numeric value (integer or float). • i[xx] Element <i>xx</i> of the local integer array. • gi[xx] Element <i>xx</i> of the global integer array.
<i>cond</i>	Condition <ul style="list-style-type: none"> • = - equal. • < - less than. • > - greater than. • <= - less or equal • >= - greater or equal. • <> - not equal.
<i>val</i>	Comparison variable All float and integer array element types can be used: <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer or float) • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array. • f[xx] - Element <i>xx</i> of the local float array. • gf[xx] - Element <i>xx</i> of the global float array. • gfa[xx][yy] - Element <i>xx/yy</i> of the global float field.
Par	Keyword for parameter file
<i>file</i>	(Drive:\Path\)\ Filename. <ul style="list-style-type: none"> • "#####" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element <i>xx</i> of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416

Example(s):

```
IfMeaVal(5<3.4E-7,Par="myfile")
    Delay(Time=5)
```

If the values on channel 5, measured with the "myfile" parameters, are less than 3.4E-7, a pause of five seconds is counted down.

```
IfMeaVal(gi[4]>12.5E-9,Par="d:\qs\dat\test1\sample5")
    Delay(Time=5)
```

If the values on the channel stored in `gi[4]`, measured with the "d:\qs\dat\test1\sample5" parameters, is greater than 12.5E-9, a pause of five seconds is counted down.

8.7.31 IfString() - Compare two text strings

Description:

If-instruction to compare two text strings. The instruction following **IfString()** is only executed, if the condition (*string1 cond string2*) is true.

Syntax:

```
IfString(string1 cond string2)
```

Parameters:

<i>string1</i>	First string to compare <ul style="list-style-type: none"> • <i>cc</i> - Direct text input • gs[xx] - Element <i>xx</i> of the global string Array.
<i>cond</i>	Condition <ul style="list-style-type: none"> • = - equal to. • < - less than. • > - greater than. • <= - less or equal. • >= - greater or equal. • <> - not equal.
<i>string2</i>	Second string to compare <ul style="list-style-type: none"> • <i>cc</i> - Direct text input • gs[xx] - Element <i>xx</i> of the global string Array.

The two strings are compared character by character, disregarding capitals. The first two characters that don't match each other in both strings determine the result of the comparison: the character that comes first in alphabet is considered to be 'smaller'. Two strings are only recognized as equal if all their characters are equal (except capital/non capital).

Example(s):

```
IfString(gs[8]="Wait")
  Delay(Time=50)
```

If the string *gs[8]* contains the text 'Wait', then wait 50 s.

```
IfString(gs[0]>gs[1])
  SetString(gs[2]=gs[0])
Else
  SetString(gs[2]=gs[1])
```

Assign the "greater" of the two strings *gs[0]* and *gs[1]* to *gs[2]*.

8.7.32 IfVar() - Testing an array or field variable

Description:

If instruction for monitoring a variable. The instruction following **IfVar()** is executed only if the condition (*var cond val*) is true.

Syntax:

```
IfVar(var cond val)
```

Parameters:

<i>var</i>	Variable to be tested. All float and integer array element types can be used: <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer or float) • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>f[xx]</i> - Element <i>xx</i> of the local float array. • <i>gf[xx]</i> - Element <i>xx</i> of the global float array. • <i>gfa[xx][yy]</i> - Element <i>xx/yy</i> of the global float field.
<i>cond</i>	Condition <ul style="list-style-type: none"> • = - equal to. • < - less than. • > - greater than. • <= - less or equal. • >= - greater or equal. • <> - not equal.
<i>val</i>	Comparison value All float and integer array element types can be used: <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer or float) • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>f[xx]</i> - Element <i>xx</i> of the local float array. • <i>gf[xx]</i> - Element <i>xx</i> of the global float array. • <i>gfa[xx][yy]</i> - Element <i>xx/yy</i> of the global float field.

Example(s):

```
IfVar(i[8]=10)
    Delay(Time=5)
```

If the variable *i[8]* has a value of 10, a pause of five seconds is counted down.

```
IfVar(gfa[gi[7]][6]<1.5E-8)
Begin
    SetDO(5=on,10;min)
    SetDO(6=on,1;min)
End
```

If the value in the array element *gfa[gi[7]][6]* is less than 1.5E-8, DO 5 is set, a pause of 10 minutes is counted down, after which DO 6 is set and another pause of 1 minute is counted down.

8.7.33**IniFileRead() - Read a value from an INI file****Description:**

Read a value from an INI file.

INI files are files with the extension '.ini', that are normally used (in Windows®) to setup configurations. They are built up using a certain format (for more information please consult your Windows® documentation):

```
...
[SectionIdentification1]
```

```

EntryA=ValueA
EntryB=ValueB
[SectionIdentification2]
EntryX=ValueX
EntryY=ValueY
...

```

You can use such files to read and write data in a program independent format.

Syntax:

```
IniFileRead(Section=SecId, Entry=EntId, Data=Var, IniFile=File)
```

Parameters:

Section	Keyword for section, where the entry is found.
<i>SecId</i>	Section identifier. <ul style="list-style-type: none"> • <i>cc</i> - Direct text input. • gs[xx] - Element <i>xx</i> of the global string array.
Entry	Keyword for entry, where the value is found.
<i>EntId</i>	Entry identifier. <ul style="list-style-type: none"> • <i>cc</i> - Direct text input. • gs[xx] - Element <i>xx</i> of the global string array.
Data	Keyword for data container.
<i>Var</i>	Variable to assign the read value to. <ul style="list-style-type: none"> • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array. • f[xx] - Element <i>xx</i> of the local float array. • gf[xx] - Element <i>xx</i> of the global float array. • gfa[xx][yy] - Element <i>xx/yy</i> of the global float field. • gs[xx] - Element <i>xx</i> of the global string array.
IniFile	Keyword for file to read from.
<i>File</i>	File name (and possibly path). <ul style="list-style-type: none"> • <i>cc</i> - Direct text input. • gs[xx] - Element <i>xx</i> of the global string array.

All values are read as standard texts. If a value isn't of the expected format, it's assumed to be empty/zero.

Example(s):

See Chapter 8.7.34 IniFileWrite() - Write a value to an INI file, 471.

8.7.34

IniFileWrite() - Write a value to an INI file

Description:

Write a value to an INI file. See also Chapter 8.7.33 IniFileRead() - Read a value from an INI file, 470.

Syntax:

```
IniFileWrite(Section=SecId, Entry=EntId, Data=Var, IniFile=File)
```

Parameters:

Section	Keyword for section, where the value is to be entered.
SecId	Section identifier. <ul style="list-style-type: none"> • <i>cc</i> - Direct text input. • gs[xx] - Element <i>xx</i> of the global string array.
Entry	Keyword for entry to assign the value to.
EntId	Entry identifier. <ul style="list-style-type: none"> • <i>cc</i> - Direct text input. • gs[xx] - Element <i>xx</i> of the global string array.
Data	Keyword for data container.
Var	Value or variable that contains the value to be written. <ul style="list-style-type: none"> • <i>nn</i> - Direct input of a numeric value (integer or float). • <i>cc</i> - Direct text input. • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array. • f[xx] - Element <i>xx</i> of the local float array. • gf[xx] - Element <i>xx</i> of the global float array. • gfa[xx][yy] - Element <i>xx/yy</i> of the global float field. • gs[xx] - Element <i>xx</i> of the global string array.
IniFile	Keyword for file to write the data to.
File	File name (and possibly path). <ul style="list-style-type: none"> • <i>cc</i> - Direct text input. • gs[xx] - Element <i>xx</i> of the global string array.

Example:

You want to run Quadstar 32-bit controlled by a supervising control unit. This unit informs Quadstar 32-bit, what measurement is to be performed and fetches the result afterwards. For this purpose, the file *remote.ini* is created, that contains in the section *WhatMea* the name of the parameter file that defines the next measurement and in the section *Result* the measured value. This value is read and interpreted by the control unit.

File *remote.ini*, after the control unit has requested to perform the measurement *PeakO2.mip*:

```

---File: remote.ini ---
[WhatMea]
ParFile=PeakO2.mip
[Result]
MeaVal=
-----

// Read name of the parameter file used for the next measurement
// into gs[0]
IniFileRead(Section="WhatMea",Entry="ParFile",Data=gs[0],_
            IniFile="c:\control\remote.ini")
// Perform measurement using the given parameter file;
// save results into gfa[0][0] (only one value)
MID(Par=gs[0], SaveGfa=0)
// Write measurement result to the matching entry in the
// ini-file
IniFileWrite(Section="Result",Entry="MeaVal",Data=gfa[0][0],_

```



```
IniFile="c:\control\remote.ini")
```

File remote.ini, after Quadstar 32-bit has performed the measurement PeakO2.mip:

```
---File: remote.ini ---
[WhatMea]
ParFile=PeakO2.mip
[Result]
MeaVal=1.23456E-6
-----
```

8.7.35 InitDisp() - Initialize the display

Description:

Initialize the display.

This instruction is required for example, if a table is displayed in Multiple Column Mode, to remove the data of previous measurements.

Syntax:

```
InitDisp()
```

Parameters:

None

Example(s):

```
InitDisp()
```

The display is cleared and initialized.

8.7.36 Loop() - Loop

Description:

Loop with counter variable.

The loop starts with the counter value from that's incremented with step every time the loop is passed. If the counter variable reaches the value to, the loop is exited. For the counter variable, only integer array elements (**i[]**, **gi[]**) are allowed. For the parameters from, to and step, also integer constants are possible.

Syntax:

```
Loop(var=from;to[;step])
```

Parameters:

<i>var</i>	Loop variable. <ul style="list-style-type: none"> i[xx] - Element xx of the local integer array. gi[xx] - Element xx of the global integer array.
<i>from</i>	Start value (-1'000'000 ... 1'000'000) <ul style="list-style-type: none"> <i>cc</i> - Numeric value (integer). i[xx] - Element xx of the local integer array. gi[xx] - Element xx of the global integer array.

to End value (-1'000'000 ... 1'000'000)

- *cc* - Numeric value (integer)
- *i[xx]* - Element *xx* of the local integer array.
- *gi[xx]* - Element *xx* of the global integer array.

step Step width (-1'000'000 ... 1'000'000) (default value = 1)

- *cc* - Numeric value (integer)
- *i[xx]* - Element *xx* of the local integer array.
- *gi[xx]* - Element *xx* of the global integer array.

Example(s):

```
Loop(i[5]=0;15)
  SetDO(i[5]=off)
```

The digital outputs of channels 0 ... 15 are switched off.

```
Loop(gi[12]=60;30;-2)
Begin
  Calculate(gfa[0][gi[12]]=gfa[1][gi[12]]/f[3]*100.0)
  Calculate(gfa[1][gi[12]]=gfa[1][gi[12]]*2)
End
```

The loop starts with *gi[12]=60* and is decremented (step = -2) until *gi[12]* has a value of equal or less than 30.

8.7.37 LoopTable() - Loop Table

Description:

Loop with value table.

LoopTable() is a special version of the **Loop()** function. It doesn't run through a series of numbers with fixed increment, but through any given number table.

The exit condition and the incremental step are missing, instead there is a table of integer values or variables.

Syntax:

```
LoopTable(var=num1[;num2][;num3][;...])
```

Parameters:

var Loop variable.

- *i[xx]* - Element *xx* of the local integer array.
- *gi[xx]* - Element *xx* of the global integer array.

num1..numX table values

- *cc* Integer value (-1'000'000 ... 1'000'000)
- *i[xx]* - Element *xx* of the local integer array.
- *gi[xx]* - Element *xx* of the global integer array.

Example(e):

```
// Valves 10, 12, 20 and 22 are opened.
LoopTable(i[34]=10;12;20;22)
  SetDO(i[34]=on)
```

```
// The DO's 7, 18, 34 etc. are turned on consecutively, in
```

```
// between the DO 12 is switched on and off
LoopTable(gi[1]=7;12;18;12;34;12;23;12;40;12;42;12;44;12;4;12)
Begin
  SetDO(gi[1]=on)
  Delay(Time=2)
  IfVar(gi[1]=12)
    SetDO(gi[1]=off)
End

// The numbers of the DO's are stored in gi[]. Switch the ones
// in position i[0,4,8,12] of gi[] on.
LoopTable(i[10]=gi[i[0]];gi[i[4]];gi[i[8]];gi[i[12]])
  SetDO(i[10]=on)
```

8.7.38 MCD() - MCD measurement

Description:

Perform a Multiple Concentration Detection (MCD) measurement.

Options:

- display the measurement data,
- store the measurement data in the field gfa[][] containing global float values,
- store the measurement data in a file containing MCD measurement data,
- store the measurement data in a file containing reference measurement data.

Syntax:

```
MCD(Par=file[;ind][,Disp=d][,Prot=p][,SaveGfa=var]_
    [,SaveRef=file[;ind][,SaveCyc=file[;ind][;m]])
```

Parameters:

Par	Keyword for parameter file.
file	(Drive:\Path\)\ Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
ind	File index - See Chapter 8.4.1.3 File names with indices, 416
Disp	Keyword for display mode.
d	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • tab - Table display. • bar - Bargraph display. • vt - Versus Time display. <p>Without an explicit entry the default mode is set to table display (tab). It can be altered with the SetDefault() instruction.</p>
Prot	Keyword for log mode

- p* Log mode.
- **[default]** - Default log mode.
 - **off** - Log function off.
 - **on** - Printer log on.
- See Chapter 8.4.3 Log output, 422.
- If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.
- SaveGfa** Keyword for storing measurement data in an array of the global data field **gfa[][]**:
- The measurement data are stored in the float array *xx* (=var) of the **gfa[xx][[]]** field. Specification of the index *yy* for **gfa[][]** is not necessary, the measured value of the 1st element is stored in **gfa[][0]**, the measured value of the 2nd element is stored in **gfa[][1]**, etc..
 - All data stored in the elements *yy* = 0 ... 63 of the float array *xx*=var of **gfa[xx][yy]** are cleared before the data are stored.
- var* Float array number (0 <= *var* <= configurable maximum).
- **cc** - Numeric value (integer).
 - **i[xx]** - Element *xx* of the local integer array.
 - **gi[xx]** - Element *xx* of the global integer array.
- SaveRef** Keyword for storing in a file containing reference data.
- file* (Drive:\Path\)
Filename.
- "#####" - Filename (e.g. "c:\qs\myfile").
 - **gs[xx]** - Element *xx* of the global character string array, if the filename is stored in **gs[]**.
- ind* File index - See Chapter 8.4.1.3 File names with indices, 416
- SaveCyc** Keyword for storing the Cycle data in a file
- file* (Drive:\Path\)
Filename.
- "#####" - Filename (e.g. "c:\qs\myfile").
 - **gs[xx]** - Element *xx* of the global character string array if the filename is stored in **gs[]**.
- ind* File index - See Chapter 8.4.1.3 File names with indices, 416
- m* Save mode
- **[append]** - Add the data to the end of the file (default).
 - **reset** - Clear the file before storing the new data.

Example(s):

```
MCD( Par="mypar" )
```

Perform an MCD measurement with the parameter file "mypar" and display the measured values as a table (default).

```
MCD( Par="c:\qs\par\mcd" , Disp=vt , SaveCyc="c:\qs\dat\mcd" )
```

Perform an MCD measurement with the parameter file "c:\qs\par\mcd", display the data in Versus Time mode and store them in the data file "c:\qs\dat\mcd".

8.7.39 Message() - Message output

Description:

Output a message in the status bar of the measurement window, the Event window (**Confirm=off**) or in a separate dialog box with an [OK] button to confirm the message (**Confirm=on**).

Syntax:

```
Message(Text=text[;strval][,Confirm=c][,Prot=p]_
        [,Delay=tim [;unit][,Save=sav][,Time=t][,Disp=d])
```

Parameters:

Text	Keyword for text.
<i>text</i>	Text. <ul style="list-style-type: none"> • "Text" - Direct input of text. • gs[xx] - Element xx of the global character string array.
<i>strval</i>	String value - See Chapter 8.4.1.4 String values, 417
Confirm	Keyword for confirmation mode. <p>With confirm, the sequence being executed can be stopped until the message (displayed in a separate dialog box) has been confirmed with the [OK] button.</p>
<i>c</i>	Confirmation mode. <ul style="list-style-type: none"> • [off] - Sequence continues without confirmation (default). • on - Sequence is stopped until the [OK] command is given.
Prot	Keyword for log mode.
<i>p</i>	Log mode. <ul style="list-style-type: none"> • [default] - Default log mode. • off - Log function off. • on - Printer log on. <p>See Chapter 8.4.3 Log output, 422.</p> <p>If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.</p>
Delay	Keyword for time delay. The message remains on the screen until the delay has been counted down.
<i>tim</i>	Time specification (1 ... 1'000'000). <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • ms - Time in milliseconds. • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
Save	Keyword for save message.

sav	Save mode. <ul style="list-style-type: none"> • [off] - Message is not stored in the Error History. • on - Message is stored in the Error History.
Time	Keyword for time information.
t	Message with or without time information. <ul style="list-style-type: none"> • [off] - Message does not contain time information. • on - Message contains time information.
Disp	Keyword for display.
d	Indication, where the message is to be displayed. <ul style="list-style-type: none"> • off - Don't display message. • [status] - Display message in status bar. • event - Display message in 'Event' window. • both - Display message in 'Event' window and status bar.

Example(s):

```
Message(Text="calibration of bottle 5 in progress")
```

Output a message, the sequence continues.

```
Message(Text=gs[i[4]],Confirm=on)
```

Output a message and stop the sequence until the operator confirms by [OK] that the sequence shall be continued.

8.7.40 MID() - MID measurement

Description:

Perform a Multiple Ion Detection (MID) measurement.

Options:

- Display the measurement data,
- store the measurement data in the field gfa[][] containing global float values,
- store the measurement data in a file containing MID measurement data,
- store the measurement data in a file containing reference measurement data.

Syntax:

```
MID(Par=file[;ind][,Disp=d][,Prot=p][,SaveGfa=var]_
    [,SaveRef=file[;ind][,SaveCyc=file[;ind][;m])
```

Parameters:

Par	Keyword for parameter file.
file	(Drive:\Path\)\ Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
ind	File index - See Chapter 8.4.1.3 File names with indices, 416
Disp	Keyword for display mode.

- d* Display mode.
- **[default]** - Default display mode.
 - **off** - Display off.
 - **tab** - Table display.
 - **bar** - Bargraph display.
 - **vt** - Versus Time display.
- Without an explicit entry the default mode is set to table display (**tab**). It can be altered with the SetDefault() instruction.
- Prot** Keyword for log mode
- p* Log mode.
- **[default]** - Default log mode.
 - **off** - Log function off.
 - **on** - Printer log on.
- See Chapter 8.4.3 Log output, 422.
- If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.
- SaveGfa** Keyword for storing measurement data in an array of the global data field **gfa[[]]**:
- The measurement data are stored in the float array *xx* (=var) of the **gfa[xx][[]]** field. Specification of the index *yy* for **gfa[[]]** is not necessary, the measured value of the 1st element is stored in **gfa[[]][0]**, the measured value of the 2nd element is stored in **gfa[[]][1]**, etc..
 - All data stored in the elements *yy* = 0 ... 63 of the float array *xx*=var of **gfa[xx][yy]** are cleared before the data are stored.
- var* Float array number (0 <= var <= configurable maximum).
- **cc** - Numeric value (integer).
 - **i[xx]** - Element *xx* of the local integer array.
 - **gi[xx]** - Element *xx* of the global integer array.
- SaveRef** Keyword for storing in a file containing reference data.
- file* (Drive:\Path\)\ Filename.
- "ffffff" - Filename (e.g. "c:\qs\myfile").
 - **gs[xx]** - Element *xx* of the global character string array, if the filename is stored in **gs[[]]**.
- ind* File index - See Chapter 8.4.1.3 File names with indices, 416
- SaveCyc** Keyword for storing the Cycle data in a file
- file* (Drive:\Path\)\ Filename
- "ffffff" - Filename (e.g. "c:\qs\myfile").
 - **gs[xx]** - Element *xx* of the global character string array if the filename is stored in **gs[[]]**.
- ind* File index - See Chapter 8.4.1.3 File names with indices, 416
- m* Save mode
- **[append]** - Add the data to the end of the file (default).
 - **reset** - Clear the file before storing the new data.

Example(s):

```
MID(Par="mypar")
```

Perform a MID measurement with the parameter file "mypar" and display the measured values as a table (default).

```
MID(Par="c:\qs\par\sample",Disp=bar,SaveCyc="c:\qs\dat\midcyc")
```

Perform a MID measurement with the parameter file "c:\qs\par\sample", display the data in Versus Time mode and store them in the data file "c:\qs\dat\midcyc".

8.7.41 MSC() - Mass scale calibration

Description:

Perform a mass scale calibration.

Syntax:

```
MSC(Par=file[;ind][,CalMode=c][,Disp=d][,Prot=p][,SaveMode=m])
```

Parameters:

Par	Keyword for parameter file.
<i>file</i>	(Drive:\Path) Filename. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
CalMode	Keyword for calibration mode.
<i>c</i>	Calibration mode. <ul style="list-style-type: none"> • [coarse] - fast and coarse search (default). • fine - Fine search of the peak top with reduced search range.
Disp	Keyword for display mode.
<i>d</i>	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>Without any explicit input, the default mode is set to on. It can be altered with the SetDefault() instruction.</p>
Prot	Keyword for log mode.
<i>p</i>	Log mode. <ul style="list-style-type: none"> • [default] - Default log mode. • off - Log function off. • on - Printer log on. <p>See Chapter 8.4.3 Log output, 422.</p> <p>If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.</p>
SaveMode	Keyword for save mode.
<i>m</i>	Save mode. <ul style="list-style-type: none"> • [append] - Add the additional basic points, i.e. mass numbers to the end of the data file (default). • reset - Clear the file before storing the data.

Example(s):

```
MSC(Par="mscal")
```

Perform a mass scale calibration using the "mscal" parameter file.

```
MSC(Par="c:\qs\par\mscal",Disp=off,SaveMode=reset)
```

Perform a mass scale calibration with the "c:\qs\par\mscal" parameter file without display. Clear the old mass numbers before storing the new adjusted mass numbers.

8.7.42**OffsetCal() - QMS Offset calibration****Description:**

Perform an offset calibration on the QMS 422/200 (not QMS 421!).

Syntax:

```
OffsetCal([Clear=onoff],[DetType=type[,SemVolt=volt]_
          ZeroMass=mass[,Resol=resol][,Disp=d][,Prot=p])
```

Parameters:

Clear Keyword for clear QMS 422 Offset-values (does not apply to QMS 200).

- onoff*
- **on** - Clear offset values of QMS 422
 - **off** - Don't clear offset values of QMS 422

DetType Keyword for detector type.

- type*
- Detector type
- **FARADAY** - Faraday
 - **SEM** - SEM
 - **CH-TRON** - Channeltron

SemVolt Keyword for SEM voltage of the channeltron (only QMS 200)

- volt*
- SEM voltage
- **common** - Common SEM voltage
 - *cc* - Direct input of a numeric value (integer or float).
 - **i[xx]** - Element *xx* of the local integer array.
 - **gi[xx]** - Element *xx* of the global integer array.

ZeroMass Keyword for zero mass.

- mass*
- *cc* - Direct input of the numeric value (Integer or Float).
 - **i[xx]** - Element *xx* of the local integer array.
 - **gi[xx]** - Element *xx* of the global integer array.
 - **f[xx]** - Element *xx* of the local float array.
 - **gf[xx]** - Element *xx* of the global float array.
 - **gfa[xx][yy]** - Element *xx/yy* of the global float field.

Resol Keyword for resolution (only QMS 200)

- resol*
- Resolution (**off**, 1 ... 255)
- **off** - Resolution switched off
 - *cc* - Direct input of the numeric value (integer or float)
 - **i[xx]** - Element *xx* of the local integer array.
 - **gi[xx]** - Element *xx* of the global integer array.

Disp Keyword for display mode.

- d** Display mode.
- **[default]** - Default display mode.
 - **off** - Display off.
 - **on** - Display as table.
- Without an explicit entry the default mode is set to on. It can be changed with the SetDefault() instruction.
- Prot** Keyword for log mode
- p** Log mode.
- **[default]** - Default log mode.
 - **off** - Log function off.
 - **on** - Printer log on.
- See Chapter 8.4.3 Log output, 422.
- If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.

Example(s):

```
OffsetCal(DetType=CH-TRON, SemVolt=common, ZeroMass=5.5,
Resol=25, Disp=off)
```

Perform an offset calibration on a QMS 200 Prisma™ (detector type channeltron, common SEM voltage, on mass number 5.5 with resolution 25) without displaying the measurement results.

```
OffsetCal(DetType=SEM, ZeroMass=f[2], Disp=on, Prot=on)
```

Perform an offset calibration on a QMS 422, on the mass number stored in f[2]. Display the measurement results and output them on the printer.

```
OffsetCal(Clear=on)
```

Set all offset calibration values of a QMS 422 to zero (delete them).

8.7.43 Process() - Data processing

Description:

This instruction is used for selecting different data (ion currents, concentrations, calculated data, etc.) and to make them available as input for further processing:

What? (**DataGfa**, **DataI**, **DataF**) -> Where? (**Disp**, **Prot**, **SaveGfa**, **SaveCyc**)

- Select data

The three keywords **DataGfa**, **DataI** and **DataF** are used for selecting the arrays (from the float data field **gfa[][]**), or for selecting the array elements from the integer arrays (**i[]**, **gi[]**) and/or from the float arrays (**f[]**, **gf[]**, **gfa[][]**).

- Processing facilities

The selected data are assembled in the order shown below and processed according to Disp, Prot (displayed, printed) or saved in the array specified under SaveGfa and/or in the file specified under SaveCyc.

Syntax:

```
Process ([,DataGfa=arr1[;arr2][;arr3]... ]_
        [,DataI=var1[;var2][;var3]... ]_
        [,DataF=var1[;var2][;var3]... ]_
        [,Disp=d][,Prot=p][,SaveGfa=var]_
        [,SaveCyc=file[;ind][;m]][,SaveRef=file [;ind]])
```

Parameters:

DataGfa Keyword for selecting float arrays from the float data field **gfa[xx][]**.

arr1;arr2;... Numbers of the arrays

- **cc** - Numeric value (integer).
- **i[xx]** - Element xx of the local integer array.
- **gi[xx]** - Element xx of the global integer array.

NOTE:

Only those elements of the selected float array are processed which have a name.

=> The gfa[] elements are named by Calculate() and SetVar() or the name is given automatically for measurement data by component and mass number.

Datal Keyword for selecting variables from the integer array.

var1;var2;... Identifier of variables (array elements).

- **i[xx]** - Element xx of the local integer array.
- **gi[xx]** - Element xx of the global integer array.

NOTE:

All selected array elements are processed. If a variable is unnamed, the designation of the array element (for example “i[4]”, “gi[10]” etc.) is used.

=> The i[] and gi[] elements are named via the Calculate() and SetVar() instructions.

DataF Keyword for selecting individual float variables:

var1;var2;... Identifier of float variables (array elements).

- **f[xx]** - Element xx of the local float array.
- **gf[xx]** - Element xx of the global float array.
- **gfa[xx][yy]** - Element xx/yy of the global float field.

NOTE:

All selected array elements are processed. If a variable is unnamed, the designation of the array element (for example, “f[3]”, “gf[5]” etc.) is used.

=> The gf[], gf[]- and gfa[] elements are named by Calculate() and SetVar(), or the name is given automatically for measurement data by component and mass number.

Disp Keyword for display mode

d Display mode

- **[default]** - Default display mode.
- **off** - Display off.
- **on** - Display on.

Without any explicit input, the default mode is set on. It can be altered with the SetDefault() instruction.

Prot	Keyword for log mode.
<i>p</i>	<p>Log mode.</p> <ul style="list-style-type: none"> • [default] - Default log mode. • off - Log function off. • on - Printer log on. <p>See Chapter 8.4.3 Log output, 422.</p> <p>If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.</p>
SaveGfa	Keyword for the array in which the selected data are to be stored.
<i>var</i>	<p>Number of the float array from the global data field gfa[[]].</p> <ul style="list-style-type: none"> • <i>cc</i> - Numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array. <p>NOTE: All data stored in the elements yy=0 ... 63 of the float array gfa[xx][yy] (xx=var) are cleared before the save operation..</p>
SaveCyc, SaveRef	Keywords for defining files in which the selected data are to be stored as cycles (SaveCyc) and/or as a reference (SaveRef).
<i>File</i>	<p>(Drive:\Path) <i>Filename</i>.</p> <ul style="list-style-type: none"> • "ffffff" <i>Filename</i> (e.g. "c:\qs\myfile"). • gs[xx] - Element <i>xx</i> of the global character string array, if the filename is stored in gs[[]].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
<i>m</i>	<p>Save mode</p> <ul style="list-style-type: none"> • [append] - Add the data to the end of the file (default). • reset - Clear the file before storing the data. <p>NOTE: In append mode no test is made whether or not the previously stored data are compatible with the new data. Since Process() creates the data dynamically, a compatibility test would have to be conducted on-line each time a data item is stored. This would require a lot of time and might adversely affect the measurement speed.</p>

Example(s):

```
Process(DataGfa=0, Disp=tab)
```

Display a table of the data stored in float array 0.

```
Process(DataGfa=4,Disp=off,SaveRef="ref3.mdr")
```

Save the data of float array 4 in the reference file "ref3.mdr" without displaying them.

```
Process(DataGfa=0;1,i[3],DataI=i[1];gi[3];gi[4],Disp=on,_,_
      SaveCyc="c:\qs\dat\mix")
```

Display data from float array 0, 1 and i[3] together with the integer variables i[1], gi[3] and gi[4]; store it all in the "c:\qs\dat\mix" Cycle data file.

8.7.44 QMSConnect() - Switch on the communication to the mass spectrometer

Description:

Switch on the communication from the PC to the selected mass spectrometer. If several mass spectrometers are online, the last used one will be switched on. If you want to switch on a different one, you have to do that by SetUnit().

Syntax:

```
QMSConnect ( )
```

Parameters:

None

8.7.45 QMSDisconnect() - Switch off the communication to the mass spectrometer

Description:

Switch off the communication from the PC to the selected mass spectrometer. If several mass spectrometers are online, the last used one will be switched off. If you want to switch off a different one, you have to select it first by SetUnit() and then you can switch it off by QMSDisconnect().

Syntax:

```
QMSDisconnect ( )
```

Parameters:

None

8.7.46 ScanAnalog() - Scan Analog

Description:

Perform a Scan Analog measurement.

Syntax:

```
ScanAnalog(Par=pFile[ ;ind]_
            [,DataGfa=arr1[ ;arr2][ ;arr3]...]_
            [,DataI=var1[ ;var2][ ;var3]...]_
            [,DataF=var1[ ;var2][ ;var3]...]_
            [,Disp=d][ ,SaveRef=file[ ;ind]]_
            [,SaveCyc=file[ ;ind][ ;m]][ ,RefFile=rFile[ ;ind]])
```

Parameters:

Par	Keyword for parameter file.
<i>pFile</i>	(Drive:\Path\) filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416

DataGfa, DataI, DataF, Keywords for additional data. Those data are stored additionally and may be displayed later in **Dispsav** as a table.

Meaning and function of **DataGfa, DataI and DataF** are explained in Chapter 8.7.43 Process() - Data processing, 482

Disp Keyword for Display mode.

d Display mode.

- **[default]** - Default display mode.
- **off** - Display off.
- **on** - Display on.

Without any explicit input, the default mode is set on. It can be altered with the SetDefault() instruction.

SaveRef Keyword for storing the reference

File (Drive:\Path\)*Filename*

- "ffffff" - Filename (e.g.. "c:\qs\myfile")
- **gs[xx]** - Element xx of the global character string array if the filename is stored in **gs[]**.

ind File index - See Chapter 8.4.1.3 File names with indices, 416

SaveCyc Keyword for storing the measurement data of a scan. The parameter *m* provides the possibility to append data to an existing file.

File (Drive:\Path\)*Filename*

- "ffffff" - Filename (e.g. "c:\qs\myfile")
- **gs[xx]** - Element xx of the global character string array if the filename is stored in **gs[]**

ind File index - See Chapter 8.4.1.3 File names with indices, 416

m Save mode

- **[append]** - Add the data to the end of the file (default).
- **reset** - Clear the file before storing the new data.

RefFile Keyword for reference file whose data are to be subtracted from the currently measured data.

rFile (Drive:\Path\)*Filename*

- "ffffff" Filename (e.g.. "c:\qs\myrfile")
- **gs[xx]** - Element xx of the global character string array if the filename is stored in **gs[]**.

NOTE:

If you declare a reference file, then the data delivered by ScanAnalog are always the difference between the currently measured and the reference. In other words, it's always the difference that is stored or displayed, not the actually measured data.

Example(s):

```
ScanAnalog(Par="anascan1 ")
```

Perform a Scan Analog measurement with the parameter file "anascan1".

```
ScanAnalog(Par="c:\qs\par\scan1 ";SaveCyc="c:\qs\dat\scadata ")
```

Perform a Scan Analog measurement with the parameter file "c:\qs\par\scan1". Saving of the data in the file "c:\qs\dat\scadata" is activated.

8.7.47

ScanBar() - Scan Bargraph

Description:

Perform a Scan Bargraph measurement.

Syntax:

```
ScanBar(Par=pFile[;ind][,Disp=d][,Prot=p]_
        [,SaveCyc=file[;ind][;m]][,RefFile=rFile[;ind]])
```

Parameters:

Par	Keyword for parameter file.
<i>pFile</i>	(Drive:\Path\) Filename. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Disp	Keyword for display mode.
<i>d</i>	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>Without any explicit input, the default mode is set on. It can be altered with the SetDefault() instruction.</p>
Prot	Keyword for log mode.
<i>p</i>	Log mode. <ul style="list-style-type: none"> • [default] - Default log mode. • off - Log function off. • on - Printer log on. <p>See Chapter 8.4.3 Log output, 422.</p> <p>If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.</p>
SaveCyc	Keyword for storing the found peaks of a scan as Cycle data in a file.
<i>File</i>	(Drive:\Path\) Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
<i>m</i>	Save mode <ul style="list-style-type: none"> • [append] - Add the new data to the end of the file. • reset - Clear the file before storing the new data.

RefFile Keyword for reference file whose data are to be subtracted from the currently measured data.

rFile (Drive:\Path\)
Filename

- "ffffff" - Filename (e.g.. "c:\qs\myrfile")
- **gs[xx]** - Element xx of the global character string array if the filename is stored in **gs[]**.

NOTE:

If you declare a reference file, then the data delivered by Scan Bargraph are always the difference between the currently measured and the reference. In other words, it's always the difference that is stored or displayed, not the actually measured data.

Example(s):

```
ScanBar ( Par="measure" )
```

Perform a Scan Bargraph measurement with the "measure" parameter file.

```
ScanBar ( Par="c:\qs\par\scan1" , Prot=on , _
          SaveCyc="c:\qs\dat\sbcbyc" )
```

Perform a Scan Bargraph measurement with the "c:\qs\par\scan1" parameter file. Log messages on the printer and save the data in file "c:\qs\dat\sbcbyc".

8.7.48

SeqExit() - Defining a sequence for sequence cancellation

Description:

Define a sequence to be executed when the current sequence is terminated or aborted:

- Each sequence can be programmed with a cancellation sequence. If a subsequence called by the **Sequence()** instruction is cancelled, then the **SeqExit()** instruction belonging to the subsequence is executed first. After that, the **SeqExit()** instruction of the main sequence is executed.
- A cancellation sequence can be locked, so it can't be cancelled itself.
- By calling **SeqExit()** again, a different sequence can be defined as the cancellation sequence. Only the last defined cancellation sequence is valid.

Syntax:

```
SeqExit ( Par=file[ ; ind] , CloseLock=OnOff )
```

Parameters:

Par Keyword for parameter file (sequence).

file (Drive:\Path\)
Filename

- "ffffff" - Filename (e.g. "c:\qs\myfile").
- **gs[xx]** - Element xx of the global character string array, if the sequence name is stored in **gs[]** .

ind File index - See Chapter 8.4.1.3 File names with indices, 416

CloseLock Keyword to lock the cancellation sequence chosen by SeqExit.

- **off** - The cancellation sequence is not locked, it can therefore be cancelled.
- **on** - The cancellation sequence is being executed completely under all circumstances, it cannot be cancelled.

Example(s):

```
SeqExit(Par="close")
```

Execute the "close"-sequence when the sequence is cancelled. "Close" may be cancelled for its part.

```
SeqExit(Par="close",CloseLock=on)
```

The sequence "close" is being executed on termination or cancellation of the sequence and cannot be cancelled.

8.7.49 Sequence() - Starting a (sub-)sequence

Description:

Call another sequence. When this subsequence has terminated, the current sequence is continued.

Syntax:

```
Sequence(Par=file [;ind])
```

Parameters:

Par Keyword for parameter file (sequence).

file (Drive:\Path\)*Filename of the sequence.*

- "ffffff" - Filename (e.g. "c:\qs\myfile").
- **gs[xx]** - Element *xx* of the global character string array, if the sequence name is stored in **gs[]**.

ind File index - See Chapter 8.4.1.3 File names with indices, 416

Example(s):

```
Sequence(Par="inlet")
```

Start the "inlet" sequence.

```
Sequence(Par="c:\par\calib")
```

Start the "c:\par\calib" sequence.

8.7.50

SetAO() - Setting an analog output channel**Description:**

Set one or all analog output channels of the control unit to a certain value.

Syntax:

```
SetAO(chan=val{ ;unit },Mode=m )
```

Parameters:

chan	AO channel (1 ... 12). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • "zzzzz" - Logical name AO#1 ... AO#12 of the output channel. • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array.
<i>val</i>	Value to be set (-10.24 ... +10.238) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of the numeric value (integer or float) • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>f[xx]</i> - Element <i>xx</i> of the local float array. • <i>gf[xx]</i> - Element <i>xx</i> of the global float array. • <i>gfa[xx][yy]</i> - Element <i>xx/yy</i> of the global float field.
<i>unit</i>	Identification of the control unit, only available in multiplex mode. <ul style="list-style-type: none"> • "zzzzz" - Logical name of the control unit. • <i>cc</i> - Direct input of the internal logical number. <p>=> Only those control units are available whose status has been set to Enable in the Parameter Setup program under [Comm] > [Node Identification]. This entry does not influence the unit selected with SetUnit().</p>
Mode	AO output mode
<i>m</i>	Mode <ul style="list-style-type: none"> • [set] - Set the AO voltage to the predefined value. • inc - Increase the AO voltage by the specified amount. If it exceeds the maximum admissible voltage, it is limited. • dec - Reduce the AO voltage by the specified amount. If it exceeds the minimum admissible voltage, it is limited.

NOTE:

The voltage effectively appearing at the analog output port is normally not the given float number in volts. This number is converted to a voltage by using the AO-Characteristic Curve ([Parameter Setup] > [Config] > [AO Characteristic Curve]).

Example(s):

```
SetAO(3=f[12])
```

Set analog output 3 to the value stored in f[12].

```
SetAO(AO#12; "Q_Unit_1"=gf[18])
```

Multiplex mode: Sets AO#12 of the unit "Q_Unit_1" to the value gf[18].

```
SetAO(i[8]=0.1,Mode=inc)
```

Increase the analog output voltage nbr i[8] by 0.1.

8.7.51 SetDefault() - Set default values

Description:

Change the default values.

Set the default values (Disp, Prot, etc.), which can be recalled by [default] in the parameter list of the individual sequencer instructions. This simplifies e.g. the logging of the current measurement, since you can set the optional prot parameter with SetDefault(prot=on). This way, prot doesn't have to be set in each instruction.

Syntax:

```
SetDefault([Disp=d][,Prot=p])
```

Parameters:

Disp Keyword for display mode.

d Display mode.

- **off** - Display off.
- **tab** - Table display.
- **bar** - Bargraph display.
- **vt** -Versus Time display.

Prot Keyword for log mode.

p Log mode.

- **off** - Log function off.
- **on** - Printer log on.

See Chapter 8.4.3 Log output, 422.

Example(s):

```
SetDefault(Prot=on).
```

Set Prot default value to **on**. Every sequencer instruction whose Prot parameter is set to default (or not specified), sends its results as a log to the printer.

8.7.52 SetDO() - Set a digital output channel

Description:

Switch one or all digital output channels of the control unit on or off. Delay the sequence continuation or display the DO's (optionally).

Syntax:

```
SetDO(num{;unit}=state[;tim[;unit]][,Disp=d])
```

Parameters:

<i>num</i>	<p>Number of the digital output channel 0 ... 31 (63, 95)</p> <ul style="list-style-type: none"> • “zzzzz” - Logical name of the output channel. (=> The logical names are defined in the Parameter Setup program under [Config] > [DI/DO Identifications].) • <i>cc</i> - Direct input of a numeric value (integer). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>all</i> - All available digital output channels.
<i>unit</i>	<p>Identification of the control unit, only in multiplex mode.</p> <ul style="list-style-type: none"> • “zzzzz” - Logical name of the control unit. • <i>cc</i> - Direct input of the internal logical number. <p>=> Only those control units are available whose status has been set to Enable in the Parameter Setup program under [Comm] > [Node Identification]. This entry does not influence the unit selected with SetUnit().</p>
<i>state</i>	<p>DO-Status</p> <ul style="list-style-type: none"> • off - Set DO status off. • on - DO status on.
<i>tim</i>	<p>Time specification (1 ... 1'000'000)</p> <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array.
<i>unit</i>	<p>Unit of time.</p> <ul style="list-style-type: none"> • ms - Time in milliseconds. • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
Disp	<p>Keyword for display</p> <p>Display the sequencer status</p> <p>See Chapter 2.7.2.2 Status display, 122</p>
<i>d</i>	<p>Display mode</p> <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p> <p>Select the channels to be displayed in the Measure program under [Sequence] > [Setup].</p>

Example(s):

```
SetDO(10;"Q_Unit_1"=on)
```

Multiplex mode: Set the digital output 10 of unit “Q_Unit_1” to on.

```
SetDO(all=off;1;min)0
```

Switch all digital output channels off, and perform a 1 minute delay.

8.7.53

SetFileInfo() - Changing the measure data documentation settings**Description:**

Changing the File Info settings (See Chapter 2.11 Setup, 148) while running a sequence. When a sequence is started, the settings made in that setup box are valid. By using the instruction **SetFileInfo()**, those settings can be changed during sequence execution.

Syntax:

```
SetFileInfo( [ ,SaveFileInfo=mode ] [ ,SaveMeaPar=mode ]_
             [ ,SaveISPar=mode ] [ ,SaveStatePar=mode ]_
             [ ,SaveZeroDat=mode ] [ ,ResetNames=mode ]_
             [ ,ResetValues=mode ] [ ,FileTitle=title ] [ ;ind ]_
             [ ,Name=num;name ] [ ,Name=... ] [ ,Value=num;val ]_
             [ ,Value=... ] )
```

Parameters:

SaveFileInfo Keyword for generally switching the file information storage on/off. If this switch is on, the items to be stored (measurement parameters, Ion Source parameters etc.) may be chosen; if it is off, no documentation is stored at all.

mode Switch the file information storage function on/off

- **off** - File information storage function disabled.
- **on** - File information storage function enabled.

SaveMeaPar Keyword for switching the measurement parameter storage on/off.

mode Switch the measurement parameter storage function on/off

- **off** - Measurement parameter storage function disabled.
- **on** - Measurement parameter storage function enabled.

SaveISPar Keyword for switching the ion source parameter storage on/off.

mode Switch the ion source parameter storage function on/off.

- **off** - Ion source parameter storage function disabled.
- **on** - Ion source parameter storage function enabled.

SaveStatePar Keyword for switching the status parameter storage on/off.

mode Switch the status parameter storage function on/off.

- **off** - Status parameter storage function disabled.
- **on** - Status parameter storage function enabled.

SaveZeroDat Keyword for switching the zero gas data storage on/off.

mode Switch the zero gas data storage function on/off.

- **off** - Zero gas data storage function disabled.
- **on** - Zero gas data storage function enabled.

ResetNames Keyword to delete additional information names.

mode Delete additional information names on/off

- **off** - Don't delete names.
- **on** - Delete all names.

ResetValues Keyword to delete additional information values.

<i>mode</i>	Delete additional information values on/off <ul style="list-style-type: none"> • off - Don't delete values. • on - Delete all values.
FileTitle	Keyword for title of the measure data file.
<i>title</i>	Title text. <ul style="list-style-type: none"> • <i>"ttttt"</i> - Direct input of the title text. • gs[xx] - Element xx of the global character string array.
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Name	Keyword for additional information name.
<i>num</i>	Index number of the name (0 ... 19). <ul style="list-style-type: none"> • cc - Direct input of the numeric value (integer) • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>name</i>	Name. <ul style="list-style-type: none"> • <i>"ttttt"</i> - Direct text input of the name. • gs[xx] - Element xx of the global character string array, if the name is stored in gs[].
Value	Keyword for additional-information value.
<i>num</i>	Index number of the value (0 ... 19). <ul style="list-style-type: none"> • cc - Direct input of the numeric value (integer) • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>val</i>	Value. <ul style="list-style-type: none"> • <i>"ttttt"</i> - Direct text input of the value. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field. • gs[xx] - Element xx of the global character string array.

Example(s):

```
SetFileInfo(SaveFileInfo=on,SaveISPar=off)
```

Store file information. No Ion Source parameters shall be stored.

```
SetFileInfo(Name=0;"Probe ID",ResetValues=on,Value=0;i[1],_
Value=gi[5];gf[5])
```

Set name nbr. 0 to "Probe ID" and delete all values. Then set value nbr. 0 to the content of i[1] and value nbr. gi[5] to the content of gf[5].

8.7.54**SetISPar() - Set Ion Source Parameters of the current unit****Description:**

Set or change parameter(s) on the ion source of the current unit. This may be useful for automatically optimizing the ion source or adapting it to varying conditions. The new values can optionally be stored on the hidden file that's assigned to this ion source; the old values will be overwritten then.

Syntax:

```
SetISPar( [Filament=filnbr] [ ,ISSet=setnbr] [ ,Emission=emicur ]_
[ ,Protection=protcur] [ ,V1=V1] [ ,V2=V2] [ ,V3=V3 ]_
[ ,V4=V4] [ ,V5=V5] [ ,V6=V6] [ ,V7=V7] [ ,V8=V8] [ ,V9=V9 ]_
[ ,RFPolarity=rfpol] ,Save=onoff )
```

Parameters:

Filament Keyword for filament to use.

filnbr nbr of the filament.

- *cc* - Direct input of a numeric value (integer).
- **i**[*xx*] - Element *xx* of the local integer array.
- **gi**[*xx*] - Element *xx* of the global integer array.

ISSet Keyword for the number of the Ion Source Set to be used

setnbr nbr of the ISSet

- *cc* - Direct input of a numeric value (integer).
- **i**[*xx*] - Element *xx* of the local integer array.
- **gi**[*xx*] - Element *xx* of the global integer array.

Emission Protection These are the Ion Source parameters that may be changed on the ion source of the currently set unit. Whether a parameter is possible/useful, what its possible range is, etc., depends very much on the ion source in use, its configuration etc.

V1
V2
V3
V4
V5
V6

Use this instruction only if you are sure, that the parameter values you intend to set are correct. Quadstar 32-bit has no possibility to check their consistency.

V7
V8
V9

This instruction can cause malfunction, if it is used without comprehensive knowledge of the meaning and relations between the parameters.

RFPolarity

Save Keyword for 'saving parameters to file yes/no'.

mode Save the new values or not.

- **off** - Use the new values for this sequence, leave file unchanged.
- **on** - Save the new values to file, overwrite the old ones.

NOTE:

Make sure that no Ion Source Sets are overwritten that you may need in further measurements.

Example(s):

```
Loop(i[0]=105;150;5)
Begin
  SetISPar(Filament=1, V1=i[0])
  ScanAnalog(Par="c:\QS50\par\temp.sap", Disp=on)
End
```

Ten Scan Analog measurements are carried out, using filament 1, while the ion source voltage V1 (IonRef) is varied from 105 volts to 150 volts, 5 volts each time. The scans will show the influence of V1 on the peak shape and height. The values for V1 are not stored to file, so it will be the same again (as before this sequence has been run), if you call the Tune Up program and there [Tune] > [Ion Source].

8.7.55

SetPar() - Setting Quadstar 32-bit parameters**Description:**

Change one or several Quadstar 32-bit parameters. The parameters correspond mainly to the display control parameters found in the program Measure under [Scan/MID/MCD] > [Setup].

NOTE:

If a keyword is not used, the last settings made in the corresponding program are taken as default settings.

Syntax:

```
SetPar([AlarmDisp=mode][,AlarmError=mode][,AlarmSound=mode]_
[,BarWidth=var][,ColorMode=mode][,Cycles=var]_
[,DispLast=mode][,DispOpt=mode][,LineType=type]_
[,Marker=mode][,MSC=mode][,NorMode=mode][,RelVal=mode]_
[,RollAxis=mode][,Scale=state][,TimeScale=tim[;unit]]_
[,XLetter=mode][,XRaster=mode][,XScale=mode]_
[,YRaster=mode][,ZeroEq=mode][,ZeroSub=mode]_
[,DDETimeout=tmt])
```

Parameters:

AlarmDisp	Keyword for graphic display of the threshold monitoring.
<i>mode</i>	Display the alarm status for all channels with enabled threshold monitoring (can be set in the Parameter Setup program under [Measure] > [Scan/MID/MCD] > [Trip]). <ul style="list-style-type: none"> • off - Alarm display off. • on - Alarm display on.
AlarmError	Keyword for threshold alarm message.
<i>mode</i>	Display threshold alarm in the Event window on/off. <ul style="list-style-type: none"> • off - Alarm error message off. • on - Alarm error message on.
AlarmSound	Keyword for audible alarm signal.
<i>mode</i>	Turn on the audible alarm warning. <ul style="list-style-type: none"> • off - Alarm sound off. • on - Alarm sound on.
BarWidth	Keyword for bar width of the bargraph.
<i>var</i>	Bar width of the bargraph in pixels (1 ... 255) <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
ColorMode	Keyword for color settings of the measure data display.
<i>mode</i>	Display mode single color/multiple color. <ul style="list-style-type: none"> • s - Single color. The color can be chosen in the Parameter Setup program under [Config] > [Screen Color]. • m - Multiple color. The measurement data are displayed in the colors defined in the channel parameter editor. <p>Each of the bargraph elements (Normal, Last) should be marked by its own grid (see Chapter 2.5.7.1 Parameter setup, 96, Chapter 2.6.8.1 Parameter setup, 113)</p>

Cycles	Keyword for number of cycles. Specifies the number of cycles to be displayed in the measurement data window before the old data are overwritten.
<i>var</i>	Number of cycles (2 ... 10000) <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer) • i[xx] - Element xx of the local integer array • gi[xx] - Element xx of the global integer array
DispLast	Keyword for Display Last Values for the bargraph display in MID and MCD mode.
<i>mode</i>	Display Last <ul style="list-style-type: none"> • off - Display Last off. • on - Display Last on.
DispOpt	Keyword for table display option.
<i>mode</i>	Display option <ul style="list-style-type: none"> • none - Optional display switched off. • integral - Display the measurement signal integrated since the start of the measurement. • <i>last</i> - Display the values of the last measurement. • MultColumn - Display the measurement data in several columns.
LineType	Keyword for line type.
<i>type</i>	Line type <ul style="list-style-type: none"> • solid - The measured values or markers are interconnected with a solid line. • dot - The measured values or markers are displayed as dots.
Marker	Keyword for marker display. Mark a measured value with a symbol. The appropriate symbol is selected automatically.
<i>mode</i>	Marker display. <ul style="list-style-type: none"> • off - Marker display off. • on - Marker display on.
MSC	Keyword for mass scale correction.
<i>mode</i>	Mass scale correction. <ul style="list-style-type: none"> • off - Mass scale correction off. The measurement is performed on the mass selected in the channel parameter editor. • on - Mass scale correction on. The mass selected in the channel parameter editor is corrected to the mass found during the mass scale calibration.
NorMode	Keyword for normalization mode.
<i>mode</i>	Normalization mode: <ul style="list-style-type: none"> • sum100 - Sum of the measured values normalized to 100%. • max100 - Largest measured value normalized to 100%. • extern - External normalization.
RelVal	Keyword for displaying relative values.
<i>mode</i>	Display relative values: <ul style="list-style-type: none"> • none - Display relative values. • sum100 - Sum of all values = 100%. • max100 - Maximum value = 100%.

RollAxis	Keyword for roll axis mode. Shifts the display to the left when the end of the X-axis is reached (unit = cycle or time).
<i>mode</i>	Roll axis mode <ul style="list-style-type: none"> • off - Roll axis off. • on - Roll axis on.
Scale	Keyword for axis mode.
<i>state</i>	Axis mode lin/log <ul style="list-style-type: none"> • lin - Linear axes. • log - Logarithmic axes.
ScanArran	Keyword for multiple scan representation.
<i>mode</i>	Axis mode vertical/horizontal. <ul style="list-style-type: none"> • Vert - Vertical, below each other. • Horiz - Horizontal, side-by-side.
TimeScale	Keyword for time scale extent
<i>tim</i>	Time specification (1 ... 1'000'000). <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
XLetter	Keyword for X-axis lettering.
<i>mode</i>	X-Axis lettering <ul style="list-style-type: none"> • off - X-axis lettering off • on - X-axis lettering on.
XRaster	Keyword for displaying the X-grid in the measurement picture
<i>mode</i>	X-grid display <ul style="list-style-type: none"> • off - X-grid display off. • on - X-grid display on.
XScale	Keyword for displaying the X-scale.
<i>mode</i>	X-scale display <ul style="list-style-type: none"> • cycle - Cycle. • timrel - Time relative. • timabs - Time absolute.
YRaster	Keyword for displaying the Y-grid in the measurement picture.
<i>mode</i>	Y-grid display. <ul style="list-style-type: none"> • off - Y-grid display off. • on - Y-grid display on.
ZeroEq	Keyword for zero equalization of the amplifier offset on the 'Zero Mass'. <ul style="list-style-type: none"> • When a QMS 422/421 is used, the zero equalization is switched on by default. When a sequence is started, a zero equalization is performed before the first measurement is started. • When a QMS 200 is used, this switch has no effect (use the instruction OffsetCal() instead).

<i>mode</i>	Zero equalization <ul style="list-style-type: none"> • off - Zero equalization on sequence start off. • on - Zero equalization on sequence start on.
ZeroSub	Keyword for zero gas subtraction. => This parameter is available only in the sequencer and (like in manual mode) always on when the sequence is started.
<i>mode</i>	Zero gas subtraction. <ul style="list-style-type: none"> • off - Zero gas subtraction off. • on - Zero gas subtraction on. The zero gas measurement signals obtained during the zero gas calibration are subtracted from the measured signals.
DDETimeOut	Keyword for DDE timeout. If data are exchanged by DDE using a network, the default setting may be too short, so the connection cannot be established.
<i>tmt</i>	Time in milliseconds (1 ... 60'000). <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.

Example(s):

```
SetPar(BarWidth=20,NorMode=sum100, Scale=lin)
```

Set the bar width to 20 pixels, the normalization mode to sum=100%, and the Scale parameter to linear (applies to bargraph display mode).

8.7.56**SetParFile() - Set value(s) in a channel parameter file****Description:**

Set or change a value in a channel parameter file. This may be useful if you run e.g. measurement series in which one parameter changes continuously.

Syntax:

```
SetParFile(ParFile=pFile[;ind],Channel=chan[,Mass][,Dwell][,FirstMass][,Width][,Speed][,Steps][,Resolution][,AmplRange][,RangeLimit][,Filter][,Offset][,DispMode][,DispRange][,DispDec])
```

Parameters:

ParFile	Keyword for parameter file.
<i>pFile</i>	(Drive:\Path\)\ Filename. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array, if the file name is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Channel	Keyword for channel number

<i>chan</i>	Number of the measure channel <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • ALL - All channels in this file
Mass	These are the channel parameters that may be changed on the file pFile; they partially exclude each other. The range of the individual parameters cannot be declared universal, they depend very much on the mass spectrometer in use, its configuration etc.
Dwell	
FirstMass	
Width	
Speed	
Steps	Use this instruction only if you are sure, that the parameter values you intend to set are correct. Quadstar 32-bit has no possibility to check their consistency.
Resolution	This instruction can cause malfunction, if it is used without comprehensive knowledge of the relations between the parameters.
AmplRange	
RangeLimit	
FilterOffset	
DispMode	
DispRange	
DispDec	
	NOTE: Copy the desired parameter file to another one and change this copy by SetParFile; this way the original file remains unchanged.

Example(s):

```
Loop(i[0]=0;90;10)
Begin
  SetParFile(ParFile="c:\as\par\temp.sap", Channel=0,_,
            FirstMass=i[0], Width=10)
  ScanAnalog(Par="c:\as\par\temp.sap", Disp=on,_,
            SaveCyc="scan";i[0])
End
```

Ten Scan Analog measurements are carried out, all with a scan width of 10 amu, starting at the masses 0, 10, 20, ..., 90. For that purpose, the parameter FirstMass in the temporary file 'temp.sap' is changed. The data are stored in the files 'scan0000' to 'scan0090'.

8.7.57**SetQMS() - Setting the QMS parameter(s)****Description:**

Set QMS parameters (ion source, total pressure gauge and SEM).

QMS 200 and QMS 422/421 have different assignment rules for the ion source parameters.

QMS 422/421

For each available ion source type, four ion source sets are stored. These sets can be assigned to the two filaments by choosing the filament (Filament) and assigning an ion source set to it (ISSet).

=> Filament 1+2 selects automatic operation, i.e. filament 2 is activated if filament 1 becomes defective. In this case, assignment of a set to a filament (1 or 2) is not possible.

QMS 200

There are totally four ion source sets. Each of them contains one complete set of ion source parameters for filament 1, 2 and 1+2.

ISSet loads the chosen set, the applying ion source parameters are automatically loaded by choosing the Filament.

=> Filament 1+2 means that both filaments are switched on concurrently. The corresponding parameters are loaded from the current ISSet.

Syntax:

```
SetQMS([CommSEM=var[;chan]][,Emission=state[;chan]_
        [,Filament=num][,ISSet=set][,ISType=typ]_
        [,TotGauge=state][,SEM=state])
```

Parameters:

CommSEM	Keyword for common SEM Voltage
<i>var</i>	Common SEM Voltage QMS 422/421: 0 ... 3500, QMS 200: 0 ... 3000 <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>chan</i>	QMU channel (0 ... 7) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Emission	Keyword for emission on/off.
<i>state</i>	Emission status. <ul style="list-style-type: none"> • off - Emission off. • on - Emission on.
<i>chan</i>	QMU channel (0 ... 7) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Filament	Keyword for filament number. (Function: See Description above)
<i>num</i>	Filament number. <ul style="list-style-type: none"> • 1 - Filament 1 • 2 - Filament 2 • 1+2 - Filament 1+2
ISSet	Keyword for ion source set. (Function: See Description above)
<i>set</i>	Ion source set (0 ... 3) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
ISType	Keyword for ion source type.
<i>type</i>	Ion source type. <ul style="list-style-type: none"> • default - Ion source type specified in the configuration • spec+ - Spec+ • spec- - Spec-
TotGauge	Keyword for turning the total pressure gauge on/off.

<i>state</i>	Status <ul style="list-style-type: none"> • off - total pressure gauge switched off. • on - total pressure gauge switched on.
SEM	Keyword for SEM status.
<i>mode</i>	SEM status. <ul style="list-style-type: none"> • off - SEM off • on - SEM on.

Example(s):

```
SetQMS(CommSEM=1800)
```

Set the common SEM voltage to 1800 V.

```
SetQMS(Emission=on;6,SEM=on)
```

Switch on the emission of QMU channel 6 and turn on the SEM.

8.7.58 SetScale() - Define a scale

Description:

Define a scale for display. Several scales may be defined and the desired one may be accessed by its identification number.

Syntax:

```
SetScale(ScaleID=scID[,Mode=m][,ScaleName=name][,Unit=u]_
        [,Min=minV][,Max=maxV][,Dec=numDecs]_
        [,DecMax=topDec])
```

Parameters:

ScaleID	Keyword for scale identification.
<i>scID</i>	Scale identification number, identifies the scale clearly. <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (0 ... 199). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Mode	Keyword for scale mode.
<i>m</i>	Scale mode (linear or logarithmic). <ul style="list-style-type: none"> • lin - linear scale. • log - logarithmic scale.
ScaleName	Keyword for name of the measured quantity.
<i>name</i>	Scale name (e.g. 'Pressure'). <ul style="list-style-type: none"> • <i>Text</i> - Direct text input. • gs[xx] - Element xx of the global character string array.
Unit	Keyword for unit of measure.
<i>u</i>	Scale unit (e.g. 'mbar'). <ul style="list-style-type: none"> • <i>Text</i> - Direct text input. • gs[xx] - Element xx of the global character string array.
Min	Keyword for start value of the linear scale.

<i>minV</i>	Scale start value (applies only to linear scale). <ul style="list-style-type: none"> • "ccc" - Direct numeric input. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.
Max	Keyword for end value of the linear scale.
<i>maxV</i>	Scale end value (applies only to linear scale) <ul style="list-style-type: none"> • "ccc" - Direct numeric input. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.
Dec	Keyword for number of decades of the logarithmic scale.
<i>numDecs</i>	Number of decades to be displayed (applies only to logarithmic scale). <ul style="list-style-type: none"> • cc - Direct input of a numeric value. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
DecMax	Keyword for highest decade of the logarithmic scale.
<i>topDec</i>	Highest decade to be displayed (applies only to logarithmic scale). <ul style="list-style-type: none"> • cc - Direct input of a numeric value. • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.

Example(s):

```
SetScale(ScaleID=0, Mode=log, ScaleName="pressure", _
        Unit="mbar", _Dec=5, DecMax=-3)
```

The scale nbr 0 shall display the pressure between 1E-3 mbar and 1E-8 mbar in logarithmic mode.

```
SetScale(ScaleID=1, Mode=lin, ScaleName="Oxygen", Unit="%", _
        Min=0, Max=25)
```

The scale nbr 1 shall display the oxygen content between 0 and 25% in linear mode.

8.7.59**SetString() - Setting a character string****Description:**

Assign a character string to an array of **gs[]**. SetString() can also be used for copying character strings from one array of **gs[]** to another array of **gs[]**.

Syntax:

```
SetString(var=string[;strval])
```

Parameters:

<i>var</i>	Variable to assign the character string to. <ul style="list-style-type: none"> • gs[xx] - Global character string array.
<i>string</i>	Character string to assign to the variable. <ul style="list-style-type: none"> • "Text" - Direct input of text. • gs[xx] - Element xx of the global character string array.
<i>strval</i>	String value - See Chapter 8.4.1.4 String values, 417

Example(s):

```
SetString(gs[12]="Please close inlet valve")
```

Define the message text.

```
SetString(gs[i[3]]="C:\qs\dat\cycle")
```

Define the cycle filename.

```
SetString(gs[i[0]]=gs[12])
```

Copy gs[12] to gs[i[0]].

8.7.60**SetUnit() - Selecting the mass spectrometer in multiplex mode****Description:**

With SetUnit() you can select a QMS by specifying its logical unit number (not to be confused with the physical node address) or its logical name. The communication to this mass spectrometer is then switched on.

The individual mass spectrometers must first be given a physical node number (node) within the range of 01 ... FF (hex).

NOTE:

PPM 422 (Plasma Process Monitor) works with Unit 0 (QMS 422 with mass scan) and Unit 1 (QMS 422 with ion energy scan) only.

Syntax:

```
SetUnit(Unit=u)
```

Parameters:

Unit	Keyword for selecting a mass spectrometer within the ArcNet network.
<i>u</i>	Unit number within the range of 0 ... 255 or logical name <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value 0 ... 255. • "zzzzz" - Logical name of the Unit.

Programming the individual nodes (define logical names and Unit numbers) is described in Chapter 3.6.3.4.2 Communication via ArcNet, 216.

Example(s):

```
SetUnit(Unit="MS1")
```

Select the mass spectrometer named MS1 of the ArcNet network.


```
SetUnit(Unit=0)
```

Select the mass spectrometer with number 0 of the ArcNet network.

8.7.61 SetVar() - Setup of variables

Description:

Assign a value to a variable.

A name (text comprising up to 12 characters) can optionally be assigned to the variable *var*. This name is subsequently used for identifying the block:

- when processing the data via the Process() instruction
- when the stored data are displayed in the Display Saved Values program.

An individual scale may be assigned to the variable as well.

Syntax:

```
SetVar([ScaleID=scID,][Name=name[;strval],]var=val)
```

Parameters:

ScaleID	Keyword for scale identification.
<i>scID</i>	Number of the scale (defined by the instruction 'SetScale', See Chapter 8.7.58 SetScale() - Define a scale, 502). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (0 ... 199). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array.
Name	Keyword for name.
<i>name</i>	Name of the variable (max. 12 characters). <ul style="list-style-type: none"> • "name" - Direct input of the variable name. • <i>gs[xx]</i> - Element <i>xx</i> of the global character string array.
strval	String value - See Chapter 8.4.1.4 String values, 417
<i>var</i>	Variable to assign a value to (and possibly a name) <ul style="list-style-type: none"> • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>f[xx]</i> - Element <i>xx</i> of the local float array. • <i>gf[xx]</i> - Element <i>xx</i> of the global float array. • <i>gfa[xx][yy]</i> - Element <i>xx/yy</i> of the global float field.
<i>val</i>	Assignment value (any value) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a number (integer or float). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>f[xx]</i> - Element <i>xx</i> of the local float array. • <i>gf[xx]</i> - Element <i>xx</i> of the global float array. • <i>gfa[xx][yy]</i> - Element <i>xx/yy</i> of the global float field.

Example(s):

```
SetVar(i[1]=1200)
```

Set element 1 of the integer array *i[]* to a value of 1200.

```
SetVar(gf[i[1]]=100.0)
```

Set the variable `gf[i[1]]` to a value of 100.

8.7.62 Sound() - Tone output

Description:

Output a tone.

The instruction `Sound()` produces a tone with a certain pitch and determined or undetermined duration using the built in PC-Speaker. A continuous tone can be stopped again by `Sound`.

Syntax:

```
Sound(Mode=soundmode [, Duration=time] , Frequency=hertz)
```

Parameters:

Mode	Keyword for sound mode.
<i>Soundmode</i>	Tonmodus. <ul style="list-style-type: none"> • On - Give out a tone for a certain time. • Ever - Give out a continuous tone. • Off - Stop the continuous tone.
Duration	Keyword for the duration of the tone (only with Mode = On).
<i>time</i>	Time specification (1 ... 1'000'000) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • ms - Time in milliseconds. • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
Frequency	Keyword for pitch in Hertz.
<i>hertz</i>	Frequency specification (1,...,1'000'000). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.

NOTE:

Make sure the value of hertz is within the range of 50 ... 5000; otherwise the tone might be inaudible.

When a continuous tone is output (**Mode = Ever**), the sequence goes on immediately. A tone with determined duration (**Mode = On**) in contrast causes the sequence to wait until the tone ends and then carry out the next instruction.

Example(s):

```
Sound(Mode=on, Duration=3, Frequency=500)
```

Produces a 500Hz tone that lasts three seconds.

```
Sound(Mode=ever, Frequency=i[1])
MID(Par="mypar")
Sound(Mode=off)
```

A tone of the pitch stored in i[1] is output and a MID measurement is performed. When the MID measurement is done, the tone is switched off.

8.7.63 TextClear() - Clear all text

Description:

Clear all text.

All text elements written to this window by TextSet() are deleted. The text window declared by TextInit() and its properties remain present.

Syntax:

```
TextClear()
```

Parameters:

None

Example(s):

See Chapter 8.7.66 TextSet() - Write a text to the text window, 508

8.7.64 TextClose() - Close text window

Description:

Close the opened text window.

The current text window declared by TextInit() is closed. Its content and properties are deleted.

Syntax:

```
TextClose()
```

Parameters:

None

Example(s):

See Chapter 8.7.66 TextSet() - Write a text to the text window, 508

8.7.65 TextInit() - Open a text window

Description:

Open a text window and declare its properties.

The lower part of the current measurement window is declared as a text window, that can be filled by any text elements generated by the sequence. In the background, a graphic can be displayed.

Syntax:

```
TextInit( [Color=BkGndColor, ]Size=WndHeight[ ,Picture=BkGndPict_
          [,XPos=x][ ,YPos=y][ ,Height=h][ ,Width=w] ] )
```

Parameters:

Color	Keyword for background color.
<i>BkGndColor</i>	Number of the color in the palette. <ul style="list-style-type: none"> • off - No background color (transparent) • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Size	Keyword for window size.
<i>WndHeight</i>	Height of the text window in pixel (the width corresponds to the measurement window) <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Picture	Keyword for background graphic.
<i>BkGndPict</i>	Name and path of the background graphic <ul style="list-style-type: none"> • "Pict.bmp" - Direct input of a name. • gs[xx] - Element xx of the global character string array.
XPos, YPos, Height, Width	Keywords for position and size of the background graphic The position of the background graphic relative to the upper left corner of the text window can be defined by the values XPos and YPos, the size by Height and Width. If nothing is defined, the graphic will fill the text window completely.
<i>x, y, h, w</i>	X- and Y- position of the graphic, and height and width. <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.

Example(s):

See instruction TextSet()

8.7.66**TextSet() - Write a text to the text window****Description:**

Write a text to the text window.

Writes a text to the text window declared by TextInit().

Syntax:

```
TextSet( ID=Ident[ ,XPos=x][ ,YPos=y][ ,Color=TextColor]_
         [,Size=TextSize],Text=OutText )
```

Parameters:

ID	Keyword for identification nbr of the text element. The identification nbr identifies the text element clearly.
<i>Ident</i>	Identification nbr of the text element. <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
XPos, YPos	Keywords for position of the text element.
<i>x, y</i>	X- and Y- position of the text element. Default: upper left corner. <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Color	Keyword for text color.
<i>TextColor</i>	Nbr of the color. Default: Quadstar 32-bit standard. <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Size	Keyword for text height.
<i>TextSize</i>	Text height in pixel. Default: 12 pixel. <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Text	Keyword for text to be output.
<i>OutText</i>	Text to be output. Default: nothing. <ul style="list-style-type: none"> • <i>"Text"</i> - Direct input of text. • gs[xx] - Element xx of the global character string array.

Example(s):

In the following example, a MID measurement is performed over 10 cycles and the current cycle nbr is displayed in the lower part of the measurement window. As a special effect, the number steps one position forward with every new cycle. At the end, the text window is erased and the measurement is performed again.

```
// Scale the measurement window and place it at the upper left
// corner
DispWindow(Disp=normal, XPos=0, YPos=0, Width=400, Height=500)

// Initialize the text window 130 pixels high, with background
// picture "marble"
TextInit(Size=130, Picture="c:\windows\marble.bmp")

// Perform 10 MID Cycles and display the cycle nbr on forward
// stepping positions
Loop(i[0]=1;10)
Begin
  // Calculate x position of cycle nbr
  Calculate(i[1] = 30*i[0])
  // Load number to be displayed
  SetString(gs[0] = "" ;i[0])
  // Measurement
  MID(Par="c:\qs\par\airdemo.mip", Disp=vt)
```

```

// Clear old text, set new text
TextClear()
TextSet(ID=0, XPos=30, YPos=50, Color=0, Size=30, _
        Text="This is Cycle Number: ")
TextSet(ID=1, XPos=i[1], YPos=100, Color=0, Size=30, _
        Text=gs[0])
End

// Remove text window, repeat measurement
TextClose()
Loop(i[0]=1;10)
    MID(Par="c:\qs32bit\par\airdemo.mip", Disp=vt)

```

8.7.67 WhileDI() - While loop for monitoring a digital input

Description:

While loop for monitoring a digital input of the control unit. The number of available channels on the...

- QMS 422/421 depends on the number of installed DI 420 boards,
- QMS 200 is limited to two. (DI 0: Extern Protection, DI 1: Extern Start)

Cancellation criteria:

- Execution of the loop continues as long as the condition `var = state` is true and the cancellation conditions are not true.
- Several cancellation conditions can be tested simultaneously: 'count down' (break) and 'time' (breaktime).

Syntax:

```

WhileDI(var{;unit}=state[,Disp=d][[,Break=cnt]_]
        [,BreakTime=tim[;unit]],BreakSeq=file[;ind])

```

Parameters:

<i>var</i>	number of the digital input channel QMS 422/421: 0 ... 31(63), QMS 200: 0 ... 1 <ul style="list-style-type: none"> • "zzzz" - Logical name of the digital input channel. (See Chapter 3.8.2.3 DI/DO identification, 237). • <i>cc</i> - Direct input of a numeric value (integer). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array.
<i>unit</i>	Identification of the control unit, only in multiplex mode. <ul style="list-style-type: none"> • "zzzz" - Logical name of the control unit. • <i>cc</i> - Direct input of internal logical number. <p>=> Only those control units are available whose status has been set to Enable in the Parameter Setup program under [Comm] > [Node Identification]. This entry does not influence the unit selected with SetUnit().</p>
<i>state</i>	DI comparison status <ul style="list-style-type: none"> • off - DI status off. • on - DI status on.

Disp	Keyword for display mode.
<i>d</i>	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>See Chapter 2.7.2.2 Status display, 122.</p> <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p> <p>Select the channels to be displayed in the Measure program under [Sequence] > [Setup].</p>
Break	Keyword for loop cancellation criterion 'count down'.
<i>cnt</i>	Loop counter (0 ... 1'000'000). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
BreakTime	Keyword for the loop cancellation criterion 'time'.
<i>tim</i>	Time specification (1 ... 1'000'000) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
BreakSeq	Keyword for sequence to be performed when the cancellation condition 'count down' or 'time' becomes true.
<i>file</i>	(Drive:\Path\)\ Filename of the cancellation sequence. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element <i>xx</i> of the global character string array, if the sequence name is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416

Example(s):

```
WhileDI(12;"Q_Unit_1"=on)
  Delay(Time=5)
```

Multiplex mode: As long as input 12 of the "Q_Unit_1" is on, wait another 5 seconds.

```
WhileDI(3=on,Break=180,BreakTime=1;h)
  MID(Par="mypar")
```

A MID measurement is performed as long as the digital input 3 is on. The loop will also be cancelled after max. 180 passes or max. 1 hour.

8.7.68

WhileDO() - While loop for monitoring a digital output**Description:**

While loop for monitoring a digital output of the control unit. The number of available channels on the...

- QMS 422/421 depends on the number of installed DO 420 boards,
- QMS 200 is limited to two. (DO 0, DO 1).

Cancellation criteria:

- Execution of the loop continues as long as the condition `var = state` is true and the cancellation conditions are not true.
- Several cancellation conditions can be tested simultaneously: 'count down' (break) and 'time' (breaktime).

Syntax:

```
WhileDO(var{;unit}=state[,Disp=d][[,Break=cnt ]_  
        [BreakTime=tim[;unit],BreakSeq=file[;ind]])
```

Parameters:

<i>var</i>	number of the digital output channel QMS 422/421: 0 ... 31(63, 95), QMS 200: 0 ... 1 <ul style="list-style-type: none"> • "zzzzz" - Logical name of the digital output channel. (See Chapter 3.8.2.3 DI/DO identification, 237). • <i>cc</i> - Direct input of a numeric value (integer). • <i>i</i>[<i>xx</i>] - Element <i>xx</i> of the local integer array. • <i>gi</i>[<i>xx</i>] - Element <i>xx</i> of the global integer array.
<i>unit</i>	Identification of the control unit, only in multiplex mode. <ul style="list-style-type: none"> • "zzzzz" - Logical name of the control unit. • <i>cc</i> - Direct input of the internal logical number. => Only those control units are available whose status has been set to Enable in the Parameter Setup program under [Comm] > [Node Identification]. This entry does not influence the unit selected with SetUnit().
<i>state</i>	DO comparison status. <ul style="list-style-type: none"> • off - DO status off. • on - DO status on.
Disp	Keyword for display mode.
<i>d</i>	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>See Chapter 2.7.2.2 Status display, 122.</p> <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p> <p>Select the channels to be displayed in the Measure program under [Sequence] > [Setup].</p>
Break	Keyword for the loop cancellation criterion 'count down'

<i>cnt</i>	Loop counter (0 ... 1'000'000) <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
BreakTime	Keyword for the loop cancellation criterion 'time'.
<i>tim</i>	Time specification (1 ... 1'000'000). <ul style="list-style-type: none"> • <i>cc</i> Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
BreakSeq	Keyword for sequence to be performed when the cancellation condition 'count down' or 'time' becomes true.
<i>file</i>	(Drive:\Path\)\ Filename the cancellation sequence. <ul style="list-style-type: none"> • "#####" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element <i>xx</i> of the global character string array, if the sequence name is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416

Example(s):

```
WhileDO("inlet4";"Q_Unit_1"=off)
  Delay(Time=5)
```

Multiplex mode: A five second pause is counted down as long as the digital output "inlet4" of Unit "Q_Unit_1" is off.

```
WhileDO(3=on,Break=180,BreakTime=1;h)
  MID(Par="mypar")
```

A MID measurement is performed as long as the digital input 3 is ON. The loop will also be cancelled after max. 180 passes or max. 1 hour.

8.7.69**WhileMeaVal() - While loop for monitoring a measured value****Description:**

While loop for monitoring a measured value on a measure channel of the QMS.

Execution of the loop continues as long as the condition *var cond val* is true and the cancellation conditions are not true.

Several cancellation conditions can be tested simultaneously: 'count down' (*break*) and 'time' (*breaktime*).

Syntax:

```
WhileMeaVal(var cond val ,Par=file[;ind][,Mode=m][,Disp=d]_
  [,Break=cnt][,BreakTime=tim[;unit]])_
  [,BreakSeq=file[;ind]])
```

Parameters:

<i>var</i>	Measure channel number (0 ... 63). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array.
<i>cond</i>	Condition <ul style="list-style-type: none"> • = - equal. • < - less than. • > - greater than. • <= - less or equal. • >= - greater or equal. • <> - not equal.
<i>val</i>	Comparison value All float and integer array elements can be inserted: <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer or float). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array. • <i>f[xx]</i> - Element <i>xx</i> of the local float array. • <i>gf[xx]</i> - Element <i>xx</i> of the global float array. • <i>gfa[xx][yy]</i> - Element <i>xx/yy</i> of the global float field.
Par	Keyword for parameter file
<i>file</i>	(Drive:\Path\)\ Filename. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • <i>gs[xx]</i> - Element <i>xx</i> of the global character string array, if the filename is stored in <i>gs[]</i>.
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Mode	Keyword for measuring mode
<i>m</i>	Measuring mode. <ul style="list-style-type: none"> • [normal] - Wait until the ion current exceeds or drops below a specified value. • integral - Wait until the integral (sum) of the determined ion currents exceeds or drops below a certain value. • delta - Wait until the difference between the momentary ion current to the previously determined ion current exceeds or drops below a certain value.
Disp	Keyword for display mode
<i>d</i>	Display mode <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on (sequencer status) <p>See Chapter 2.7.2.2 Status display, 122.</p> <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p>
Break	Keyword for the loop cancellation criterion 'count down'.
<i>cnt</i>	Loop counter (0 ... 1'000'000). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • <i>i[xx]</i> - Element <i>xx</i> of the local integer array. • <i>gi[xx]</i> - Element <i>xx</i> of the global integer array.

BreakTime	Keyword for the loop cancellation criterion 'time'.
<i>tim</i>	Time specification (1 ... 1'000'000) <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • [s] - Time in seconds (default). • min - Time in minutes. • h - Time in hours. • d - Time in days.
BreakSeq	Keyword for sequence to be performed when the cancellation condition 'count down' or 'time' becomes true.
<i>file</i>	(Drive:\Path\)\ Filename of the cancellation sequence. <ul style="list-style-type: none"> • "#####" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the sequence name is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416

Example(s):

```
WhileMeaVal(i[2]<9E-10,par="d:\qs\par\myfile",Mode=integral)
    Delay(Time=5)
```

As long as the integral (sum of measured values on channel i[2], determined with the parameter file "d:\qs\par\myfile") is less than 9E-10, wait another five seconds.

8.7.70**WhileString() - While loop for monitoring a string****Description:**

While loop for monitoring a string.

Execution of the loop continues as long as the condition string1 cond string2 is true and the cancellation conditions are not true.

Several cancellation conditions can be tested simultaneously: 'count down' (break) and 'time' (breaktime).

Syntax:

```
WhileString(string1 cond string2[,Disp=d][,Break=cnt]_
    [,BreakTime=tim[;unit]][,BreakSeq=file[;ind]])
```

Parameters:

<i>string1</i>	First string to compare; text or element of the global string array gs[] <ul style="list-style-type: none"> • <i>cc</i> - Direct text input • gs[xx] - Element <i>xx</i> of the global string Array.
<i>cond</i>	Condition <ul style="list-style-type: none"> • = - equal to. • < - less than. • > - greater than. • <= - less or equal. • >= - greater or equal. • <> - not equal.
<i>string2</i>	Second string to compare; text or element of the global string array gs[] <ul style="list-style-type: none"> • <i>cc</i> - Direct text input • gs[xx] - Element <i>xx</i> of the global string Array.
Disp	Keyword for display mode.
<i>d</i>	Display mode <ul style="list-style-type: none"> • [default] Default display mode. • off - Display off. • on - Display on. <p>See Chapter 2.7.2.2 Status display, ¶ 122.</p> <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p>
Break	Keyword for the loop cancellation criterion 'count down'.
<i>cnt</i>	Loop counter (0 ... 1'000'000). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
BreakTime	Keyword for loop cancellation criterion 'time'.
<i>tim</i>	Time specification (1 ... 1'000'000). <ul style="list-style-type: none"> • <i>cc</i> - Direct input of a numeric value (integer). • i[xx] - Element <i>xx</i> of the local integer array. • gi[xx] - Element <i>xx</i> of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • [s] - Time in seconds (default) • min - Time in minutes • h - Time in hours • d - Time in days
BreakSeq	Keyword for sequence to be performed when the cancellation condition 'count down' or 'time' becomes true.
<i>file</i>	(<i>Drive:\Path</i>) Filename of the cancellation sequence. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element <i>xx</i> of the global character string array, if the sequence name is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, ¶ 416

The two strings are compared character by character, disregarding capitals. The first two characters that don't match each other in both strings determine the result of the comparison: the character that comes first in alphabet is considered to be 'smaller'. Two

strings are only recognized as equal if all their characters are equal (except capital/non capital).

Examples:

```
WhileString(gs[2]<>"exit")
  MID(Par="midpar")
```

While gs[2] does not contain the text "exit", repeat the MID measurement.

```
WhileString(gs[i[2]]>gs[1],Disp=on,Break=10,BreakTime=60;s,
  BreakSeq="c:\qs\par\seq")
```

```
Begin
  SetDO(5=off)
  Delay(Time=5)
  SetDO(6=on)
  Delay(Time=1)
  SetDO(6=off)
  SetDO(5=on)
  IfMeaVal(3>1E-5,par="testpar")
  SetVar(i[2]=0)
End
```

While gs[i[2]]>gs[1], the gas inlet is to be purged. After max. 10 passes or max. 60 seconds, the while loop shall be cancelled in any case and the "c:\qs\par\seq" sequence executed.

8.7.71

WhileVar() - While loop for monitoring a variable

Description:

While loop for monitoring a variable.

Execution of the loop continues as long as the condition var cond val is true and the cancellation conditions are not true.

Several cancellation conditions can be tested simultaneously: 'count down' (break) and 'time' (breaktime).

Syntax:

```
WhileVar(var cond val[,Disp=d][,Break=cnt]_
  [,BreakTime=tim[;unit]][,BreakSeq=file[;ind]])
```

Parameters:

<i>var</i>	Variable to be tested All float and integer array elements can be inserted: <ul style="list-style-type: none"> • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.
<i>cond</i>	Condition <ul style="list-style-type: none"> • = - equal • < - less than • > - greater than • <= - less or equal • >= - greater or equal • <> - not equal
<i>val</i>	Comparison value All float and integer array elements can be inserted: <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer or float). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.
Disp	Keyword for display mode.
<i>d</i>	Display mode <ul style="list-style-type: none"> • [default] Default display mode. • off - Display off. • on - Display on. <p>See Chapter 2.7.2.2 Status display, 122.</p> <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p>
Break	Keyword for the loop cancellation criterion 'count down'.
<i>cnt</i>	Loop counter (0 ... 1'000'000). <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
BreakTime	Keyword for loop cancellation criterion 'time'.
<i>tim</i>	Time specification (1 ... 1'000'000). <ul style="list-style-type: none"> • cc - Direct input of a numeric value (integer). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
<i>unit</i>	Unit of time <ul style="list-style-type: none"> • [s] - Time in seconds (default) • min - Time in minutes • h - Time in hours • d - Time in days
BreakSeq	Keyword for sequence to be performed when the cancellation condition 'count down' or 'time' becomes true.

- file* (Drive:\Path) Filename of the cancellation sequence.
- "ffffff" - Filename (e.g. "c:\qs\myfile").
 - **gs[xx]** - Element xx of the global character string array, if the sequence name is stored in **gs[]**.
- ind* File index - See Chapter 8.4.1.3 File names with indices, 416

Example(s):

```
WhileVar (gf[21]<45.5)
Begin
  MID(Par="midpar", SaveGfa=0)
  Calculate(gf[21]=gfa[0][0]*srt(2)/10)
End
```

As long as the value of gf[21] is less than 45.5, the MID measurement is performed repeatedly and gf[21] is used for a calculation.

```
WhileVar (i[2]=3, Disp=on, Break=60, BreakTime=10;min, _
  BreakSeq="c:\qs\par\seq")
Begin
  SetDO(5=off)
  Delay(Time=5)
  SetDO(6=on)
  Delay(Time=1)
  SetDO(6=off)
  SetDO(5=on)
  IfMeaVal(3>1E-5, par="testpar")
  SetVar(i[2]=0)
End
```

As long as the variable i[2] has a value of 3, the gas inlet is to be purged. After max. 60 passes or max. 10 minutes, the while loop shall be cancelled in any case and the "c:\qs\par\seq" sequence executed.

8.7.72 ZeroGas() - Zero gas measurement

Description:

Perform a zero gas measurement.

If contaminations of an ultra-pure gas in the lower ppm range (parts per million) need to be determined, the residual gas background in the analysis chamber must be subtracted, otherwise the result will be incorrect.

For additional information on zero gas measurement refer to Chapter 2.8.2 Background measurement (Zero Gas), 137

Syntax:

```
ZeroGas(Par=file[;ind][,Disp=d][,Prot=p][,SaveMode=m])
```

Parameters:

Par	Keyword for parameter file.
<i>file</i>	(Drive:\Path\) Filename. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Disp	Keyword for display mode.
<i>d</i>	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>Without any explicit declaration, the default display mode is on. The default settings can be altered by SetDefault().</p>
Prot	Keyword for log mode.
<i>p</i>	Log mode. <ul style="list-style-type: none"> • [default] - Default log mode. • off - Log function off. • on - Printer log on. <p>See Chapter 8.4.3 Log output, 422.</p> <p>If nothing else is specified, the default log mode is off. It can be altered with the SetDefault() instruction.</p>
SaveMode	Keyword for save mode.
<i>m</i>	Save mode. <ul style="list-style-type: none"> • reset - Clear the data file before storing the data. • [append] - Append the new data to the end of the file.

Example(s):

```
ZeroGas ( Par="zerocal" )
```

Perform a zero gas measurement with the "zerocal" parameter file.

```
ZeroGas ( Par="c:\qs\par\zero1" , Prot=on , SaveMode=reset )
```

Perform a zero gas measurement with the "c:\qs\par\zero1" parameter file, printer log on, and clear the zero gas file before saving the measurement data.

8.8 Process Monitor PPM 422

The sequencer instructions described on the following pages, which all begin with 'PPM...', are exclusively available for the Plasma Process Monitor PPM 422 .

The following standard instructions are replaced by corresponding PPM 422 specific instructions:

Standard instruction	Replaced by PPM instruction
ScanAnalog()	PPMScanMassAnalog()
MID()	PPMMID()

8.8.1 PPMGetCorr() - Get correction value CORR

Description:

Load CORR value and assign it to a variable (e.g. for a calculation).

Syntax:

```
PPMGetCorr ( SaveCorr=var )
```

Parameters:

SaveCorr Keyword for the correction value CORR

var Variable in which the correction value is to be stored.

- **f[xx]** - Element xx of the local float array.
- **gf[xx]** - Element xx of the global float array.
- **gfa[xx][yy]** - Element xx/yy of the global float field.

Example(s):

```
PPMGetCorr ( SaveCorr=gf [ 3 ] )
```

Store the correction value CORR in the variable gf[3].

8.8.2 PPMMID() - PPM 422 MID measurement

Description:

Perform a Multiple Ion Detection (MID) measurement (the data are collected from the M-Chassis).

Options:

- Display the measurement data
- Store the measurement data in the field **gfa[[]]** containing the global float values
- Store the measurement data in a file containing MID measurement data
- Store the measurement data in a file containing reference measurement data

Syntax:

```
PPMMID ( Par=file[ ;ind][ ,Disp=d][ ,Prot=p][ ,SaveGfa=var]_
        [ ,SaveRef=file[ ;ind][ ,SaveCyc=file[ ;ind][ ;m] ] )
```

Parameters:

Par	Keyword for parameter file
<i>file</i>	(Drive:\Path\) Filename <ul style="list-style-type: none"> • “<i>ffffff</i>” - Filename (e.g. “c:\qs\myfile”) • gs[<i>xx</i>] - Element <i>xx</i> of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
Disp	Keyword for display mode.
<i>d</i>	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • tab - Display as table. • bar - Display as bargraph. • vt - Display as Versus Time. <p>If nothing is entered, the mode for default is set to Display as table (tab). The default mode can be changed with the SetDefault() instruction.</p>
Prot	Keyword for log mode.
<i>p</i>	Log mode. <ul style="list-style-type: none"> • [default] - Default log mode. • off - Log output off. • on - Log output to printer. <p>Also see Chapter 8.4.3 Log output, 422.</p> <p>If nothing else is entered, the mode for default is off. The default mode can be changed with the SetDefault() instruction.</p>
SaveGfa	Keyword for storing measurement data in an array of the global data field gfa [][]. <ul style="list-style-type: none"> • The measurement data are stored in the float array <i>xx</i> (=var) of the field gfa[<i>xx</i>][]. Specification of the index <i>yy</i> for gfa[][] is not necessary, the measured values of the activated channels are stored consecutively in gfa[][0], gfa[][1], gfa[][2], etc., irrespective of their channel number. • All data stored in the elements <i>yy</i> = 0 ... 63 of the float array <i>xx</i>=var of gfa[<i>xx</i>][<i>yy</i>] are cleared before the data are saved.
<i>var</i>	Float array number ($0 \leq var \leq$ configurable maximum). <ul style="list-style-type: none"> • <i>cc</i> - numeric value (integer). • i[<i>xx</i>] - Element <i>xx</i> of the local integer array. • gi[<i>xx</i>] - Element <i>xx</i> of the global integer array.
SaveRef	Keyword for storing into a file containing reference data.
<i>file</i>	(Drive:\Path\) Filename. <ul style="list-style-type: none"> • “<i>ffffff</i>” - Filename (e.g. “c:\qs\myfile”). • gs[<i>xx</i>] - Element <i>xx</i> of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
SaveCyc	Keyword for storing the cycle data in a file.

<i>file</i>	(Drive:\Path\) Filename. <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile"). • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
<i>m</i>	Save mode <ul style="list-style-type: none"> • [append] - Add the new data to the end of the file (default). • reset - Clear the file before storing the data.

Example(s):

```
PPMID(Par="mypar")
```

Perform an MID measurement with the parameter file "mypar" and display the measured values as a table (default).

```
PPMMID(Par="c:\qs\par\sample", Disp=bar, _
        SaveCyc="c:\qs\dat\midcyc")
```

Perform an MID measurement with the parameter file "c:\qs\par\sample", display the data in Bar mode and store them in the data file "c:\qs\dat\midcyc".

8.8.3 PPMNeutrals() - Neutral particle analysis

Description:

PPMNeutrals() sets the mode for neutral particle analysis.

Syntax:

```
PPMNeutrals(Mode=n)
```

Parameters:

Mode	Keyword for neutral particle mode.
<i>n</i>	Neutral particle mode <ul style="list-style-type: none"> • coarse - coarse • fine - fine

Example(s):

```
PPMNeutrals(Mode=coarse)
```

8.8.4 PPMScanEnergy() - PPM 422 Energy Scan

Description:

This instruction performs an energy scan.

A scan is programmed on the E-chassis. The data are diverted to the E-chassis, so that they can be read in synchronously with the E-scan. Start energy and energy width are declared in this instruction, the other measurement parameters are taken from the Scan Analog parameter file.

Syntax:

```
PPMScanEnergy(Par=file[;ind],EStart=esta,EWidth=ewid,[,Disp=d])_
```

```
[ ,SaveRef=file[ ;ind]][ ,SaveCyc=file[ ;ind][ ;m]]
```

Parameters:

Par	Keyword for Scan Analog parameter file (*.SAP)
file	(Drive:\Path\) Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character array if the filename is stored in gs[].
ind	File index - See Chapter 8.4.1.3 File names with indices, 416
EStart	Keyword for energy start value.
esta	Start value in the range of -12 ... 512 eV <ul style="list-style-type: none"> • cc - Direct input of the numeric value (float). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.
EWidth	Keyword for energy width
ewid	EWidth within the range of -12 ... 512 eV <ul style="list-style-type: none"> • cc - Direct input of the numeric value (float). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array.
Disp	Keyword for display mode.
d	Display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>If nothing else is entered, the mode for default is set to on. The basic default setting can be changed by the SetDefault() instruction.</p>
SaveRef	Keyword for storing the measurement data as a reference file.
File	(Drive:\Path\) Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array if the filename is stored in gs[].
ind	File index - See Chapter 8.4.1.3 File names with indices, 416
SaveCyc	Keyword for storing the measurement data as a data file.
file	(Drive:\Path\) Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
ind	File index - See Chapter 8.4.1.3 File names with indices, 416
m	Save mode <ul style="list-style-type: none"> • [append] - Add the data to the end of the file (default). • reset - Clear the file before storing the new data.

NOTE:

The other measurement parameters are taken from the Scan Analog parameter file (*.SAP), e.g. mass = First Mass, measurement speed = Speed etc.

EStart is calculated as:

$$EStart = Vestart - SHIFT - CORR,$$

where the energy scan control voltage **Vestart** must be within the range 0 ... 512. If **Estart** < 0, **SHIFT** and **CORR** have to be set accurately. The sign of the effectively measured ion energy depends on the signs of the ions and **Ve**. Also see Chapter 5.9.1.6 Relations (Ion energy and **Ve**, **CORR** and **SHIFT**), 363.

	Ve+	Ve-
Ions Pos	-12 < E < 512	0 < E < 512
Ions Neg	-512 < E < 0	-12 < E < 512

Example(s):

```
PPMScanEnergy(Par="c:\qs\par\sample",EStart=5.0,EWidth=75.0,
               Disp=off,SaveCyc="c:\qs\dat\ppmcyc")
```

Load the measurement parameters from the file sample.sap and perform an energy scan from 5.0 eV ... 80.0 eV. The display is switched off and the data are stored in the file ppmcyc.sac:

8.8.5**PPMScanMassAnalog() - PPM 422 Scan Mass Analog****Description:**

Perform a mass scan and save the corresponding energy voltage.

Syntax:

```
PPMScanMassAnalog(Par=file[;ind][,Energy=ef][,Disp=d]_
                  [,SaveRef=file[;ind]]_
                  [,SaveCyc=file[;ind][;m]))
```

Parameters:

Par	Keyword for parameter file.
file	(Drive:\Path\)\ Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
ind	File index - See Chapter 8.4.1.3 File names with indices, 416
Energy	Keyword for energy value.
ef	Energy value within the range of -12 ... 512 eV <ul style="list-style-type: none"> • cc - Direct input of the numeric value (float). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.

Disp	Keyword for display mode.
<i>d</i>	display mode. <ul style="list-style-type: none"> • [default] - Default display mode. • off - Display off. • on - Display on. <p>If nothing else is entered, the mode for default is set to on. The basic default setting can be changed by the SetDefault() instruction.</p>
SaveRef	Keyword for storing the measurement data as a reference file.
<i>File</i>	(<i>Drive:\Path\</i>) Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
SaveCyc	Keyword for storing the measurement data as a data file.
<i>file</i>	(<i>Drive:\Path\</i>) Filename <ul style="list-style-type: none"> • "ffffff" - Filename (e.g. "c:\qs\myfile") • gs[xx] - Element xx of the global character string array, if the filename is stored in gs[].
<i>ind</i>	File index - See Chapter 8.4.1.3 File names with indices, 416
<i>m</i>	Save mode <ul style="list-style-type: none"> • [append] - Add the data to the end of the file (default). • reset - Clear the file before storing the new data.

Example(s):

```
PPMScanMassAnalog(Par="anascan1",Energy=200)
```

Perform a PPM 422 Scan Analog measurement with the parameter file "anascan1" and an energy of 200eV.

8.8.6 PPMSetEnergy() - Entering the ion energy

Description:

PPMSetEnergy() sets the ion energy of the E-chassis, considering the CORR and SHIFT values.

Syntax:

```
PPMSetEnergy( Energy=ef )
```

Parameters:

Energy	Keyword for energy value.
<i>ef</i>	Energy value within the range of -512 ... 512 eV <ul style="list-style-type: none"> • <i>cc</i> - Direct input of the numeric value (float). • i[xx] - Element xx of the local integer array. • gi[xx] - Element xx of the global integer array. • f[xx] - Element xx of the local float array. • gf[xx] - Element xx of the global float array. • gfa[xx][yy] - Element xx/yy of the global float field.

Example(s):

```
PPMSetEnergy(Energy=82.5)
```

Set the ion energy to 82.5 V.

```
SetVar(gf[0] = 5.0)
WhileVar(gf[0]<= 150.0)
Begin
  PPMSetEnergy(Energy=gf[0])
  Delay(5)
  ...
  //Get Data (from M-Chassis)
  MID(...)
  ...
  ...
  Calculate(gf[0]=gf[0]+0.5)
End
```

Perform a MID measurement with an energy of 5.0 ... 150.0 V and a step width of 0.5 V.

8.8.7 PPMSetEPolarity() - Polarity of the energy voltage

Description:

PPMSetEPolarity() sets the polarity of the energy voltage.

Syntax:

```
PPMSetEPolarity(Polarity=p)
```

Parameters:

Polarity Keyword for polarity of the Scan Energy voltage.

p Polarity (0+/-)

- **pos** - Ve+
- **neg** - Ve-

Example(s):

```
PPMSetEPolarity(Polarity=neg)
```

8.8.8 PPMSetIS() - Entering the ion source parameters

Description:

This instruction enables you to:

- Select the ion source sets for the E-chassis (IS-Set-E) and the M-chassis (IS-Set-M) separately
- Modify ion source voltages.
- Parameter Shift: modify all shift voltages simultaneously

Syntax:**For Plasma
Sampling Mode
Ions +/-Ions -**

```
PPMSetIS([IsSetM=n][,IsSetE=e][,Itro=mvari5][,Foc1=mvari8]_
[,Entry=mvarf9][,Centr=evarf2][,Foc2=evarf3]_
[,Faxis=evarf4][,Case=evari5][,Dfli=evari6]_
[,Dflo=evari7][,Mirr=evari8][,Infl=evarf9]_
[,Shift=volt])
```

**For Plasma
Sampling Mode
Neutrals Coarse**

```
PPMSetIS([IsSetM=n][,IsSetE=e][,CathC=mvarf2][,IonRef=mvari1]_
[,ExtrC=mvari5][,Foc1=mvari8][,Centr=evarf2]_
[,Foc2=evarf3],Faxis=evarf4][,Case=evari5]_
[,Dfli=evari6][,Dflo=evari7][,Mirr=evari8]_
[,Infl=evarf9][,Shift=volt])
```

**For Plasma
Sampling Mode
Neutrals Fine**

```
PPMSetIS([IsSetM=n][,IsSetE=e][,CathC=mvarf2][,IonRef=mvari1]_
[,ExtrF=mvari5][,Foc1=mvari8][,Centr=evarf2]_
[,Foc2=evarf3],Faxis=evarf4][,Case=evari5]_
[,Dfli=evari6][,Dflo=evari7][,Mirr=evari8]_
[,Infl=evarf9][,Shift=volt])
```

NOTE:

If you modify voltages of the ion source sets, the modifications remain local, i.e. they are modified on the mass spectrometers, but they are not stored to the corresponding IS-Sets.

Please note that some of the parameters are mutually exclusive:

- ExtrC - ExtrF is inhibited.
- ExtrF - ExtrC is inhibited.
- CathC, ExtrC or ExtrF - Itro and Entry are inhibited.

Parameters:

IsSetM	Load the ion source set for the M-chassis.
<i>n</i>	Number of the selected IS-Set 0 ... 15 of the group IS-M <ul style="list-style-type: none"> • <i>cc</i> - Direct input of the numeric value (float). • <i>i[xx]</i> - Element xx of the local integer array. • <i>gi[xx]</i> - Element xx of the global integer array.
IsSetE	Load the ion source set for the E-chassis.
<i>e</i>	Number of the desired IS set 0 ... 15 of the group IS-E <ul style="list-style-type: none"> • <i>cc</i> - Direct input of the numeric value (float). • <i>i[xx]</i> - Element xx of the local integer array. • <i>gi[xx]</i> - Element xx of the global integer array. <p>Keywords for the IS voltages of the sets IS-E and IS-M with corresponding validity range and step width:</p>
IonRef	<i>mvari1</i> - IS-M V1 : 0 ... 150 V, Step 1 V
CathC	<i>mvari2</i> - IS-M V2 : 0 ... 125 V, Step 0.5 V
Itro	<i>mvari5</i> - IS-M V5 : 0 ... 450 V, Step 2 V
ExtrC	<i>mvari5</i> - IS-M V5 : 0 ... 450 V, Step 2 V
ExtrF	<i>mvari5</i> - IS-M V5 : 0 ... 4.500 V, Step 0.02 V

Foc1	<i>mvari8</i> - IS-M V8 : -125 ... +125 V, Step 1 V
Entry	<i>mvari9</i> - IS-M V9 : 0 ... 60.00 V, Step 0.25 V
Centr	<i>evarf2</i> - IS-E V2 : 0.0 ... 125.0 V, Step 0.5 V
Foc2	<i>evarf3</i> - IS-E V3 : -30.00 ... +30.00 V, Step 0.25 V
Faxis	<i>evarf4</i> - IS-E V4 : 0.00 ... 60.00 V, Step 0.25 V
Case	<i>evari5</i> - IS-E V5 : 0 ... 450 V, Step 2 V
Dfli	<i>evari6</i> - IS-E V6 : 0 ... 450 V, Step 2 V
Dflo	<i>evari7</i> - IS-E V7 : 0 ... 250 V, Step 1 V
Mirr	<i>evari8</i> - IS-E V8 : -125 ... +125 V, Step 1 V
Infl	<i>evarf9</i> - IS-E V9 : 0 ... 60 V, Step 0.25 V

mvarfx, evarfxx Float variables for M-chassis or E-chassis

- *cc* - Direct input of the numeric value (float).
- **i[xx]** - Element *xx* of the local integer array.
- **gi[xx]** - Element *xx* of the global integer array.
- **f[xx]** - Element *xx* of the local float array.
- **gf[xx]** - Element *xx* of the global float array.
- **gfa[xx][yy]** - Element *xx/yy* of the global float field.

mvarix, evarix Integer variables for M-chassis or E-chassis.

- *cc* - Direct input of the numeric value (integer).
- **i[xx]** - Element *xx* of the local integer array.
- **gi[xx]** - Element *xx* of the global integer array.

Shift Keyword for simultaneously modifying the shift-voltages within the range of 0 ... 10 V (step width : 2 V).

The following voltages (called 'shift-voltages') can be modified by Shift:

IS-Set-E : V2-E...V8-E (Centr, Foc2, Faxis, Case, Dfli, Dflo, Mirr)

and for Plasma Sampling Mode 'Ion+/-' also the IS-Set-M voltages

V5-M, V8-M, and V9-M (ExtrC/ExtrF, Foc1, Entry)

volt Voltage to be added to the shift-voltages.

- *cc* - Numeric value (integer).
- **i[xx]** - Element *xx* of the local integer array.
- **gi[xx]** - Element *xx* of the global integer array.

Example(s):

```
PPMSetIS(IsSetE=0, IsSetM=2)
```

Load the IS-Set-E=0 and IS-Set-M=2

```
SetVar(gf[0] = 10.0)
PPMSetIS(IsSetM=12, IsSetE=12, Mirr=40, SHIFT=4)
WhileVar(gf[0] <= 120.0)
Begin
  ...
  PPMSetIS(Centr=gf[0])
  Calculate(gf[0]=gf[0]+0.5)
  ...
  ...
End
```

Load the IS-Set-E=12, IS-Set-M=12, add 4 V to the shift-voltages, set the IS voltage V8-E to 40 volts and increase the voltage V2-E step by step from 10.00 V ... 120.00 V by adding 0.5 V each time.

8.8.9 PPMSetPSM() - Entering the PPM 422 Plasma Sampling Mode (PSM)

Description:

Choose PPMSetPSM() to set the PPM 422 Plasma Sampling Mode.

Syntax:

```
PPMSetPSM(Mode=m)
```

Parameters:

Mode	Keyword for Plasma Sampling Mode (PSM)
<i>m</i>	Mode <ul style="list-style-type: none"> • IonPos - Ions+ • IonNeg - Ions- • Neutral - Neutrals

NOTE:

Please note that before a PPMSetPSM() instruction always the PPMSetIS() instruction must be called in order to load the corresponding ion source sets (E-Set and M-Set).

Example(s):

```
PPMSetPSM(Mode=IonPos)
PPMSetPSM(Mode=Neutral)
```

9. Access Rights Control

In the administration program Quadstar 32-bit Access Rights Control (subsequently referred to as ACCESSC) you can define access rights for different users. You can for example lock certain functions for certain users by not enabling the corresponding menu entries for them. You may as well set up a separate work environment for each user (separate directories for measure data, measure parameters, equipment in use etc.). The maximum number of users is limited to 20.

Basically we must distinguish between the User name, User password and the ACCESSC password (administrator password).

- The User name is defined for each Quadstar 32-bit user in the ACCESSC control program and enables him to use the functions of the Quadstar 32-bit programs he is authorized for.
- The User password can be defined by the Quadstar 32-bit user. It must be entered when a Quadstar 32-bit program is called for the first time. As long as the user has not defined a personal password, the default password <ENTER> (ENTER key) is used.
- The ACCESSC password is defined in the ACCESSC control program. It prevents the utilization of ACCESSC by unauthorized persons. This password protection is switched off after a new software installation so that ACCESSC can initially be called up without a password.

9.1 Login

9.1.1 Login to a Quadstar 32-bit program

If Quadstar 32-bit is password protected, you must enter your user name and user password to log into a Quadstar 32-bit program:




Fig. 9-1

User name and user password must be entered only when the first Quadstar 32-bit program is started, subsequent Quadstar 32-bit programs of the same installation recognize the user automatically. After the successful login, your work environment is setup. If you don't want to use the default password <ENTER> (ENTER key), you can define your own personal password by Change Password:



Fig. 9-2

If the password check is switched off (control field [Quadstar 32-bit Login enabled/Access Rights used] in ACCESSC not checked), but the Multiuser option is switched on (control field [Quadstar 32-bit multiuser enabled] in ACCESSC checked), only the user name needs to be selected from the list field.




Fig. 9-3

In this way, you can use your personal work environment without any password check.

9.1.1.1 Forgotten user password

If a user has forgotten his User password (for the Quadstar 32-bit programs) it can be reset in ACCESSC by [User] > [Delete Password]. See Chapter 9.2.4.2.7 Delete user Password, (539).

The user can then log into a Quadstar 32-bit program with the user name and the default user password <ENTER> (ENTER key). He can subsequently define his personal user password with the Change Password command.

9.1.1.2 Logging into ACCESSC

When ACCESSC is password protected ([Protection] > [ACCESSC Login Setup] > [Password enabled]), the administrator password must be entered when the ACCESSC program is started:



Fig. 9-4

The ACCESSC password is defined under [Protection] > [ACCESSC Login Setup] .

9.1.1.3 Forgotten ACCESSC password

If the ACCESSC password is lost, the Q_ACCESS.DAT file must be copied from the installation diskette into the execution file path of Quadstar 32-bit. Subsequently the ACCESSC control program can be started without a password like after the initial software installation. The previously defined user passwords will not be affected.

9.1.2 Safeguarding the password protection

Access privileges can be changed at any time using the ACCESSC control program, if it is not password protected. Therefore it is important, that ACCESSC is protected by another password, and the installation and backup diskettes are stored in a safe place.

If you have set up an extensive password protection (for multiusers), you should make backup copies of the Q_ACCESS.DAT, QUADSTAR.INI and Q_PSW.DAT files.

9.2 Control of the access privileges




Fig. 9-5

9.2.1 Selecting the Quadstar 32-bit program

Choose the program, for which you want to determine access rights:




Fig. 9-6

9.2.2 Defining the access privileges



Fig. 9-7

For displaying the access privileges of other users, select their names from the drop-down list boxes. The user privileges are enabled by checking the corresponding control boxes in the column below the user name.

9.2.3 Switch the password protection on/off

By the control fields in the lower left corner, you can switch on/off the password protection, the automatic start of sequences and the different work environments:

- QUADSTAR 32-bit Login enabled / Access Rights used
- Run sequence by application start enabled
- QUADSTAR 32-bit Multiuser enabled

Fig. 9-8

Quadstar 32-bit Login enabled/Access Rights used

- Checked: Quadstar 32-bit password protection enabled.
=> When the first Quadstar 32-bit program is started, input of the user name and user password is requested. The granted access privileges and the specific work environment are then effective.
- Unchecked: Quadstar 32-bit password protection disabled.
=> No input of the user name or user password is requested.
=> The access privileges are not restricted, the work environment is the same for all users.

Run sequence by application start enabled

- Checked: A sequence can be executed automatically, when the Measure program is started.

Quadstar 32-bit Multiuser enabled

- Checked: If the password protection is switched off, a separate work environment for each user can although be supplied by using this option.

9.2.4 Main menu



Fig. 9-9

File

Display the version number, store the program settings and quit the program.

User

Modify the user password (change, delete etc.) and the user-specific work environment.

Protection

Password protection for the ACCESSC program.

Help

Displays the help contents. See Chapter 1.9 Hints concerning the work with Quadstar 32-bit, 56.

9.2.4.1 File menu



Fig. 9-10

About Access Rights Control...

Display the current version number of the ACCESSSC program.

Save

Save all settings.

Print...

Print access rights of one or all users.

Exit

Quit the Quadstar 32-bit ACCESSSC program.

9.2.4.1.1 About Access Rights Control



Fig. 9-11

Display the version number of the Quadstar 32-bit ACCESSSC program:

9.2.4.2 User menu




Fig. 9-12

Change Name...

Change an existing user name.

Insert...

Insert a new user name.

Delete...

Delete a user name.

Copy...

Copy user privileges.

Path Setup...

Define the paths for the user's working environment.

Set Predefined Access...

Set predefined access privileges.

Delete Password...

Delete a user password.

9.2.4.2.1 Change name

Modify an existing user name. Under Select Name choose the user name to be changed. Now enter under Edit Name a new name and confirm by [OK]; the old name is replaced, but the already defined user privileges remain unchanged.




Fig. 9-13

9.2.4.2.2 Insert

By Insert New User you can add more user names.




Fig. 9-14

NOTE:

The default password is <ENTER> (ENTER key). When a Quadstar 32-bit program is started, the user can define his own password by using Change Password.

9.2.4.2.3 Delete

The user name and the corresponding access privileges are deleted from the list.




Fig. 9-15

9.2.4.2.4 Copy

The access privileges of the user selected in the Source list box are copied to the user selected in the Destination list box by clicking the [OK] button:




Fig. 9-16

9.2.4.2.5 Define the work environment (Path Setup)




Fig. 9-17

Execution File Path

Directory of the Quadstar 32-bit program files (.EXE).

Auxiliary File Path

Directory for auxiliary files (the information about the equipment in use is stored here, for example).

Parameter File Path

Directory for the measure parameter files.

Data File Path

Directory for the measure data files.

[Select Path]

Choose the drive and directory.

9.2.4.2.6 Set Predefined Access




Fig. 9-18

Assign predefined user privileges:

- DEFAULT - Adopt the factory settings.
- ENABLE_ALL - Enable all access privileges
- DISABLE_ALL - Disable all access privileges

NOTE:

The selection applies to all Quadstar 32-bit programs.

9.2.4.2.7 Delete user Password



Fig. 9-19

Choose Delete Password to delete a user password. This is useful e.g., if a user has forgotten his password.

9.2.5 Protection menu



Fig. 9-20

9.2.5.1 Administrator password (ACCESSSC Login Setup...)

Password protection of the Quadstar 32-bit ACCESSSC-program:



Fig. 9-21

New Password

Enter a new ACCESSC password. The old password is deleted.

Login ACCESSC


- Password disabled - Deactivate the password protection
- Password enabled - Activate the password protection.
=> The password protection is not switched on/off unless Set is activated.

10. File Handling, Data Formats, Error Messages, Troubleshooting

10.1 File Manager

The file manager supports the handling of all parameter and measurement data files. It can be used for opening files, creating new files, deleting old files, and it also provides an easy method for file copying. Some files relate to others in such a way that they can only be used together as a database (e.g. MCD-measurement: File.MCC, File.MCO and File.MCP belong together). Therefore, Quadstar 32-bit files should always be manipulated jointly by using the Quadstar 32-bit file manager, since this way all necessary files will be covered.

Procedure

Quadstar 32-bit calls the file manager whenever a file must be accessed that's not set yet. For the following descriptions, the assumption is that the file manager dialog box is displayed on the screen as shown in Fig. 10-1,  542:

Loading a file

- 1 Use the mouse to select the desired path under «Look in:».
- 2 Position the cursor on the File name text box.
- 3 Type in the file name or select a file name in the File Name list box.
- 4 Press [OK].

Creating a new file

- 1 Use the mouse to select the desired path under «Look in:».
- 2 Position the cursor on the File name text box.
- 3 Type in the name of the file to be created.
- 4 Press [OK].
- 5 Choose Yes to answer the question: "File not existing, create ?"

Copying a file

- 1 Use the mouse to select the desired path under «Look in:».
- 2 Position the cursor on the File name text box.
- 3 Type in the source file name or use the mouse to select a file name in the File Name list box.
- 4 Press the Copy button.
- 5 In the «to» text box enter the name of the destination file.
- 6 Press [OK].

- 7 Choose the new file in the File list box and press Open.

Deleting a file

- 1 Use the mouse to select the desired path under «Look in:».
- 2 Position the cursor on the File name text box.
- 3 Type in the file name or use the mouse to select a file name in the Files list box.
- 4 Press [Delete].
- 5 Confirm the Quadstar 32-bit prompt with Yes or No.

Renaming a file

- 1 Use the mouse to select the desired path under «Look in:».
- 2 Position the cursor on the File name text box.
- 3 Type in the source file name or use the mouse to select a file name in the Files list box.
- 4 Choose Rename.
- 5 Position the cursor on the «to» text box and change the name.
- 6 Press [OK].

Dialog box of the file manager




Fig. 10-1

Control elements

[OPEN]

Open the selected file.

[Cancel]

Quit the file manager without activity.

[Copy]

Copy the contents of a file into an existing file or create a new file with the same contents.

[Rename]

Rename a file.

[Delete]

Delete a file.

NOTE:

Please note that a parameter set, e.g. for gas-specific calibration, can comprise several files. For this reason, Quadstar 32-bit files should always be manipulated with the Quadstar 32-bit file manager so that all related files are handled appropriately.

Filename

Enter the name of the file to be opened.

The file name can be entered directly or be chosen from a list of files. This list shows only files of the type that's chosen in the field «Files of type».

NOTE:

To display a list of files with a specific extension, enter the asterisk (*), a dot and the three character extension.

Example: If you want to view all files of a directory that have the extension .SAC (containing Scan Analog Cycles), enter *.sac.


Look in

Browse the directory in which the file to be opened is located.

Files of type

In this list Quadstar 32-bit displays only the file types that are relevant to the called program. If several file types exist, mark the type of file to be opened.

File Title of '...'

In this field the file title appears - if there is one - as soon as the file is selected (not opened yet). The file title can only be displayed here, but not edited. Enter the file title (within the measure data documentation) before the measurement starts. See Chapter 2.1.5 Documenting the measurement data (File Info),  69).

10.2 Path and file names

Path and file names may be chosen almost freely, with the following restrictions:

- They may contain regular letters (capital or small) and figures in any order.
- They may contain the special characters ! # \$ % & () - \ { } ~ as well.
- They are to contain no other/special characters (such as space, Å, ö, é, %, @ etc.).
- A file name can consist of maximum 255 characters including three characters extension (e.g. a_long_filename.abc).
- Path and file name together can consist of maximum 255 characters.

10.3 Quadstar 32-bit files

Quadstar 32-bit uses universal data formats for storing the measurement data.

10.3.1 Overview

The Quadstar 32-bit files are allocated to the following sub directories:

- Programs and setup utilities. (EXE Group)
- Auxiliary files that are unique to the application. (AUXI Group)
- Parameter files whose names can be selected via the file manager. (PAR Group)
- Data files whose names can be selected via the file manager. (DAT Group)

The following tables contain a choice of the used files and file types:

Programs

Index	Name	Group	Designation
n/a	Accessc.exe	EXE	
n/a	ComServ.exe	EXE	Communication Server for ARCNET server mode
n/a	DispSav.exe	EXE	Presentation of stored data
n/a	Measure.exe	EXE	Measurement program
n/a	Parset.exe	EXE	Parameter setup program
n/a	Service.exe	EXE	Service
n/a	TuneUp.exe	EXE	Tuning
n/a	Utility.exe	EXE	Auxiliary functions

Additional programs

Index	Name	Group	Designation
n/a	PrtExe.exe	EXE	Printer executer for handling delayed line printing
n/a	PrtServ.exe	EXE	Printer server for spooling line printing

Help files

Index	Name	Group	Designation
n/a	Quadstar 32-bit.chm	EXE	Quadstar 32-bit Online HTML-Help

Tools

Index	Name	Group	Designation
n/a	Exchange.exe	EXE	Exchange of pathnames in sequences during installation
n/a	Exchange.pif	EXE	Console settings for running EXCHANGE.EXE

Dynamic Link Libraries (DLLs)

Index	Name	Group	Designation
n/a	General.dll	EXE	General Quadstar 32-bit functions
n/a	SeqEdit.dll	EXE	Quadstar 32-bit sequence editor
n/a	Resource.dll	EXE	Quadstar 32-bit resources
n/a	UserCtrl.dll	EXE	Quadstar 32-bit specific controls
n/a	Ss32d25.dll	EXE	Tools

Quadstar 32-bit settings

Index	Name	Group	Designation
220	Q_Option.dat	EXE	Contains authorization information (encrypted)
n/a	Quadstar.ini	EXE	Quadstar 32-bit initialization file

Access settings

Index	Name	Group	Designation
203	Q_Access.dat	EXE	Password for accessing ACCESSSC program
200	Q_PSW.dat	EXE	Access rights, passwords and user names for ACCESSSC program

AUXI (device dependent)

Index	Name	Group	Designation
131	Q_AIChar.par	AUXI	Analog Input characteristic curve settings
132	Q_AOChar.par	AUXI	Analog Output characteristic curve settings
320	Q_calfac.lib	AUXI	Calibration factors for MCD
101	Q_config.par	AUXI	QMS hardware settings
130	Q_IOD.par	AUXI	Digital I/O settings
120	Q_is.par	AUXI	Ion source parameters (QMS422/421)
123	Q_is200.par	AUXI	Ion source parameters (QMS200)
121	Q_ischp.par	AUXI	Ion source channel parameters
201	Q_MasCal.dat	AUXI	Adjusted mass numbers
107	Q_Setup1.par	AUXI	Setup parameters (multiplex)
202	Q_Zero.dat	AUXI	Zero gas measurement data

AUXI (device independent)

Index	Name	Group	Designation
104	Q_ColTab.par	AUXI	Screen color settings
109	Q_CPECfg.par	AUXI	Channel parameter editor settings
210	Q_EH.dat	AUXI	Error history (generated at runtime)
211	Q_EH_AI.dat	AUXI	Error history - supplementary information (generated at runtime)
n/a	Q_Fragm.wri	AUXI	Key fragment ions table
122	Q_leak.par	AUXI	Leak test parameters
n/a	Q_peasrc.par	AUXI	Peak search parameters
105	Q_Print.par	AUXI	Printer settings
102	Q_QMSCom.par	AUXI	QMS communication settings
106	Q_SeqOpt.par	AUXI	Sequence editor settings
100	Q_Setup.par	AUXI	Setup parameters
300	Q_Spec.LIB	AUXI	Spectra library: data file
301	Q_Spec_h.LIB	AUXI	Spectra library: header file
133	Q_NodeID.PAR	AUXI	Node settings for ARCNET mode
108	Q_FIn.par	AUXI	File info settings
110	Q_pso.par	AUXI	Peak shape optimizer parameters

AUXI (only PPM)

Index	Name	Group	Designation
10002	P_setup.par	AUXI	PPM Setup parameter
10000	P_is.par	AUXI	PPM Ion source parameters

DAT (general)

Index	Name	Group	Designation
700	*.mdc	DAT	MID/MCD cycle data files
600	*.mdr	DAT	MID/MCD reference data files
710	*.sac	DAT	Scan Analog data files
610	*.sar	DAT	Scan Analog reference files
720	*.sbc	DAT	Scan Bargraph cycle data files
721	*.sbh	DAT	Scan Bargraph cycle data header files
n/a	*.sbr	DAT	Scan Bargraph reference files

DAT (grafic editor settings)

Index	Name	Group	Designation
602	*.amr	DAT	MID/MCD reference files
612	*.aar	DAT	Scan Analog reference files
622	*.abr	DAT	Scan Bargraph reference files
702	*.amd	DAT	MID/MCD cycle data files
712	*.aad	DAT	Scan Analog data files
723	*.abd	DAT	Scan Bargraph cycle data files

DAT (additional file information)

Index	Name	Group	Designation
601	*.imr	DAT	MID/MCD reference files
611	*.iar	DAT	Scan Analog reference files
621	*.ibr	DAT	Scan Bargraph reference files
701	*.imd	DAT	MID/MCD cycle data files
711	*.iad	DAT	Scan Analog data files
722	*.ibd	DAT	Scan Bargraph cycle data files

PAR

Index	Name	Group	Designation
531	*.gcc	PAR	Gas specific calibration: channel parameters
530	*.gcp	PAR	Gas specific calibration: matrix parameters
521	*.mcc	PAR	MCD channel parameters
522	*.mco	PAR	MCD concentration parameters
520	*.mcp	PAR	MCD matrix parameters
500	*.mip	PAR	MID parameters
503	*.msp	PAR	Mass scale calibration parameters
501	*.sap	PAR	Scan Analog parameters
502	*.sbp	PAR	Scan Bargraph parameters
n/a	*.seq	PAR	Sequencer text files (source code)
510	*.sqe	PAR	Executable sequencer files (compiled code)

* = any valid file name (for example AIRDEMO)

10.4 File formats

For detailed information on structure and use of Quadstar 32-bit files please contact your nearest service office.

10.5 Data export ASCII

In the Display Saved Values program you can convert binary data to the ASCII format. Choose this option if Quadstar 32-bit measurement data are to be processed in other programs such as Microsoft® EXCEL® etc.

Please note that the parameters:

- Column Separator (Separator for dividing columns)
- Line Separator (Separator for dividing lines)
- TimeFraction (Parameter for outputting the fraction of time)

can be chosen in the Parameter Setup program under [Setup] > [ASCII Conversion]. If the default values are chosen, the data will be exported in a format that can be read by EXCEL.

10.5.1 Scan Analog data

Scan Analog measurement data can be converted to the ASCII format in the Display Saved Values program under [Scan] > [Analog Data] > [File] > [Convert to ASCII]. Only 1 cycle per ASCII file can be exported.

NOTE:

such files can become very large, since there are up to 64 measurement values per mass (QMS 200: 32) to be stored.

The following example shows the first cycle of the file airdemo.sac, in the range of 16 ... 21 amu (here only continued until 16.5 amu):

```

ASCII SCAN                                airdemo.sac
ANALOG DATA:
DATE:          14.03.96      TIME:      14:46:36

Number of stored                        21
SCA cycles
Converted SCA cycle                    1

INFO SCAN:                                1
Data Titel      Ion Counter
Data Unit       cps
Scan Titel      Mass
Scan Unit       amu
Comment
First Mass      0
Scan Width      60
Values/Mass     64
Zoom Start      16
Zoom End        21

      Cycle      CycleDate      CycleTime
      1          14.03.96      14:46:37:42
ScanData      1
      16          0
      16.02      0
      16.03      0
      16.05      0
      16.06      0
      16.08      0
      16.09      0
      16.11      0.649044
      16.13      1.66857
      16.14      2.92892
      16.16      4.33886
      16.17      5.72936
      16.19      6.90862
      16.2        7.69983
      16.22      7.97851
      16.23      7.69983
      16.25      6.90862
      16.27      5.72936
      16.28      4.33886
      16.3        2.92892
      16.31      1.66857
      16.33      0.675122
      16.34      1.61E-15
      16.36      0
      16.38      0
      16.39      0
      16.41      0
      16.42      0
      16.44      0
      16.46      0
      16.48      0
      16.49      0

```

Fig. 10-2

10.5.2 Scan Bargraph data

Scan Bargraph measurement data can be converted to the ASCII format with [Scan] > [Bargraph Cycles] > [File] > [Convert to ASCII].

The following example shows peaks 18, 28 and 32 exported from 12 Scan Bargraph spectra.

```

ASCII SCAN BARGRAPH CYCLES :airdemo.sbc
DATE : 24.04.96      TIME : 14:25:07
CONVERTED CYCLES :      12

Number of stored SCB cycl      100
Printed start SCB cycle      1
Printed end SCB cycle      12
Min Total Ion Current      3.36588E-005[A]
Max Total Ion Current      3.36606E-005[A]

First mass      10
Scan width      35

Enabled m/e      Statistic Info
18      100 [%]
28      100 [%]
32      100 [%]

```

Cycle	Date	Time	RelTime[s]	TIC[A]	m/e	Intensity[A]
1	24.04.92	14:25:10:13	3.13	3.37E-05	18	1.25E-06
2	24.04.92	14:25:14:14	7.14	3.37E-05	18	1.25E-06
3	24.04.92	14:25:18:10	11.1	3.37E-05	18	1.25E-06
4	24.04.92	14:25:22:11	15.11	3.37E-05	18	1.25E-06
5	24.04.92	14:25:26:17	19.17	3.37E-05	18	1.25E-06
6	24.04.92	14:25:30:24	23.24	3.37E-05	18	1.25E-06
7	24.04.92	14:25:34:13	27.13	3.37E-05	18	1.25E-06
8	24.04.92	14:25:38:14	31.14	3.37E-05	18	1.25E-06
9	24.04.92	14:25:42:21	35.21	3.37E-05	18	1.25E-06
10	24.04.92	14:25:46:22	39.22	3.37E-05	18	1.25E-06
11	24.04.92	14:25:49:79	42.79	3.37E-05	18	1.25E-06
12	24.04.92	14:25:54:18	47.18	3.37E-05	18	1.25E-06

Fig. 10-3

10.5.3

MID/MCD Data

MID/MCD measurement data can be converted to the ASCII format with [Process] > [Cycles] > [File] > [Convert to ASCII].

The following example shows the intensities of the masses 14, 16, 18, 28, 32, 40, 44 exported from the data block 0 (MID-Cycle 50 ... 72) of the airdemo.mdc file.

```

ASCII MID CYCLES :          airdemo.mdc
DATE :      24.04.96   TIME :   13:10:03
CONVERTED CYCLES :          21

Number of stored cycles      200
Printed start cycle          50
Printed end cycle            70
Number of stored datablocks  1

Datablock 0 Ion Current      [A]
'0/0' 14  min: 8.32E-06  max: 8.32E-06
'0/1' 16  min: 2.49E-06  max: 2.49E-06
'0/3' 18  min: 1.25E-06  max: 1.25E-06
'0/5' 28  min: 1.00E-05  max: 1.00E-05
'0/7' 32  min: 8.32E-06  max: 8.32E-06
'0/8' 40  min: 1.67E-06  max: 1.67E-06
'0/9' 44  min: 6.31E-07  max: 6.32E-07

Cycle Date      Time RelTime[s]      '0/0'      '0/1'      '0/3'      '0/5'      '0/7'      '0/8'      '0/9'
50 24.04.92 13:10:17:99 14.99 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
51 24.04.92 13:10:18:32 15.32 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
52 24.04.92 13:10:18:60 15.6 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
53 24.04.92 13:10:18:92 15.92 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
54 24.04.92 13:10:19:25 16.25 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
55 24.04.92 13:10:19:53 16.53 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
56 24.04.92 13:10:19:86 16.86 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
57 24.04.92 13:10:20:13 17.13 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
58 24.04.92 13:10:20:41 17.41 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
59 24.04.92 13:10:20:68 17.68 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
60 24.04.92 13:10:21:10 18.01 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
61 24.04.92 13:10:21:29 18.29 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
62 24.04.92 13:10:21:62 18.62 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
63 24.04.92 13:10:21:89 18.89 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
64 24.04.92 13:10:22:22 19.22 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
65 24.04.92 13:10:22:55 19.55 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
66 24.04.92 13:10:22:82 19.82 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.31E-07
67 24.04.92 13:10:23:15 20.15 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
68 24.04.92 13:10:23:43 20.43 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07
69 24.04.92 13:10:23:70 20.7 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.31E-07
70 24.04.92 13:10:23:98 20.98 8.32E-06 2.49E-06 1.25E-06 1.00E-05 8.32E-06 1.67E-06 6.32E-07

```

Fig. 10-4

10.6 Error messages/warnings

10.6.1 General

Some of the Quadstar 32-bit messages that are displayed are based on the error messages and warnings transmitted by QMS mass spectrometer; they contain only a brief comment and are not explained here. For detailed information please refer to the QMS 422/421 or QMS 200 operating instructions:

QMS Error 'XX'

See operating instructions of the corresponding mass spectrometer under "Troubleshooting".

QMS Warning 'XX'

See operating instructions of the corresponding mass spectrometer under "Warnings".

10.6.2 Quadstar 32-bit Error Messages/Warnings

Error/Warning	Description
0.4 <= mass deviation < 0.5	Mass scale calibration was not possible, the deviation exceeds the tolerance limit. => Selected mass number is incorrect or inaccurate mass calibration on the RF generator.
ADC error	QMS error 29 (service message).
additional data doesn't match existing file	The additional data to be stored with the scan analog data doesn't match the already stored data.
adjust mass too high	Peak Adjust error; the selected mass is too far to the right of the peak.
adjust mass too low	Peak Adjust error; the selected mass is too far to the left of the peak.
adjust range and mass too low	Peak Adjust error; the electrometer Range is too low and the mass number too far to the left of the peak.
adjust range or threshold too high	Peak Adjust error; the electrometer Range or the threshold is too high.
adjust range too low	Peak Adjust error; the electrometer Range is too low.
adjust range too low and mass too high	Peak Adjust error; the electrometer Range is too low and the mass number too far to the right of the peak.
alarm at cycle: 'Cycle', ident: 'Ident', value: 'Value'	Alarm (Level A , B) for cycle: 'Cycle', Identification: 'Ident', measured value: 'Value'.
AO-DAC error	QMS error 28 (service message).

Error/Warning	Description
ArcNet driver initialisation error	The ArcNet driver for Windows could not be initialized correctly, so Quadstar 32-bit cannot communicate via ArcNet. Restart Windows, update your Quadstar 32-bit installation or contact our service office.
ArcNet parameters of another installation used	Warning: If several Quadstar 32-bit-installations use the ArcNet server, the ArcNet-protocol-parameters of the first installation are used.
ArcNet server not started or client not connected	Error on starting or announcing to ArcNet server. Check the amount of free resources or restart Windows®.
argument of lgn() <= 0 at sequence 'SeqName' line 'LineNbr'	The argument of the lgn() function (logarithm to base e) in 'LineNbr' of 'SeqName' is 0.
argument of log()<=0 at sequence 'SeqName' line 'LineNbr'	The argument of the log() function (logarithm to base 10) in 'LineNbr' of 'SeqName' is 0.
argument of srt() <0 at sequence 'SeqName' line 'LineNbr'	The argument of the srt() function (square root calculation) in 'LineNbr' of the 'SeqName' is 0.
argument out of range at sequence 'SeqName' line 'LineNbr'	The argument in 'LineNbr' of 'SeqName' is outside the limits.
ASCII line buffer overflow	Overflow in buffer for ASCII-Lineprinting. Check whether your sequence prints correctly.
at least one user has to be defined to save data	At least one user must be defined when the ACCESSC program is terminated.
auto not possible	Automatic operation is inhibited (filament 1+2) in degas mode (QMS Warning 20).
autorange algorithm couldn't detect measure value	The (amplifier) measure range can not be determined automatically, for example because the measured values vary too much.
bad file handle	Windows® system error, access to an opened file is denied.
baseline too high at mass: 'Mass-Nbr', offset: 'Offset'%, range:'1E-XX'	Quadstar 32-bit warning. The electrometer offset 'Offset' determined for the electrometer range '1E-XX' (measured on Zero-Mass 'MassNbr') amounts to 'Offset'%, i.e. it is >1.99%. (=> The selected Zero-Mass is unfavorable or the hardware must be realigned)

Error/Warning	Description
baseline too low at mass: 'Mass-Nbr', offset: 'Offset%', range:'1E-XX'	The electrometer offset 'Offset' determined for the electrometer range '1E-XX' (measured on Zero-Mass 'MassNbr') amounts to 'Offset%', i.e. it is <1.99%). (=> The selected Zero-Mass is unfavorable or the hardware must be realigned)
baudrate of QMS not correct	The QMS uses a bad transmission rate.
buffer not o.k.	QMS warning 22. Error in the QMS measurement data buffer management.
calibration factor of internal standard not defined	Calibration factor of the internal standard is not defined. (=> Define the internal standard in the selected calibration matrix)
calibration matrix undefined	Parameter setup error in the selected calibration matrix. => Verify that all selected mass numbers are assigned to a component.
CAN error	QMS error message 40. Malfunction of the ClosedAreaNetwork controller.
can't load printer install	Dialog box for the printer installation cannot be loaded.
Can't open communication Line for online printing	The defined printer port for outputting the sequence log cannot be used. => Make sure the printer port defined in "Sequence Printing Setup" is actually installed. Correct the port specification if necessary.
can't save bitmap(DIB) to .BMP file	The picture to be stored can not be written to file, e.g. because there is not enough disk space.
CD error	QMS error message 35. SEM high voltage on HV 421 overloaded or faulty.
channel doesn't exist in parameter file at sequence 'SeqName' line 'LineNbr'	The specified channel in the given parameter file doesn't exist or is switched OFF.
channel not enabled	Selected QMS channel is not in the ENABLE state.
CLA 421 AMPE overrange	Chopper Lock-In amplifier (electrometer) overrange. => Adapt the channel parameter Amp.Range.

Error/Warning	Description
CLA 421 AMP1 overrange	Chopper Lock-In amplifier 1 overrange. => Adapt the channel parameter Amplifier 1.
CLA 421 AMP2 overrange	Chopper Lock-In amplifier 2 overrange. => Adapt the channel parameter Amplifier 2.
CLA 421 AMP3 overrange	Chopper Lock-In amplifier 3 overrange. => Adapt the channel parameter Amplifier 3.
client identification mismatch in Arc-Net server	Client could not be identified by ArcNet server. System error message. Restart Windows® or call our service office.
client not disconnected	Client could not be cancelled by ArcNet server. System error message. Restart Windows® or call our service office.
cold cathode off	Sequencer test warning. => The cold cathode gauge should be switched on by SetQMS(Pen=on) before the first measurement is started.
communication error from QC to CS	QMS error message 42. No response from the console.
communication error on LAN	QMS error message 44. Faulty communication via LocalAreaNetwork interface.
communication error on RS232C	QMS error message 45. Faulty communication via RS232C interface.
communication error QMS control DSP	QMS error message 43. Faulty communication via dual port RAM.
component 'CompName' already defined	The selected component name 'CompName' is already defined in the current matrix.
concentration not defined	The concentration of a component is not defined.
corrupted identification file (Q_OPTION.DAT)	The Quadstar 32-bit identification file cannot be read. => Copy the original file from the installation diskette.
corrupted password file (Q_PSW.DAT)	The Q_PSW.DAT file defines passwords and access rights. The file is automatically cleared and must subsequently be recreated by the ACCESSC application.

Error/Warning	Description
cross device link	Input error, e.g. different disk drive or directories in rename function.
data assignment error in ArcNet server communication	The synchronous communication to ArcNet server is disturbed. Restart Windows® or call our service office.
DDE Error at sequence 'SeqName' line 'LineNbr'	The data exchange via DDE was not possible on line 'LineNbr' of the sequence 'SeqName'.
degas not possible	QMS warning 19. The degas function is inhibited in 'Auto'-mode (Filament 1+2).
device already open	The selected interface cannot be opened because it claims to be open already.
device can't print bit maps	The selected printer does not support graphic output.
device not open	Access to the selected interface is denied because it has not been opened.
different number of users found in (Quadstar.INI) and (Q_PSW.DAT)	The two files Quadstar.INI and Q_PSW.DAT do not have the same number of users. => This error message is displayed if a different number of users is found in these two files when ACCESSC is started => With the subsequent prompt "Load the default access rights for all users ?" the user is able to automatically adapt the Q_PSW.DAT file to the Quadstar.INI file. In this case, however, all access rights of all users are set to DEFAULT. If this automatic adaptation is not wanted, the ACCESSC program is cancelled.
different revision number in identification file 'FileName'	Update Error, the Quadstar 32-bit program (EXE-File) does not have the same revision number as the identification file (Q_OPTION.DAT).
different software versions (ArcNet server and application)	Quadstar 32-bit-application and ArcNet server cannot work together because their software versions are different. Restart Windows® or check your installation(s). Different Quadstar 32-bit-versions can not be run simultaneously on the same PC.

Error/Warning	Description
different usernames found in (Quadstar.INI) and (Q_PSW.DAT)	The two files Quadstar.INI and Q_PSW.DAT must have the same user names. => If the same number of users but different user names are found when ACCESSC is started, the user names in the Q_PSW.DAT file are updated according to the specifications in the Quadstar.INI file. Changes are displayed. => User names can only be specified in the ACCESSC application.
directory 'DirName' not empty, not possible to remove	The directory 'DirName' is not empty so it can't be removed.
division by zero at sequence 'SeqName' line 'LineNbr'	Division by zero occurred in line 'LineNbr' of the sequence 'SeqName'.
drive change error	The selected disk drive is invalid.
DSP dual port RAM defect	QMS error message 25. Access to signal processor defective.
DSP EPROM defect	QMS error message 24. The checksum of the EPROM signal processor is invalid.
emission error	QMS error message 39. The emission current does not correspond to the specified nominal value.
emission error on QME #'QMENbr'	The emission of the QME 125, number 'QMENbr' has failed (Multiplex operation via the QMU 112).
emission not possible (external)	The emission cannot be switched on or off by Quadstar 32-bit (Sequencer SetQMS()) because in Parameter Setup the option "Emission External" is selected. See Chapter 3.8.1 QMS Configuration, 234.
emission off	The measurement cannot be started; turn on the emission.
empty frame received	Communication error; an empty frame has been received from the QMS .
empty frame sent	Error in ArcNet driver; the frame to be transmitted is empty.
energy parameters don't match existing file 'FileName'	The scan comparison cannot be performed; the energy range of the reference scan does not match the energy range of the current scan (only PPM).
error in deleting printer reservation	System error. Restart Windows® or contact our service office.

Error/Warning	Description
error in printer reservation	System error. Restart Windows® or contact our service office.
error sending ASCII line to printer server	System error. Restart Windows® or contact our service office.
external normalization factor is zero	The external normalization factor is zero. => Perform a gas-specific calibration. Make sure that in the calibration matrix the determination of the external normalization factor is activated.
fatal communication error 'Number'	Error message 'Number' of the window communications driver.
filament 1 defect	QMS error message 37. Filament 1 defective.
filament 2 defect	QMS error message 38. Filament 2 defective.
file mismatch found, file(s) might be corrupted	Files of a compound database don't match each other or are damaged. Save it all first and try it again.
'FileName' isn't of the expected file extension	Due to an incorrect file name extension the selected 'FileName' is not accepted.
'FileName' isn't of the expected file type	The type of the selected 'FileName' does not correspond to the required file type.
file 'FileName' exceeds maximal size	The file 'FileName' is larger than allowed.
file 'FileName' exceeds maximal size, file truncated	The file 'FileName' is larger than allowed and is therefore shortened.
file 'FileName' not existing at sequence 'SeqName' line 'LineNbr'	The file 'FileName' on the specified line 'LineNbr' of the sequence 'SeqName' does not exist.
file already existing 'FileName'	The file 'FileName' cannot be created because it exists already.
file check impossible at 'SeqName' line 'LineNbr', variable file name	Warning in the test for file existence. => Test not possible; the file specified in line 'LineNbr' of the sequence 'SeqName' is a file with a variable name.
file close error 'FileName'	The file 'FileName' cannot be closed.
file manager error	Error when working with the file manager.
file name too long at sequence 'SeqName' line 'LineNbr'	The file name specified in the line 'LineNbr' of the sequence 'SeqName' contains more than 255 characters.

Error/Warning	Description
file not existing 'FileName'	Access to the file 'FileName' is denied because it does not exist.
file permission denied 'FileName'	Access to the file 'FileName' is not allowed.
file positioning error 'FileName'	Positioning error when writing or reading the file 'FileName'.
file (Quadstar.INI/entry lastuser) not found	No user name can be loaded. The 'lastuser' entry in the Quadstar.INI file is missing => The application is cancelled. => To remedy the fault load the Quadstar.INI file from the installation diskette. All previously defined users and access rights will be lost.
file too small 'FileName'	The selected file 'FileName' is too small for storing the desired number of bytes. => e.g. parameter file with no channel enabled
floating point overflow at sequence 'SeqName' line 'LineNbr'	Floating point value in the line 'LineNbr' of the sequence 'SeqName' is out of range.
font manager error	The fonts required for the sequencer cannot be loaded. => Possible causes: - Wrong Windows® version. => Incorrect Windows® installation. => Insufficient working storage or resources etc.
formula must contain at least one capital X	The entered formula is invalid. A correct input must contain a capital X as the variable.
formular dimension too small	The selected size of the form to be printed is so small that no reasonable output is possible.
gfa[][] overflow at sequence 'SeqName' line 'LineNbr'	In line 'LineNbr' of the sequence 'SeqName' an array containing more than 64 elements was assigned to the global float field.
handshaking error (<ENQ> received after <ENQ>)	Contention error in the communication with the QMS (SECS protocol).
handshaking error (first received char no <ENQ>)	Protocol error when receiving a frame from the QMS.
handshaking error (no <EOT> received after <ENQ>)	Protocol error when establishing a connection to the QMS before transmitting a frame.
hardware COMx not available (or locked by another device)	The selected communication line COMx is not available or already occupied by another device.
hardware not present at sequence 'SeqName' line 'LineNbr'	The hardware option specified in line 'LineNbr' of the sequence 'SeqName' is not installed.

Error/Warning	Description
illegal 'Option' configuration	Invalid hardware configuration.
illegal file name syntax	Invalid syntax in a file name.
illegal file name syntax at sequence 'SeqName' line 'LineNbr'	Invalid syntax in a file name in the line 'LineNbr' of the sequence 'SeqName'.
illegal input format	The input format is invalid.
illegal Option specified	Input of a configuration that is not valid in this place for a given option.
illegal parameter combination	The chosen combination of parameters is not allowed.
incomplete definition	Definition is incomplete (e.g. file name is missing etc.).
incorrect file name	Entered file name is not correct.
incorrect frame length	Protocol error when receiving a frame from the QMS => The specified frame length does not match the transmitted number of bytes.
incorrect password verification	Each password change must be confirmed with a password verification. The password verification was invalid. The old password is still in effect.
input incomplete, not all path names defined	The data defined in the ACCESSC application are incomplete and cannot be stored. => After this error message the user is automatically prompted to define the missing paths (working environment of a user) in the Set Path dialog box.
input must be between 'Low' and 'High'	The specified value is outside the 'Low' and 'High' limits.
integrator error	QMS error message 47. Integrator value < measured background value.
internal standard not defined	No Internal standard has been defined in the current calibration matrix.
invalid argument 'FileName'	Windows® system error, e.g. invalid file positioning etc.
invalid AO identifier	Entered (logical) name of the analog output is not defined.
invalid DI identifier	Entered (logical) name of the digital input is not defined.
invalid DO identifier	Entered (logical) name of the digital output is not defined.
invalid file name identifier	This file type is invalid in this context.

Error/Warning	Description
invalid line	Type of the entered line does not exist.
invalid name	Entered name does not exist.
invalid password	Wrong password has been entered (User-name was correct).
invalid QMS type for PPM	A connected QMS is not suitable for Plasma Process Monitoring (only PPM).
invalid state	Entered status does not exist (at this point).
invalid unit identifier	Entered (logical) name or node number of the Mass Spectrometer Controller does not exist.
invalid unit identifier at sequence 'SeqName' line 'LineNbr'	The mass spectrometer called in the sequence 'SeqName' on line 'LineNbr' is not defined.
invalid username	During the log-in procedure of an application an unknown user name has been entered.
ion source error	QMS error message 36. IS 420 voltage V1 ... V9 is outside the set point tolerance.
ion source set not loaded at sequence 'SeqName' line 'LineNbr'	The access to the ion source parameters on the specified line 'LineNbr' of the sequence 'SeqName' is not possible. The ion source parameters have to be loaded first.
less than 2 cycles stored 'FileName'	Less than 2 cycles are stored in the selected file 'FileName', the desired display/action is not possible.
less than 2 masses defined	Less than 2 mass numbers are defined, the desired display/action is not possible.
line end expected	Syntax error, the selected sequencer line contains invalid entries such as comments, etc. behind the correct line code.
line from type '*' expected	Sequencer symbol '*' for end of comments missing.
line from type 'End' expected	Sequencer line containing the end-of-block instruction ('End' line) missing.
line too long	Entered sequencer line too long.
long string in short packet mode	Error in ArcNet driver; the message to be transmitted is too long for the selected "Short Packet Mode".
mass 'MassNbr' already defined	Selected mass number 'MassNbr' is already defined.

Error/Warning	Description
mass deviation >= 0.5	The mass deviation is greater than 0.5 amu and can no longer be corrected by Quadstar 32-bit. => Selected mass number incorrect or inaccurate mass calibration on the RF generator
mass 'MassNbr' out of adjusted range	The 'MassNbr' mass number entered is outside the range defined in the Measure program under [Calibration] > [Mass Scale].
matrix full	The format (max. 64 mass numbers × 64 components) of the analysis matrix is exceeded.
matrix singular	The matrix is incorrectly defined, the equation system cannot be solved.
matrix contains undefined calibration factors	The analysis matrix makes reference to non-existing calibration factors. => Set the undefined calibration factors of the analysis matrix to zero or calibrate them.
measured value of internal standard is zero, calibration impossible	The value measured on the mass number selected for the Internal Standard is zero; the calibration factors cannot be determined.
measure data dropout, measurement too fast	In measure mode Dynamic Measurement, the QMS-internal data buffer grew completely full. The QMS was stopped, the buffer content was transferred to the PC if possible and then the QMS was restarted. Either the measurement is too fast or the transmission to the PC is too slow.
measure value OVERRANGE	The measured value exceeds the measure range.
memory allocation error	System error during the temporary allocation of working storage by Quadstar 32-bit. => Possibly insufficient working storage available. Terminate all applications that are no longer needed.
memory allocation error in ArcNet server	Not enough memory for ArcNet server. Close all applications that are not absolutely necessary.
memory re-allocation error	System error during the access to temporarily allocated working storage by Quadstar 32-bit. => Possibly insufficient working storage available. Terminate all applications that are no longer needed.
missing scale identifier	There is no ID for the scale defined.
missing text identifier	There is no ID for the text window defined.
missing unit identifier	There is no Mass Spectrometer defined.

Error/Warning	Description
monitor-DAC error	QMS error 26 (service message).
name too long	The entered name contains more than 255 characters.
name too long at sequence 'SeqName' line 'LineNbr'	The name specified in the line 'LineNbr' of the sequence 'SeqName' is longer than 255 characters.
negative Calibration factor determined at 'Comp'/Mass', previous factor kept	A negative calibration factor was calculated for the component 'Comp' on 'Mass'. This factor will not be saved to the calibration factor library and the old calibration factor will be retained.
new password not stored	It was not possible to store the newly entered user password.
no channel on state 'ENABLED'	Parameter setup error; all measurement channels of the selected parameter file are in the OFF or SKIP state or all selected QMU channels in the SEM/Emission control field are switched OFF.
no component name defined	Parameter setup error, no component name defined.
no data available from QMS	No data available from the QMS.
no data selected	Display is only possible if data (data blocks) have been selected.
no error history file found	Printing error, the error history file to be printed does not exist.
no high voltage supply	QMS warning 17. Configuration error; the selected hardware (HV 420/421) does not exist or is not declared.
no ion source supply	QMS warning 18. Configuration error; the selected hardware (IS 420) does not exist or is not declared.
no line from type '/'* allowed	Input of the sequencer symbol '/'* for beginning of comments is not allowed in this place.
no line from type 'Begin' allowed	Input of a sequencer line with an instruction block begin ('Begin' line) is not allowed.
no line from type 'Else' allowed	Input of the sequencer instruction 'Else' is not allowed.
no line from type 'End' allowed	Input of a sequencer line with an instruction block end ('End' line) is not allowed.
no mass number defined	Parameter setup error, no mass numbers defined.
no more channels available	No additional QMS channels are available.

Error/Warning	Description
no more scale entries possible at sequence 'SeqName' line 'LineNbr'	No more entries for the scale are allowed.
no more text entries possible at sequence 'SeqName' line 'LineNbr'	In the specified text window there is no more space for further entries..
no more users allowed	No additional access rights for users can be defined because all 20 possible users have been defined.
no space left on Device	Storing of the data has been cancelled or cannot be started because no space is available on the selected disk drive.
no such file or directory '\Path\File-name'	The file selected in '\Path\FileName' does not exist or is located in a different path.
not enough memory	A memory area requested by Quadstar 32-bit is not available.
not enough memory for measure display	The auxiliary window for displaying the measurement data cannot be opened because insufficient memory is available.
not enough memory for print screen	The working storage required for outputting the selected picture is not available on the system.
not enough memory, reduce rank or select less data	Not enough working storage available on the system. Try to get more storage by closing other applications or reduce the rank of the matrix or select less data.
not enough resources, please close other applications	The computer does not have enough system resources to perform a function. Please close other applications on your computer in order to get more resources.
nothing selected; please check groups to be printed	The print job cannot be executed, because there is no item selected. Check the groups to be printed and print again.
not possible to create the file 'File-Name'	The desired file 'FileName' cannot be created.
not possible to open the file 'File-Name'	The selected file 'FileName' cannot be opened.
no unit enabled	The communication to the Mass Spectrometer cannot be opened because there is no unit set on state enabled.
no unit selected to get the last client identifier	System error. Restart Windows® or contact our service office.

Error/Warning	Description
no user privilege	The menu option or the execution of a command is intentionally inhibited for this user. => The access privileges are granted or inhibited with the ACCESSC program. => If the access privilege is denied, the corresponding options/dialog boxes are automatically inhibited in the program (without message "no user privilege").
No user privilege (Can't identify the QMS xxx type, Prog-Nr check failed)	The type of QMS belonging to this software could not be determined. Contact our service office.
No user privilege (QMS 200/QMS 421 instead of QMS 422 software installed)	Installed Quadstar 32-bit software is set up for use with a QMS421 or QMS200. Install a Quadstar 32-bit software for QMS422. If no effect, contact our service office.
No user privilege (Software key-number out of range)	The read software key number is bad. Restart Windows® or reinstall Quadstar 32-bit. If no effect, contact our service office.
No user privilege (Unknown software type)	The software type could not be figured out. Restart Windows® or reinstall Quadstar 32-bit. If no effect, contact our service office.
no valid scale ID at sequence 'SeqName' line 'LineNbr'.	There is no scale with the specified identifier.
number of calibration factors > 'NbrOfFact'	No more than 'NbrOfFact' calibration factors can be stored for a component.
old version of sequence file 'FileName', please recompile	The sequence 'FileName' has been compiled with a version that is older than the current Quadstar 32-bit version. Please recompile this file in the Sequence Editor.
only one filament available	Filament 2 cannot be selected, because only one filament is available (QMS Warning 21).
parameter 'Group' : 'ParName' changed from 'OldValue' to 'NewValue'	The 'ParName' parameter in the 'Group' group was changed from 'OldValue' to 'NewValue'.
parameter 'Group' : 'ParName' of channel 'Qms-Chan' changed from 'OldValue' to 'NewValue'	The 'ParName' parameter from the 'QmsChan' QMS channel in the 'Group' group was changed from 'OldValue' to 'NewValue'.
parameter lost	QMS error message 41. Value in the NOVRAM lost (battery is discharged or EPROM has been changed).
parameter out of range	The parameter is not within the valid range.

Error/Warning	Description
parameter corrected to current configuration	Warning, Quadstar 32-bit had to adjust the parameter to the changed configuration.
path 'PathName' not existing at sequence 'SeqName' line 'LineNbr'	The path 'PathName' on the specified line 'LineNbr' of the sequence 'SeqName' does not exist.
path name too long	The entered path name is too long and is not accepted.
printer client initialization error	Check whether WINDOWS has still enough free system resources. If no effect, call our service office.
printer currently not available	The printer is busy at the moment. Wait until it gets ready or stop the other printing task and send your task again.
printer not installed or printer error	Printer not installed or general printer error.
printer server not connected	System error. Restart Windows® or contact our service office.
printer server not disconnected	System error. Restart Windows® or contact our service office.
print screen already active	Printing already in progress. => Close Print Screen or wait until the selected picture area has been transmitted to the print manager (or has been printed).
print screen not available	The Print Screen option is blocked at the moment. => The current sequence uses the same printer port as the standard printer and therefore the printer is unavailable for Print Screen.
Process data doesn't match existing file 'FileName' at sequence 'SeqName' line 'LineNbr'	The data created by the Process()-instruction at sequence 'SeqName' on line 'LineNbr' does not match the file 'FileName' (different number of data components). => Correct the data processing.
Process() table too large at sequence 'SeqName' line 'LineNbr'	The memory required by the Process () instruction programmed in the line 'LineNbr' of the sequence 'SeqName' is not available => Simplify the Process() instruction.
QMS control EPROM defect	QMS error message 20. Checksum of the QMS EPROM is incorrect.
QMS control idle error	QMS error message 33. Operating system of the QMS overloaded.

Error/Warning	Description
QMS control buffer RAM defect	QMS error message 23. Working storage defective.
QMS control dual port RAM defect	QMS error message 22. Access to QMS controller defective.
QMS control NOVRAM defect	QMS error message 21, Parameter memory defective.
QMS control stack overflow	QMS error message 17. Internal instruction stack overflow.
QMS control watchdog error	QMS error message 19. Watchdog timer overflow.
QMS error bit'x' set	See QMS 422/421 or QMS 200 operating instructions under "Trouble shooting".
QMS error detected, no data	No measurement data because a QMS hardware error has been detected.
QMS hardware not present	The hardware option is not installed in the QMS or cannot be recognized.
QMS offline	Communication with the QMS cannot be established, the QMS is operating in console mode.
QMS off, not available or wrong node ID configured	Communication with the QMS cannot be established (connection line cut off, QMS not switched on, wrong node ID configured on the QMS, ...).
QMS temperature overrange	QMS error 31 (service message). Rise in temperature inside of the QMS 200. => Cooling insufficient or temperature sensor inside of the QMS 200 faulty. Check if QMS is overheated !
QMS temperature underrange	QMS error 30 (service message). Temperature sensor inside of the QMS 200 faulty or environment frosty.
QMS warning bit'x' set	See QMS/QMI or QMS 200 User's guide, "Warnings".
Quadstar 32-bit installation overflow in ArcNet server	Too many Quadstar 422-Installations. You may run up to eight installations of the same software version in eight different EXE-paths with ArcNet server.
Quadstar 32-bit instance overflow in ArcNet server	ArcNet server allows to run up to 16 Quadstar 32-bit-applications simultaneously that can have communication with QMS's.
read error 'FileName'	Error while reading the file 'FileName'.

Error/Warning	Description
receive checksum error	Transmission error; the checksum of a received frame is incorrect.
reduce border size, print job aborted	The page margins are too large, the page cannot be printed. Reduce the margin size(s) and print again.
reference file not compatible	The reference file doesn't match the measure file, a comparison is not possible. Often caused by measuring with different settings for First Mass, Width, Steps, Mean Value etc.
reset error	QMS error message 18. Reset circuit defective.
residual disk space on drive 'Drive-Name' is 'Number' MB	The available memory on the disk drive 'DriveName' is only 'Number' MBytes.
resolution-DAC error	QMS error 27 (service message).
result too large	Windows® system error. Error in the specification of a file or the path.
RF error	QMS error message 46. Error in the RF stage QMH 400/410.
RF-ID-Error	QMS error message 32. RF print cannot be identified. => Missing print or bad connection.
ringbuffer overflow	Memory error in the QMS data buffer. There is not enough space in the buffer for the number of measured cycles requested (in Display after Measure). => Reduce the number of cycles.
RXD/TXD communication buffer flush error	Error in the Windows® communication driver in an attempt to clear the I/O buffers.
saved zero gas SEM voltage unequal to SEM voltage for mass(es) 'Mass1', 'Mass2',....	The actual SEM voltage does not agree with the SEM voltage the (stored) intensities of the zero gas have been measured with.
scale ID out of range at sequence 'SeqName' line 'LineNbr'	The given scale identifier is out of the allowed range.
SEM error	QMS error message 34. Error in the SEM high voltage.
SEM voltage off	The SEM voltage is not switched on.
send/receive time out detected	Timeout error during data transmission. => The QMS did not respond or did not transmit the required number of bytes within a frame.
sequence editor already open	The sequence editor is already open.

Error/Warning	Description
sequence file 'FileName' not compiled	The 'FileName' sequence file selected has not been compiled.
source and destination file name equal at sequence 'SeqName' line 'LineNbr'	Source and destination file name are equal, therefore the file manipulation cannot be done.
source and destination pathname equal at sequence 'SeqName' line 'LineNbr'.	Source and destination path name are equal, therefore the file manipulation cannot be done.
source and destination pathname not equal at sequence 'SeqName' line 'LineNbr'. Taking source path	Source and destination path name are not equal, therefore the source path is used for the file manipulation.
source or destination file name not valid at sequence 'SeqName' line 'LineNbr'	The file manipulation cannot be done, because one of the file names is not valid.
source path and destination path equal	Source and destination path are equal.
source path and destination path unequal	Error in the path specification.
string too long for selected packet mode	ArcNet driver error; the frame to be transmitted is too long.
syntax error	Syntax error in the entered line or word.
the selected unit is not member of the installation	The selected unit can't be used by the calling Quadstar 32-bit-program because it is already occupied by another Quadstar 32-bit installation. Make sure the defined groups of units/nodes of the different installations do not overlap each other.
too many clients announced	too many Quadstar 32-bit applications running that may use the printer. A maximum of 20 applications is allowed.
too many clocks or timers	A required timer cannot be started because too many timers are in operation (used by other applications).
too many enabled channels in file 'FileName', data not saved	The saving of scan measurement data is permitted in single channel mode only. In the parameter file 'FileName' more than one channel is in state enabled.
too many entries	The sequencer instruction contains more entries than allowed. => Make sure the arguments of the instruction are correct.

Error/Warning	Description
too many interlocked sequences at sequence 'SeqName' line 'LineNbr'	Too many nested sequences in line 'LineNbr' of the sequence 'SeqName'. => Reduce the number of the Sequence() instructions.
too many interlocks	Too many nestings and loops. => Reduce the number of loops.
too many interlocks at sequence 'SeqName' line 'LineNbr'	Too many nested loops, etc. in the 'LineNbr' of the sequence 'SeqName'. => Reduce the number of loops.
too many sequencers started	Too many sequencers running, the printer is limited to cooperate with a maximum of 6 sequencers.
transmission denied (receiver busy)	The attempt to establish a connection with the QMS is being denied due to an overload condition.
transmission wrong (<NAK> received)	A frame transmitted by Quadstar 32-bit is not accepted by the QMS due to transmission errors.
transmission to QMS failed	ArcNet driver error; the QMS signals 'not ready' for a data transmission.
unbalanced number of brackets	The number of opening brackets does not match the number of closing brackets.
undefined data type received	A data frame transmitted by the QMS contains a measurement data identification number that is unknown to Quadstar 32-bit.
unexpected type	The processed variable has an unexpected type.
unit is already used by another instance	The selected unit is already in use of another Quadstar 32-bit application. End the other application or restart Windows®.
unknown client identifier in ArcNet server	ArcNet server is unable to identify the calling Quadstar 32-bit-application. Restart Windows® or contact our service office.
unknown line type	The line contains an unknown sequencer instruction.

Error/Warning	Description
user authorization failure	<p>Attempt to use a program (or program element) that is password protected with Quadstar 32-bit ACCESSC.</p> <p>=> If you have forgotten the user password (for the applications PARSET; MEASURE, DISPSAV, SERVICE, TUNEUP, UTILITY), it can be reset with the ACCESSC program. You can then log in with "Username" and 'Enter' and define a new user password.</p> <p>=> If you have forgotten the ACCESSC Password (for the ACCESSC application), you must copy the Q_ACCESS.DAT file from the installation diskette into the EXE path of ACCESSC before ACCESSC is started. You can then start ACCESSC without a password and define a new ACCESSC password.</p>
username already existing, user not accepted	<p>The newly defined user exists already.</p> <p>=> Define a different user name.</p>
value out of range	The entered value is not compatible with the format, i.e. it is outside the valid limits.
variable index out of range at sequence 'SeqName' line 'LineNbr'	The index for an array selected in line 'LineNbr' of the sequence 'SeqName' is larger than the entered limit (for this index).
window creation error	The desired window can not be opened due to a memory problem.
windows communication driver error	Receive or transmit error in the Windows® communications driver.
write error 'FileName'	Error when writing to the file 'FileName'.
wrong block ID received	A frame transmitted by QMS contains an unknown frame identification number.
wrong communication line selected	The actual sequence instruction is not allowed with the currently opened communication line.
wrong configuration	During the test of a sequence an illegal access to non-existing hardware has been detected.
wrong data type received from QMS	The data frame transmitted by QMS contains measurement data of a wrong type.
wrong version of identification file (Q_OPTION.DAT)	<p>The Quadstar 32-bit identification file does not correspond to the current version.</p> <p>=> Copy the original file from the installation diskette.</p>

10.7 Troubleshooting


Problem	Possible causes	Solution
The measurement produces very small intensities, although the parameters are set reasonably.	<ul style="list-style-type: none"> Emission and/or SEM switched off. Zero Gas Subtraction ON, the currently measured spectrum has been measured as Zero Gas by mistake. 	<ul style="list-style-type: none"> Switch Emission and/or SEM on (CTRL+S). Make a new Zero Gas Measurement by using a real zero gas. Check in program Dispsav (Auxiliary : Zero Gas), whether the Zero Gas values are still valid; if Not, delete them and save again (File : Save). Switch the Zero Gas Subtraction off (Parameters : Setup).
A certain measurement produces various intensities, depending on the measure time (Speed/Dwell).	<ul style="list-style-type: none"> Pause is not set correctly. 	<ul style="list-style-type: none"> Set the parameter Pause in Channel Parameter Editor under Amplifier to a safe value (normally 1.0).
A certain measurement produces under various conditions always the same intensities.	<ul style="list-style-type: none"> Simulation mode (SIM-symbol in status bar) is active. 	<ul style="list-style-type: none"> Switch off the simulation mode under Parameter Setup: Setup : General : Simulation Mode : OFF.
The measure curve of the QMS 422/200 seems to be shifted against the baseline.	<ul style="list-style-type: none"> Offset of the measure amplifier not calibrated or calibration invalid. 	<ul style="list-style-type: none"> Recalibrate the measure amplifier. See Chapter 2.8.3 QMS Offset, ¶ 140.
The QMS does not measure on the specified masses, but on slightly different masses.	<ul style="list-style-type: none"> Mass Scale Correction is switched on. 	<ul style="list-style-type: none"> Switch off the Mass Scale Correction. See Chapter 2.8.1 Mass scale calibration (Mass Scale), ¶ 134 or Setup of the current measure mode.
Emission/SEM can't be switched on in Quadstar 32-bit.	<ul style="list-style-type: none"> Emission External/HV External option is active. Protection (SEM+Fil)-parameter is set incorrectly. 	<ul style="list-style-type: none"> Switch off Emission External/HV External. See Chapter 3.8.1 QMS Configuration, ¶ 234. If you are sure, that emission and SEM are not controlled or protected externally, then switch over to INTERN. See Chapter 3.8.1 QMS Configuration, ¶ 234; otherwise check the vacuum control unit.
The Analog outputs of the QMS (AO's) apparently don't work.	<ul style="list-style-type: none"> AO-Characteristic-Curve has still default values; it doesn't match the required range. 	<ul style="list-style-type: none"> Set up the AO-Characteristic Curve (program Parameter Setup : Config : AO Characteristic Curve) in such a way, that the range (e.g. 0 ... 100%) is assigned to the AO range (e.g. - 5.12V ... +5.12V) adequately.
The printer could not print the lines of the sequence. Where can I find the 'lost' information?	<ul style="list-style-type: none"> Paper jam, out of paper, connection to the printer interrupted etc. 	<ul style="list-style-type: none"> Force the remaining lines to be output by using the Print-Manager. See Chapter 8.4.3 Log output, ¶ 422

Problem	Possible causes	Solution
If I print a screen hardcopy by CTRL+SHIFT+F12 while a measurement is in progress, the error message Send/Receive Timeout Detected appears. How can I avoid that?	<ul style="list-style-type: none"> The printing requires too much CPU time, the interface to the QMS cannot be served any more and therefore reports a timeout error. 	<ul style="list-style-type: none"> Windows® problem, unsolvable. Choose the print region as small as possible. Store the graphic in a file first and print it after the measurement is done. Try other equipment (printer, PC).
I want to print on a network printer from Quadstar 32-bit, but instead of that, the printer on LPT1 or none at all is printing.	<ul style="list-style-type: none"> Quadstar 32-bit identifies printers exclusively by ports (LPT1 ... 3, COM1 ... 4). Therefore, the network path can not be found. 	<ul style="list-style-type: none"> Define a virtual printer port for the used network path.
The printed items (parameters, matrix elements etc.) are not equal to the ones I see on the screen. What's wrong?	<ul style="list-style-type: none"> Some items can only be printed correctly by using files as data base. 	<ul style="list-style-type: none"> Save the corresponding files before printing them.
Quadstar 32-bit shows everything on the screen in monochrome color or supplies no more than two colors under Config : Screen Color to choose from, although my graphics system supplies much more colors. What's wrong?	<ul style="list-style-type: none"> In Parameter Setup program under Config : Screen Color the color palette Black&White is selected. 	<ul style="list-style-type: none"> Select a different color palette.

11. Residual Gas Analysis Program for Prisma™

The RGA Program supplied with the mass spectrometer Prisma™ allows for quick and easy measurement of residual gases in vacuum chambers:

- The program leads you, step by step, to the desired results.
- For all common measuring tasks there are ready-to-use set up measurements.
- The RGA Program is easy to be used via dialog boxes and leads rapidly to the desired results (overview spectra, measuring tables, gas concentrations, partial pressures etc.).
- Auxiliary programs allow for fine tuning of the Prisma™ mass spectrometer.

This RGA Program was developed as a Quadstar 32-bit sequence for residual gas analysis when operating the Prisma™. It consists of a main sequence and a number of sub-sequences. (A sequence is a program, that has been written with the sequencer, i.e. in the macro programming language of Quadstar 32-bit; see Chapter 8. Process Control,  409).

11.1 Overview of the operating modes

11.1.1 Main selection box

In the main selection box you can choose the type of measurement and start some auxiliary programs for fine tuning of the mass spectrometer:




Fig. 11-1

Scan

Measure a spectrum across a given mass range:

- [Analog] - Display an analog spectrum.
- [Bargraph] - Search for peaks in the spectrum and display them.
- [Stair] - Measure only on integer mass numbers.

Currents

Measure the ion currents on predefined mass numbers:

- [Table] - The intensities are displayed as a table.
- [Bargraph] - The intensities are displayed as a bar chart.
- [Vs. Time] - The intensities are displayed Versus Time.

Concentrations

Determine the concentrations of the major components of a typical residual gas (Ar, CO₂, H₂, H₂O, He, C_xH_y, O₂ and N₂/CO):

- [Table] - The concentrations are displayed as a table.
- [Bargraph] - The concentrations are displayed as a bar chart.
- [Vs. Time] - The concentrations are displayed Versus Time.

Part. Press.

Determine the partial pressures of the major components of a typical residual gas (Ar, CO₂, H₂, H₂O, He, C_xH_y, O₂ and N₂/CO):

- [Table] - The partial pressures are displayed as a table.
- [Bargraph] - The partial pressures are displayed as a bar chart.
- [Vs. Time] - The partial pressures are displayed Versus Time.

Auxiliaries

Auxiliary programs for adjustment of the mass spectrometer:

- [Mass Scale] - Adjust the mass scale.
- [Backgr.] - Determine the residual gas background.
- [Leak] - Perform a helium leak test.
- [Resolution] - Change the mass resolution (smaller number = sharper discrimination of peaks).
- [Cal. SEM] - Optimize the SEM voltage (in such a way, that the gain is 2000).
- [Cal. Offset] - Calibrate the amplifier offset.
- [Cal. Sens] - Determine the mass spectrometer sensitivity for partial pressure measurements.

11.1.1.1

Scan

Scan offers the measurement types Analog, Bargraph and Stair to choose from:




Fig. 11-2

11.1.1.2

Currents

Currents determines ion currents on predefined mass numbers. The results can be displayed as table, bar chart or Versus Time:




Fig. 11-3

11.1.1.3 Concentr.

Concentr. determines concentrations of gas components. The results can be displayed as table, bar chart or Versus Time:




Fig. 11-4

11.1.1.4 Part. Press.

Part. Press. determines partial pressures of gas components. The results can be displayed as table, bar chart or Versus Time:




Fig. 11-5

11.1.1.5 Auxiliaries

The auxiliary programs serve for adjusting the mass spectrometer:




Fig. 11-6

Mass Scale
Compensate mass scale deviations.

Backgr.

Determine residual gas background.

Leak

Perform a helium leak test.

Resolution

Change the resolution.

Cal. SEM

Optimize the SEM voltage.

Cal. Offset

Compensate the amplifier offset.

Cal. Sens.

Determine the mass spectrometer sensitivity (for partial pressure measurements).

11.2 Working with the RGA program

11.2.1 Start of the RGA program

The RGA program is delivered in compiled form and is automatically installed during the installation of Quadstar 32-bit for a Prisma™:




Fig. 11-7

Double clicking on the icon RGA Program starts the program and shows the dialog box Introduction.

NOTE:

If Quadstar 32-bit reports an error instead, e.g. because the installation path is not equal to the standard path, you should start the Parset program and load a sequence by [Sequence] > [Editor]. Then choose [File] > [Compile All] and click on [Compile All], so all sequences - i.e. the RGA-program as well- are recompiled. Now you can start the RGA-program as described above.

If it still doesn't work, then you should study the sections 'Installation/Configuring Quadstar 32-bit' in Chapter 1. Introduction, 17.

11.2.2 Introduction

First you have to state a few properties of your Prisma™:




Fig. 11-8

Mass Range

Mass range of your Prisma™ detector.

Detector Type

Type of your Prisma™ detector.

If the settings are correct, click on [Accept selection and proceed]. The program now performs an automatic amplifier offset correction which takes a few seconds.

11.2.3 Switching on the emission

If the emission is switched off, you are asked to switch it on:




Fig. 11-9

Click [Yes, start]





Fig. 11-10

	⚠ CAUTION
	<p>Make sure the total pressure in the vacuum chamber is low enough before activating the button [Yes] in the second dialog box. The ion source may be damaged if the pressure is higher than 10E-4 mbar.</p>

11.2.4 Display of an overview spectrum

In the following dialog box you can choose whether an overview spectrum is to be measured and displayed or not; if not, you can proceed to the main selection box directly:




Fig. 11-11

[Start Spectrum]

Choose [Start Spectrum] to display the spectrum from mass 0 ... 50.




Fig. 11-12

When the measurement is finished, a dialog box containing an [OK] button appears. If you want to go on to the main selection box, click on the [OK] button.

[EXIT SEQUENCE]

Exit the RGA program.

[Proceed to Selection Boxes]

Continue to the main selection box.

11.2.5 Switching on the Channeltron

If you have chosen Channeltron-/Faraday detector in the Introduction and the SEM voltage has not been switched on yet, the following dialog box appears:





Fig. 11-13

If the total pressure in the vacuum chamber is below $1 \times 10E-5$ mbar, click on 'Yes'. The SEM voltage of the channeltron is then switched on (the SEM symbol turns to yellow). If the total pressure is higher than $1 \times 10E-5$ mbar, don't turn on the SEM voltage; click on 'No' to proceed. Your Prisma™ will work like one equipped with a Faraday detector in this case.

	⚠ CAUTION
	Switching on the SEM voltage at a total pressure $>10E-5$ mbar may cause damage to the Channeltron.

11.2.5.1 Calibrating the Channeltron

After the SEM voltage has been switched on, Quadstar 32-bit determines the high voltage necessary for a gain of 2000 of the output current compared to the current measured with the Faraday cup. This calibration takes about 40 seconds. The following message is displayed during this time:




Fig. 11-14

The determined high voltage is displayed in the status bar afterwards:



Fig. 11-15

11.3 Main Selection Box

You have now reached the Main Selection Box. Here you can start the different types of measurements and auxiliary programs you have learned about previously:




Fig. 11-16

11.3.1 Terminate RGA program [EXIT SEQUENCE]

Most dialog boxes contain the [EXIT SEQUENCE] button. This allows for exiting the program. Quadstar 32-bit asks however, whether the program is really to be terminated:




Fig. 11-17

- Choose [Yes], if you want to terminate the program. The settings you have made already (e.g. emission switched on) will remain.

- Choose [No], if you want to continue the program; in that case, the next dialog box will be opened.

11.3.2 Back to previous screen

Within the RGA program you can step back to the previous screen at any time by pressing <SHIFT> and <F4> simultaneously or by selecting the name of the previous screen in the menu File.

11.3.3 Measurement programs

11.3.3.1 Scan measurements

[Analog]

Analog spectrum across the given mass range.

[Bargraph]

Peak search across the given mass range

[Stair]

The intensities on all integer masses of the given mass range are measured.

11.3.3.1.1 Enter measure and display parameters

If you click on a scan measurement (Analog, Bargraph or Stair), the following dialog box appears:

Measurement Speed	First Mass	Width
Medium	0	50
Max. Decade	Number of Decades	Y-Scale
E-8 A	6	Log
Save Data	Detector	
OFF	Faraday	

EXIT O.K.

Press Shift + F4 to return to Main Selection Menu

Fig. 11-18

Here you can enter the most important properties of the measurement and the data display:

Measurement Speed

Speed with which the samples are taken.

First Mass

Start mass, at which the measurement begins.

Width

Mass range from the start mass onward.

Max. Decade

Upper limit of the display.

Number of Decades

Number of decades to be displayed (when logarithmic display is chosen).

Y-Scale

Display of the Y-Scale (linear or logarithmic).

Save Data

Store the measured data on/off.

Detector

Operation mode of the detector (Faraday or Channeltron).

11.3.3.1.2 Storage of measured data

If you have chosen ON under Save Data in the dialog box above, the measured data will be stored.




Fig. 11-19

For that, you have to specify the name and path of the file to contain the data. You may enter them manually or choose them by the file manager [File Manager].




Fig. 11-20

When you have chosen a file and clicked on the [OPEN] button, the file name and path will be inserted into the input field of the dialog box Data File Selection:





Fig. 11-21

	⚠ CAUTION
	<p>If you choose an already existing file, the newly measured data must be compatible with the already stored ones (i.e., at least First Mass and Width must be equal). If this is not granted, the file may be destroyed.</p>

11.3.3.1.3 Analog

[Analog]

The complete mass range is measured through continuously and the results are displayed as an analog curve.




Fig. 11-22

11.3.3.1.4 Bargraph

[Bargraph]

The given mass range is searched for peaks. The maximum values of the found peaks are displayed over their masses.




Fig. 11-23

11.3.3.1.5 Stair

[Stair]

On all integer masses of the specified mass range the intensities are measured and displayed.




Fig. 11-24

11.3.3.2 Currents

In these measurements, the ion currents are determined on all given mass numbers and displayed as table, bar chart or Versus Time.

[Table]

Table display.

[Bargraph]

Bar chart display.

[Vs. Time]

Display Versus Time.

When you have chosen a display mode, a dialog box appears, where you can define a few things like measurement parameters, display properties and data storage if necessary. See Chapter 2.4 Scan menu, 76.

11.3.3.2.1 Table

[Table]

The measured intensities are displayed as a table. The entries on the left of the table mean:

Nbr

Number of the measured value

Type

Type of the measured value (here: Cur., i.e. ion current)

Ident

Identification of the measured value (here: mass number)

Unit

Unit of measure (here: A, i.e. Ampere)

[M1] QUADSTAR 32-bit Measurement - [Scan Bargraph < scsf.sbp >]

Measurement Number: 5				MID	MID	MID	MID
Nbr	Type	Ident	Unit	midc.mip 17:18:39	midc.mip 17:18:31	midc.mip 17:18:34	midc.mip 17:18:36
0	Cur.	2.00	A	2.977E-10	2.985E-10	2.962E-10	2.964E-10
1	Cur.	4.00	A	5.775E-14	5.333E-14	6.168E-14	6.470E-14
2	Cur.	14.00	A	9.374E-12	1.085E-11	1.038E-11	9.959E-12
3	Cur.	16.00	A	1.072E-10	1.100E-10	1.085E-10	1.060E-10
4	Cur.	18.00	A	1.333E-08	1.378E-08	1.362E-08	1.350E-08
5	Cur.	28.00	A	2.303E-10	2.483E-10	2.428E-10	2.355E-10
6	Cur.	32.00	A	2.160E-11	2.403E-11	2.320E-11	2.301E-11
7	Cur.	40.00	A	1.433E-09	1.428E-09	1.429E-09	1.433E-09
8	Cur.	44.00	A	1.096E-10	1.141E-10	1.118E-10	1.094E-10
9	Cur.	55.00	A	1.350E-12	1.498E-12	1.469E-12	1.455E-12
10	Cur.	57.00	A	7.154E-13	6.625E-13	7.349E-13	7.049E-13
11							
12							

Press SHIFT + F4 to return to previous screen

Fig. 11-25

On the right of the table the new measure data are displayed alternatingly.

11.3.3.2.2 Bargraph

[Bargraph]

The measured ion currents are displayed as bar chart:




Fig. 11-26

11.3.3.2.3 Versus Time

[Vs. Time]

The measured ion currents are displayed as a curve Versus Time. On the right of the screen the current values are displayed numerically as well:




Fig. 11-27

11.3.3.3 Concentr.

In these measurements, the concentrations of certain gases in the vacuum chamber are determined and displayed as table, bar chart or Versus Time.

[Table]

Table display.

[Bargraph]

Bar chart display.

[Vs. Time]

Display Versus Time.

When you have chosen a display mode, a dialog box appears, where you can define a few things like measurement parameters, display properties and data storage if necessary. See Chapter 2.4 Scan menu, 76.

11.3.3.3.1 Table

[Table]

The concentrations are displayed as a table:

Measurement Number: 13				MCD	MCD	MCD	MCD
Nbr	Type	Ident	Unit	mdc.mcp 17:50:31	mdc.mcp 17:50:22	mdc.mcp 17:50:25	mdc.mcp 17:50:28
0	Conc.	Ar		69.501%	69.476%	69.442%	69.401%
1	Conc.	CO2		4165.234ppm	4095.518ppm	4021.052ppm	4126.750ppm
2	Conc.	H2		3.317%	3.270%	3.300%	3.345%
3	Conc.	H2O		25.583%	25.693%	25.697%	25.658%
4	Conc.	CxHy		699.674ppm	537.784ppm	523.342ppm	585.370ppm
5	Conc.	N2 / CO		1.005%	9962.084ppm	9956.739ppm	1.011%
6	Conc.	O2		1032.170ppm	951.417ppm	1043.714ppm	1050.766ppm
7	Conc.	He		48.763ppm	60.304ppm	60.768ppm	85.036ppm
8							
9							
10							
11							
12							

Fig. 11-28

11.3.3.3.2 Bargraph

[Bargraph]

The concentrations are displayed as a bar chart:




Fig. 11-29

11.3.3.3 Versus Time

[Vs. Time]

The concentrations are displayed as a curve Versus Time. On the right of the screen the current values are displayed numerically as well:




Fig. 11-30

11.3.3.4 Part. Press. (Partial Pressures)

In these measurements, the partial pressures of certain gases in the vacuum chamber are determined and displayed as table, bar chart or Versus Time. To make that possible, the auxiliary program 'Cal. Sens.' must be carried out first and the current total pressure be given.

[Table]

Table display.

[Bargraph]

Bar chart display.

[Vs. Time]

Display Versus Time.

When you have chosen a display mode, a dialog box appears, where you can define a few things like measurement parameters, display properties and data storage if necessary. See Chapter 2.4 Scan menu, 76.

11.3.3.4.1 Table

[Table]

The partial pressures are displayed as a table:

Measurement Number: 62				Process	Process	Process	Process
Nbr	Type	Ident	Unit	18:04:50	18:04:57	18:05:04	18:05:11
0	Part.Pres	TOTAL	mbar	6.247E-09	6.262E-09	6.255E-09	6.261E-09
1	Part.Pres	Ar	mbar	4.809E-09	4.834E-09	4.840E-09	4.842E-09
2	Part.Pres	CO2	mbar	1.936E-11	1.879E-11	1.897E-11	2.045E-11
3	Part.Pres	H2	mbar	2.029E-10	1.826E-10	1.757E-10	1.782E-10
4	Part.Pres	H2O	mbar	1.153E-09	1.162E-09	1.157E-09	1.158E-09
5	Part.Pres	CxHy	mbar	3.561E-12	3.151E-12	3.659E-12	3.989E-12
6	Part.Pres	N2 / CO	mbar	5.396E-11	5.416E-11	5.106E-11	5.184E-11
7	Part.Pres	O2	mbar	5.251E-12	5.537E-12	5.913E-12	5.491E-12
8	Part.Pres	He	mbar	6.796E-13	1.749E-12	1.790E-12	1.784E-12
9							
10							
11							
12							

Fig. 11-31

11.3.3.4.2 Bargraph

[Bargraph]

The partial pressures are displayed as a bar chart:




Fig. 11-32

11.3.3.4.3 Versus Time

[Vs. Time]

The partial pressures are displayed as a curve Versus Time. On the right of the screen the current values are displayed numerically as well:




Fig. 11-33

11.3.3.5 Auxiliaries

You can choose among the following auxiliary programs:

[Mass Scale]

Calibrate the mass scale.

[Backgr.]

Determine the residual gas background.

[Leak]

Perform a helium leak test.

[Resolution]

Change the mass resolution (smaller number = sharper discrimination of peaks).

[Cal. SEM]

Calibrate the SEM voltage (in such a way, that the gain is 2000).

[Cal. Offset]

Calibrate the amplifier offset.

[Cal. Sens.]

Determine the mass spectrometer sensitivity for partial pressure measurements.

11.3.3.5.1 Calibrate the mass scale

[Mass Scale]

For accurate measurements of ion currents, concentrations and partial pressures it is necessary to always measure at the peak maximum. Therefore a calibration of the mass scale must be done which corrects deviations up to ± 0.5 amu from the nominal mass.




Fig. 11-34

11.3.3.5.2 Determine the residual gas background

[Backgr.]

In every vacuum chamber there are fairly constant gas parts, that do not come from the inlet gas sample; they falsify the measurement results. Those parts can be measured by pumping the vacuum chamber to its base pressure. In the following measurements those parts are subtracted in order to get correct measurement results.

You can, if desired, merely delete the values of the old background measurement. This may be indicated if you expect that the old values are not accurate any more (e.g. because you are using channeltron instead of faraday detector) but you can't or don't want to remeasure the background.




Fig. 11-35

[Reset old background values]

Delete old, inaccurate values

[Determine actual background values]

Measure new background values

Choose the type of detector, with which the background is to be measured:





Fig. 11-36

	⚠ CAUTION
	If you choose a different detector later or if the SEM voltage changes, the previous background values will become invalid.

You are then asked, whether the system is at its base pressure (the lowest possible pressure in the vacuum chamber):




Fig. 11-37

[YES]

The background is remeasured and the old values are overwritten.




Fig. 11-38

[NO]

Pump your system to its base pressure and restart the measurement afterwards:




Fig. 11-39

11.3.3.5.3 Perform a helium leak test

[Leak]

To find leaks in the vacuum system, the intensity of helium is measured. If you apply helium near the leak, Quadstar 32-bit reacts optically and acoustically to the higher intensity: the measurement curve shows a higher level and the signal tone rises.




Fig. 11-40

Choose the length of the time axis in such a way, that the whole search for the leak is recorded. After that, the leak test starts:




Fig. 11-41

11.3.3.5.4 Change the mass resolution

[Resolution]

This auxiliary program allows for change of the peak discrimination. That may be necessary, if peaks standing close to each other can not be discriminated (=> decrease resolution value) or if you need maximum sensitivity (=> increase resolution value).




Fig. 11-42

11.3.3.5.5 Calibrate the SEM voltage

[Cal. SEM]

To give the SEM of a Prisma™ with channeltron detector optimal working conditions, its voltage is set to a value which provides a gain of about 2000. The calibration process and the set voltage are displayed in the status bar.



Fig. 11-43

11.3.3.5.6 Calibrate the amplifier offset

[Cal. Offset]

On request, the Prisma™ determines by itself all necessary correction values to compensate offsets of the measure amplifier under different conditions.




Fig. 11-44

11.3.3.5.7 Determine the mass spectrometer sensitivity

[Cal. Sens.]

Determine the mass spectrometer sensitivity for partial pressure measurements.

In order to perform partial pressure measurements with the Prisma™, its sensitivity must be determined first. For that purpose, declare the current total pressure as it is measured by a total pressure gauge (e.g. PKR 250/251). Quadstar 32-bit calculates the sensitivity of the Prisma™ for the individual gases then.




Fig. 11-45

If you have entered the pressure value, press [OK]:




Fig. 11-46

NOTE:

The RGA program was developed using the Quadstar 32-bit Process Control (Sequencer). The source code of this sequence(s) is supplied and can be adapted to new tasks. Because the source code is write-protected, a copy must be made first and one has to work with the copy. This way you can create new applications without extensive programming work.

Writing sequences is explained in detail in Chapter 8. Process Control, 409.

12. Sputter Process Monitor Program

The routines delivered together with the SPM 200 Sputter Process Monitor allow for an easy measurement of process gases and residual gases in sputter process chambers.

- The parameter settings needed for solving the common problems are integrated already
- The SPM 200 Program is easily operated via dialog-fields and delivers results quickly
- Additional Auxiliary Programs enable a fine-tuning of the SPM 200

The SPM 200 Program is a Quadstar 32-bit sequence developed for the operation of the SPM 200. It consists out of a main sequence and several subroutines which are edited within the Process Control Module of Quadstar 32-bit. The different measurement modes are based on specific parameter files edited within the Parameter Setup program of Quadstar 32-bit.

12.1 Summary of the different operation modes

12.1.1 The Measurement Selection Box

The Measurement Selection box is the main selection box within the SPM 200 program. From this box different measurement modes can be selected and started.




Fig. 12-1

Argon Sputter Gas

As Argon is one of the most widely used sputter gases, several routines are included to measure the impurities within that gas.

Measured components: Ar, CO₂, H₂, H₂O, N₂ and O₂

Processes: e.g. Al, Ti, W- sputtering

For more information see Chapter 12.3 Impurity determinations, 616

- [High speed measurement]
Approx. cycle time: 2 s
For very fast determination of the trace impurities in Ar-sputter gas.
- [Medium speed measurement]
Approx. cycle time: 12 s
For medium speed determination of the trace impurities in Ar-sputter gas.
- [Low speed measurement]
Approx. cycle time: 25 s
For low speed determination of the trace impurities in Ar-sputter gas.

Others

Several routines are included for the determination of impurities in various gas mixtures and pure Nitrogen.

For more information see Chapter 12.3 Impurity determinations, 616.

- [Argon and Nitrogen]
Processes: e.g. TiN
Approx. cycle time: 2 s
- [Argon, Nitrogen and Oxygen]
Approx. cycle time: 2 s
- [Nitrogen]
Approx. cycle time: 10 s

Auxiliary Routines

The routines within this field include parameter settings mostly used during troubleshooting at the process chamber or for defining the status of the plant.

- [Residual Gas Analysis]
Measured components: Ar, CH_x (Hydrocarbons), CO₂, H₂, H₂O, N₂ and O₂
Approx. cycle time: 2 s
For more information see Chapter 12.4 Residual Gas Measurements, 619.
- [Helium Leaktest]
Performs a Helium leaktest while continuously measuring and displaying the ion current on mass 4 amu
- [Scan Analog]
A scan within the mass range 0 ... 50 amu is displayed.
Scan rate: 2 s/amu
- [Offset Determination]
The offset of the SPM 200 electrometer amplifier is determined for different ranges and measurement speeds. The values obtained are included in all following measurements. The measurement takes about 10 s.
- [Mass Scale Adjust]
The exact peak positions of several mass numbers are defined. Peak positions of mass numbers not included in this measurement will be inter-/extrapolated. The measurement takes about 20 s.

NOTE:

The [Mass Scale Adjust] routine has to be run at or close to the base pressure of the process chamber. At pressures above 1×10^{-4} mbar the peak positions of low mass peaks can not be defined precisely any more. It is recommended to run this routine before the first measurements are performed, routinely about once a week.

12.2 General functions of the SPM 200 sequence

12.2.1 Starting the SPM 200 sequence

The SPM 200 sequence is delivered in a version compiled already together with the Quadstar 32-bit Measure program. It will be installed automatically with a separate icon during the installation of Quadstar 32-bit:



Fig. 12-2

After double-clicking the SPM 200 Program icon, the SPM 200 sequence is started and the dialog window Introduction will be displayed.




Fig. 12-3

12.2.2 Leaving the SPM 200 sequence

Most of the windows within the SPM 200 sequence contain the button [EXIT SEQUENCE]. This allows for leaving the sequence in a defined way. The actuation of the [EXIT SEQUENCE] button has to be confirmed.




Fig. 12-4

If the button [NO] is pressed, the sequence will continue with the next dialog window.

If the [YES] button is pressed, the sequence will be stopped and the Measure program will be closed.

12.2.3 Jumping backwards within the SPM 200 sequence

Within the SPM 200 sequence it is possible always to jump backwards to a previous screen by simultaneously pressing the <SHIFT> and the <F4> key. Alternatively select <File> from the menu and select the text of the last window displayed.

12.2.4 Switching on the filament and the emission

After starting the SPM 200 sequence the emission status is checked. If the emission is not running already, the following windows will be displayed.




Fig. 12-5





Fig. 12-6

⚠ CAUTION	
	<p>Before switching on the Emission, check the total pressure in the process chamber, the SPM 200 is connected to. The total pressure has to be below 1×10^{-2} mbar to start the emission without any danger for the filament.</p> <p>If the SPM 200 is separated by an isolation valve from the process chamber, the SPM 200 pumping station has to be running for at least 30 minutes before switching on the emission to allow for an evacuation of the volume between the ion source and the isolation valve.</p>

12.2.5 Switching on the Channeltron Voltage

After starting the SPM 200 sequence the Channeltron Voltage status is checked. If the Channeltron Voltage is not switched on already, the following window will be displayed.





Fig. 12-7

⚠ CAUTION	
	<p>Before switching on the Channeltron Voltage, check the total pressure in the process chamber, the SPM 200 is connected to. The total pressure has to be below 1×10^{-2} mbar to start the Channeltron Voltage without any danger for the channeltron.</p> <p>If the SPM 200 is separated by an isolation valve from the process chamber, the SPM 200 pumping station has to be running for at least 30 minutes before switching on the Channeltron Voltage to allow for an evacuation of the volume between the ion source and the isolation valve.</p>

12.2.6 Calibration of the Channeltron Voltage

The Channeltron Voltage is adjusted to a value, that gives an amplification of 2.000 relative to the Faraday-data. This adjustment allows for comparable results independent on the aging of the Channeltron.

The amplification adjustment is performed by a comparison of the ion currents of measurements in the Integral Scan mode, as long as the SPM 200 is measuring the residual gas composition only. If Argon is detected, the adjustment will be performed by measurements on mass 36 amu. The adjustment of the Channeltron Voltage takes about 40s. During this time the following information is displayed on the screen.




Fig. 12-8

The resulting Channeltron Voltage is displayed after completion of the adjustment routine.



Fig. 12-9

12.2.7 Display of an Overview Spectrum

An overview scan can be started immediately after the determination of the Channeltron High Voltage is finished.




Fig. 12-10

An analog scan will be started after pressing the Start Spectrum button.




Fig. 12-11

12.2.8 Saving of measurement data

Measurement data of Impurity Determinations, Residual Gas Measurements and Analog Scans can be saved to the hard disk or the floppy disk drive.




Fig. 12-12

12.2.8.1 Saving data from Analog Scans

The following window will be displayed only, if data from Analog Scans have to be stored. Storage of 1, 3, 10 or all scans is enabled. If data from more than one scan are stored, it will be possible to display these data in a 3-dimensional graph.

The measurement and display of Analog Scans will be repeated independently on the selection until the measurement is stopped manually.




Fig. 12-13

12.2.8.2 Assigning data file names

A name has to be assigned to each data file before the storage can start. The file name can either be defined manually by typing it into the edit field or by selecting a file with the file manager.




Fig. 12-14

12.2.8.3 File Manager

The File Manager allows for the definition of the file name, the directory and the drive.




Fig. 12-15

NOTE:

The functions of the File Manager are described more extensively in Chapter 10.1 File Manager, 541.

The file name selected within the File Manager will be displayed in the Data File Selection window, after it has been confirmed by pressing the [Open] button.





Fig. 12-16

	⚠ CAUTION
	Measurement data that shall be appended to an already existing file have to be compatible with the data already stored. Compatibility has to be tested with regard to names and numbers of gas components, mass numbers and analog input data. Any incompatibility will destroy the data files!

12.3 Impurity determinations

The impurities in sputter gases are determined in a concentration mode as %, ppm or ppb. The differing sensitivities for the gases Ar, CO₂, N₂, O₂, H₂O and H₂ are taken into account.

12.3.1 Selection of different Display Modes

Data obtained during the determination of impurities in a sputter gas can be displayed either in Tabular form or in a Versus Time mode. The display mode can be selected from the following window.




Fig. 12-17

Choosing the display mode Tabular results in a multicolumn display showing the component names, the number of measurements, the measurement parameter file name, the time of the measurements and the results obtained.

Display of data from Impurity Determinations in Tabular form

Measurement Number: 4				MCD argf.mcp 16:10:58	MCD argf.mcp 16:11:01	MCD argf.mcp 16:11:04	MCD argf.mcp 16:11:06
Nbr	Type	Ident	Unit				
0	Conc.	Ar		99.964%	99.964%	99.964%	99.964%
1	Conc.	CO ₂		9.476ppm	9.603ppm	9.788ppm	9.552ppm
2	Conc.	H ₂		126.879ppm	126.862ppm	127.214ppm	127.965ppm
3	Conc.	H ₂ O		104.833ppm	104.711ppm	105.432ppm	105.631ppm
4	Conc.	N ₂ / CO		108.471ppm	108.134ppm	108.172ppm	108.720ppm
5	Conc.	O ₂		3.741ppm	3.792ppm	3.550ppm	3.664ppm
6							
7							
8							
9							
10							
11							
12							

Fig. 12-18

Choosing the display mode Versus Time results in a display in graphical mode showing the measured data in a concentration/time graph. The data of the last measured cycle are displayed in tabular form additionally.

The length of the time axis can be adapted to the current application.




Fig. 12-19

Display of data from Impurity Determinations in Versus Time mode.




Fig. 12-20

12.3.1.1 Adapting the concentration scale for Versus Time displays

The concentration scale is set from 100ppb to 100% as a standard. If a different scale should be selected, proceed as follows:

Procedure

- 1 Start the program PARSET
- 2 Select MEASURE
- 3 Select MCD
- 4 Select the desired file from the file manager, the analysis matrix of this routine will be displayed.

- 5 Select EDIT
- 6 Select CONCENTRATION PARAMETERS
- 7 Define the Disp. F.S.R.
- 8 Define the Disp. Decades
- 9 Close the windows and save the parameter changes.

12.4 Residual Gas Measurements

Two different measurement modes are possible for the determination of the residual gas composition within the process chamber:

- The partial pressures of the residual gases can be determined combining the data from the mass spectrometer and the data from a total pressure gauge, if the analog output voltage of a total pressure gauge is connected to the Analog Input AI#0 of the SPM 200.

The analog output characteristics of the total pressure gauge has to be defined in the corresponding sequence RESGAS2.SEQ, line 27. If the pressure/output characteristic cannot be given in a mathematical formula, the characteristic has to be defined as an AI-Characteristic curve.

Procedure

Proceed as follows:

- 1 Start the Parameter Setup program
 - 2 Select Config
 - 3 Select AI Characteristic Curve
 - 4 Define as many as possible Basic Points for the definition of the output characteristic
 - 5 Check the characteristic curve using the Diagram button
 - 6 Close the edit field by pressing the [OK]- button
 - 7 Change line 27 in he sequence RESGAS2.SQE to a comment-line
- The partial pressures of the residual gases will be determined via the absolute sensitivity of the SPM 200, as long as no signals from a total pressure gauge are available.

The following window gives the opportunity for the selection between these two modes.




Fig. 12-21

The Residual Gas composition can be determined correctly only, if the pressure in the process chamber is below 5×10^{-5} mbar. If the pressure is above this value, overflow of ion current measurements may occur and meaningless results would be obtained.




Fig. 12-22

12.4.1 Selection of different Display Modes

Data obtained from Residual Gas measurements can be displayed either in Tabular form or in a Versus Time mode. The display mode can be selected from the following window.




Fig. 12-23

Choosing the display mode Tabular results in a multicolumn display showing the component names, the number of measurements, the measurement parameter file name, the time of the measurements and the results obtained. This display is similar to the window for results from impurity determinations.

Choosing the display mode Versus Time results in a display in graphical mode showing the measured data in a concentration/time graph. The data of the last measured cycle are displayed in tabular form additionally.

The length of the time axis can be adapted to the current application:




Fig. 12-24

Display of data from Residual Gas determinations in Versus Time mode:




Fig. 12-25

12.4.1.1 Adapting the pressure scale for Residual Gas measurements

The pressure scale for Residual Gas measurements is set from 1×10^{-11} to 1×10^{-4} mbar as a standard. If a different scale should be selected, proceed as follows:

- 1 Start the program MEASURE
- 2 Select SEQUENCE
- 3 Select SETUP
- 4 Define the Maximum Display Value and the Number of Decades
- 5 Close the SETUP window again

Index

A

Abbreviations and terms	23
Access privileges	
See Password protection	
Access Rights Control	
See Password protection	
Add a title to a data file in sequencer mode	
See Sequencer	
Adjust	175
Input of the threshold values	177
Interpretation of the comments/error messages	176
AI 421	
See Analog input	
AI/AO Manual	387
AIRDEMO, Description of example sequence	430
AIRDEMO2, Description of the sample sequence	431
Alarm	
Entering the tone frequency and duration	227
Alarm indications	
See MID/MCD measurement type	87
Alarm thresholds	
See Switching functions	
Analog input	
Activating the channels for a measurement	186
Configuration	238
Default values for characteristic curve	72
Displaying the Characteristic	241
Enter Formula for Characteristic	244
Entering the channel designations	240
Limits for entering the characteristic curve	239
Manual, QMS 200	386
Manual, QMS 422/421	386
Measurement of analog input signals	72
Parameter setup for...	242
Analog output	
Configuration	240
Displaying the Characteristic	241
Enter Formula for Characteristic	244
Entering the channel designations	240
Manual, QMS 200	386
Manual, QMS 422/421	386
Maximum output voltages	240
Parameter setup for..	242
Analyzer, Selection of	236
AO 421	240
ArcNet	211
ASCII	232

Audible alarm in sequencer operation	
See Sequencer, Error messages	
AUTO and AUTO D, Amplifier Mode	64
Axis roll mode, defining in sequencer operation	
See Sequencer	

B

Background calibration	137
Background measurement	137
Background subtraction	91
Bar width for bargraph display in sequencer operation	
See Sequencer	
Bargraph of saved MID/MCD-measure data	293
Binary to ASCII conversion	
Parameters	232

C

Calibration	
Background measurement	137
Background measurement, Purpose of the background measurement	137
Background subtraction at MID/MCD	91
Gas specific sensitivity	142
Gas specific sensitivity, Parameter setup	199
Mass scale, Comments and Error messages	135
Mass scale, Display of the Characteristic curve	304
Mass scale, Gas-specific parameter setup	199
Mass scale, Input / optimize the parameters	199
Mass scale, Performing the calibration	135
Mass spectrometer, Calibration factor library	142
Mass spectrometer, Calibration matrix	201
Mass spectrometer, Defining the external normalization	205
Mass spectrometer, Defining the internal standard	205
Mass spectrometer, Extending the calibration matrix	202
QMS offset	140
Calibration Background	
Displaying the measurement data	304
Capture Screens	58
CD-SEM, Declaration of the signal source	236
Change password	
See Password protection, Login	
Channel parameter editor	158
Channel parameters	204
Overview	159
Choose the measuring range	64
CH-TRON, Declaration of the signal source	236
CLA 421	
Declaring the	236
Operation in leak test mode	382
Clipboard, screen save	59
Color selection for screen output	245
Common SEM Voltage	223
Communication	
See Interface to the QMS	
Computer, System requirements	22

Configuration of the QMS, changing the...	234
Copy, Copying of channel parameters	170
CP 400, Declaring the...	236
Creating a measurement data file	541
Cycles, Entering the number in sequencer operation	
See Sequencer	

D

Data block, definition of...	295
Data exchange with other programs (DDE)	424
Data file Handling	
See File Manager	
DDE	424
Timeout	499
Decades	
Displayed number of...	66
Decimal positions for table output, Input of...	227
Degas	344
QMA 125, Display the emission current in Degas mode	345
QMA 125, Input of the degas time	345
QMA 125, Start/Stop	345
QMA 200/4XX	344
QMA 200/4XX, Input of degas time	344
QMA 200/4XX, Input of emission current	344
QMA 4XX, Maximum filament current	345
DI 420	
See DI/DO	
DI/DO	
Behavior in simulation mode	237
Manual, QMS 200	146, 208, 341, 386
Manual, QMS 422/421	145, 207, 341, 386
Status indication in sequencer mode	122
Disp. Decades	66
Disp. F.S.R.	66
Display	
Documentation of measure data	270
Display after Measurement	64
Display Mode	80, 85
Display range	
Configuring the...	65
Display Saved Values	
MID (Multiple Ion Detection), Cycle measurement data	286
Process (MID- / MCD measurement data)	286
Scan Analog measurement data, Magnifier function	261
Scan Analog measurement data, Manual elimination of freak values	302
Scan Analog measurement data, Menu	253
Scan Analog measurement data, Numeric display of measurement data	253
Scan Analog measurement data, Storing as Reference Scan	254
Scan Bargraph measurement data, Intensity versus Time (IVT)	274
Scan Bargraph measurement data, Magnifying a screen section	282
Scan Bargraph measurement data, Numeric display of measurement data	275
Scan Bargraph measurement data, Total Ion Current (TIC)	274
Displaying the stored measurement data	
See Display Saved Values	
DO 420	
Installing multiple boards	187

See DI/DO	
Documentation of measure data	69
Display	270
setup	69
Dynamic Measurement method	64

E

Emission	
On/Off switching	223
Emission current input/change	
See Ion Source	
Emission External	236
Emission on/off	344, 389
Entering the I/O channel identifications	
See DI/DO (Identifications of the...)	
Entering the logical names of the I/O channels	
See DI/DO (Identifications of the...)	
Entering the unit of measure for total pressure measuring	
See PE 420/PI 420, TPR/PKR	
Error History	
Define maximum number of lines	227
Displaying the error messages	304, 315
for MCD Bargraph	114
for MCD Table	110
for MCD versus Time	117
for MID Bargraph	97
for MID Table	93
for MID versus Time	100
On/off switching of the recording function	227
Printing the message file	315
Error messages	553
Error, QMS...	
See Quadstar 32-bit Error messages	
Export Spectra Library	406
Extensions	
Installation of...	408
External emission control	236
External HV-supply (SEM)	236
External Standard	
See Calibration	134

F

F.S.R.	66
FARADAY, Declaration of the signal source...	236
Filament	
QMA 125, Selection of...	324
QMA 200/4XX, Selection of...	321
Selecting the Intern/Extern control	236
Filament protection (automatic) for the QMA 200	323
Filament protection Extern (QMA 200/QMA 4XX)	323
File Info	69
File information storage function on/off in sequencer mode	
See Sequencer	
File manager	

Copying a file	541
Creating a new file	541
Deleting a file	542
Loading a file	541
Renaming a file	542
File names	544
File names with indices (Sequencer)	416
File title	
See Documentation, Measurement data	
Files, Overview of Quadstar 32-bit files	544
FIR Filter	
Activating the filter mode (SCAN-F, PEAK-F)	160
Scan Analog under SCAN-F	62
Font selection	173
Function	
Tune Up	319

G

Gas specific sensitivity	142
Gas-specific calibration	143
General Setup, Parameter input	226
Getting started...	34

H

Help Table	100, 117
HV-external supply (SEM)	236

I

Indiced File names (Sequencer)	416
Insert screens from the clipboard	59
Installation	
Extensions	408
Integral in sequencer operation	
See Sequencer, Display Last Values	
Integral Scan	81
Integration of stored MID/MCD measurement data...	291
Intensity vs. Time representation of stored MID/MCD measurement data...	294
Interface to QMS	
Activation	221
Automatic adaptation of the baud rate	221
Communication type selection Via RS232C	210
Configuration of multiple serial interfaces under Windows	210
Configuration of the ArcNet board	212
Connection cable PC - QMS	210
Deactivation	222
For communication via RS232C	215
Installation of the ArcNet board	212
Installation of the ArcNet board, Troubleshooting	212
Parameter setup	214
Selecting the transmission mode via ArcNet	214
Setting the node address	217
Transmission types	214

Ion source parameter storage function on/off in sequencer mode	
See Sequencer	
Ion Source parameters	320
Copying/exchanging channel parameters	328, 360
Copying/exchanging IS sets	327, 359
Displaying reference measurement data	370
Displaying Sample measurement data	344, 371
Displaying Scan Analog measurement data	370
Increment/Decrement Single and Common	362
Input of the IS voltages	322
input of the protection	322
Loading / Storing defaults	329
Loading the default values	329
Optimizing the IS Voltages	329
QMA 125	323
QMA 200/4XX	320
QMA 200/4XX, Assignment of IS sets to the filaments	322
QMA 200/4XX, Displaying reference measurement data	370
QMA 200/4XX, Displaying Sample measurement data	371
QMA 200/4XX, Displaying Scan Analog measurement data	370
QMA 200/4XX, Entering the emission current	322
QMA 200/4XX, Entering the IS voltages	322
QMA 200/4XX, Entering the protection for normal operation	323
QMA 200/4XX, Filament selection	322
QMA 4XX, Optimizing the IS voltages	329
Ion source, Selection of the type...	236
IVT (Intensity versus Time)	274

K

Key	
See Software key	

L

Leak test	379
Parameter setup	380
QMS Channel parameters	381
Lin/Log changeover in sequencer operation	
See Sequencer	
Lock-In amplifier CLA 421, Declaring the...	236
Log In	532
Log output	231

M

Main features of Quadstar 32-bit	19
Manual menu	207
Manual mode	
Digital I/O	145, 207, 341
Manual mode for hardware extensions	236
Manual mode for hardware options	386
Analog I/O	387
Digital I/O	386
Mass range, Selection of	236

Mass scale calibration	
Purpose of the calibration	134
Mass scale correction in sequencer operation	
See Sequencer or calibration, Mass scale	
Mass spectrometer	
Declare type	236
Matrix solution method	103
MCD measurement type	103
Activating the display of concentration	107, 184
Activating the matrix intensity display	184
Analysis matrix	180
Analysis matrix, Expansion of...	181
Analysis matrix, Maximum size	181
Calibration factor library	103
Calibration matrix	201
Display the measurement data of the last measured cycle	118
Dynamic measurement, Performance test	64
Input of threshold values	184
Insert Spectra from Library	185
Matrix solution method	103
Matrix solution method, Literature	103
MCD Bargraph, Activating the Error-History	114
MCD Bargraph, Parameter setup	113
MCD Table, Activating the Error-History	110
MCD Table, Activating the temporary storage	110
MCD Table, Column width for Multiple Column display	110
MCD Table, Multiple column display	110
MCD Table, Setup	109
MCD Versus Time, Activating the Error-History	117
MCD Versus Time, Setup	116
Normalization	205
Parameter input	181
Save measurement data, Input number of cycles	106
Save measurement data, Reference data	106
Storing the measurement data	105
Switching functions	186
Measure and Display	63
Measure data	
Setup documentation of...	148
Measure data documentation	69
Measure parameter input	156
Measurement data file	
MCD/MID measurement data in ASCII	552
Scan Analog measurement data to ASCII	549
Scan Bargraph measurement data in ASCII	550
Measurement parameter files	
Creating a new parameter file	157
Select directory in which to save	228
Storing parameter file after input	156
Measurement parameter storage function on/off in sequencer mode	
See Sequencer	
Measurement parameters	
Channel parameter input	183
Copying of...	170
Editor, Configuration example	166
Editor, Creating a new configuration	167, 172
Editor, Creating a new register card....	167
Editor, Deleting a register card....	168

Editor, Deleting the configuration	167, 173
Editor, Renaming a register card	168
Editor, Saving the configuration	167, 172
Editor, Selection of a configuration	167
Modifying the channel state	165
Overview	159
Print	154
Select directory	182, 193
Select the directory for storing	538
Selecting a parameter file	156
Test measurement for parameters (Measure)	174
MID measurement type	87
Display the measurement data of the last cycle	101
Dynamic measurement, Performance test	64
MID Bargraph, Activating the Error-History	97
MID Bargraph, Activating the temporary storage	97
MID Bargraph, Graphic/audible alarm indication	97
MID Bargraph, Setup	96
MID Table, Activating the Error-History	93
MID Table, Activating the temporary storage	93
MID Table, Column width for Multiple Column display	93
MID Table, Multiple column display	94
MID Table, Parameter setup	92
MID Versus Time, Activating the Error-History	100
MID Versus Time, Activating the temporary storage	100
MID Versus Time, Graphic/audible alarm indication	100
MID Versus Time, Parameter setup	99
Save measurement data, Input number of cycles	89
Save measurement data, Reference data	89
Save measurement data, Selecting the save mode	89
Switching functions	95, 111
Monochrome	
Color selection	245
Multiple Concentration Detection	
See MCD measuring type	
Multiple Ion Detection	
See MID measuring type	
Multiple Quadstar 32-bit Installation	27
Multiple scan representation, Vertical/horizontal changeover in sequencer operation	
See Sequencer	
Multiplex operation	
Data collection possibilities	212
Installation of the hub	212
Unit-specific files	213

N

Names	
Path, files	544
NIST, Export of the Spectra Library	406
Normalization mode in sequencer operation	
See Sequencer	

O

Offset alignment via Zero Mass	90
--------------------------------	----

Optimize peak shape. 333

P

Parallel multiplex operation with several mass spectrometers	219
Parameter editor	
Control elements	158
Parameter file handling	
See File Manager	
Password protection	
Access privileges, Selecting the Quadstar 32-bit application....	534
Automatic start-up of a sequence	535
Changing the User name	537
Copying the Access privileges	538
Defining the...	534
Deleting the User name	538
Deleting the User password	539
Dialog box	534
Entering the Access privileges	534
Entering the User name	537
Forgotten ACCESSC password	533
Forgotten user password	532
Login	532
Login, Into ACCESSC	533
Login, To a Quadstar 32-bit program	532
Multiuser mode, activation of	535
Number of users	531
Safeguarding the ACCESSC data...	533
Setup of predefined Access privileges	539
Storing the parameters (user-specific)...	537
Path names	544
PE 420	
Entering the unit of measure	226
Peak processor	
See Scan Bargraph measuring mode	
Peak search	178
Input of the search criteria.	179
Input of the search parameters	179
Peak Shape Optimizer	333
Start measurement	335
Performance test.	64
PI 420	
Entering the unit of measure	226
PPM 422	
CORR parameter, input of...	372
Functional diagram Ion-optical potentials	350
Ion source potentials	362
Parameter overview	350
Parameter setup of an E-Scan	354
Parameter setup of an M-Scan	355
Print	
General Settings	154
Ion Source Settings.	325
Measurement parameters	169
Print Screens.	59
Print sequence log.	231, 422
Programming the sequencer	

See Sequencer	
Protection input/change	322

Q

QMS	
Hardware configuration, Controlling the filament and the SEM supply	236
Hardware configuration, Declaration of the installed analyzer family	236
Hardware configuration, Declaration of the ion source type	236
Hardware configuration, Declaration of the mass range	236
Hardware configuration, Declaration of the signal source type	236
Manual selection of a control unit in the ArcNet network	218
Setting the node addresses of several control units in the ArcNet network	218
Type declaration	236
QMS 200	
Analog I/O (Manual)	387
Digital I/O (Manual)	386
Temperature monitoring	337
QMS 422/421	
Analog I/O (Manual)	387
Digital I/O (Manual)	386
QMS error	
See Quadstar 32-bit Error messages	
QMS offset	140
QMS warnings	
See Quadstar 32-bit Error messages	
QMU 112	
Channel display in MID Bargraph	98
Declaring the...	236
Display of the channel for MID-Table	302
Displaying the channel in MID-Versus Time	287
Entering the common SEM voltage	225
MCD measurement type	108
Operation in leak test mode	382
Switching SEM On/Off	225
Quadstar 32-bit	
Computer system requirements	22
Expansion stages	20
Installation	26
Installation in a directory, non standard	27
Installation on a hard disk	27
Main features	19
Multiple Installation	27
Process Control	21
Product description	18
Quantitative Analysis	20
Scope of delivery	21
Standard	19
Update	27
Quadstar 32-bit error messages	
List of error messages	553

R

Reference data	
Load MCD measurement data	106

Load MID measurement data	89
Save MCD reference data	106
Save MID reference data	89
Storing of Scan Analog measurement data	254
Relations between Ion energy and Ve, CORR and SHIFT	363
Residual gas analysis program	575
RF generator of the QMS 200	
Tuning	337
RGA Program	575
RS232C	
See Interface to the QMS	

S

Sample files, types of	22
Save difference	81, 86
Save Screens	59
Save, temporary for MID/MCD	91
Scan Analog measurement type	77
Integral Scan	81
Save difference	81
Scan Analog Setup	80
Setup, Display Mode	80
Setup, Scan Arrangement	80
Temporary storage of measurement data	81
Scan Arrangement	80, 85
Scan Bargraph	
PEAK-F	86
PEAK-L	86
STAIR	86
Scan Bargraph measurement type	
Display Mode	85
Parameter setup	85
Save difference	86
Scan Arrangement	85
Screens, insert from the clipboard	59
Screenshots	59
SEM	
Common SEM, Entering the voltage	223
Common SEM, Influence on an analysis	138
On/Off switching	223
Selecting the Intern/Extern control	236
Sequence	
Example, AIRDEMO	430
Example, AIRDEMO2	431
See Sequencer	
Sequencer	
Add a title to a data file	494
Axis roll mode (RollAxis parameter)	498
Bar width for MID/MCD (BarWidth parameter)	496
Changeover (ScanArran parameter)	498
Color mode single/multiple (ColorMode parameter)	496
Cycles, Entering the number (Cycles parameter)	497
DI/DO status, Displaying the	123
Display Last in MID/MCD Bargraph Mode (DispLast parameter)	497
Display Last Values, Activating the Table display (DispOpt parameter)	497
Display Last Values, Integral (DispOpt parameter)	497

Display Range for Int/Float values.....	93
Display relative values (RelVal parameter).....	497
Error messages, Audible alarm (AlarmSound parameter).....	496
File information storage function, Switching on/off.....	493
File names with indices.....	416
Ion source parameter storage function on/off.....	493
Lin/Log changeover (Scale parameter).....	498
Line type solid/dot (LineType parameter).....	497
Marker display (Marker parameter).....	497
Mass scale correction (MSC parameter).....	497
Measurement parameter storage function on/off.....	493
Measuring with.....	120
Multiple scan representation, Vertical/horizontal.....	498
Normalization mode (NormMode parameter).....	497
Parameter file, displaying/editing.....	198
Print a sequence.....	193
Program a sequence.....	190
Programming, Input.....	190
Reset Additional Info names.....	494
Reset Additional Info values.....	494
Scan Arrangement.....	127
Set Additional Info name.....	494
Set Additional Info value.....	494
Starting a sequence.....	122
Status parameter storage function on/off.....	493
String values.....	417
Threshold monitoring, Alarm message (AlarmError parameter).....	496
Threshold monitoring, Display of... (AlarmDisp parameter).....	496
Time scale extent (Parameter TimeScale).....	498
Translate all available sequences.....	193
Variable types.....	197, 411
X-axis lettering.....	498
X-Raster in measurement screen.....	498
Y-Raster in measurement screen.....	498
Zero equalization (Zero Equalize).....	498
Zero gas data storage function on/off.....	493
Sequencer (Display of measurement data).....	123
Gas-specific calibration.....	131
Mass scale calibration.....	132
MCD Bargraph.....	129
MCD Table.....	129
MCD Versus Time.....	130
MID Bargraph.....	126
MID Table.....	125
MID Versus Time.....	126
Scan Analog.....	127
Scan Bargraph.....	128
SetPar() instruction.....	123
Zero gas measurement.....	131
Sequencer (Editing of sequences).....	
Compilation of sequences.....	192
Font selection.....	192
Input of array limits.....	197
Memory requirements for variable arrays.....	198
Variable types.....	197
Sequencer (Syntax).....	
*/ (End of comment).....	428
/* (Begin of comment).....	428

Auxiliary instructions, Blank line	427
Auxiliary instructions, Block begin	426
Auxiliary instructions, Block end	426
Auxiliary instructions, Blocks, Definition of	426
Control structures, If-Else instructions	419
Control structures, Loop instructions	420
Control structures, Quadstar 32-bit dialog box	431
Control structures, While instructions	421
Default values	415
Else	427
Single line comment	428
Syntax of a sequence	418
Syntax of a sequencer instruction	414
Variables, Global float array gf[.]	411
Variables, Global float data field gfa[.][]	411
Variables, Global integer array gi[.]	411
Variables, Global string array gs[.]	411
Variables, Integer array i[.]	411
Sequencer instructions	
Begin	426
Calculate()	439
DDEClose()	441
DDEExecute()	441
DDEInit()	442
DDEPoke()	442
DDERequest()	443
Delay()	444
Dialog()	445
DispEventWindow()	448
DispWindow()	449
End	426
Error()	450
Execute()	451
Exit()	451
ExtractString()	452
FileCopy()	453
FileDelete()	453
FileMove()	454
FileRename()	454
FormatDate()	455
FormatTime()	456
GetDate()	457
GetErrorCode()	459
GetSysInfo()	462
GetTime()	457
GSC()	463
IfBit()	464
IfDI()	465
IfDO()	466
IfMeaVal()	467
IfString()	469
IfVar()	469
IniFileRead()	470
IniFileWrite()	471
InitDisp()	473
Loop()	473
LoopDialog()	448
LoopTable()	474

MCD()	475
Message()	477
MID()	478
MSC()	480
OffsetCal()	481
PPMGetCorr()	521
PPMMID()	521
PPMNeutrals()	523
PPMScanEnergy()	523
PPMScanMassAnalog()	525
PPMSetEnergy()	526
PPMSetEPolarity()	527
PPMSetIS()	527
PPMSetPSM()	530
Process()	482
QMSConnect()	485
QMSDisconnect()	485
ScanAnalog()	485
ScanBar()	487
SeqExit()	488
Sequence()	489
SetAO()	490
SetDefault()	491
SetDO()	491
SetFileInfo()	493
SetISPar()	494
SetPar()	496
SetParFile()	499
SetQMS()	500
SetScale()	502
SetString()	503
SetUnit()	504
SetVar()	505
Sound()	506
TextClear()	507
TextClose()	507
TextInit()	507
TextSet()	508
WhileDI()	510
WhileDO()	512
WhileMeaVal()	513
WhileString()	515
WhileVar()	517
ZeroGas()	519
Set Unit, choose a QMS Unit	218
Setup of Quadstar 32-bit	26
Show/Hide supplementary graphic elements..	253
Simulation	
Description	228
Display in the status bar	228
External mode	228
Internal mode	228
on/off	228
Software key, display	153
Software options, display	153
Software update	27
Spectra library	393
Copying a spectrum	401

Displaying a stored spectrum	395
Editing a spectrum	404
Entering a new spectrum	403
Export to Text File	406
Go to a spectrum	403
Load a comparison spectrum	396
Normalizing a spectrum	400
Numeric display of spectra	396
Open Library	394
Overlapping of spectra	398
Printing the spectra	405
Removing a spectrum	402
Searching a spectrum by formula	402
Sorting the library	403
Storing the spectra	398
Subtraction of spectra	399
Standard deviation of measurement data	291
Starting the Quadstar 32-bit programs	57
State, input of channel state	165
Statistical information	
Calculation of the data	290
on stored measurement data	289
Status bar	
Description of	65
Output of messages	65
Symbol for simulation	228
Status parameter storage function on/off in sequencer mode	
See Sequencer	
Storing the spectra in the library	398
String values (Sequencer)	417
Subtract reference	91
Subtraction of spectra	
See Spectra library	
Switching functions	186
A and B with hysteresis	187
Behavior of	188
Graphic switching symbols	188
Logical combinations, type of	186
Outputting the state to a digital output channel	187
System-Information, displaying the	153

T

Temperature monitoring of the QMS 200	337
Temporary save for MID/MCD	91
Threshold monitoring, Alarm message	
See Sequencer	
Threshold monitoring, Displaying in sequencer	
See Sequencer	
Threshold values	
Display for Cycle measurement data	302
Thresholds	
See Switching functions	
TIC (Total Ion Current)	274
Title	
See Documentation, Measurement data	
Total pressure measurement	

See PE 420/PI 420, TPR/PKR	
Total pressure monitoring	73, 382
Activation as a signal source for a QMS channel	73
Automatic ON/OFF switching of the PE 420/PKR 250 option	73
PE 420, Entering the unit of measure...	384
PI 420, Entering the unit of measure...	384
PKR 250/251, Entering the unit of measure...	384
Setup parameters	383
TPR 250, Entering the unit of measure...	384
TPR/PKR	
Entering the unit of measure (QMS 200)	226
Trip parameter	
See Switching functions	
Troubleshooting	573
Tune Data (QMS 200), Displaying the	337
Tune Mass Scale (QMS 200)	338
Tune Up	
Function	372
Measurement	335
Tuning the QMS 200 RF generator	337

U

Update to new version	
See Quadstar 32-bit, Update	

W

Warnings	
See Quadstar 32-bit Error messages	

X

X-axis lettering in sequencer operation	498
X-Raster in measurement screen in sequencer operation	498
X-scale in the measurement screen in sequencer operation	498

Y

Y-axis lettering	68
Y-Raster in the measurement screen in sequencer operation	498

Z

Zero equalization (Zero Equalize) in sequencer operation	498
Zero Gas	137
Zero gas data storage function on/off in sequencer mode	493
Zero Mass	90