



Quadrupole Mass Spectrometer System QMG 422



Validity

This document applies to QMG 422 systems with QMS 422 and QMI 422 control units equipped with the modules listed on →  65, some of which come with their own operating instructions.

It is valid for firmware numbers:

DSP / PRG No. **BG 509 732 -..**

QMS / PRG No. **BG 509 733 -..**

CS 422 / PRG No. **BG 509 734 -..**

The above numbers can be read out with *config-TEST* or the Quadstar 422 software. The letter (A...Z) at the end of the number represents the modification index which indicates the firmware level. This operating manual remains valid as long as only the index changes. In most cases the function is enhanced but also additional functions may be included that are not described in this edition.

We reserve the right to make engineering changes without notice.

Contents

Validity	2
1 Safety	5
1.1 Symbols used	5
1.2 Intended use	6
1.3 Safety information	6
1.4 Liability and warranty	6
1.5 Courses	6
2 Overview	7
2.1 QMS 422 control unit	8
2.2 QMI 422 control unit	8
2.3 EP 422 Electrometer preamplifier	8
2.4 Family 400 components	9
2.5 Family 125 components	9
2.6 Options	10
3 Technical data	11
3.1 General	11
3.2 QMS 422 control unit	11
3.3 CS 422 operator console	11
3.4 QMI 422 control unit	11
3.5 QC 422 quadruple controller	12
3.6 IS 420 Ion source supply	16
3.7 EP 422 Electrometer	17
3.8 CP 400 Ion counter preamplifier	17
3.9 HV 420 High voltage supply	18
3.10 HV 421 High voltage supply	18
3.11 DI 420 Digital Input	19
3.12 DO 420A Digital Output	19
3.13 AI 421 Analog Input	20
3.14 OH 421 Optical Hub	20
4 Installation	21
4.1 QMS/QMI 422	21
4.2 Overall system	21
4.3 EP 422	22
4.4 CP 400	22
4.5 Cabling with QMA 400	23
4.6 Cabling with QMA 125	24
4.7 RS-232-C interface	25
4.8 LAN interface	26
4.9 Installing/removing options	28
5 Description	31
5.1 Operator console	31
5.2 Functions	33
5.3 Parameter list	35
6 Operation	46
6.1 Initial start up	46
6.2 Filament protection	47
6.3 Degas	47
6.4 Offset correction	47
6.5 Ion counter	48
6.6 Extern input	48
6.7 Simulation	48
6.8 Measurement cycle	49
6.9 Electrometer modes	50
6.10 Mass scan modes	50
6.11 Integral spectrum	52
6.12 Adjust	52
6.13 Analog outputs	55
6.14 Switching functions	58

7 Troubleshooting	60
7.1 General	60
7.2 Warnings	60
7.3 Error messages	60
7.4 Measurement signal problems	62
7.5 General problems	63
7.6 Service interventions	64
8 Maintenance	64
9 Decommissioning	64
10 Spare parts and accessories	65
 Appendix	 66
A: Default parameter values	66
B: Literature	68
C: Index	69
 Declaration of Contamination	 70
Declaration of Conformity	71

1 Safety

1.1 Symbols used



DANGER

Information on preventing any kind of physical injury.



WARNING

Information on preventing extensive equipment and environmental damage.



Caution

Information on correct handling or use. Disregard can lead to malfunctions or minor equipment damage.



Skilled personnel

Instructions marked with this symbol may only be carried out by persons who have suitable technical training and the necessary experience to do it safely.

italic-ITALIC:ITALIC *Function-PARAMETERNAME:PARAMETERVALUE*
Example: *mass-FIRST:12* (starting mass 12 of the mass scan)

1.2 Intended use

The QMG 422 is a mass spectrometer designed for gas analysis in the high vacuum range. It may be used only for this purpose. The instructions in this user's guide and in those of the accessories must be conscientiously followed.



DANGER

The QMG 422 is not intended to produce measurement results on which the safety of persons or large assets depend. For such applications the safety must be ensured by additional measures.

1.3 Safety information

Adhere to the applicable regulations and take the necessary precautions for the process media used.

When returning products that have been exposed to the vacuum for maintenance or repair, enclose a declaration of contamination.

Adhere to the forwarding regulations and prescriptions of the countries and forwarding agencies concerned.

Before handling any used instruments or components, find out whether they are contaminated. Adhere to the relevant regulations and take the necessary precautions when handling contaminated parts.

Pass on the safety information to other users.

1.4 Liability and warranty

Pfeiffer Vacuum assumes no liability and the warranty becomes null and void if the custodian or third parties

- disregard the information in this document
- use the product in a non-conforming manner
- make any kind of changes (modifications, alterations etc.) to the product
- use the product with accessories not listed in the corresponding product documentation

The custodian assumes the responsibility in conjunction with the process media used.

1.5 Courses



Courses

Pfeiffer Vacuum offers application, operating and maintenance courses for the best use of this product. Please inquire with your local Pfeiffer Vacuum partner.


2 Overview

Family 400

The QMG 422 comprises two equipment families:

In the family 400 the analyzers QMA 400 and QMA 430 (8 mm rod diameter) and QMA 410 (16 mm rod diameter) and the HF generators QMH 400/410 with mass ranges of 128 to 2048 amu are used.

400 This symbol refers to information that is applicable only to Series 400 components.



For an explanation of the abbreviations see Sections 2.3 to 2.7


For ion counting the CP 400 ion counter preamplifier rather than the EP422/2 electrometer is used.

With a QMI 422 rather than QMS 422 no system bus, IS 420 ion source supply, HV 420/421 high voltage supply and options are available.

Family 125


The family 125 uses the QME 125 mass filter electronics with the mass ranges 100 and 200 amu and QMA 125 analyzers with 6 mm rod diameter.

125 This symbol refers to information that is applicable only to Series 125 equipment.



With QMS 422 rather than QMI 422 the system bus and consequently the options are also available for this family.

2.1 QMS 422 control unit




The basic unit comprises the power supply, QC 422 quadrupole controller and the system bus.

For family 400 components the IS 420 and HV 420 or HV 421 are installed.

Input/output modules can be installed in either family.

2.2 QMI 422 control unit




The QMI 422 comprises the power supply and QC 422 quadrupole controller and allows computer operation of family 125 equipment.

QMH/QMA 400 can be operated, but without ion source and SEM supply.

The only option available is the AO 421 analog output or IC 421 ion counter.

2.3 EP 422 Electrometer preamplifier




The EP 422 amplifies the very small ion current signals of the analyzer to voltage levels that are suitable for further processing. It is installed directly on the analyzers in order to minimize parasitic noise.

- Compact, simple installation on QMA
- Low-noise, low-drift, little vibration sensitivity
- Fast response and quick recovery from overdriving

400 On analyzers with 90° off-axis SEM, two EP 422 can be connected. This allows simple changeover from Faraday to SEM mode.

2.4 Family 400 components

QMH 400/410




The HF generator produces the high-frequency voltage required for mass separation (→ [1]).

QMH type	Range	QMA type	Rod ø
QMH 400-1	128 amu	QMA 410	16 mm
QMH 400-5	300 amu	QMA 430	8 mm
QMH 400-5	512 amu	QMA 400	8 mm
QMH 410-1	1024 amu	QMA 400	8 mm
QMH 410-2	2048 amu	QMA 400	8 mm
QMH 410-3	340 amu	QMA 410	16 mm
QMH 402	16 amu	QMA 410	16 mm

Note: In the following QMH 400 or QMA 400 always refers to all types if nothing else is specified.

QMA 400/410/430




The 400 analyzer comprises the ion source, mass filter, ion collector and housing (→ [2]).

Ion collector types:

- SEM 217: 90° off-axis with integrated Faraday
- SEM 218: 90° off-axis with integrated Faraday and separate conversion dynode CD

On request: Faraday collector only


IS 420



The ion source supply is installed in the QMS 422, and supplies the ion source with the necessary operating voltages.

- Programmable potentials, short-circuit-proof.
- Polarity reversible for positive and negative ions
- Normal mode/degas mode
- Suitable for all ion source type of the QMA 400

HV 420 / HV 421




The high-voltage supply is installed in the QMS 422 and supplies the SEM with the necessary high voltage.

HV 420: For positive ions with SEM 217

HV 421: For positive and negative ions with SEM 217 or for positive ions with the separate conversion dynode of the SEM 218.

2.5 Family 125 components

QME 125



The mass filter electronics QME 125 → [3] comprises the ion source supply, high frequency generator and SEM high voltage supply.

QME 125-1: Mass range 100

QME 125-2: Mass range 200

A special version with 6 m cable between QME 125-1 and QMA 125 is available.

QMA 125




The QMA 125 analyzer → [4] comprises the ion source, mass filter, ion collector and housing.

3 Ion collector types are available:

Faraday, Faraday/Channeltron, 90° SEM


2.6 Options

CS 422




The operator console of the CS 422 is installed in the QMS 422 and allows manual operation.

AO 421 / IC 421



The analog output AO 421 and the ion counter IC 421 are mounted on one board which is directly installed in the QC 422 (without system bus).


CP 400



The ion counter preamplifier comprises the pulse coupling, amplifier and pulse height discriminator with adjustable threshold.

It is installed directly on the QMA with 90° off-axis SEM and is connected to the IC 421.

Input/Output




AI 421: 16 channel analog input module

DI 420: 32 bit digital input module

DO 420A: 32 bit digital output module


Vacuum measurement



PI 420: Dual Pirani module for coarse and fine vacuum (→ [8]).


PE 420: Penning module for high vacuum (→ [9]).

OPA 200



OPA 200 network controller board for the ARCNET® local area network (LAN). It is installed in the PC (→ [6]).


OH 421



5-Port optical hub (star distributor) for the ARCNET® local area network (LAN)

Up to 255 nodes can be cascaded

OHA 200



5 or 10 port optical hub (star distributor) for the ARCNET® local area network (LAN) (→ [7]).

3 Technical data

3.1 General

This information applies to all components unless specified otherwise.

Ambient conditions

Temperature	Storage: $-40\text{ }^{\circ}\text{C}$... $+65\text{ }^{\circ}\text{C}$ / Operation: $+5\text{ }^{\circ}\text{C}$... $+40\text{ }^{\circ}\text{C}$
Relative humidity	$\leq 80\%$ up to $+31\text{ }^{\circ}\text{C}$, decreasing linearly to 50% at $+40\text{ }^{\circ}\text{C}$
Use	indoors, altitude up to 2000 m
Type of protection	IP 30: protection $>2.5\text{ mm}$ against particles no protection against water


Standards

Safety	EN 61010-1: Protection class 1, pollution degree 2, overvoltage category II
EMC	EN 50081-2, EN 50082-2

3.2 QMS 422 control unit

Power: 90 ... 265 VAC, 47 ... 63 Hz, $300\text{ W}_{\text{max}}$

Dimensions:



Weight: 9.6 kg with QC 422 (without additional modules)

Number of slots: Total 17, used by QC 422: 3

3.3 CS 422 operator console

Matching to QMS 422


Backlit LCD display, 4 lines of 40 characters each, 5 status LEDs, membrane keyboard

Weight: 0.75 kg

3.4 QMI 422 control unit

Power: 90 ... 265 VAC, 47 ... 63 Hz, $200\text{ W}_{\text{max}}$

Dimensions:



Weight: 6.5 kg with QC 422

3.5 QC 422 quadruple controller

Slots	3 (with and without AO 421 or IC 421)
Number per unit	1
Weight	0.67 kg without / 0.9 kg with AO/IC421
Number of measurement channels	64
Operating modes	<i>MONO</i> / <i>MULTI</i> channel
Measurement cycles	1 ... 10'000, or <i>REPEAT</i>
Channel switching time	1.5 ... 3 ms (with min. <i>PAUSE</i>)

Mass scan modes

mass-MODE	Purpose
<i>SCAN-N</i>	Analog scan normal
<i>SCAN-F</i>	Analog scan with FIR filter for measured value
<i>STAIR</i>	Scan Bargraph
<i>SAMPLE</i>	Single mass and MID (Multiple Ion Detection)
<i>PEAK-L</i>	Peak search with level criterion
<i>PEAK-F</i>	Peak search with FIR FILTER

Mass scale resolution

SPEED	STEPS per mass 1)	
	<i>FIX</i> -Range	<i>AUTO</i> -Range
0.5... 1 ms/u	16/u	---
2 ...5 ms/u	32/u	---
10...20 ms/u	64/u ²⁾	16/u
50...100 ms/u	64/u ²⁾	32/u
0.2...60s/u	64/u ²⁾	64/u ²⁾

¹⁾ See STEPS → 43 ²⁾ 32 at mass range 2048

Measurement speeds

mass-MODE	EP 422 or ext. input		Ion counter
	<i>FIX</i> -Range	<i>AUTO</i> -Range	<i>AUTO</i> -Range
<i>SAMPLE</i>	0.5 ms ... 60 s	0.5 ms ... 60 s	1 ms ... 60 s
<i>STAIR</i>	0.5 ms/u ... 60 s/u	2 ms/u ... 60 s/u	2 ms/u ... 60 s/u
<i>SCAN</i>	0.5 ms/u ... 60 s/u	10 ms/u ... 60 s/u	20 ms/u ... 60 s/u
<i>PEAK</i>	0.5 ms/u ... 60 s/u	10 ms/u ... 60 s/u	20 ms/u ... 60 s/u

Detectors

detect-TYPE	
<i>FARAD</i>	Faraday collector, EP 422
<i>SEM</i>	SEM (type configurable), EP 422
<i>ION-CNT</i>	Ion counter, CP 400 / IC 421
<i>EXTERN</i>	External analog input of the QC 422
<i>PIRANI</i>	Pirani
<i>PENNING</i>	Cold cathode
<i>A-INPUT</i>	Analog signal via AI 421 module

Measurement ranges and resolution


Detector type	Meas. ranges	Modes	Resolution
<i>FARAD, SEM</i>	$10^{-12} \dots 10^{-5}$ A fsd	<i>FIX-</i> and <i>AUTO-</i> Range	16 bit * (per range)
<i>EXTERN</i>	GAIN 1: ± 10.240 V GAIN 10: ± 1.024 V	<i>FIX</i> -Range	16 bit *
<i>ION-CNT</i>	10-2...108 cps, meaning full use up to 106...107 cps	<i>AUTO</i> -Range	in mass-MODE: <i>SAMPLE</i> : 1/DWELL <i>STAIR</i> : 2u/SPEED <i>SCAN</i> : STEPS/SPEED

*) Further increased by averaging

Analog filter

Type	Two-stage lowpass, effective for electrometer and external input
Filter time constant	automatic or selectable in eight steps: 18, 85, 400 μ s / 1.7, 8, 40, 180, 800 ms

Filter step response



τ_{63} : Filter time constant

Settling time to $\pm 1\%$:

$$t_s \approx 4 \times \tau_{63}$$

Digital filter

NORMAL (N)	Low pass (average value)
FIR (F)	Finite Impulse Response
AVERAGE	Average formed across several measurement cycles



Ion sources

400	Also refer to → 16
Types	Axial, cross beam, grid, sputter process monitor, Spec+, Spec-
Parameter sets	4 per ion source
Potentials	V1 ... V9
125	Also refer to → [3]
Types	Axial, cross beam, grid, sputter process monitor
Emission	Standard: 0.07 ... 2 mA; Degas: 0.7 ... 20 mA / 500 V


Switching functions

<i>trip-TYPE: ABS</i>	2 absolute switching functions per channel
<i>trip-TYPE: HYST</i>	1 hysteresis switching function per channel
Reaction time with DO 420A	<1 ms after measurement is completed

RS-232-C interface


Detailed description	Also refer to →  [5]
Measured data buffer	256 kB
Protocol	ASCII or binary protocol (according to SECS-1 standard) 1 start, 8 data, 1 stop bit, no parity
Connector	9 pin D-sub, see (→  25).
Baud rate	300 ^{*)} , 1200, 2400, 4800, 9600, 19200 baud
Cable length	≤15 m, shielded for baud rate 19200 baud >15 m, shielded at reduced baud rate
*) Only in conjunction with ASCII protocol	

LAN interface

Type	ARCNET® with fiber optics
Connection	JIS F07 / TOSLINK
Type of fiber	PCF 200/230 or 200/300 or APF 980/1000 µm
Distance	→  26
Baud rate	2.5 Mbit/s
Wavelength	800 nm
Length of fiber-optic conductor	0 ... 1000 m
Transmission distance	3000 m _{max} (cascaded)

Connections on QC 422

Connector ctrl



Pin

5	ELM OUT+
9	ELM OUT-
4	EXT IN-
8	EXT IN+
3	RUN IN
7	--
2	EXT PROT
6	GND D
1	GND D

Analog output of the electrometer signal
0 ... ± 10.24 V, R_{out} : $2 \times 200 \Omega$, short-circuit proof
Analog input in place of electrometer signal.
 ± 10.24 V, common mode signal: ± 2 V_{max}, R_{in} = 2×50 k Ω
see detect-*TYPE:EXTERN* 45
Input for ext. measurement cycle start
TTL, int. pull-up 4.7 k Ω to +5 V.
see cycle-*TRIG* 44
Input for filament and SEM protection
TTL, int. pull-up 4.7 k Ω to +5 V.
see config-*CTRL/SEM+FIL* 43

sync

Trigger signal for recording instrument
TTL, R_{out} > 1 k Ω

elm

see above
connector *ctrl*: ELM OUT

mon (Monitor)

Analog output of the measured value
after digital processing.
0 ... ± 10.24 V
 R_{out} : $2 \times 200 \Omega$, short-circuit proof

scan

Analog mass number signal
0 ... ± 10.24 V
 R_{out} : $2 \times 200 \Omega$ short-circuit proof
Laboratory sockets, $\varnothing 2$ mm

QMH/QME

Connector for QMH 400/410
or QME 125

RS 232

Serial interface → 25

CP (Option)

Connection for CP 400

key
(not used)


AO (Option)

12 Analog outputs
0 ... ± 10.24 V, 1 mA_{max}
12-Bit, monotonic
Gain: $\pm 1.7\%$, Offset: ± 21 mV_{max}
 R_{out} : $2 \times 200 \Omega$, short-circuit proof
Rise time: 10 μ s to 99%

LAN

Fiberoptic network interface for high
transmission rates and long distances


TOSLINK TODX 296, duplex
Typ JIS F07



T Transmit
R Receive

View towards QC

D-Sub 25 pin male



View towards QC

Pin

1	CH 1+	14	CH 1-
2	CH 2+	15	CH 2-
3	CH 3+	16	CH 3-
4	CH 4+	17	CH 4-
5	CH 5+	18	CH 5-
6	CH 6+	19	CH 6-
7	CH 7+	20	CH 7-
8	CH 8+	21	CH 8-
9	CH 9+	22	CH 9-
10	CH 10+	23	CH 10-
11	CH 11+	24	CH 11-
12	CH 12+	25	CH 12-
13	open		




For EMC reasons only shielded cables may be used on D-Sub connectors. The shield must be connected to chassis ground. The opposite end may not have a different ground potential.

3.6 IS 420 Ion source supply


Slots	5
No. of IS 420 per unit	max. 1
Supply	5 V / 0.6 A; ± 24 V / 2 A (2.4 A with Degas)
Fuses	F1, F2 → 65
Filament supply	0...10 V / 5 A _{max} / 50 W _{max} with Fil1+2: 1.4...2 V on Fil2
Filament modes	1 / 2 / 1+2 (1 in operation, 2 preheated)
Protection	0 ... 5 A Resolution 10 mA
Emission normal	0 ... 2 mA Resolution 10 μ A
Emission Degas	0 ... 20 mA Resolution 0.1 mA
Signal SPEC SRC ON	23 V / 70 mA $R_i = 110 \Omega$
Ion source cable	Connector PEEK +260 °C Cable SIR -25...+180 °C
Weight	1.45 kg

With ion src-TYPE:SPEC \pm and EMISS:OFF SPEC SRC ON becomes active. A relay for changing over the ion source lines can be actuated.

Connector towards IS 420	1	QMA GND	9	Filament common
	2	SPEC SRC RET	10	V4, Field axis
	3	V6, Deflection inner	11	V0, Ref.Gnd
	4	V3, Focus	12	Screen
	5	V9, Wehnelt	13	V8, Reserve
	6	V5, Extraction	14	V1, Ionref
	7	Filament +	15	SPEC SRC ON
	8	Filament - / Cathode	16	V7, Deflection outer



	Electrode name	Ref. direction	Range [V]	Increm. [V]	Current [mA _{max}]	Degas potential to V0 [V]
V1	IONREF	V1-V0	0 ... 150	1	± 2	+550
V2	CATH	V1-V2	0 ... 125	0.5	+2	+7
V3	FOCUS	V1-V3	-30 ... +30	0.25	± 2	+550
V4	F-AXIS	V1-V4	0 ... 60	0.25	± 0.5	0
V5	EXTRACT	V1-V5	0 ... 450	2	± 0.1	0
V6	DEF-I	V1-V6	0 ... 450	2	± 0.1	0
V7	DEF-O	V1-V7	0 ... 250	1	± 0.1	0
V8	---	V1-V8	-125 ... +125	1	± 0.1	0
V9	WEHNELT	V2-V9	0 ... 60	0.25	± 0.1	+7






The external voltage source for V0 must be reliably limited to 2 mA_{max} and isolated for 750 V.

The relay connected to SPEC SRC ON and SPEC SRC RET incl. lines must be isolated for 750 V and be protected against accidental contact.

3.7 EP 422 Electrometer

Interface to	QMH 400/410, QME 125
Voltage supply	±16 VDC, ± 0.2 V / 10 mA _{max} / ripple 10 mV _{max}
Output	±10 V / 2 mA _{max}
Input impedance	100 kΩ
Input connector	Type TNC
Output connector	D-Sub 9-pin
Temperature	Operation: 0 ... 50°C, Storage: -40 ... +70°C
Weight	150 g


Dimension



Measurement range	Sensitivity	Tolerance at 25 °C	Rise time 10...90 %	Offset at 25 °C
±10 ⁻⁵ A	10 ⁻⁶ A / V	±1%	50 μs	±0.5 mV
±10 ⁻⁷ A	10 ⁻⁸ A / V	±1%	90 μs	±0.5 mV
±10 ⁻⁹ A	10 ⁻¹⁰ A / V	±2%	1.9 ms	±2 mV
±10 ⁻¹¹ A	10 ⁻¹² A / V	±2%	2.6 ms	-50 ... +150 mV
Drift	Offset doubles with each 10°C of temperature increase			
Noise	typ. 2×10 ⁻¹³ A _{pp} unfiltered			

3.8 CP 400 Ion counter preamplifier

Input	Installed directly on SEM feed-throughs of the QMA Pulse width 10 ns _{typ} / pulse height 1...5 mV Impedance 50 Ω / double pulse resolution ≤20 ns Protection against arcs in SEM with HV 420/421, QME 125			
High voltage	SHV connector HV+ and HV- HV+: 6.7 kV to GND HV-: 6 kV to GND SEM: 3.5 kV between HV+ and HV-			
Output	ECL level complementary			
Discriminator threshold	Control voltage LEVEL+ to LEVEL- 0.1 ... 1 V corresponds to pulse height 1 ... 5 mV; common mode ±0.5 V _{max}			
QC connector	Pin 1	QMA-GND	Pin 5	V+ (+12...15V / 0.12 A)
D-Sub 15 male	2	Identification	6	LEVEL-
	3	OUT-	7	LEVEL+
	4	OUT+	8	V- (-12...15 V / 0.05 A)
			9 ... 15	not connected
Weight	0.5 kg			



3.9 HV 420 High voltage supply

Slots	2	
Number per unit	max. 1	not simultaneously with HV 421
Supply	5 V / 0.4 A ; +24 V / 0.15 A ; -24 V / 0.05 A	
Fuse F1	→ 65	
SEM voltage HV-	0 ... -3500 V	Resolution 1 V, ripple <10 mV _{pp}
Admissible load	15 MΩ	Current limitation <1 mA
Internal resistance	620 kΩ	
Settling time	0.8 s to 0.1%	
HV connector	SHV	Signal ground from QMA via HV cable
HV test connection	1 V pro 1 kV	R _i = 2 kΩ
Potential isolation	0.5 V _{max}	Between chassis and QMA-GND
Weight	0.42 kg	



3.10 HV 421 High voltage supply

Slots	3	
Number per unit	max. 1	not simultaneously with HV 420
Supply	5 V / 0.2 A ; +24 V / 0.2 A ; -24 V / 0.2 A	
HV connector	SHV	
SEM voltage HV-/HV+	0 (-750 ¹)...-3500 V	Resolution 1 V, ripple <10 mV _{pp}
Admissible load	17 MΩ	Current limitation 0.8 mA
Settling time	<0.7 s	to 1%
CD voltage to GND	-6300 V	Adjustable -4.7 ... -6.4 kV
Admissible load	100 MΩ	Ripple <10 mV _{pp}
Bias voltage HV-	+3.1 / -3.1 kV	Current limitation < 0.5 mA
Test terminals	1 V per 1 kV	For SPEC- / SPEC+ adjustable 2.4...3.2 kV
Weight	1.7 kg	CD Test also for HV-
¹⁾ In SPEC+ mode		

3.11 DI 420

Digital Input



Slots	1
Number per unit	max. 2
Supply	5 V / 0.45 A
Number of inputs	32
Input signals	24 VDC $\pm 25\%$ / 10 mA, low true
Switching threshold	10.6 ... 16 V _{typ}
Insulation	30 V _{eff} / 60 VDC to GND and between input groups
Protection	+35 V / -30 V continuous; 100 V max. 1 s
Connector	2, 16 inputs each, 32-pin DIN 41612 type C/2
Weight	0.24 kg



SOURCE0..3 / 4...7 / 8...11 / 12...15 are internally connected, common plus poles of the 4 insulated input groups. SOURCE15 has no pin of its own.

3.12 DO 420A Digital Output


Slots	1
Number per unit	max. 3
Supply	5V / 0.4 A
Outputs	32, open collector Darlington drivers
Output signal	max. 28 VDC / 100 mA, $U_{sat} < 1.7 \text{ V}$ @ 100 mA, low true
Protection	Zener diode 33V/1W and parallel diode 1A
Insulation	30 V_{eff} / 60 VDC to GND and between groups
Connector	2, 16 outputs each, 32-pin DIN 41612 type C/2
Weight	0.24 kg




RET0..4 / 5...9 / 10...15 are internally connected, common negative poles of the 3 insulated outputs groups. RET15 has no pin of its own.

3.13 AI 421 Analog Input

Slots	1
Number per unit	max. 1
Supply	5V / 1 A
Inputs	16, differential, $I_{in} < \pm 150 \text{ nA @ } 70^\circ \text{C}$
Measurement range	-10.24 ... +10.235 VDC, linear -10.0 ... +10.0 V
Resolution	12 bit monotone
Accuracy	$\pm 0.1\%$ FSR
Measurement interval	40 ms for all 16 inputs
Protection	70 V to AGND, max. 8 inputs simultaneously
Insulation	30 V_{eff} / 60 VDC between AGND and GND
Connector	2, 8 inputs each, 32-pin DIN 41612 type C/2
Weight	0.3 kg



Example: Connection to input AIN3



Analog ground AGND must be connected to an admissible potential (see above "Insulation"), preferably to ground.

The cable must be shielded. Ground the shield, but not to AGND.

Twisted-pair conductors for each input provide the best signal-to-noise ratio. Open inputs produce unpredictable values.

3.14 OH 421 Optical Hub

Slots	1
Number per unit	any
Voltage supply	+5 V $\pm 2.5\%$, 200 mA typical
Connection points	5
Optical interfaces	→ 14
Weight	0.17 kg

4 Installation

4.1 QMS/QMI 422



DANGER

Before you connect the equipment make sure that the line voltage corresponds to the specifications on the nameplate.

A 3-conductor power cable with protective ground must be used.

The power outlet must have a protective ground contact.

Extensions without protective ground conductor are inadmissible.

To ensure continuity of the protective ground, always connect the power cable before all other cables. Conversely, unplug all other cables before the power cable.

Do not yet switch on the equipment!



WARNING

In rack installations the temperature inside the rack must not exceed 40°C. Ensure adequate air circulation.

The air filters inside the unit should be periodically checked and serviced (refer to 64).

In desktop installation the air should be able to enter through the lateral inlets and exit through the rear panel slots without obstruction.

4.2 Overall system

Install peripheral components such as the analyzer, QMH 400/410, QME 125 etc. in accordance with the information in the respective user's guides.

All components involved must be grounded to a single point. Utilization of a single power distributor is recommended. The only exception is the computer, but only if it is connected to the QC 422 by means of a fiber-optic link.



Skilled personnel

Make sure that the QMA, the vacuum chamber and the entire equipment is always connected to the protective ground.

Hazardous voltages up to 600 V are present on the QMA. If this unit can be touched by the user when the vacuum system is open, additional protection is required, e.g.:

Mech. protection against contact

Forced disconnection of the QMS/QMI 422 line voltage by means of a door contact

The electrode system of the QMA must not be subjected to hazardous external voltages (from contact, arcing, plasma, ion or electron beams, etc.). If such danger hazards in the vacuum system appropriate protection measures must be taken there (arrangement, shielding, grounding, etc.) that reliably preclude such influences. In addition the QMS/QMI 422 must have a permanent ground connection (no plug!). On the QMS 422 the ground terminal is located behind the power inlet, on the QMI 422 there is an M4 thread on the rear panel. Prepare this ground connection from yellow/green stranded copper wire:

2.5 mm² if mech. protected (according to DIN VDE 110 T540)

4.0 mm² if unprotected

Also refer to the standards applicable to your system.

**Skilled personnel**


When the QMA is in operation, hazardous voltages up to 600 VDC are present. Under unfavorable conditions other built-in components in the vacuum chamber (e.g. gauge heads) can be subjected to this voltage. If as a result such components become dangerous to touch (also take into consideration the lines and the connected equipment!), they must be arranged or protected in such a way that no contact, no arcs, and no charge carrier flow can occur.

4.3 EP 422

Connect the EP 422 to the corresponding connector on the QMA. Position it in such a way that it does not touch the surrounding connectors and firmly tighten the knurled nut.

Connect the control cable to the **EP** connector of the QME 125 (→  24) or the **ep1** or **ep2** connector (see p.23) of the QMH 400. Lock the connector with the slide.

For optimum signal stability the EP 422 must be protected from vibrations, temperature fluctuations, high temperature, humidity and strong magnetic alternating fields.

The Teflon cable (max. 200 °C →  65) allows remote operation if the temperature on the QMA is too high. The EP 422 must be mounted outside the hot area (M3 threads on the housing). The cable must not be subjected to vibrations. Increased noise levels must be expected.

125

The same applies for the QME 125-1 with 6 m cable length, however the maximum cable temperature is 70 °C.

4.4 CP 400

**DANGER**


The CP 400 may not be operated with a high voltage supply that can deliver hazardous voltages or currents.

Switch the unit off and detach all cables before you open the cover.

Operation of the equipment with the cover removed is not allowed.



The inside of the CP 400 may not be touched or contaminated. Finger smudges can cause noise pulses or even arcing.

Remove the SEM connector plate of the QMA (→  [2], [4]).

Remove the 6 screws A and the cover of the CP 400.

Unfasten the 3 screws C so that the full cross-section of the sockets is exposed. Caution, do not lose the screws!


Carefully slide the CP 400 over the feedthroughs; the latter must not be stressed because they can break!

Fix the CP 400 with the 4 screws B.

Tighten the 3 screws C.

Fasten the cover with all 6 screws A including washers.

For conversion to EP 422 operation perform the above steps in reverse order.



4.5 Cabling with QMA 400

Cable the system in accordance with the following tables
Short TNC and short SHV = short circuit plugs

1. Faraday cup

Configuration config-SYSTEM-DETECT:FARAD

Module	Connector	Connection	l [m]	Module	Connector	Comments
QC422	QMH	control cable QMH	3	QMH	(QC)	Optional extension 7m
IS420	QMA	Ion source	3	QMA	IS	or 10m
QMH	RF+ RF-	Radio frequency	0.7	QMA	RF A RF B	Polarity see test report
QMH	FA	Field axis	0.7	QMA	FA	
EP422	Input	Meas. signal	---	QMA	EP(FARAD)	
	---	control cable	0.8	QMH	ep1/farad	

2. SEM 217, HV 420 or HV 421

In addition to 1. Configuration config-SYSTEM-DETECT:SEM

Module	Connector	Connection	l [m]	Module	Connector	Comments
HV42x	HV-	High voltage	3	QMA	HV-	or 10 m
HV421	HV+	short SHV	---	---	---	
	CD	remains open				
---	---	short HV	---	QMA	HV+	
EP422	Input	Meas. signal	---	QMA	EP(SEM)	
	---	control cable	0.8	QMH	ep2/sem	
		short TNC	---	QMA	EP(FARAD)	if only 1 EP

3. SEM 218 (CD-SEM)

In addition to 1. Configuration config-SYSTEM-DETECT:CD-SEM

Module	Connector	Connection	l [m]	Module	Connector	Comment
HV421	CD	High voltage CD	3	QMA	CD	or 10 m
	HV-	High voltage	3	QMA	HV-	or 10 m
	HV +	short SHV	---	---	---	
EP422	Input	Meas. signal	---	QMA	EP(SEM)	
	---	control cable	0.8	QMH	ep2/sem	
		short TNC	---	QMA	EP(FARAD)	if only 1 EP

4. SEM 217, HV 420 and ion counter

In addition to 1. Configuration config-SYSTEM-DETECT:SEM

Module	Connector	Connection	l [m]	Module	Connector	Comments
HV420	HV-	High voltage	3	CP400	HV-	or 10 m
---	---	short SHV	---	CP400	HV+	
		short TNC	---	QMA	EP(FARAD)	if no EP
QC422	CP	control cable	3	CP400	QC	or 10 m

5. SEM 218, HV 421 and ion counter

In addition to 1. Configuration config-SYSTEM-DETECT:H-SEM

Module	Connector	Connection	l [m]	Module	Connector	Comments
HV421	CD	remains open				
	HV-	High voltage	3	CP	HV-	or 10 m
	HV+	High voltage	3	CP	HV+	or 10 m
		short TNC	---	QMA	EP(FARAD)	If no EP
QC422	CP	control cable	3	CP400	QC	or 10 m

4.6 Cabling with QMA 125



The **polarity** switch on the QME 125 must be in the "+" position.

1. Faraday cup

Configuration config-SYSTEM-DETECT:FARAD

Module	Connector	Connection	l [m]	Module	Connector	Comments
QC422	QME	control cable QME	3	QME	QMS	or 10/20 m
QME	QMA	Analyzer	0.2	QMA	QME	ev. 6 m
EP422	Input	Meas. signal	---	QMA	EP	ev. 6 m
	---	control cable	0.8	QME	EP	

2. Channeltron

In addition to 1. Configuration config-SYSTEM-DETECT:CH-TRON

Module	Connector	Connection	l [m]	Module	Connector	Comments
QME	HV -	High voltage	0.3	QMA	HV -	ev. 6 m

3. 90° SEM

In addition to 1. Configuration config-SYSTEM-DETECT:SEM


Module	Connector	Connection	l [m]	Module	Connector	Comments
QME	HV -	High voltage	0.3	QMA	HV -	ev. 6 m
EP422	Input	Meas. signal short TNC	---	QMA	EP(SEM) EP	in SEM mode
EP422	Input	Meas. signal short TNC	---	QMA	EP EP(SEM)	in Faraday mode
	---	control cable	0.8	QME	EP	

4. 90° SEM and ion counter

In addition to 1. Configuration config-SYSTEM-DETECT:SEM

Module	Connector	Connection	l [m]	Module	Connector	Comments
QME	HV -	High voltage	0.3	CP	HV -	ev. 6 m
		short SHV	---	CP	HV +	
QC422	CP	Ion counter	3	CP	QC	
EP422	Input	Meas. signal	---	QMA	EP	
		short TNC	---	QMA	EP	without EP422
d.o	---	control cable	0.8	QME	EP	

4.7 RS-232-C interface



4.8 LAN interface

Configuring the transmission distance

On the QC 422 check the setting of jumper X20 and correct it, if necessary.
Installation/removal of QC 422 (→ 10).

X20 determines the **Transmission power**, it is determined by the **receiver of the remote station** and the type of fiber-optic (FO) conductor..

The factory default setting **medium** provides a broad compatibility range when old and new FO modules are mixed.

Procedure:

Determine the length and type of the FO conductor (glass fiber PCF or plastic fiber APF).

Determine the FO module types to be interconnected:

Modules with serial number:W.... or higher ("W" is increased annually) are equipped with the new FO module type (applies to all component types)

In case of doubt (e.g. after a module has been replaced) open the unit and read off the FO module type, that is, No. TODX 29? on the FO connector.



This is necessary because the receivers of earlier LAN modules can be overdriven in short distances, with the new modules this is no longer the case.

Set jumper X20 on the **transmitter side** as shown in the table and the diagram:

Setting X20	New FO module TODX 296	Old FO module TODX 294	Old PC interface with SMA connectors
short	0...500 (0...1.4) m	0...150 (0...1.2) m	not allowed
medium *)	0...750 (0...2.2) m	0...400 (0...2.0) m	0...50 m
long	0...1000 (0...3.0) m	400...700 (0...3.0) m	0...300 m

These values apply to PCF glass fibers, the values in (...) for APF plastic fibers.

*) Factory setting




In case of transmission problems test the adjacent setting in order to compensate possible atypical attenuation of the FO conductor.

Cabling PC-QMG

Install the PC interface board into the PC according to its Operating manual.

Remove the protective caps from the fiber optic connectors and establish the fiber optic link.



Do not kink the fiber optic conductor. The minimum bending radius is 15 mm


OH 421

Installation/removal → 29, there is no address setting.

Settings:

Jumper	Function	Purpose
X12: IN *)	Central hub	First hub directly linked to the PC FO connectors X1 ... X5 are peer-to-peer
X12: OUT	Expansion hub	Connected to other hub Connect X1 always in the direction toward the PC! FO connectors X2 ... X5 are peer-to-peer
X6...X10	Transmitter distance setting for connectors X1 ... X5 For FO lengths refer to table on 26	


*) OH 421 No. BG 442 455-T without X12 can only be used as central hub.



Network



Unused FO connectors should always be closed off with a dummy plug to prevent disturbance by parasitic light.



QC 422, QC 421, QMS 200, OH 421, OHA 200, OH 200 and OPA 200 are compatible.

Note the technical specifications concerning FO conductor lengths and total transmission distance.

Set up all transmitters according the table on 26.

The 1st hub must be configured as the central hub, all others as expansion hubs.

4.9 Installing/removing options

Options are factory installed if they have been ordered together with the system. They can also be installed in the field at any time.



Skilled personnel

Work on open equipment may only be performed by specialists. Switch off the unit before any manipulations on the equipment. Wait 10 s and detach all cables (power cable last). For commission perform these steps in reverse order.



WARNING


Work may only be performed on ESD protected benches while observing appropriate working methods. The modules should always be stored in antistatic bags. Defects caused by the disregard of this warning will void the warranty.

CS 422

1. Detach the power plug on the QMS 422, wait 10 s
2. Remove the cover panel
3. Unfasten 6 screws of the front panel.
4. Detach the electrical connections of the front panel:
 - Ground connection
 - Flat-pin terminal on the power switch (note the pin assignment)
 - LED connection of the bus board (connector J20)
5. Establish the electrical connections on the CS 422:
 - Ground connection (sequence: head of screw, lock washer, plain washer, cable lug, plain washer)
 - Wiring to the power switch, same pin assignment as before.
 - Flat cable to bus board (connector J20)
6. Fasten the front panel
7. Mount the cover panel




Disconnect all cables from the unit, unplug the power cable



Bus modules



Disconnect all cables from the unit, unplug the power cable



To prevent damage to the connectors the module to be installed must be accurately pushed into the circuit board guides. Firmly tighten the screws. Loose screws cause malfunctions.



DANGER

As hazardous voltages are present inside the unit empty slots must be closed off with blanking plates (→ 65).

Never connect or detach cables while the equipment is switched on.

Never install or remove modules when the equipment is switched on. After power off wait 10 s before you touch or move any modules.


Always tighten the screws firmly!

Before installation check the module address according to the specifications of the individual modules.

AO 421 / IC 421


The analog output AO 421 and the ion counter IC 421 are mounted on one board.

- Disconnect the power cable and wait 10 sec.
- Install/remove the QC 422 as described above
- Remove the connector cover(s) ❶ according to the option to be installed.
- Install the AO/IC 421 as shown in the illustration. Install the hexagon pins ❷ without washers and secure them with Loctite.




IS 420

Check address setting:
Address 175 500(octal)
Jumpers A3,4,5,7,10 inserted on
IC 420



HV 420

Check address setting:
Address 175 400(octal)
Jumpers A1,2,3,4,5,6,7,10
inserted



HV 421

The address 175 404 is fixed. Jumpers: X6 all OUT, X7 IN, X8:2-4 connected.
Jumper plugs are inserted on the SHV terminals, depending on the operating mode
of the SEM (→ 23).

Setting the CD voltage

The CD or bias voltage can be set with trimmer R4. Measure the value with a DVM
on **CD-test**. Before making any adjustments switch off the unit and remove the HV
421. One counterclockwise rotation reduces the CD voltage by approx. 150 V or
the bias voltage by approx. 80 V. Reinstall, re-measure and correct the setting, if
necessary. Do not change the settings of any other trimmers!

DI 420

Up to 2 DI 420 can be installed.
Address setting: Jumpers A6, 9, 10 inserted
DI 420 #1 174 600(octal), rotary switch position 0
DI 420 #2 174 604(octal), rotary switch position 1
For connecting the inputs see technical specifications 19

DO 420A

Up to 3 DO 420A can be installed. (DO 420 is not suited)
Address setting: Jumpers A9,10 inserted
DO 420 #1 174 700(octal), rotary switch position 0
DO 420 #2 174 704(octal), rotary switch position 1
DO 420 #3 174 710(octal), rotary switch position 2
For connecting the outputs see technical specifications 19

AI 421

One AI 421 can be installed (AI 420 is not suited)
Jumper setting:
Jumpers: A6, 7, 9, 10, B, C, D, E inserted
Address 174 400(octal), rotary switch position 0
For connecting the inputs see technical specifications 20

PI 420 / PE 420

PI 420: Address: 174 000, see also → 8
PE 420: Address: 174 100, see also → 9


5 Description


5.1 Operator console

In the basic version of the QMS 422 and on the QMI 422 the front panel contains only a power switch and the corresponding LED.

The CS 422 operator console is available as an option. It can easily be retrofitted and is highly recommended for learning purposes, e.g. for software development. Even if you do not have a CS 422 you should carefully read the following information. This will enhance your understanding also in computer mode.


Contrast setting




With the   keys you can optimize the display contrast.

Parameter display

After a function has been selected, e.g. (*mass* key) the function (SETUP MASS), its parameters (here *MODE*, *SPEED...RESOL*) and parameter values and the measured value of the selected channel are displayed.




Parameter menus depend on the parameter state; in the picture below a changeover to *MODE: SAMPLE* has been made



In computer mode the parameters are displayed for approx. 30 s when a key is pressed. Input operations are disabled. Subsequently the measured value display reappears automatically.


Measured value display

After the start (**run** key) of the measurement cycle the measured value of the momentary measurement channel is displayed. In multichannel mode (*cycle-MODE: MULTI*) the progression of the channels can be seen (possibly with gaps in fast processes)



In computer mode the display is updated for monitoring purposes.

Softkeys




Choose the parameter to be entered with one of the 6 softkeys. After the key has been pressed the parameter value flashes.

If the parameter has only two values (e.g. ON/OFF or $\times 1/\times 10$) you can change it by simply pressing the soft key again.

Parameter values are entered or changed via the numeric keypad.


Number pad



- ... numeric keys
- change sign
- decimal point
- exponential entry
- cancel and return to measured value display
- Accept new value (Enter)
- Change parameters in small increments ...
- ... and in large increments. When the keys are pressed continuously the entire value range is scanned. The new value becomes effective without

Function groups

The operation is subdivided into four **function groups**:




Their function is described in the next Chapter.

Each group comprises several function keys for calling a function (e.g. *mass* in the *channels* group). Each function contains up to 6 parameters.

5.2 Functions


Parameters for operating the mass filter and for ion detection can be stored in up to 64 channels. When the measurement is performed the channels are processed sequentially and cyclically.

Channels group



Control keys for all parameters of a measurement channel

select



Choose a channel for parameter input.

All values entered under **channels** relate to this **selected** channel.

detect

Determines the signal source (detector, e.g. Faraday or SEM).

mass

Mass scan parameters such as mass number, speed, etc.

amplif

Measurement amplifier parameters, measurement ranges, Autorange/Fixrange....

aux

Enable or skip the channel during the measurement operation and copy parameter sets to different channels.


output

Parameters for analog output of measured values, e.g. linear or logarithmic.

trip

Parameters of the switching functions of each measurement channel.

General group



Keys for general settings such as configuration, initialization, maintenance and service, and processing of error messages.

di/do

Operation of DI-/DO-Bits.



config

Input and display of system and equipment configuration.

error

Error messages are displayed here in detail.

ion source group



Parameters of the selected ion source set. The set is chosen under *operation-ion src* and displayed here.

emiss

400 Emission current and filament protection


v1...v6

400 Ion source voltages V1...V6.

v7...

400 Ion source voltages V7...V9.

operation group



Control keys for operation of:

- SEM
- ion source
- measurement process

sem hv

sem hv defines the global SEM high voltage. It is effective in all channels for which no individual voltage is specified.

sem

Switch SEM high voltage on/off

The *sem* LED is on when the high voltage is switched on and flashes when it is inhibited by the EXT-PROT signal on the **ctrl** connector.

ion src

400 Ion source mode: Filament selection, degas, etc.

The ion source parameters are channel-independent and are selected as a complete set. This allows fast parameter change.

The ion source parameters of the set belonging to the active filament can be reviewed in the *ion source* group.

125 Control of ion source supply for QME 125

filam

Switch filament on/off

The *filam* LED is on when the emission is on and flashes when it is inhibited by the EXT-PROT signal on the **ctrl** connector.

cycle

Choose measurement cycle mode or offset or adjustment measurement. Input of the corresponding parameters.

run/halt

START / STOP of the measurement defined under *cycle*.


The *mono* or *multi* LED is on while the measurement cycle is running and flashes while waiting for external triggering by the RUN-IN signal on the **ctrl** connector.

The *halt* LED is on when the measurement cycle is stopped and flashes if the measurement cycle has been stopped by external triggering.

5.3 Parameter list

The parameters below are listed in alphabetic order by name.

Notation:



Additional information on the utilization, the advantages and disadvantages of individual settings can be found in Chapter 6.

ADJ-TYP (cycle)

<i>ADJ-TYP</i>	Type of search with <i>cycle-FUNCT:ADJUST</i>
<i>COARSE</i>	Coarse search → 52
<i>FINE</i>	Fine search → 54.

AI-CH (detect)

<i>AI-CH</i>	Only with <i>detect-TYPE:A-INPUT</i>
0 ... 15	AI 421 channel number to be measured

AO-CH (output)

<i>AO-CH</i>	Not with <i>detect-TYPE:PIRANI, PENNING, A-INPUT</i>
1 ... 12	Output channel of AO 421 or IC 421 for the measured value of the selected channel.
<i>NONE</i>	No analog output assigned

In *halt* condition the outputs are set to 0 V, except when they are seized by computer outputs.

AO-MODE (output)

Format selection for analog output to AO / IC 421 and **mon** (→ 55).

Electrometer in Fix-Range:

<i>AO-MODE</i>	<i>amplif-MODE:FIX</i> and <i>detect-TYPE:FARAD,SEM,EXTERN</i>
<i>LIN</i>	Linear output in selected measurement <i>RANGE</i>
<i>LOG 3D</i>	Logarithmic, 3 decades, 3 1/3 V / dec. within <i>RANGE</i>

Electrometer in Auto-Range:

<i>AO-MODE</i>	<i>amplif-MODE:AUTO, AUTO-D</i> and <i>detect-TYPE:FARAD,SEM</i>
<i>LIN</i>	Linear output within range selected with <i>O-RNG</i>
<i>LOG 3D</i>	Logarithmic, 3 decades, 3 1/3 V / dec. within range selected with <i>O-RNG</i>
<i>LOG 8D*</i>	Logarithmic, 8 decades, 1.25 V / dec. across all ranges

*) With *mass-MODE:SCAN* or *PEAK* and *SPEED* < 100 ms/u, with *STAIR* < 10ms/u a changeover to 3 decades occurs automatically

Ion counter:





<i>AO-MODE</i>	<i>LOG-DEC</i>	<i>detect-TYPE:ION-CNT</i>
<i>LIN</i>	----	Linear output within the range selected with <i>O-RANGE</i>
<i>LOG</i>	3 DEC	Logarithmic, 3 decades, 3 1/3 V / dec. within the range selected with <i>O-RANGE</i> .
<i>LOG</i>	10 DEC*	Logarithmic, 10 decades, 1 V / dec. 10 ⁻¹ ...10 ⁸ cps

*) With *mass-MODE:SCAN* or *PEAK* and *SPEED* < 50 ms/u a changeover to 3 decades occurs automatically.

AVERAGE (mass)

Moving average across measurement cycles → 51

<i>AVERAGE</i>	Only with <i>mass-MODE:SAMPLE</i> and <i>detect-TYPE:FARAD,SEM,ION-CNT,EXTERN</i>
1	No averaging across measurement cycles
2, 4, 8... ...512, 1024	Number of measurement cycles for forming average

BAUD (config-CTRL)	<i>BAUD</i>	Only with <i>config-CTRL-MODE:ASCII, BIN and MODEM</i>
	300, 1200, 2400 4800, 9600, 19200	Baud rate of the RS232 interface, can always be set on the CS 422 (300 baud only with ASCII protocol)
BEGIN (cycle)	<i>BEGIN</i>	With <i>cycle-FUNCT:CYCLE</i> and <i>cycle-MODE:MULTI</i>
	0 ... 63	Start channel of the measurement cycle with <i>cycle-MODE:MULTI</i>
BIT (di/do-DIG-OUT)	See DIG-OUT	
CALIB (amplif)	<i>CALIB</i>	Calibration factor for measured value.
	$\pm 1E^{-10} \dots \pm 9,99E^{+10}$	The raw measured value is multiplied times <i>CALIB</i> Can be used for either normalizing a measured value, e.g. from 7.7 V to 10 V (100%) or for conversion from [A] to [mbar]. In the following cases multiplication takes place only with the mantissa of <i>CALIB</i> : <ul style="list-style-type: none"> For the computer interface with <i>amplif-MODE:FIX</i> or <i>detect-TYPE:EXTERN</i> and <i>mass-MODE:SCAN</i> or <i>PEAK</i> Always for the analog signals on the AO/IC 421 and <i>mon</i>
CATH (v1...v6)	→ V2	
CLEAR (di/do-DIG-OUT)	→ DIG-OUT	
CLEAR (error)	Deletes all pending error messages	
CLEAR (cycle-FUNCT-OFFSET)	Sets all offset values to zero and consequently disables offset correction	
COPY TO CH (aux)	<i>COPY TO CH</i>	Copies parameters of the selected channel to another channel
	0 .. 63	Target channel for copying process
	SURE ?	Confirm copy function by pressing 
COPY TO ALL (aux)	<i>COPY TO ALL</i>	Copies the parameters of the selected channel to the channels <i>cycle-BEGIN...-END</i> .
	SURE ?	Confirm copy function by pressing 
COPY TO SET (ion src)	<i>COPY TO SET</i>	Copies the IQ set activated under <i>ion src-FIL1</i> or <i>FIL2</i>
	SET 0 .. SET 3	Target set for the copying process (only with <i>ion src-MODE:NORMAL</i>)
CP-LEV (amplif)	<i>CP-LEV</i>	Only with <i>detect-TYPE:ION-CNT</i>
	0.10 .. 1.00 V	Response threshold of the CP 400 →  17
CS 422 (config)	→ under TEST	
CTRL (config)	See BAUD, MODE, NODE or SEM+FIL	
CTRL (ion src)	<i>CTRL</i>	Only with <i>ion src-MODE:DEGAS</i>
	STOP	Switch Degas off
	START	Switch Degas on
	SURE ?	Confirm Degas activation with 
	RUN	Degas switched on

CYCLES (cycle)	<i>CYCLES</i>	With <i>cycle-FUNCT:CYCLE</i>
	<i>REPEAT (0)</i> <i>1 ... 10'000</i>	The measurement cycle is repeated endlessly. Number of measurement cycles to be executed
D-EMIS (ion src)	400 <i>D-EMIS</i>	With <i>ion src-MODE:DEGAS</i>
	<i>0.0 ... 20.0 mA</i>	Emission current in Degas mode.
D-PROT (ion src)	400 <i>D-PROT</i>	With <i>ion src-MODE:DEGAS</i>
	<i>0.00 ... 5.00 A</i>	Maximum filament current in Degas mode
D-TIME (ion src)	<i>D-TIME</i>	With <i>ion src-MODE:DEGAS</i>
	<i>MANUAL (0)</i> <i>1 ... 99 min</i>	Degas runs until stop command is given Degas duration. The remaining time is displayed.
DEF-I (v1...v6)	Deflection inside, see V6	
DETECT (config-SYSTEM)	<i>DETECT</i>	Specification of the existing signal source (ion collector)
	<i>FARAD</i>	Faraday collector
	<i>SEM</i>	90° SEM
	400 <i>CD-SEM</i>	90° SEM with conversion dynode
	400 <i>H-SEM</i>	High SEM, only with <i>config-SYSTEM-OPTION:CP</i>
	125 <i>CH-TRON</i>	Channeltron/Faraday combination
DIG-IN (di/do)	Status indication of the DI 420 input bits, not dependent on measurement channel. <i>NEXT</i> switches to the next 32 bits.	
DIG-OUT (di/do)	Display / manual operation of the DO 420A output bits, not dependent on measurement channel.	
	<i>BIT</i>	Choose DO bit to be operated
	<i>SET</i>	Set DO bit
	<i>CLEAR</i>	Clear DO bits
	<i>NEXT</i>	Advance to next 32 DO bits
DISP-T (config)	See <i>TEST-CS 422</i>	
DO-A, DO-B (trip)	Assignment of a switching function A or B to any bit of a DO 420A. If several switching functions are assigned to the same bit they are combined in an AND function → Section 6.14.	
	<i>DO-A, DO-B</i>	With <i>mass-MODE:SAMPLE</i> or <i>detect-TYPE:PIRANI, PENNING; A-INPUT</i>
	<i>OFF</i>	No assignment, output remains high impedance
	<i>0 ... 95</i>	Assignment of the switching functions to the DO 420A output bit
DSP (config)	The DO 420A outputs can also be operated manually or via interfaces. See <i>TEST</i>	

DWELL (mass)

Measurement time on mass number *MASS* with *mass-MODE:SAMPLE*

Electrometer or Extern:

<i>DWELL</i>	<i>detect-TYPE:FARAD, SEM or EXTERN</i>
0.5, 1, 2, 5, 10, 20, 50 ms 0.1, 0.2, 0.5, 1, 2, 5, 10, 20, 60 s	The measured value is determined by averaging across the <i>DWELL</i> time.

Ion counter:

<i>DWELL</i>	<i>detect-TYPE:ION-CNT</i>
1 ms ... 60 s	Counting rate = pulse count / <i>DWELL</i>

E-PROT (emiss)

400 <i>E-PROT</i>	Filament protection
0.00 ... 5.00 A	When the set current threshold is exceeded the filament is switched off and an error message is output.

EMI-CUR (ion src)

125 <i>EMI-CUR</i>	Switch between displaying emission current and electrometer value.
OFF	Display electrometer value
ON	Display emission current

The measurement cycle must be stopped (halt).

EMISS (emiss)

400 <i>EMISS</i>	Emission current set point
OFF	Emission switched off
0.01 ... 2.00 mA	Emission current

END (cycle)

<i>END</i>	With <i>cycle-FUNCT:CYCLE</i> and <i>cycle-MODE:MULTI</i>
0 ... 63	Ending channel of the measurement cycle with <i>cycle-MODE:MULTI</i>

EPROM-T (config)

See *TEST-CS 422*, *TEST-DSP* and *TEST-QMS*

EXTRACT (v1...v6)

See *V5*

F-AXIS (v1...v6)

See *V4*

FIL1, FIL2 (ion src)

Assignment of an ion source parameter set to filament 1 or 2. Applies also to the assignment with *SPEC±* (without filament).

<i>FIL1, FIL2</i>	With <i>ion src-MODE:NORMAL</i>
<i>SET 0 ... SET 3</i>	Assignment of the ion source parameter set

FILAM (ion src)

Filament selection for ion sources containing two filaments

400 <i>FILAM</i>	Not with <i>config-SYSTEM-IS-TYP:AXIAL</i>
1	Filament 1
2	Filament 2
1+2	Filament 1 in operation, filament 2 is pre-heated. If filament 1 is defective filament 2 is automatically activated. This results in a brief fading of the emission.
125 <i>FILAM</i>	filament select switch on QME 125 → [3] set to remote !
1	Filament 1 in normal operation / 1+2 with DEGAS
2	Filament 2

FILTER (amplif)

Time constant of the analog filter for the electrometer signal

<i>FILTER</i>	With detect-TYPE:FARAD, SEM and EXTERN			
18 μ s ... 800 ms	Manual setting for special requirements. Choose a filter value that is appropriate to the measuring speed.			
AUTO	The analog filter is automatically set as follows:			
	<i>SPEED</i> [ms/u] <i>DWELL</i> [ms]	<i>FILTER</i>	<i>SPEED</i> [s/u] <i>DWELL</i> [s]	<i>FILTER</i>
	0.5 / 1	18 μ s	0.2 / 0.5	8 ms*
	2 / 5	85 μ s	1 / 2	40 ms**
	10 / 20	400 μ s	5 / 10	180 ms
	50 / 100	1.7 ms	20 / 60	800 ms

*) Minimum value in RANGE 10^{-11} **) Minimum value in range RANGE 10^{-12} With *amplif-MODE:AUTO* or *AUTO-D* an optimum filter is automatically used.With *mass-MODE:SAMPLE* (MID mode) a fast filter results in faster settling times which means that *PAUSE* can be shortened.Signals above the maximum range of ± 10.24 V (e.g. noise) are clipped. In this case the subsequent processing (e.g. averaging) may possibly be incorrect. In critical cases the electrometer signal must be analyzed with an oscilloscope or with Quadstar 422 in *mass-MODE:SCAN-N* so that it can be optimized.With *mass-MODE:SCAN-F* the FIR filter provides for additional filtering.

FILTER (config)

See *TEST-SERVICE*

FIRST (mass)

Starting mass number of the mass scan

<i>FIRST</i>	With <i>mass-MODE: SCAN, STAIR</i> and <i>PEAK</i>
0.00 ... max. 2047.99	The maximum value depends on the mass range

The mass number is displayed as a decimal value; internally steps of $\frac{1}{64}$ u are used. At high *SPEED* the resolution decreases to $\frac{1}{32}$ or $\frac{1}{16}$ (\rightarrow 6.12).

FOCUS (v1...v6)

See V3

F.S.+ / F.S.– (config)

See *TEST-SERVICE*



FUNCT (cycle)

<i>FUNCT</i>	Measurement cycle mode, can only be changed in <i>halt</i> condition
<i>CYCLE</i>	Normal measurement operation
<i>ADJUST</i>	Adjustment to peak top with <i>mass-MODE:SAMPLE</i> \rightarrow 6.12
<i>OFFSET</i>	Offset correction of the EP 422 \rightarrow 47

GAIN (amplif)

<i>GAIN</i>	With <i>detect-TYPE:EXTERN</i>
$\times 1$, $\times 10$, $\times -1$, $\times -10$	Post-amplification factor for the Extern signal

INIT (config)

<i>RESET</i>	Load standard parameters (according to <i>config-SYSTEM</i>) \rightarrow 66
<i>FACTORY</i>	Load standard parameters set by pressing 
<i>SURE ?</i>	Confirm with  , the old parameters are irretrievably lost!

IONREF (v1...v6)

See V1

IS-TYP (config)

IS-TYP	Specify the ion source installed in the QMA.
AXIAL	Axial ion source
CB	Cross-beam ion source
GRID	Grid ion source
SPM	Sputter process monitor ion source
<input type="checkbox"/> 400 SPEC+	Special ion source positive ions
<input type="checkbox"/> 400 SPEC-	Special ion source negative ions

KEY-T (config)

See TEST-CS 422

LEVEL-A, -B (trip)

Threshold values of the switching functions

LEVEL-A LEVEL-B	TYPE	With mass-MODE: SAMPLE or detect-TYPE: PIRANI, PENNING, A-INPUT
1×10^{-24}	ABS	Threshold value of the switching function A or B
..... $9,99 \times 10^{+24}$	HYST	Upper (A) and lower (B) threshold
	OFF	Switching function off

If with TYPE:HYST : LEVEL-A < 1.1×LEVEL-B this minimum hysteresis is automatically set.

LOG-DEC (output)

See AO-MODE: ion counter → 35

MASS (cycle)

With cycle-FUNCT:ADJUST and cycle-MODE:MONO, See below MASS(mass)

MASS (mass)

In SAMPLE (MID) mode measurement takes place on this mass number during the measurement time DWELL and the average value of the measurement signal is formed.

The measurement resolution is up to 24 bits (mantissa)

MASS	With mass-MODE: SAMPLE
0.00 ... max. 2047.99	The maximum value depends on the mass range

MASS-R (config-SYSTEM)

Configuration specification of the existing measurement rage (HF generator)

MASS-R	QME-Typ	MASS-R	QMH-Typ	MASS-R	QMH-Typ
100	QME 125-1	128	QMH 400-1	1024	QMH 410-1
200	QME 125-2	512	QMH 400-5	2048	QMH 410-2
		300	d.o.+QMA 430	340	QMH 410-3

MODE (amplif)

Operating mode of the electrometer amplifier

MODE	With detect-TYPE: FARAD and SEM
AUTO	Automatic changeover across all measurement ranges, very universal
AUTO-D	Automatic changeover down to the lower search limit RANGE-L.
FIX	Manual range selection for fastest measurements

MODE (config-CTRL)


MODE	Controlling interface
CS 422	Console CS 422
ASCII	RS-232-C in ASCII format
BIN	RS-232-C in binary format
MODEM	RS-232-C with modem in binary format
LAN	Arcnet interface

Each interface can switch to itself and thereby interrupt others.

MODE (config)

See TEST-SERVICE

MODE (cycle)	<i>MODE</i>	Measurement cycle mode. The cycle is started/stopped with <i>run/halt</i>
	<i>MONO</i>	Single channel measurement in selected channel
	<i>MULTI</i>	Measurement of the channels between <i>BEGIN</i> and <i>END</i> . Channels that are in <i>aux-STATE:SKIP</i> state will be skipped.
MODE (ion src)	<i>MODE</i>	Ion source mode
	<i>NORMAL</i>	Normal operation with the parameters defined in the ion source set.
	<i>DEGAS</i>	Degas mode. The necessary parameters are entered directly.
MODE (mass)	Mass scan mode, for details refer to 50	
	<i>MODE</i>	Not with <i>detect-TYPE:PIRANI</i> , <i>PENNING</i> and <i>A-INPUT</i>
	<i>SCAN-N</i>	Normal spectrum from the start mass <i>FIRST</i> across the scan width <i>WIDTH</i> at the speed set with <i>SPEED</i> .
	<i>SCAN-F</i>	Same, with FIR filter.
	<i>STAIR</i>	Spectrum with integer mass jumps
	<i>SAMPLE</i>	Measurement on mass <i>MASS</i> with averaging across <i>DWELL</i> time.
	<i>PEAK-L</i>	Peak search (Level criterion) from <i>FIRST</i> via <i>WIDTH</i> with the speed <i>SPEED</i> . Significant data reduction because only the intensities and mass number of detected peaks are output.
	<i>PEAK-F</i>	Same, with FIR filter.
MONITOR (output)	Format of the measured value at the analog output <i>mon</i>	
	<i>MONITOR</i>	Not with <i>detect-TYPE:PIRANI</i> , <i>PENNING</i> and <i>A-INPUT</i>
	<i>LIN / LOG</i>	See <i>AO-MODE</i>
	<i>RNG-CODE</i>	Range-Code: E-12=1V.. E-5=8V (only for test purposes)
NEG (config)	See <i>TEST-SERVICE</i>	
NEXT (di/do)	Display next 32 bits with <i>DIG-IN</i> and <i>DIG-OUT</i>	
NEXT (error)	Next error message, if more than one exists.	
NODE (config-CTRL)	<i>NODE</i>	Only with <i>config-CTRL-MODE:LAN</i>
	1 ... 255	ARCNET node address
OFFSET (config)	See <i>TEST-SERVICE</i>	
OFFSET (cycle-FUNCT)	Offset correction, see <i>FUNCT(cycle)</i>	
OPTION (config-SYSTEM)	<i>OPTION</i>	Configuration input for CP 400 ion counter preamplifier
	<i>NO</i>	No CP 400
	<i>CP</i>	CP 400 exists
O-RNG (output)	Output-Range of AO 421 and <i>mon</i> , see <i>AO-MODE</i>	
	<i>O-RNG</i>	with <i>detect-TYPE:ION-CNT</i> or <i>amplif-MODE:AUTO</i> , <i>AUTO-D</i>
	<i>E-1 ... E+8</i>	In ion counter mode
	<i>E-5 ... E-12</i>	In electrometer mode
PAUSE (amplif)	<p>Measurement pause during channel change (→ 49).</p> <p>Not with <i>detect-TYPE:PIRANI</i>, <i>PENNING</i> and <i>A-INPUT</i>.</p> <p>The actual pause time is displayed if it can be calculated. <i>AUTO</i> is displayed with <i>amplif-MODE:AUTO</i> or <i>AUTO-D</i>. However, <i>P-CAL</i> is still effective.</p> <p>Press the <i>PAUSE</i> soft key to get to the <i>P-CAL</i> submenu.</p>	
P-CAL (amplif-PAUSE)	<i>P-CAL</i>	
	0.0 ... 9.9	Pause time calibration factor → 49

PE-CTRL (detect)	<table> <tr> <td>PE-CTRL</td><td>With <i>detect-TYPE:PENNING</i></td></tr> <tr> <td>OFF</td><td>Penning switched off</td></tr> <tr> <td>ON</td><td>Penning switched on (wait for ignition, pressure dependent)</td></tr> </table>	PE-CTRL	With <i>detect-TYPE:PENNING</i>	OFF	Penning switched off	ON	Penning switched on (wait for ignition, pressure dependent)						
PE-CTRL	With <i>detect-TYPE:PENNING</i>												
OFF	Penning switched off												
ON	Penning switched on (wait for ignition, pressure dependent)												
PI-CH (detect)	<table> <tr> <td>PI-CH</td><td>With <i>detect-TYPE:PIRANI</i></td></tr> <tr> <td>0 / 1</td><td>Pirani channel to be measured</td></tr> </table>	PI-CH	With <i>detect-TYPE:PIRANI</i>	0 / 1	Pirani channel to be measured								
PI-CH	With <i>detect-TYPE:PIRANI</i>												
0 / 1	Pirani channel to be measured												
POS (config)	See <i>TEST-SERVICE</i>												
PRG-NR (config)	See <i>TEST-CS 422, -DSP and -QMS</i>												
QMA (config-SYSTEM)	<table> <tr> <td>QMA</td><td>Configuration input of the QMA type based on which the unit recognizes the family</td></tr> <tr> <td>125 <input type="text" value="125"/></td><td>QMA with 6 mm rod system</td></tr> <tr> <td>400 <input type="text" value="400"/></td><td>QMA with 8 mm rod system</td></tr> <tr> <td>410 <input type="text" value="400"/></td><td>QMA with 16 mm rod system</td></tr> <tr> <td>430 <input type="text" value="400"/></td><td>QMA with 8 mm rod system (stainless steel)</td></tr> </table>	QMA	Configuration input of the QMA type based on which the unit recognizes the family	125 <input type="text" value="125"/>	QMA with 6 mm rod system	400 <input type="text" value="400"/>	QMA with 8 mm rod system	410 <input type="text" value="400"/>	QMA with 16 mm rod system	430 <input type="text" value="400"/>	QMA with 8 mm rod system (stainless steel)		
QMA	Configuration input of the QMA type based on which the unit recognizes the family												
125 <input type="text" value="125"/>	QMA with 6 mm rod system												
400 <input type="text" value="400"/>	QMA with 8 mm rod system												
410 <input type="text" value="400"/>	QMA with 16 mm rod system												
430 <input type="text" value="400"/>	QMA with 8 mm rod system (stainless steel)												
QMS (config)	See <i>TEST</i>												
QMS-HW (config)	The unit detects its modules automatically, as far as possible, and displays them.												
RAM-T (config)	See <i>TEST-CS 422, -DSP and -QMS</i>												
RANGE (amplif)	<table> <tr> <td>RANGE</td><td>With <i>amplif-MODE:FIX</i> and <i>detect-TYPE:FARAD</i> or <i>SEM</i></td></tr> <tr> <td>E-12 ... E-5</td><td>Manual electrometer range selection</td></tr> </table>	RANGE	With <i>amplif-MODE:FIX</i> and <i>detect-TYPE:FARAD</i> or <i>SEM</i>	E-12 ... E-5	Manual electrometer range selection								
RANGE	With <i>amplif-MODE:FIX</i> and <i>detect-TYPE:FARAD</i> or <i>SEM</i>												
E-12 ... E-5	Manual electrometer range selection												
RANGE-L (amplif)	<table> <tr> <td>RANGE-L</td><td>With <i>amplif-MODE:AUTO-D</i> and <i>detect-TYPE:FARAD</i> or <i>SEM</i></td></tr> <tr> <td>E-12 ... E-5</td><td>Lower search limit with <i>AUTO-D</i></td></tr> </table>	RANGE-L	With <i>amplif-MODE:AUTO-D</i> and <i>detect-TYPE:FARAD</i> or <i>SEM</i>	E-12 ... E-5	Lower search limit with <i>AUTO-D</i>								
RANGE-L	With <i>amplif-MODE:AUTO-D</i> and <i>detect-TYPE:FARAD</i> or <i>SEM</i>												
E-12 ... E-5	Lower search limit with <i>AUTO-D</i>												
RESET (config)	See <i>INIT</i>												
RESOL (mass)	<div> <input type="text" value="400"/> Setting of the mass peak separation (resolution) </div> <table> <tr> <td>RESOL</td><td>With <i>detect-TYPE:FARAD, SEM, ION-CNT, EXTERN</i></td></tr> <tr> <td>OFF (0)</td><td>Integral mass spectrum (DC OFF)</td></tr> <tr> <td>1 ... 255</td><td>Mass peak separation. The peak width is approximately proportional to the set number. Unit resolution at 20...30 (with QMH 400-1: ≈100)</td></tr> </table> <p>Often unit resolution suffices, that is, resolution of the adjacent peaks. Decreasing the mass peak separation (larger number!) causes wider peaks and higher sensitivity.</p> <table> <tr> <td>125 RESOL</td><td>Changeover Spectrum/Integral</td></tr> <tr> <td>OFF</td><td>Integral mass spectrum (DC OFF)</td></tr> <tr> <td>ON</td><td>Normal mass spectrum (DC ON)</td></tr> </table>	RESOL	With <i>detect-TYPE:FARAD, SEM, ION-CNT, EXTERN</i>	OFF (0)	Integral mass spectrum (DC OFF)	1 ... 255	Mass peak separation. The peak width is approximately proportional to the set number. Unit resolution at 20...30 (with QMH 400-1: ≈100)	125 RESOL	Changeover Spectrum/Integral	OFF	Integral mass spectrum (DC OFF)	ON	Normal mass spectrum (DC ON)
RESOL	With <i>detect-TYPE:FARAD, SEM, ION-CNT, EXTERN</i>												
OFF (0)	Integral mass spectrum (DC OFF)												
1 ... 255	Mass peak separation. The peak width is approximately proportional to the set number. Unit resolution at 20...30 (with QMH 400-1: ≈100)												
125 RESOL	Changeover Spectrum/Integral												
OFF	Integral mass spectrum (DC OFF)												
ON	Normal mass spectrum (DC ON)												
RETURN (....)	Return from a submenu to the preceding menu												
SELF (config)	See <i>TEST-SERVICE</i>												
SELF/CH (config)	See <i>TEST-SERVICE</i>												
SEM (detect)	<p>Channel related SEM high voltage</p> <table> <tr> <td>SEM</td><td>With <i>detect-TYPE:SEM, ION-CNT, EXTERN</i></td></tr> <tr> <td>SEM-HV (0)</td><td>The global value entered with <i>sem hv-SEM-VOLTAGE</i> is applicable.</td></tr> <tr> <td>1 ... 3500 V</td><td>Individual SEM high voltage for the selected measurement channel.</td></tr> </table> <p>The individual SEM high voltage leads to long settling times and makes sense only in special cases. With High SEM (<i>config-SYSTEM-DETECT:H-SEM</i>) the minimum value is 750 V (→  18).</p>	SEM	With <i>detect-TYPE:SEM, ION-CNT, EXTERN</i>	SEM-HV (0)	The global value entered with <i>sem hv-SEM-VOLTAGE</i> is applicable.	1 ... 3500 V	Individual SEM high voltage for the selected measurement channel.						
SEM	With <i>detect-TYPE:SEM, ION-CNT, EXTERN</i>												
SEM-HV (0)	The global value entered with <i>sem hv-SEM-VOLTAGE</i> is applicable.												
1 ... 3500 V	Individual SEM high voltage for the selected measurement channel.												

SEM+FIL (config-CTRL)

<i>SEM+FIL</i>	Control of filament and SEM supply
<i>INTERN</i>	Control with CS 422 or interface.
<i>EXTERN</i>	Control with EXT-PROT signal on ctrl connector: <ul style="list-style-type: none"> • Contact open: SEM+FIL switched off. • Contact closed: SEM+FIL switched on.
<i>EXT-PROT</i>	Switch-off with EXT-PROT signal on ctrl connector: <ul style="list-style-type: none"> • Contact open: Switch-off and inhibition of switch-on. • Contact closed: Enables switching on SEM+FIL via CS 422 or interface.

SEM-VOLTAGE (sem hv)

Global SEM high voltage. It is valid for all measurement channels for which no individual setting has been defined with detect-SEM:SEM-HV.

<i>SEM-VOLTAGE</i>	Not with <i>config-SYSTEM-DETECT:FARAD</i>
<i>0 ... 3500 V</i>	Global SEM high voltage

With High SEM (*config-SYSTEM-DETECT:H-SEM*) the minimum value is 750 V
→ 18

SERVICE (config)

See TEST

SET (di/do)

See DIG-OUT

SET (ion src)

See FIL1, FIL2

SIMUL (config)

<i>SIMUL</i>	Simulation spectrum for test purposes (→ 48).
<i>OFF</i>	Simulation switched off.
<i>INTERN</i>	Simulation via QC internal measurement path.
<i>EXTERN</i>	Simulation via QC external connection. Only for factory use, additional hardware required.

If no error message exists the warning SIMULATION is displayed.

SPEED (mass)

Speed for mass scan

<i>SPEED</i>	Not with <i>mass-MODE:SAMPLE</i> or <i>detect-TYPE:PIRANI, PENNING, A-INPUT</i>
<i>0.5, 1, 2, 5, 10, 20, 50 ms/u</i> <i>0.1, 0.2, 0.5, 1, 2, 5, 10, 20, 60 s/u</i>	With <i>detect-TYPE:FARAD, SEM, EXTERN</i> and <i>amplif-MODE:FIX</i>
<i>10 ms/u ... 60 s/u</i>	With <i>detect-TYPE:FARAD, SEM, EXTERN</i> and <i>amplif-MODE:AUTO, AUTO-D</i>
<i>2 ms/u ... 60 s/u</i>	With <i>detect-TYPE:ION-CNT</i> and <i>mass-MODE:STAIR</i>
<i>20 ms/u ... 60 s/u</i>	With <i>detect-TYPE:ION-CNT</i> and <i>mass-MODE:SCAN</i> and <i>PEAK</i>

STATE (aux)

<i>STATE</i>	Enable or skip a channel in multichannel mode
<i>SKIP</i>	Skip channel.
<i>ENABLE</i>	Measure channel.

STEPS

Reduces the number of measured values/u transmitted via the interface with *mass-MODE:SCAN* to ½ or ¼; this parameter can only be operated via the interface.

SYSTEM (config)

From the system configuration defined here the unit determined the possible operating modes and parameter sets.

Expansion with bus modules is detected automatically by the unit.

See sub menus *QMA, MASS-R, IS-TYP, DETECT, OPTION*.

TEST (config)

Test and alignment programs for service purposes.

The individual tests are initiated with the soft keys, endless tests are terminated with del. During a test the word BUSY is displayed, the result is shown in the status line for approximately ten seconds.

- CS 422** Test of the CS 422 console with:
- DISP-T** Endless test of the LC display. After the test has been canceled the Display-RAM test result is displayed.
- EPROM-T** After the EPROM test the result is displayed and the checksum is displayed at the soft key.
- KEY-T** Endless test of the keyboard. Consecutively press all keys to display the corresponding value.
- PRG-NR** The program number of the installed firmware (program version) is displayed.
- RAM-T** After the RAM test the result is displayed.
- DSP-...** Test of the signal processor: EPROM-T, PRG-NR, RAM-T as above
- QMS-...** Test of the system controller: EPROM-T, PRG-NR, RAM-T as above
- SERVICE** Test programs only for factory use

THRESH (mass)

With mass-MODE:PEAK-L and PEAK-F or with cycle-FUNCT:ADJUST

THRESH (...ADJUST)

Minimum intensity at which a peak is detected by the peak processor and adjust algorithm.

THRESH	<i>amplif-MODE:FIX and detect-TYPE:FARAD, SEM or EXTERN</i>
0.01, 0.03, 0,1... 30 % f.s.d	With Fixrange in % of the full scale deflection
THRESH	<i>amplif-MODE: AUTO or AUTO-D and detect-TYPE:FARAD, SEM</i>
1E-15, 1E-14.... 1E- 8	With Autorange in [A]
THRESH	<i>detect-TYPE:ION-CNT</i>
1E0, 1E1.... 1E7	In ion counting mode in counts per seconds [cps]

TRIG (cycle)

TRIG	Selection of measurement cycle start
INTERN	Start/Stop is performed via CS 422 or interface.
EXT-AUTO	Start on positive slope of the ext. start signal RUN-IN (see p 11). The cycle runs until terminated with <i>halt</i> or the specified number of measurement cycles has been attained.
EXT-NORM	Start on positive edge of RUN-IN. The cycle runs as long as RUN-IN is high, or until it is terminated with <i>halt</i> or the number of measurement cycles specified with <i>CYCLES</i> has been attained.
EXT-SNGL	Start on positive edge of RUN-IN. The unit must first be armed with <i>run</i> . The cycle runs until it is terminated with <i>halt</i> or the number of measurement cycles specified with <i>CYCLES</i> has been attained.

TYPE (detect)

TYPE	Selection of signal source, depends on the configuration
<i>FARAD</i>	Electrometer signal from Faraday collector
<i>SEM</i>	Electrometer signal with SEM
<i>ION-CNT</i>	Ion counter
<i>EXTERN</i>	External analog signal in place of EP 422 signal. Filter and processing functions of the QC 422 are used.
<i>PIRANI</i>	Total pressure measurement with Pirani module
<i>PENNING</i>	Total pressure measurement with Penning module
<i>A-INPUT</i>	Analog signals on AI 421. Filter and processing functions of the QC 422 are not used.

TYPE (ion src)

400 Changeover to special ion sources. The electrode names are replaced by "V1...V9" and all potentials are made accessible.

TYPE	With <i>ion src-MODE:NORMAL</i>
<i>xyz</i>	Normal ion source according to <i>config-SYSTEM:IS-TYP</i>
<i>SPEC+</i>	Special ion source, detection of positive ions
<i>SPEC-</i>	Special ion source, detection of negative ions. The potentials of the IS 420 and the bias voltage of the HV 421 with <i>config-SYSTEM-DETECT:H-SEM</i> are inverted.

With *SPEC±* and *EMISS* = 0 the SPEC-SRC-ON signal on the QMA connector of the IS 420 is active. In this way an external relay for changing over the ion source electrodes can be controlled.

TYPE (trip)

TYPE	Mode of switching functions (→ 6.14)
<i>OFF</i>	Switching function not active. The DO bit is available for other applications.
<i>ABS</i>	A and B are independent switching functions with one threshold value each.
<i>HYST</i>	A and B form a switching function with hysteresis. status changes when the upper or lower threshold value is exceeded.

V1... V9

400 The V... designations appears with *ion src-TYPE :SPEC±*, with standard ion sources the electrode names are displayed.
See tech. Data 16

WEHNELT (ion src)

See V9

WIDTH (mass)

Mass scan width of the measurement channel

WIDTH	Not with <i>mass-MODE:SAMPLE</i> or <i>detect-TYPE:PIRANI, PENNING, A-INPUT</i>
-2047 ... +2047	The maximum value depends on the mass range

Negative *WIDTH* results in a backward scan. In this way small peaks that are 1 mass above a very large peak can be measured more effectively.

6 Operation

The following description is applicable to units equipped with CS 422 operator console. They apply analogously also without CS 422.

6.1 Initial start up



DANGER

Before you switch on the power make sure that all components have been installed correctly (see Chapter Installation) and that the installation conforms to the technical data → 11 subseq.

The main power switch is located in the lower left-hand section of the front panel. After power on the unit performs a self-test and after a few seconds responds with a beep. Press any key to activate the measurement or parameter display.



In complete (factory aligned) systems the values that have been determined as optimal are stored in the controller. Do not modify these, go directly to 6.4.

After a change of the ion source setting, a filament change, replacement of components, etc. the following steps should be performed.

Configuration

Configure the unit as follows if the system you are putting into service has not been factory aligned.

- In the function group *general* press the *config* function.
- Choose *SYSTEM* with the corresponding softkey.
- Press *QMA* softkey and enter your *QMA* type with by pressing .
- Press *MASS-R*, enter the mass range by pressing .
- Under *DETECT* enter the ion collector type of your analyzer and under *IS-TYP* enter the existing ion source type.
- If you use the CP 400 set *OPTION* to *CP* with .
- Confirm the configuration input by pressing the *RETURN* softkey.

QMH 400/410

Please refer to the QMH 400/410 user's guide (→ [1]).



Complete factory supplied systems have already been optimally aligned. Do not change any settings without valid reason. Optimize **tune** if **best hit** does not light up.

QME 125

Please refer to the user's guide of the QME 125 → [3] and QMA 125 → [4], but do not change any factory settings.

Set the **filament select** switch to **remote** and **polarity** to **"+"** (→ [3]).

6.2 Filament protection

- 400** Optimize the filament current cut-off point for optimum protection of the filament.
- Switch off the emission: press *filam*, the *filam* lamp goes out
 - Ensure that the pressure is $\leq 10^{-4}$ mbar.
 - Switch on *filam*: *filam* light turns on.
 - If *Emission error* *) is displayed continue with f).
 - If this is not the case, switch off *filam*, reduce *emiss-E-PROT* with  by 0.1 A and turn *filam* on again; repeat until *Emission error* appears.
 - Switch off *filam*, increase *emiss-E-PROT* with  by 0.1 A and turn *filam* on again; repeat until *Emission error* no longer appears.

*) In case of ****ERROR**** press the *error* key to read out the type of error.



If the switching threshold is abnormally high a fault exists in the ion source, the vacuum system, or the unit. Investigate the cause in order to prevent destruction of the filament.

If a new filament has been installed re-optimize after a few hours.

The *D-PROT* cut-off point for DEGAS is set analogously.

6.3 Degas

Please refer to the user's guide of the QMA and the ion sources. Switch to Degas only if the conditions specified there are fulfilled.

- 400** Optimize the filament protection for Degas mode as described above, however with *ion src-D-PROT* rather than *emiss-E-PROT*.

Switch on Degas with *ion src-CTRL:START* and confirm with .

- 125** →  [4] and [3]

6.4 Offset correction

The EP 422 is an amplifier for very small currents. Its zero (offset) must be corrected occasionally. This is intentionally not performed automatically to prevent periodic dead times in data acquisition.


Perform an offset correction when the system is put into service for the first time, after the system (particularly the EP 422) has attained a stable temperature.

Subsequent repetition is advisable if deviations from the zero line occur.

Zero line shifts are often caused by ion and electron currents; the offset correction does not compensate these.

Choose *cycle-FUNCT:OFFSET* and start the offset measurement with *run*.

The following appears:



Perform the measurement for SEM and Faraday if both EP 422 are installed. If only one EP 422 exists choose DETECT according to the current operating mode.

The mass number should be selected in such a way that no ion current occurs. With *CLEAR* you can disable the offset correction; the offset values are set to zero. The values of the ranges 10^{-5} , 10^{-7} , 10^{-9} , 10^{-11} A are not displayed, they are available at the interfaces.

400 For offset measurement the SEM voltage is switched off automatically, *RESOL* = 1, *F.A.* = 0 V and *EMISS* = 10 μ A are set. If you switch *filam* off, measurement takes place without emission.

125 For offset measurement the SEM voltage and the emission are switched off automatically.

6.5 Ion counter

With the IC 421 and CP 400 ion counter the measured value is displayed as a quasi-logarithmic bar across ten decades. The range is always selected automatically; for linear analog output the display range is selected with *O-RNG*.

The discriminator threshold of the CP 400 preamplifier is set with *amplif-CP-LEV*. Recommended value: 0.1...0.3 V with *SEM-VOLTAGE* : 2500 V.

The ion polarity is selected by choosing the ion source type *ion src-TYPE:SPEC+* for positive ions or *SPEC-* for negative ions.

6.6 Extern input

With *detect-TYPE:EXTERN* you can capture analog measured values of a different unit in place of the EP 422. Connect the signal to the EXT IN pins of the **ctrl** connector on the QC 422 (see p.11)

amplif-GAIN determines the gain ($\times \pm 1/\times \pm 10$).

6.7 Simulation

For experiments and tests a simulated spectrum according to the following diagram is available with *config-SIMUL:INTERN*. Use it to familiarize yourself with the operating procedures. For this purpose you do not need a vacuum system and also no equipment such as QMH or QMA.


The spectrum is generated as an HF generator control signal **scan** and inserted into the QC 422 electrometer signal path. It can be used in all *mass modes*.



For normal measures the simulation must be switched off: *config-SIMUL:OFF*.

Only the detector types *FARAD* and *SEM* are admissible.

config-SIMUL:EXTERN is reserved for factory tests.



The peak intensities for *RANGE* 10^{-5} , 10^{-7} , 10^{-9} and 10^{-11} are identical. The intensity is amplified by a factor of 10 for the intermediate ranges. The spectrum is repeated periodically from mass 64.

6.8 Measurement cycle

With *cycle-MODE* choose single channel (*MONO*) or multichannel (*MULTI*) mode.

The measurement cycle is started and stopped with *run/halt* and its state is indicated with the *mono*, *multi* and *halt* LEDs.

The number of measurement cycles is chosen with *cycle-CYCLES*: 1...10000 or repeating (0).

You can choose external control with *cycle-TRIG*.

HALT

The measurement cycle is stopped. The unit measures in *mass-MODE:SAMPLE* (even if a different *MODE* has been entered) in the selected channel on the mass defined with *MASS* (or *FIRST*) at the speed selected with *DWELL* (or *SPEED*).

The switching functions are OFF, the signals on *elm*, *mon* and *AO* are available.

MONO

Single channel mode: enter *cycle-FUNCT:CYCLE* and *cycle-MODE:MONO*.

The unit measures in the measurement channel chosen with *select*. If the channel is changed in the *run* state the measurement is cancelled and the new channel is started.

MONO is suitable for measurement tasks in manual mode. A separate channel is programmed for each task. This means that the parameters are continually available. You can quickly change the measurement task by changing the channel.

MULTI

Multichannel mode: enter *cycle-FUNCT:CYCLE* and *cycle-MODE:MULTI*

The 64 channels can be programmed with any parameters. The channels located between *cycle-BEGIN* and *cycle-END* are processed sequentially if they are not set to *STATE:SKIP*.

It is advantageous to first optimize each channel involved in *MONO* mode.

With the *aux-COPY* function load additional channels with the parameters of the first optimized channel. Subsequently you only need to adapt a few parameters of the individual channels.

To achieve the shortest measuring time the channels with identical detector type, electrometer range and SEM voltage should directly follow each other.

The measurement cycle time *CYCLE-TIME* consists of the measurement and pause times of the involved channels. It is measured by the built-in clock.

The clock starts with *run*, stops with *halt* and is displayed by pressing *cycle*. With *cycles-CYCLES:1* (no. of cycles = 1) you measure the time for one cycle.


PAUSE

When the channel is changed in *MULTI* channel mode the data acquisition must pause until the new measured value is stable. The pause time is set automatically by the QC 422 based on *RANGE*, *FILTER*, *MASS*, *SEM* voltage and *detect-TYPE*.

The pause time is displayed under the softkey *amplif-PAUSE*. If you want to achieve shorter cycle times or greater accuracy you can optimize it with *amplif-P-CAL*. The minimum value is 1 ms with *P-CAL:0.0*.

You can reduce *P-CAL* in each channel until its measured value deviates inadmissibly. The preceding channel should not have a measured value that is almost identical, otherwise there is practically no transient response and the value of *P-CAL* would be too small. Change e.g. *MASS* of the preceding channel by ½ u to determine whether or not its measured value is without influence on the one of the selected channel.

With the oscilloscope (triggering on the falling slope of **sync**) the transient response during the pause can be observed at **elm** and (see also 55).



6.9 Electrometer modes

AUTO

With *amplif-MODE:AUTO* the electrometer range is set automatically across all decades. This results in a huge dynamic response of over 10 decades or 200 dB. Use *AUTO* whenever possible. In this way you achieve the best resolution of the measured value and no overdriving of the amplifier can occur.

AUTO-D

With *amplif-MODE:AUTO-D* (Auto down) the range is limited in the downward direction. This is usefully for noisy measurement signals and can lead to faster measurements. You can define the available dynamic response with *RANGE-L*.

FIX

With *amplif-MODE:FIX* choose the measurement range with *RANGE* manually. This allows fastest measurements with a limited dynamic response.

With *Scan-SPEED* < 10 ms/u there is only Fixrange.

For the most accurate measurements (e.g. isotope ratios) Fixrange is recommended because the mutual tolerances of the measurement ranges are eliminated or can be calibrated.


In the following diagram of the measurement signals the raw signal *elm* (→ 55) is represented always in *FIX-RANGE* because it is difficult to follow in *AUTO-RANGE*.

In all operating modes the measured values (except on *elm*) are multiplied times *CALIB* before they are output.

6.10 Mass scan modes

SCAN-N

The *mass-MODE:SCAN-N* (SCAN-Normal) mode is used for recording an analog spectrum across the range defined with *FIRST* and *WIDTH*.




The number of steps per mass depends on *SPEED* and the mass range (→ 12). With *SCAN-N* the average value of the mass signal is output with each mass step.

Example: With *SPEED* 100 ms/u and mass scale resolution $\frac{1}{64}$ u there is an integration time per step of $100 \text{ ms/u} \times \frac{1}{64} \text{ u} = 1.56 \text{ ms}$

With *SCAN-N* you obtain a direct image of the measured values captured by the measuring amplifier or the ion counter. This mode is particularly suitable for analyzing raw data, e.g. for optimizing parameter values.

SCAN-F


With *SCAN-F* the measured values are additionally subjected to an FIR filter algorithm (Finite Impulse Response).



The FIR filter largely eliminates noise and interference so that also very small peaks can be detected against the background. Statistical intensity fluctuations which on account of the 90° SEM arrangement frequently account for the major portion of the noise are particularly well suppressed by the FIR filter. For this reason you should always use *SCAN-F*, except in the few special cases where raw data are actually required.

SAMPLE

With *mass-MODE:SAMPLE* the measurement is performed on the constant mass number *MASS*. In most cases it will be set to a peak top (*ADJUST* → 6.12



AVERAGE

With *AVERAGE* > 1 a moving average (M) is formed across the number (n) of measurement cycles since *RUN*. Beginning with the first measurement cycle it supplies a value that becomes more stable with increasing number of cycles. In this way *DWELL* can be shortened without significant impairment of the filter effect (faster settling time).

The following recursive formulas apply:

- $n < AVERAGE$: $M_{new} = M_{old} + (M_{new} - M_{old}) / n$
- $n \geq AVERAGE$: $M_{new} = M_{old} + (M_{new} - M_{old}) / AVERAGE$

Time constant of the averaging: $\tau \approx AVERAGE \times \text{cycle time}$


PEAK

The *mass-MODE:PEAK* (peak processing) is an intelligent data reduction process which searches the spectrum for peaks in real time mode. Instead of 64 measured values/u only the **intensity and mass number** of detected peaks are output on the computer interface.

Mass scan is same as with *SCAN*

The marker at **mon** means that a peak of the displayed height has been detected

The value at **AO** remains until a new peak is detected



Peak Processing runs with all *SPEED* settings. The peak search extends across the range defined with *FIRST* and *WIDTH*. The peak criteria of *ADJ-TYP:COARSE* apply (→ Section 6.12).

There are two methods:


- With **PEAK-L** (Level) the peak processing algorithm is applied to the normal spectrum (SCAN-N).
- With **PEAK-F** the peak processing algorithm is applied to the measured values processed with the FIR filter. This is advantageous because parasitic signals have largely been removed from the measured values so that a very low **THRESH** can be used.

STAIR

With **mass-MODE:STAIR** integer mass jumps across the range **FIRST... WIDTH** are performed. A bargraph spectrum is created.

After each mass jump the average value across approx. half the dwell time is formed.


Example: With **SPEED** 100 ms/u the averaging time is ≈ 50 ms



The start mass of each channel must be on a peak maximum. See **ADJUST** (\rightarrow 52). If the peak maximums are not hit, large measuring errors are unavoidable. For this reason you should limit **WIDTH** per channel to approx. 10% of the mass range. In this way you can compensate deviations of the mass scale by correcting the corresponding starting mass **FIRST**.

6.11 Integral spectrum

With **mass-RESOL:OFF** an integral spectrum is created that can be used, e.g. for total pressure measurement.



6.12 Adjust

With **cycle-FUNCT:ADJUST** you can automatically optimize the mass number **MASS** to the peak maximum in **SAMPLE** (or **STAIR**) mode.

The measurement channel must be set to **aux-STATE:ENABLE**.

This possibility is advantageously used, for example, to optimize the system after turn on and particularly after several parameters have been changed.

Adjust COARSE

With **ADJ-TYP:COARSE** a range of $\pm \frac{1}{2} u$ around the mass number **MASS** is normally searched for a peak. The search range will possibly be enlarged by $\pm \frac{1}{4} u$. If possible use **amplif-MODE AUTO** for **ADJUST**, it will be easier to obtain a result.


Peak criteria:

Four criteria must be met for a peak to be detected:

- $I_{\max} > 2 I_{\min}$
- $I_{\text{end}} < 0.5 I_{\max}$
- $I_{\max} > I_{\text{tresh}}$
- $\Delta m_{50\%} \geq \frac{1}{8} u^*)$ at $\frac{1}{2} I_{\max}$
 $^*) \frac{1}{4} u$ with **mass-MODE:PEAK**

Time:

$$t_{\text{Adjust}} \approx 0.5 \dots 1.25 \text{ DWELL}$$




Start:

Stop measuring cycle: *run/halt:halt*
 Choose *mass-MODE:SAMPLE*
 Choose *cycle-FUNCT:ADJUST*
 Choose *cycle-MODE:MONO* or *MULTI*
 Choose *cycle-ADJ-TYP:COARSE*
 Start *ADJUST* measurement: *run/halt:run*

Procedure: Signal mon

With successful Adjust the mass number *MASS* of the measured channel is updated with the new value. If unsuccessful it remains unchanged.



Status message: After expiration a status message is displayed:

Symbols					Status code
ADJUST STOP STATUS CH 14: →					000001
SETUP CYCLE **					**
FUNCT	MODE	ADJ-TYP	THRESH	RANGE	MASS
ADJUST	MONO	COARSE	0.01%	1E-05	13.46

The 3 main parameters *THRESH*, *RANGE*, *MASS* can be entered here directly. You do not have to switch back to the *channels* group.

Symbols:


- OK** The adjust was successful
- Increase *MASS* slightly
- ← Decrease *MASS* slightly
- ↑ Increase Intensity or lower *THRESH*
- ↓ Decrease Intensity
- I** Peak too narrow (e.g. parasitic pulse or poor peak shape). Repeat *ADJUST*. If unsuccessful: investigate peak shape.

Status code:

	Peak width	Intensity		Mass number <i>MASS</i>		
	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
Bit=1:	too narrow	< <i>THRESH</i>	too high	too low *)	too high	too low
Symbol:	I	↑	↓	→	←	→

*) and intensity not dropped back to ½

Examples:



*) or preferably use *amplif-MODE:AUTO*

Adjust FINE


A peak maximum within the range of $\pm 1/8$ u around the mass number *MASS* is searched. Also in this case *amplif-MODE:AUTO* is recommended.

Peak criteria:

- a) $I_{\max} > I_1 + 1,5\%$
- b) $I_2 < I_{\max} - 1,5\%$
- c) $I_{\max} > \text{THRESH}$
- d) No overdriving
- e) $I_3 < I_{\min} + 1,5\%$

time:

$$t_{\text{Adjust}} \approx 16 \text{ DWELL}$$




Start: If you are not sure that a peak is located within the searched range, first perform an ADJUST COARSE. Start as described under COARSE, however with *cycle-ADJ-TYP:FINE*

Procedure:

Signal *mon*

If Adjust was successful the *MASS* of the measured channel is updated with the new value, if it was unsuccessful it remains unchanged




Status message: As in ADJUST-COARSE, however, without information on mass number and peak width.

Bits 1, 2 and 5 are always zero.

	---	Intensity		Mass number <i>MASS</i>		
	Bit 5	Bit 4	Bit 3	Bit 2	Bit 1	Bit 0
Bit=1:	---	<THRESH	too high	---	---	no Peak
Symbol:		↑	↓			→ ←

Examples of unsuccessful fine searches:



6.13 Analog outputs

On the QC 422 there are some signal connections (see pin assignment 15). Nowadays measurement data are generally acquired with computers. For investigating possible measuring problems (noise, transient response, etc.) and for special cases the analog signals can be very useful.

Connectors

elm

Analog filtered electrometer signal. It can be readily evaluated only in *amp-MODE:FIX* (Fixrange). With *Autorange* it becomes difficult to follow. **elm** is highly suitable for assessing the quality of the raw measured values. The calibration factor **CALIB** has no influence on **elm**.

mon

Monitor sequentially supplies the measurement signals of all measurement channels after they have been processed by the signal processor. It also serves as the analog output of the ion counter and can be used in linear or logarithmic format.

AO

The AO/IC 421 option can output up to twelve analog measured values in linear or logarithmic format.

You can assign one or several measurement channels (*output-AO-CH*) to each of the twelve AO channels. If several measured values are assigned to the same AO channel they will be output sequentially.

Via the computer interface also data from the PC software can be output on the AO 421.

sync

Indicates the start of the measurement and is suitable for triggering an oscilloscope.

scan

Proportional to the momentary mass number
0...10.24 V for QMH 400/410; 0 ... 10.00 V for QME 125

The behavior of the above signals in the various operating modes is described beginning on 50.

Output formats

The measured values at **mon** and **AO** have the following formats:

LIN Linear

LOG 3D Logarithmic across 3 decades

LOG 8D Logarithmic across 8 decades

LOG Logarithmic across 3 or 10 decades

The possible choices (operating mode dependent) can be found in the following tables, the scaling from the corresponding formulas and diagrams.

The formulas and graphics apply to positive measurement signals.

In the negative range the characteristics are mirror imaged at the zero according to the formulas a) ... e) and n) ... q)

Note: $\log = \log_{10}$ $U_o = U_{\text{output}}$ at **mon** or **AO**

Mass units: U : [V] I : [A] Counting rate: [cps] counts per second, s^{-1}

Electrometer operation

detect-TYPE:FARAD, SEM


mass-MODE	output-AO-MODE, MONITOR	amplif-MODE	Output format		Formulas
			Decades	V/Decade	
SCAN-N SCAN-F	LIN	FIX,AUTO AUTO-D	1	10	a)
STAIR PEAK-L	LOG 3D	FIX,AUTO ¹⁾ AUTO-D ¹⁾	3	3.333	b), c), d)
PEAK-F ADJ-COARSE ³⁾	LOG 8D	AUTO ²⁾ AUTO-D ²⁾	8	1.25	e)
SAMPLE ADJ-FINE ³⁾	LIN	FIX,AUTO AUTO-D	1	10	a)
	LOG 3D	FIX	3	3.333	b), c), d)
	LOG 8D	AUTO AUTO-D	8	1.25	e)

¹⁾ only for SCAN-SPEED 10...50 ms/u, with STAIR 2... 5 ms/u


²⁾ only for SCAN-SPEED ≥ 100ms/u, with STAIR ≥ 10 ms/u

³⁾ only at **mon** connector


LIN a) for U₀: 0...10 V: $I = 0.1 \times \text{RANGE} \times U_0$



LOG 3D b) for U₀: 1V ... 10V $I = 10^{-3} \times \text{RANGE} \times 10^{0.3 U_0}$
 c) 0.25V ... 1V $I = 1.333 \times 10^{-3} \times \text{RANGE} \times (U_0 + 0.5)$
 d) 0V ... 0.25V $I = 4 \times 10^{-3} \times \text{RANGE} \times U_0$
 Applicable with *amplif-MODE:AUTO*: RANGE = O-RNG.



LOG 8D e) for U₀: 0...10 V: $I = 10^{(0.8 U_0 - 13)}$




Ion counting operation

detect-TYPE:ION-CNT

mass-MODE	output-AO-MODE, MONITOR	output-LOG-DEC	Output format		Formulas
			Decades	V/Decade	
SCAN, STAIR	LIN	---	1	10	g)
SAMPLE	LOG	3 DEC	3	3.333	h), i), k)
PEAK ADJUST ³⁾		10 DEC (>20 ms/u)	10	1	l), m)

LIN

g) for U_o : 0 ... 10 V: Rate = $0.1 \times O-RNG \times U_o$




LOG 3 DEC

h) for U_o : 1V ... 10V Rate = $10^{-3} \times O-RNG \times 10^{0.3 U_o}$


i) 0.25V ... 1V Rate = $1.333 \times 10^{-3} \times O-RNG \times (U_o + 0.5)$

k) 0V ... 0.25V Rate = $4 \times 10^{-3} \times O-RNG \times U_o$



LOG 10 DEC

l) for U_o : 0V...10V Rate = $10^{(U_o - 2)}$




Extern input

detect-TYPE:EXTERN

mass-MODE	output-AO-MODE, MONITOR	Output format		
		Decades	V/Decade	Formulas
SCAN, STAIR	LIN	1	10	n)
SAMPLE, PEAK	LOG 3D	3	3.333	o), p), q)

LIN

n) for U_o : 0V ... 10V $U_{Extern} = (1 / GAIN) \times U_o$




LOG 3D

o) for U_o : 1V ... 10V $U_{Extern} = (0.01 / GAIN) \times 10^{0.3 U_o}$

p) 0.25V ... 1V $U_{Extern} = (0.01333 / GAIN) \times (U_o + 0.5)$


q) 0V ... 0.25V $U_{Extern} = (0.04 / GAIN) \times U_o$



Recording

y/t Recorder

For recording analog measured values (spectra or versus time) the following possibilities exist:




Paper feed should match SPEED, e.g. 1 mm/s for 1 s/u.

Sensitivity: 10 V for full scale deflection

Mirror imaged spectra can be inverted with *amplif-CALIB: -1*.

Instead of **mon** you can also use one or several **AO** outputs (→ 15).

y/t - Oscilloscope



The same applies as for the y/t recorder.

Additional tips:

Choose a **WIDTH** that is somewhat larger (e.g. 14) than the range to be represented (e.g. 10 u) and optimize for minimum flicker!

sync triggers the oscilloscope with each start of the selected channel. Trigger manually to the positive slope, for **PAUSE** to the negative slope.

Digital storage scopes are advantageous

x/y Recording

The **scan** output supplies the mass number signal (0 ... 10.24 V for the full mass range). x/y recording is only suitable for special cases.

6.14 Switching functions

With the TRIP switching functions measured values can be monitored in *mass-MODE: SAMPLE* or with *detect-TYPE: PI, PE, AI*.

With *halt* the state of the switching functions is OFF.

Each measuring channel has two switching functions: TRIP A and TRIP B.

These can be assigned to the output bits of the DO 420A modules as desired and also be interrogated via the computer interface.

If several switching functions are assigned to the same DO bit they are logically combined with an AND function. There is no warning if the DO bits are already assigned.


The DO 420A can also be controlled from the computer interface. Simultaneous assignment of switching functions is not advisable.

Vacuum relay

Choose *trip-TYPE: ABS*

Enter *trip-LEVEL-A*

ABS-TRIP A switches ON when threshold A is exceeded




With *trip-DO-A:xx* assign a DO bit if DO 420A exists

Overpressure relay

Choose *trip-TYPE: ABS*

Enter *trip-LEVEL-B*

ABS-TRIP B switches ON when threshold B is exceeded



With *trip-DO-B:xy* assign a DO bit if DO 420A exists.

Window comparator


The window comparator is only possible in conjunction with the DO 420A module

Choose *trip-TYPE* : ABS

With *trip-LEVEL-A* enter the upper threshold

With *trip-LEVEL-B* enter the lower threshold

With *trip-DO-A:xx*, *trip-DO-B:xx* assign A and B to the same DO bit.



Hysteresis function


The hysteresis prevents fluttering with unsteady signals. The minimum hysteresis is 10%

Choose *trip-TYPE* : HYST


With *trip-LEVEL-A* enter the lower threshold

With *trip-LEVEL-B* enter the upper threshold

TRIP A switches ON when the signal drops below LEVEL-A and switches OFF when it exceeds LEVEL-B



TRIP B works inversely to A



With *trip-DO-A* : xx, *trip-DO-B* : xy assign the DO bits if DO 420A exists

Window and hysteresis

This requires two measurement channels and one DO 420A.

With *select:x* choose the first measurement channel x

Choose *trip-TYPE* : HYST

With *trip-LEVEL-A* enter the lower threshold

With *trip-LEVEL-B* enter the upper threshold


With *trip-DO-A* assign a DO bit

With *aux-COPY TO CH:x+1* copy channel x to channel x+1

With *select:x+1* choose the measuring channel x+1

With *trip-LEVEL-A* enter the lower threshold

With *trip-LEVEL-B* enter the upper threshold



Many other combinations are feasible.

7 Troubleshooting

7.1 General



Skilled personnel

Work on an **open unit** (as specifically instructed in some parts of this Chapter) may only be performed by **skilled personnel**.






Corresponding warnings are not given separately on each occasion!

The relevant safety instructions given in the corresponding Chapters must be conscientiously followed.

Protection against electrostatic discharges (ESD) is absolutely essential, otherwise the Pfeiffer Vacuum warranty becomes null and void.

7.2 Warnings


Warnings related to operator action are displayed for approx. 10 sec.

No	Warning	Meaning, comments
1	** OP ERROR **	Operator error, illegal entry
2	** ↑ ONLY ↓ **	Parameter change with    
3	** > MAX **	Input value too large
4	** < MIN **	Input value too small
5	** SYNTAX ! **	Incorrect format
6	** REMOTE **	Unit set to computer operation
7	** CH SKIP **	Selected channel is in <i>aux-STATE:SKIP</i>
8	** EXTERN ! **	Control by external signal
9	** HARDWARE **	Necessary hardware does not exist
10	** ENTER ONLY **	Confirmation only possible with 
	** SIMULATION **	Simulation mode activated
No	Warning	Meaning, comments
17	** NO HV **	No HV 420 or HV 421 exists
18	** NO IS **	No IS 420 exists
19	** CANNOT DEG **	<i>DEGAS</i> not possible with <i>FILAM:1+2</i>
20	** NO AUTO **	<i>FILAM:1+2</i> not selectable with <i>DEGAS</i>
21	** ONLY F1 **	Only one filament available
22	** BUFFER **	Buffer management not OK

Other operator information is displayed in suitable locations.

7.3 Error messages

For many error types ****ERROR !!**** is displayed. To obtain detailed information press the *error* key.



NEXT displays the next error messages if more than one exists.

CLEAR deletes all messages (unresolved errors reappear immediately)

RETURN jumps to the preceding display information



Find out whether or not error messages are reproducible.

For this purpose switch all involved components OFF and ON again.

Restart the computer and the software.

For sporadic errors or errors that are difficult to reproduce → 63

The following tables help you to take appropriate action in response to reproducible error messages or faults.

This information and methods apply to the most probable cases, however, exceptions are feasible.

7.3.1 ERROR table

No	Description	Possible cause / Test methods	Correction
2	Communication CS to QC	CS 422, its cable or QC 422 defective	Replace
3	CS 422 stack overflow	CS 422 defective	Replace
4	CS 422 idle error, op. syst. overloaded		
7	CS 422 watchdog error		
14	CS 422 Display RAM	run testprogram: <i>config-TEST-CS422:EPROM-T or RAM-T</i>	Replace CS 422 if not ok
15	CS 422 EPROM checksum incorr.		
16	CS 422 RAM		
17	QMS controller stack overflow	QC 422 defective	Replace
18	Reset error		
19	QMS controller watchdog error		
20	QMS controller EPROM checksum	run testprogram: <i>config-TEST-QMS:EPROM-T</i>	Replace QC 422 if not ok
21	QMS controller NOVRAM	run testprogram: <i>config-TEST-QMS:RAM-T</i>	Replace QC 422 if not ok
22	QMS controller dual port RAM		
23	QMS controller buffer RAM		
24	DSP-EPROM checksum	run testprogram: <i>config-TEST-DSP:EPROM-T or RAM-T</i>	Replace QC 422 if not ok
25	DSP dual port RAM		
26	Monitor-DAC error	Only for factory use with special instruments	
27	Resolution-DAC error		
28	AO-DAC error		
29	ADC error		
33	QMS controller idle error (Op. sys.)	QC 422 defective	Replace QC 422
34	SEM error	SEM high voltage overloaded or defective Cabling/jumper error Faulty insulation Arcing SEM should have 18 MΩ CD should have > 100 MΩ Measure 1V per 1kV at test socket Measure actual high voltage	→ 23 and 24 Check by detaching cable Switch off unit for 2 min. Ohm meter Insulation tester DVM High voltage probe
35	CD error	CD Voltage HV 421 overloaded or defective	See Error 34
36	Ion source error	V1...V9 overloaded or defective Switch off the unit, detach IS cable: if error disappears if error persists	Check insulation of Cable /QMA F1, F2 (on IS 420) or IS 420 defective
37	Filament 1 defective	Test filament 1 or 2 ($\approx 1 \Omega$) and cable	QMA 400 → 2, IS 420 pin assignment → 16
38	Filament 2 defective		
39	Emission error (Emission \neq set point)	Pressure too high Filament burnt out (QMA 125), Filament transport protection not removed Switch <i>emiss:OFF</i> -... <i>ON</i> , if unsuccessful Adjust PROTECTION Cable interrupted or insulation fault Wiring in QMA interrupted or Insulation fault Wrong settings	$p < 10^{-4}$ mbar Check with ohm meter, replace, if defective; see QMA 125 → 4 Remove → 2, [4] → 6.2 or QME 125 → 3 Measure QMA 400 → 2 QMA 125 → 4 Test report and → 2, [4]
40	CAN error	QC 422 defective (CAN for CS 422)	Replace QC 422
41	Parameter lost	QC 422 NOVRAM defective	Replace QC 422
42	Communication from QC to CS	CS 422, its cable or QC 422 defective	Replace
43	Communication QMS-controller / DSP	Quadrupole controller QC 422 defective	Replace QC 422
44	Communication LAN	LAN communication not o.k.	Check connection, settings, parameters → 26, → 6 etc.
45	Communication RS-232-C	RS-232 communication not o.k.	Check connection, settings, parameters
46	RF error	QMH 4x0 in heat-up phase Error message from QMH 4x0	Heat-up time approx. 10 min. QMH 4x0 → 1

7.4 Measurement signal problems

Problem	Possible cause / Test methods	Correction
No measurement signal with EP 422 or CP 400	Try simulation	Replace QC 422 if not ok
	Cabling not o.k.	Check, see p. 23 u. 24
	Wrong detector selection	Set <i>detect-TYPE</i> correctly
	No emission	Switch on <i>filam</i>
	Emission too low	Adjust: see test report
	No SEM high voltage	Switch on <i>sem</i> .; see also error 34. QME 125 fuse F1 → [3]
	SEM high voltage too low	Increase: <i>detect-SEM</i> channel dependent <i>sem hv-SEM-VOLTAGE</i> channel independent See also Error 34 for HV 420
	No high voltage with HV 421	See Error 34
	Field axis voltage too low	See test report v1..v6-F-AXIS (possibly V4) QME125 → [3]
	Wrong ion source parameters	See test report, QMA400 → [2], QMA125 → [4]
	Resolution too high: try integral spectrum with <i>mass-RESOL:OFF</i> .	If integral spectrum exists adjust resolution. QMH 4x0 → [1], QME 125 → [3]
No measured value with EP 422	EP in wrong connector of QMH	Check → [23]
	SCAN-N with 10 ms/u, range 10^{-9} A Disconnect EP 422 from QMA, touch input with screwdriver	50 (60) Hz signal should appear Check signal with an oscilloscope, if not o.k. replace EP 422
No measured value with CP 400	Threshold too high	Decrease <i>amplif-CP-LEV</i> → [48]
Electrometer signal negative	EP1 and EP2 cable mixed up	Check, see p. 23
	POLARITY switch on QME 125 set to "-"	Set POLARITY to "+"
Electrometer offset strongly mass dependent	Loose ground connection (EP input or below QMA connector plate)	Tighten / correct
	Open shielding below QMA connector plate	QMA 400 → [2], QMA 125 → [4]
Electrometer signal: - not zero between peaks - negative / small peaks missing	Offset not aligned	Perform offset correction → [47]
Offset in range 10^{-12} very high	Temperature of the EP 422 too high	Decrease
	Bad insulation collector to flange (good: $>>1$ GΩ)	Correct insulation fault → QMA 400 → [2], QMA 125 → [4]
	Moisture in electrometer or on analyzer connector	Dry with warm air (no over 60 °C) → [64]
Electrometer signal sensitive to vibrations	Knurled nut loose on EP input Shield below connector plate of QMA not correctly installed	Tighten Correct shielding QMA 400 → [2], QMA 125 → [4]
High noise signals with EP 422	Analyze signal (<i>amplif-MODE:FIX</i> , <i>SCAN-N</i> , high <i>SPEED</i>) with PC or <i>elm</i> signal with oscilloscope	Remedy noise or choose slower <i>SPEED</i> , <i>DWELL</i> and/or increase <i>amplif-FILTER</i>
High counting rate with CP 400 also besides peaks	Corona or arcing in CP 400 or in HV 420 in HV 421 in high voltage cables in QMA or SEM	Open CP 400 and dry with war air (<50°C), Remove dust. Replace HV 421 Replace cable QMA 400 → [2], QMA 125 → [4]
High counting rate with CP 400 also without high voltage	Poor ground connection, shielding open, coupling of parasitic signals, e.g. with isolated system set-up.	Correct Establish shielding or decouple
Measurement signal: -Limited to values <10 V -Jumps to 10.24 V	With <i>SEM</i> and <i>FIX-range</i> 10^{-11} and 10^{-12} A	EP 422 overdriven Use Autorange or <i>RANGE</i> 10^{-10} and higher SEM voltage

Problem	Possible cause / Test methods	Correction
Unsatisfactory peak shape, poor sensitivity	Small emission (0.1 mA) cannot be set on QME 125	Replace the insulators in the ion source, see QMA 125 → [4]
	Ion source insulation in analyzer bad (good: >100 MΩ)	Replace the insulators in the ion source, see QMA 400 → [2], QMA 125 → [4]
Unsatisfactory peak shape, poor sensitivity	Ion source or rod system in analyzer contaminated or defective	Clean QMA 400 → [2], QMA 125 → [4]
	SEM voltage too low	See Error 34
	SEM contaminated or defective	Replace SEM, QMA 400 → [2], QMA 125 → [4]
Peaks become wider/narrower with increasing mass number	Incorrect setting of <i>resolution coarse</i>	Adjust QMH 4x0 → [1], QME 125 → [3]

7.5 General problems

Problem	Possible cause / Test methods	Correction
Fans not running, no indication on CS 422 or <i>power-LED</i>	Line voltage missing or too low	Check line voltage
	Power cable defective	Replace power cable
	Short circuit in external unit such as QMH 4x0 or QME 125	Switch off / unplug ext. equipment / switch on again. If unsuccessful replace the defective unit.
	Short circuit on bus or in wiring Defective power supply	Trace Replace power supply (manufacturer's warranty becomes void if the power supply is opened)
Fans running, display on CS 422 is blank	Contrast strongly out of adjustment	Adjust contrast with [Grid Icon]
	CS 422 or its cable defective	Replace
	QC 422 not correctly installed or defective	→ [29] or replace
CS 422 keys dead	Control via interface	Manual control: <i>config-CTRL-MODE:CS 422</i>
QMH 4x0 connector does not fit into QC 422	QMH cable with old locking device	Use adapter → [65]
Sporadic error messages	EMC problems	Correct the ground connection → [21] Detach QMG cable from noise sources Identify noise source and eliminate noise Use LAN (fiber optics)
	Line voltage dips	Check supply voltage quality
No ArcNet communication	Check status LEDs of the OPA 200, check FO connection	OPA 200 → [6]
	FO connectors contaminated	Clean, e.g. with alcohol
	HUB has no power	Check
	Stray light	Mount caps on all unused FO connectors!
	Wrong settings	Check jumper settings on all LAN units QC 422, OH 421 → [26] and → [6], [7]
	Wrong node addresses	Correct
	QC 422 defective	Replace QC 422
	Defective LAN port device	Reduce system to minimum and then rebuild it in steps

7.6 Service interventions

If you are unable to remedy a fault or if you are not allowed to do so due to the lack of skilled personnel, please contact the responsible service location.

If you need advice or if you want to return the equipment for repair, please supply a comprehensive description of the error together with:

- Description of fault, e.g. hard copies, recorder charts and text
- Application conditions and operating modes under which the error occurs
- In case of sporadic errors all observations that could help to reproduce the error
- Type, series, software and firmware numbers of all components involved.

Products that have been exposed to vacuum conditions must always be accompanied by a completed declaration of contamination (→ 69).

8 Maintenance

The QMS 422 and QMI 422 have lateral ventilation inlets. Their filters are to be cleaned before the air circulation becomes obstructed. The cleaning interval depends on the local dust evolution. Dry dust can easily be removed with the aid of a vacuum cleaner.



Skilled personnel

If necessary remove the filters and wash them in a mild soap solution. Dry them well before you reinstall them!

Defective filters should be replaced (→ 65).

The installed fans should be checked semi-annually. Replace them if they are not running smoothly or are overly noisy.

Dusty circuit boards can be cleaned with compressed air (max. 2 bar). Make sure that no components get damaged or bent.

Moisture (condensation) in the EP 422 can lead to unstable behavior (offset fluctuations). Open the EP 422 and dry it with a hair dryer (max. 60 °C).

Refer to the maintenance instructions in all the user's guides of the components that form part of the system. See list of literature on 68.

9 Decommissioning

Please contact your Pfeiffer Vacuum service location on instructions of how to dispose of your system.

10 Spare parts and accessories

	Ordering number
QMS 422 Basic unit without QC 422	PT 444 580-T
CS 422 Operator console	PT 444 650-T
1 Filter mat for QMS 422	B 5099 154 FD
1 Ventilator for QMS 422	B 5099 130 CD
1 Blanking plate, 4 subunits (20 mm wide)	BG 544 775-T
Power supply for QMS 422 and QMI 422 (90...265 V)	B 5181 214 QY
QMI 422 Control unit with QC 422	PTM D27 280
1 Air filter mat for QMI 422	B 5099 154 FB
EP 422 Electrometer	BG 444 570-T
Input cable TNC/TNC, 200°C, l = 0.5 m	B 4564 401 EB
Input cable TNC/TNC, 70°C, l = 6 m, low-noise	B 4564 401 E2
TNC short circuit plug (fits QMA)	B 4728 138 BC
QC 422 Quadruple controller without options	BG 444 590 -T
AO 421 Analog output (incl. connector)	BG 442 328-T
IC 421 Ion counter and analog output	BG 442 320-T
Ion counter preamplifier CP 400	BG 442 210-T
Cable CP 400-QC 422: 3m	BG 448 134-T
Cable CP 400-QC 422: 10m	BG 448 199-T
Ion source supply IS 420	BG 512 900-T
Ion source cable 3 m	BG 548 082-T
Ion source cable 10 m	BG 548 083-T
Fuses F1, F2 2.5 A slow	B 4666 444
HV 420 SEM high voltage supply	BG 546 040-T
Fuse F1 0.2 AT	B 4666 422
HV 421 SEM high voltage supply	BG 442 250-T
HV cable 3m	BG 541 978-T
HV cable 10m	BG 541 979-T
AI 421 Analog input (incl. connector)	BG 442 240 -T
DI 420 Digital input	BG 512 830 -T
1 Connector housing	BG 531 194-T
1 Multipoint connector, solder version	B 4717 306 DL
DO 420A Digital output	BG 512 842 -T
Connector, see DI 420	
PI 420 Pirani module	BG 512 715 -T
PE 420 Cold cathode module	BG 512 726 -T
OH 421 Optical hub	BG 442 465-T
OPA 200 Optical PC Arcnet interface	B 5278 503 KT
OHA 200 Optical hub 5-port	PT 442 510-T
OHA 200 Optical hub 10-port	BG 442 520-T
Fiber-optic conductor PCF 10m	B 5159 615 2H
Fiber-optic conductor PCF 20m	B 5159 615 2K
Fiber-optic conductor PCF 50m	B 5159 615 2Q
Fiber-optic conductor APF 1 m	B 5159 615 2C
Fiber-optic conductor APF 3 m	B 5159 615 2D
Other lengths on request	
Adapter for QMH 400 / QME125 with sliding lock	B 4720 786 CD

Appendix

A: Default parameter values

For activating the default parameters see Parameter *INIT*.

Channels

Parameter	Function	Default value	Parameter	Function	Default value
AI-CH	detect	0	MODE	amplif	FIX
AO-CH	output	1	MODE	mass	SCAN-N
AO-MODE	output	LIN	MONITOR	output	LIN
AVERAGE	mass	1	O-RNG	output	E-1
CALIB	amplif	1.000 E0	P-CAL	amplif	1.0
COPY TO XX	aux	0	PE-CTRL	detect	OFF
CP-LEV	amplif	0.00 V	PI-CH	detect	1
DO-A	trip	OFF	RANGE	amplif	E-5
DO-B	trip	OFF	RANGE-L	amplif	E-5
DWELL	mass	1 s/u	RESOL	mass	25
FILTER	amplif	AUTO	SEM	detect	SEM-HV
FIRST	mass	14.00	SPEED	mass	1 s/u
GAIN	amplif	×1	STATE	aux	ENABLE
LEVEL-A	trip	1.00 E-6	THRESH	mass	0.3 %F.S.
LEVEL-B	trip	1.00 E-5	TYPE	detect	SEM
LOG-DEC	output	3 DEC	TYPE	trip	ABS
MASS	mass	14.00	WIDTH	mass	+16

General

Parameter	Function	Default value	Parameter	Function	Default value
BAUD	config	2400 Bit/s	MODE	config	CS 422
DETECT	config	SEM	NODE	config	176
IS-TYP	config	CB	OPTION	config	NO
MASS-R	config	512	QMA	config	400
			SEM+FIL	config	INTERN

Ion source










400

Parameter	Function	Axial	CB	Grid	SPM	Spec+/-
E-PROT	emiss	4.40 A	4.40 A	4.40 A	3.50 A	0.00 A
EMISS	emiss	1.00 mA	1.00 mA	1.00 mA	0.50 mA	OFF
V1	v1...v6	90 V	90 V	90 V	40 V	0 V
V2	v1...v6	70.0 V	70.0 V	70.0 V	40.0 V	0 V
V3	v1...v6	+20.0 V	+20.0 V	0.0 V	0,0 V	0 V
V4	v1...v6	15 V	15 V	15 V	3 V	0 V
V5	v1...v6	0 V	250 V	0 V	0 V	0 V
V6	v1...v6	300 V	300 V	300 V	110 V	0 V
V7	v7...	0 V	0 V	0 V	0 V	0 V
V8	v7...	0 V	0 V	0 V	0 V	0 V
V9	v7...	30 V	0 V	0 V	30 V	0 V

Operation

Parameter	Function	Default value	Parameter	Function	Default value
<i>ADJ-TYP</i>	<i>cycle</i>	COARSE	<i>FIL1</i>	<i>ion src</i>	SET 0
<i>BEGIN</i>	<i>cycle</i>	0	<i>FIL2</i>	<i>ion src</i>	SET 0
<i>COPY TO SET</i>	<i>ion src</i>	SET 0	<i>FUNCT</i>	<i>cycle</i>	CYCLE
<i>CYCLES</i>	<i>cycle</i>	REPEAT	<i>MODE</i>	<i>cycle</i>	MONO
<i>D-EMIS</i>	<i>ion src</i>	10.0 mA	<i>MODE</i>	<i>ion src</i>	NORMAL
<i>D-PROT</i>	<i>ion src</i>	4.00 A	<i>SEM-VOLTAGE</i>	<i>sem hv</i>	1500 V
<i>D-TIME</i>	<i>ion src</i>	10 min	<i>TRIG</i>	<i>cycle</i>	INTERN
<i>END</i>	<i>cycle</i>	63	<i>TYPE</i>	<i>ion src</i>	CB
<i>FILAM</i>	<i>ion src</i>	1			

B: Literature

-  [1] Operating instructions RF Generator QMH 400 / 410
BG 805 982 BE
-  [2] Operating instructions QMA 400 / 410 / 430 Analyzer
BK 805 983 BE
-  [3] Operating instructions Quadrupole electronics QME 125
BG 800 987 BE
-  [4] Operating instructions Analyzer QMA 125
BG 805 988 BE
-  [5] Communication protocol Quadrupole Controller QC 422
BG 805 986 BE
-  [6] User's guide Network Controller Board OPA 200
vina02d3
-  [7] Operating manual Optical Hub OHA 200
-  [8] Operating instructions Pirani module PI 420
-  [9] Operating instructions Penning Module PE 420

Ordering source
Pfeiffer Vacuum GmbH Emmeliusstrasse 33 D-35614 Asslar

C: Index

Note:

Information on the parameters can be found in the alphabetical list beginning on page 35

Accessories	65	Group	Parameter	
ADJUST		channels	display	31
coarse	52	general	Parameters	
fine	54	ion source	alphabetic list	35
AI 421	10; 20; 30	operation	Pause	49
Analog filter	13	HV 420	PE 420	
Analog measured		HV 421	Penning module	10; 30
value output	55	IC 421	PEAK	
Analog signal		Installation	maximum	52
Extern	48	Integral spectrum	processing	51
AO Analog		Ion counter	PI 420	
outputs	55	Ion polarity	Pirani module	10; 30
AO 421	10; 29	IS 420	Programming	68
ARCNET		LAN interface	Protection	47
Interface	14; 26	Literature	QC 422	12
network	27	Maintenance	QMA 125	9
Average	51	Mass scan	QMA 400/410/430	9
Bus module		modes	QME 125	9
installation	29	Measured value	QMH 400/410	9
Cabling	23	display	QMI 422	8
CD Voltage		recording	QMS 422	8
adjustment	30	Measured value output	Recording	58
Configuration	46	analog	Resolution	
Connections		lin, log	Mass scale	12
analog signals	55	Measurement cycle	Measured value	13
on IS 420	16	halt, mono, multi	RS-232-C	14; 25
on QC422	15	Measurement signal	Safety symbols and	
CP 400	10; 17; 22	problems	information	5
CS 422	10; 28	mon	SAMPLE	51
DI 420	10; 19; 30	Monitor output	scan	
DO 420A	10; 19; 30	Network	Mass number signal	55
Electrometer		Offset	SCAN-F	50
Auto, Fix	50	correction	SCAN-N	50
elm		OH 421	Service	
Electrometer output	55	OHA 200	interventions	64
EP 422	17; 22	Operation	Simulation	48
Error		Options	Spare parts	65
messages	60	installing/	STAIR	52
Extern Input	48	removing	Switching	
Fault isolation	60	Overview	functions	58
Filament		Ordering numbers	sync	
protection	47	Overview	Synchronization signal	55
FIR-Filter	50	series 125/400 family	Technical data	11
First time			TRIP	58
operation	46		Troubleshooting	60
			Vacuum	
			measurement	10
			Validity	2
			Zero line	47

Declaration of Contamination

The service, repair, and/or disposal of vacuum equipment and components will only be carried out if a correctly completed declaration has been submitted. Non-completion will result in delay.
This declaration may only be completed (in block letters) and signed by authorized and qualified staff.

1 Description of product Type _____ Part number _____ Serial number _____	2 Reason for return _____ _____ _____
---	---

3 Operating fluid(s) used (Must be drained before shipping.) _____ _____

4 Process related contamination of product <table style="width: 100%;"> <tr> <td style="width: 40%;">toxic</td> <td style="width: 20%;">no <input type="checkbox"/> 1)</td> <td style="width: 20%;">yes <input type="checkbox"/></td> </tr> <tr> <td>caustic</td> <td>no <input type="checkbox"/> 1)</td> <td>yes <input type="checkbox"/></td> </tr> <tr> <td>biological hazard</td> <td>no <input type="checkbox"/></td> <td>yes <input type="checkbox"/> 2)</td> </tr> <tr> <td>explosive</td> <td>no <input type="checkbox"/></td> <td>yes <input type="checkbox"/> 2)</td> </tr> <tr> <td>radioactive</td> <td>no <input type="checkbox"/></td> <td>yes <input type="checkbox"/> 2)</td> </tr> <tr> <td>other harmful substances</td> <td>no <input type="checkbox"/> 1)</td> <td>yes <input type="checkbox"/></td> </tr> </table>	toxic	no <input type="checkbox"/> 1)	yes <input type="checkbox"/>	caustic	no <input type="checkbox"/> 1)	yes <input type="checkbox"/>	biological hazard	no <input type="checkbox"/>	yes <input type="checkbox"/> 2)	explosive	no <input type="checkbox"/>	yes <input type="checkbox"/> 2)	radioactive	no <input type="checkbox"/>	yes <input type="checkbox"/> 2)	other harmful substances	no <input type="checkbox"/> 1)	yes <input type="checkbox"/>	 2) Products thus contaminated will not be accepted without written evidence of decontamination!
toxic	no <input type="checkbox"/> 1)	yes <input type="checkbox"/>																	
caustic	no <input type="checkbox"/> 1)	yes <input type="checkbox"/>																	
biological hazard	no <input type="checkbox"/>	yes <input type="checkbox"/> 2)																	
explosive	no <input type="checkbox"/>	yes <input type="checkbox"/> 2)																	
radioactive	no <input type="checkbox"/>	yes <input type="checkbox"/> 2)																	
other harmful substances	no <input type="checkbox"/> 1)	yes <input type="checkbox"/>																	

The product is free of any substances which are damaging to health
 yes ☐

1) or not containing any amount of hazardous residues that exceed the permissible exposure limits

5 Harmful substances, gases and/or by-products Please list all substances, gases, and by-products which the product may have come into contact with:			
Trade/product name Manufacturer	Chemical name (preferably with formula)	Precautions associated with substance	Action in case of human contact

6 Legally binding declaration: We hereby declare that the information on this form is complete and accurate and that we will assume any further costs that may arise. The contaminated product will be dispatched in accordance with the applicable regulations.	
Organization/company _____ Address _____ Phone _____ E-mail _____ Name _____	Post code, place _____ Fax _____
Date and legally binding signature _____ <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"></div> <div style="width: 45%; text-align: center;"> Company stamp _____ </div> </div>	

Original: German BG 805 981 BD (0201)



bg805981be

PFEIFFER  **VACUUM**

*Berliner Strasse 43
D-35614 Asslar
Deutschland
Tel +49 (0) 6441 802-0
Fax +49 (0) 6441 802-202
info@pfeiffer-vacuum.de*

www.pfeiffer-vacuum.de