



**Elastomer seal, FKM,
DN 16 CF**




Elastomer seal



Elastomer seal, FKM, DN 16 CF

- Reusable

Dimensions



Technical Data

Technical Data	Elastomer seal, FKM, DN 16 CF
Temperature range min.	-15 °C
Temperature range max.	150 °C
Package	2 pcs.
O-ring material	FKM
Temperature range	-15 – 150 °C
Connection flange	DN 16 CF
Pressure range	1E-8 hPa – to ambient pressur · 10
Temperature range max.	150 °C 302 °F 423.15 K
Temperature range min.	-15 °C 5 °F 258.15 K

Order number

Order number	Elastomer seal, FKM, DN 16 CF
Elastomer seal, FKM, DN 16 CF	402DFL016-S2

Your Success. Our Passion.

We give our best for you every day –
worldwide!

Are you looking for an optimum vacuum solution?

Talk to us:

Pfeiffer Vacuum Components & Solutions GmbH Germany

T +49 551 99963-0

info-cs@pfeiffer-vacuum.com

Or scan the barcode, to visit our web page:



<https://webportal.pfeiffer-vacuum.com/global/en/contact>



Errors and/or changes excepted. - 11/1/2024

Follow Us On Social Media
#pfeiffervacuum



www.pfeiffer-vacuum.com

PFEIFFER  **VACUUM**